

K071895

NOV 14 2007

STERIS®



**510(k) Summary
For
Verify® 275F 3 Indicator**

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060
Phone: (216) 354-2600
Fax No: (216) 639-4459

Contact: Jack Scoville
Fellow, Regulatory Affairs

Telephone: 440 392 7330
Fax No: 440 397 9198

Submission Date: July 06, 2007

**STERIS ABBREVIATED 510(k) PREMARKET NOTIFICATION
VERIFY® 275F 3 INDICATOR**

1. Device Name

Indicator: Verify® 275F 3 Indicator.

Common Name: Chemical Indicator.

Classification Name: Physical/chemical sterilization process indicator (21 CFR 880.2800 (b), Product Code JOJ).

2. Predicate Devices

- 3M SteriGage Chemical Indicators (K894399)
- DANA SteriScan Indicators (K012195)
- SteriTec Integrgraph (Cardinal Steam Integrators¹) (K960441)
- STERIS Verify Integrators (K902958, K002937)
- EZ Test Steam Biological Indicator (K963841)

3. Device Description

The proposed Verify® 275F 3 Indicator consist of:

- A 22 mm x 143 mm polypropylene strip with one 12 mm chemical indicator ink spot.

The indicator ink spot is located on one end of the strip, adjacent to a reference block exhibiting the endpoint color. The indicator ink on the proposed Verify® 275F 3 Indicator changes from yellow to blue/purple color when the steam sterilization cycle is complete.

The Verify® 275F 3 Indicator can be used to monitor 3 minute 275°F (135°C) pre-vacuum steam sterilization cycles.

4. Intended Use

The Verify® 275F 3 Indicator is a chemical indicator intended for use by health care providers to accompany products being sterilized through a sterilization procedure. The Verify® 275F 3 Indicator is an emulating indicator intended for use in steam sterilization. The indicator changes color from yellow to blue/purple when exposed to 275°F (135°C) pre-vacuum sterilization cycles for 3 minutes. The performance of

¹ Cardinal is a private label brand produced by Steritec under K960441.

**STERIS ABBREVIATED 510(k) PREMARKET NOTIFICATION
VERIFY® 275F 3 INDICATOR**

the Verify® 275F 3 Indicator meets the requirements of ANSI/AAMI / ISO 11140-1:2005 for emulating [Class 6] steam indicators.

5. Description of Safety and Substantial Equivalence

The proposed and predicate devices are all single use indicators for use in steam sterilization cycles. The differences between the proposed Verify® 275F 3 Indicator and the predicate device are limited to differences in design, material, and parameters of the sterilization cycles these indicators are designed to monitor. These differences do not raise any new issues of safety and efficacy.

6. Performance Testing

Performance testing was conducted to verify that the proposed Verify® 275F 3 Indicator meet the requirements for emulating [Class 6] indicators as defined in ANSI/AAMI ISO 11140-1:2005 using a resistometer to ANSI/AAMI ST44.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

NOV 14 2007

Mr. John Scoville
Fellow, Regulatory Affairs Sterilization Technology
STERIS Corporation
5960 Heisley Road
Mentor, Ohio 44060

Re: K071895
Trade/Device Name: Verify[®] 275F 3 Indicator
Regulation Number: 880.2800
Regulation Name: Sterilization process Indicator
Regulatory Class: II
Product Code: JOJ
Dated: October 25, 2007
Received: October 26, 2007

Dear Mr. Scoville:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

STERIS ABBREVIATED 510(k) PREMARKET NOTIFICATION
VERIFY® 275F 3 INDICATOR

Indications for Use

510(k) Number (if known): K071895

Device Name: Verify® 275F 3 Indicator

Indications for Use:

The Verify® 275F 3 Indicator is a chemical indicator intended for use by health care providers to accompany products being sterilized through a sterilization procedure. The Verify® 275F 3 Indicator is an emulating indicator intended for use in steam sterilization. The indicator changes color from yellow to blue/purple when exposed to 275°F (135°C) pre-vacuum sterilization cycles for 3 minutes as indicated in the following table:

MODEL	TEMPERATURE	STERILIZATION TYPE	TIME
Verify 275F 3	275°F (135°C)	Pre-vacuum steam	3 minutes

Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use X
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Shela H. Murphy
Biotechnology, General Use
Medical Device
K071895

Page 1 of 1