1031161

JUL 2 8 2008

### 510 (k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

1. Manufacturer's Name, Address, Telephone, and Contact Person, Date of Preparation:

Manufacturer: Dade Behring Marburg GmbH Emil-von-Behring Str. 76 35041 Marburg, Germany

**Contact Information:** 

Helen M. Lee Siemens Healthcare Diagnostics Inc. 500 GBC Drive, M/S 514 Newark, DE 19714-6101 302.631.8706 302.631.6299 (fax)

Date of Preparation: July 17, 2008

2. Name of Products: Dimension Vista® System IGG Flex® reagent cartridge

Dimension Vista® System Protein 1 Calibrator Dimension Vista® System Protein 3 Control

3. Identification of the Legally marketed Device:

Beckman Coulter IMMAGE® Immunochemistry System Urine Immunoglobulin G (k951635) N Antisera to Human Immunoglobulins (IgG, IgA, and IgM) (k042735) Dimension Vista® System Protein 1 Calibrator (k073561) Dimension Vista® System Protein 3 Control (k073561)

#### 4. Device Descriptions:

Dimension Vista® System Immunoglobulin G Flex® reagent cartridge

Proteins contained in human body fluids from immune complexes in an immunochemical reaction with specific antibodies. These complexes scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the respective protein in the sample. The result is evaluated by comparison with a standard of known concentration.

#### Dimension Vista® System Protein 1 Calibrator

PROT1 CAL is a multi-analyte, liquid human serum based product containing:

a\_-Acid Glycoprotein

Immunoglobulin G

 $\alpha_1$ -Antitrypsin

Immunoglobulin G Subclass 1

 $\beta_2$ -Microglobulin

Immunoglobutin G Subclass 2

C3 Complement

Immunoglobulin G Subclass 3

C4 Complement

Immunoglobulin G Subclass 4

Ceruloplasmin

Immunoglobulin M

Haptoglobin

Prealbumin

Hemopexin Homocysteine Retinol Binding Protein soluble Transferrin Receptor

Immunoglobulin A

Transferrin

Immunoglobulin E

Dimension Vista® System Protein 3 Control

PROT3 CON is a multi-analyte, lyophilized, polygeline and rabbit albumin based product containing:

a\_-Microglobulin

Immunoglobulin G

Albumin

#### 5. Device Intended Uses:

# Dimension Vista Immunoglobulin G Flex reagent cartridge

The IGG method is an *in vitro* diagnostic test for the quantitative measurement of immunoglobulin G in human serum, heparinized plasma, cerebrospinal fluid (CSF) and urine on the Dimension Vista System. Measurements of IgG aid in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents.

#### Dimension Vista® System Protein 1 Calibrator

PROT1 CAL is an in vitro diagnostic product for the calibration of the Dimension Vista System for:

α<sub>4</sub>-Acid Glycoprotein (A1AG)

Immunoglobulin G (IGG) [serum/plasma]

Immunoglobulin G (IGG-C) [cerebrospinal fluid]

Immunoglobulin G (IGG-U) [urine]

α<sub>1</sub>-Antitrypsin (A1AT)

Immunoglobulin G Subclass 1 (IGG1)

 $\beta_{a}$ -Microglobulin (B2MIC)

Immunoglobulin G Subclass 2 (IGG2)

C3 Complement (C3)

Immunoglobulin G Subclass 3 (IGG3)

C4 Complement (C4)
Ceruloplasmin (CER)
Haptoglobin (HAPT)
Hemopexin (HPX)
Homocysteine (HCYS)
Immunoglobulin A (IGA)
Immunoglobulin E (IGE)

Immunoglobulin G Subclass 4 (IGG4) Immunoglobulin M (IGM) Prealbumin (PREALB) Retinol Binding Protein (RBP) soluble Transferrin Receptor (STFR) Transferrin (TRF)

#### Dimension Vista® System Protein 3 Control

PROT3 CON is an assayed quality control for the assessment of precision and analytical bias on the Dimension Vista System in quantitative determination of:

α<sub>4</sub>-Microglobulin (A1MIC) Immunoglobulin G (IGG-C)\* Immunoglobulin G (IGG-U) \*\* Microalbumin (MALB)
specialty Albumin (sALB) \*

\* For cerebrospinal fluid (CSF)

### 6. Medical Device to which equivalence is claimed and comparison information:

The Beckman Coulter IMMAGE® Immunochemistry System Urine Immunoglobulin G was determined to be substantially equivalent in 510(k) Premarket Notification k951635 and is the predicate for adding the urine sample matrix to the previously cleared Dimension Vista® IGG assay.

The N Antisera to Human Immunoglobulins (IgG, IgA, and IgM) on the BN Systems was determined to be substantially equivalent in 510(k) Premarket Notification k042735 and is the predicate for the addition of a pre-reaction step to the previously cleared IGG-C (CSF). The operating principle and reagent composition have not changed. For your convenience, the sections with changes have been highlighted in the Draft Instructions for Use included in this submission.

Dimension Vista® Protein 1 Calibrator, modified to include IGG-U (IGG in urine), is substantially equivalent in intended use to the current Dimension Vista® Protein 1 Calibrator (k073561). The modified Dimension Vista® Protein 1 Calibrator, like the current product, is intended to be used for the calibration of human protein assays on the Dimension Vista® System.

Dimension Vista® System Protein 3 Control modified to include IGG-U (IGG in urine), is substantially equivalent in Intended Use to the current Dimension Vista® System Protein 3 Control (K073561). The modified Dimension Vista® System Protein 3 Control like the current product is intended to be used as an assayed intralaboratory quality control on the Dimension Vista® System.

<sup>\*\*</sup> For urine

#### 7. Conclusion:

The studies included in this submission demonstrate correlation to and equivalent performance between the predicate Beckman Coulter IMMAGE® Immunochemistry System Urine Immunoglobulin G for urine sample matrix and the predicate N Antisera to Human Immunoglobulins (IgG IgA, and IgM) on the BN ProSpec® System for the addition of a pre-reaction step to the Dimension Vista® IGG-C assay.

The modified Dimension Vista® System IGG assay, modified Dimension Vista® System Protein 1 Calibrator and modified Dimension Vista® System Protein 3 Control are substantially equivalent to the legally marketed devices based upon the information above.







Food and Drug Administration 2098 Gaither Road Rockville MD 20850

Siemens Healthcare Diagnostics, Inc. c/o Ms. Helen M. Lee Regulatory Affairs and Compliance Manager 500 GBC Drive, MS 514 Newark, DE 19714-6101

JUL 2 8 2008

Re: k081161

Trade/Device Name: Dimension Vista® System Immunoglobulin G Flex® Reagent Cartridge

Dimension Vista® System Protein 1 Calibrator Dimension Vista® System Protein 3 Control

Regulation Number: 21 CFR 866,5510

Regulation Name: Immunoglobulins A, G, M, D and E immunological test system

Regulatory Class: Class II Product Code: CFN, JIX, JJY

Dated: July 2, 2008 Received: July 3, 2008

Dear Ms. Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding

of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <a href="http://www.fda.gov/cdrh/industry/support/index.html">http://www.fda.gov/cdrh/industry/support/index.html</a>.

Sincerely yours,

Maria M. Chan, Ph.D.

maria m chan

Acting Division Director

Division of Immunology and Hematology Devices Office of In Vitro Diagnostic Device Evaluation and Safety

Center for Devices and Radiological Health

Enclosure

## **Indications for Use Statement**

510(k) Number (if known): 1081161

Device Name:		
Dimension Vista® System Protein 3 Control		
Indications for Use:		
PROT3 CON is an assayed, low level intralaboratory quality control for assessment of		
precision and analytical bias on the Dimension Vista System in the determination of α <sub>1</sub> -Microglobulin (A1MIC), specialty Albumin (sALB)*, Immunoglobulin G (IGG -C)*, Immunoglobulin G (IGG-U)** and Microalbumin (MALB).		
* For Cerebrospinal fluid (CSF)  ** For urine		
Prescription Use X Over-The-Counter-Use		
(Per 21 CFR 801 Subpart D) (21 CFR 801)		
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)		
Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)		
Division Sign-Off		
Division Sign-Off		

Office of In Vitro Diagnostic Device Evaluation and Safety

## **Indications for Use Statement**

510(k) Number (if known):   C8   16		
Device Name:		
Dimension Vista® System Protein 1 Calibrator		
Indications for Use: PROT1 CAL is an <i>in vitro</i> diagnostic product system for: α <sub>1</sub> -Acid Glycoprotein (A1AG)	for the calibration of the Dimension Vista <sup>®</sup> Immunoglobulin G (IGG) [serum/plasma], Immunoglobulin G (IGG-C) [cerebrospinal	
α <sub>1</sub> -Antitrypsin (A1AT) β <sub>2</sub> -Microglobulin (B2MIC) C3 Complement (C3) C4 Complement (C4) Ceruloplasmin (CER) Haptoglobin (HAPT) Hemopexin (HPX) Homocysteine (HCYS) Immunoglobulin A (IGA) Immunoglobulin E (IGE)	fluid] and Immunoglobulin G (IGG-U)  [urine]  Immunoglobulin G Subclass 1 (IGG1)  Immunoglobulin G Subclass 2 (IGG2)  Immunoglobulin G Subclass 3 (IGG3)  Immunoglobulin G Subclass 4 (IGG4)  Immunoglobulin M (IGM)  Prealbumin (PREALB)  Retinol Binding Protein (RBP)  soluble Transferrin Receptor (STFR)  Transferrin (TRF)	
Prescription Use X (Per 21 CER 801 Subport D)	Over-The-Counter-Use	
(Per 21 CFR 801 Subpart D)	(21 CFR 801)	
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)		
Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)		
Mana m chan Division Sign-Off		
Office of In Vitro Diagnostic		

# **Indications for Use Statement**

510(k) Number (if known): K081161		
Device Name:		
Dimension Vista® System Immunoglobulin G Flex® Reagent cartridge		
Indications for Use: The IGG method is an <i>in vitro</i> diagnostic test for the quantitative measurement of immunoglobulin G in human serum, heparinized plasma, cerebrospinal fluid (CSF) and urine on the Dimension Vista® System. Measurements of IgG aid in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents.		
Prescription Use X Over-The-Counter-Use (Per 21 CFR 801 Subpart D) (21 CFR 801)		
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)		
Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)		
Mara In chan Division Sign-Off		
Office of <b>in Vitro Diagnostic</b> Device Evaluation and Safety		

510(K) KO&1161