

K 083060

**II. SUMMARY AND CERTIFICATION**

**A. 510(k) Summary**

**Submitter:** SterilMed, Inc.  
**Contact Person:** Dennis Toussaint  
 11400 73<sup>rd</sup> Avenue North  
 Maple Grove, MN 55369  
 Ph: 763-488-3410  
 Fax: 763-488-2051

**Date Prepared:** October 10, 2008  
**Trade Name:** Reprocessed Harmonic Scalpel  
**Classification Name:** Scalpel, Ultrasonic, Reprocessed  
**Classification Number:** Unclassified  
**Product Code:** NLQ

JAN - 9 2009

<b>Predicate Devices:</b>	The reprocessed harmonic scalpels are substantially equivalent to Ethicon Harmonic WAVE™ harmonic scalpels.
<b>Device Description:</b>	<p>SterilMed reprocessed harmonic scalpels are used in combination with a hand piece, generator and torque wrench and are intended to be used in soft tissue surgery for simultaneous cutting and coagulation of vessels and tissue. The instrument has a scissor handle with hand control capabilities consisting of MIN and MAX buttons. The handle housing has an integrated mechanism for limiting the force that can be applied when closing the distal mechanism. The instrument has an 18 cm shaft length, 8.5 mm shaft diameter, active blade length of 18 mm, and utilizes a straight blade and clamp arm.</p> <p>Note: Only the harmonic scalpel is the subject of this submission, the reusable hand piece, generator, and any other related equipment are not included in the scope of this submission.</p>
<b>Intended Use:</b>	<p>The reprocessed harmonic scalpels are indicated to be used for cutting of soft tissue and providing hemostasis when control of bleeding and minimal thermal injury is desired.</p> <p>The instrument can be used as an adjunct to, or a substitute for, electrosurgery, lasers, and steel scalpels in abdominal, pediatric, gynecologic, exposure to orthopedic structures (such as spine and joint space) and other open and endoscopic procedures.</p>
<b>Functional and Safety Testing:</b>	Representative samples of reprocessed harmonic scalpels were tested to demonstrate appropriate functional characteristics. Process validation testing was performed to validate the cleaning and sterilization procedures as well as device packaging. In addition, the manufacturing process includes visual and validated functional testing of all products produced.
<b>Conclusion:</b>	<p>The reprocessed harmonic scalpels are substantially equivalent to Ethicon Harmonic WAVE™ harmonic scalpels.</p> <p>This conclusion is based upon the devices' similarities in functional design (principle of operation), materials, indications for use and methods of construction.</p>



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

SterilMed, Inc.  
% Mr. Dennis Toussaint  
Director of Regulatory Affairs  
11400 73<sup>rd</sup> Avenue North  
Maple Grove, Minnesota 55369

JAN - 9 2009

Re: K083060  
Trade/Device Name: Reprocessed Harmonic Scalpels  
Regulatory Class: Unclassified  
Product Code: NLQ, LFL  
Dated: October 10, 2008  
Received: October 14, 2008

Dear Mr. Toussaint:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K083060

Device Name: Reprocessed Harmonic Scalpels

### Indications for Use:

The reprocessed harmonic scalpels are indicated to be used for cutting of soft tissue and providing hemostasis when control of bleeding and minimal thermal injury is desired.

The instrument can be used as an adjunct to or a substitute for electrosurgery, lasers, and steel scalpels in abdominal, pediatric, gynecologic, exposure to orthopedic structures (such as spine and joint space) and other open and endoscopic procedures.


Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
Division of General, Restorative,  
and Neurological Devices

510(k) Number \_\_\_\_\_

K083060

**Devices included in this Premarket Notification Submission – 510(k) K083060**

Manufacturer	Model#
Ethicon	WAVE18S