

K090751

## Attachment IV.2 – 510(k) Summary (revised)

### General Information

OCT 28 2009

**Owner's Name:** Apollo Spine  
**Address:** 307 Placentia Ave, Suite 111  
Newport Beach CA 92663  
**Telephone Number:** (949) 645-7746  
**Fax Number:** (949) 645-7749  
**Contact Person:** Kamran Aflatoon

**Subject Device Name:** Comet Anterior Cervical Plate System  
**Trade Name:** Comet Anterior Cervical Plate System  
**Common/Usual Name:** Anterior Cervical Plate System  
**Classification Name:** KWQ – Spinal Intervertebral Body Fixation Orthosis  
21 CFR 888.3060; Class II

**Predicate Device Name:** Comet Anterior Cervical Plate System  
**Trade Name:** Comet Anterior Cervical Plate System  
**Common/Usual Name:** Anterior Cervical Plate System  
**Classification Name:** KWQ – Spinal Intervertebral Body Fixation Orthosis  
21 CFR 888.3060; Class II  
**Premarket Notification:** K082504, SE date January 27, 2009

**Predicate Device Name:** Synthes Anterior Cervical Plate System  
**Trade Name:** Synthes Anterior Cervical Plate System  
**Common/Usual Name:** Anterior Cervical Plate System  
**Classification Name:** KWQ – Spinal Intervertebral Body Fixation Orthosis  
21 CFR 888.3060; Class II  
**Premarket Notification:** K926453, SE date Oct. 12, 1993

### Device Description

Multi-level anterior cervical plates and fixation screws are being added to the previously cleared single-level Comet Anterior Cervical Plate System. The multi-level device additions consist of dual-level plate lengths ranging from 24mm through 48mm and triple-level plate lengths ranging from 36mm through 72mm. Similarly to the single-level devices, the multi-level devices contain barbed plate fixation pins that are permanently affixed to the plate. Plate fixation screws are used to secure the plate to each vertebrae. An additional plate fixation screw diameter is being added in a 3.3mm diameter in 14mm, 16mm and 18mm lengths.

### Indications for Use

The indications for use for the multi-level device additions are the same as those for the previously cleared single-level devices of the Comet Anterior Cervical Plate System. The indications are as follows:

The Comet Anterior Cervical Plate System is intended for anterior fixation, immobilization and stabilization of adjacent cervical vertebral bodies as an adjunct to intervertebral fusion by autogenous and/or allogenic bone graft.

The Comet Anterior Cervical Plate System is indicated for anterior cervical fixation (levels C2-C7) for the following indications: degenerative disc disease (defined as neck pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis, trauma (i.e., fracture or dislocation), spinal stenosis, deformities or curvatures (i.e., scoliosis, kyphosis, and/or lordosis), tumor, pseudoarthrosis, and failed previous fusion.

### **Performance Testing**

Performance data demonstrated that the multi-level device additions to the Comet Anterior Cervical Plate System are substantially equivalent to the predicate single-level devices of the Comet Anterior Cervical Plate System and/or met pre-determined acceptance criteria. The risks associated with use of the new devices were found acceptable when evaluated by FMEA.

Bench tests performed in accordance with FDA's May 2004 *Guidance for Industry and Staff Spinal System 510(k)s* and ASTM F1717-04 *Standard Test Methods for Spinal Implant Constructs in a Vertebroctomy Model* included assessments of static and dynamic axial compression bending and static torsion.

No biocompatibility testing was conducted; all materials used in the manufacture of the Comet Anterior Cervical Plate System device have been previously cleared for similar devices.

### **Conclusion**

The Comet Anterior Cervical Plate System multi-level device additions meet all the pre-determined acceptance criteria of the testing performed to confirm safety and effectiveness; the Comet Anterior Cervical Plate System multi-level device additions are substantially equivalent to the predicate device, the single-level devices of the Comet Anterior Cervical Plate System.



Apollo Spine  
% STD MED, Inc.  
Ms. Christie Santagate  
75 Mill Street  
Stoughton, Massachusetts 02072

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

OCT 23 2009

Re: K090751  
Trade/Device Name: Comet Anterior Cervical Plate System  
Regulation Number: 21 CFR 888.3060  
Regulation Name: Spinal intervertebral body fixation orthosis  
Regulatory Class: II  
Product Code: KWQ  
Dated: October 19, 2009  
Received: October 22, 2009

Dear Ms. Santagate:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Christie Santagate

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic,  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

