

510(k) SUMMARY

K091146

This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of 21 CFR 807.92.

Name, Address, Phone and Fax number of the Applicant

Accuray Incorporated
1310 Chesapeake Terrace
Sunnyvale, California 94089
Ph: (408) 716-4600
Fax: (408) 789-4264

MAY - 1 2009

Contact Person

Anne Schlagenhaft

Date Prepared

April 15, 2009

Device Name

Trade Name: CyberKnife® Robotic Radiosurgery System
Classification Name: Medical charged particle radiotherapy device

Device Description

The CyberKnife Robotic Radiosurgery System is a computer controlled medical system for planning and performing minimally invasive stereotactic radiosurgery and precision radiotherapy using a treatment radiation generator, linear accelerator, manipulator (robot), and a target locating subsystem to accurately deliver radiation to the treatment target. The CyberKnife System uses skull tracking, fiducial tracking, Xsight™ Spine Tracking, Xsight™ Lung Tracking, and Synchrony™ Respiratory Tracking for dynamic positioning and pointing of the linear accelerator.

Intended Use

The CyberKnife Robotic Radiosurgery System is intended to provide treatment planning and image-guided stereotactic radiosurgery and precision radiotherapy for tumors, lesions and conditions anywhere in the body when radiation treatment is indicated.

Substantial Equivalence

The CyberKnife Robotic Radiosurgery System is substantially equivalent to the predicate device. The intended use, principles of operation, technological characteristics and labeling are the same or equivalent.



MAY - 1 2009

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Anne Schlagenhaft
Senior Regulatory Affairs Associate
Accuray, Inc.
1310 Chesapeake Terrace
SUNNYVALE CA 94089

Re: K091146

Trade/Device Name: CyberKnife® Robotic Radiosurgery System
Regulation Number: 21 CFR 892.5050
Regulation Name: Medical charged-particle radiation therapy system
Regulatory Class: II
Product Code: MUJ
Dated: April 17, 2009
Received: April 20, 2009

Dear Dr. Reiss:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter.

21 CFR 876.xxx	(Gastroenterology/Renal/Urology)	(240) 276-0115
21 CFR 884.xxx	(Obstetrics/Gynecology)	(240) 276-0115
21 CFR 892.xxx	(Radiology)	(240) 276-0120
Other		(240) 276-0100

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please contact the CDRH/Office of Surveillance and Biometrics/Division of Postmarket Surveillance at 240-276-3464. For more information regarding the reporting of adverse events, please go to <http://www.fda.gov/cdrh/mdr/>.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Janine M. Morris
Acting Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

K091146

Device Name: CyberKnife® Robotic Radiosurgery System

Indications For Use:

The CyberKnife® Robotic Radiosurgery System is indicated for treatment planning and image guided stereotactic radiosurgery and precision radiotherapy for lesions, tumors and conditions anywhere in the body when radiation treatment is indicated.

Prescription Use
 (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
 (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Ronald M. Zhang
(Division Sign-Off)
Division of Reproductive, Abdominal and
Radiological Devices
510(k) Number K091146

Page 1 of _____