

510(K) Summary

This is a 510(K) summary in accordance with CFR807.82(c).

K091981

A. Submitter Information:

Submitter: Lightmed Corporation
Address: NO.1-1, Lane1, Pao-An St. Sec. 3,
Shulin City, Taipei Hsien 23861, Taiwan
Owner/Operator Number: Mr. Gary Lee, President / 9040850
Contact person: Anita Chen, Regulatory Affair
TEL: +886-2-2688-1726
FAX: +886-2-2676-4920

SEP 28 2009

B. Device Name:

Product Name: LightLas 577 Medical Optical Pumped Semiconductor Laser

Common name: Ophthalmic Laser, Surgical Laser

Classification name:

86 HQF, Laser, Ophthalmic

Regulation Number:

21 CFR 886.4390, Ophthalmic Laser

Regulatory Class: II

Performance standards: 21 CFR 1040.10 & 1040.11

C. Predicate Device Names:

- (4) Device Names: LightLas 561 (K063297)
- (5) Device Names: Novus Varia (K022181)
- (6) Device Names: Family of IRIDEX IQ Laser System(K071687)

D. Device Description:

The LightLas 577 Laser beam has a wavelength of 577nm, which is in the visible spectrum and is a Yellow light. A red aiming beam is used to position the treatment Yellow beam prior to delivery.

The LightLas 577 is a solid state, Optical Pumped Semiconductor surgical laser. It's an instrument used in the photothermolysis (photocoagulation) of soft tissue at an emission wavelength of 577nm.

LightLas LaserLink : Laser delivery adapter that may be coupled to each of the above Selecta models and connected to a currently cleared LightLas 577 retinal photocoagulator to allow use of the slit lamp to deliver 577 nm continuous wave laser energy for retinal photocoagulation.

Compatible delivery devices include: slit lamps, slit lamp adapters/ attachments, laser indirect ophthalmoscopes(LIO) and endoprobe.

The intended use has not changed from the predicate devices.

E. Intended Use:

The LightLas 577 Medical Optical Pumped Semiconductor Laser is intended for use in ophthalmic surgical procedures. A complete list is contained in the Indications for Use Statement.

Ophthalmology:

- Retinal Photocoagulation
- Pan Retinal Photocoagulation
- Endophotocoagulation
- Macular Photocoagulation
- Laser Trabeculoplasty (to treat Open Angle Glaucoma)

The intended use has not changed from the predicate devices.

F. Technological Characteristics summary & Substantial Equivalence

The LightLas 577 Medical Optical Pumped Semiconductor Laser has the same indications for use as the Irides-IQ 577 (K071687). They have similar functional elements such as treatment wavelengths, pulse rates, treatment power, spot size and cooling system. Control systems such as the door interlock, and the safety systems and displays are constantly monitored in these systems for user intervention during a

procedure or maintenance.

G. Performance Data Summary:

The appropriate testing including safety, performance and functional testing to determine substantial equivalence of the LightLas 577 Medical Optical Pumped Semiconductor Laser System.

H. Conclusion

The LightLas 577 Medical Optical Pumped Semiconductor Laser is substantially equivalent to predicate devices currently legally marketed for the indication of retinal photocoagulation, laser trabeculoplasty, the treatment of vascular and pigmented skin lesions, and other laser treatments.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration
10903 New Hampshire Avenue
Document Mail Center - WO66-G609
Silver Spring, MD 20993-0002

Lightmed Corporation
% Ms. Anita Chen
No. 1-1, Lane 1, Pao-An St. Sec. 3
Shulin City, Taipei 23861
Taiwan

SEP 28 2009

Re: K091981

Trade/Device Name: LightLas 577 Medical Optical Pumped Semiconductor Laser
Regulation Number: 21 CFR 886.4390
Regulation Name: Ophthalmic laser
Regulatory Class: Class II
Product Code: HQF
Dated: June 4, 2009
Received: July 1, 2009

Dear Ms. Chen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 - Ms. Anita Chen

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/cdrh/mdr/> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,


FOR Mark N. Melkerson
Director
Division of Surgical, Orthopedic
and Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use Statement as Requested by FDA

510(K) Number (if Known): K091981

Product / Trade Name: LightLas 577 Medical Optical Pumped Semiconductor Laser

Indications for Use:

Ophthalmology:

- Retinal Photocoagulation
- Pan Retinal Photocoagulation
- Endophotocoagulation
- Macular Treatments
- Laser Trabeculoplasty

(7) The intended use has not changed from the predicate devices (K071687, K063297, K022181)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDR, Office of Device Evaluation (ODE)

Prescription Use: X OR Over-The-Counter Use: _____
(Per 21 CFR 801.109)

Neil R. P. Odeh
(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

510(k) Number K091981