

510(k) Summary

MAY 12 2010

Introduction

This summary is intended to comply with requirements of the SMDA and 21CFR§807.92. FDA may make this summary available to the public within 30 days following a finding of substantial equivalence.

510(k) Applicant

Applied Physiology Pty Ltd
119 Willoughby Rd., Crows Nest
2065 Sydney, NSW, Australia
Tel: +61 (2) 9956 3838 / Fax: +61 (2) 9439 2157

510(k) Correspondent

Robert N. Clark, President and Senior Consultant
Medical Device Regulatory Advisors, Inc.
13605 West 7th Ave., Golden, CO 80401 USA
Tel: 303-463-0900 / Fax: 303-558-3833

Date Prepared

April 30, 2010

Trade Name of Device

Navigator™

Common Name of Device

Computer, diagnostic, pre-programmed, single-function

Classification Name

Single function pre-programmed diagnostic computer

510(k) Classification

Class II 21 CFR 870.1435

Predicate Devices

The Navigator is substantially equivalent to the following predicate device(s):

- *K072735 Pulsion PiCCO₂*
- *K023960 LiDCOplus Hemodynamic Monitor*

Device Description

The Navigator is a physiologically integrated system designed to assist clinicians in the management of the systemic circulatory state in the critically ill patient.

Indications for Use:

The Navigator™ is indicated for the acquisition, processing, and display of hemodynamic parameters, in order to assist the clinician in achieving and maintaining a prescribed target hemodynamic stability.

Mean Systemic Filling Pressure (P_{ms}) and Heart Efficiency (E_h) parameters are calculated by the Navigator from Cardiac Output (CO), Mean Arterial Pressure (MAP) and Right Atrial Pressure (RAP) data, obtained directly from patient monitoring equipment, and other data obtained by manual clinician entry. The Navigator system provides a visual indication of the patient's circulatory status in relation to predetermined goals.

Intended Use

The Navigator™ is valuable in any environment where resuscitation, stabilization and optimization of hemodynamic and oxygen metabolism is required. It provides clinicians with monitoring and support information that assists with management of the circulatory state of critically ill patients. The device provides clinicians with a graphical display of monitoring and support information as a visual aid in determining a patient's circulatory state. Navigator™ is used in conjunction with standard monitoring in the following environments including:

- *Intensive care*
- *OR/amb. anesthesia*
- *High dependency and step-down units*
- *Emergency departments*
- *Coronary care*
- *Acute care*

This includes a broad range of patients (between the ages of 15-90 and weights of 40-150 kilograms) with unstable circulations presenting to the intensive care unit (ICU) or critical care units and those undergoing or suffering from:

- *Pre-operative post open heart surgery.*
- *Pre-operative post major surgery.*
- *Sedated patients.*
- *Drug overdose*
- *Septic shock.*
- *Renal failure*
- *Major burns.*
- *Major trauma.*
- *Hypovolemic shock*
- *Cardiogenic shock*

The device is applicable to critically ill patients requiring circulatory support in whom MAP, RAP and CO are being monitored regularly.

Risk Management

This device has been designed to either completely eliminate or mitigate known health hazards associated with the use of the device. Health hazard risk reduction has been accomplished by rigorous application of a risk management program according to standard ISO 14971:2007.

Non-Clinical Testing

Non-clinical testing was performed in order to validate the design against the company's specific design requirements, and to assure conformance with the following voluntary standards:

IEC 60601-1: Medical Electrical Equipment – Part 1: General Requirements for Safety (1988), including Amendment 1 (1991) & Amendment 2 (1995).

BS EN 60601-1-2:2002, incorporating amendment A1:2006. This standard is identical with EMC standard IEC 60601-1-2:2001, incorporating amendment 1:2004: Medical electrical equipment – Part 1-2: General requirements for safety - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests (Edition 2:2001 with Amendment 1:2004).

Rigorous software verification and validation testing was successfully performed by an independent software testing service. Testing included verification and validation to approved Navigator System Requirements, Operational Requirements, and Data Processing Requirements.

The combination of Electrical Safety testing, EMC testing, and Software Verification and Validation provide a high level of confidence that both the Navigator hardware and software are at least as safe and effective as the predicate devices.

Clinical Testing

Applied Physiology Pty Ltd conducted multi-center, open, randomized, controlled human clinical trials in order to collect safety and efficacy data supporting substantial equivalence.

Methodology:

Following surgery, eligible patients were to be randomized on admission to ICU to receive care guided by Navigator™ or conventional care while CO was being monitored. All patients were to be connected to the Navigator™. The screen of the patients in the control group was to have the top section blank, the right hand side was to display actual values of MAP, CO and Right Atrial Pressure (RAP) as slaved from the bedside monitor, along with the patient's screening and randomization number and initials. The arm of the study to which the patient was randomized (control or Navigator™), was also to be shown.

Number of subjects (planned and analyzed):

Sufficient patients were to be consented to allow 100 patients to complete the study (50 in each treatment arm). A total of 112 patients were enrolled into the study and formed the intent to treat (ITT) population. Of these, 107 patients completed the study as planned. A total of 105 patients formed the modified intent to treat (MITT) population (57 patients in the Navigator™ arm and 48 patients in the control arm) these being the patients connected to Navigator™ with three hours or more of data.

Efficacy Results:

The primary endpoint was the average distance (mean standardized) to the central point of the cardiovascular zone over the period the patient was connected to the Navigator™.

The study demonstrated non-inferiority to Navigator over conventional care. There was no statistically significant difference between the two treatment arms with regards to the mean standardized distance to the central point of the cardiovascular zone.

Safety Results:

All 112 patients enrolled in the study were included in the safety evaluation. Adverse events were collected from the time that the patient was connected to Navigator™ until the follow up

visit. Serious adverse events were collected from the time the patient was randomized until the follow-up visit.

Four hundred and thirteen (413) adverse events involving 99 patients (52 (88.1%) patients in the Navigator™ arm and 47 (88.7%) patients on the conventional care arm) were reported during the study.

The most frequently reported adverse events were anaemia, hyperglycaemia, decreased haemoglobin, hypotension and abnormal blood glucose. None of the adverse events experienced by patients on the Navigator™ arm or the conventional care arm of the study were considered to have a causal relationship to the device.

One patient was withdrawn from the study due to an adverse event (cardiac tamponade) considered to be serious but not related to the device). Two patients died during the study due to events considered to be not related to the device.

Ten (10) device incidents occurred during the study, but none were associated with any clinical consequences. All occurrences were recognizable by the operator or physician, and necessary corrections were made in the production device. These incidents fell into the following categories:

- *Four occurrences of Navigator screen buttons or settings that did not function as intended. These occurrences were recognizable by the operator.*
- *Four occurrences of data communication problems with peripheral equipment or failure to display data on peripheral equipment.*
- *One occurrence where a Swan-Ganz catheter appeared to be incorrectly positioned on x-ray and the operator was confident of MAP, CO and CI readings.*
- *One occurrence where a discrepancy was reported with MAP display when an aortic balloon pump was inflated.*

Thirty-five (35) serious adverse events involving 24 patients (12 on the Navigator™ arm and 12 on the conventional care arm) were reported during the study. No serious adverse events were considered to be device related.

Review of adverse events per treatment group, suggested no trends in the presentation (type, frequency or severity) of adverse events per study arm.

Conclusion

The study demonstrated non-inferiority to Navigator over conventional care. There was no statistically significant difference between the two treatment arms with regards to the mean distance to the central point of the cardiovascular zone.

The incidence of adverse events and serious adverse events was comparable across the two study arms, with no trends in the type, frequency or severity of event. No adverse or serious adverse events were considered to be related to Navigator™.

Device incidents which occurred during the study were not associated with any clinical consequences and did not impact on patient safety.

Data from this study indicate that the use of Navigator™ is safe when compared to conventional ICU care.

Substantial Equivalence

Applied Medical Technology Pty Ltd believes that the Navigator is safe and effective when used as instructed by knowledgeable and trained personnel. Navigator performs at least as safely and effectively as conventional care using predicate devices, and is therefore substantially equivalent to the marketed predicate device(s).



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

MAY 12 2010

Applied Physiology Pty., Ltd.

~~c/o Mr. Robert N. Clark~~

President and Senior Consultant
Medical Device Regulatory Advisors, Inc.
4251 Kipling Street, Suite 565
Wheat Ridge, CO 80033-2899

Re: K092219

Device Name: Navigator™ Clinical Guidance System

Regulation Number: 21 CFR 870.1435

Regulation Name: Single function pre-programmed diagnostic computer

Regulatory Class: Class II (Two)

Product Code: DXG

Dated: April 30, 2010

Received: May 7, 2010

Dear Mr. Clark

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K092219

Device Name: Navigator™

Indications for Use:

The Navigator™ is indicated for the acquisition, processing, and display of hemodynamic parameters, in order to assist the clinician in achieving and maintaining a prescribed target hemodynamic stability.

Mean Systemic Filling Pressure (P_{ms}) and Heart Efficiency (E_h) parameters are calculated by the Navigator from Cardiac Output (CO), Mean Arterial Pressure (MAP) and Right Atrial Pressure (RAP) data, obtained directly from patient monitoring equipment, and other data obtained by manual clinician entry. The Navigator system provides a visual indication of the patient's circulatory status in relation to predetermined goals.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Cardiovascular Devices
510(k) Number K092219