

K092541

510(k) Summary

Interface Bone Void Filler

**510(k) Summary**

NOV - 9 2010

**BioStructures, LLC**

**Interface Bone Void Filler**

October 1, 2010

**ADMINISTRATIVE INFORMATION**

Manufacturer Name: BioStructures, LLC  
3700 Campus Drive, Suite 204  
Newport Beach, CA 92660  
Telephone: +1 (949) 553-1717  
Fax: +1 (949) 553-0407

Official Contact: Russell Cook, CEO

Representative/Consultant: Kevin A. Thomas, PhD  
Floyd G. Larson  
PaxMed International, LLC  
11234 El Camino Real, Suite 200  
San Diego, CA 92130  
Telephone: +1 (858) 792-1235  
Fax: +1 (858) 792-1236  
Email: kthomas@paxmed.com  
flarson@paxmed.com

**DEVICE NAME AND CLASSIFICATION**

Trade/Proprietary Name: Interface Bone Void Filler

Common Name: Filler, bone void, calcium compound  
Classification Regulations: Resorbable calcium salt bone void filler device  
21 CFR 888.3045  
Class II

Product Code: MQV  
Classification Panel: Orthopaedic and Rehabilitation Devices Panel  
Reviewing Branch: Restorative Devices Branch

## INTENDED USE

Interface Bone Void Filler is indicated for bony voids or gaps that are not intrinsic to the stability of the bony structures. Interface Bone Void Filler is indicated to be gently packed into bony voids or gaps of the skeletal system (the extremities and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The product provides a bone void filler that resorbs and is replaced with bone during the healing process.

## DEVICE DESCRIPTION

Interface Bone Void Filler is a synthetic bioactive bone graft for use in the repair of osseous defects. It is supplied as irregular synthetic granules of bioactive glass, sized from 200 microns to 420 microns. When implanted in living tissue, the material undergoes a time dependent surface modification. The surface reaction results in the formation of a calcium phosphate layer, which is equivalent in composition and structure to the hydroxyapatite found in bone mineral. The biological apatite layer of the granules provides an osteoconductive scaffold for the generation of new osseous tissue. New bone infiltrates around the granules allowing the repair of the defect as the granules are absorbed. The elemental composition of Interface Bone Void Filler granules is Si, Ca, Na, and P. Interface Bone Void Filler conforms to ASTM specification F1538 for 45S5 bioactive glass.

## EQUIVALENCE TO MARKETED DEVICE

BioStructures, LLC demonstrated that for the purposes of FDA's regulation of medical devices, Interface Bone Void Filler is substantially equivalent in indications and design principles to NovaBone Resorbable Bone Graft Substitute from NovaBone Products, LLC, cleared under K021336.

The subject device and the predicate device have the same intended use and have the same technological characteristics. The subject device and the predicate device are synthetic granules of bioactive glass with elemental composition of Si, Ca, Na, and P. The subject device and the predicate device have similar particle sizes, and are packaged in similar materials and sterilized using similar methods.

To demonstrate equivalence detailed side-by-side material characterization was performed including chemical composition, physical properties and performance characteristics. Chemical composition was analyzed by scanning electron microscopy with energy dispersive x-ray analysis (SEM/EDXA). Trace elemental analysis was performed by inductively coupled plasma/optical emission spectroscopy (ICP/OES). Crystallinity was analyzed by Fourier Transform Infrared Spectroscopy (FT-IR) and X-ray Diffraction (XRD). Physical properties were evaluated by scanning electron microscopy and particle size was determined by laser diffraction. Dissolution testing was performed by monitoring the concentration of calcium in

media by a calcium-specific electrode in an appropriate dissolution apparatus. The analytical characterization demonstrated equivalent chemical composition, physical properties and performance characteristics for the Interface Bone Void Filler and the NovaBone Resorbable Bone Graft Substitute devices.

The radiographic, morphometric and histologic performance of the subject Interface Bone Void Filler device were compared to that of the predicate NovaBone device in an animal model. The results of the study demonstrated that the performance of the subject Interface Bone Void Filler device was equivalent to that of the predicate NovaBone device.

Overall, Interface Bone Void Filler has the following similarities to the predicate devices:

- has the same intended use,
- uses the same operating principle,
- incorporates the same basic design,
- incorporates the same or very similar materials, and
- has similar packaging and is sterilized using the same processes.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

BioStructures, LLS  
% Mr. Russell Cook, CEO  
3700 Campus, Drive Suite 204  
Newport Beach, California 92660

Re: K092541

NOV - 9 2010

Trade/Device Name: Interface Bone Void Filler  
Regulation Number: 21 CFR 888.3045  
Regulation Name: Resorbable calcium salt bone void filler device  
Regulatory Class: Class II  
Product Code: MQV  
Dated: October 1, 2010  
Received: October 4, 2010

Dear Mr. Cook:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson

Director

Division of Surgical, Orthopedic,  
and Restorative Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

**Indications for Use**

510(k) Number: K092541

NOV - 9 2010

Device Name: Interface Bone Void Filler

Indications for Use:


Interface Bone Void Filler is indicated for bony voids or gaps that are not intrinsic to the stability of the bony structures. Interface Bone Void Filler is indicated to be gently packed into bony voids or gaps of the skeletal system (the extremities and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The product provides a bone void filler that resorbs and is replaced with bone during the healing process.

Prescription Use  X  AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of  1

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number  K092541