

**Section 5: 510(k) Summary**

OCT 22 2009

The following information is provided as required by 21 CFR § 807.87 for ENKO Ltd.'s premarket notification. In response to the Safe Medical Devices Act of 1990, the following is a summary of the safety and effectiveness information upon which the substantial equivalence determination is based.

The safety and effectiveness of the DEBP is based upon a determination of the substantial equivalence as well as the safety and effectiveness of its predicate device(s).

**Applicant:** ENKO Ltd.  
10006 sokak, No:64  
A.O.S.B. Cigli - Izmir 35620  
TURKIYE  
Tel: 0090 232 3767806  
Fax: 0090 232 3767792  
[www.enkoelektronik.com](http://www.enkoelektronik.com)  
Attn: Sinan Kazazoglu

**Contact:** Calley Herzog  
Biologics Consulting Group, Inc.  
13417 Quivas St.  
Westminster, CO 80234  
Ph. 720-883-3633  
Fax. 720-293-0014

**Date of Submission:** 07/01/2009

**Proprietary Name:** DEBP

**Common Name:** Powered Breast Pump

**Regulatory Class:** Class II

**Product Codes:** HGX

**Predicate Device(s):** Medela swing Breastpump, By Medela Inc. (K053052) & Hollister Inc. Expresse and Premier Powered Breast Pump (Lactaline Personal) (K973501)

**Device Description:** The Double Electric Breast Pump (DEBP) is a powered breast pump. The pumping can be performed on one breast or on both breasts at the same time. The Lansinoh

DEBP can be powered by 6 AA batteries or an AC adaptor provided with the pump. The pumping system consists of a diaphragm-type vacuum pump which is driven by a microcontroller controlled DC electric motor. The user interface consists of a front panel keypad and LCD display. The user is able to control cycle speed and vacuum level. The Lansinoh DEBP is capable of providing vacuum levels from 50 to 250 mmHg, with cycling rates up to 1.85 cycles per second.

**Intended Use:** The DEBP is intended to express and collect the mother's milk of a nursing woman for the purpose of feeding the collected milk to a baby.

**Performance Testing:** Suction curves are provided to illustrate the performance of the DEBP. Additionally a backflow test was conducted to ensure satisfactory performance of the pump in the unlikely event that milk were to backflow into the pump unit.

The DEBP will be tested to meet:

- IEC 60601-1, "Medical Electrical Equipment, General Requirements for Safety"
- IEC 60601-1-2:2007, "Medical electrical equipment. General requirements for basic safety and essential performance. Collateral standard. Electromagnetic compatibility."

**Substantial Equivalence:** The Lansinoh DEBP is substantially equivalent to the predicate devices in intended use, technological characteristics and device design. The table below provides a comparison of the DEBP to the predicate devices.

	<b>New Device</b>	<b>Predicate Device</b>	<b>Predicate Device</b>
<b>Manufacturer</b>	ENKO Ltd.	Ameda/Hollister	Medela
<b>Device Name</b>	<b>DEBP</b>	<b>Expresse &amp; Premeier Breast Pumps</b>	<b>Medela Swing</b>
<b>510(k) #</b>	not yet assigned	K973501	K053052
<b>Intended Use</b>	The DEBP is intended to express and collect the mother's milk of a nursing woman for the purpose of feeding the collected milk to a baby.	The Lactaline Personal Breast Pumps are intended to express and collect the mother's milk of a nursing woman for the purpose of feeding the collected milk to a baby.	The Swing Breast Pump is a powered breast pump to be used by lactating women to express and collect milk from their breasts.
<b>Pumping Suction</b>	50 - 250 mmHg	<100 - 265 mbar (<75 - 199 mmHg)	0 - 250 mmHg

<b>Stimulation/Let Down Phase</b>			
<b>Suction Levels</b>	50 – 150 mmHg	n/a	0 - 250 mmHg
<b>Cycles per Second</b>	1.85 (fixed)	n/a	up to 2.17
<b>Expression Phase</b>			
<b>Suction Levels</b>	50 – 250 mmHg	<100 - 360 mbar	0 - 250 mmHg
<b>Cycles per Second</b>	0.51 – 1.0	0.5 – 1.0	up to 2.17
<b>Suction Settings</b>	8	3	11
<b>Power Supply</b>	a) 6 AA alkaline batteries b) AC Adapter	a) rechargeable NiCd batteries b) 6 AA alkaline batteries c) AC Adapter d) 12 V adapter for use in motor vehicle	AC Adapter
<b>Pumping Option</b>	Single or Double	Single or Double	Single
<b>Back Flow Protection</b>	Yes	Yes	No
<b>Let Down Function</b>	Yes	No	Yes
<b>Cycling/Suction Control Mechanism</b>	Microcontroller	Microprocessor	Microprocessor



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room – WO66-G609  
Silver Spring, MD 20993-0002

OCT 22 2009

ENKO Ltd.  
c/o Mr. Ned Devine  
Sr. Staff Engineer  
Underwriters Laboratories, Inc.  
333 Pfingsten Road  
NORTHBROOK IL 60062-2096

Re: K092783  
Trade/Device Name: DEBP  
Regulation Number: 21 CFR §884.5160  
Regulation Name: Powered breast pump  
Regulatory Class: II  
Product Code: HGX  
Dated: October 5, 2009  
Received: October 7, 2009

Dear Mr. Devine:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris". The signature is fluid and cursive, with the first name being the most prominent.

Janine M. Morris  
Acting Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

### Section 4: Indications for Use Statement

510(k) Number: ~~To be assigned~~ K092783

Device Name: DEBP


Indications for Use:

The DEBP is intended to express and collect the mother's milk of a nursing woman for the purpose of feeding the collected milk to a baby.

Prescription Use \_\_\_\_\_ AND/OR Over-The-Counter Use  X   
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal, and  
Radiological Devices  
510(k) Number  K092783