

SOMANETICS VITAL SYNC SYSTEM 510(K) PREMARKET NOTIFICATION

510(k) Summary

Date of Submission: April 21, 2010 APR 29 2010

Device Trade Name: Vital Sync™ System

Device Common Name: Cardiac monitor (without alarms)

Device Classification Name: Monitor, Physiological, Patient (without arrhythmia detection or alarms)

Regulation and Code: 21 CFR 870.2300, Product Code MWI

Submitted by: Somanetics Corporation
2600 Troy Center Drive
Troy, MI 48084
Phone: 248-244-1400
Fax: 248-244-0978

Contact Person: Ronald A. Widman
Vice President, Medical Affairs
248-244-1449

Predicate Device: CentraView, ICU DataSystems, K033283

Device Description: The Vital Sync™ System is a bedside data management system that receives historic digital data produced by primary external devices through device specific cables, accepts manual data entry, and displays and stores this information for review and archiving by healthcare professionals.

Accessories: Serial Port Concentrator, 4 Port, Model 5000-SPC4
Serial Port Concentrator, 8 Port, Model 5000-SPC8
USB Cable, 5 meters, Model 5000-USB5
USB Cable, 1 meter, Model 5000-USB1
USB Extension Cable, 5 m, Model 5000-USBE5
Keyboard, Model 5000-KE
Mouse, Model 5000-MO
Roll Stand with Single Mount, Model 5000-RS
Roll Stand with Dual Mount, Model 5000-RS2
Desk Stand, Model 5000-DS
Serial Port Cable, 14 ft (4.27 m), Model SPAC-14
Serial Port Cable, 7 ft (2.13 m), Model SPAC-07
Serial Port Adapters, Models MA022-MA033
Printer w/ Isolation Transformer, Model 5000-PIK

SOMANETICS VITAL SYNC SYSTEM 510(K) PREMARKET NOTIFICATION

Advanced Clinical Tools Pkg, Model 5000
Advanced Clinical Tools Pkg, Model 5000-ACT

Intended Use: The Vital Sync™ System is intended for display and recording of multiple physiological parameters of adult, pediatric and neonatal patients. It is not intended for alarm notification, nor is it intended to control any of the independent bedside devices it is connected to. A listing of supported devices and displayed parameters is attached.

Technological Characteristics:..... Intended use, operating principle, performance claims, and technological characteristics of the device, including design, materials, chemical composition and energy source are identical to the predicate device.

Performance Data: Bench testing and verification and validation activities were performed to establish the performance, reliability and functionality of the Vital Sync System. Clinical testing was not required to establish substantial equivalence. Hazard analysis established the safety and system level testing and validation needed to demonstrate substantial equivalence. Testing modes included error handling, system faults, power cycling, consistency, recall, user interface, performance and maintenance procedures. Features tested include automated data capture, manual data entry, data formatting, report generation, trending, data backup and security. Additionally, the communication interface for each of the supported devices was stress-tested to ensure safety and compatibility with the Vital Sync System. Pass/fail criteria were established based on the published specifications of both the predicate and the current device. The results demonstrate substantial equivalence with the predicate device.

Conclusion Drawn from the Testing:..... The conclusion drawn is that the device and the revised indications for use are substantially equivalent to the predicate device and do not raise new questions of safety and effectiveness.

SOMANETICS VITAL SYNC SYSTEM 510(K) PREMARKET NOTIFICATION

Vital Sync System Displayed Parameters

Vital Sync Displayed Parameters

Amplitude Integrated EEG Left	End Diastolic Volume
Amplitude Integrated EEG Right	End Diastolic Volume Index
Anesthetic Agent	Stat End Diastolic Volume Index
Air Temperature	Stat End Diastolic Volume
Air Temperature Setting	Exhaled Minute Volume
Airway Temperature	End Inspiratory Pressure
Arterial Base Excess	End Systolic Volume
Arterial Bicarbonate	End Systolic Volume Index
Arterial pH	End Tidal CO2
Arterial Temperature	Environement Temp
Average Heart Rate	Esophageal Temperature
Axillary Temp	Exhaled Tidal Volume
Bi-level Positive Airway Pressure	Exhalation Time
BIS	Expired Positive Airway Pressure
Bladder Temperature	Inspired Fraction of Oxygen
Blood Temperature/Pulm. Artery Temperature	Infusate Temp
Body Surface Area	Flow Rate 1
Bolus Cardiac Output	Flow Rate 2
Bolus Cardiac Index	Flow Rate 3
Brain PO2	Flow Velocity
Brain Temp	Gastric pCO2
Calculated SO2	Heart Rate
Cardiac Index	Heater Output Percent
Cardiac Output	Heater Output Percent Setting
Stat Cardiac Index	Hematocrit
Stat Cardiac Output	Hemoglobin
Cardioplegia Line Pressure	High Inspired Pressure Setting
Central Venous Pressure	Humidity
Cerebral Perfusion Pressure	Humidity Setting
Cerebral Blood Flow	Variation of Contractility Index
Ch1 rSO2	Inspired:Expired Ratio
Ch2 rSO2	Inspired CO2
Ch3 rSO2	Inspired O2 Setting
Ch4 rSO2	Insp Pos Air Pressure
Contractility Index	Inspiratory Pressure
Contnuous Cardiac Index	Inspiratory Resistance
Continuous Cardiac Output	Inspiratory Tidal Volume
Control Temp	Inspiratory Time
Core Temperature	Intracranial Pressure
Coronary Sinus Pressure	Diastolic Arterial Blood Pressure
Continuous Positive Airway Pressure	Mean Arterial Blood Pressure
Delta Pressure	Systolic Arterial Blood Pressure
Dynamic Compliance	% Leak in Tidal Volume (Insp/Exp)
Ejection Fraction	Left Atrial Pressure
Emboli 1	Left Cardiac Work
Emboli 2	Left Cardiac Work Index
Emboli 3	Left Stroke Work

SOMANETICS VITAL SYNC SYSTEM 510(K) PREMARKET NOTIFICATION**Vital Sync System Displayed Parameters (cont'd)****Vital Sync System Displayed Parameters
(cont'd)**

Left Stroke Work Index	Pump Flow
Left Ventricular Ejection Time	Venous pCO2
Line Pressure	Venous pO2
Mean Airway pressure	Rate Pressure Product
Mattress Temperature	Rectal Temperature
Mattress Temperature Setting	Right Ventricular Ejection fraction
Respiration Rate	Stat Right Ventricular Ejection fraction
Ventilation Mode	Arterial SO2
Myocardial Temperature	Set Point Temperature
Nasopharyngeal Temperature	Skin Temperature
Nitric Oxide	Skin Temperature Setting
Nitric Dioxide	Pulse SO2 1
Diastolic Cuff Blood Pressure	Pulse SO2 2
Mean Cuff Blood Pressure	Spectral Edge Frequency
Systolic Cuff Blood Pressure	Spontaneous Respiration rate
Oxygen Consumption	Static Compliance
Oxygen Extraction Index	ST Interval
Oral Temp	Stat Stroke Volume Index
Arterial pCO2	Stat Stroke Volume
Arterial pO2	Venous SO2
Peak Flow	Oxygenator Sweep
Plateau Time	Systolic Time Ratio
Positive Pressure Duration	Thoracic Fluid Index
Potassium	Tidal Volume Setting
Pressure Control	Transcutaneous pCO2
Pressure Limit	Transcutaneous pO2
Pressure Sensitivity	Tympanic Temperature
Pressure Support	Venous pH
Pulmonary Capillary Wedge Pressure	Venous Temperature
Pulmonary Arterial Diastolic pressure	Volume Flow 1
Pulmonary Arterial Mean Pressure	Volume Flow 2
Pulmonary Arterial Systolic Pressure	Volume Flow 3
Pulse Amplitude	Water Temperature
Pulse Rate	Weight

SOMANETICS VITAL SYNC SYSTEM 510(K) PREMARKET NOTIFICATION**Vital Sync System Supported Devices (Listed by Company)****Atom Medical**

Infa Warmer V505

Baxter Healthcare

AS 50

Colleague IP

FloGuard 6201

FloGuard 6301

Bird

VIP Gold/Sterling Ventilator

Cardiotronic/Osypka

Aesculon Noninvasive Cardiac Output

Cincinnati Sub-zero

Blanketrol II

Blanketrol III

**Covidien/Tyco/Malinckrodt/Nellcor/Puritan
Bennett**

PB 840 Ventilator

N200 Pulse Oximeter

N295 Pulse Oximeter

N395 Pulse Oximeter

N595 Pulse Oximeter

InfantStar 500 Ventilator

InfantStar 950 Ventilator

Datex/Ohmeda

S/5 Patient Monitor

Draeger

Babylog 8000 Ventilator

Babylog 8000SC Ventilator

Evita Ventilator

Evita 2 Ventilator

Evita 4 Ventilator

Infinity Patient Monitor

SC7000

SC8000

SC9000XL

Edwards Life Sciences

Vigilance Hemodynamic Monitor

Vigilance II Hemodynamic Monitor

Vigileo Hemodynamic Monitor

General Electric

Dash 2000 Patient Monitor

Dash 3000 Patient Monitor

Dash 4000 Patient Monitor

Solar 8000i Patient Monitor

Solar 8000M Patient Monitor

Maquet/Siemens

Servo 300 Ventilator

Servo i Ventilator

Maquet Perfusion Pump System

Masimo

SET Radical-7 Pulse Oximeter

SET Radical-9 Pulse Oximeter

Mennen Medical

Horizon 2000 Patient Monitor

Phillips/Agilent/Hewlett Packard

CMS 2001 Patient Monitor

V24 Patient Monitor

V26 Patient Monitor

MP5 Patient Monitor

MP20 Patient Monitor

MP40 Patient Monitor

MP50 Patient Monitor

MP60 Patient Monitor

MP70 Patient Monitor

MP80 Patient Monitor

MP90 Patient Monitor

Somanetics Corporation

INVOS 5100B Cerebral/Somatic Oximeter

INVOS 5100C Cerebral/Somatic Oximeter

Sorin Biomedica

Sorin S3 Perfusion Pump System

Sorin C5 Perfusion Pump System

Sorin S5 Perfusion Pump System

Sorin SIII Encore Perfusion Pump System

Viasys

Viasys Avea Ventilator



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room –WO66-G609
Silver Spring, MD 20993-0002

Somanetics Corporation
c/o Mr. Ronald A. Widman
Vice President, Medical Affairs
2600 Troy Center Drive
Troy, MI 48084

APR 29 2010

Re: K093422
Device Name: Vital Sync™ System
Regulation Number: 21 CFR 870.2300
Regulation Name: Patient Physiological Monitor (without arrhythmia detection or alarms)
Regulatory Class: Class II (Two)
Product Code: MWI
Dated: April 21, 2010
Received: April 22, 2010

Dear Mr. Widman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,


For Bram D. Zuckerman, M.D.
Director

Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K093422

Device Name: Somanetics Vital Sync™ System and Accessories

Indications For Use:

The Vital Sync™ System is intended for display and recording of multiple physiological parameters of adult, pediatric and neonatal patients. It is not intended for alarm notification, nor is it intended to control any of the independent bedside devices it is connected to.

Vital Sync System displayed parameters are listed on the following pages.

Prescription Use X
(Part 21 CFR 801 subpart D)

OR

Over-The-Counter Use
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K093422

Vital Sync Displayed Parameters

Amplitude Integrated EEG Left	End Diastolic Volume
Amplitude Integrated EEG Right	End Diastolic Volume Index
Anesthetic Agent	Stat End Diastolic Volume Index
Air Temperature	Stat End Diastolic Volume
Air Temperature Setting	Exhaled Minute Volume
Airway Temperature	End Inspiratory Pressure
Arterial Base Excess	End Sysstolic Volume
Arterial Bicarbonate	End Sysstolic Volume Index
Arterial pH	End Tidal CO2
Arterial Temperature	Environement Temp
Average Heart Rate	Esophageal Temperature
Axillary Temp	Exhaled Tidal Volume
Bi-level Positive Airway Pressure	Exhalation Time
BIS	Expired Positive Airway Pressure
Bladder Temperature	Inspired Fraction of Oxygen
Blood Temperature/Pulm. Artery Temperature	Infusate Temp
Body Surface Area	Flow Rate 1
Bolus Cardiac Output	Flow Rate 2
Bolus Cardiac Index	Flow Rate 3
Brain PO2	Flow Velocity
Brain Temp	Gastric pCO2
Calculated SO2	Heart Rate
Cardiac Index	Heater Output Percent
Cardiac Output	Heater Output Percent Setting
Stat Cardiac Index	Hematocrit
Stat Cardiac Output	Hemoglobin
Cardioplegia Line Pressure	High Inspired Pressure Setting
Central Venous Pressure	Humidity
Cerebral Perfusion Pressure	Humidity Setting
Cerebral Blood Flow	Variation of Contractility Index
Ch1 rSO2	Inspired:Expired Ratio
Ch2 rSO2	Inspired CO2
Ch3 rSO2	Inspired O2 Setting
Ch4 rSO2	Insp Pos Air Pressure
Contractility Index	Inspiratory Pressure
Contuous Cardiac Index	Inspiratory Resistance
Continous Cardiac Output	Inspiratory Tidal Volume
Control Temp	Inspiratory Time
Core Temperature	Intracranial Pressure
Coronary Sinus Pressure	Diastolic Arterial Blood Pressure
Continuous Positive Airway Pressure	Mean Arterial Blood Pressure
Delta Pressure	Systolic Arterial Blood Pressure
Dynamic Compliance	% Leak in Tidal Volume (Insp/Exp)
Ejection Fraction	Left Atrial Pressure
Emboli 1	Left Cardiac Work
Emboli 2	Left Cardiac Work Index
Emboli 3	Left Stroke Work

**Vital Sync System Displayed Parameters
(cont'd)**

Left Stroke Work Index	Pump Flow
Left Ventricular Ejection Time	Venous pCO2
Line Pressure	Venous pO2
Mean Airway pressure	Rate Pressure Product
Mattress Temperature	Rectal Temperature
Mattress Temperature Setting	Right Ventricular Ejection fraction
Respiration Rate	Stat Right Ventricular Ejection fraction
Ventilation Mode	Arterial SO2
Myocardial Temperature	Set Point Temperature
Nasopharyngeal Temperature	Skin Temperature
Nitric Oxide	Skin Temperature Setting
Nitric Dioxide	Pulse SO2 1
Diastolic Cuff Blood Pressure	Pulse SO2 2
Mean Cuff Blood Pressure	Spectral Edge Frequency
Systolic Cuff Blood Pressure	Spontaneous Respiration rate
Oxygen Consumption	Static Compliance
Oxygen Extraction Index	ST Interval
Oral Temp	Stat Stroke Volume Index
Arterial pCO2	Stat Stroke Volume
Arterial pO2	Venous SO2
Peak Flow	Oxygenator Sweep
Plateau Time	Systolic Time Ratio
Positive Pressure Duration	Thoracic Fluid Index
Potassium	Tidal Volume Setting
Pressure Control	Transcutaneous pCO2
Pressure Limit	Transcutaneous pO2
Pressure Sensitivity	Tympanic Temperature
Pressure Support	Venous pH
Pulmonary Capillary Wedge Pressure	Venous Temperature
Pulmonary Arterial Diastolic pressure	Volume Flow 1
Pulmonary Arterial Mean Pressure	Volume Flow 2
Pulmonary Arterial Systolic Pressure	Volume Flow 3
Pulse Amplitude	Water Temperature
Pulse Rate	Weight