

**510(k) Summary**

Page 1 of 2

7-Dec-09

EasyGlide Ltd.  
30 Ha'Ella St.  
Kfar Truman, 73150  
Israel

Tel - 011-972-528565644  
Fax – 011-972-776201003

*K093779*

*pg 1 of 2*

**Official Contact:** Izhak Fabian - CEO

**Proprietary or Trade Name:** ClearPath UGI

**MAR - 9 2010**

**Common/Usual Name:** Gastroscope accessory

**Classification Name/Code:** FDS - endoscope and accessories  
CFR 876.1500

**Device:** ClearPath UGI

**Predicate Devices:** K091305 – Easy Glide – ClearPath  
K060907 – Karl Storz video gastroscope

**Device Description:**

The ClearPath UGI's main purpose is to improve procedure reliability by improving procedural visualization during endoscopic procedures.

The ClearPath UGI utilizes a suction / irrigation head for the purpose of irrigating or cleaning the upper digestive tract.

The ClearPath UGI is composed of two major units:

- The disposable Irrigator and
- the reusable Control cabinet.

The disposable Irrigator is attached to a standard endoscope by means of silicone bands and does not hinder the endoscope's functionality nor affects the procedure sequence in any way.

The Control cabinet supplies the water flow and vacuum control using a peristaltic pump. The standard medical facility vacuum system is the vacuum source.

**Indications for Use:**

The ClearPath-UGI System is intended for irrigating or cleaning the upper digestive tract and evacuating the irrigation fluid, blood and bile in the upper GI tract during endoscopic procedures.

It is for use only by trained medical personnel located in hospitals, clinics, and doctors' offices.

*K093979*

*pg 2 of 2*

**Patient Population:**

Individuals undergoing procedures where a gastroscope is used.

**Environment of Use:**

Hospitals, clinics, and doctors' offices.

**Summary of substantial equivalence:**

We have performed testing to demonstrate that the operating performance of the ClearPath UGI is equivalent to the predicate for flow rate, irrigation and vacuum pressures.

	<b>Proposed ClearPath UGI</b>	<b>ClearPath EasyGlide</b>
<b>510(k) number</b>		K091305
<b>Air / Water pressure specifications</b>	Up to 30 Psi	Up to 30 Psi
<b>Flow rate</b>	Up to 300 ml/min	Up to 300 ml/min
<b>Suction specifications</b>	Standard wall suction Approximately 0.5 Bar	Standard wall suction Approximately 0.5 Bar

The proposed ClearPath UGI is only a modification of the predicate ClearPath in that the disposable tubing set is shorter in length and the irrigation tip has been redesigned. No modification have been made to the controller or software.

The ClearPath UGI is viewed as substantially equivalent to the predicate devices because:

**Indications –**

- Identical to predicate – K060907 – Karl Storz

**Technology –**

- Identical technology used – K091305 – EasyGlide ClearPath

**Materials –**

- The materials in patient contact are identical to the predicate device, ClearPath K091305.

**Environment of Use –**

- Identical to predicate – K091305 – EasyGlide ClearPath

**Performance specifications –**

- Irrigation pressures, flow rate and suction pressures are identical to the predicate – K091305 EasyGlide ClearPath



DEPARTMENT OF HEALTH & HUMAN SERVICES

MAR 9 2010

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Mail Center - WO66-G609  
Silver Spring, MD 20993-0002

EasyGlide Ltd.  
% Mr. Paul Dryden  
President  
ProMedic, Incorporated  
24301 Woodsage Drive  
BONITA SPRINGS FL 34134

Re: K093779  
Trade/Device Name: ClearPath UGI  
Regulation Number: 21 CFR §876.1500  
Regulation Name: Endoscope and accessories  
Regulatory Class: II  
Product Code: FDS  
Dated: December 8, 2009  
Received: December 9, 2009

Dear Mr. Dryden:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related

adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris", is written over the typed name and title.

Janine M. Morris  
Acting Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

**Indications for Use Statement**

**510(k) Number:** K093779 (To be assigned)

**Device Name:** ClearPath UGI

**Indications for Use:**

The ClearPath-UGI System is intended for irrigating or cleaning the upper digestive tract and evacuating the irrigation fluid, blood and bile in the upper GI tract during endoscopic procedures.

It is for use only by trained medical personnel located in hospitals, clinics, and doctors' offices.

**Prescription Use XX**  
(Part 21 CFR 801 Subpart D)

or

**Over-the-counter use** \_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



---

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K093779