510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY DEVICE ONLY TEMPLATE

A. 510(k) Number:

K030328

B. Analyte:

D-Dimer

C. Type of Test:

Quantitative, enzyme-linked fluorescent Immunoassay

D. Applicant:

bioMerieux, Inc.

E. Proprietary and Established Names:

VIDAS D-DIMER NEW (DD2) ASSAY

F. Regulatory Information:

1. Regulation section:

21 CRD 864.7320

2. Classification:

Class II

3. Product Code:

DAP

4. Panel:

81 Hematology

G. Intended Use:

1. Indication(s) for use:

The VIDAS®D-Dimer New is an automated, quantitative test for use on the VIDAS analyzer for the immunoenzymatic determination of cross-linked fibrin degradation products (FbDP) containing the D-dimer domain in citrated human plasma using the Enzyme Linked Fluorescent Assay (ELFA) technique.

2. Special condition for use statement(s):

The VIDAS®D-Dimer New is indicated for use in conjunction with a clinical pretest probability (PTP) assessment model to exclude deep vein thrombosis (DVT) in outpatients suspected of DVT.

3. Special instrument Requirements:

The VIDAS®D-Dimer New is intended for use on the VIDAS analyzer

H. Device Description:

The VIDAS®D-Dimer New (DD2) is an automated, quantitative test for d-dimer, intended for use on the VIDAS analyzer (K891385). Fibrin degradation products (FbDP) in human plasma are determined using the enzyme-linked fluorescent immunoassay (ELFA) technique. The instrument controls all assay steps and assay temperatures. A pipette tip like disposable device, the Solid Phase Receptacle (SPR), serves as the solid phase as well as a pipettor for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed DD2 Reagent Strips.

I. Substantial Equivalence Information:

- 1. <u>Predicate device name(s):</u> VIDAS D-Dimer (DD) New Assay
- 2. Predicate K number(s): K020810
- 3. Comparison with predicate:

Similarities								
Item	Device	Predicate						
Sample requirements	Citrated plasma	Same						
Differences								
Item	Device	Predicate						
Indications for use	To exclude DVT in	Aid in diagnosis of DVT						
	conjunction with a PTP	and PE						

J. Standard/Guidance Document Referenced (if applicable):

K. Test Principle:

enzyme-linked fluorescent immunoassay (ELFA)

L. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
 - a. Precision/Reproducibility:

			Within-run	Total Precision
			Precision	
Plasma	N	Conc.	CV	CV
		(ng FEU/ml)	(%)	(%)
Level 1	80	264	5.0	5.7
Level 2	80	549	3.9	5.8
Level 3	80	7283	5.3	7.1

- b. Linearity/assay reportable range:
 - 45 10,000 ng FEU/ml
- c. Traceability (controls, calibrators, or method):
- *d.* Detection limit:45 ng FEU/ml
- e. Analytical specificity:
- f. Assay cut-off: 500 ng FEU/ml
- 2. Comparison studies:
 - a. Method comparison with predicate device:
 - b. Matrix comparison:
- 3. Clinical studies:
 - a. Clinical sensitivity: 100% (95% CI, 95.0-100)
 - *b. Clinical specificity:* 33% (95% CI, 27.0-39.1)
 - *c. Other clinical supportive data (when a and b are not applicable):* Negative predictive value 100% (95% CI, 95.3-100)
- 4. Clinical cut-off: 500 ng FEU/ml
- 5. Expected values/Reference range:

<500 ng FEU/ml

Patients	N	%Clinical	%Clinical	%Negative
		Sensitivity	Specificity	Predictive Value
		(95% CI)	(95% CI)	(95% CI)
Suspected DVT	295	100.0	39.7	100.0
With Low PTP		(18/18)	(110/277)	(110/110)
		(81.5-100.0)	(33.8-45.7)	(96.7-100.0)
Suspected DVT	189	100.0	26.7	100.0
With Moderate		(17/17)	(46/172)	(46/46)
PTP		(80.5-100.0)	(20.3-34.0)	(92.3-100.0)
Suspected	71	100.0	18.0	100.0
DVT with		(21/21)	(8/50)	(8/8)
High PTP		(83.9-100.0)	(7.2-29.1)	(83.1-100.0)

M. Conclusion:

Data has demonstrated that this device is substantially equivalent to a legally marketed device.