510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY DEVICE ONLY TEMPLATE

A. 510(k) Number:

K032458

B. Analyte:

Calibrator – estradiol, follicle stimulating hormone (FSH), leuteinizing hormone (LH), progesterone, prolactin, triiodothyronine (T3)

C. Type of Test: Calibrator

D. Applicant:

Abbott Laboratories

E. Proprietary and Established Names:

Abbott ARCHITECT[®] Estradiol Calibrators Abbott AxSYM[®] Estradiol Standard Calibrators Abbott AxSYM[®] Estradiol Master Calibrators Abbott ARCHITECT[®] FSH Calibrators Abbott ARCHITECT[®] LH Calibrators Abbott ARCHITECT[®] Progesterone Calibrators Abbott ARCHITECT[®] Prolactin Calibrators Abbott AXSYM[®] Total T3 Standard Calibrators Abbott AxSYM[®] Total T3 Master Calibrators

F. Regulatory Information:

- <u>Regulation section</u>: 21 CFR § 862.1150, Calibrator
- 2. <u>Classification</u>: Class II
- 3. <u>Product Code:</u> JIT, Secondary Calibrators
- 4. <u>Panel:</u> Clinical Chemistry (75)

G. Intended Use:

1. Indication(s) for use:

Abbott ARCHITECT Estradiol Calibrators are devices intended for use in the ARCHITECT Estradiol assay test system to establish points of reference that are used in the quantitative determination of estradiol in human specimens.

Abbott AxSYM Estradiol Standard Calibrators and AxSYM Estradiol Master Calibrators are devices intended for use in the AxSYM Estradiol assay test system to establish points of reference that are used in the quantitative determination of estradiol in human specimens. Abbott ARCHITECT FSH Calibrators are devices intended for use in the ARCHITECT FSH assay test system to establish points of reference that are used in the quantitative determination of follicle-stimulating hormone (FSH) in human specimens.

ARCHITECT LH Calibrators are devices intended for use in the ARCHITECT LH assay test system to establish points of reference that are used in the quantitative determination of luteinizing hormone (LH) in human specimens.

ARCHITECT Progesterone Calibrators are devices intended for use in the ARCHITECT Progesterone assay test system to establish points of reference that are used in the quantitative determination of progesterone in human specimens.

ARCHITECT Prolactin Calibrators are devices intended for use in the ARCHITECT Prolactin assay test system to establish points of reference that are used in the quantitative determination of prolactin in human specimens.

AxSYM Total T3 Standard Calibrators and AxSYM Total T3 Master Calibrators are devices intended for use in the AxSYM Total T3 assay test system to establish points of reference that are used in the quantitative determination of total T3 in human specimens.

2. <u>Special condition for use statement(s)</u>:

For professional use.

3. Special instrument Requirements:

The ARCHITECT test system or the AxSYM test system is required.

H. Device Description:

Calibrator components and validation procedures are described below.

I. Substantial Equivalence Information:

- Predicate device name(s): IMx[®] Estradiol Calibrators Abbott FSH Calibrators AxSYM[®] LH Standard Calibrators AxSYM[®] Progesterone Standard Calibrators Abbott Prolactin Calibrators Ciba-Corning ACS:180 Prolactin Assay Calibrator AxSYM[®] Total T3 Standard and Master Calibrators
- 2. <u>Predicate K number(s):</u> K951629 K890135

K935611 K955025 K896162 K934517

3. <u>Comparison with predicate:</u>

All calibrators are identical to their predicate in intended use.

Similarities <u>ARCHITECT[®] Estradiol Calibrators</u>		
Item	Device	Predicate
matrix	Tris buffer based	Tris buffer based
traceability	Standardized to an internal	Standardized to an internal
	reference standard –	reference standard –
	manufactured	manufactured
	gravimetrically using a	gravimetrically using a
	stock solution of estradiol	stock solution of estradiol
	(not less that 97% pure by	(not less that 97% pure by
	HPLC) at each	HPLC) at each
	concentration level	concentration level
Differences		
Item	Device	Predicate
number	2 levels (0 and 1600 pg/mL)	6 levels (0, 50, 250, 750,
		1500, 3000 pg/mL

Similarities AxSYM [®] Estradiol Standard Calibrators and AxSYM [®] Estradiol			
Master Calibrators			
Item	Device	Predicate	
matrix	Tris buffer based	Tris buffer based	
	Differences		
Item	Device	Predicate	
traceability	Standardized to an internal reference standard for AxSYM Estradiol – standardized to correlate with GC/MS	Standardized to an internal reference standard – manufactured gravimetrically using a stock solution of estradiol (not less that 97% pure by HPLC) at each concentration level	
concentration	6 levels (0, 50, 100, 200, 500, and 1000 pg/mL) (Master contains two of those – 0 and 200 pg/mL – to establish points of reference)	6 levels (0, 50, 250, 750, 1500, 3000 pg/mL	

Similarities <u>ARCHITECT[®] FSH Calibrators</u>		
Item	Device	Predicate
matrix	Bovine serum	Bovine serum
traceability	Manufactured by addition	Manufactured by
	of FSH to a target	gravimetric addition of FSH
	concentration – referenced	to a target concentration –
	against WHO FSH 2 nd	referenced against WHO
	International Reference	FSH 2 nd International
	Preparation (IRP) 78/549	Reference Preparation (IRP)
		78/549
Differences		
Item	Device	Predicate
number	2 levels - 0, $100 mIU/mL$	6 levels – 0, 1, 10, 25, 100,
		250 mIU/mL

Similarities <u>ARCHITECT[®] LH Calibrators</u>		
Item	Device	Predicate
matrix	calf serum	calf serum
traceability	Manufactured	Manufactured
	gravimetrically-referenced	gravimetrically-referenced
	to WHO LH Human,	to WHO LH Human,
	Pituitary 2 nd International	Pituitary 2 nd International
	Standard 80/552 at each	Standard 80/552 at each
	concentration	concentration
Differences		
Item	Device	Predicate
number	2 levels - 2, $100 mIU/mL$	6 levels – 0, 2, 10, 25, 100,
		250 mIU/mL

Simil	Similarities <u>ARCHITECT[®] Progesterone Calibrators</u>		
Item	Device	Predicate	
	Same intended use		
	Differences		
Item	Device	Predicate	
number	2 levels – 0.7, 40 ng/mL	6 levels – 0, 0.7, 2, 7,	
		20, 40 ng/mL	
matrix	Human serum	Serum for level 0 and Tris	
		buffer based for the	
		remaining levels	

Similarities <u>ARCHITECT[®] Prolactin Calibrators</u>		
Item	Device	Predicate
traceability	Referenced to WHO 3 rd	Referenced to WHO 3 rd
	International Standard	International Standard
	84/500 for Prolactin	84/500 for Prolactin
matrix	Tris buffer based	Tris buffer based
Differences		
Item	Device	Predicate
number	2 levels - 5, 30 ng/mL	6 levels – 0, 5, 10, 30,
		80, 200 ng/mL

Similarities AxSYM [®] Total T3 Standard and Master Calibrators			
Item	Device	Predicate	
matrix	Bovine serum	Bovine serum	
number	6 levels – 0, 0.5, 1, 2, 4, 8	6 levels – 0, 0.5, 1, 2, 4, 8	
	ng/mL (only 0 and 1 ng/mL	ng/mL (only 0 and 1 ng/mL	
	for Master)	for Master)	
	Differences		
Item	Device	Predicate	
traceability	Manufactured using USP	Internal reference standards	
	Grade L-triiodothyronine,	are manufactured using L-	
	sodium salt and signal	triiodothyronine Sodium	
	matched to internal	(HPLC purity 95.0 –	
	reference standards which	101.0%). Primary and	
	are traceable to the USP	secondary standards are	
	Reference Standard L-	manufactured	
	triiodothyronine (free acid)	gravimetrically using this	
	at each concentration level.	reference standard. All list	
	The stock solution	material is tested against	
	concentration is determined	these primary and	
	by HPLC	secondary calibrators.	

J. Standard/Guidance Document Referenced (if applicable):

Abbreviated 510(k) Submissions for In Vitro Diagnostic Calibrators

K. Test Principle:

Not applicable as this submission is for calibrators (for two different test systems).

L. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
 - *a. Precision/Reproducibility:* Not applicable. See traceability below.
 - *b. Linearity/assay reportable range:* Not applicable.

c. Traceability (controls, calibrators, or method):

Calibrator validation and stability studies are summarized for each of the calibrators. The sponsor specifies the concentrations of materials evaluated in the studies, the frequency of testing, the method for testing the materials, environmental conditions of storage, and acceptance criteria for the study. Real time studies are being performed to support the stability claims in the labeling.

Abbott ARCHITECT[®] Estradiol Calibrators are standardized to an internal reference standard and manufactured gravimetrically using a stock solution of estradiol (not less that 97% pure by HPLC) at each concentration level.

Abbott AxSYM[®] Estradiol Standard Calibrators and Abbott AxSYM[®] Estradiol Master Calibrators are standardized to an internal reference standard for AxSYM Estradiol which correlates with GC/MS.

Abbott ARCHITECT[®] FSH Calibrators are manufactured by addition of follicle stimulating hormone to a target concentration and referenced against WHO FSH 2nd International Reference Preparation (IRP) 78/549.

Abbott ARCHITECT[®] LH Calibrators are manufactured gravimetrically and referenced to WHO Leuteinizing Hormone Human, Pituitary 2nd International Standard 80/552 at each concentration.

Abbott ARCHITECT[®] Progesterone Calibrators are manufactured by spiking stripped normal male human serum with calculated amounts of a progesterone stock solution. This stock solution is manufactured by adding Progesterone USP gravimetrically to a target concentration which is then assigned to the stock solution. The levels are tested by RLU matching to Progesterone Reference Calibrators. The concentration is adjusted, if necessary, by adding stripped normal male human serum or progesterone stock solution.

Abbott ARCHITECT[®] Prolactin Calibrators are referenced to WHO 3rd International Standard 84/500 for Prolactin.

Abbott AxSYM[®] Total T3 Standard Calibrators and Abbott AxSYM[®] Total T3 Master Calibrators are manufactured using USP Grade L-triiodothyronine, sodium salt and signal matched to internal reference standards which are traceable to the USP Reference Standard L-triiodothyronine (free acid) at each concentration level. The stock solution concentration is determined by HPLC.

d. Detection limit:

Not applicable.

e. Analytical specificity: Not applicable.

- *f.* Assay cut-off: Not applicable.
- 2. Comparison studies:
 - *a. Method comparison with predicate device:* Not applicable.
 - *b. Matrix comparison:* Not applicable. The sponsor summarizes the methods used for analysis and validation of the calibrator material and matrix.
- 3. Clinical studies:
 - *a. Clinical sensitivity:* Not applicable.
 - *b. Clinical specificity:* Not applicable.
 - *c. Other clinical supportive data (when a and b are not applicable):* Not applicable.
- 4. <u>Clinical cut-off:</u> Not applicable.
- 5. <u>Expected values/Reference range:</u> Not applicable.

M. Conclusion:

I recommend that the Abbott ARCHITECT[®] Estradiol Calibrators, Abbott AxSYM[®] Estradiol Standard Calibrators, Abbott AxSYM[®] Estradiol Master Calibrators, Abbott ARCHITECT[®] FSH Calibrators, Abbott ARCHITECT[®] LH Calibrators, Abbott ARCHITECT[®] Progesterone Calibrators, Abbott ARCHITECT[®] Prolactin Calibrators, Abbott AxSYM[®] Total T3 Standard Calibrators, and Abbott AxSYM[®] Total T3 Master Calibrators are substantially equivalent to the legally marketed predicate devices.