510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY DEVICE ONLY TEMPLATE

- **A. 510(k) Number:** k033860
- **B.** Purpose of the Submission: New device
- C. Analyte: Hydroxybutyric Acid
- **D. Type of Test:** Calibrator
- E. Applicant: Wako Chemicals, USA, Inc.
- F. Proprietary and Established Names: Wako Total Ketone Body Calibrators

G. Regulatory Information:

- 1. <u>Regulation section:</u> 21 CFR 862.1150, Calibrators, Secondary
- 2. <u>Classification:</u> II
- 3. <u>Product Code:</u> JIT
- 4. <u>Panel:</u> Chemistry (75)

H. Intended Use:

- 1. <u>Intended use(s):</u> Refer to Indications for use.
- 2. Indication(s) for use:

The Wako Total Ketone Body Calibrators are intended to be used with the Wako Total Ketone Bodies and the Wako Autokit 3-HB test kits to establish points of reference that are used in the determination of values in the measurement of total ketone bodies and 3-HB in human serum or plasma.

(Reviewer Note: 3-HB is a common abbreviation for 3 -Hydroxybutyrate.)

- 3. <u>Special condition for use statement(s):</u> The calibrators are for Rx in vitro diagnostic use.
- 4. <u>Special instrument requirements:</u> Wako 30R®

I. Device Description:

The products consist of two levels of liquid ready to use calibrators. They consist of a matrix of Goods Buffer, polyethylene glycol and sodium azide and are prepared by adding 3-Hydroxybutyrate.

J. Substantial Equivalence Information:

- 1. <u>Predicate device name(s):</u> Stanbio Standard, 1 mM Sodium D-3-Hydroxybutyrate (Formerly GDS)
- 2. <u>Predicate K number(s):</u> k910108
- 3. <u>Comparison with predicate:</u> Both are liquid based calibrators intended to provide reference points for ketone assays. They are intended to calibrate different assays made by different manufacturers.

K. Standard/Guidance Document Referenced (if applicable):

The sponsor did not reference any standards in their submission.

L. Test Principle: Not applicable

M. Performance Characteristics (if/when applicable):

- 1. <u>Analytical performance:</u>
 - *a. Precision/Reproducibility:* Not applicable.
 - *b. Linearity/assay reportable range:* Not applicable.
 - c. Traceability (controls, calibrators, or method):

There are two levels of calibrators: 40 µmol/L and 300µmol/L of 3-Hydroxybutyrate.

Calibrators are prepared by adding 3-Hydroxybuturate to buffer. The grade of chemicals utilized was not specified.

The value assignment and traceability is provided through use of the molar extinction coefficient (K-factor) of NADH.

Results from Real Time Stability studies are summarized. The calibrators were evaluated for up to 15 months when stored at 10 degrees C and were both levels were within 1% of the expected value at all measurement times.

- *d. Detection limit:* Not applicable.
- *e. Analytical specificity:* Not applicable.
- *f.* Assay cut-off: Not applicable.
- 2. <u>Comparison studies:</u>
 - *a. Method comparison with predicate device:* Not applicable.
 - *b. Matrix comparison:* Not applicable.
- 3. Clinical studies:
 - *a. Clinical sensitivity:* Not applicable. Clinical studies are not typically submitted for this device type.
 - *b. Clinical specificity:* Not applicable. Clinical studies are not typically submitted for this device type.
 - c. Other clinical supportive data (when a and b are not applicable):
- 4. <u>Clinical cut-off:</u> Not applicable.
- 5. <u>Expected values/Reference range:</u> Not applicable.

N. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.