

510(k) SUBMISSION TEMPLATE

A. 510(k) Number: K040407

B. Analyte: Antibodies to Protective Antigen (PA) protein, *Bacillus anthracis*

C. Type of Test: ELISA

D. Applicant: Immunetics, Inc.

E. Proprietary and Established Names of the Product: QuickELISA™ Anthrax-PA Kit

F. Regulatory Information:

1. Regulation section: Unclassified
2. Classification: Unclassified
3. Product code: **NRL - ENZYME LINKED IMMUNOABSORBENT ASSAY, ANTIBODY, B. ANTHRACIS**
4. Panel: 83

G. Intended Use and Indication for Use:

1. Intended use(s): The Immunetics® QuickELISA™ Anthrax-PA Kit is intended for use in the qualitative detection of antibodies to the Protective Antigen (PA) protein of *B. anthracis* in human serum. The assay should be used only on serum samples from individuals with a clinical history, signs or symptoms consistent with anthrax infection as an aid in the diagnosis of anthrax, or from recipients of anthrax vaccine.
2. Indication(s) for use: The assay should be used only on serum samples from individuals with a clinical history, signs or symptoms consistent with anthrax infection as an aid in the diagnosis of anthrax, or from recipients of anthrax vaccine.
3. Special conditions for use statement(s):
 - The diagnosis of anthrax infection must be made based on history, signs, symptoms, exposure likelihood, and when possible other laboratory data, in addition to the presence of antibodies to *B. anthracis* PA. Negative results in ELISA should not be used in isolation from other evidence to exclude anthrax.
 - The assay has not been evaluated in individuals without clinical signs or symptoms who were presumed exposed and who subsequently developed cutaneous or inhalation anthrax. There is no information available to interpret the meaning of a positive or negative test result for such individuals.
 - The assay has not been evaluated with specimens from patients infected by the gastrointestinal route; expected results with such infections are unknown.
 - The minimum level of anti-PA antibodies that confers protection following vaccination is not known. The QuickELISA measures total antibody and the

relationship between this value and protective immunity has not been established.

- The affinity and/or avidity of anti-PA IgG and IgM for the rPA antigen have not been determined with this assay.
- This assay is for Rx use.

4. Special instrument requirements: ELISA reader, dual wavelength; ELISA plate washer; microplate shaker

H. Device Description:

The QuickELISA assay kit includes necessary reagents to perform laboratory testing, using immunochemical reactions in an ELISA plate format. Conjugates are manufactured from a recombinant protein product produced by an avirulent plasmidless *B. anthracis* strain. Kit reagents include positive controls (sera from rabbits immunized with recombinant PA), Conjugate A (streptavidin-rPA), Conjugate B (rPA-horseradish peroxidase), a negative control (normal human serum), microwell strips (precoated with biotinylated bovine serum albumin, sample diluent, TMB substrate, wash buffer, and stop reagent. An ELISA reader (dual wavelength), ELISA plate washer and microplate shaker are needed to perform the testing. The cutoff was determined based on preclinical testing to be above the highest negative sample and validated to detect approximately 300 ng/mL of PA-specific antibody (pool of vaccinee sera).

I. Substantial Equivalence Information (if known):

1. Predicate device name(s): preamendments device, US Army Biological Laboratories reagents for modified agar diffusion assay
2. Predicate 510(k) number(s): NA
3. Comparison with predicate:

Similarities		
Item	Device	Predicate
Intended Use	Detection of anti-PA antibodies	Same
Specimen Type	Serum	Same
Negative Control	Normal human serum	Same
Differences		
Item	Device	Predicate
Antigen	Recombinant PA (Protective Antigen)	Purified PA from Sterne strain culture
Assay Endpoint	Spectrophotometric measure (Net OD-450)	Visualized band of identity
Positive Control	Pooled sera from rabbits immunized with rPA	Serum dilution from human vaccinee
Testing Time	Overnight (18-24 h)	35 min (with agitation); 80 min (without agitation)

J. Standard/Guidance Document Referenced (if applicable): NA

K. Test Principle: Immunochemical using ELISA microwell plate format.

L. Performance Characteristics:

1. Analytical performance:
 - a. *Precision(Repeatability/Reproducibility):*

<i>Reproducibility Panel Sample</i>	<i>Agitated Incubations: Intra-Assay</i>		<i>Agitated Inter-Assay</i>	<i>Agitated Inter-Lot</i>	<i>Agitated Inter-Site</i>	<i>Stationary Intra-Assay</i>	<i>Stationary vs. Agitated</i>
	<i>Net OD-450</i>	<i>CV %</i>	<i>CV %</i>	<i>CV %</i>	<i>CV %</i>	<i>CV %</i>	<i>CV %</i>
PC	2.233	4.4	2.0	11.2	35.4	10.8	15.1
Low PC	0.860	4.1	1.5	10.2	13.8	13.2	11.6
NC	0.018	7.2	6.2	5.4	18.1	10.8	27.9
HiP1	3.809	6.2	0.4	2.2	3.9	6.6	1.3
HiP2	3.830	5.6	0.5	0.9	3.5	7.6	1.2
HiP3	3.806	5.2	2.0	1.7	1.9	6.9	0.8
LoP1	2.225	6.1	2.5	22.6	15.3	3.2	19.0
LoP2	0.609	7.9	3.0	26.3	4.0	4.5	27.1
LoP3	0.442	8.3	6.4	20.5	10.0	7.9	22.9
N1	0.019	9.5	21.5	19.2	7.9	10.9	12.8
N2	0.021	8.0	23.1	14.9	23.9	9.9	12.8
N3	0.021	5.9	9.5	18.5	19.8	3.0	37.3
N4	0.021	12.2	6.4	13.6	20.9	2.3	33.5
N5	0.023	14.3	8.6	24.2	20.2	6.1	32.4
REF-1	0.357	8.6	8.4	15.7	13.7	7.6	15.4
MeP1	2.575	6.7	7.8	16.0	15.9	5.0	5.2

- b. *Linearity/Assay measuring (reportable) range:* NA
 - c. *Calibrators, and Controls:* Pooled sera from rabbits immunized with rPA
 - d. *Traceability of the assay:* NA
 - e. *Analytical limits at low levels:* approx. 300 ng of PA-specific antibody
 - f. *Analytical specificity:*
Sera from 225 individuals with other antibody responses associated with various disease conditions, and 583 sera from healthy blood donors were tested with QuickELISA. Positive results were obtained for 5/583 (0.9%) blood donors and 1/225 (0.4%) of patients with other disease conditions.
 - g. *Analytical characterization of cut-off:*
Titrations of pooled human sera from AVA vaccinees (standardized using radial immunodiffusion and nephelometry with US National Reference Preparation for Specific Human Serum Proteins, and mass value using USFDA reference serum, 1983 *H. influenzae* type b)
2. Comparison studies using clinical specimens:
 - a. *Method comparison:* NA
 - b. *Matrix description and comparison:* NA
3. Clinical Studies :
 - a. *Clinical sensitivity and clinical specificity:* NA – anthrax infection would not be expected in a US population.

