510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY ONLY TEMPLATE

A. 510(k) Number:

k043242

B. Purpose for Submission:

This is a new device for an Acetaminophen assay configured with other analytes which have been previously cleared (k012745) (see device description section).

C. Measurand:

Acetaminophen/paracetamol (APAP) Amphetamines (AMP), barbiturates (BAR), benzodiazepines (BZO), cocaine (COC), methamphetamines (mAMP), opiates (OPI), phencyclidine (PCP), THC and tricyclic antidepressants (TCA) in urine.

D. Type of Test:

Qualitative, Fluorescence Immunoassay

E. Applicant:

Biosite Incorporated

F. Proprietary and Established Names:

Triage® TOX Drug Screen (available in different configurations)

G. Regulatory Information:

1. Regulation section:

862.3030 - Acetaminophen Test System

862.3100 - Amphetamine Test System

862.3150 - Barbiturate Test System

862.3170 - Benzodiazepine Test System

862.3250 - Cocaine and Cocaine Metabolite Test System,

862.3610 - Methamphetamine Test System

862.3650 - Opiate Test System

862.3100 - Amphetamine test system (Phencyclidine)

862.3870 - Cannabinoid Test System and

862.3910 - Tricyclic Antidepressant Drug Test System

2. Classification:

Class II

3. Product code:

LDP, DKZ, DIS, JXM, JXO, LAF, DJG, LCM, LDJ and LFG; respectively

4. Panel:

91, Toxicology

H. Intended Use:

1. <u>Intended use(s):</u>

See indications for use.

<u>Indication(s) for use:</u>

The Triage Tox Drug Screen is a Fluoresence immunoassay intended to be used with the Triage MeterPlus for the point-of-care qualitative determination of the presence of major metabolites above the threshold concentrations of up to 8 distinct drug classes, including assays for acetaminophen/paracetamol, amphetamines, barbiturates, benzodiazepines, cocaine, methamphetamines, opiates, phencyclidine, THC and tricyclic antidepressants in urine. The acetaminophen/paracetamol assay will yield positive results when acetaminophen/paracetamol is ingested at or above therapeutic doses.

The threshold concentrations are provided below:

Acetaminophen/Paracetamol	APAP	5 μg/mL
Amphetamines	AMP	1000 ng/mL
Methamphetamines	mAMP	1000 ng/mL
Barbiturates	BAR	300 ng/mL
Benzodiazepines	BZO	300 ng/mL
Cocaine	COC	300 ng/mL
Opiates	OPI	300 ng/mL
Phencyclidine	PCP	25 ng/mL
THC	THC	50 ng/mL
Tricyclic Antidepressants	TCA	1000 ng/mL

2. Special conditions for use statement(s):

This test provides only preliminary test results. Clinical consideration and professional judgment must be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result. In order to obtain a confirmed analytical result, a more specific alternate chemical method is needed. Gas Chromatography/Mass Spectroscopy (GC/MS) is the preferred confirmatory method.

A quantitative serum acetaminophen/paracetamol measurement is the common confirmatory method for preliminary positive acetaminophen/paracetamol results.

For Prescription and Point-of-Care use. The Point-of-Care classification was address in a previously cleared submission (k012745)

4. Special instrument requirements:

Triage® MeterPlus; k973547

I. Device Description:

The Triage TOX Drug Screen is a competitive fluorescence immunoassay used for the qualitative determination of major metabolites of acetaminophen/paracetamol (APAP), amphetamines (AMP), methamphetamines (mAMP), barbiturates (BAR), benzodiazepines (BZO), cocaine (COC), opiates (OPI), phencyclidine (PCP), THC and tricyclic antidepressants (TCA). All analytes were previously cleared (k012745) except for the acetaminophen analyte. It is a single use test device and is used in conjunction with the Triage ® MeterPlus. The device contains murine monoclonal antibodies and drug metabolites labeled with a fluorescent dye or immobilized on the solid phase and stabilizers. The testing device is inserted into and read by the Triage MeterPlus. Threshold concentrations are used to separate a negative result from a presumptive positive result.

No human source materials are used during the manufacturing of this product.

J. Substantial Equivalence Information:

1. Predicate device name(s):

Triage® TOX Drug Screen

2. Predicate 510(k) number(s):

k012745

3. Comparison with predicate:

Similarities				
Item	Device	Predicate		
Intended use	Qualitative determination	Qualitative determination		
	of major metabolites in	of major metabolites in		
	human urine.	human urine.		
Matrix	Human Urine	HumanUrine		
Test Principle	Fluorescence	Fluorescence		
	immunoassay	immunoassay		
	Differences			
Item	Device	Predicate		
Anaylets	Acetaminophen,	Amphetamines,		
	Amphetamines,	Methamphetamines,		
	Methamphetamines,	Barbiturates,		
	Barbiturates,	Benzodiazepines,		
	Benzodiazepines,	Cocaine, Opiates,		
	Cocaine, Opiates,	Phencyclidine, THC and		
	Phencyclidine, THC and	Tricyclic Antidepressants		
	Tricyclic Antidepressants			

K. Standard/Guidance Document Referenced (if applicable):

None stated.

L. Test Principle:

The Triage Tox Drug Screen is a competitive fluorescence immunoassay used for the qualitative determination of parent compound and major metabolites of drugs of abuse in urine specimens. A urine sample is added to the device which then passes through a filter. The sample moves by capillary action into a reaction chamber and is allowed to react with murine monoclonal fluorescent antibody conjugates or mix fluorescent drug conjugates. After an incubation period, the reaction mixture flows through the device detection lane. The presence of drug or drug metabolite in the urine specimen prevents binding of the fluorescent conjugates to the solid phase on the detection zone. Unbound fluorescent conjugate is washed from the detection lane by excess urine. The concentration of the drug metabolite in the urine specimen is inversely related to the fluorescence bound to the detection zone. Threshold concentrations are used to separate a negative result from a presumptive positive result.

M. Performance Characteristics (if/when applicable):

See k012745 for amphetamines, barbiturates, benzodiazepines, cocaine, methamphetamines, opiates, phencyclidine, THC and tricyclic antidepressants. Sponsor states that performance studies have been conducted for the addition of acetaminophen through its internal Design Control process, the performance characteristics are exactly the same and data are on file at the Biosite.

1. Analytical performance:

a. Precision/Reproducibility:

The threshold for acetaminophen of 5 ng/ml was challenged by testing specimens containing the drug or drug metabolite spiked into drug-free urine at concentrations within 25% and 50% of the threshold. Each specimen was tested using the Triage® Tox Drug Screen. The data paralleled the expected agreement based on the coefficient of variation of the assay.

b. Linearity/assay reportable range:

Not applicable. The assay is intended for qualitative use.

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Users are instructed to run external controls with each new lot or shipment of test materials, or every 30 days, and as otherwise required by federal, state or local guidelines. Commercial control materials are required but are not specifically identified in the labeling.

Each device has internal process controls which run with every sample. The controls monitor various procedural and functional parameters such as sufficient sample added to the test device and that the unbound fluorescent label has been washed sufficiently from the detection zone. During testing if the results are within the limits the results are reported. If the results are not within the limits the results are not reported and an error message will be displayed.

d. Detection limit:

Cutoff validation and sensitivity were evaluated. Four contrived samples were prepared by adding acetaminophen to drug-free urine at concentrations within 25% and 50% of the cutoff concentration. Samples were measured 10 times using the Triage TOX Drug Screen.

Table IV.1: APAP Cutoff Validation

	0.5x	0.75x	1.25x	1.5x
POS	0	0	9	10
NEG	10	10	1	0
Concordance	100%	100%	90%	100%

e. Analytical specificity:

Substances that may cause interference were added to samples containing acetaminophen within 25% of the cutoff concentration. Cross-reactivity substances were added to drug-free urine. Samples were measured in duplicate. The substances did not interfere or have a cross-reaction with the acetaminophen assay at the concentrations listed below:

Table IV.2.1: Interfering Substances

Substance	Concentration	Substance	Concentration
Acetone	5 mg/mL	Hemoglobin	1.2 mg/mL
Acetylsalicylic Acid	1 mg/mL	Human Serum Albumin	5 mg/mL
Ascorbic Acid	15 mg/ml	Ibupprofen	1 mg/mL
Bilirubin	$2.5 \mu g/mL$	Ketamine	25 μg/mL
Caffeine	0.125 mg/mL	Oxalic Acid	10 mg/mL
Creatinine	2.5 mg/mL	Riboflavin	75μg/mL
Dextrose	20 mg/mL	Scopalamine	62.5 μg/mL
Ethanol	5 mg mL	Sodium Chloride	30 mg/mL
Fluoxetine	0.5 mg/mL	Urea	30 mg/mL
Gamma Globulin	5 mg/mL		

Table IV.2.2: Cross-Reacting Substances

Substance	Concentration	Substance	Concentration
Acetopromazine	100 μg/mL	Methoxyphenamine	100 μg/mL
Benzphentamine	100 μg/mL	Methylphenidate	100 μg/mL
Benztropine Methane	100 μg/mL	Naloxone	80 μg/mL
Bupropin	100 μg/mL	Naproxen	100 μg/mL
Butyrophenine	100 μg/mL	Norpsuedoephedrine	100 μg/mL
Cimetidine	100 μg/mL	Phenethylamine	100 μg/mL
Clonidine	100 μg/mL	Phenmetrazine	100 μg/mL
Cotinine	100 μg/mL	Phenylephrine	100 μg/mL
Dextromethorphan	100 μg/mL	Phenylhydantoin, d/l-5-	100 μg/mL
		(p-hydroxyphenyl)-5-	
Dextrorphan	100 μg/mL	Phenylpropanolamine	100 μg/mL
Diphenhydramine	100 μg/mL	Promrthazine	100 μg/mL
Dopamine	100 μg/mL	Propanolol, d/l	100 μg/mL
Epinephrine, 1-	100 μg/mL	Propoxyphene	100 μg/mL
Fenfluramine	20 μg/mL	Pseudoephedrine, d-	100 μg/mL
Glutethimide	100 μg/mL	Quinacrine	100 μg/mL
KetorolacTromethane	100 μg/mL	Ranitidine	100 μg/mL
Levorphanol	50 μg/mL	Thioridazine	100 μg/mL
Meperidine	100 μg/mL	Tramadol	100 μg/mL
Mesoridazine	100 μg/mL	Tyamine	60 μg/mL
Methadone d/l	100 μg/mL	Tranylcypromine	100 μg/mL
Methaqualone	100 μg/mL	Zolpidem	100 μg/mL

Acetaminophen was added at 75% and 125% of the cutoff concentration to urine samples with specific gravities ranging from 1.006 - 1.026 and to urine samples with the pH ranging from 4.5 - 7.5. Six (6) replicates of each sample were measured using the Triage TOX Drug Screen. There is no effect of specimen specific gravity or pH on the acetaminophen assay results.

f. Assay cut-off:

Characterization of how the device performs analytically around the claimed cutoff concentration appears in precision and sensitivity sections, above.

The cutoff concentration of this qualitative assay was arbitrarily chosen by the sponsor. There is no known clinical significance regarding the cutoff concentration. It indicates only that the analyte was or was not present at that level.

2. Comparison studies:

See k012745 for amphetamines, barbiturates, benzodiazepines, cocaine, methamphetamines, opiates, phencyclidine, THC and tricyclic antidepressants.

a. Method comparison with predicate device:

Because the candidate was compared to a reference method, GC/MS. It was not compared to a predicate device.

102 specimens with known acetaminophen concentration by GC/MS were measured using the Triage TOX Drug Screen. 8 specimens were within 25% below the cutoff concentration and 12 specimens were with 25% above the cutoff concentration. The cutoff concentration for both methods was 5 μ g/mL.

The overall agreement of results from both methods is 97.1% (95% confidence interval 93.8%-100%).

b. Matrix comparison:

Not applicable. The assay is intended for only one sample matrix.

3. Clinical studies:

a. Clinical Sensitivity:

Not applicable.

b. Clinical specificity:

Not applicable.

c. Other clinical supportive data (when a. and b. are not applicable):

4. Clinical cut-off:

Not applicable.

5. Expected values/Reference range:

Not applicable.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.