510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY ONLY TEMPLATE

A. 510(k) Number:

k060027

B. Purpose for Submission:

Clearance of new device

C. Measurand:

Controls and Calibrators for dehydroepiandrosterone sulfate (DHEA-S) assay

D. Type of Test:

Not applicable. This submission is for clearance of controls and calibrators.

E. Applicant:

BIOKIT S.A.

F. Proprietary and Established Names:

ARCHITECT DHEA-S CALIBRATORS (A-F) AND CONTROLS (LOW, MEDIUM AND HIGH)

G. Regulatory Information:

Product Code	Classification	Regulation Section	Panel
Single (Specified)	<u>Class I</u>	<u>21 CFR 862.1660,</u>	75 CLINICAL
Analyte Controls		Quality control	CHEMISTRY (CH)
(Assayed And		material (assayed	
Unassayed) (JJX)		and unassayed).	
Calibrator,	Class II	21 CFR 862.1150,	75 CLINICAL
Secondary (JIT)		Calibrator.	CHEMISTRY (CH)

H. Intended Use:

1. Intended use(s):

Calibrators

The ARCHITECT[®] DHEA-S Calibrators are for the calibration of the ARCHITECT i System when used for the quantitative determination of DHEA-S in human serum and plasma.

Controls

The ARCHITECT[®] DHEA-S Controls are for the verification of the accuracy and precision of the ARCHITECT i System when used for the quantitative determination of DHEA-S in human serum and plasma.

2. Indication(s) for use:

Calibrators

The ARCHITECT[®] DHEA-S Calibrators are for the calibration of the ARCHITECT i System when used for the quantitative determination of DHEA-S in human serum and plasma.

Controls

The ARCHITECT[®] DHEA-S Controls are for the verification of the accuracy and precision of the ARCHITECT i System when used for the quantitative determination of DHEA-S in human serum and plasma.

3. <u>Special conditions for use statement(s):</u>

For Prescription use only.

4. Special instrument requirements:

The ARCHITECT[®] DHEA-S Calibrators and Controls are designed to be used for the calibration of the ARCHITECT[®] DHEA-S on the ARCHITECT[®] *i*2000 or *i*2000sr instrument platform. The instrument system was 510k cleared through submission number K983212.

I. Device Description:

Calibrator A contains human serum nonreactive for HBsAg, anti-HIV-1/HIV-2, HIV NAT, anti-HBc, HCV NAT, and anti-HCV. Calibrators B-F contain purified synthetic DHEA-S in human serum nonreactive for HBsAg, anti-HIV-1/HIV-2, HIV NAT, anti-HBc, HCV NAT, and anti-HCV. Preservative: Sodium azide. Each DHEA-S Calibrator kit contains 6 bottles of Calibrators (a 2.0 ml fill volume per bottle) with the following DHEA-S concentrations: 0, 5, 12, 60, 300, and 1500 µg/dL.

The controls contain purified synthetic DHEA-S in human serum nonreactive for HBsAg, anti-HIV-1/HIV-2, HIV NAT, anti-HBc, HCV NAT, and anti-HCV. Preservative: Sodium azide. Each DHEA-S control kit contains 3 bottles of controls (a 4.0 ml fill volume per bottle) with the following target concentrations (low, Medium, and High): 10, 100, and 1000 μ g/dL.

J. Substantial Equivalence Information:

Predicate	K935806	
Describe the item being compared		
DPC IMMULITE ® DHEA-SO4 Calibrators and Controls		

Similarities - Calibrators			
Item	Device	Predicate	
Intended Use	The ARCHITECT® DHEA-S Calibrators are for the calibration of the ARCHITECT <i>i</i> System when used for the quantitative determination of DHEA -S in human serum and plasma.	DHEA-SO4 For the quantitative determination of dehydroepiandrosterone sulfate in serum.	
System Methodology	Chemiluminescent Microparticle Immunoassay (CMIA)	Chemiluminescent a solid phase enzyme immunoassay	
Assay Protocols	Competitive assay	Competitive assay	
Matrix	DHEA-S (synthetic) in human serum with preservative. (SodiumAzide)	Lyophilized DHEA -SO4 in human serum with preservative	

Differences - Calibrators			
Item	Device	Predicate	
Platform	ARCHITECT <i>i</i> System	Immulite 2000 Analyzer	
Calibration	0, 5, 12, 60, 300, and 1500	15, 1,000, ug/dI	
Range/Levels	μg/dL	15-1,000 μg/dL	

Similarities - Controls			
Item	Device	Predicate	
Intended Use	The ARCHITECT® DHEA-S Calibrators are for the calibration of the ARCHITECT <i>i</i> System when used for the quantitative determination of DHEA -S in human serum and plasma.	DHEA-SO4 For the quantitative determination of dehydroepiandrosterone sulfate in serum.	
System Methodology	Chemiluminescent Microparticle Immunoassay (CMIA)	Chemiluminescent a solid phase enzyme immunoassay	
Assay Protocols	Competitive assay	Competitive assay	
Matrix	DHEA-S (synthetic) in human serum with preservative. (SodiumAzide)	Lyophilized DHEA -SO4 in human serum with preservative	

Differences - Controls			
Item	Device	Predicate	
Platform	ARCHITECT <i>i</i> System	Immulite 2000 Analyzer	
	Single constituent, DHEA -S	Multi constituent, Human serum	
Matrix	(synthetic) in human serum	based tri-level control containing	
	with preservative. (Sodium	over 25 constituents commonly	
	Azide)	measured by immunoassay.	
Levels	10, 100, and 1000 µg/dL	Three levels, see package insert for	
		measured values	
Assay sample	Serum and plasma	Serum	
type	Serum and plasma	Seruili	

K. Standard/Guidance Document Referenced (if applicable):

STANDARDS

Title and Reference Number

Interference Testing in Clinical Chemistry; Approved Guideline (CLSI EP 7-A)

Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline (CLSI EP5-A)

Document Title	Office	Division	Web Page
Guidance for Industry - Abbreviated 510(k) Submissions for In Vitro Diagnostic Calibrators; Final	OIVD		http://www.fda.gov/cdrh/ode/calibrator.html
Points to Consider Guidance Document on Assayed and Unassayed Quality Control Material; Draft	OIVD		http://www.fda.gov/cdrh/ode/99.html

L. Test Principle:

Not applicable. This submission is for calibrators and controls.

M. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
 - a. Precision/Reproducibility:

Not applicable

b. Linearity/assay reportable range:

Not applicable

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Primary calibrators and controls are prepared gravimetrically by adding Synthetic DHEA-S obtained from a commercial source to a human serum-based matrix. The secondary calibrators and controls are also prepared gravimetrically by adding Synthetic DHEA-S obtained from a commercial source to a human serum-based matrix. Relative light unit (RLU) testing is performed and RLU values are compared to the corresponding Primary Calibrator or Control. The sponsor's acceptance criteria is that RLU variation must be within +/- 1.5%.

Target values for secondary calibrators are 0, 5, 12, 60, 300, and 1500 μ g/dL DHEA-S. Target values for secondary controls are 10, 100, and 1000 μ g/dL DHEA-S.

The DHEA-S Calibrators and Controls are stable for 6 months if stored at 2 to 8 °C. Protocols and acceptance criteria for stability testing were described and found to be acceptable.

d. Detection limit:

Not applicable

e. Analytical specificity:

Not applicable

f. Assay cut-off:

Not applicable

2. Comparison studies:

a. Method comparison with predicate device:

Not applicable

b. Matrix comparison:

Not applicable

- 3. <u>Clinical studies</u>:
 - a. Clinical Sensitivity:

Not applicable

b. Clinical specificity:

Not applicable

c. Other clinical supportive data (when a. and b. are not applicable):

Not applicable

4. Clinical cut-off:

Not applicable

5. Expected values/Reference range:

Not applicable

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.