510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY ONLY TEMPLATE

A. 510(k) Number:

k060854

B. Purpose for Submission: New Device

C. Measurand:

Cholesterol HDL Cholesterol LDL Cholesterol Triglycerides HDL Cholesterol Calibrator LDL Cholesterol Calibrator

D. Type of Test:

Quantitative immunoturbidimetric and enzymatic photometric assays

E. Applicant:

HORIBA ABX

F. Proprietary and Established Names:

ABX PENTRA 400: Lipoproteins

G. Regulatory Information:

<u>Regulation section:</u> 21CFR Sec. 862.1150 - Calibrator. 21CFR Sec. 862.1660 - Quality control material (assayed and unassayed). 21CFR Sec. 862.1705 - Triglyceride test system. 21CFR Sec. 862.1175 - Cholesterol (total) test system. 21CFR Sec. 862.1475 - Lipoprotein test system.

2. <u>Classification:</u>

Class II, Class I reserved, and Class I meets the limitations to exemptions 21 CFR §862.9 (c) (4)

3. <u>Product code:</u>

JIT - Calibrator, Secondary JJY - Multi-Analyte Controls, All Kinds (Assayed and Unassayed) CDT - Lipase Hydrolysis/Glycerol Kinase Enzyme CHH - Enzymatic Esterase--Oxidase, Cholesterol LBR - LDL & VLDL Precipitation, HDL MRR – Low Density Lipoprotein 4. <u>Panel:</u> Chemistry (75)

H. Intended Use:

- 1. <u>Intended use(s):</u> See indications for use below.
- 2. Indication(s) for use:

ABX PENTRA Cholesterol CP reagent with associated calibrators and controls are for quantitative in vitro diagnostic determination of cholesterol in human serum and plasma based on an enzymatic photometric test (Trinder's reaction).

ABX PENTRA HDL Direct CP reagent with associated calibrators and controls are for quantitative in vitro diagnostic determination of High Density Lipoprotein Cholesterol (HDL-C) in human serum and plasma based on an enzymatic assay with accelerator selective detergent methodology.

ABX PENTRA LDL Direct CP reagent with associated calibrators and controls are for quantitative in vitro diagnostic determination of Low Density Lipoprotein Cholesterol (LDL-C) in human serum and plasma based on an enzymatic colorimetric assay.

Cholesterol measurements are used in the diagnosis and treatment of disorders involving excess cholesterol in the blood and lipid and lipoprotein metabolism disorders.

ABX PENTRA Triglycerides CP reagent with associated calibrators and controls are for quantitative in vitro diagnostic determination of triglycerides in human serum and plasma based on an enzymatic colorimetric assay.

Measurements obtained by this device are used in the diagnosis and treatment of patients with diabetes mellitus, nephrosis, liver obstruction, other diseases involving lipid metabolism, or various endocrine disorders.

The ABX PENTRA HDL Cal is a calibrator for use in the calibration of quantitative Horiba ABX PENTRA HDL Direct CP method on Horiba ABX clinical chemistry analyzers.

The ABX PENTRA LDL Cal is a calibrator for use in the calibration of quantitative Horiba ABX PENTRA LDL Direct CP method on Horiba ABX clinical chemistry analyzers.

- 3. <u>Special conditions for use statement(s):</u> Prescription Use
- 4. <u>Special instrument requirements:</u> ABX PENTRA 400 (k052007)

I. Device Description:

ABX Pentra Cholesterol CP is ready-to-use.

Reagent	Good's buffer pH 6.7 Phenol 4-Aminoantipyrine Cholesterol esterase (CHE) Cholesterol oxidase (CHO) Peroxidase (POD) Sodium azide	$\begin{array}{l} 50 \text{ mmol/L} \\ 5 \text{ mmol/L} \\ 0.3 \text{ mmol/L} \\ \geq 200 \text{ U/L} \\ \geq 50 \text{ U/L} \\ \geq 3 \text{ kU/L} \\ 0.95 \text{ g/L} \end{array}$
ABX Pentra HI	DL Direct CP is ready-to-use.	
Reagent 1	Good's Buffer Cholesterol oxidase Peroxidase N,N-bis(4-sulphobutyl)- m- toluidine- disodium (DSBmT) Accelerator Preservative	< 1000 U/L < 1300 ppg U/L < 1 mM < 1 mM < 0.06 %
Reagent 2	Good's Buffer Cholesterol esterase 4- Aminoantipyrine (4-AAP) Detergent Restrainer Preservative Ascorbic acid oxidase	< 1500 U/L < 1mM < 2 % < 0.15% < 0.06 % < 3000 U/L
ABX Pentra LI	DL Direct CP is ready-to-use.	
Reagent 1	MES Buffer (pH 6.3) Detergent 1 Cholesterol Esterase Cholesterol Oxidase Peroxidase 4-aminoantipyrine Ascorbic Acid Oxidase Preservative	< 1.0 % < 1500 U/L < 1500 U/L < 1300 ppg U/L < 0.1 % < 3000 U/L
Reagent 2	MES Buffer (pH 6.3) Detergent 2 N,N-bis(4-sulfobutyl)-toluidine, disodium (DsBmT) Preservative	< 1.0 % < 1.0 mM

ABX Pentra Triglycerides CP is ready-to-use.

Reagent	Pipes free acid	50 mmol/L
_	Sodium hydroxide	3.36 g/L
	Triton X-100	1 ml/L
	Magnesium salt	14.8 mmol/L
	p-chlorophenol	2.7 mmol/L
	ATP	3.15 mmol/L
	Sodium azide	7.99 mmol/L
	Potassium ferrocyanide	10 µmol/L
	4-aminoantipyrine	0.31 mmol/L
	Lipoprotein lipase	\geq 2000 U/L
	Glycerokinase	\geq 500 U/L
	Glycerol phosphate Oxidase	\geq 4000 U/L
	Peroxidase	\geq 500 U/L

ABX Pentra HDL Cal is a lyophilized calibrator. It is a preparation of lyophilized human serum containing lipoproteins from the various lipoprotein classes including high-density lipoproteins.

ABX Pentra LDL Cal is a lyophilized calibrator. It is a preparation of lyophilized human serum containing lipoproteins from the various lipoprotein classes including low-density lipoproteins. • The kit is composed of 2 vials of calibrator (lyophilized for 1 mL)

J. Substantial Equivalence Information:

- Predicate device name(s): Roche Cholesterol Reagent Ultra N-Geneous HDL Cholesterol Reagent N-Geneous LDL Cholesterol Reagent Roche Reagent for Triglycerides Roche Apolipoprotein Standard Ultra N-Geneous HDL Cholesterol Reagent Calibrator N-Geneous LDL Cholesterol Reagent Calibrator Liquichek Lipids Control
- 2. <u>Predicate 510(k) number(s):</u> k941573, k021316, k971573, k893973, k021316, k971573, respectively
- 3. Comparison with predicate:

	Predicate device (k941573):	Device :
Device Name	Cholesterol	ABX Pentra Cholesterol CP
Manufactured by	Roche, USA	HORIBA ABX, France
Instrument	COBAS MIRA chemistry	ABX PENTRA 400
	system	
Analytes	Cholesterol	Cholesterol

	Predicate device (k941573):	Device :
Device Name	Cholesterol	ABX Pentra Cholesterol CP
Method :	Enzymatic photometric test	Enzymatic photometric test
Specimen :	Serum	Serum
	Plasma	Plasma
Certification :	Traceability to the National	Traceability to the National
	Reference System for	Reference System for
	Cholesterol	Cholesterol
Component reagent	Single-reagent bottle, ready to	Mono-reagent cassette, ready to
matrices	USE: DEACENT : 4	use: DEACENT : Cood's huffer
	Aminoantinyrina Cholesterol	REAGENT . Good S bullet, Phonol 4 Aminoantinyring
	Esterase Cholesterol Oxidase	Cholesterol Esterase, Cholesterol
	Perovidase Phenol huffers	Oxidase Perovidase Sodium
	stabilizers fillers	azide
Format	Liquid	Liquid
Labels	-	Horiba ABX specific label
Packaging	Single-reagent bottles	Mono-reagent cassette :
	REAGENT : 15 x 12 ml	REAGENT : 99 mL
	or 4 x 12 mL	
Controls	Commercially available quality	Recommended quality control
	control (not included): 2 levels	material (not included):
	of serum-based controls with	ABX Pentra N Control
	known cholesterol values	(Normal control)
		ABX Pentra P Control
Calibratara	Commonoially available	(Pathologic control)
Camprators	collibrator (not included)	material (not included):
	canorator (not included)	ABX Pentra Multical
Performance data :		
Number of tests	-	380 tests
Sample volume	3 µL/test	3 µL/test
Detection limit	2 m a/dI	4
Detection minit	2 mg/dL	4 mg/dL
Accuracy and	CV Total < 2.4%	CV Total < 3 01%
Precision		
	2 mg/dL - 525 mg/dL	2.55 mg/dL - 583.26 mg/dL
Measuring range		
	525 mg/dL (with automatic post-	580 mg/dL
Upper linearity limit	dilution : 1050 mg/dL)	
	Calibration required monthly	8 days
Calibration stability	Canoration required montiny	o days
	Until the expiration date when	24 months at 2-8°C
Closed reagent	stored at 2-8°C	
stability		
	30 days at 2-8°C	on-board stability (refrigerated
		area): 48 days

	Predicate device (k941573):	Device :
Device Name	Cholesterol	ABX Pentra Cholesterol CP
Open Reagent stability		
	Predicate device (k021316):	Device :
Device Name	Ultra N-geneous® HDL Cholesterol Reagent	ABX Pentra HDL Direct CP
Manufactured by	Genzyme, USA	HORIBA ABX, France
Instrument	Automated clinical chemistry analyzer	ABX PENTRA 400
Analytes	High Density Lipoprotein Cholesterol (HDL-C)	High Density Lipoprotein Cholesterol (HDL-C)
Method :	Accelerator Selective detergent Methodology	Accelerator Selective detergent Methodology
Specimen :	Serum Plasma	Serum Plasma
Component reagent matrices	Single reagent bottles, ready to use: REAGENT 1 : Buffer, Cholesterol oxidase, Peroxidase, N,N-bis(4-sulphobutyl)-m- toluidine-disodium (DSBmT), Accelerator, Preservative, Ascorbic oxidase REAGENT 2 : Buffer, Cholesterol esterase, 4- Aminoantipyrine (4-AAP), Detergent, Preservative	Bi-reagent cassette, ready to use REAGENT 1 : Good's buffer, Cholesterol oxidase, Peroxidase, N,N-bis(4-sulphobutyl)-m- toluidine-disodium (DSBmT), Accelerator, Preservative REAGENT 2 : Good's buffer, Cholesterol esterase, 4- Aminoantipyrine (4-AAP), Detergent, Restrainer, Preservative, Ascorbic acid oxidase
Format	Liquid	Liquid
Labels	-	Horiba ABX specific label
Packaging	Single-reagent bottles REAGENT 1 : 1 x 250 mL REAGENT 2 : 1 x 80 mL	Bi-reagent cassette : REAGENT 1 : 62 mL REAGENT 2 : 21 mL
Controls	Recommended quality control material (not included): Normal and near the concentrations for decision- making controls	Recommended quality control material (not included): ABX Pentra N Control (Normal control) ABX Pentra P Control (Pathologic control)
Calibrators	Recommended calibration material (not included): Ultra N_geneous® HDL Cholesterol Calibrator	Recommended calibration material (not included): ABX Pentra HDL Cal
Performance data :		
Number of tests	-	240 tests
Sample volume	3 µL/test	3 µL/test
Detection limit	-	1.16 mg/dL

	Predicate device (k941573):	Device :
Device Name	Cholesterol	ABX Pentra Cholesterol CP
Accuracy and Precision	CV Total < 1.5%	CV Total < 3.52%
	33.6 mg/dL – 133 mg/dL	5.4 mg/dL - 151.9 mg/dL
Measuring range	200 mg/dL	154.8 mg/dL
Upper linearity limit		C C
	-	
Calibration stability		14 days
Closed reagent stability	Until the expiration date when stored at 2-8°C	22 months at 2-8°C
	4 weeks at 2-8°C	on-board stability (refrigerated area): 31 days
Open Reagent stability		

	Predicate device (k971573):	Device :
Device Name	N-geneous® LDL	ABX Pentra LDL Direct CP
	Cholesterol Reagent	
Manufactured by	Genzyme, USA	HORIBA ABX, France
Instrument	Automated clinical chemistry	ABX PENTRA 400
	analyzer	
Analytes	Low Density Lipoprotein	Low Density Lipoprotein
	Cholesterol (LDL-C)	Cholesterol (LDL-C)
Method :	Enzymatic colorimetric assay	Enzymatic colorimetric assay
Specimen :	Serum	Serum
	Plasma	Plasma
Component reagent	Single reagent bottles, ready to	Bi-reagent cassette, ready to use
matrices	use:	
	REAGENT 1 : MES buffer,	REAGENT 1 : MES buffer,
	Detergent 1, Cholesterol	Detergent 1, Cholesterol
	Esterase, Cholesterol Oxidase,	Esterase, Cholesterol Oxidase,
	Peroxidase, 4-Aminoantipyrine,	Peroxidase, 4-Aminoantipyrine,
	Ascorbic Acid Oxidase,	Ascorbic Acid Oxidase,
	Preservative	Preservative
	REAGENT 2 : MES buffer,	REAGENT 2 : MES buffer,
	Detergent 2, N,N-bis(4-	Detergent 2, N,N-bis(4-
	sulfobutyl)-m-toluidine-	sulfobutyl)-toluidine-disodium
	disodium (DSBmT),	(DSBmT), Preservative,
	Preservative	
Format	Liquid	Liquid
Labels	-	Horiba ABX specific label
Packaging	Single-reagent bottles :	Bi-reagent cassette :
	REAGENT 1 : 1 x 30 mL	REAGENT 1 : 28 mL
	REAGENT 2 : 1 x 10 mL	REAGENT 2 : 10 mL

	Predicate device (k971573):	Device :
Device Name	N-geneous [®] LDL	ABX Pentra LDL Direct CP
	Cholesterol Reagent	
	5	
Controls	Commercially available quality control material (not included)	Recommended quality control material (not included): ABX Pentra N Control (Normal control) ABX Pentra P Control (Pathologic control)
Calibrators	Recommended calibration	Recommended calibration
	material (not included):	material (not included):
	N_geneous® LDL Cholesterol	ABX Pentra LDL Cal
	Calibrator	
Performance data :		
Number of tests	-	100 tests
Sample volume	3 µL/test	2.4 µL/test
Detection limit	-	1.55 mg/dL
Analytical sensitivity	0.278 mg/dL	-
Accuracy and Precision	CV Total < 2.27%	CV Total < 6.39%
	6.6 mg/dL – 992 mg/dL	1.35 mg/dL – 369.39 mg/dL
Measuring range		
	992 mg/dL	387 mg/dL
Upper linearity limit		
	-	12 days
Calibration stability	Until the expiration date when	18 months at 2 8°C
Closed reagent	stored at 2-8°C	18 months at 2-8 C
stability	stored at 2-6 C	
~~~~y	4 weeks at 2-8°C	on-board stability (refrigerated
		area): 97 days
Open Reagent stability		· · ·

	Predicate device (k893973):	Device :
Device Name	Triglycerides	ABX Pentra Triglycerides
		СР
Manufactured by	Roche, USA	HORIBA ABX, France
Instrument	COBAS chemistry system	ABX PENTRA 400
Analytes	Triglycerides	Triglycerides
Method :	Enzymatic colorimetric assay	Enzymatic colorimetric assay
Specimen :	Serum	Serum
	Plasma	Plasma
<b>Component reagent</b>	Single reagent bottles, the	Mono-reagent cassette, ready to
matrices	combined reagents contain the	use

	Predicate device (k893973):	Device :
Device Name	Triglycerides	ABX Pentra Triglycerides
		СР
	following components:	REAGENT : Pipes free acid,
	ATP, 4-Aminoantipyrine,	Sodium hydroxide, Triton X-
	Glycerol kinase, Glycerol	100, Magnesium salt, p-
	phosphate oxidase, Lipases,	chlorophenol, ATP, Sodium
	Peroxidase, Sodium N-Ethyl-N-	azide, Potassium ferrocyanide,
	(3-Sulfopropyl)-m-anisidine	4-Aminoantipyrine, Lipoprotein
	(ESPAS), buffers, stabilizers,	lipase, Glycerokinase, Glycerol
	fillers	phosphate Oxidase, Peroxidase
Format	Liquid	Liquid
Labels	-	Horiba ABX specific label
Packaging	Kit composed of single-reagent	Mono-reagent cassette :
0.0	bottles	REAGENT : 99 mL
	REAGENT 1 : 2 x 80 mL	
	REAGENT 2 : 2 x 40 ml	
Controls	Commercially available quality	Recommended quality control
	control (not included): 2 levels	material (not included):
		ABX Pentra N Control
		(Normal control)
		ABX Pentra P Control
		(Pathologic control)
Calibrators	Commercially available	Recommended calibration
	calibrator (not included)	material (not included):
		ABX Pentra Multical
Additional Reagents	-	Cleaning solution (not included):
		ABX Pentra Clean-Chem CP /
		ABX Pentra Clean-Chem 99 CP
Performance data :		
Number of tests	-	327 tests
Sample volume	4 µL/test	3 µL/test
Detection limit		7 ma/dI
Detection mint	-	/ mg/dL
A courses and	CV Total < 2.6%	CV Total < 2.83%
Accuracy and Provision	C V 10tal < 5.070	C V 10tal < 2.8570
	10  mg/dI 662 mg/dI	3.1  mg/dI $1/3/ mg/dI$
Measuring range		$5.1 \operatorname{Ing/ull} = 1454 \operatorname{Ing/ull}$
Wiedsur ing Fange	900 mg/dL (with automatic post-	1470  mg/dI with an
Upper linearity limit	dilution : 1800 mg/dL)	automatique post-dilution · 5580
opper mieunty mint		mg/dL
	_	
Calibration stability		14 days
	Until the expiration date when	
<b>Closed reagent</b>	stored at 2-8°C	16 months at 2-8°C
stability		
	Combined reagent is stable:	
	1 year at 2-8°C	on-board stability (refrigerated
Open Reagent	60 days at room temperature	area): 48 days

	Predicate device (k893973):	Device :
Device Name	Triglycerides	ABX Pentra Triglycerides
		СР
stability	(15-25°C)	

Device NameUltra N-geneous® HDL Cholesterol CalibratorABX Pentra HDL CalManufactured byGenzyme, USAHORIBA ABX, FranceInstrumentChemistry systemsABX PENTRA 400Method :Calibration of Genzyme Ultra N geneous HDL Cholesterol assayCalibration of HORIBA ABX HDL Cholesterol measurement methodComponent matricesVials (lyophilized)Vials (lyophilized)		Predicate device (k021316):	Device :
Cholesterol CalibratorManufactured byGenzyme, USAHORIBA ABX, FranceInstrumentChemistry systemsABX PENTRA 400Method :Calibration of Genzyme Ultra N- geneous HDL Cholesterol assayCalibration of HORIBA ABXComponent matricesVials (lyophilized)Vials (lyophilized)	Device Name	Ultra N-geneous® HDL	ABX Pentra HDL Cal
Manufactured byGenzyme, USAHORIBA ABX, FranceInstrumentChemistry systemsABX PENTRA 400Method :Calibration of Genzyme Ultra N geneous HDL Cholesterol assayCalibration of HORIBA ABX HDL Cholesterol measurement methodComponent matricesVials (lyophilized)Vials (lyophilized)		<b>Cholesterol Calibrator</b>	
InstrumentChemistry systemsABX PENTRA 400Method :Calibration of Genzyme Ultra N geneous HDL Cholesterol assayCalibration of HORIBA ABX HDL Cholesterol measurement methodComponent matricesVials (lyophilized)Vials (lyophilized)	Manufactured by	Genzyme, USA	HORIBA ABX, France
Method :     Calibration of Genzyme Ultra N- geneous HDL Cholesterol assay     Calibration of HORIBA ABX HDL Cholesterol measurement method       Component matrices     Vials (lyophilized)     Vials (lyophilized)	Instrument	Chemistry systems	ABX PENTRA 400
geneous HDL Cholesterol assay   HDL Cholesterol measurement method     Component matrices   Vials (lyophilized)   Vials (lyophilized)	Method :	Calibration of Genzyme Ultra N-	Calibration of HORIBA ABX
Component matrices       Vials (lyophilized)       method		geneous HDL Cholesterol assay	HDL Cholesterol measurement
Component matricesVials (lyophilized)Vials (lyophilized)			method
	<b>Component matrices</b>	Vials (lyophilized)	Vials (lyophilized)
Human serum based Human serum based, with		Human serum based	Human serum based, with
preservative			preservative
HDL Cholesterol HDL Cholesterol		HDL Cholesterol	HDL Cholesterol
Calibrated molecules	Calibrated molecules		
Format Lyophilized Lyophilized	Format	Lyophilized	Lyophilized
Labels - Horiba ABX specific label	Labels	-	Horiba ABX specific label
PackagingKit composed of :Kit composed of :	Packaging	Kit composed of :	Kit composed of :
1 x 1 mL vial 2 x 1 mL vial		1 x 1 mL vial	2 x 1 mL vial
Performance data :	Performance data :		
Calibration value- Assigned by procedures- Assigned by procedures	Calibration value	- Assigned by procedures	- Assigned by procedures
traceable to the National traceable to the National		traceable to the National	traceable to the National
Reference System for Reference System for		Reference System for	Reference System for
Cholesterol (NRS/CHOL) Cholesterol (NRS/CHOL)		Cholesterol (NRS/CHOL)	Cholesterol (NRS/CHOL)
I ne concentration of		-	- The concentration of
component is fot specific			component is lot specific
The concentration of The concentration is given in		The concentration of	The concentration is given in
component is given on the vial the enclosed appex		component is given on the vial	the enclosed anney
label		label	the enclosed annex
		laber	
<b>Closed stability</b> Up to the expiration date at 2- 24 months at 2-8°C	Closed stability	Up to the expiration date at 2-	24 months at 2-8°C
8°C	chosed stability	8°C	2 • montais at 2 • • •
<b>Open stability</b> 14 days at 2-8°C 14 days at 2-8°C	<b>Open stability</b>	14 days at 2-8°C	14 days at 2-8°C
4 weeks at $-70^{\circ}$ C 4 weeks at $-70^{\circ}$ C	- ⁻ V	4 weeks at -70°C	4 weeks at -70°C

	Predicate device (k971573):	Device :
<b>Device Name</b>	N-geneous [™] LDL	ABX Pentra LDL Cal
	Cholesterol Calibrator	
Manufactured by	Genzyme, USA	HORIBA ABX, France
Instrument	Chemistry systems	ABX PENTRA 400
Method :	Calibration of Genzyme N-	Calibration of HORIBA ABX

	Predicate device (k971573):	Device :
Device Name	N-geneous [™] LDL	ABX Pentra LDL Cal
	<b>Cholesterol Calibrator</b>	
	geneous LDL Cholesterol assay	LDL Cholesterol measurement
		method
<b>Component matrices</b>	Vials (lyophilized)	Vials (lyophilized)
	Human serum based	Human serum based, with
		preservative
Calibrated maleaulas	LDL Cholesterol	LDL Cholesterol
Calibrated molecules		
Format	Lyophilized	Lyophilized
Labels	-	Horiba ABX specific label
Packaging	Kit composed of :	Kit composed of :
	1 x 1 mL vial	2 x 1 mL vial
Performance data :		
Calibration value	A seize ad her and a drives	A saise ad her end as dunna
Cambration value	- Assigned by procedures	- Assigned by procedures
	Pafarance System for	Pafaranaa System for
	Cholesterol (NRS/CHOL)	Cholesterol (NRS/CHOL)
	cholesteror (NRS/CHOL)	cholesteror (NRS/CHOL)
	_	- The concentration of
		component is lot specific
		component is for specific
	- The concentration of	- The concentration is given in
	component is given on the vial	the enclosed annex
	label	
<b>Closed stability</b>	Up to the expiration date at 2-	24 months at 2-8°C
	8°C	
Open stability	2 weeks at 2-8°C	2 weeks at 2-8°C

# K. Standard/Guidance Document Referenced (if applicable):

CLSI - Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline-Second Edition - EP05-A2

CLSI - Evaluation of the Linearity of Quantitative Analytical Methods; Approved Guideline - EP06-A

CLSI - Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Edition - EP09-A2

# L. Test Principle:

# **ABX PENTRA Cholesterol CP**

"CHOD-PAP": enzymatic photometric test.

Determination of cholesterol after enzymatic hydrolysis and oxidation (3,4). The colorimetric indicator is quinoneimine which is generated from 4-aminoantipyrine and phenol by hydrogen peroxide under the catalytic action of peroxidase (Trinder's reaction)

# **ABX PENTRA HDL Direct CP**

The assay is a homogeneous method for directly measuring HDL-C levels in serum or plasma without the need for any off-line pretreatment or centrifugation steps.

The method is in a two reagent format and depends on the properties of a unique detergent, as illustrated. This method is based on accelerating the reaction of cholesterol oxidase (CO) with non-HDL unesterified cholesterol and dissolving HDL selectively using a specific detergent.

In the first reagent, non-HDL unesterified cholesterol is subject to an enzyme reaction and the peroxide generated is consumed by a peroxidase reaction with DSBmT yielding a colorless product.

The second reagent consists of a detergent capable of solubilizing HDL specifically, cholesterol esterase (CE) and chromagenic coupler to develop color for the quantitative determination of HDL-C.

# ABX PENTRA LDL Direct CP

The assay is an homogeneous method for directly measuring LDL-C levels in serum or plasma, without the need for any off-line pretreatment or centrifugation steps.

The method is in a two reagent format and depends on the properties of a unique detergent. This detergent (Reagent 1) solubilizes only the non LDL lipoprotein particles. The cholesterol released is consumed by cholesterol esterase and cholesterol oxidase in a non color forming reaction. A second detergent (Reagent 2) solubilizes the remaining LDL particles and a chromogenic coupler allows for color formation. The enzyme reaction with LDL-C in the presence of the coupler produces color that is proportional to the amount of LDL cholesterol present in the sample.

# **ABX PENTRA Triglycerides CP**

Enzymatic determination of triglycerides according to the following reactions:



(DHAP = Dihydroxyacetone phosphate, 4-AAP = 4-aminoantipyrine)

#### M. Performance Characteristics (if/when applicable):

- 1. <u>Analytical performance:</u>
  - a. Precision/Reproducibility: ABX PENTRA Cholesterol CP

	Mean value mg/dL	CV %
Normal control	113	0.82
Pathological control	186	0.74
Specimen 1	117	1.21
Specimen 2	191	0.53
Specimen 3	389	0.62

Two control materials and three specimens were tested 20times in a single run. The results are presented in the table below.

In accordance with CLSI guideline EP-5A, two control materials and two specimens were tested in duplicate for 20 days, two times a day (N=80). The results are presented in the table below.

	Mean value mg/dL	CV %
Normal control	109	2.96
Pathological control	183	2.34
Specimen 1	170	2.80
Specimen 2	250	3.01

# ABX PENTRA HDL Direct CP

Two control materials and three specimens were tested 20times in a single run. The results are presented in the table below.

	Mean value mg/dL	CV %
Normal control	35.82	1.29
Pathological control	81.72	0.79
Specimen 1	27.94	1.32
Specimen 2	48.59	1.91
Specimen 3	97.39	0.62

In accordance with CLSI guideline EP-5A, two control materials and two specimens were tested in duplicate for 20 days, two times a day (N=80). The results are presented in the table below.

	Mean value mg/dL	CV %
Normal control	35.83	2.88
Pathological control	80.35	3.06
Specimen 1	47.07	3.52
Specimen 2	80.16	2.69

# **ABX PENTRA LDL Direct CP**

Two control materials and three specimens were tested 20times in a single run. The results are presented in the table below.

	Mean value mg/dL	CV %
Normal control	61.26	1.01
Pathological control	75.08	2.82
Specimen 1	111.26	0.91
Specimen 2	141.45	1.00
Specimen 3	191.16	0.63

In accordance with CLSI guideline EP-5A, two control materials and two specimens were tested in duplicate for 20 days, two times a day (N=80). The results are presented in the table below.

	Mean value mg/dl	CV %
Normal control	60.64	5.59
Pathological control	74.27	6.39
Specimen 1	156.58	3.94
Specimen 2	191.62	4.04

#### **ABX PENTRA Triglycerides CP**

Two control materials and three specimens were tested 20times in a single run. The results are presented in the table below.

	Mean value mg/dL	CV %
Normal control	126	2.52
Pathological control	214	0.82
Specimen 1	60	2.83
Specimen 2	108	1.84
Specimen 3	232	1.00

In accordance with CLSI guideline EP-5A, two control materials and two specimens were tested in duplicate for 20 days, two times a day (N=80). The results are presented in the table below.

	Mean value mg/dL	CV %
Normal control	128	1.91
Pathological control	216	1.70
Specimen 1	132	1.57
Specimen 2	243	1.37

b. Linearity/assay reportable range:

The reagent linearity is determined according to the recommendations found in the CLSI (NCCLS), EP6-A protocol. As a result of the testing, the following represent the linear range for each assay.

#### **ABX PENTRA Cholesterol CP**

Low linearity: 4 mg/dL High linearity: 580 mg/dL

# **ABX PENTRA HDL Direct CP**

Low linearity: 1.16 mg/dL High linearity: 154.8 mg/dL.

# **ABX PENTRA LDL Direct CP**

Low linearity: 1.55 mg/dL High linearity: 387 mg/dL

# **ABX PENTRA Triglycerides CP**

Low linearity: 7 mg/dL High linearity: 1470 mg/dL, with automatic post-dilution: 5880 mg/dL

c. Traceability, Stability, Expected values (controls, calibrators, or methods): Calibrators for Cholesterol and Triglyceride, and controls for Cholesterol, HDL, LDL and Triglyceride see k052007.

#### ABX PENTRA and HDL, LDL Calibrators and L/H Controls

Real time stability at 2-8°C was investigated for Calibrator, High Controls and Low Controls.

Real time stability study on reconstituted material was performed at regular intervals over a period providing data to support claims at 2-8°C (All) and -20°C (Lipoproteins), and -70°C (HDL).

Traceability

#### ABX PENTRA Cholesterol CP Certification

Traceability to the National Reference System for Cholesterol was established by performing a direct comparison with the cholesterol reference method using human specimens that cover the National Cholesterol Education Program (NCEP) medical decision points. The ability to meet the NCEP's performance criteria for accuracy was demonstrated through testing and certification by the Cholesterol Reference Method Laboratory Network (CRMLN).

The results of the direct comparison and precision studies are as follows:

Among-run %CV : 1.5% Average %Bias : -2.5% Total Error : 5.6%

d. Detection limit:

# **ABX PENTRA Cholesterol CP:**

Minimum Detection Limit (MDL) is calculated from 30 measurements of saline. Formula: MDL = mean of measurements + 4.65 SD (mean of measurement = 0 when negative)

From the MDL value, the minimum detection limit assigned for this method is 4 mg/dL.

# **ABX PENTRA HDL Direct CP**

From the MDL value, the test range low assigned for this method is 1.16 mg/dL.

# **ABX PENTRA LDL Direct CP**

From the MDL value, the test range low assigned for this method is 1.55 mg/dL.

# **ABX PENTRA Triglycerides CP**

From the MDL value, the test range low assigned for this method is 7 mg/dL.

e. Analytical specificity:

# **ABX PENTRA Cholesterol CP**

Hemoglobin: No significant influence is observed up to 336 mg/dL Triglycerides: No significant influence is observed up to 612.5 mg/dL (as Intralipid®, representative of lipemia) Total Bilirubin: No significant influence is observed up to 20.5 mg/dL Direct Bilirubin: No significant influence is observed up to 6.8 mg/dL Other limitations are given by way of literature reference known to affect this methodology.

# **ABX PENTRA HDL Direct CP**

Hemoglobin: No significant influence is observed up to 479 mg/dl. Triglycerides: No significant influence is observed up to 612.5 mg/dL. (as Intralipid®, representative of lipemia) Total Bilirubin: No significant influence is observed up to 11.7 mg/dl.

Direct Bilirubin: No significant influence is observed up to 28.1 mg/dl.

# **ABX PENTRA LDL Direct CP**

Hemoglobin: No significant influence is observed up to 336 mg/dl. Triglycerides: No significant influence is observed up to 612.5 mg/dL. (as Intralipid®, representative of lipemia)

Total Bilirubin: No significant influence is observed up to 29.2 mg/dl. Direct Bilirubin: No significant influence is observed up to 10.8 mg/dl.

# **ABX PENTRA Triglycerides CP**

Hemoglobin: No significant influence is observed up to 500 mg/dl Total Bilirubin: No significant influence is observed up to 22.5 mg/dl Direct Bilirubin: No significant influence is observed up to 22.5 mg/dl

f. Assay cut-off: Not Applicable

# 2. Comparison studies:

a. Method comparison with predicate device:

# **ABX PENTRA Cholesterol CP**

135 patient samples (serum) are correlated with a commercial reagent taken as reference according to the recommendations found in the CLSI (NCCLS), EP9-A2 protocol. Values ranged from 2.55 to 583.26 mg/dL.

The equation for the regression line obtained is:

Y = 0.95 x + 1.90 with a correlation coefficient  $r^2 = 0.9943$ .

# **ABX PENTRA HDL Direct CP**

121 patient samples are correlated with a commercial reagent taken as reference according to the recommendations found in the CLSI (NCCLS), EP9-A2 protocol. Values ranged from 3.5 to 155.3 mg/dL The equation for the regression line obtained is:

Y = 0.91 x + 1.98 with a correlation coefficient  $r^2 = 0.9768$ 

# **ABX PENTRA LDL Direct CP**

122 patient samples are correlated with a commercial reagent taken as reference according to the recommendations found in the CLSI (NCCLS), EP9-A2 protocol. Values ranged from 0.3 to 385.6 mg/dL The equation for the regression line obtained is:

Y = 0.96 x - 0.21, with a correlation coefficient  $r^2 = 0.9963$ 

# **ABX PENTRA Triglycerides CP**

135 patient samples are correlated with a commercial reagent taken as reference according to the recommendations found in the CLSI (NCCLS), EP9-A2 protocol. Values ranged from 8.2 to 274.2 mg/dL The equation for the regression line obtained is:

Y = 0.99 x + 0.20, with a correlation coefficient  $r^2 = 0.9994$ 

# b. Matrix comparison:

# **ABX PENTRA Cholesterol CP**

To demonstrate equivalence of Cholesterol results in serum and Plasma Heparin-Lithium samples, comparison study was performed. 40 samples (103-373 mg/dL) were evaluated on Pentra 400 analyzer using ABX Pentra Cholesterol CP reagent.

The equation for the regression line obtained is:

Y = 0.9903 x + 1.8467 with a correlation coefficient r = 0.996

# **ABX PENTRA HDL Direct CP**

To demonstrate equivalence of HDL Cholesterol results in serum and Plasma Heparin-Lithium samples, comparison study was performed. 43 samples (35-118 mg/dL) were evaluated on Pentra 400 analyzer using ABX Pentra HDL Direct CP reagent.

The equation for the regression line obtained is:

Y = 0.9467 x + 1.5252 with a correlation coefficient r = 0.996

# **ABX PENTRA LDL Direct CP**

To demonstrate equivalence of LDL Cholesterol results in serum and Plasma Heparin-Lithium samples, comparison study was performed. 41 samples (53-362 mg/dL) were evaluated on Pentra 400 analyzer using ABX Pentra LDL Direct CP reagent.

The equation for the regression line obtained is:

Y = 0.9681 x + 0.2352 with a correlation coefficient r = 0.999

# ABX PENTRA Triglycerides CP

To demonstrate equivalence of Triglycerides results in serum and Plasma Heparin-Lithium samples, comparison study was performed. 41 samples (47-728 mg/dL) were evaluated on Pentra 400 analyzer using ABX Pentra Triglycerides CP reagent. The equation for the regression line obtained is:

Y = 0.9694 x - 1.7717 with a correlation coefficient r = 0.999

- 3. <u>Clinical studies</u>:
  - *a. Clinical Sensitivity:* Not Applicable
  - *b. Clinical specificity:* Not Applicable
  - c. Other clinical supportive data (when a. and b. are not applicable): Not Applicable
- 4. <u>Clinical cut-off:</u> Not Applicable
- 5. Expected values/Reference range: ABX PENTRA Cholesterol CP Desirable: < 200 mg/dL Borderline high risk: 200 - 239 mg/dL High risk: > 240 mg/dL Cholesterol Education Program (NCEP) medical decision points

# ABX PENTRA HDL Direct CP

Men: 30 - 70 mg/dL Women: 30 - 85 mg/dL (literature) According to the NCEP, HDL values greater than or equal to 40mg/dL are considered desirable, and values greater than or equal to 60mg/dL are considered to offer some protection against coronary heart disease. Values below 40 mg/dL are considered to be a significant independent risk factor for coronary heart disease.

# ABX PENTRA LDL Direct CP

The following NCEP cut-points for patient classification are used for the prevention and management of coronary heart disease.

LDL CholesterolClassification< 130 mg/dL</td>Desirable130-159 mg/dLBorderline High Risk160 mg/dLHigh Risk

# ABX PENTRA Triglycerides CP

The following NCEP cut- points been classified according to the risk of developing cardiovascular diseases: Normal: < 150 mg/dL Low risk: 150 - 200 mg/dL High: 200 - 500 mg/dL Extremely high: > 500 mg/dL

# N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

# **O.** Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.