# 510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY ONLY TEMPLATE

#### **A.** 510(k) Number:

k072078

#### **B.** Purpose for Submission:

Addition of the urine matrix to the assay procedure

#### C. Measurand:

Beta-2-microglobulin

## **D.** Type of Test:

Latex particle enhanced immunoturbidimetric assay

# E. Applicant:

Biokit S.A.

# F. Proprietary and Established Names:

Quantia Beta-2 Microglobulin

# **G. Regulatory Information:**

<b>Product Code</b>	Classification	<b>Regulation Section</b>	Panel
System, Test, Beta-	Class II	21 CFR 866.5630,	82 IMMUNOLOGY
2-Microglobulin		Beta-2-	(IM)
Immunological		microglobulin	
(JZG)		immunological test	
		system.	

#### H. Intended Use:

#### 1. Intended use(s):

The Quantia Beta-2 Microglobulin is intended as a latex particle enhanced immunoturbidimetric assay for the *in vitro* quantitative determination of beta-2-microglobulin concentration in human serum, plasma (EDTA) or urine on the AEROSET ® Instrument as an aid in the diagnosis of active rheumatoid arthritis and kidney disease.

The Quantia Beta-2 Microglobulin is intended to be used with the already cleared Quantia PROTEINS Control (k050596) and the Beta-2 Microglobulin Standard (k050613).

## 2. Indication(s) for use:

Same as Intended Use

# 3. Special conditions for use statement(s):

For prescription use only

## 4. Special instrument requirements:

The reagents are for use on the Abbott AEROSET® instrument

## I. Device Description:

The Quantia Beta-2 Microglobulin kit contains 4 bottles of Reagent 1 (R1) (6 mL each) and 4 bottles of Reagent 2 (R2) (3 mL each). R1 buffer is sodium dihydrogen

phosphate dihydrate with polyethylene glycol and preservative (sodium azide). R2 is a suspension of polystyrene latex particles of uniform size coated with the IgG fraction of rabbit anti-human Beta-2-microglobulin specific serum with preservative (sodium azide).

## J. Substantial Equivalence Information:

#### Predicate: IL Test Beta-2-Microglobulin K943686

#### Similarities:

The Quantia Beta-2 Microglobulin and the IL Test Beta-2-Microglobulin are both manufactured by Biokit and are both intended for the quantitative *in vitro* diagnostic determination of beta-2-microglubulin. They also use the same methodology: Particle Enhanced Immunoturbidimetry.

The Quantia Beta-2 Microglobulin and the IL Test Beta-2- Microglobulin, have the same composition: Latex Reagent Suspension of polystyrene latex particles coated with rabbit IgG anti-human Beta-2 Microglobulin in a buffer containing bovine serum albumin and < 0.1 % w/w sodium azide; and

Reaction Buffer Phosphate buffer 40mM containing bovine serum albumin and sodium azide < 0.1 % w/w.

#### Differences:

Specimen type: both assays can use serum and urine as samples. They differ in the plasma types. The device can use EDTA and the predicate can use EDTA and sodium heparin.

# K. Standard/Guidance Document Referenced (if applicable):

#### **STANDARDS**

#### **Title and Reference Number**

CLSI Interference Testing in Clinical Chemistry; Approved Guideline (EP 7-A)

CLSI Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline (EP09-A2)

CLSI Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline (EP5-A)

# L. Test Principle:

When a sample containing Beta-2 Microglobulin is mixed with the reagent, a clear agglutination occurs which can be measured by turbidimetry. Results are expressed in mg/L of Beta-2-microglobulin based on the 1<sup>st</sup> WHO International Standard (B2M) established in 1985.

## M. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
  - a. Precision/Reproducibility: CLSI EP5-A was followed.

Samples/Runs	Mean (mg/L)	CV (%) Within Run	CV (%) Total
2/40	0.066	4.2	5.4
2/40	0.094	1.7	3.1
2/40	0.204	1.5	2.6
2/40	0.302	1.6	2.3

## b. Linearity/assay reportable range:

Linearity was assessed according to CLSI EP6-A. The overall reportable range for serum, plasma and urine is 0.025 to 96 mg/L. The assay was linear from 0.025 to 1.6 mg/L (the reportable range for urine samples) with the automatic rerun capability (Dilution Protocol 2); from 0.25 to 16 mg/L without the automatic rerun capability; and from16 to 96 mg/L with the automatic rerun capability (Dilution Protocol 1).

#### Prozone

The manufacturer was asked to run a sample higher than 114 mg/L to demonstrate the instrument give a result of >16 mg/L. They tested samples from 80.4 to 200.9 mg/L. In all cases the instrument gave a result of > 96 (serum linearity upper limit) and triggered the Dilution Protocol 1. The assay did not demonstrate prozone effect with specimens up to 200 mg/L.

c. Traceability, Stability, Expected values (controls, calibrators, or methods): The assay calibrators are standardized against WHO reference material B2M. Stability of the reagents was established at 19 months by testing 4 different lots. Reagents should be stored at 2-8°C.

#### d. Detection limit:

The limit of quantitation (LOQ) was defined as the minimum quantity of analyte that can be measured with a within-run CV below 20% and an error below 20%. The LOQ for the assay was determined to be 0.25 mg/L without the automatic rerun capability and 0.025 mg/L with the automatic rerun capability for serum, plasma and urine. This was established by running serial dilutions of the 4 mg/L calibrator. Through 0.025 mg/L the %CV ranged from 1.0 to 5.5% and error ranged from 3.8 to 9.3%.

The limit of detection (LOD) was defined as the mean reported value + 2SD for a sample free of analyte. The LOD was determined to be 0.046 mg/L without the rerun capability and 0.025 mg/L (LOQ) with the automatic rerun capability.

## e. Analytical specificity:

CLSI guideline EP7-A was followed.

Substance tested	Concentration (mg/dL)	Outcome
Conjugated bilirubin	20.9	No interference
Ascorbic acid	20	Interference below 10%
Hemoglobin	23.6	Interference below 10%
Protein (IgG)	100	No interference

Urine pH showed no significant positive or negative influence on the result.

## f. Assay cut-off:

See Expected values/Reference range

# 2. <u>Comparison studies:</u>

# a. Method comparison with predicate device:

The comparison studies were performed using 110 urine samples with values ranging from 0.01 to 18.85 mg/L. The required specifications were: Slope 1.0  $\pm$  0.20; and correlation coefficient  $r \ge 0.950$ . The following results were obtained:

#### **AEROSET versus ILab 900**

Parameter	Outcome
Slope	1.088 (95% CI: 1.061 to 1.127)
Intercept	-0.001 (95% CI: -0.003 to 0.002)
Range (mg/L)	0.01- 18.85
Mean x (mg/L)	2.47
Mean y (mg/L)	2.467
R	0.9894

# b. Matrix comparison:

Urine was the only matrix compared.

#### 3. Clinical studies:

a. Clinical Sensitivity:

Not determined

b. Clinical specificity:

Not determined

c. Other clinical supportive data (when a. and b. are not applicable):

Not applicable

# 4. Clinical cut-off:

Not determined

# 5. Expected values/Reference range:

Concentrations of Beta-2-microglobulin in urine from healthy subjects averaged 0.098 mg/L with an upper normal limit of 0.32 mg/L (literature).

## N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

## O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantially equivalent decision.