510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY AND INSTRUMENT COMBINATION TEMPLATE

A. 510(k) Number:

k103788

B. Purpose for Submission:

New device

C. Measurand:

Sodium, potassium, chloride, and glucose

D. Type of Test:

Quantitative, hexokinase method for glucose and potentiometric method for sodium, potassium, and chloride

E. Applicant:

HORIBA ABX

F. Proprietary and Established Names:

PENTRA C200,,I.S.E. Module, and ABX PENTRA Glucose HK CP

G. Regulatory Information:

Product Code	Classification	Regulation Section	Panel
CFR	Class II	21 CFR §862.1345: Glucose test	Chemistry
		system.	(75)
JGS	Class II	21 CFR §862.1665: Sodium test	Chemistry
		system.	(75)
CEM	Class II	21 CFR §862.1600: Potassium test	Chemistry
		system.	(75)
CGZ	Class II	21 CFR §862.1170: Chloride test	Chemistry
		system.	(75)
JJE	Class I	21 CFR §862.2160: Discrete	Chemistry
		Photometric Chemistry Analyzer for	(75)
		Clinical Use.	

H. Intended Use:

1. Intended use(s):

See indications for use below

2. Indication(s) for use:

The PENTRA C200 is a discrete photometric benchtop chemistry analyzer for use in clinical laboratories. It is not intended for use in Point of Care Settings. It duplicates manual analytical procedures by performing various steps such as pipetting, mixing, heating and measuring color intensity. The PENTRA C200 is intended for quantitative measurements of a variety of analytes: Glucose, Sodium, Potassium, and Chloride.

ABX Pentra Glucose HK CP reagent with associated calibrators and controls are for quantitative in vitro diagnostic determination of glucose in serum and plasma using glucose hexokinase method by colorimetry. Glucose measurements are used on the diagnosis and treatment of carbohydrate metabolism disorders including diabetes mellitus, neonatal hypoglycemia, and idiopathic hypoglycemia, and of pancreatic islet cell carcinoma.

The I.S.E. (Ion Selective Electrode) module is intended for the quantitative determination of Sodium, Potassium and Chloride in serum, plasma, and urine by potentiometry using ion selective electrode with associated reference solution, calibrators and controls. Measurement of Sodium, Potassium and Chloride are used in diagnosis and treatment diseases involving electrolyte imbalance.

3. Special conditions for use statement(s):

For prescription use only It is not intended for use in Point of Care settings

4. Special instrument requirements:

PENTRA C200 Clinical Chemistry Analyzer

I. Device Description:

The PENTRA C200 is a bench-top clinical chemistry analyzer using two measuring principles: absorbance and ion selective electrodes.

The instrument may be summarized as follows:

- Multi-parametric (up to 15 simultaneous tests + 3 ISE tests)
- On routine or Stat
- 90 (without ISE) to 360 tests / hour (with ISE) (analytical cycle of 40 seconds)
- random access working on primary tubes or sample cups

- Reagent cassettes are compact and ready-to-use
- On-board bar-code readers are used to identify newly loaded reagent cassettes and samples for patient identification.

ISE module - PENTRA C200 includes an optional Ion Selective Electrode module which contains a sodium electrode, a potassium electrode, a chloride electrode, a reference electrode, two calibrators and two controls. The direct application is intended for serum and plasma, indirect application is intended for urine.

Glucose reagent contains a bi-reagent cassette, ready to use, reagent 1 (R1) and reagent 2 (R2). R1 reagent contains NAD, ATP, buffer, and sodium azide. R2 contains hexokinase, G-6-PDH, magnesium sulphate, and sodium azide.

J. Substantial Equivalence Information:

1. Predicate device name(s):

HORIBA ABX PENTRA 400

2. Predicate 510(k) number(s):

k052007

3. Comparison with predicate:

For instrument and ISE: Similarities and Differences:

	Candidate device	Predicate device (k052007)
	PENTRA C200	ABX PENTRA 400
Intended Use	The PENTRA C200 is intended for quantitative measurements for Glucose, Sodium, Potassium, and Chloride.	Same
Instrument Type	Bench-top	Same
Touch Screen Interface	Yes	Same
I.S.E module:		
ISE parameters	Na, K, Cl	Same
Methods	Direct & indirect	Same
Sample volume in μL	30 μL for 3 parameters	60 μL for 3 parameters
Sodium Electrode	Glass membrane selective to	Same

	NT .	
	Na+ ions	_
Potassium Electrode	Plastic membrane selective	Same
	to K+ ions	
Chloride Electrode	Plastic membrane selective	Same
	to Cl- ions	
Type of samples	Serum, plasma, urine	Same
Throughput		
Without I.S.E.	Up to 360 tests/hour	Up to 420 tests/hour
With I.S.E	270 tests/hour	180 tests/hour
Parameters on board	15 cassettes + 3 ISE	52 mono or twin cassettes +
		3 ISE
		0 102
Sample capacity	Up to 15 sample tube or	Up to 60 (6 racks of 10
Sample capacity	cups (handling of STAT	samples)
	available)	samples)
Sample volume	2 to 45µL	2 to 290 mJ
Sample volume	•	2 to 380 μL
Dilution of patient sample	Yes	Yes
Fibrin detecting device	Yes	Yes
D		
Reagents:	C	T · · · 1
Type of reagents	Same	Liquid
Disposable/Washable	Same	Disposable cuvettes
cuvettes		
D: :1 C	C	
Principles of	Same	Spectrophotometry:
Measurement		Colorimetry and
		Turbidimetry:
		parallel bi-chromatic
		measurement of light
		absorbance
		(tungsten halogen lamp)
		Potentiometry : Direct
		(Serum or
		Plasma) and Indirect (Urine)
Photometer type: light		
source	Same	Halogen
Type of measurement	Same	Mono, Bi chromatic
Dimensions (HxWxD)	75.5 x 72.5 x 57 cm	100 x 65 x 57 cm
Weight	95 kg	120 kg

For glucose assay: Similarities and Differences:

	Candidate device	Predicate device (k052007)
Intended Use	ABX Pentra Glucose HK CP reagent with associated calibrators and controls are for quantitative in vitro diagnostic determination of glucose in serum and plasma using glucose hexokinase method by colorimetry.	Same
Instrument	ABX PENTRA C200	ABX PENTRA 400
Analyte	Glucose	Glucose
Method:	Same	Enzymatic method using hexokinase coupled with glucose-6-phosphate dehydrogenase
Sample type:	Serum Plasma (Lithium heparin, Fluoride oxalate)	Serum Plasma
Reagent component	Same	Bi-reagent cassette, ready to use REAGENT 1: NAD, ATP, Buffer, Sodium azide REAGENT 2: hexokinase, G-6- PDH, Magnesium sulphate, Sodium azide
Format	Same	Liquid
Packaging	R1-R2: (30 mL/10 mL) R1-R2: (50 mL/50 mL) R1-R2: (80 mL/10 mL)	Bi-reagent cassette : R1 : 46 mL R2 : 12 mL
Calibrators	Same	The ABX PENTRA Multical
Controls	Same	The ABX PENTRA N control and P control

K. Standard/Guidance Document Referenced (if applicable):

C28-A3: Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory

EN13640: Stability testing of in vitro diagnostic reagents, August 2002

EP05-A2: Evaluation of Precision Performance of Quantitative Measurement Methods

EP06-A: Evaluation of the Linearity of Quantitative Measurement

EP07-A2: Interference Testing in Clinical Chemistry

EP09-A2: Method Comparison and Bias Estimation Using Patient Samples

EP17-A: Protocols for Determination of Limits of Detection and Limits of

Quantitation

L. Test Principle:

There are two different measuring principles employed, potentiometry and photometry/colorimetry.

Potentiometry for sodium potassium and chloride: Electrical potential measured between the reference electrode flowed by a reference solution and the specific selective electrode flowed by the sample. Slopes of the electrodes are determined with two standard solutions of known concentrations and stored by the instrument.

Colorimetry for glucose: Glucose measurement employed the hexokinase method coupled with glucose-6-phosphate dehydrogenase. The amount of absorbance change is proportional to the concentration of the glucose being present in the sample.

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. Precision/Reproducibility:

A precision study was performed to assess the with-run precision and total precision of the PENTRA C200 system by replicate measurements of controls materials and human samples (Serum, lithium heparin plasma, and urine). The precision study was performed in duplicates, twice a day, for twenty days. The precision results were summarized below:

For glucose (serum): within-run precision (N=20)

	Mean value	CV%
	mg/dL	
Serum control 1	91	0.76
Serum control 2	244	0.75
Serum sample 1	40	1.81
Serum sample 2	88	0.51
Serum sample 3	314	0.61

For glucose (serum): total precision (N=80)

	Mean value	CV%
	mg/dL	
Serum control 1	93	1.99
Serum control 2	246	1.60
Serum sample 1	41	1.81
Serum sample 2	86	1.58
Serum sample 3	304	1.40

For sodium (serum and plasma): within-run precision (N=20)

	Mean value mmol/L	CV%
Serum control 1	137.57	0.13
Serum control 2	152.62	0.15
Serum sample 1	105.40	0.57
Serum sample 2	144.87	0.18
Serum sample 3	176.81	0.14
Plasma sample 1	117.12	0.18
Plasma sample 2	143.50	0.15
Plasma sample 3	176.71	0.13

For sodium (urine): within-run precision (N=20)

	Mean value	CV%
	mmol/L	
Urine control 1	101.04	2.22
Urine control 2	179.11	1.28
Urine sample 1	77.16	2.46
Urine sample 2	114.73	1.50
Urine sample 3	231.60	0.79

For sodium (serum): total precision (N=80)

	Mean value	CV%
	mmol/L	
Serum control 1	130.02	0.88
Serum control 2	148.49	1.10
Serum sample 1	107.85	0.80
Serum sample 2	138.93	0.62
Serum sample 3	151.11	0.78

For sodium (urine): total precision (N=80)

	Mean value	CV%
	mmol/L	
Urine control 1	76.96	4.39
Urine control 2	158.07	2.88
Urine sample 1	119.65	4.91
Urine sample 2	141.62	3.14

For potassium (serum and plasma): within-run precision (N=20)

	Mean value	CV%
	mmol/L	
Serum control 1	3.58	0.25
Serum control 2	6.29	0.35
Serum sample 1	3.56	0.85
Serum sample 2	4.17	0.55
Serum sample 3	5.15	0.85
Plasma sample 1	2.37	0.57
Plasma sample 2	3.96	0.70
Plasma sample 3	7.16	0.75

For potassium (urine): within-run precision (N=20)

	Mean value	CV%
	mmol/L	
Urine control 1	32.70	0.88
Urine control 2	68.96	0.81
Urine sample 1	26.90	1.03
Urine sample 2	99.56	0.91
Urine sample 3	207.59	0.42

For potassium (serum): total precision (N=80)

	Mean value CV%	
	mmol/L	
Serum control 1	3.56	0.87
Serum control 2	6.43	1.07
Serum sample 1	4.02	0.86
Serum sample 2	4.65	0.73
Serum sample 3	4.79	0.85

For potassium (urine): total precision (N=80)

	Mean value CV%	
	mmol/L	
Urine control 1	28.75	1.72
Urine control 2	62.84	1.64
Urine sample 1	31.15	1.54
Urine sample 2	106.05	1.67
Urine sample 3	61.37	2.87

For chloride (serum and plasma): within-run precision (N=20)

	Mean value	CV%
	mmol/L	
Serum control 1	89.54	0.26
Serum control 2	116.25	0.23
Serum sample 1	81.02	0.96
Serum sample 2	109.52	0.51
Serum sample 3	145.23	0.76
Plasma sample 1	83.76	0.74
Plasma sample 2	110.78	0.63
Plasma sample 3	143.57	0.55

For chloride (urine): within-run precision (N=20)

	Mean value CV%	
	mmol/L	
Urine control 1	98.64	2.74
Urine control 2	176.08	1.71
Urine sample 1	81.86	2.35
Urine sample 2	148.30	2.13
Urine sample 3	174.52	1.26

For chloride (serum): total precision (N=80)

	Mean value CV%	
	mmol/L	
Urine control 1	85.88	1.20
Urine control 2	113.36	1.55
Serum sample 1	83.26	0.99
Serum sample 2	104.72	0.96
Serum sample 3	115.32	0.88

For chloride (urine): total precision (N=80)

	Mean value CV%	
	mmol/L	
Urine control 1	98.26	4.59
Urine control 2	172.06	1.56
Urine sample 1	116.68	4.02
Urine sample 2	169.16	3.12

b. Linearity/assay reportable range:

A linearity study was performed according to the recommendations found in the CLSI EP6-A guideline. Serum or urine samples containing the target analytes (glucose, sodium, potassium, and chloride) in different concentrations were prepared. Each sample was run in four replicates per level. The observed values were plotted against the expected values and an appropriate line fitted by standard linear regression. Results were summarized in the tables below:

Table 1: Serum samples

	Sodium	Potassium	Chloride	Glucose
	(serum)	(serum)	(serum)	(serum)
Correlation (r)	0.9997	0.9996	0.9995	0.9999
Slope	1.0273	0.9630	1.0210	0.9160
Intercept	- 2.753	0.1613	-2.1902	2.7394
Range tested	79.4-204.1	1.01-10.13	67.3-173.0	1.5-933
	(mmol/L)	(mmol/L)	(mmol/L)	mg/dL
Levels tested	9	10	10	13

Table 2: Urine samples

	Sodium	Potassium	Chloride
	(urine)	(urine)	(urine)
Correlation (r)	0.9996	0.9996	0.9915
Slope	0.9979	0.9982	0.9585
Intercept	0.3203	-1.0224	12.2347
Range	57.6-329.9	21.48-	41.87-366.32
	mmol/L	261.53	mmol/L
		(mmol/L)	
Levels tested	10	10	10

Results of the study support the sponsor claims for the following measuring ranges:

Analyte	Measuring range
Glucose (serum)	5-900 mg/dL
Sodium (serum)	90 – 190 mmol/L
Sodium (urine)	60 – 280 mmol/L
Potassium (serum)	2-9.5 mmol/L
Potassium (urine)	25 – 250 mmol/L
Chloride (serum)	70 – 170 mmol/L
Chloride (urine)	70 – 280 mmol/L

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Traceability of all the analytes:

Analyte	Traceability
Sodium	Gravimetric method using standard
	materials
Potassium	Gravimetric method using standard
	materials
Chloride	Gravimetric method using standard
	materials
Glucose	Reference material: NIST SRM 965a

The ABX Pentra Multical is the calibrator used for calibration for the glucose assay and has been previously cleared (k052007). The ABX Pentra Standard 1 and 2 are calibrators used for calibration for the sodium, potassium, and chloride assays and have been previously cleared (k052007).

The glucose assay has a calibration stability of 20 days. The glucose reagent has a shelf-life of 36 months (when stored at 2 to 8°C) and an on-board (refrigerated) stability of 39 days. Protocols and acceptance criteria has been provided and were found to be adequate.

d. Detection limit:

See linearity study in M.2.b. above.

In addition, the sponsor performed a limit of detection study for the glucose assay according to a modified protocol based on the CLSI EP17-A guideline. Limit of blank (LoB) was determined by testing 90 times of a blank sample on 3 different instruments. Limit of detection (LoD) was determined by testing four low serum samples 20 times on one instrument. Limit of quantitation

(LoQ) was determined by testing three low samples 5 times on three instruments over 5 days. LoQ is defined as the lowest concentration where total precision is < 10 %CV. The sponsor determined that for the glucose assay, the LoB is 0.72 mg/dL, LoD is 1.08 mg/dL, and LoQ is 4.9 mg/dL.

The sponsor claims for the following measuring ranges:

Analyte	Measuring range
Glucose (serum)	5-900 mg/dL
Sodium (serum)	90 – 190 mmol/L
Sodium (urine)	60 – 280 mmol/L
Potassium (serum)	2-9.5 mmol/L
Potassium (urine)	25 – 250 mmol/L
Chloride (serum)	70 – 170 mmol/L
Chloride (urine)	70 – 280 mmol/L

e. Analytical specificity:

A study of common interfering substances was conducted and evaluated following recommendations by CLSI EP7A2. Two sample pools were prepared at the following target concentrations:

Analyte	Target Concentrations	
	Low level	High level
Glucose (serum)	90 mg/dL	234 mg/dL
Sodium (serum)	120 mmol/L	150 mmol/L
Sodium (urine)	80 mmol/L	150 mmol/L
Potassium (serum)	2 mmol/L	6 mmol/L
Potassium (urine)	60 mmol/L	150 mmol/L
Chloride (serum)	90 mmol/L	110 mmol/L
Chloride (urine)	100 mmol/L	200 mmol/L

Various concentrations of interfering substances were spiked into the serum and urine sample pools. The sponsor states that interferences are considered to be non-significant if the bias between the spiked and non-spiked samples is within \pm 10%. Results are summarized in the following:

I. Glucose (serum):

Hemoglobin: No significant interference observed up to 603 mg/dL Triglycerides: No significant interference observed up to 10,500 mg/dL Total Bilirubin: No significant interference observed up to 24.4 mg/dL Direct Bilirubin: No significant interference observed up to 37.7 mg/dL

Protein: No significant interference observed up to 12 g/dL Intralipid: No significant interference observed up to 640 mg/dL Acetylsalicyclic acid: No significant interference observed up to 65 mg/dL Bicarbonate: No significant interference observed up to 40 mmol/L

II. Sodium, potassium, and chloride (serum)

Hemoglobin: No significant interference observed up to 200 mg/dL The sponsor included a limitation in the labeling that hemolyzed samples should not be used.

Triglycerides: No significant interference observed up to 740 mg/dL Total Bilirubin: No significant interference observed up to 19.9 mg/dL Urea: No significant interference observed up to 258 mg/dL Protein: No significant interference observed up to 12 g/dL Acetylsalicyclic acid: No significant interference observed up to 65 mg/dL L-Glutathione reduced: No significant interference observed up to 92.2 mg/dL

Methyl Dopa: No significant interference observed up to 16.9 mg/dL Lithium: No significant interference observed up to 118 mg/dL Bicarbonate: No significant interference observed up to 50 mmol/L

For sodium (serum):

Ammonium Nitrate: No significant influence observed up to 40 mmol/L Ammonium Bromide: No significant influence observed up to 37.5 mmol/L

Probenecid: No significant influence observed up to 2100 µmol/L

For potassium (serum):

Ammonium Nitrate: No significant influence observed up to 40 mmol/L Ammonium Bromide: No significant influence observed up to 37.5 mmol/L.

Probenecid: Significant interference observed from 500 µmol/L

For chloride (serum):

Probenecid: No significant influence observed up to $1200~\mu mol/L$ Ammonium Nitrate: Significant interference observed from 0.2~mmol/L Ammonium Bromide: Significant interference observed from 1.25~mmol/L

III. Sodium, potassium, and chloride (urine)

Hemoglobin: No significant interference observed up to 1000 mg/dL Total Bilirubin: No significant interference observed up to 8.8 mg/dL Urea: No significant interference observed up to 3600 mg/dL

Protein: No significant interference observed up to 0.2 g/dL Ascorbic Acid: No significant interference observed up to 60 mg/dL Boric acid: No significant interference observed up to 867 mg/dL

Based on the significant interference substances identified above, the sponsor states the following in their limitations section in the labeling of the sodium, potassium and chloride:

- 1. Do not use hemolyzed samples.
- 2. Do not use samples containing Probenecid acid, ammonium nitrate and ammonium bromide.

f. Assay cut-off:

Not Applicable

2. Comparison studies:

a. Method comparison with predicate device:

Method comparison studies were conducted following CLSI EP9-A2. This study was conducted for glucose (serum), chloride (serum, plasma, and urine), potassium (serum, plasma and urine), and sodium (serum, plasma, and urine). Plasma samples were collected from lithium heparin plasma samples. The method comparison studies were performed on the ABX Pentra C200 and compared to the predicate device (ABX Pentra 400). Left over clinical laboratory serum, plasma, and urine samples covering the measuring range were analyzed in duplicate and a singlet set of data were used for linear regression analysis. To cover the high and low ends of the measuring range, some specimens were spiked with stock solutions and others were diluted with saline (glucose samples) and DI water (sodium, potassium, and chloride samples). Deming regression results were summarized in the table below:

Summary of test results using Deming Regressions are as follows:

Analyte	Sample type	Slope	Intercept	r	N	Sample range tested
Glucose (mg/dL)	serum	0.98	0.248	0.9920	103	7.02-821.52 mg/dL
Sodium (mmol/L)	serum	0.96	6.42	0.9909	129	103.4-166.5 mmol/L
	plasma	1.05	-5.32	0.9989	132	93.4-165.9 mmol/L
	urine	1.01	-2.20	0.9946	101	67.64-273.61 mmol/L
Potassium (mmol/L)	serum	1.01	-0.06	0.9989	122	2.52-8.32 mmol/L
,	plasma	1.01	-0.09	0.9992	125	2.24-9.47 mmol/L
	urine	1.02	-0.27	0.9987	159	25.22-226.70 mmol/L
Chloride (mmol/L)	serum	0.96	3.74	0.9980	170	74.37-166.92 mmol/L
	plasma	1.04	-4.17	0.9984	131	71.15-144.56 mmol/L
	urine	1.04	-5.63	0.9933	112	73.05-262.46 mmol/L

b. Matrix comparison:

Anticoagulated plasma and serum were evaluated to demonstrate equivalency of these matrices for measurement of glucose, sodium, potassium, and chloride. The following anticoagulants were tested in the study: lithium heparin and fluoride oxalate for the glucose assay and lithium heparin for the sodium, potassium and chloride assays. Matched samples for serum and plasma were collected for each comparison. All samples were assayed in duplicate over several days on the PENTRA C200. Results of the study are summarized in the following tables:

Glucose serum vs. Lithium heparin plasma

Number of samples	40
Slope	1.0460
Intercept	-0.36158
R	0.999
Sample range tested	47.52-821.9 mg/dL

Glucose serum vs fluoride oxalate plasma

Number of samples	49
Slope	0.9963
Intercept	-0.8242
R	0.999
Sample range tested	59.04-787.50 mg/dL

Sodium serum vs lithium heparin plasma

Number of samples	26
Slope	1.0162
Intercept	-2.5289
R	0.975
Sample range tested	128.20-141.30 mmol/L

Potassium serum vs lithium heparin plasma

1 0 000001001111 0 0 1 001111 1 0 1 1 1 1	man pulling properties
Number of samples	26
Slope	0.9829
Intercept	0.0307
R	0.999
Sample range tested	3.33-9.06 mmol/L

Chloride serum/lithium heparin plasma

Number of samples	26
Slope	1.0166
Intercept	-1.8639
R	0.993
Sample range tested	93.55-110.95 mmol/L

The sponsor determined that lithium heparin is an acceptable anticoagulant for sodium, potassium, chloride and glucose. In addition, fluoride oxalate is also acceptable for glucose measurement.

3. Clinical studies:

- a. Clinical Sensitivity:Not Applicable
- b. Clinical specificity: Not Applicable
- c. Other clinical supportive data (when a. and b. are not applicable): Not Applicable
- 4. Clinical cut-off:

Not Applicable

5. Expected values/Reference range:

Glucose⁽¹⁾

Serum: 70 - 115 mg/dL

Chloride (2),(3)

Serum:101 – 110 mmol/L Urine:110-250 mmol/24h

Potassium (2), (3)

Serum:3.7 – 5.5 mmol/L Urine:25 – 125 mmol/24h

Sodium (2), (3)

Serum: 136 – 145 mmol/L Urine: 40-220 mmol/24h

- 1. Thomas L., Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results. 1st ed and Frankfurt: TH Books Verlagsgesellschaft; 1998 p.132
- 2. Tietz, Fundamentals of Clinical Chemistry, 4th ed., p.983-990
- 3. Tietz, Fundamentals of Clinical Chemistry, 5th ed., p. 970-1009

N. Instrument Name:

PENTRA C200

O. System Descriptions:

1. Modes of Operation:

Routine or Stat random access

2. Software:

FDA	has rev	iewed appl	licant's Hazard Analysis and software development
proce	esses for	this line of	of product types:
Yes	X	or No	

3. Specimen Identification:

Real Patient sample bar-code ID

4. Specimen Sampling and Handling:

Primary tubes or sample cups

5 Calibration:

ISE calibration is stable for 8 hours and needs to be calibrated at the beginning of each day and every 8 hours.

6. Quality Control:

At least two levels of controls should be assayed daily and after a calibration. In the labeling the sponsor states that the user should follow the state, local and state guidelines for testing quality control materials.

P. Other Supportive Instrument Performance Characteristics Data Not Covered In The "Performance Characteristics" Section above:

A dilution study was provided to show that the instrument can automatically diluted sample with a dilution factor of 3. Six serum samples (spiked) with glucose concentrations range from 900 mg/dL to 2700 mg/dL were manually diluted X3 and assayed on the Pentra C200 instrument. Results were compared with the automatic dilution protocol of the Pentra C200 instrument. Results showed that the instrument can automatically dilute the serum samples accurately.

An ambient temperature study on the Pentra C200 was performed to demonstrate that ambient temperature between 15 to 30°C will not have a significant impact on the performance (ISE and glucose) of the Pentra C200. The study protocol and acceptance criteria was found to be adequate.

Q. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

R. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.