## SPECIAL 510(k): Device Modification OIR Review Memorandum

To: Princeton BioMeditech Corporation RE: K132465

This 510(k) submission contains information/data on modifications made to the SUBMITTER'S own Class II, Class III or Class I devices requiring 510(k). The following items are present and acceptable (delete/add items as necessary):

1. The name and 510(k) number of the SUBMITTER'S previously cleared device:

BioSign Flu A+B

510(k) number: K083746

- 2. Submitter's statement that the **INDICATION/INTENDED USE** of the modified device as described in its labeling **HAS NOT CHANGED** along with the proposed labeling which includes instructions for use, package labeling, and, if available, advertisements or promotional materials (labeling changes are permitted as long as they do not affect the intended use).
- 3. Description of the device **MODIFICATION(S)**:

The modification presented in this special 510(k) consisted of expanded reactivity table to include reactivity information for one H7N9 influenza A virus (A/Anhui/1/2013), four (4) H3N2v viruses (A/Indiana/10/2011, A/Indiana/08/2011, A/Minnesota/10/2011, A/Minnesota/10/2011 X-203), and an influenza B virus (B/Texas/39/2006). The firm tested the ability of the BioSign Flu A+B test to detect the six aforementioned viruses. The viruses used were obtained from the Centers for Disease Control and Prevention as non-infectious beta-propiolactone inactivated virus. A LoD study was performed with each of the viruses using the same procedure employed in the original submission. Each titered virus was diluted until the minimal visual signal intensity appeared on the test line. This was defined as the lowest reacting level of the virus. Each virus was then tested in triplicate at that dilution. All virus strains tested were detected in 3 out of 3 tests at the lowest reacting level. The empirically determined LoD's for each virus are listed below:

- A/Anhui/1/2013 (H7N9) 7.94 x 10<sup>6</sup> EID<sub>50</sub>/mL
- A/Indiana/10/2011 (H3N2v) 2.34 x 10<sup>3</sup> TCID<sub>50</sub>/mL
- A/Indiana/08/2011 (H3N2v) 2.87 x 10<sup>6</sup> TCID<sub>50</sub>/mL
- A/Minnesota/10/2011 (H3N2v) 2.13 x 10<sup>6</sup> TCID<sub>50</sub>/mL
- A/Minnesota/10/2011 X-203 (H3N2v) 2.28 x 10<sup>3</sup> TCID<sub>50</sub>/mL
- B/Texas/39/2006 2.34 x 10<sup>4</sup> TCID<sub>50</sub>/mL

The BioSign Flu A+B test and Status Flu A&B test package inserts have been updated to include the additional analytical reactivity information. Status Flu A&B is the name of the same device being sold by LifeSign LLC under agreement with Princeton BioMeditech Corporation.

4. The FUNDAMENTAL SCIENTIFIC TECHNOLOGY of the modified device has not changed.

## 5. Comparison Information

## **Similarities**

	Modified Device	Predicate Device		
Features	BioSign Flu A+B test	BioSign Flu A+B test		
	The BioSign Flu A+B test is an in	The BioSign Flu A+B test is an in		
	vitro rapid qualitative test that detects	vitro rapid qualitative test that detects		
	influenza type A and type B	influenza type A and type B antigens		
	nucleoprotein antigens directly from	directly from nasal swab,		
	nasal swab, nasopharyngeal swab,	nasopharyngeal swab, and		
	and nasopharyngeal aspirate/wash	nasopharyngeal aspirate/wash		
	specimens obtained from patients	specimens of patients with signs and		
	with signs and symptoms of	symptoms of respiratory infection. It		
	respiratory infection. It is intended to	is intended to aid in the rapid		
	aid in the rapid differential diagnosis	differential diagnosis of influenza A		
	of influenza A and B viral infections.	and B viral infections. Negative test		
	Negative test results are presumptive	results are presumptive and it is		
	and it is recommended these results	recommended these results be		
	be confirmed by viral culture.	confirmed by viral culture. Negative		
	Negative results do not preclude influenza virus infection and should	results do not preclude influenza virus infection and should not be		
	not be used as the sole basis for	used as the sole basis for treatment		
	treatment or other patient	or other patient management		
	management decisions. The test is	decisions. The test is intended for		
	intended for professional and	professional and laboratory use.		
	laboratory use. Performance	Performance characteristics for		
	characteristics for influenza were	influenza were established during the		
	established during the 2007-2009	2007-2009 influenza seasons when		
late a de d l le e	influenza seasons when influenza A	influenza A viruses New		
Intended Use	viruses A/New Caledonia/20/99	Caledonia/20/99 (H1N1), Solomon		
	(H1N1), A/Solomon Islands/3/2006	Islands/3/2006 (H1N1),		
	(H1N1), A/Brisbane/59/2007 (H1N1),	Brisbane/59/2007 (H1N1),		
	A/California/07/2009 (H1N1),	California/07/2009 (H1N1),		
	A/Wisconsin/67/2005 (H3N2),	A/Wisconsin/67/2005 (H3N2),		
	A/Brisbane/10/2007 (H3N2), and	A/Brisbane/10/2007 (H3N2), and		
	influenza B viruses B/Ohio/01/2005,	influenza B viruses Ohio/01/2005,		
	B/Florida/4/2006,	Florida/4/2006, Brisbane/60/2008		
	B/Brisbane/60/2008 were the	were the predominant influenza		
	predominant influenza viruses in circulation according to the Flu	viruses in circulation according to the Flu Activity & Surveillance report by		
	Activity & Surveillance report by	CDC. Performance characteristics		
	CDC. Performance characteristics	may vary against other emerging		
	may vary against other emerging	influenza viruses. If infection with a		
	influenza viruses. If infection with a	novel Influenza virus is suspected		
	novel Influenza virus is suspected	based on current clinical and		
	based on current clinical and	epidemiological screening criteria		
	epidemiological screening criteria	recommended by public health		
	recommended by public health	authorities, specimens should be		
	authorities, specimens should be	collected with appropriate infection		
	collected with appropriate infection	control precautions for novel virulent		
	control precautions for novel virulent	Influenza viruses and sent to state or		
	Influenza viruses and sent to state or	local health department for testing.		
	local health department for testing.	Viral culture should not be attempted		

	Viral culture should not be attempted	in these cases unless a BSL+3		
	in these cases unless a BSL+3	facility is available to receive and		
	facility is available to receive and	culture specimens.		
	culture specimens.			
Sample Type		Nasal swab		
	Same as predicate device	Nasopharyngeal swab		
		Nasopharyngeal aspirate/wash		
Analytical Principle	Samo as pradicata davias	Solid phase chromatographic		
	Same as predicate device	immunoassay		
Extraction	Same as predicate device	Incubated 1 minute in extraction		
		reagent		
Read Result Time	Same as predicate device	10 Minutes		
Test Line	Same as predicate device	Colloidal gold		
Internal Control	Same as predicate device	Reddish-purple line		
Control Samples		Positive Control Swab: Influenza		
		A and B antigens (non-infective		
		recombinant nucleoprotein)		
	Same as predicate device	Negative Control Swab:		
		Inactivated Group B		
		Streptococcus antigen (non-		
		infective)		

#### **Differences**

The package insert has been updated to include detection of the A/Anhui/1/2013, A/Indiana/10/2011, A/Indiana/08/2011, A/Minnesota/10/2011, A/Minnesota/10/2011 X-203, and B/Texas/39/2006 viruses at the following limits of detection:

- A/Anhui/1/2013 (H7N9) 7.94 x 10<sup>6</sup> EID<sub>50</sub>/mL
- A/Indiana/10/2011 (H3N2v) 2.34 x 10<sup>3</sup> TCID<sub>50</sub>/mL
- A/Indiana/08/2011 (H3N2v) 2.87 x 10<sup>6</sup> TCID<sub>50</sub>/mL
- A/Minnesota/10/2011 (H3N2v) 2.13 x 10<sup>6</sup> TCID<sub>50</sub>/mL
- A/Minnesota/10/2011 X-203 (H3N2v) 2.28 x 10<sup>3</sup> TCID<sub>50</sub>/mL
- B/Texas/39/2006 2.34 x 10<sup>4</sup> TCID<sub>50</sub>/mL

Although this test has been shown to detect these H7N9, H3N2v, and type B/Texas/39/2006 viruses cultured from positive human respiratory specimens, the performance characteristics of this device with clinical specimens that are positive for these influenza viruses have not been established.

Two minor grammatical changes were made within the Intended Use, but those changes do not alter the meaning of the IU.

#### 6. **Design Control Activities Summary**:

Analytical reactivity testing was conducted for the H7N9 virus, four H3N2v viruses, and the influenza B virus using identical methods employed in the original submission for the unmodified device.

The risk analysis method used to assess the impact of the modification, adding additional viruses to the analytical sensitivity section of the package insert, was Failure Modes and Effects Analysis (FMEA). Based on the result of the risk analysis, the verification activities required and acceptance criteria were

identified. Since the change is adding detection levels of additional strains without changing anything in the test device, including fundamental scientific technology or indications for use, no risk is involved for this change except as listed below:

Change	Hazard	Resolution of Risk	Testing Performed	Test Method	Acceptance Criteria	Acceptance Criteria Met?
Addition of new virus strains to package insert	Non- detection of the virus strain added	Confirm Analytical Sensitivity for all additional strains	Analytical Sensitivity Testing conducted for each of the added strains	Tested in triplicate for each dilution of each additional strain	Positive Results at 10 minutes for each virus at the determined analytical sensitivity	* Yes
	Misinterpreta tion of test use: test used for detection of the additional strains from human specimen	Labeling: Limitation of test added as a footnote below the inclusivity table for all additional strains	n/a	n/a	n/a	n/a

<sup>\*</sup> This was indicated after all experiments were completed.

A "Declaration of Conformity" statement was submitted for the manufacturing facility and validation activities and signed by the Regulatory Affairs manager and the Quality Assurance manager. The statements indicate that:

- 1. The manufacturing facility is in conformance with design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review.
- 2. The validation activities, as required by the risk analysis, for the modification were performed by the designated individuals and the results demonstrated that the predetermined acceptance criteria were met.

In conclusion, based on the results of the analytical reactivity testing the modified labeling is truthful and accurate. The changes do not affect the performance of the test and it is therefore substantially equivalent to the current cleared test.

# 7. Truthful and Accurate Statement, a 510(k) Summary or Statement and the Indications for Use Enclosure.

The labeling for this modified subject device has been reviewed to verify that the indication/intended use for the device is unaffected by the modification. In addition, the submitter's description of the particular modification(s) and the comparative information between the modified and unmodified devices demonstrate that the fundamental scientific technology has not changed. The submitter has provided the design control information as specified in The New 510(k) Paradigm and on this basis, I recommend the device be determined substantially equivalent to the previously cleared device.