

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY TEMPLATE**

A. 510(k) Number:

k180427

B. Purpose for Submission:

New device

C. Measurand:

Fentanyl

D. Type of Test:

Qualitative, homogeneous enzyme immunoassay

E. Applicant:

ARK Diagnostics, Inc.

F. Proprietary and Established Names:

ARK Fentanyl Assay

G. Regulatory Information:

Product Code	Classification	Regulation Section	Panel
DJG	Class II	21 CFR 862.3650, Opiate test system	Toxicology (91)

H. Intended Use:

1. Intended use(s):

Refer to Indications for Use.

2. Indication(s) for use:

The ARK Fentanyl Assay is an immunoassay intended for the qualitative determination of fentanyl in human urine at a cutoff concentration of 1.0

ng/mL. The assay is intended for use in laboratories with automated clinical chemistry analyzers. This in vitro diagnostic device is for prescription use only.

The ARK Fentanyl Assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/tandem Mass Spectrometry (LC-MS/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be exercised with any drug test result, particularly when the preliminary test result is positive.

3. Special conditions for use statement(s):

For prescription use only.

4. Special instrument requirements:

Performance data was collected on the Beckman Coulter AU680 analyzer.

I. Device Description:

The ARK Fentanyl Assay is a homogeneous enzyme immunoassay used for the detection of fentanyl in human urine. The ARK Fentanyl Assay consists of reagents R1 rabbit anti-fentanyl polyclonal antibody with substrate and R2 fentanyl derivative labeled with bacterial recombinant G6PDH enzyme. The ARK Fentanyl Assay test system includes separately provided kits for the ARK Fentanyl Assay, ARK Fentanyl Calibrator and ARK Fentanyl Control.

J. Substantial Equivalence Information:

1. Predicate device name(s):

Immunoassay SEFRIA Fentanyl Urine Enzyme Immunoassay

2. Predicate 510(k) number(s):

k161216

3. Comparison with predicate:

Similarities		
Item	Candidate Device	Predicate device – k161216
Intended Use	Same	For the determination of fentanyl in human urine at a cutoff concentration of 1.0 ng/mL.
Methodology	Same	Homogenous enzyme immunoassay (EIA)
Sample Matrix	Same	Urine
User Environment	Same	Clinical laboratories; Prescription use only
Platform Required	Same	Automated clinical chemistry analyzer
Reagents Form	Same	Liquid – Ready to use
Reagent Materials	Same	Two (2) reagent system: Antibody/substrate reagent and enzyme labeled conjugate Sodium azide preservative
Stability	Same	2-8°C until expiration date
Antibody	Same	Rabbit antibodies to fentanyl
Cutoff	Same	1.0 ng/mL

Differences		
Item	Candidate Device	Predicate device – k161216
Detection Wavelength	Absorbance change measured spectrophotometrically at 340 nm.	Absorbance change measured spectrophotometrically at 570 nm.

K. Standard/Guidance Document Referenced (if applicable):

- CLSI EP07-A2: Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition
- CLSI EP12-A2: User Protocol For Evaluation Of Qualitative Test Performance; Approved Guideline - Second Edition

L. Test Principle:

The ARK Fentanyl Assay is a homogeneous enzyme immunoassay technique used for the analysis of fentanyl in human urine. The assay is based on competition between fentanyl in the patient specimen and fentanyl labeled with recombinant glucose-6-phosphate dehydrogenase (rG6PDH) for antibody binding sites. As the latter binds antibody, enzyme activity decreases. In the presence of fentanyl from

the patient specimen, enzyme activity increases and is proportional to the fentanyl concentration. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH in the presence of glucose-6-phosphate (G6P), resulting in an absorbance change that is measured spectrophotometrically.

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. Precision/Reproducibility:

A precision study data was conducted over 20 non-consecutive days, with two runs per day in replicates of four, on a single Beckman Coulter AU680 automated clinical chemistry analyzer. Samples were prepared in a human urine pool negative for fentanyl.

Results are summarized below:

Human Urine (ng/mL)	Relative % Cutoff	# of Results	Results
0.00	-100	160	160 Negative
0.25	-75	160	160 Negative
0.50	-50	160	160 Negative
0.75	-25	160	160 Negative
1.00	Cutoff	160	97 Negative; 63 Positive
1.25	25	160	160 Positive
1.50	50	160	160 Positive
1.75	75	160	160 Positive
2.00	100	160	160 Positive

b. Linearity/assay reportable range:

Not applicable. This device is intended for qualitative use only.

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

ARK Fentanyl Calibrators are traceable to a commercially available certified standard solution, and are prepared by volumetric dilution, with concentrations verified by high-performance liquid chromatography.

d. Detection limit:

Not applicable.

e. Analytical specificity:

The sponsor evaluated the cross-reactivity of the following metabolites and structural analogs of Fentanyl by spiking them into drug-free, negative human urine to determine the minimum concentration that would give a positive result approximately equivalent to the 1.0 ng/mL fentanyl cutoff. These concentrations were used to determine the percent cross-reactivity according to the formula:

$$\text{Percent cross-reactivity} = (\text{Cutoff concentration} / \text{lowest concentration of cross-reactant causing a positive result}) \times 100$$

Cross-reactivity:

Norfentanyl (major metabolite)

Compound	Concentration tested (ng/mL)	Percent Cross reactivity
Norfentanyl (major metabolite)*	2.5	10
	300	0.33

*Cross-reactivity of norfentanyl was observed to be non-linear

Other Metabolites and Structural Analogs of Fentanyl

Compound	Lowest concentration that produces a response approximately equivalent to the cutoff (ng/mL)	Percent Cross reactivity
Acetyl fentanyl	1.20	83.33
Acrylfentanyl	1.20	83.33
ω -1-Hydroxyfentanyl	1.20	83.33
Isobutyryl fentanyl	1.50	66.67
Ocfentanil	1.50	66.67
Butyryl fentanyl	1.60	62.50
Furanyl fentanyl	1.75	57.14
Valeryl fentanyl	2.50	40.00
(\pm) β -hydroxythiofentanyl	2.80	35.71

Compound	Lowest concentration that produces a response approximately equivalent to the cutoff (ng/mL)	Percent Cross reactivity
4-Fluoro-isobutyryl fentanyl	3.00	33.33
Para-fluorobutyryl fentanyl (p-FBF)	3.00	33.33
Para-fluoro fentanyl	3.00	33.33
(±)-3-cis-methyl fentanyl	5.00	20.00
Despropionyl fentanyl (4- ANPP)	75.00	1.33
Carfentanil	500	0.20
Sufentanil	625	0.16
Norcarfentanil	5,000	<0.02
Acetyl norfentanyl	10,000	0.01
Remifentanil	10,000	<0.01
Alfentanil	100,000	<0.001

The sponsor also provided results from a study evaluating whether other structurally similar compounds cross-reacted with the candidate device. None of the compounds below caused a positive result at the concentrations indicated.

Compound	Concentration Tested (ng/mL)
6-Acetyl morphine	10,000
Amphetamine	100,000
Buprenorphine	100,000
Buprenorphine glucuronide	50,000
Codeine	100,000
Dextromethorphan	100,000
Dihydrocodeine	100,000
EDDP	100,000
EMDP	50,000
Fluoxetine	50,000
Heroin	30,000

Compound	Concentration Tested (ng/mL)
Hydrocodone	100,000
Hydromorphone	100,000
Ketamine	100,000
Levorphanol	50,000
Meperidine	100,000
Methadone	100,000
Morphine	100,000
Morphine-3-glucuronide	50,000
Naloxone	50,000
Naltrexone	50,000
Norbuprenorphine	50,000
Norcodeine	50,000
Norketamine	100,000
Normeperidine	100,000
Normorphine	50,000
Noroxycodone	100,000
Oxycodone	100,000
Oxymorphone	50,000
Pentazocine (Talwin)	10,000
Pipamperone	100,000
Risperidone	2,000
Tapentadol	50,000
Thioridazine	50,000
Tilidine	50,000
Tramadol	100,000
Tramadol-O-Desmethyl	100,000
Tramadol-N-Desmethyl	100,000
Trazodone	10,000

The sponsor performed a study to evaluate whether common endogenous compounds cause negative or positive interference. These compounds were spiked, at the concentrations listed below, into urine samples containing fentanyl at concentrations corresponding to -50% and +50% of the cutoff. None of the compounds tested caused a positive result at -50% of the cutoff or a negative result at +50% of the cutoff.

Compound	Concentration Tested (mg/dL)	Candidate Device Result	
		0.5 ng/mL Fentanyl	1.5 ng/mL Fentanyl
Acetone	1000	NEG	POS
Ascorbic Acid	560	NEG	POS
Bilirubin	2.0	NEG	POS
Creatinine	500	NEG	POS
Ethanol	1000	NEG	POS
Galactose	10	NEG	POS
Gamma globulin	500	NEG	POS
Glucose	3000	NEG	POS
Hemoglobin	500	NEG	POS
Human Albumin	500	NEG	POS
Oxalic Acid	100	NEG	POS
Riboflavin	7.5	NEG	POS
Sodium Chloride	4000	NEG	POS
Urea	2000	NEG	POS

The sponsor performed a study to evaluate whether low or high pH causes negative or positive interference with the candidate device. Samples containing fentanyl concentrations corresponding to -50% and +50% of the cutoff were adjusted to pH levels ranging from 3.0 to 11.0. None of the pH levels tested caused a positive result at -50% of the cutoff or a negative result at +50% of the cutoff.

Compound Tested	Candidate Device Result	
	0.5 ng/mL Fentanyl (POS/NEG)	1.5 ng/mL Fentanyl (POS/NEG)
Urine pH 3.0	NEG	POS
Urine pH 4.0	NEG	POS
Urine pH 5.0	NEG	POS
Urine pH 6.0	NEG	POS
Urine pH 7.0	NEG	POS
Urine pH 8.0	NEG	POS
Urine pH 9.0	NEG	POS
Urine pH 10.0	NEG	POS
Urine pH 11.0	NEG	POS

The sponsor performed a study to evaluate whether low or high specific gravity causes negative or positive interference. Samples containing

fentanyl concentrations corresponding to -50% and +50% of the cutoff were adjusted to specific gravity levels ranging from 1.001 to 1.030. None of the specific gravity levels tested caused a positive result at -50% of the cutoff or a negative result at +50% of the cutoff.

Specific Gravity Tested	Candidate Device Result	
	0.5 ng/ml Fentanyl N=3 (POS/NEG)	1.5 ng/ml Fentanyl N=3 (POS/NEG)
1.001	NEG	POS
1.002	NEG	POS
1.004	NEG	POS
1.010	NEG	POS
1.018	NEG	POS
1.024	NEG	POS
1.026	NEG	POS
1.028	NEG	POS
1.030	NEG	POS

The sponsor evaluated boric acid and determined that it causes false negative results. The use of boric acid as a preservative is contraindicated in the device labeling.

Compound	Concentration Tested	Candidate Device Result	
		0.5 ng/ml Fentanyl Mean of N=3 (ng/mL)	1.5 ng/ml Fentanyl Mean of N=3 (ng/mL)
Boric Acid	1% w/v	NEG	NEG

The sponsor performed a study to evaluate whether structurally unrelated compounds cause negative or positive interference with the candidate device. These compounds at the concentrations listed in the table below were spiked into urine samples containing fentanyl concentrations corresponding to -50% and +50% of the cutoff. None of the compounds tested caused a positive result at -50% of the cutoff or a negative result at +50% of the cutoff.

Compound	Concentration Tested (µg/mL)	Candidate Device Result	
		0.5 ng/mL Fentanyl (POS/NEG)	1.5 ng/mL Fentanyl (POS/NEG)
Acetaminophen	500	NEG	POS
Acetylsalicylic acid	1000	NEG	POS
Albuterol	100	NEG	POS
Amitriptyline	35	NEG	POS
Amobarbital	100	NEG	POS
Benzoylcegonine	100	NEG	POS
Bupropion	50	NEG	POS
Caffeine	100	NEG	POS
Carbamazepine	100	NEG	POS
Chlorpromazine	50	NEG	POS
Clomipramine	50	NEG	POS
Cyclobenzaprine	10	NEG	POS
Desipramine	50	NEG	POS
Doxepin	50	NEG	POS
Ecgonine	100	NEG	POS
Ephedrine	100	NEG	POS
Fluphenazine	100	NEG	POS
Ibuprofen	500	NEG	POS
Imipramine	30	NEG	POS
Lidocaine	50	NEG	POS
Maprotiline	50	NEG	POS
Methapyrilene	10	NEG	POS
Methaqualone	50	NEG	POS
Metronidazole	300	NEG	POS
Nicotine	10	NEG	POS
Nortriptyline	25	NEG	POS
Oxazepam	100	NEG	POS
Phencyclidine	100	NEG	POS
Phenobarbital	100	NEG	POS
Propoxyphene	50	NEG	POS
Ranitidine	100	NEG	POS
Secobarbital	100	NEG	POS
Valproic acid	250	NEG	POS
Venlafaxine	100	NEG	POS

f. Assay cut-off:

Characterization of how the device performs analytically around the claimed cutoff concentration of 1 ng/mL fentanyl is described in the

precision section, M.1.a. above.

2. Comparison studies:

a. *Method comparison with predicate device:*

100 confirmed fentanyl-positive and 50 confirmed fentanyl-negative unaltered clinical urine specimens were analyzed with the ARK Fentanyl Assay in single replicates on a Beckman Coulter AU680 analyzer and compared to results obtained by LC-MS/MS. Results are summarized in the following table:

Candidate Device Result	Low Negative Less than 50% below the Cutoff (< 0.5 ng/mL by LC-MS/MS)	Near Cutoff Negative Between 50% below the Cutoff and the Cutoff (0.5 – 0.9 ng/mL by LC-MS/MS)	Near Cutoff Positive Between the Cutoff and 50% above the Cutoff (1.0 – 1.5 ng/mL by LC-MS/MS)	High Positive Greater than 50% above the Cutoff (> 1.5 ng/mL by LC-MS/MS)
Positive	1*	20	12	64
Negative	50	3	0	0

Discordant results:

*Norfentanyl was detected in this discordant sample and contributed to the positive result obtained by the candidate device for this sample.

Sample Number	Immunoassay result	Fentanyl (ng/mL by LC-MS/MS)	Sample Number	Immunoassay result	Fentanyl (ng/mL by LC-MS/MS)
051	Positive	0.7	062	Positive	0.6
052*	Positive	0.4	063	Positive	0.9
053	Positive	0.9	064	Positive	0.8
054	Positive	0.9	065	Positive	0.5
055	Positive	0.6	066	Positive	0.7
056	Positive	0.6	069	Positive	0.5
057	Positive	0.9	070	Positive	0.6
058	Positive	0.5	072	Positive	0.6

Sample Number	Immunoassay result	Fentanyl (ng/mL by LC-MS/MS)	Sample Number	Immunoassay result	Fentanyl (ng/mL by LC-MS/MS)
059	Positive	0.9	073	Positive	0.8
060	Positive	0.5	074	Positive	0.6
061	Positive	0.9			

b. Matrix comparison:

Not applicable. The device is intended to be used with human urine samples only.

3. Clinical studies:

a. Clinical Sensitivity:

Not applicable.

b. Clinical specificity:

Not applicable.

c. Other clinical supportive data (when a. and b. are not applicable):

Not applicable.

4. Clinical cut-off:

Not applicable.

5. Expected values/Reference range:

Not applicable.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Parts 801 and 809, as applicable.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.