



510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY

I Background Information:

A 510(k) Number

K230956

B Applicant

BD Integrated Diagnostic Solutions

C Proprietary and Established Names

BD Respiratory Viral Panel (BD RVP) for BD MAX System; BD Respiratory Viral Panel-SCV2 (BD RVP-SCV2) for BD MAX System

D Regulatory Information

Product Code(s)	Classification	Regulation Section	Panel
QOF	Class II	21 CFR 866.3981 - Device To Detect And Identify Nucleic Acid Targets In Respiratory Specimens From Microbial Agents That Cause The SARS-Cov-2 Respiratory Infection And Other Microbial Agents When In A Multi-Target Test	MI - Microbiology
QQX	Class II	21 CFR 866.3981 - Device to detect and identify nucleic acid targets in respiratory specimens from microbial agents that cause the SARS-CoV-2 respiratory infection and other microbial agents when in a multi-target test	MI - Microbiology

II Submission/Device Overview:

A Purpose for Submission:

BD Respiratory Viral Panel (BD RVP)

The purpose of this submission is to demonstrate that the BD Respiratory Viral Panel (SARS-CoV-2, Influenza A, Influenza B and RSV) for use on the BD MAX System, is substantially equivalent to the BioFire Respiratory Panel 2.1 (RP2.1) (DEN200031) and to obtain clearance for the BD Respiratory Viral Panel (SARS-CoV-2, Influenza A, Influenza B and RSV) Assay.

BD Respiratory Viral Panel-SCV2 (BD RVP-SCV2)

The purpose of this submission is to demonstrate that the BD Respiratory Viral Panel-SCV2 (SARS-CoV-2) for use on the BD MAX System, is substantially equivalent to the BioFire Respiratory Panel 2.1 (RP2.1) (DEN200031) and to obtain clearance for the BD Respiratory Viral Panel-SCV2 (SARS-CoV-2) Assay.

B Measurand:

BD Respiratory Viral Panel (BD RVP)

The BD Respiratory Viral Panel detects and identifies nucleic acids from the following pathogens: Severe Acute Respiratory Syndrome (SARS)-Coronavirus-2 (SARS-CoV-2), Influenza A, Influenza B, and Respiratory Syncytial Virus (RSV).

For each specimen, the amplification and detection of SARS-CoV-2, influenza A, influenza B, RSV and the internal control, human RNase P, occurs in a single reaction with the use of multiplexed primers and probes. The multiplexed primers and probes target RNA from the nucleocapsid phosphoprotein gene (N1 and N2 regions) of SARS-CoV-2, matrix (M1) gene of influenza A, matrix (M1) gene and hemagglutinin (HA) gene of influenza B, N gene and M gene from RSV and the RNase P gene from the human genome.

A total of eighteen (18) primers are involved in the DNA amplification process and nine (9) molecular probes are involved in the detection process of the BD Respiratory Viral Panel for BD MAX System.

BD Respiratory Viral Panel-SCV2 (BD RVP-SCV2)

The BD Respiratory Viral Panel-SCV2 detects and identifies nucleic acids from the following pathogens: Severe Acute Respiratory Syndrome (SARS)-Coronavirus-2 (SARS-CoV-2).

The BD Respiratory Viral Panel-SCV2 is produced using the same manufacturing processes and assay reagents as the Respiratory Viral Panel, but requires a different Assay Definition File (ADF) which, via software-mediated masking, only reports results for the SARS-CoV-2 target.

C Type of Test:

Multiplex nucleic acid assay for use with the BD MAX System for the qualitative detection of viral pathogens in individuals with signs and symptoms of respiratory tract infection consistent with COVID-19, Influenza A, Influenza B, or RSV.

III Intended Use/Indications for Use:

A Intended Use(s):

BD Respiratory Viral Panel for BD MAX System:

BD Respiratory Viral Panel for BD MAX System is an automated multiplexed real-time reverse transcriptase polymerase chain reaction (RT-PCR) test intended for the simultaneous, qualitative detection and differentiation of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), influenza A, influenza B, and/or respiratory syncytial virus (RSV) nucleic acid in nasopharyngeal swab (NPS) and anterior nasal swab (ANS) specimens from individuals with signs and symptoms of respiratory tract infection. Clinical signs and symptoms of respiratory tract infection due to SARS-CoV-2, influenza, and RSV can be similar.

BD Respiratory Viral Panel for BD MAX System is intended for use as an aid in the differential diagnosis of SARS-CoV-2, influenza A, influenza B, and/or RSV infection if used in conjunction with other clinical and epidemiological information, and laboratory findings. SARS-CoV-2, influenza A, influenza B, and RSV viral nucleic acid are generally detectable in NPS and ANS specimens during the acute phase of infection.

Positive results do not rule out co-infection with other organisms. The agent(s) detected by the BD Respiratory Viral Panel for BD MAX System may not be the definitive cause of disease.

Negative results do not preclude SARS-CoV-2, influenza A, influenza B, and/or RSV infection.

The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

BD Respiratory Viral Panel-SCV2 for BD MAX System:

BD Respiratory Viral Panel-SCV2 for BD MAX System is an automated multiplexed real-time reverse transcriptase polymerase chain reaction (RT-PCR) test intended for the simultaneous, qualitative detection of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) viral nucleic acid in nasopharyngeal swab (NPS) and anterior nasal swab (ANS) specimens from individuals with signs and symptoms of respiratory tract infection. SARS-CoV-2 viral RNA is generally detectable in NPS and ANS specimens during the acute phase of infection.

The BD Respiratory Viral Panel-SCV2 for BD MAX System is intended for use as an aid in the diagnosis of SARS-CoV-2 infection if used in conjunction with other clinical and epidemiological information, and laboratory findings.

Positive results do not rule out co-infection with other organisms. The agent detected by the BD Respiratory Viral Panel-SCV2 for BD MAX System may not be the definitive cause of disease.

Negative results do not preclude SARS-CoV-2 infection.

The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.

B Indication(s) for Use:

See Intended Use(s)

C Special Conditions for Use Statement(s):

Rx - For Prescription Use Only

For in vitro diagnostic use only

D Special Instrument Requirements:

The BD Respiratory Viral Panel Assay and BD Respiratory Viral Panel-SCV2 Assay are performed on the BD MAX System.

IV Device/System Characteristics:

A Device Description:

BD Respiratory Viral Panel (BD RVP) and BD Respiratory Viral Panel-SCV2 (BD RVP-SCV2)

The BD Respiratory Viral Panel (BD RVP) and BD Respiratory Viral Panel-SCV2 (BD RVP-SCV2) along with the BD MAX System comprise an instrument with associated hardware and accessories, disposable microfluidic cartridges, master mixes, unitized reagent strips, and extraction reagents. The instrument automates sample preparation including target lysis, Total Nucleic Acid (TNA) extraction and concentration, reagent rehydration, target nucleic acid amplification and detection using real-time PCR. The assay includes a Sample Processing Control (SPC) that is present in the Extraction Tube. The SPC monitors RNA extraction steps, thermal cycling steps, reagent integrity and the presence of inhibitory substances. The BD MAX System software automatically interprets test results.

For the BD Respiratory Viral Panel for BD MAX System the analytes reported are SARS-CoV-2, Influenza A, Influenza B and RSV. For the BD Respiratory Viral Panel-SCV2 for BD MAX System, the analyte is SARS-CoV-2. The difference in analyte(s) reported is determined by an assay definition file (ADF).

B Principle of Operation:

Test Principle

The BD Respiratory Viral Panel for BD MAX System and BD Respiratory Viral Panel-SCV2 for BD MAX System assays are designed for use with a nasopharyngeal or anterior nasal swabs collected in BD Universal Viral Transport System (UVT) or Copan Universal Transport Media System (UTM). Once collected, the UVT/UTM patient sample is vortexed and 750 ul is transferred to the BD Molecular RVP Sample Buffer Tube (SBT) provided with the BD

Respiratory Viral Panel for BD MAX System. placed in the BD MAX System. For all sample types the SBTs are vortexed and then loaded into the BD MAX system along with the Unitized Reagent Strips, Master Mix, Extraction Tubes, and PCR Cartridges. No further operator intervention is necessary.

The BD RVP Unitized Reagent Strip contains a combination of lytic and extraction reagents designed to perform cell lysis and TNA extraction. Nucleic acids released from the target organisms are captured on magnetic affinity beads. The beads, together with the bound nucleic acids, are washed and the nucleic acids are eluted by a combination of heat and pH variation. Eluted TNA is added to neutralization buffer, mixed, and transferred to BD Respiratory Viral Panel master mix for rehydration. After reconstitution, the BD MAX System dispenses a fixed volume of RT-PCR-ready solution containing extracted nucleic acids into the PCR Cartridge. Microvalves on the cartridge are sealed by the system prior to initiating PCR in order to contain the amplification mixture and thus prevent evaporation and contamination.

The amplified cDNA targets are detected using hydrolysis (TaqMan) probes, labeled at one end with a fluorescent reporter dye (fluorophore), and at the other end, with a quencher moiety. Probes labeled with different fluorophores are used to detect the target analytes in different optical channels of the BD MAX System. When the probes are in their native state, the fluorescence of the fluorophore is quenched due to its proximity to the quencher. However, in the presence of target cDNA, the probes hybridize to their complementary sequences and are hydrolyzed by the 5'–3' exonuclease activity of the DNA polymerase as it synthesizes the nascent strand along the cDNA template. As a result, the fluorophores are separated from the quencher molecules and fluorescence is emitted. The amount of fluorescence detected in the optical channels is directly proportional to the quantity of the corresponding probe that is hydrolyzed. The BD MAX System monitors these signals at each cycle of the PCR and interprets the data at the end of the reaction to provide qualitative test results for each analyte.

Identification of SARS-CoV-2, influenza A, influenza B, RSV, and RNase P occurs by the use of target-specific primers and fluorescent-labeled probes that hybridize to conserved regions in the viral genomes (**Table 1**).

Table 1. Assay Primer and Probe Targets

Analyte	Gene Targeted	Instrument Channel
SARS-CoV-2	Nucleocapsid gene (N1 and N2 regions)	FAM
Influenza A Virus	Matrix (M1) Gene	Cy5
Influenza B Virus	Matrix (M1) & Hemagglutinin (HA) Genes	Cy5.5
Respiratory Syncytial Virus A/B	N & M Genes	VIC
Extraction Control	Human RNase P gene	ROX

C Instrument Description Information:

1. Instrument Name:

BD MAX System, software version 5.14A

2. Specimen Identification:

Specimen identification can be entered either *via* barcode scanning or by manual entry.

3. Specimen Sampling and Handling:

Use of the BD Respiratory Viral Panel or the BD Respiratory Viral Panel-SCV2 requires either a nasopharyngeal swab (NPS) or anterior nasal swab (ANS) to be collected according to standard procedures and expressed into BD Universal Transport System (UVT) or Copan Universal Transport Media System (UTM). **Note:** the BD UVT and the Copan UTM are a similar device with identical chemical formulations. See Section B. Comparison Studies/Matrix Comparison.

750 uL of patient sample is transferred to the sample buffer tube (SBT) provided with the BD Respiratory Viral Panel for BD MAX System kit, which contains 4.5% Triton X-100 Reduced, and placed in the BD MAX sample rack. Unitized Reagent Strips are placed in the sample rack and securely seated. Foil-sealed dried Extraction Tubes and Master Mix Tubes are snapped into the appropriate positions on each Unitized Reagent Strip. The sample rack is placed in the BD MAX instrument along with a BD PCR cartridge. The BD MAX System automates running of the assay and result reporting.

The BD Respiratory Viral Panel Assay Extraction Tube and Unitized Reagent Strip (URS) contain a combination of lytic and extraction reagents designed to perform cell lysis and Total Nucleic Acid (TNA) extraction. 600 µL is transferred from the sample buffer tube into the URS strip for cell lysis. Following cell lysis, the released nucleic acids are captured by magnetic affinity beads. The beads with the bound nucleic acids are washed and the TNA is eluted, and neutralization buffer added. 12.5 µL of eluted TNA is then transferred to BD Respiratory Viral Panel Master Mix containing RT-polymerase, dNTPs, primers and probes. This final rehydrated Master Mix is transferred to a BD PCR cartridge for the initiation of reverse transcriptase PCR mediated conversion of RNA to cDNA and subsequent real-time PCR.

4. Calibration:

The BD MAX system preventative maintenance is performed twice per year by a BD Field Service Engineer.

5. Quality Control:

Control Materials

a) External Controls

External Control materials (assay run controls) are not provided as part of the BD Respiratory Viral Panel for BD MAX System. External Positive and Negative Controls are not used by the BD MAX System software for the purpose of sample test result interpretation. External Controls are treated as if they were patient samples. BD recommends the use of Microbiologics controls (**Table 2**) which were utilized during assay validation.

Laboratories must establish the number, type and frequency of testing control materials according to guidelines or requirements of local, provincial, state, federal, and/or country

regulations or accreditation organizations in order to monitor the effectiveness of the entire analytical process.

Table 2. Recommended External Controls

External Controls	Catalog #
Microbiologics Helix Elite Synthetic Standard SARS-CoV-2 Synthetic RNA (N gene Targets)	HE0060S
Microbiologics Helix Elite Inactivated SARS-CoV-2 Whole Virus (Pellet)	HE0065N
Microbiologics Helix Elite Inactivated SARS-CoV-2 Whole Virus (Swab)	HE0066NS
Microbiologics Helix Elite Inactivated Standard Inactivated influenza A/B and Respiratory Syncytial Virus	HE0044N
Microbiologics Helix Elite Inactivated Standard Negative Cellularity Control (Pellet)	HE0058N
Microbiologics Helix Elite Flu/RSV/SARS-CoV-2 Control Panel (Inactivated Swab)	8246
Microbiologics Helix Elite Inactivated Standard Negative Cellularity Control (Swab)	HE0067NS

b) Extraction and Internal Amplification Control

The human RNase P gene is present in all appropriately collected patient samples. It is co-amplified with SARS-CoV-2, influenza A, influenza B and RSV gene targets (if present) and will serve as both an endogenous nucleic acid extraction control and internal amplification control. In the event that SARS-CoV-2, influenza A, influenza B, RSV are negative, an RNase P result must be positive for the SARS-CoV-2, influenza A, influenza B, RSV results to be valid negative results. When either SARS-CoV-2, influenza A, influenza B and/or RSV target results are positive, RNase P result is ignored. An Unresolved (UNR) result is indicative of specimen-associated inhibition or reagent failure.

V Substantial Equivalence Information:

A Predicate Device Name(s):

BioFire Respiratory Panel 2.1 (RP2.1)

B Predicate 510(k) Number(s):

DEN200031

C Comparison with Predicate(s):

Table 3. BD Respiratory Viral Panel for BD MAX System Substantial Equivalence Comparison

Device & Predicate	<u>K230956</u>	<u>DEN200031</u>
--------------------	----------------	------------------

Device(s):		
Device Trade Name	BD Respiratory Viral Panel for BD MAX System	BioFire Respiratory Panel 2.1 (RP2.1)
Regulation Number	21 CFR 866.3981	Same
Regulation Name	Multi-Target Respiratory Specimen Nucleic Acid Test Including SARS-Cov-2 And Other Microbial Agents	Same
Product Code	QOF	Same
Intended Use/Indications For Use	<p>BD Respiratory Viral Panel for BD MAX System is an automated multiplexed real-time reverse transcriptase polymerase chain reaction (RT- PCR) test intended for the simultaneous, qualitative detection and differentiation of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), influenza A, influenza B, and/or respiratory syncytial virus (RSV) nucleic acid in nasopharyngeal swab (NPS) and anterior nasal swab (ANS) specimens from individuals with signs and symptoms of respiratory tract infection. Clinical signs and symptoms of respiratory tract infection due to SARS-CoV-2, influenza, and RSV can be similar.</p> <p>BD Respiratory Viral Panel for BD MAX System is intended for use as an aid in the differential diagnosis of SARS-CoV-2, influenza A, influenza B, and/or RSV infection if used in conjunction with other clinical and epidemiological information, and laboratory findings. SARS-CoV- 2, influenza A, influenza B, and RSV viral nucleic acid are generally detectable in NPS and ANS specimens during the acute phase of infection.</p> <p>Positive results do not rule out co-infection with other organisms. The agent(s) detected by the BD Respiratory Viral Panel for BD MAX System may not be the definitive cause of disease.</p> <p>Negative results do not preclude SARS-CoV-2, influenza A, influenza B, and/or RSV infection.</p>	<p>The BIOFIRE Respiratory Panel 2.1 (RP2.1) is a PCR based multiplexed nucleic acid test intended for use with the BIOFIRE FilmArray 2.0 or BIOFIRE FilmArray Torch systems for the simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) obtained from individuals suspected of respiratory tract infections, including COVID-19.</p> <p>The following organism types and subtypes are identified using the BIOFIRE RP2.1:</p> <p>Adenovirus, Coronavirus 229E Coronavirus HKU1 Coronavirus NL63 Coronavirus OC43 Severe Acute Respiratory Syndrome Coronavirus-2 (SARS-CoV-2) Human Metapneumovirus Human Rhinovirus/Enterovirus Influenza A, including subtypes H1, H1-2009, and H3 Influenza B Parainfluenza Virus 1 Parainfluenza Virus 2 Parainfluenza Virus 3 Parainfluenza Virus 4 Respiratory Syncytial Virus Bordetella parapertussis (IS1001) Bordetella pertussis (ptxP) Chlamydia pneumoniae, and Mycoplasma pneumoniae Nucleic acids from the</p>

	<p>The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.</p>	<p>respiratory viral and bacterial organisms identified by this test are generally detectable in NPS specimens during the acute phase of infection. The detection and identification of specific viral and bacterial nucleic acids from individuals exhibiting signs and/or symptoms of respiratory infection is indicative of the presence of the identified microorganism and aids in the diagnosis of respiratory infection if used in conjunction with other clinical and epidemiological information. The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions. Negative results in the setting of a respiratory illness may be due to infection with pathogens that are not detected by this test, or lower respiratory tract infection that may not be detected by an NPS specimen. Positive results do not rule out coinfection with other organisms. The agent(s) detected by the BIOFIRE RP2.1 may not be the definite cause of disease. Additional laboratory testing (e.g., bacterial and viral culture, immune-fluorescence, and radiography) may be necessary when evaluating a patient with possible respiratory tract infection.</p>
Condition for use	<p>For prescription use For in vitro diagnostic use only.</p>	<p>Same</p>
Sample Types	<p>Nasopharyngeal swab specimen Nasal swab specimen</p>	<p>Nasopharyngeal swab specimen</p>
Patient Population	<p>Individuals suspected of respiratory viral infection, including COVID-19.</p>	<p>Individuals suspected of respiratory tract infections, including COVID-19</p>
Analyte Targets	<p>The following organism types are identified using the BD Respiratory Viral Panel for BD MAX System:</p> <p>Severe Acute Respiratory Syndrome Coronavirus (SARS-CoV-2)</p>	<p>The following organism types and subtypes are identified using the BioFire RP2.1:</p> <ul style="list-style-type: none"> • Adenovirus • Coronavirus 229E • Coronavirus HKU1, • Coronavirus NL63

	Influenza A Influenza B and Respiratory Syncytial Virus (RSV)	<ul style="list-style-type: none"> • Coronavirus OC43 • Severe Acute Respiratory Syndrome Coronavirus (SARS-CoV-2) • Human Metapneumovirus • Human Rhinovirus/Enterovirus • Influenza A, including subtypes H1, H1-2009, and H3, • Influenza B, • Parainfluenza Virus 1, • Parainfluenza Virus 2, • Parainfluenza Virus 3, • Parainfluenza Virus 4, • Respiratory Syncytial Virus, • Bordetella parapertussis (IS1001), • Bordetella pertussis (ptxP), • Chlamydia pneumoniae, and • Mycoplasma pneumoniae
Sample Preparation Procedure	Automated by BD MAX System	Automated by BioFire FilmArray 2.0 or BioFire FilmArray Torch systems
Amplification Technology	Real-Time PCR	Nested multiplex RT-PCR
Analyte	RNA	RNA/DNA
Detection Chemistry	Paired reporter and quencher fluorescence labeled probes (TaqMan Technology) using fluorescence resonance energy transfer	Two Step Nested multiplex PCR: -Reverse transcription, followed by a multiplexed first stage PCR reaction (PCR1). Multiple simultaneous second stage PCR reactions (PCR2) to amplify sequences within the PCR1 products using fluorescence double stranded binding dye. Endpoint melting curve data to detect target specific amplicons
Control used	<ol style="list-style-type: none"> 1. The RNA Internal Control (RNase P) 2. External Positive and negative controls 	Two process controls: <ol style="list-style-type: none"> 1. RNA Process Control (IC) 2. PCR2 Control (A positive result indicates that PCR2 was successful)
Result Analysis	Based on PCR cycle threshold analysis	Endpoint melting curve data to detect target-specific amplicons
Test Interpretation	Automated test interpretation and	Same

	report generation. User cannot access raw data.	
Time to Result	About 2 hours	About 45 min

Table 4. BD Respiratory Viral Panel-SCV2 for BD MAX System Substantial Equivalence Comparison

Device & Predicate Device(s):	<u>K230956</u>	<u>DEN200031</u>
Device Trade Name	BD Respiratory Viral Panel-SCV2 for BD MAX System	BioFire Respiratory Panel 2.1 (RP2.1)
Regulation Number	21 CFR 866.3981	Same
Regulation Name	Multi-Target Respiratory Specimen Nucleic Acid Test Including Sars-Cov-2 And Other Microbial Agents	Same
Product Code	QOF	QOF
Intended Use/Indications For Use	<p>BD Respiratory Viral Panel-SCV2 for BD MAX System is an automated multiplexed real-time reverse transcriptase polymerase chain reaction (RT-PCR) test intended for the simultaneous, qualitative detection of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) viral nucleic acid in nasopharyngeal swab (NPS) and anterior nasal swab (ANS) specimens from individuals with signs and symptoms of respiratory tract infection. SARS-CoV-2 viral RNA is generally detectable in NPS and ANS specimens during the acute phase of infection.</p> <p>The BD Respiratory Viral Panel-SCV2 for BD MAX System is intended for use as an aid in the diagnosis of SARS-CoV-2 infection if used in conjunction with other clinical and epidemiological information, and laboratory findings.</p> <p>Positive results do not rule out co-infection with other organisms. The agent detected by the BD Respiratory Viral Panel-SCV2 for BD MAX System may not be the definitive cause of disease.</p> <p>Negative results do not preclude SARS-CoV-2 infection.</p>	<p>The BIOFIRE Respiratory Panel 2.1 (RP2.1) is a PCR based multiplexed nucleic acid test intended for use with the BIOFIRE FilmArray 2.0 or BIOFIRE FilmArray Torch systems for the simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) obtained from individuals suspected of respiratory tract infections, including COVID-19.</p> <p>The following organism types and subtypes are identified using the BIOFIRE RP2.1:</p> <ul style="list-style-type: none"> • Adenovirus • Coronavirus 229E • Coronavirus HKU1 • Coronavirus NL63 • Coronavirus OC43 • Severe Acute Respiratory Syndrome Coronavirus (SARS-CoV-2) • Human Metapneumovirus • Human Rhinovirus/Enterovirus • Influenza A, including subtypes H1, H1-2009, and H3 • Influenza B • Parainfluenza Virus 1

	<p>The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.</p>	<ul style="list-style-type: none"> • Parainfluenza Virus 2 • Parainfluenza Virus 3 • Parainfluenza Virus 4 • Respiratory Syncytial Virus • Bordetella parapertussis (IS1001) • Bordetella pertussis (ptxP) • Chlamydia pneumoniae, and • Mycoplasma pneumoniae <p>Nucleic acids from the respiratory viral and bacterial organisms identified by this test are generally detectable in NPS specimens during the acute phase of infection. The detection and identification of specific viral and bacterial nucleic acids from individuals exhibiting signs and/or symptoms of respiratory infection is indicative of the presence of the identified microorganism and aids in the diagnosis of respiratory infection if used in conjunction with other clinical and epidemiological information. The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions.</p> <p>Negative results in the setting of a respiratory illness may be due to infection with pathogens that are not detected by this test, or lower respiratory tract infection that may not be detected by an NPS specimen. Positive results do not rule out coinfection with other organisms. The agent(s) detected by the BIOFIRE RP2.1 may not be the definite cause of disease. Additional laboratory testing (e.g., bacterial and viral culture, immunofluorescence, and radiography) may be necessary when evaluating a patient with possible respiratory tract infection.</p>
Condition for use	For prescription use	Same

	For in vitro diagnostic use only.	
Sample Types	Nasopharyngeal swab specimen Nasal swab specimen	Nasopharyngeal swab specimen
Patient Population	Individuals suspected of COVID-19 by their healthcare provider	Individuals suspected of respiratory tract infections, including COVID-19
Analyte Targets	SARS-CoV-2	<p>The following organism types and subtypes are identified using the BioFire RP2.1:</p> <ul style="list-style-type: none"> • Adenovirus, • Coronavirus 229E, • Coronavirus HKU1, • Coronavirus NL63, • Coronavirus OC43, • Severe Acute Respiratory Syndrome Coronavirus (SARS-CoV-2), • Human Metapneumovirus, • Human Rhinovirus/Enterovirus, • Influenza A, including subtypes H1, H1-2009, and H3, • Influenza B, • Parainfluenza Virus 1, • Parainfluenza Virus 2, • Parainfluenza Virus 3, • Parainfluenza Virus 4, • Respiratory Syncytial Virus, • Bordetella parapertussis (IS1001), • Bordetella pertussis (ptxP), • Chlamydia pneumoniae, and • Mycoplasma pneumoniae
Sample Preparation Procedure	Automated by BD MAX System	Automated by BioFire FilmArray 2.0 or BioFire FilmArray Torch systems
Amplification Technology	Real-Time PCR	Nested multiplex RT-PCR
Analyte	RNA	RNA/DNA
Detection Chemistry	Paired reporter and quencher fluorescence labeled probes (TaqMan Technology) using fluorescence resonance energy transfer	Two Step Nested multiplex PCR: -Reverse transcription, followed by a multiplexed first stage PCR reaction (PCR1). Multiple simultaneous second stage PCR reactions (PCR2) to amplify sequences within the PCR1 products using double stranded binding dye. Endpoint melting curve data to detect target specific amplicons.

Control used	1. The RNA Internal Control (RNase P) 2. External Positive and negative controls	Two process controls: 1. RNA Process Control (IC) 2. PCR2 Control (A positive result indicates that PCR2 was successful)
Result Analysis	Based on PCR cycle threshold analysis	Endpoint melting curve data to detect target-specific amplicons
Test Interpretation	Automated test interpretation and report generation. User cannot access raw data.	Same
Time to Result	About 2 hours	About 45 min

VI Standards/Guidance Documents Referenced:

Quality System

- ISO 13485:2016-Medical devices -- Quality management systems
- 21 CFR 820: Quality System Regulations (QSR)
- IVDR 2017/746: REGULATION (EU) 2017/746 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU

Regulatory

- The 510(k) Program: Evaluating Substantial Equivalence in Premarket Notifications [510(k)], Guidance for Industry and Food and Drug Administration Staff (July 2014)
- Electronic Submission Template for Medical Device 510(k) Submissions, Guidance for Industry and Food and Drug Administration Staff (September 2022)

Risk Management

- ISO 14971: 2019 - Medical devices — Application of risk management to medical devices

Labeling

- ISO 15223-1: 2021 (“Medical devices -Symbols to be used with medical device labels, labelling, and information to be supplied -Part 1: General requirements”)
- ISO 20417:2021 - Medical devices — Information to be supplied by the manufacturer
- 21 CFR 809 Subpart B –Labeling before the device is shipped in interstate commerce.
- Unique Device Identification: Policy Regarding Compliance Dates for Class I and Unclassified Devices, Direct Marking, and Global Unique Device Identification

Database Requirements for Certain Devices Guidance for Industry and Food and Drug Administration Staff (July 2022)

Human Factor

- IEC 62366-1: 2015 Ed. 1.1 (2020) – Medical Devices - Part 1 Application of usability engineering to medical devices - Edition 1.1
- Applying Human Factor and Usability Engineering to Medical Devices: Guidance for Industry and Food and Drug Administration Staff (February 2016)

Design / Performance

- ISO DTS 5798: Quality Practice for detection of Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) by nucleic acid amplification methods
- CLSI EP5-A3: Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition.
- CLSI EP12-A2: User Protocol for Evaluation of Qualitative test Performance: Approved Guideline-Second Edition.
- CLSI EP07-A2: Interference Testing in Clinical Chemistry: Approved Guideline-Second Edition
- CLSI EP17-A2: Protocols for Determination of Limits of Detection and Limits of Quantitation: Approved Guideline-Second Edition.
- CLSI MM13-A: Collection, Transport, Preparation, and Storage of Specimens for Molecular Methods; Approved Guideline.
- CLSI MM17-A: Verification and Validation of Multiplex Nucleic Acid Assays.
- MDCG 2021-21 August 2021: Guidance on performance evaluation of SARS-CoV-2 in vitro diagnostic medical devices
- SARS-CoV-2 Common Specifications: Annex XIII of Commission Implementing Regulation (EU) 2022/1107 of 4 July 2022 laying down common specifications for certain class D in vitro diagnostic medical devices in accordance with Regulation (EU) 2017/746 of the European Parliament and of the Council (Text with EEA relevance)
- Class II Special Controls Guidance Document: Respiratory Viral Panel Multiplex Nucleic Assay. October 9, 2009.
- Instructions and requirements for Emergency Use Listing (EUL) submission: In vitro diagnostics detecting SARS-CoV-2 nucleic acid and rapid diagnostics tests detecting SARS-CoV-2 antigens (WHO PQDx_347 v4)
- Policy for Coronavirus Disease-2019 Tests During the Public Health Emergency (Revised) - Immediately in Effect Guidance for Clinical Laboratories, Commercial Manufacturers, and Food and Drug Administration Staff

Clinical

- ISO 14155:2020 - Clinical investigation of medical devices for human subjects — Good clinical practice
- ISO 20916:2019 - In vitro diagnostic medical devices — Clinical performance studies using specimens from human subjects — Good study practice
- 21 CFR Part 50 - Protection of Human Subjects – eCFR
- 21 CFR Part 54 – Financial Disclosure by Clinical Investigators

- 21 CFR Part 56 - Institutional Review Boards - eCFR
- 21 CFR Parts 812 - Investigational Device Exemptions – eCFR
- Declaration of Helsinki,
- Good Clinical Practice (ICH E6)
- Good Clinical Practice (GCP)

VII Performance Characteristics

A Analytical Performance:

Note: The BD Respiratory Viral Panel – SCV2 for BD MAX System (Panel detects SARS-CoV-2, Influenza A, Influenza B, RSV) is comprised of the same formulation, composition, and principle of operation as the BD Respiratory Viral Panel for BD MAX System; however, the BD Respiratory Viral Panel-SCV2 for BD MAX System utilizes a different Assay Definition File (ADF), which masks the results for influenza A/B and RSV and only reports the detection of SARS-CoV-2. As such, the supporting analytical and clinical validation data generated for the BD Respiratory Viral Panel for BD MAX System also supports validation of the BD Respiratory Viral Panel – SCV2 for BD MAX System.

1. Precision/Reproducibility:

a. Within-Laboratory Precision

Within-laboratory precision was evaluated for the BD Respiratory Viral Panel at one (1) site with one (1) reagent lot. Testing was performed over twelve (12) days, with two (2) operators performing two (2) runs per day for a total of forty-eight (48) runs. The viral materials used to generate the positive panel members were contrived in simulated nasopharyngeal matrix (See: B. Comparison Studies/Matrix Comparison, below) and included SARS-CoV-2, Influenza A, Influenza B, and RSV. Each panel member was tested in three (3) replicates. The following target concentrations were used for each target organism contained in each panel member:

- Moderate positive (MP): 3x LoD
- Low Positive (LP): 2x LoD
- True Negative (TN): No target

The results are shown in **Table 5**.

Table 5. Overall Precision Study Results Using One Lot of the BD Respiratory Viral Panel (Percent Agreement with Expected Results)

Sample Concentration	SARS-CoV-2 Percent Agreement (N), 95% CI	Flu A Percent Agreement (N), 95% CI	Flu B Percent Agreement (N), 95% CI	RSV Percent Agreement (N), 95% CI
	100%	100%	100%	100%

Moderate Positive (3x LoD)	(144/144) 97.4-100	(144/144) 97.4-100	(144/144) 97.4-100	(144/144) 97.4-100
Low Positive (2x LoD)	100% (144/144) 97.4-100	97.2% (140/144) 93.1-98.9	99.3% (143/144) 96.2-99.9	99.3% (143/144) 96.2-99.9
True Negative^a	100% (288/288) 98.7-100	100% (288/288) 98.7-100	100% (288/288) 98.7-100	100% (288/288) 98.7-100

^a For the True Negative category, the reported agreement indicates percent of negative results.

b. Reproducibility

For the Site-to-Site reproducibility study, three (3) sites (two external and one internal) were provided the same panels as described above for the Precision study. Each site performed testing on five (5) distinct days (consecutive or not), wherein each day, one (1) panel was tested by two (2) technologists. Each panel member was tested in three (3) replicates.

The qualitative reproducibility is presented below in **Table 6** by target analyte. Ct, internal criterion used to determine a final assay result, was selected as an additional means of assessing assay reproducibility. Overall mean Ct values with variance components (SD and %CV) are shown in **Table 7**.

Table 6. Site-to Site Reproducibility Study using (1) Lot of the BD Respiratory Viral Panel (Percent Agreement with Expected Results)

Sample Concentration	SARS-CoV-2 Percent Agreement (N), 95% CI	Flu A Percent Agreement (N), 95% CI	Flu B Percent Agreement (N), 95% CI	RSV Percent Agreement (N), 95% CI
Moderate Positive (3x LoD)	100% (90/90) 95.9-100	97.8% (88/90) 92.3-99.4	100% (90/90) 95.9-100	100% (90/90) 95.9-100
Low Positive (2x LoD)	100% (90/90) 95.9-100	96.7% (87/90) 90.7-98.9	100% (90/90) 95.9-100	100% (90/90) 95.9-100
True Negative^a	100% (180/180) 97.9-100	100% (180/180) 97.9-100	100% (180/180) 97.9-100	100% (180/180) 97.9-100

^a For the True Negative category, the reported agreement indicates percent of negative results.

Table 7. Site-to-Site Reproducibility Across Sites, Days, Runs, and Replicates – Ct Values

Target	Level	N	Mean Ct	Within Run		Between Run		Between Day		Between Site		Total	
				SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)
CoV-2	LP	90	33.3	0.74	2.2	0.00	0.0	0.00	0.0	0.54	1.6	0.92	2.8
CoV-2	MP	90	33.0	0.45	1.4	0.08	0.2	0.00	0.0	0.74	2.2	0.87	2.6

Target	Level	N	Mean Ct	Within Run		Between Run		Between Day		Between Site		Total	
				SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)
Flu A	LP	87	34.9	1.31	3.8	0.36	1.0	0.00	0.0	0.51	1.5	1.45	4.2
Flu A	MP	88	33.5	1.03	3.1	0.37	1.1	0.00	0.0	0.20	0.6	1.11	3.3
Flu B	LP	90	33.6	1.25	3.7	0.00	0.0	0.29	0.9	0.20	0.6	1.29	3.9
Flu B	MP	90	33.0	0.67	2.0	0.00	0.0	0.20	0.6	0.15	0.5	0.72	2.2
RSV	LP	90	32.4	1.32	4.1	0.00	0.0	0.00	0.0	0.00	0.0	1.32	4.1
RSV	MP	90	31.9	0.92	2.9	0.00	0.0	0.00	0.0	0.13	0.4	0.92	2.9

For the Lot-to-Lot reproducibility study, one (1) internal site was provided the same panels as described for the Precision study above. Three (3) reagent lots were tested across five (5) distinct days (consecutive or not) using one (1) BD MAX, wherein each day, two (2) panels were tested by two (2) technologists. Each panel member was tested in three (3) replicates.

The qualitative reproducibility is presented below in **Table 8** by target analyte. Ct, internal criterion used to determine a final assay result, was selected as an additional means of assessing assay reproducibility. Overall mean Ct values with variance components (SD and %CV) are shown in **Table 9**.

Table 8. Lot-to-Lot Reproducibility Study Results using Three (3) Lots of the BD Respiratory Viral Panel (Percent Agreement with Expected Results)

Sample Concentration	SARS-CoV-2 (N), 95% CI	Flu A (N), 95% CI	Flu B (N), 95% CI	RSV (N), 95% CI
Moderate Positive (3x LoD)	99.4% (179/180) 96.9-99.9	100% (180/180) 97.9-100	98.9% (178/180) 96.0-99.7	100% (180/180) 97.9-100
Low Positive (2x LoD)	100% (180/180) 97.9-100	97.8% (176/180) 94.4-99.1	100% (180/180) 97.9-100	100% (180/180) 97.9-100
True Negative ^a	100% (360/360) 98.9-100	100% (360/360) 98.9-100	100% (360/360) 98.9-100	100% (360/360) 98.9-100

^a – For the True Negative category, the reported agreement indicates percent of negative results.

Table 9. Lot-to-Lot Reproducibility across Operators, Days, Runs, and Replicates – Ct Values

Target	Level	N	Mean Ct	Lot		Day		Operator		Run		Within Run (Repeatability)		Total	
				SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
CoV-2	MP	179	33.7	0.26	0.8	0.06	0.2	0.10	0.3	0.00	0.0	0.61	1.8	0.67	2.0
CoV-2	LP	180	33.9	0.24	0.7	0.18	0.5	0.12	0.3	0.00	0.0	0.64	1.9	0.72	2.1
Flu A	MP	180	33.2	0.26	0.8	0.00	0.0	0.00	0.0	0.00	0.0	1.05	3.2	1.08	3.3
Flu A	LP	176	34.3	0.44	1.3	0.44	1.3	0.00	0.0	0.23	0.7	1.40	4.1	1.55	4.5
Flu B	MP	178	33.3	0.30	0.9	0.36	1.1	0.00	0.0	0.16	0.5	1.30	3.9	1.39	4.2
Flu B	LP	180	34.1	0.00	0.0	0.17	0.5	0.00	0.0	0.22	0.6	1.20	3.5	1.23	3.6

Target	Level	N	Mean Ct	Lot		Day		Operator		Run		Within Run (Repeatability)		Total	
				SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
RSV	MP	180	31.9	0.58	1.8	0.00	0.0	0.00	0.0	0.00	0.0	1.12	3.5	1.26	4.0
RSV	LP	180	32.7	0.60	1.8	0.00	0.0	0.00	0.0	0.05	0.1	0.87	2.7	1.06	3.2

2. Linearity:

Not applicable. This is a qualitative assay.

3. Analytical Specificity/Interference:

Analytical Reactivity (Inclusivity)

a. Wet-Testing

SARS-CoV-2, FluA, FluB, and RSV

This study was performed to determine the analytical reactivity of the BD Respiratory Viral Panel BD MAX System assay to detect clinically relevant strains, serotypes, or subtypes of the target species (*i.e.*, SARS-CoV-2, Flu A, Flu B, and RSV). The inclusivity panel was prepared by spiking various heat-inactivated SARS-CoV-2 strains and purified cultures of Flu A/Flu B/RSV virus serotypes/subtypes/strains encompassing temporal and geographical diversity into negative clinical NP swab VTM/UTM matrix at a concentration of ~3x LoD and testing in triplicate.

Strains that did not yield 100% reactivity at 3x LoD were prepared at higher concentrations and tested until the minimum concentration that achieved 100% reactivity was reached. The strains evaluated and the lowest concentration that achieved 100% reactivity are shown in **Table 10**.

Table 10. Inclusivity Specifications and Results

Subtype	Type	X LoD	Result				
			SARS-CoV-2	Flu A	FluB	RSV	
SARS-CoV-2	SARS-CoV-2						
	Hong Kong/VM200001061/2020	3	3/3	0/3	0/3	0/3	
	Italy-INMI1	3	3/3	0/3	0/3	0/3	
	Alpha (B.1.1.7) USA/CA CDC 5574/2020	3	3/3	0/3	0/3	0/3	
	Alpha (B.1.1.7) England/204820464/2020	3	3/3	0/3	0/3	0/3	
	Beta (B.1.351) South Africa/KRISP-K005325/2020	3	3/3	0/3	0/3	0/3	
	Kappa (B.1.617.1) USA/CA-Stanford-15 S02/2021	3	3/3	0/3	0/3	0/3	
	Gamma (P1) Japan/TY7-503/2021	3	3/3	0/3	0/3	0/3	
	Delta (B.1.617.2) USA/PHC658/2021	3	3/3	0/3	0/3	0/3	
Iota (B.1.526_2021) NY-Wadsworth-21025952-01/2021	3	3/3	0/3	0/3	0/3		

	Zeta (P2_2021) NY-Wadsworth-21006055-01/2021	3	3/3	0/3	0/3	0/3
	Omicron (USA/GA-EHC-2811C/2021)	3	3/3	0/3	0/3	0/3
H1N1	Influenza A					
	A/H1N1(pdm09)/Bangladesh/3002/2015	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)/Iowa/53/2015	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)/Michigan/272/2017	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)/Michigan/45/2015	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)/St. Petersburg/61/2015	3	1/3 ^a	3/3	0/3	0/3
	A/H1N1(pdm09)/Wisconsin/505/2018	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)/Wisconsin/588/2019	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)AVR/Louisiana/08/2013	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)AVR/Maryland/08/2013	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)AVR/New York/18/2009	3	0/3	3/3	0/3	0/3
	A/H1N1(pdm09)AVR/North Carolina/4/2014	3	0/3	3/3	0/3	0/3
	A/H1N1/Guangdong-Maonan/SWL 1536/19	3	0/3	3/3	0/3	0/3
A/H1N1(pdm09)/Idaho/07/2018	3	0/3	3/3	0/3	0/3	
H3N2	A/H3N2/Alaska/232/2015	3	0/3	3/3	0/3	0/3
	A/H3N2/Arizona/45/2018	3	0/3	3/3	0/3	0/3
	A/H3N2/California/02/2014	3	0/3	3/3	0/3	0/3
	A/H3N2/Hong Kong/2671/19	3	0/3	3/3	0/3	0/3
	A/H3N2/Norway/466/14	3	0/3	3/3	0/3	0/3
	A/H3N2/Perth/16/09	3	0/3	3/3	0/3	0/3
	A/H3N2/Singapore/INFIMH-16-0019/2016	3	0/3	3/3	0/3	0/3
	A/H3N2/South Australia/55/14	3	0/3	3/3	0/3	0/3
	A/H3N2/Stockholm/6/14	3	0/3	3/3	0/3	0/3
	A/H3N2/Texas/71/2017	3	0/3	3/3	0/3	0/3
	A/H3N2/Victoria/361/11	3	0/3	3/3	0/3	0/3
	A/H3N2/Wisconsin/04/2018	3	0/3	3/3	0/3	0/3
	H5N1	A/H5N1/common magpie/Hong Kong/645/2006	3	0/3	3/3	0/3
H5N2	A/H5N2/pheasant/New Jersey/1355/1998	3	0/3	3/3	0/3	0/3
H7N2	A/H7N2/turkey/Virginia/4529/2002	3	0/3	3/3	0/3	0/3
H7N7	A/H7N7/mallard/Netherlands/12/2000	3	0/3	1/3	0/3	0/3
		6	0/3	3/3	0/3	0/3
H7N9	A/H7N9/Anhui/1/2013	3	0/3	1/3	0/3	0/3
		6	0/3	3/3	0/3	0/3
H9N2	A/H9N2/chicken/Hong Kong/G9/1997	3	0/3	2/3	0/3	0/3
		6	0/3	3/3	0/3	0/3
Victoria	Influenza B					
	B/Colorado/6/2017	3	0/3	0/3	3/3	0/3
	B/Hawaii/01/2018	3	0/3	0/3	3/3	0/3
	B/Hong Kong/286/2017	3	0/3	0/3	3/3	0/3
	B/Missouri/12/2018	3	0/3	0/3	3/3	0/3
B/Nevada/3/2011	3	0/3	0/3	3/3	0/3	
Yamagata	B/Guangdong-Liwan/1133/2014	3	0/3	0/3	3/3	0/3
	B/Indiana/17/2017	3	0/3	0/3	3/3	0/3
	B/Oklahoma/10/2018	3	0/3	0/3	3/3	0/3
	B/Utah/9/2014	3	0/3	0/3	3/3	0/3
	B/Wisconsin/10/2016	3	0/3	0/3	3/3	0/3
RSV A	Respiratory Syncytial Virus (RSV)					
	(RSV-A) 12/2014	3	0/3	0/3	0/3	3/3
	(RSV-A) 2/2015	3	0/3	0/3	0/3	3/3
	(RSV-A) 4/2015	3	0/3	1/3 ^a	0/3	3/3
RSV B	(RSV-B) 12/2014	3	1/3 ^b	0/3	0/3	3/3
	(RSV-B) 3/2015	3	0/3	0/3	0/3	3/3

^a-Positive result for SARS-CoV-2 amplification during the testing of strain A/H1N1(pdm09)/St. Petersburg/61/2015. Curve analysis shows clear amplification of the influenza A target and RNaseP endogenous control, no amplification of the influenza B and RSV targets, and weak/delayed amplification of the SARS-CoV-2 target. The curve is suggestive of a viral titer below the limit of detection for SARS-CoV-2 target.

^b-Amplification of targets other than RSV. Curve analysis shows clear amplification of the RNaseP endogenous control and RSV targets, no amplification of the influenza A or B targets, and weak/delayed amplification of the SARS-CoV-2 target. No root cause could be identified.

The results from this study demonstrate that the BD Respiratory Viral Panel BD MAX System is capable of detecting multiple clinically relevant strains of SARS-CoV-2, Flu A, Flu B, and RSV.

b. *In-silico*

The inclusivity of the BD MAX System BD Respiratory Viral Panel SARS-CoV-2/Flu A/B/RSV assay was evaluated using *in silico* analysis of the forward primers, reverse primers, and probes for the SARS-CoV-2, Flu A, Flu B and RSV target systems in relation to sequences available in the NCBI GenBank and GISAID gene databases.

Database Retrieval and Alignments of RSV, FluA and FluB Primer/Probe Sets

Sequences were retrieved from the NCBI GenBank nucleotide database for each species. The sequence coordinates of the RNA design regions were identified using blastn alignments between a reference sequence and each species dataset. These regions were extracted from all sequences and subsequently pairwise aligned using R Bioconductor, requiring global alignments of each oligo against local alignments to the subject. Ambiguous bases were called as a match when the paired nucleotide matched one of the corresponding bases, outlined in standard IUPAC code conventions. (Dixon, Bielka, & Cantor, 1986) Primers were grouped by sets and in the case of RSV, where there are two redundant detection systems; results were compared across primer-sets to determine overall system capability to detect each individual isolate.

NCBI accession numbers were used to retrieve associated meta-data from the GenBank files for each, and the results were used to identify the lineage, type and/or host of each sequenced isolate when available

Based on *in silico* analysis of all sequences available as of February 16, 2023 in GISAID and NCBI databases, the BD MAX System BD Respiratory Viral Panel SARS-CoV-2/Flu A/B/RSV assay is predicted to detect $\geq 99.8\%$ Flu A, $\sim 100\%$ for Flu B, and $\geq 99.4\%$ for RSV A/B.

Database Retrieval and Alignments of SARS-CoV-2 Primer/Probe Sets

All available full-length SARS-CoV-2 genome sequences were retrieved from the EpiCoV database at <https://www.gisaid.org> as of February 16, 2023. Sequences meeting internal quality control criteria (n = 9,330,351), and the nucleocapsid N gene sequence was extracted by comparison to reference sequence “SARS-CoV- 2/human/USA/WA-CDC-WA1/2020” (NCBI Accession MN985325). To evaluate the potential impact of sequence variation on the

BD MAX SARS-CoV2 N1 and N2 PCR primer sets, all included sequences were compared through alignment with the BD MAX SARS-CoV-2 primers and probes.

Based on the *in-silico* analysis of GISAID and NCBI sequences available up to February 16, 2023 for SARS-CoV-2, the BD MAX System BD Respiratory Viral Panel SARS-CoV-2/Flu A/B/RSV assay, is predicted to detect $\geq 99.4\%$ of the sequenced isolates (n = 9,330,351).

Cross-Reactivity/Microbial Interference

a. Exclusivity In-silico Assessment

The aim of this study is to perform an *in-silico* search for nucleic acid sequences from species other than the intended targets that match the assay primers and probes. This analysis serves to identify potential cross reactivity concerns due to primer and probe interactions with genetic material from species other than the intended targets. The analyses for primer/probe interactions for SARS-CoV-2, Flu A/B, and RSV were performed separately.

No clinically relevant cross-reactivity was identified for SARS-CoV-2, Influenza A, Influenza B, or RSV.

b. Wet-Testing

The purpose of this study was to evaluate the analytical specificity (cross-reactivity) and microbial interference of the BD Respiratory Viral Panel for BD MAX System and BD Respiratory Viral Panel – SCV2 for BD MAX System. Analytical specificity or cross-reactivity is defined as the ability to generate negative results for the analytes (SARS-CoV-2, Flu A/B, RSV, and RNase P) in the presence of viruses, yeasts, fungi, and bacteria that are phylogenetically-related or likely to be found in respiratory tract clinical specimens, and potential microbial interference is assessing whether the candidate device can generate positive results for the analytes in the presence of viruses, yeasts, fungi, and bacteria that are phylogenetically-related or likely to be found in respiratory tract clinical specimens. Three (3) replicates in simulated nasopharyngeal matrix (See: B. Comparison Studies/Matrix Comparison, below) with and without the target analytes were tested for each potential cross-reactant/ microbial interferent listed in **Table 11** below and must give a negative SARS-Cov-2, Flu A, Flu B, and RSV result in order to be deemed non-cross-reactive, and a positive SARS-CoV-2, Flu A, Flu B, and RSV result in order to be deemed non-microbial interferent. No cross-reactivity or microbial interference was observed at the concentrations tested.

Table 11. Cross-reactivity and Microbial Interference Study Microorganisms Evaluated by Wet-Testing

Target	Concentration of Target Tested in SBT	Positive Testing (assessing potential microbial interference)	Negative Testing (assessing potential cross-
--------	---------------------------------------	---	--

						reactivity)
		SARS-CoV-2	Flu A	Flu B	RSV	Negative
Adenovirus – Type 1	1.00E+05 ¹ TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Adenovirus – Type 4	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Adenovirus – Type 7	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
<i>Aspergillus flavus</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Aspergillus fumigatus</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Aspergillus niger</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Aspergillus terreus</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Bordetella parapertussis	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Bordetella pertussis	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Candida albicans</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Chlamydia pneumonia</i>	1.00E+06 IFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Corynebacterium diphtheriae</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Cytomegalovirus	1.00E+05 copies/mL	3/3	3/3	3/3	3/3	3/3
Enterovirus B (Echovirus 6)	1.00E+05 units/mL	3/3	3/3	3/3	3/3	3/3
Enterovirus C (Coxsackievirus A16)	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Enterovirus D	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Epstein Barr virus	1.00E+05 copies/mL	3/3	3/3	3/3	3/3	3/3
<i>Escherichia coli</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Fusobacterium necrophorum</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Haemophilus influenzae</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Herpes simplex virus Type 1	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Herpes simplex virus Type 2	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Human coronavirus 229E	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Human coronavirus HKU1	1.00E+05 GC/mL	3/3	3/3	3/3	3/3	3/3
Human coronavirus NL63	1.00E+07 copies/mL	3/3	3/3	3/3	3/3	3/3
Human coronavirus OC43	1.00E+05 GE/mL	3/3	3/3	3/3	3/3	3/3

Human Metapneumovirus (hMPV)	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
<i>Lactobacillus acidophilus</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Legionella pneumophila</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Measles	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
MERS-coronavirus	1.00E+07 copies/mL	3/3	3/3	3/3	3/3	3/3
<i>Moraxella catarrhalis</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Mumps	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
<i>Mycoplasma genitalium</i>	1.82E+05 cells/mL	3/3	3/3	3/3	3/3	3/3
<i>Mycobacterium tuberculosis</i>	1.00E+06 copies/mL	3/3	3/3	3/3	3/3	3/3
<i>Mycoplasma pneumoniae</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Neisseria meningitidis</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Neisseria gonorrhoeae</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Parainfluenza virus 1	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Parainfluenza virus 2	2.12E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Parainfluenza virus 3	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
Parainfluenza virus 4	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
<i>Pneumocystis jirovecii</i> (PJP)	1.00E+06 cells/mL	3/3	3/3	3/3	3/3	3/3
Pooled human expressed nasopharyngeal swab matrix	n/a	3/3	3/3	3/3	3/3	3/3
<i>Pseudomonas aeruginosa</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Rhinovirus	1.00E+05 TCID ₅₀ /mL	3/3	3/3	3/3	3/3	3/3
SARS-coronavirus	1.00E+05 GE /mL	3/3	3/3	3/3	3/3	3/3
<i>Staphylococcus aureus</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Staphylococcus epidermis</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Streptococcus pneumoniae</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Streptococcus pyogenes</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
<i>Streptococcus salivarius</i>	1.00E+06 CFU/mL	3/3	3/3	3/3	3/3	3/3
Varicella-zoster virus	1.00E+07 copies/mL	3/3	3/3	3/3	3/3	3/3

¹CFU = Colony Forming Units; GE = Genome Equivalents; IFU = Inclusion Forming Units; TCID₅₀ = Median Tissue Culture Infectious Dose

Interfering Substances

This study evaluated the performance of the BD Respiratory Viral Panel for BD MAX System and BD Respiratory Viral Panel – SCV2 for BD MAX System in the presence of medications, over the counter products, and other potentially interfering substances found in a clinical respiratory specimen. Assay results were evaluated to determine if the presence of potentially interfering substances in target analyte negative or target analyte positive samples had an effect on assay performance. The study utilized simulated nasopharyngeal matrix containing mucin and human DNA to establish the interference claims (See: B. Comparison Studies/Matrix Comparison, for validation of equivalency between matrices).

Each interferent was tested with three (3) positive replicates containing SARS-CoV-2, Influenza A, Influenza B, and RSV and three negative replicates containing only the interferent with the simulated nasopharyngeal matrix for a total of six (6) samples per interferent level. Three (3) positive replicates containing SARS-CoV-2, Influenza A, Influenza B, and RSV were run without interferent to demonstrate passing acceptance criteria for positive samples. Whole blood (Human) at 2% v/v was found to interfere with the assay of both SARS-CoV-2 and Influenza B. When the amount of blood was titrated downward to 0.2% v/v, it no longer interfered with either assay. **Table 12.**

Table 12. Interfering Substance Results

Substance	Active Ingredient	Concentration Tested	Positive Testing (Positive/Total)				Negative Testing (Negative / Total)	Result
			SARS-CoV-2	Influenza A	Influenza B	RSV		
None	N/A	N/A	3/3	3/3	3/3	3/3	N/A	N/A
Oral anesthetic and analgesic	Benzocaine	0.8 mg/mL	3/3	3/3	3/3	3/3	3/3	NI
	Menthol		3/3	3/3	3/3	3/3	3/3	NI
Biologicals	Purified Mucin	60 µg/mL	3/3	3/3	3/3	3/3	3/3	NI
	Whole Blood (human)	2% v/v	1/3	3/3	1/3	3/3	3/3	I
		0.2% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Leukocytes	2% v/v	3/3	3/3	3/3	3/3	3/3	NI
Nasal Sprays / Drops	Zinc	1 mg/mL	3/3	3/3	3/3	3/3	3/3	NI
	Phenylephrine	5% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Oxymetazoline	5% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Sodium Chloride with preservatives	5% v/v	3/3	3/3	3/3	3/3	3/3	NI
Corticosteroids	Beclomethasone	17% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Dexamethasone	17% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Flunisolide	17% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Triamcinolone	17% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Budesonide	17% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Mometasone	17% v/v	3/3	3/3	3/3	3/3	3/3	NI
Nasal Gel - Homeopathic Allergy Relief	Luffa operculata							
	Sulfur							

	Galphimia glauca	5% v/v	3/3	3/3	3/3	3/3	3/3	NI
	Histaminum hydrochloricum							
Antiviral Drug	Zanamivir	3.3 mg/mL	3/3	3/3	3/3	3/3	3/3	NI
Antibiotic	Mupirocin	10 mg/mL	3/3	3/3	3/3	3/3	3/3	NI
Antibacterial	Tobramycin	4 µg/mL	3/3	3/3	3/3	3/3	3/3	NI
NI = Non-Interferent, I = Interferent								

FluMist Quadrivalent

In a supplementary interfering substances study, the FluMist Quadrivalent vaccine (MedImmune, LLC) was assayed for its impact on the performance of the BD Respiratory Viral Panel for BD MAX System and the BD Respiratory Viral Panel – SCV2 for BD MAX System. FluMist vaccine contains live Influenza A and B strains and was found to interfere with the BD Respiratory Viral Panel at concentrations greater than 6.67E-12% in nasopharyngeal specimens. **Tables 13 & 14.**

Table 13. Interfering Substance Positive Testing Results

Substance	Active Ingredient	Concentration Tested	Positive Testing (Positive/Total)				Result
			SARS- CoV-2	Influenza A	Influenza B	RSV	
None	N/A	N/A	3/3	3/3	3/3	3/3	N/A
FluMist Quadrivalent	Live, attenuated Flu A and Flu B strains	6.67%	3/3	3/3	3/3	3/3	NI
		6.67E-04%	3/3	3/3	3/3	3/3	NI
		6.67E-08%	3/3	3/3	3/3	3/3	NI
		6.67E-12%	3/3	3/3	3/3	3/3	NI
NI = Non-Interferent, I = Interferent							

Table 14. Interfering Substance Negative Testing Results

Substance	Active Ingredient	Concentration Tested	Negative Testing (Negative/Total)				Result
			SARS- CoV-2	Influenza A	Influenza B	RSV	
FluMist Quadrivalent	Live, attenuated Flu A and Flu B strains	6.67%	3/3	0/3	0/3	3/3	I
		6.67E-04%	3/3	0/3	0/3	3/3	I
		6.67E-08%	3/3	2/3	2/3	3/3	I
		6.67E-12%	3/3	3/3	3/3	3/3	NI
NI = Non-Interferent, I = Interferent							

Competitive Interference

The purpose of this study was to demonstrate the ability of the BD Respiratory Viral Panel and BD Respiratory Viral Panel – SCV2 for the BD Max system assay to detect the targeted analytes in cases of simulated mixed infection, when three (3) viral targets are low and one (1) target is high (competitive interference).

One (1) representative strain of each targeted organism (SARS-CoV-2, Flu A, Flu B and RSV) was tested. Co-infection panels were made by spiking one (1) target organism at high

concentration (1.0E+06 TCID₅₀/mL or cp/mL) and three (3) other target organisms at low concentration (2x LoD) in simulated nasopharyngeal matrix with heat-inactivated virus (See: B. Comparison Studies/Matrix Comparison, for validation of equivalency between matrices). The study was performed with 20 replicates per condition. Acceptance criteria ≥95%. None of the analytes present at a very high concentration interfered with the detection of low levels of the other three analytes. See **Table 15**.

Table 15. Mixed Infection Results for the BD Respiratory Viral Panel for the BD MAX System

Low (2X LoD)/ High (1E+6 cp/mL)	Positivity			
	SARS-CoV-2	Flu A	Flu B	RSV
High (SCV2)	100% (20/20)	100% (20/20)	100% (20/20)	100% (20/20)
High (Flu A)	100% (20/20)	100% (20/20)	100% (20/20)	100% (20/20)
High (Flu B)	100% (20/20)	95% (19/20)	100% (20/20)	100% (20/20)
High (RSV)	95% (19/20)	95% (19/20)	95% (19/20)	100% (20/20)
All Low	100% (20/20)	95% (19/20)	100% (20/20)	100% (20/20)

4. Assay Reportable Range:

Not applicable; this is a qualitative assay.

5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):

Traceability

This study was conducted in order to establish the limit of detection (LoD) of the SARS-CoV-2 target using the WHO international standard [Heat inactivated, England/02/2020, NIBSC code: 20/146] with the BD Respiratory Viral Panel and BD Respiratory Viral Panel-SCV2 for BD MAX System. The data generated from this study will be used to support the conformity of the assay to common specifications in the EU for SARS-CoV-2 assays. The LoD of the WHO international standard was determined in both clinical nasal matrix and clinical nasopharyngeal matrix.

The LoD confirmation was performed with 20 replicates for each matrix type used, thus aligning with CLSI 17-A2 Evaluation of the *Detection Capability for Clinical Laboratory Measurement Procedures* and WHO PQDx_347. The LoD of SARS-CoV-2 [England/02/2020, NIBSC code: 20/146] was determined to be 726 IU/mL. Results are shown in **Table 16**.

Table 16. WHO International Standard LoD Results

Level	Concentration (IU/mL)	Positivity		SARS-CoV-2 Ct		RNaseP Ct	
		Nasopharyngeal	Nasal	Nasopharyngeal	Nasal	Nasopharyngeal	Nasal
		100%	100%				

3x	726	(20/20)	(20/20)	35.0	35.1	22.2	23.7
1x	242	50% (10/20)	65% (13/20)	37.0	37.2	22.1	23.9
0.3x	80.7	30% (6/20)	35% (7/20)	36.8	37.3	22.0	23.8

a) Specimen Stability

Specimen stability studies were performed to support the following stability claims:

Table 17. BD Respiratory Viral Panel Assay for BD MAX System Specimen Stability

Specimen Stability	Temperature	Duration
In UVT/UTM	25 ± 2 °C	12 hours
	2–8 °C	72 hours
In BD Molecular RVP Sample Buffer Tube	25 ± 2 °C	24 hours
	2–8 °C	120 hours

Freeze/thaw specimen stability studies were performed and indicate that SARS-CoV-2, Influenza A, Influenza B, and RSV are stable for up to two (2) freeze / thaw cycles in either clinical nasopharyngeal matrix or clinical nasal matrix expressed in UVT/UTM (neat) or post transfer of UVT/UTM matrix into sample buffer tubes (nested). Once thawed, samples are stable for up to 12 hours at 25 ± 2 °C. Specimens were constructed using natural clinical nasopharyngeal or anterior nasal matrix spiked at 3X LoD for each analyte.

b) Reagent Kit Stability

To evaluate reagent kit stability, the reagents were refrigerated at (5±3°C) and stored at “room temperature” (RT 25 ± 2°C) in parallel. At defined time points during the study, the Master Mix (MM) and Extraction Tube (Tx) foil pouches are being 1) unsealed, 2) opened, 3) closed and then, 4) refrigerated and stored at RT in parallel to simulate in-use conditions across the targeted temperature range of 2-25°C. Reagent tubes are stable for up to 31 days at 2-25 °C after initial opening and re-sealing of the pouch.

c) Reagent Kit Shipping Stability

The study also assessed reagent kit shipping stability. The purpose is to verify the stability of the reagents throughout their stored shelf-life at 2-25°C after being exposed to conditions which could reasonably be encountered during normal shipping conditions. Simulated shipping stability was performed by subjecting reagents to shipping conditions. This is performed in an environmental chamber that adjusts both temperature and humidity to mimic shipping conditions. After shipping exposure, reagents are tested for a post-exposure baseline and subsequently included in the real-time stability arm of the study.

6. Detection Limit:

SARS-CoV-2 LoD Determination in Clinical NP and NS Matrix in UVT/UTM

The purpose of this study was to establish the Limit of Detection (LoD) for SARS-CoV-2 in clinical nasopharyngeal (NP) and nasal (NS) matrices for the BD Respiratory Viral Panel for BD MAX System.

LoD studies determine the lowest detectable concentration of virus at which approximately 95% of all (true positive) replicates test positive. LoD was estimated using probit analysis. Confirmation of the estimated LoD with heat-inactivated SARS-CoV-2 (2019-nCoV/USA-WA1/2020) was performed with one (1) reagent lot in replicates of 20 prepared in clinical nasopharyngeal and clinical nasal matrix. Refer to **Table 18** for the BD Respiratory Viral Panel-SCV2 LoD for SARS-CoV-2 in clinical nasopharyngeal matrix in UVT/UTM. The LoD for SARS-CoV-2 in NP matrix in UVT/UTM was determined to be 700 copies/mL.

Table 18. LoD determination for SARS-CoV-2 in Clinical NP Matrix

Concentration			Positivity	Mean Ct	
IU/mL	Copies/mL	TCID ₅₀ /mL		SARS-CoV-2	RNaseP
112	323	4.96E+00	90% (18/20)	37.1	23.4
153	442	6.79E+00	90% (18/20)	35.4	24.2
194	561	8.62E+00	85% (17/20)	35.4	22.5
242	700	1.08E+01	95% (19/20)	34.2	23.3

The LoD determination for SARS-CoV-2 in clinical anterior nasal matrix in UVT/UTM was carried out using the data from the NPS study as a starting point. Results are shown in **Table 19**. Results for NP matrix, show the LoD for SARS-CoV-2 in ANS Matrix in UVT/UTM is 700 copies/mL.

Table 19. LoD determination for SARS-CoV-2 in Clinical ANS Matrix

Concentration			Positivity	Mean Ct	
IU/mL	Copies/mL	TCID ₅₀ /mL		SARS-CoV-2	RNaseP
153	442	6.79E+00	80% (16/20)	36.3	25.5
242	700	1.08E+01	100% (20/20)	35.9	25.9

FluA, FluB, and RSV LoD Determination in Clinical NPS and ANS Matrix in UVT/UTM

The LoDs of live strains of Influenza A [H1N1/Brisbane] & [H3N2/Kansas/14/17], Influenza B [Victoria/Colorado] & [Yamagata/Phuket/3073/13], RSV A/B were determined in both NP and

ANS clinical matrices in UVT/UTM. Twenty (20) replicates of each dilution were carried out. **Table 20** shows the compiled LoD determinations for Flu A, FluB, and RSV A/B in clinical NP matrix in UVT/UTM.

Clinical Nasopharyngeal Matrix

Table 20. Analytical LoDs for Influenza A, Influenza B, and RSB A/B in NP Matrix in UVT/UTM

Strain ID	Confirmed Limit of Detection	
Influenza A/H1N1/Brisbane	500 copies/mL	5.59E-03 TCID50/mL
Influenza A/H3N2/Kansas	500 copies/mL	2.75E-01 TCID50/mL
Influenza B/Colorado	33 copies/mL	6.80E-03 TCID50/mL
Influenza B/Phuket/3073/13	100 copies/mL	2.90E-02 TCID50/mL
RSV A	143 copies/mL	3.09E-02 TCID50/mL
RSV B	687 copies/mL	1.68E-02 TCID50/mL

Clinical Anterior Nasal Matrix

Twenty replicates of each dilution were carried out. **Table 21** shows the compiled LoD determinations for Flu A, FluB, and RSV A/B in clinical ANS matrix in UVT/UTM.

Table 21 Analytical LoDs for Influenza A, Influenza B, and RSB A/B in ANS Matrix in UVT/UTM

Strain ID	Confirmed Limit of Detection	
Influenza A/H1N1/Brisbane	500 copies/mL	5.59E-03 TCID50/mL
Influenza A/H3N2/Kansas	500 copies/mL	2.75E-01 TCID50/mL
Influenza B/Colorado	33 copies/mL	6.80E-03 TCID50/mL
Influenza B/Phuket/3073/13	33 copies/mL	9.56E-03 TCID50/mL
RSV A	143 copies/mL	3.09E-02 TCID50/mL
RSV B	229 copies/mL	5.61E-03 TCID50/mL

7. Assay Cut-Off

This series of experiments was conducted in order to determine assay thresholds for the SARS-CoV-2, Flu A, Flu B, and RSV target assays by utilizing both contrived positive and negative nasopharyngeal swab (NPS) and simulated nasopharyngeal (SIM) specimens in UVT/UTM. NP specimens were representative for both NPS and anterior nasal swab (ANS) specimen types as NPS specimens are considered worst case. NPS specimens provide the most amount of background material while simulated matrix represents a specimen with the least amount of background material. Specimens were tested over six (6) master mix (MM) reagent lots across fifteen (15) instruments over three (3) test days.

The assay cut-off for the BD Respiratory Viral Panel for BD MAX was established based on PCR-based metrics taken together by the BD MAX software algorithm to make the qualitative decision whether a curve is to be considered positive or negative. These metrics (e.g., Ct, RFU endpoints, signal-to-noise ratios) are initially set by default parameters defined by the instrument. As the product undergoes product development, the data are supplemented, and the algorithm is adjusted (“trained”) using viral cultures spiked into clinical background matrices at levels surrounding the limit of detection and expected clinical range.

8. Accuracy (Instrument):
Not Applicable

9. Carry-Over:

Carryover / Cross-Contamination

This study was conducted to investigate within-run and between-run carryover while processing samples with high viral load of SARS-CoV-2 in the BD Respiratory Viral Panel for BD MAX System. High positive samples contained heat inactivated SARS-CoV-2 spiked into pooled nasal swab matrix at a concentration of $\geq 1.94E+07$ copies/mL. The negative samples consisted of simulated nasopharyngeal (NP) matrix without any target analyte. Twelve (12) replicates of the high positive panel member and 12 replicates of the negative panel member were tested in nine (9) runs by alternating negative and positive samples, using three BD MAX Systems. A total of 108 positive and 108 negative samples were tested. Of the 108 negative samples tested, one (1) false positive result was obtained (0.93%, 95% CI: 0.16–5.06%).

B Comparison Studies:

1. Method Comparison with Predicate Device:

Not Applicable

2. Matrix Comparison:

Matrix Equivalency

The purpose of this study was to demonstrate the functional equivalency among clinical nasopharyngeal swab (NPS) specimens in either BD Universal Viral Transport (UVT/UTM) media, clinical nasal swab specimens expressed in either BD Universal Viral Transport (UVT/UTM) media and simulated nasopharyngeal/nasal matrix (contains mucin and human cells/DNA) when tested with the BD Respiratory Viral Panel for BD MAX System assay. The study utilized contrived low (2x LoD) and moderate (5x LoD) positive samples and negative samples. The viral samples employed were: heat-inactivated SARS-CoV-2 (USA-WA/2020); Influenza A/H1N1/Brisbane; Influenza B/Yamagata/Phuket/3073/13, and RSV/A. The results are shown in

Table 22 and support that the BD MAX System BD Respiratory Viral Panel SARS-CoV-2/Flu A/B/RSV assay can be used with the evaluated viral transport media types, and the simulated nasopharyngeal/nasal matrix (contains mucin and human cells/DNA) is similar to the natural nasopharyngeal/nasal swabs in UVT/UTM matrix with regard to the candidate assay performance.

Table 22 Matrix Equivalency Results

Conc.	Matrix	Positivity				Ct				
		SARS-CoV-2	Influenza A	Influenza B	RSV	SARS-CoV-2	Influenza A	Influenza B	RSV	RNase P
Negative	Clinical NP Matrix	0%	0%	0%	0%	n/a	n/a	n/a	n/a	22.1
		(0/15)	(0/15)	(0/15)	(0/15)					
	Clinical Nasal Matrix	0%	0%	0%	0%	n/a	n/a	n/a	n/a	24.6
		(0/15)	(0/15)	(0/15)	(0/15)					
	Simulated NP/NS Matrix	0%	0%	0%	0%	n/a	n/a	n/a	n/a	23.2
		(0/15)	(0/15)	(0/15)	(0/15)					
2X LOD	Clinical NP Matrix	100%	100%	97%	100%	33.9	33.2	33.2	32.4	22.6
		(30/30)	(30/30)	(29/30)	(30/30)					
	Clinical Nasal Matrix	100%	100%	100%	100%	34.1	33.1	33.8	32.6	24.3
		(30/30)	(30/30)	(30/30)	(30/30)					
	Simulated NP/NS Matrix	100%	97%	100%	100%	33.7	32.6	34.0	32.5	23.4
		(30/30)	(29/30) ^a	(30/30)	(30/30)					
5X LOD	Clinical NP Matrix	100%	100%	100%	100%	31.7	31.5	32.3	30.8	22.5
		(15/15)	(15/15)	(15/15)	(15/15)					
	Clinical Nasal Matrix	100%	100%	100%	100%	32.1	31.5	32.6	31.0	24.6
		(15/15)	(15/15)	(15/15)	(15/15)					
	Simulated NP/NS Matrix	100%	100%	100%	100%	31.8	31.9	32.4	31.3	24.0
		(15/15)	(15/15)	(15/15)	(15/15)					

^a There was one non-reportable result during the study

C Clinical Studies:

Prospective Clinical Study

The clinical performance of the BD Respiratory Viral Panel for BD MAX System was established in a multi-center prospective study conducted with paired nasopharyngeal swab (NPS) and anterior nasal swab (ANS) specimens in UVT/UTM, collected following receipt of informed consent from study participants. All specimens were prospectively collected (i.e., all comers between two time points that met the clinical study inclusion criteria) from patients with signs and symptoms of respiratory tract infections during the 2022 respiratory illness season. Subject matched, paired NPS and ANS specimens were collected from six (6) geographically diverse clinical sites, four (4) in the U.S. and two (2) in Europe and were

tested with the BD Respiratory Viral Panel for BD MAX System at four (4) U.S. testing sites.

The BD Respiratory Viral Panel for BD MAX System was evaluated for SARS-CoV-2 performance by comparing the candidate device testing results to a composite comparator algorithm (CCA) consisting of three (3) highly sensitive U.S. FDA EUA SARS-CoV-2 molecular tests. A final CCA result was assigned when two of the three composite comparator assays were in concordance. The comparator method utilized to establish performance for the Flu A, Flu B, and RSV targets was a U.S. FDA-cleared molecular Flu A/B/RSV assay. All comparator testing was performed in accordance with the respective package inserts.

The BD Respiratory Viral Panel for BD MAX System was evaluated with both Category I specimens (i.e., prospectively collected, tested fresh) and Category II specimens (i.e., prospectively collected, tested frozen/thawed).

A total of 2005 paired NPS and ANS specimens were enrolled for the prospective clinical study between January 2022 and August 2022. Category I specimens were collected from April 2022 to August 2022, while Category II specimens were collected from January 2022 to April 2022.

Of these 2005 subjects, 360 were excluded as non-compliant (protocol deviations, *etc.*). This left a total of 1645 Compliant Subjects. A demographic summary of the compliant subjects is shown in **Table 23**.

For the SARS-CoV-2 evaluation, an additional 100 NP specimens were excluded (protocol deviations, *etc.*) leaving 1545 NPS specimens for performance evaluation. Similarly, an additional 84 NS specimens were excluded (protocol deviations, *etc.*) leaving 1561 ANS specimens for evaluation.

Of the 1545 valid NPS specimens, 1021 (66.1%) were tested fresh (Category I specimens) and 524 (33.9%) were tested frozen/thawed (Category II specimens) with the BD Respiratory Viral Panel for BD MAX System. Of the 1561 valid ANS specimens, 1021 (65.4%) were tested fresh and 540 (34.6%) were tested frozen/thawed with the BD Respiratory Viral Panel for BD MAX System.

For the Flu A, Flu B, and RSV evaluations, an additional 83 NS specimens were excluded (protocol deviations, *etc.*) leaving a total of 1562 NPS specimens. For the NS specimens, an additional 81 specimens were excluded (protocol deviations, *etc.*) leaving 1564 ANS specimens for evaluation.

Of the 1562 valid NPS specimens, 1022 (65.4%) were tested fresh and 540 (34.6%) were tested frozen/thawed with the BD Respiratory Viral Panel for BD MAX. Of the 1564 valid NS specimens, 1025 (65.5%) were tested fresh and 539 (34.5%) were tested frozen/thawed with the BD Respiratory Viral Panel for BD MAX System.

Table 23. Prospective Specimens: Demographic Summary of Compliant Subjects

Demographics	Characteristics	Total (N=1645)
Gender	Female	61.8% (1016/1645)
	Male	38.2% (629/1645)
Age Group	0 - 5 years	1.2% (20/1645)
	6 - 21 years	10.3% (170/1645)
	22 - 59 years	58.2% (957/1645)
	> 59 years	30.3% (498/1645)
Ethnicity	Hispanic / Latino	47.1% (775/1645)
	Non-Hispanic / Latino	49.5% (814/1645)
	Not Reported	3.4% (56/1645)
Race	American Indian or Alaska Native	0.3% (5/1645)
	Asian	0.3% (5/1645)
	Black or African American	13.3% (219/1645)
	Native Hawaiian or Other Pacific Islander	0.0% (0/1645)
	White	84.8% (1395/1645)
	Mixed Race	0.9% (15/1645)
	Not Reported	0.1% (1/1645)
	Unknown	0.3% (5/1645)
Patient Population	Outpatient	95.9% (1578/1645)
	Hospitalized	3.0% (49/1645)
	Emergency	1.0% (17/1645)
	Unknown	0.1% (1/1645)
Immuno-compromised	Yes	3.2% (52/1645)
	No	95.2% (1566/1645)
	Unknown	1.6% (27/1645)

A summary of the BD Respiratory Viral Panel for BD MAX System assay prospective clinical study performance is provided in **Tables 24 & 25**. Positive Percent Agreement (PPA) was calculated as $100\% \times (TP / (TP + FN))$. True positive (TP) indicates that both the BD Respiratory Viral Panel for BD MAX System assay and the comparator method had a positive result for the specific analyte, and false negative (FN) indicates that the BD Respiratory Viral Panel for BD MAX System was negative while the comparator result was positive. Negative Percent Agreement (NPA) was calculated as $100\% \times (TN / (TN + FP))$. True negative (TN) indicates that both the BD Respiratory Viral Panel for BD MAX System assay and the comparator method had negative results, and false positive (FP) indicates that the BD Respiratory Viral Panel for BD MAX System assay was positive while the comparator result was negative.

Table 24. BD Respiratory Viral Panel for BD MAX System Performance with Prospective NPS Specimens (Fresh, Frozen, and Fresh/Frozen Combined)

Analyte		Positive Percent Agreement			Negative Percent Agreement		
		TP/ (TP+FN)	%	95% CI	TN/ (TN+FP)	%	95% CI
SARS-CoV-2 ^a	Fresh	370/372	99.5	98.1-99.9	641/649	98.8	97.6-99.4
	Frozen	147/151	97.4	93.4-99.0	358/373	96.0	93.5-97.5
	Combined	517/523	98.9	97.5-99.5	999/1022	97.7	96.6-98.5
Flu A ^b	Fresh	58/60	96.7	88.6-99.1	955/962	99.3	98.5-99.6
	Frozen	1/1	100	20.7-100	539/539	100	99.3-100

	Combined	59/61	96.7	88.8-99.1	1494/1501	99.5	99.0-99.8
Flu B ^c	Fresh	0/0	0	NC	1021/1022	99.9	99.4-100
	Frozen	0/0	0	NC	540/540	100	99.3-100
	Combined	0/0	0	NC	1561/1562	99.9	99.6-100
RSV	Fresh	11/11	100	74.1-100	1011/1011	100	99.6-100
	Frozen	1/1	100	20.7-100	539/539	100	99.3-100
	Combined	12/12	100	75.8-100	1550/1550	100	99.8-100

TP – true positive; FN – false negative; TN – true negative; FP – false positive; NC – not calculable

^a SARS-CoV-2 was detected in 3/6 FN specimens with all three composite comparator methods. SARS-CoV-2 was detected in 15/23 FP specimens with one of the three composite comparator methods.

^b Flu A was detected in both FN specimens when tested with an independent molecular method. Flu A was detected in 3/7 FP specimens when tested with an independent molecular method. Flu A was Equivocal in 1/7 FP specimens when tested with an independent method.

^c Flu B was not detected in the single FP specimen when tested with an independent molecular method.

Table 25. BD Respiratory Viral Panel for BD MAX System Performance with Prospective ANS Specimens (Fresh, Frozen, and Fresh/Frozen Combined)

Analyte		Positive Percent Agreement			Negative Percent Agreement		
		TP/ (TP+FN)	%	95% CI	TN/ (TN+FP)	%	95% CI
SARS-CoV-2 ^a	Fresh	340/344	98.8	97.0-99.5	665/677	98.2	96.9-99.0
	Frozen	138/142	97.2	93.0-98.9	385/398	96.7	94.5-98.1
	Combined	478/486	98.4	96.8-99.2	1050/1075	97.7	96.6-98.4
Flu A ^b	Fresh	61/63	96.8	89.1-99.1	958/962	99.6	98.9-99.8
	Frozen	1/1	100	20.7-100	538/538	100	99.3-100
	Combined	62/64	96.9	89.3	1496/1500	99.7	99.3-99.9
Flu B	Fresh	0/0	0	NC	1025/1025	100	99.3-100
	Frozen	0/0	0	NC	539/539	100	99.3-100
	Combined	0/0	0	NC	1564/1564	100	99.8-100
RSV ^c	Fresh	11/11	100	74.1-100	1013/1014	99.9	99.4-100
	Frozen	0/1	0	0-79.3	538/538	100	99.3-100
	Combined	11/12	91.7	64.6-98.5	1551/1552	99.9	99.6-100

TP – true positive; FN – false negative; TN – true negative; FP – false positive; NC – not calculable

^a SARS-CoV-2 was detected in 2/8 FN specimens with all three composite comparator methods. SARS-CoV-2 was detected in 17/25 FP specimens with one of the three composite comparator methods.

^b Flu A was not detected in both FN specimens when tested with an independent molecular method. Flu A was detected in 1/4 FP specimens when tested with an independent molecular method.

^c RSV was detected in the single FN specimen when tested with an independent molecular method. RSV was detected in the single FP specimen when tested with an independent molecular method.

Non-Reportable Results Rates

Table 26 and Table 27 provide the non-reportable results rates for nasopharyngeal (NP) and anterior nasal swab (ANS) specimens, respectively.

Table 26. Prospective NP Specimens: SARS-CoV-2/Flu/RSV UNR, IND, INC and Total Non-reportable Rate

Sample Type (NP)	Unresolved (UNR)		Indeterminate (IND)		Incomplete (INC)		Total Non-Reportable (UNR+IND+INC)	
	Initial (95% CI)	Valid Repeat (95% CI)	Initial (95% CI)	Valid Repeat (95% CI)	Initial (95% CI)	Valid Repeat (95% CI)	Initial (95% CI)	Valid Repeat (95% CI)

Fresh	0.0% (0/1022) (0.0%, 0.4%)	0.0% (0/1022) (0.0%, 0.4%)	0.0% (0/1022) (0.0%, 0.4%)	0.0% (0/1022) (0.0%, 0.4%)	1.3% (13/1022) (0.7%, 2.2%)	0.0% (0/1022) (0.0%, 0.4%)	1.3% (13/1022) (0.7%, 2.2%)	0.0% (0/1022) (0.0%, 0.4%)
Frozen	0.2% (1/541) (0.0%, 1.0%)	0.2% (1/541) (0.0%, 1.0%)	0.0% (0/541) (0.0%, 0.7%)	0.0% (0/541) (0.0%, 0.7%)	0.0% (0/541) (0.0%, 0.7%)	0.0% (0/541) (0.0%, 0.7%)	0.2% (1/541) (0.0%, 1.0%)	0.2% (1/541) (0.0%, 1.0%)
Total	0.1% (1/1563) (0.0%, 0.4%)	0.1% (1/1563) (0.0%, 0.4%)	0.0% (0/1563) (0.0%, 0.2%)	0.0% (0/1563) (0.0%, 0.2%)	0.8% (13/1563) (0.5%, 1.4%)	0.0% (0/1563) (0.0%, 0.2%)	0.9% (14/1563) (0.5%, 1.5%)	0.1% (1/1563) (0.0%, 0.4%)

Table 27. Prospective ANS Specimens: SARS-CoV-2/Flu/RSV UNR, IND, INC and Total Non-reportable Rate

Sample Type (ANS)	Unresolved (UNR)		Indeterminate (IND)		Incomplete (INC)		Total Non-Reportable (UNR+IND+INC)	
	Initial (95% CI)	Valid Repeat (95% CI)	Initial (95% CI)	Valid Repeat (95% CI)	Initial (95% CI)	Valid Repeat (95% CI)	Initial (95% CI)	Valid Repeat (95% CI)
Fresh	0.1% (1/1027) (0.0%, 0.5%)	0.1% (1/1027) (0.0%, 0.5%)	0.0% (0/1027) (0.0%, 0.4%)	0.0% (0/1027) (0.0%, 0.4%)	1.3% (13/1027) (0.7%, 2.2%)	0.0% (0/1027) (0.0%, 0.4%)	1.4% (14/1027) (0.8%, 2.3%)	0.1% (1/1027) (0.0%, 0.5%)
Frozen	0.2% (1/541) (0.0%, 1.0%)	0.2% (1/541) (0.0%, 1.0%)	0.4% (2/541) (0.1%, 1.3%)	0.0% (0/541) (0.0%, 0.7%)	0.0% (0/541) (0.0%, 0.7%)	0.0% (0/541) (0.0%, 0.7%)	0.6% (3/541) (0.2%, 1.6%)	0.2% (1/541) (0.0%, 1.0%)
Total	0.1% (2/1568) (0.0%, 0.5%)	0.1% (2/1568) (0.0%, 0.5%)	0.1% (2/1568) (0.0%, 0.5%)	0.0% (0/1568) (0.0%, 0.2%)	0.8% (13/1568) (0.5%, 1.4%)	0.0% (0/1568) (0.0%, 0.2%)	1.1% (17/1568) (0.7%, 1.7%)	0.1% (2/1568) (0.0%, 0.5%)

Co-infections Detected by the BD Respiratory Viral Panel for BD MAX Assay in the Prospectively Collected Specimens

The BD Respiratory Viral Panel for BD MAX System detected a total of 4 of the 5 co-infections (3 ANS, 1 NP) in which SARS-CoV-2 and Flu A were detected by the comparators. Only Flu A was detected by the candidate in the remaining sample (ANS).

The candidate assay detected 2 additional samples with SARS-CoV-2 and Flu A co-infections (1 ANS, 1 NP) and 1 additional sample with SARS-CoV-2 and Flu B co-infection (NP). However, the comparators detected only SARS-CoV-2 in all three samples.

Retrospective Clinical Study

Flu B and RSV were of low prevalence and were not encountered in sufficiently large numbers during the prospective clinical study to adequately demonstrate assay performance. To supplement the results of the prospective clinical study, an evaluation of 20 RSV standard of care (SoC) positive and 14 Flu B SoC positive Category III ANS specimens was undertaken. These specimens, in addition to negative specimens collected under a BD collection protocol, in UVT/UTM, were collected between February 2021 and February 2023

and were randomized and de-identified. Five (5) of the 20 RSV SoC positive samples provided by the source laboratory, that were not confirmed as positive by the comparator, are included in the total N=187. They were counted as part of the Negative Percent Agreement (NPA). One (1) specimen generated an Unresolved (UNR) assay result for both Flu B and RSV and was excluded from the performance analysis. A summary of the demographic information of the retrospective ANS samples is provided in **Table 28**.

Table 28. Retrospective ANS Specimens. Demographic Summary of Compliant Subjects

Demographics	Characteristics	Total (N=187)
Gender	Female	50.8% (95/187)
	Male	49.2% (92/187)
Age Group	0 - 5 years	8.6% (16/187)
	6 - 21 years	16.0% (30/187)
	22 - 59 years	62.0% (116/187)
	> 59 years	13.4% (25/187)

A summary of the BD Respiratory Viral Panel for BD MAX System assay retrospective ANS specimen clinical study performance is provided in **Table 29**.

Table 29. BD Respiratory Viral Panel for BD MAX System Flu B and RSV Performance with Retrospective ANS Specimens

Analyte ¹		Positive Percent Agreement			Negative Percent Agreement		
		TP/ (TP+FN)	%	95% CI	TN/ (TN+FP)	%	95% CI
Flu B	Retrospective	12/12	100%	75.7%-100%	172/174	98.9%	95.9%-99.7%
RSV	Retrospective	15/15	100%	79.6%-100%	170/171	99.4%	96.8%-99.9%

TP – true positive; FN – false negative; TN – true negative; FP – false positive

¹One specimen generated an Unresolved (UNR) assay result for both Flu B and RSV. An UNR is generated when the target result and the Sample Processing Control (SPC) are both negative and no instrument error was recorded.

In addition, clinical performance characteristics of the BD Respiratory Viral Panel for BD MAX System were also assessed from a total of 240 frozen retrospective nasopharyngeal swabs (NPS) in UVT/UTM obtained from two (2) external sources with historical positive or negative results for either influenza B or RSV. The specimens were collected as part of routine patient care between December 2019 and January 2022. All the specimens were tested in a blinded and randomized fashion with the BD Respiratory Viral Panel for BD MAX System at least three different testing sites and FDA-cleared high sensitivity RT-PCR assay as comparator at one reference testing site. **Table 30** summarizes the demographic information (gender and age group) by percentage for compliant subjects for retrospective NPS specimens.

Table 30. Retrospective NPS Specimens: Demographic Summary of Compliant Subjects

Demographics	Characteristics	Total (N=240)
Gender	Female	37.5% (90/240)
	Male	37.5% (90/240)
	Unknown	25.0% (60/240)
Age Group	0 - 5 years	26.7% (64/240)
	6 - 21 years	20.8% (50/240)

22 - 59 years	32.5% (78/240)
> 59 years	20.0% (48/240)

A summary of the BD Respiratory Viral Panel for BD MAX System assay retrospective NPS clinical study performance is provided in **Table 31**.

Table 31. BD Respiratory Viral Panel for BD MAX System Flu B and RSV Performance with Retrospective NPS Specimens

Analyte		Positive Percent Agreement			Negative Percent Agreement		
		TP/ (TP+FN)	%	95% CI	TN/ (TN+FP)	%	95% CI
Flu B	Retrospective	58/58	100%	93.8%-100%	180/182	98.9%	96.1%-99.7%
RSV	Retrospective	62/63	98.4%	91.5%-99.7%	177/177	100%	97.9%-100%

TP – true positive; FN – false negative; TN – true negative; FP – false positive

1. Other Clinical Supportive Data (When 1. and 2. Are Not Applicable):
Not applicable.

D Clinical Cut-Off:
Not applicable.

E Expected Values/Reference Range:

The BD MAX RVP SARS-CoV-2/Flu A/B/RSV assay prospective clinical study included a total of 1562 prospectively collected NPS specimens and 1566 prospectively collected ANS specimens. The stratification of positivity rates by analyte (as determined by the BD MAX RVP SARS-CoV-2/Flu A/B/RSV test) and location for the prospective samples is shown in **Table 32**. Sample collection was carried out at 8 geographic locations. For the CTMD site, the same principal investigator had oversight of three clinics (two (2) in Florida and one (1) in North Carolina). Data from those three (3) clinics were aggregated.

Table 32. Prospective Specimens¹: Positivity Rates based on BD Max RVP

Target	Collection Site ²	NPS	ANS
SARS-CoV-2	BUL	60.9% (28/46)	59.6% (28/47)
	CMR	18.9% (62/328)	17.1% (56/328)
	CTMD	39.9% (311/780)	37.2% (291/783)
	MAC	38.0% (101/266)	35.5% (94/265)
	ROM	28.9% (37/128)	25.8% (33/128)
	VHMC	14.3% (2/14)	13.3% (2/15)
	Overall		34.6% (541/1562)
Flu A	BUL	0.0% (0/46)	0.0% (0/47)
	CMR	3.4% (11/328)	3.0% (10/328)
	CTMD	6.7% (52/780)	6.6% (52/783)
	MAC	0.0% (0/266)	0.0% (0/265)
	ROM	0.0% (0/128)	0.0% (0/128)
	VHMC	21.4% (3/14)	26.7% (4/15)

Target	Collection Site ²	NPS	ANS
	Overall	4.2% (66/1562)	4.2% (66/1566)
Flu B	BUL	0.0% (0/46)	0.0% (0/47)
	CMR	0.0% (0/328)	0.0% (0/328)
	CTMD	0.1% (1/780)	0.0% (0/783)
	MAC	0.0% (0/266)	0.0% (0/265)
	ROM	0.0% (0/128)	0.0% (0/128)
	VHMC	0.0% (0/14)	0.0% (0/15)
	Overall	0.1% (1/1562)	0.0% (0/1566)
RSV	BUL	0.0% (0/46)	0.0% (0/47)
	CMR	0.9% (3/328)	0.9% (3/328)
	CTMD	1.0% (8/780)	1.1% (9/783)
	MAC	0.0% (0/266)	0.0% (0/265)
	ROM	0.8% (1/128)	0.0% (0/128)
	VHMC	0.0% (0/14)	0.0% (0/15)
	Overall	0.8% (12/1562)	0.8% (12/1566)

¹Calculations are based on specimens that are compliant at the subject level, and compliant and reportable for MAX, regardless of the comparator compliancy and results.

²BUL (Sofia, Bulgaria), CMR (Providence, RI), CTMD (Palm Springs, FL / Ashville, NC / Coconut Creek, FL), MAC (Savannah, GA), ROM (Bucuresti, Romania), VHMC (Phoenix, AZ)

F Other Supportive Instrument Performance Characteristics Data:

Not Applicable

VIII Proposed Labeling:

The labeling supports the finding of substantial equivalence for this device.

IX Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.