

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY AND INSTRUMENT**

I Background Information:

A 510(k) Number

K233606

B Applicant

Biokit S.A.

C Proprietary and Established Names

ADVIA Centaur EBV-VCA IgM

D Regulatory Information

Product Code(s)	Classification	Regulation Section	Panel
LSE	Class I, reserved	21 CFR 866.3235 - Epstein-Barr Virus Serological Reagents	MI - Microbiology

II Submission/Device Overview:

A Purpose for Submission:

To obtain clearance of a new device.

B Measurand:

IgM antibodies to the viral capsid antigen (VCA) of Epstein Barr virus (EBV).

C Type of Test:

Fully automated qualitative chemiluminescent immunoassay.

III Intended Use/Indications for Use:

A Intended Use(s):

The ADVIA Centaur EBV-VCA IgM (EBVM) assay is for in vitro diagnostic use in the qualitative detection of IgM antibodies to the viral capsid antigen (VCA) of the Epstein-Barr virus (EBV) in human pediatric (2-21 years old) and adult serum and plasma (EDTA and lithium heparin) using the ADVIA Centaur XP system. When used in conjunction with other EBV markers, this assay is intended for use as an aid in the diagnosis of Epstein-Barr virus infection, such as infectious mononucleosis.

B Indication(s) for Use:

Not applicable

C Special Conditions for Use Statement(s):

Rx - For Prescription Use Only

D Special Instrument Requirements

The assay is intended for use on the fully automated ADVIA Centaur XP manufactured by Siemens Healthcare Diagnostics Inc.

IV Device/System Characteristics:**A Device Description:**

The ADVIA Centaur EBV-VCA IgM (EBVM) is a qualitative, serological, two-step sandwich chemiluminescent immunoassay performed on the fully automated ADVIA Centaur XP system intended for use in conjunction with other EBV markers as an aid in the diagnosis of Epstein-Barr virus infection, such as infectious mononucleosis.

The ADVIA Centaur EBVM assay kit consists of the ReadyPack primary reagent pack containing Lite Reagent, Solid Phase, and Ancillary Well Reagent, The ReadyPack ancillary reagent pack containing Ancillary Reagent (diluent), and low and high calibrators. External positive and negative controls are required to perform the assay and are available separately.

B Principle of Operation:

The ADVIA Centaur EBV-VCA IgM assay is a fully automated two-step sandwich immunoassay using acridinium ester chemiluminescent technology. The specimen is incubated with the Ancillary Well Reagent and the Solid Phase, which contains an EBV-VCA IgM specific antigen, Viral Capsid Antigen (VCA) p18 (biotinylated peptide) coated onto streptavidin magnetic microparticles. Antigen-antibody complexes will form if anti EBV-VCA IgM antibody is present in the specimen. The Lite Reagent contains monoclonal anti-human IgM labeled with acridinium ester and is used to detect EBV-VCA IgM in the specimen.

Index values for samples are determined automatically by software as ratio of the sample relative light units (RLU) to the cutoff RLU and reported in Index Values and as Nonreactive or Reactive:

- Nonreactive: < 1.00 Index. These samples are considered negative.
- Reactive: ≥ 1.00 Index. These samples are considered positive.

C Instrument Description Information:**1. Instrument Name:**

ADVIA Centaur XP

2. Specimen Identification:

Serum, EDTA plasma, and lithium heparin plasma are the validated specimen types for the ADVIA Centaur EBV-VCA IgM assay.

3. Specimen Sampling and Handling:

Samples in the primary collection device (serum stored on the clot, plasma stored on packed red cells, and samples processed and stored in gel-barrier blood collection tubes) are stable for up to 7 days at 2 - 8°C. Separated samples are stable for up to 4 hours at room temperature, up to 7 days at 2 - 8°C, and up to one month at $\leq -20^{\circ}\text{C}$ with no more than one freeze/thaw cycle. The assay requires 20 μL of sample for a single determination. The instrument performs all handling procedures automatically.

4. Calibration:

The ADVIA Centaur EBVM assay kit includes two calibrators (high and low). Calibration is performed as follows:

- At the end of the 28-day calibration interval.
- When changing lot numbers of primary reagent packs.
- When indicated by quality control results.
- After major maintenance or service, as indicated by quality control results

Before initiating calibration on each new lot of reagent, the assay master curve values need to be entered by scanning the master curve card.

5. Quality Control:

External positive and negative controls are available separately. The controls should be used at least once during each day that samples are analyzed. In addition, quality control is performed:

- Following a valid calibration
- With use of a new lot of reagent
- When troubleshooting test results that do not match clinical conditions or symptoms

Acceptable performance is achieved when the analyte values obtained are within the expected control interval for the system.

V Substantial Equivalence Information:

A Predicate Device Name(s):

Diasorin Liaison EBV IgM Assay

B Predicate 510(k) Number(s):

K040120

C Comparison with Predicate(s):

Device & Predicate Device(s):	<u>K233606</u>	<u>K040120</u>
Device Trade Name	ADVIA Centaur EBV-VCA IgM	LIAISON EBV IgM
Intended Use	The ADVIA Centaur EBV-VCA IgM (EBVM) assay is for in vitro diagnostic use in the qualitative detection of IgM antibodies to the viral capsid antigen (VCA) of the Epstein-Barr virus (EBV) in human pediatric (2-21 years old) and adult serum and plasma (EDTA and lithium heparin) using the ADVIA Centaur XP system. When used in conjunction with other EBV markers, this assay is intended for use as an aid in the diagnosis of Epstein-Barr virus infection, such as infectious mononucleosis.	The LIAISON EBV IgM assay uses chemiluminescent immunoassay (CLIA) technology on the LIAISON Analyzer family for the qualitative determination of specific IgM antibodies to Epstein-Barr virus (EBV) viral capsid antigen (VCA) p18 synthetic peptide in human serum. When performed in conjunction with other EBV markers, this assay can be used as an aid in the clinical laboratory diagnosis of Epstein-Barr Viral Syndrome in patients with signs and symptoms of EBV infection such as infectious mononucleosis.
General Device Characteristic Similarities		
Measurand	Same	IgM antibodies to the viral capsid antigen (VCA) of Epstein Barr virus (EBV).
Technology	Same	Chemiluminescent immunoassay (CLIA)
Antigen	Same	p18
General Device Characteristic Differences		
Sample Type	Human serum and plasma (EDTA and lithium heparin)	Human serum
Target Population	Adults and pediatrics	Adults
	As an aid in the diagnosis of Epstein-Barr virus infection, such as infectious mononucleosis	As an aid in the clinical laboratory diagnosis of Epstein-Barr Viral Syndrome in patients with signs and symptoms of EBV infection such as infectious mononucleosis

VI Standards/Guidance Documents Referenced:

- CLSI EP05-A3 Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - Third Edition (Reaffirmed: September 2019).
- CLSI EP09c Measurement Procedure Comparison and Bias Estimation Using Patient Samples - 3rd Edition (2018).
- CLSI EP07 Interference Testing in Clinical Chemistry; Approved Guideline - 3rd Edition (2018).
- CLSI EP37 Supplemental Tables for Interference Testing in Clinical Chemistry - 1st Edition (2018).
- CLSI EP25-A (Replaces EP25-P) Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline (2009).
- CLSI EP12-A2 User Protocol For Evaluation of Qualitative Test Performance; Approved Guideline - 2nd Edition (2008).
- CLSI I/LA21-A2 Clinical Evaluation of Immunoassays; Approved Guideline - 2nd Edition.

VII Performance Characteristics (if/when applicable):

A Analytical Performance:

1. Precision/Reproducibility:

- 1) Precision study was performed using three lots of the ADVIA Centaur EBV-VCA IgM with two runs per day for 20 days measuring duplicates of each external control and ten contrived samples with levels of analyte ranging from negative to high positive. The within laboratory precision study results for one representative reagent lot are shown in **Table 1**.

Table 1. Precision study results. SD: standard deviation of mean value (Index); %CV: coefficient of variation.												
Sample ID	Mean Value (Index)	N	Repeatability		Between-Run		Between-Day		Between-Lot		Total Precision	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
Control 1	0.28	240	0.012	N/A	0.014	N/A	0.010	N/A	0.011	N/A	0.024	N/A
Control 2	3.02	240	0.071	2.3%	0.049	1.6%	0.088	2.9%	0.214	7.1%	0.247	8.2%
Serum A	0.75	240	0.019	2.5%	0.018	2.4%	0.017	2.3%	0.081	10.9%	0.087	11.6%
Serum B	0.96	240	0.024	2.5%	0.023	2.4%	0.023	2.4%	0.076	8.0%	0.086	9.0%
Serum C	1.39	240	0.034	2.4%	0.030	2.2%	0.041	2.9%	0.119	8.5%	0.133	9.6%
Serum D	3.16	240	0.075	2.4%	0.053	1.7%	0.101	3.2%	0.173	5.5%	0.221	7.0%
Serum E	7.22	240	0.217	3.0%	0.140	1.9%	0.269	3.7%	0.531	7.4%	0.649	9.0%
Plasma, EDTA A	0.74	240	0.019	2.6%	0.026	3.4%	0.023	3.0%	0.048	6.4%	0.062	8.3%
Plasma, EDTA B	1.00	240	0.022	2.2%	0.031	3.1%	0.027	2.7%	0.061	6.1%	0.077	7.7%
Plasma, EDTA C	1.40	240	0.029	2.1%	0.045	3.2%	0.041	2.9%	0.118	8.4%	0.136	9.7%

Plasma, EDTA D	3.35	240	0.081	2.4%	0.062	1.8%	0.111	3.3%	0.308	9.2%	0.342	10.2%
Plasma, EDTA E	7.63	240	0.218	2.9%	0.367	4.8%	0.320	4.2%	0.466	6.1%	0.709	9.3%

- 2) Reproducibility of the ADVIA Centaur EBV-VCA IgM assay was assessed at three external testing sites using one reagent lot and three ADVIA Centaur XP instruments (one per testing site) evaluating external controls and eight samples with levels of analyte ranging from high negative to positive over five days with two runs per day and triplicate measurements. Study results are presented in **Table 2**.

Table 2. Reproducibility study results.											
		Within-Run Repeatability		Between- Run		Between- Day		Between- Site		Reproducibility	
Sample	Mean (Index)	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	CV (%)
Serum A	0.85	0.023	2.7%	0.028	3.3%	0.004	0.5%	0.017	2.0%	0.040	4.7%
Serum B	1.00	0.024	2.4%	0.020	2.1%	0.008	0.8%	0.011	1.1%	0.034	3.5%
Serum C	1.57	0.038	2.4%	0.022	1.4%	0.014	0.9%	0.028	1.8%	0.054	3.4%
Serum D	3.30	0.099	3.0%	0.025	0.8%	0.009	0.3%	0.081	2.5%	0.131	4.0%
Plasma EDTA A	0.75	0.023	3.0%	0.026	3.4%	0.000	0.0%	0.013	1.8%	0.037	4.9%
Plasma EDTA B	1.00	0.035	3.5%	0.023	2.3%	0.000	0.0%	0.009	0.9%	0.043	4.3%
Plasma EDTA C	1.49	0.036	2.4%	0.038	2.5%	0.000	0.0%	0.016	1.1%	0.055	3.7%
Plasma EDTA D	3.72	0.103	2.8%	0.088	2.4%	0.054	1.4%	0.106	2.8%	0.180	4.9%
Control 1	0.24	0.012	N/A	0.004	N/A	0.006	N/A	0.009	N/A	0.017	N/A
Control 2	3.20	0.248	7.8%	0.000	0.0%	0.075	2.3%	0.145	4.5%	0.297	9.3%

2. Linearity

Not Applicable

3. Analytical Specificity/Interference:

1) Cross-reactivity

Biokit evaluated the influence of potentially cross-reacting antibodies to infections other than EBV (as well as autoimmune disorders) on the ADVIA Centaur EBV-VCA IgM. A total of 353 samples confirmed to be EBV IgM-negative on a comparative assay and positive for antibodies for 32 different diseases/conditions were tested in singlicate with the candidate device on ADVIA Centaur XP system. The study results are presented in **Table 3**.

Table 3. Results of cross-reactivity study.

Clinical Category	Number Tested	Candidate Device Results		Comparative Assay Results		
		-	+	-	Equivocal	+
Cytomegalovirus (CMV) IgG	10	10	0	10	0	0
Cytomegalovirus (CMV) IgM	10	7	3	7	0	3
Parvovirus B19 IgM	10	10	0	9	0	1
<i>Toxoplasma gondii</i> IgG	10	10	0	10	0	0
<i>Toxoplasma gondii</i> IgM	28	22	6	24	0	4
Rubella IgG	20	18	2	19	1	0
Rubella IgM	11	11	0	10	1	0
Hepatitis B Virus (HBV) IgM	11	11	0	11	0	0
Hepatitis A Virus (HAV) IgM	10	7	3	8	0	2
Hepatitis C Virus (HCV)	11	11	0	10	0	1
Human Immunodeficiency Virus (HIV)	10	9	1	9	0	1
Herpes Simplex Virus (HSV1) IgG	11	11	0	11	0	0
Herpes Simplex Virus (HSV1) IgM	10	10	0	10	0	0
Herpes Simplex Virus (HSV2) IgG	12	12	0	12	0	0
Herpes Simplex Virus (HSV2) IgM	11	9	2	10	0	1
<i>Treponema pallidum</i> (Syphilis)	11	11	0	11	0	0
Varicella Zoster Virus (VZV) IgM	10	9	1	10	0	0
Measles virus IgM	10	10	0	10	0	0
Mumps virus IgM	10	10	0	10	0	0
<i>Borrelia burgdorferi</i> IgM (Lyme IgM)	11	11	0	11	0	0
Influenza virus IgM	10	10	0	10	0	0
<i>Mycoplasma pneumoniae</i> IgM	18	15	3	16	0	2
Antinuclear antibodies (ANA)	10	10	0	10	0	0
Rheumatic Factors (RF)	11	10	1	11	0	0
Human Anti-mouse antibodies (HAMA)	10	10	0	10	0	0
Human Herpes Virus (HHV6)	14	14	0	14	0	0
Systemic Lupus Erythematosus (SLE)	10	10	0	10	0	0
Flu vaccinated patients	10	10	0	10	0	0
Elevated IgG	10	10	0	10	0	0
Elevated IgM	11	11	0	11	0	0
Epstein Barr Virus (EBV) Viral Capsid Antigen (VCA) IgG	10	10	0	10	0	0
Epstein Barr Virus (EBV) Epstein-Barr Nuclear Antigen (EBNA) IgG	10	9	1	9	0	1
Total	371	348	23	353	2	16

2) Interference

Biokit evaluated the effect of the presence of commonly found endogenous substances on the performance of the ADVIA Centaur EBV-VCA IgM using negative, low positive, and high positive serum, EDTA plasma, and lithium heparin plasma contrived samples. No interference (false positive or false negative) was observed at the tested concentrations reported in **Table 4**.

Table 4. Interference study results.	
Substance	Substance Test Concentration
Hemoglobin	1000 mg/dL
Bilirubin, conjugated	40 mg/dL
Bilirubin, unconjugated	40 mg/dL
Intralipid	1500 mg/dL
Biotin	3510 ng/mL
Cholesterol	502 mg/dL
Protein (hyperproteinemic)	15 g/dL
Protein (hypoproteinemic)	3 g/dL

3) Class Specificity

Provided data supports that ADVIA Centaur EBV-VCA IgM assay selectively detects EBV IgM antibodies.

4. Assay Reportable Range

Not Applicable

5. Stability:

1) Reagent Stability

The data support storage conditions for the kit reagents as listed in **Table 5**.

Table 5. Reagent stability.		
Reagent	Storage Condition	Stability
Primary and ancillary reagent packs	Unopened at 2 – 8°C	12 months
	Onboard	28 days
Calibrators	Unopened at 2 – 8°C	12 months
	Opened at 2 – 8°C	60 days
	At room temperature	8 hours
Washes	Unopened at 2 – 25°C	12 months
	Onboard	1 month
Controls	Opened at 2 – 8°C	60 days

The data also support calibration interval stability of 28 days.

2) Sample Stability

The data support storage conditions for samples collected in serum, lithium heparin plasma, and EDTA plasma listed in **Table 6**.

Table 6. Summary of serum, EDTA plasma, and lithium heparin plasma stability.		
Storage Condition	Stability	
After centrifugation at 2 – 8°C	7 days	
Separated samples at room temperature	4 hours	
Separated samples at 2 – 8°C	7 days	
Separated samples at -20°C	1 month	
Freeze/thaw cycles	1	

6. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):
Not applicable.

7. Detection Limit:
Not applicable.

8. Assay Cut-Off:

A study to establish the ADVIA Centaur EBV-VCA IgM cut-off was performed using 231 clinical human leftover remnant and intended use serum samples. The cut-off value was determined by performing concordance and Receiver Operator Characteristic (ROC) analysis, using predicate (DiaSorin Liaison EBV IgM) results as the standard. The estimated cut-off value for ADVIA Centaur EBV-VCA IgM assay is 1.0 Index value.

9. Carry-Over:

The ADVIA Centaur EBV-VCA IgM is not susceptible to carry-over.

B Comparison Studies:

1. Method Comparison with Predicate Device:

Biokit conducted a multisite study to evaluate clinical performance of the ADVIA Centaur EBV-VCA IgM on the ADVIA Centaur XP system by determining Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) with the predicate device, DiaSorin Liaison EBV IgM. A total of 1,428 leftover samples were collected over a contiguous time period from individuals for whom an EBV test was ordered (population 1). Of these, 188 were from unclassified serostatus individuals. Two hundred and two (202) samples with a known positive result (population 2) were evaluated as well. Of these 4 were from unclassified serostatus individuals. The study results otherwise showed that these populations included individuals with acute infection, past infection, or no serologic evidence of EBV infection. Testing was conducted at three sites in the U.S. representative of testing environments where the candidate device will be used. The results are presented in **Tables 7 and 8**:

Table 7. Results of the method comparison study evaluating population 1.				
ADVIA Centaur EBV-VCA IgM	LIAISON EBV IgM			
	Negative	Equivocal	Positive	Total
Nonreactive	1294	2	26	1322
Reactive	41	3	62	106
Total	1335	5	88	1428
	NPA = 96.71% (1294/1338) 95% CI: 95.61% - 97.54%		PPA = 68.89% (62/90) 95% CI: 58.72% - 77.51%	

Out of 90 samples that were positive on the reference assay, 53 were primary acute, and of those, 48 samples were positive on the ADVIA Centaur EBVM assay. The PPA in primary acute patients is therefore 90.57% with 95% CI of 79.75% - 95.90%. Primary acute in this study is defined as the presence of either EBV IgM or heterophile antibodies, and absence of EBNA IgG.

Table 8. Results of the method comparison study evaluating population 2.				
ADVIA Centaur EBV-VCA IgM	LIAISON EBV IgM			
	Negative	Equivocal	Positive	Total
Nonreactive	0	0	0	0
Reactive	0	0	202	202
Total	0	0	202	202
			PPA = 100.00% (202/202) 95% CI: 98.13% - 100.00%	

The PPA and NPA in pediatric subjects are presented in **Tables 9** and **10**. For population 1, 84 samples were from unclassified serostatus individuals. For population 2, three samples were from unclassified serostatus individuals.

Table 9. Results of the method comparison study evaluating the pediatric population in population 1.				
ADVIA Centaur EBV-VCA IgM	LIAISON EBV IgM			
	Negative	Equivocal	Positive	Total
Nonreactive	402	0	9	411
Reactive	19	1	48	68
Total	421	1	57	479
	NPA = 95.26% (402/422) 95% CI: 92.79% - 96.91%		PPA = 84.21% (48/57) 95% CI: 72.64% - 91.47%	

Out of 57 pediatric subject samples that were positive on the reference assay, 46 were primary acute, and of those, 42 samples were positive on the ADVIA Centaur EBVM assay. The PPA in primary acute pediatric patients is therefore 91.30% with 95% CI of 79.68% -

96.57%. Primary acute in this study is defined as the presence of either EBV IgM or heterophile antibodies, and absence of EBNA IgG.

Table 10. Results of the method comparison study evaluating the pediatric population in population 2.				
ADVIA Centaur EBV-VCA IgM	LIAISON EBV IgM			
	Negative	Equivocal	Positive	Total
Nonreactive	0	0	0	0
Reactive	0	0	155	155
Total	0	0	155	155
			PPA = 100.00% (155/155) 95% CI: 97.58% - 100.00%	

2. Matrix Comparison:

Biokit performed this study to assess the performance of different sample matrices with the ADVIA Centaur EBV-VCA IgM assay. The study evaluated 70 sets of matched samples collected in serum separator tube (SST), EDTA plasma, and lithium heparin plasma measured in singlicate using one lot of reagents. The data supports comparable performance of samples collected in serum and EDTA and lithium heparin plasma (**Table 11**).

Table 11. Results of regression analysis performed for matrix comparison study.				
Tube (y) vs. Serum (x)	Regression Equation	Sample Interval	N^a	r^b
Plasma, EDTA	y=1.00x - 0.03 Index	0.12-12.35 Index	70	1.00
Plasma, lithium heparin	y=1.00x - 0.05 Index	0.10-11.49 Index	70	1.00

^a Number of samples tested.
^b Correlation coefficient.

C Clinical Studies:

See Section B. Comparison Studies for clinical study results.

D Clinical Cut-Off:

Not Applicable

E Expected Values/Reference Range:

The ADVIA Centaur EBV-VCA IgM results for population 1 (N=1,428 samples) for all sites combined by age, group, and sex are summarized in **Table 12**.

Table 12. Expected values established using the ADVIA Centaur XP system and confirmed by assay comparison.						
Age Group (Years)	Sex	Reactive		Nonreactive		Total
		N ^a	% ^b	N	%	
2-21	Male	32	14.1%	195	85.9%	227
	Female	36	14.3%	216	85.7%	252
	Overall	68	14.2%	411	85.8%	479
13-21	Male	25	18.5%	110	81.5%	135
	Female	32	18.7%	139	81.3%	171
	Overall	57	18.6%	249	81.4%	306
>21	Male	11	3.3%	324	96.7%	335
	Female	27	4.4%	587	95.6%	614
	Overall	38	4.0%	911	96.0%	949
Total	Male	43	7.7%	519	92.3%	562
	Female	63	7.3%	803	92.7%	866
	Overall	106	7.4%	1322	92.6%	1428

^a Number of samples tested.

^b Percentages are for the numbers of reactive and nonreactive in a given row.

VIII Proposed Labeling:

The labeling supports the finding of substantial equivalence for this device.

IX Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.