TOPCON TRC-NW200 510(K) SUMMARY

Name of Device:

Topcon Model TRC-NW200

Non-Mydriatic Retinal Camera

Common / Generic Device Name:

Retinal Camera

Classification Name:

Camera, Ophthalmic, AC-Powered

Regulation Number:

886-1120

Product Code:

HKI

Submitter:

Topcon Corporation

Address:

37 West Century Road, Paramus, NJ 07652

Telephone Number:

(201) 261-9450

Facsimile Number:

(201) 387-2710 Donald H. Winfield

Contact Person:

Date Prepared:

August 12, 2004

Predicated Devices:

Manufacturer: 1)

Canon, Inc.

Trade Name:

Canon

Model Name:

CR6-45NM

Classification Name:

86 HKI

510(k) Number:

K 980246

2) Manufacturer: Canon, Inc.

Trade Name:

Canon

Model Name:

CR-DGi

Classification Name:

86 HKI

510(k) Number:

K 031629

Manufacturer: 3)

Nidek, Inc.

Trade Name:

Nidek NM-1000

Model Name: Classification Name:

86 HKI

510(k) Number:

K 014274

Description of Device:

The Topcon TRC-NW200 is a non-mydriatic ophthalmic camera that captures and displays images of the eye. The images are captured by using halogen illumination with an infrared filter for observation and a xenon flash bulb for photography. There is no laser illumination. Alignment dots and focusing bars assist the operator in properly positioning the camera. The operator presses a button on the top of the joystick to capture the image. The image is then displayed on the built-in LCD display or, if desired, transferred to a personal computer.

Intended Use:

Topcon Model TRC-NW200 Non-Mydriatic Camera is intended for use in capturing images of the retina and the anterior segment of the eye and presenting the data to the eyecare professional, without the use of a mydriatic.

Substantial Equivalence:

The TRC-NW200 and the predicated devices have the same intended use and indications. All of these devices are non-mydriatic ophthalmic cameras that capture and display images of the eye to aid the health care professional in diagnosing and/or monitoring diseases of the eye that may be observed and photographed. The NW200 and the predicated devices acquire these images according to the same principles of operation and technical characteristics. The question of safety and effectiveness is the same for the NW200 and the predicated devices. Therefore, the TRC-NW200 is substantially equivalent to legally marketed devices intended to capture images of the retina and the anterior segment of the eye and presenting this data to the eyecare professional without the use of a mydriatic.



NOV - 5 2004

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Topcon Medical Systems, Inc. c/o Donald Winfield Vice President 37 West Century Road Paramus, NJ 07652

Re: K041367

Trade/Device Name: Topcon Non -Mydriatic Retinal Camera, Model TRC-NW200

Regulation Number: 21 CFR 886.1120 Regulation Name: Ophthalmic Camera

Regulatory Class: Class II

Product Code: HKI Dated: May 21, 2004 Received: May 24, 2004

Dear Mr. Winfield:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

A. Ralph Rosenthal, M.D.

Director

Division of Ophthalmic and Ear,

A Palpi forentbal

Nose and Throat Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if Known): <u>K041367</u>

Device Name:	Non-Mydriatic Retinal	Camera, Model T	<u>RC-NW200</u>	
Indication for Use:				
The Topcon Model TRC-NW200 Non-Mydriatic Camera is intended for use in capturing images of the retina and the anterior segment of the eye and presenting the data to the eyecare professional, without the use of a mydriatic.				
Prescription Us		AND/OR	Over The Country H-	
(Part 21 CFR 80		AND/OK	Over-The-Counter Use(21 CFR 807 Subpart C)	_
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)				
	Concurrence of CDF	RH, Office of Devi	ce Evaluation (ODE)	
Dexender 11/05/2004				
(Division Sign-Off) Division of Ophthalmic Ear,				
Nose and Throat Devises				