



**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY**

I Background Information:

A 510(k) Number

K260065

B Applicant

Guangzhou Wondfo Biotech Co., Ltd.

C Proprietary and Established Names

SAFElife™ T-Dip Multi-Drug Urine Test Panel
SAFElife™ T-Dip Multi-Drug Urine Test Panel Dx

D Regulatory Information

Product Code(s)	Classification	Regulation Section	Panel
NFT	Class II	21 CFR 862.3100	Toxicology (91)
NGL	Class II	21 CFR 862.3650	Toxicology (91)
PTH	Class II	21 CFR 862.3150	Toxicology (91)
NFV	Class II	21 CFR 862.3170	Toxicology (91)
NFY	Class II	21 CFR 862.3250	Toxicology (91)
PTG	Class II	21 CFR 862.3620	Toxicology (91)
NGG	Class II	21 CFR 862.3610	Toxicology (91)
LCM	Class II	Unclassified	Toxicology (91)
QBF	Class II	21 CFR 862.3700	Toxicology (91)
QAW	Class II	21 CFR 862.3910	Toxicology (91)
NFW	Class II	21 CFR 862.3870	Toxicology (91)

II Submission/Device Overview:

A Purpose for Submission:

Addition of new analytes to a previously cleared device

B Measurand:

6-Monoacetylmorphine
Amphetamine
Secobarbital
Buprenorphine
Oxazepam
Benzoylcegonine
2-ethylidene-1,5 dimethyl-3,3-diphenylpyrrolidine (EDDP)
Fentanyl
3,4 Methylenedioxymethamphetamine
D(+)-Methamphetamine
Morphine
Methadone
Norfentanyl
Oxycodone
Phencyclidine
d-Propoxyphene
Nortriptyline
11-nor- Δ^9 -THC-9-COOH
Tramadol

C Type of Test:

Qualitative, lateral flow immunoassay

III Intended Use/Indications for Use:

A Intended Use(s):

See Indications for Use below.

B Indication(s) for Use:

SAFElife™ T Dip Multi Drug Urine Test Panel

SAFElife™ T Dip Multi Drug Urine Test Panel is a competitive binding, lateral flow immunochromatographic assay for qualitative and simultaneous detection of 6-Monoacetylmorphine, Amphetamine, Buprenorphine, Secobarbital, Oxazepam, Cocaine, 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP), Fentanyl, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Opiates, Methadone, Norfentanyl, Oxycodone, Phencyclidine, Propoxyphene, Nortriptyline, Cannabinoids and Tramadol in human urine at the cutoff concentrations of:

Drug (Identifier)	Cutoff Level
6-Monoacetylmorphine (6-MAM)	10 ng/mL
Amphetamine (AMP)	1000 ng/mL or 500 ng/mL
Buprenorphine (BUP)	10 ng/mL

Secobarbital (BAR)	300 ng/mL
Oxazepam (BZO)	300 ng/mL
Cocaine (COC)	300 ng/mL or 150 ng/mL
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
Fentanyl (FTY)	1 ng/mL
Methamphetamine (MET/mAMP)	1000 ng/mL or 500 ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP)	300 ng/mL
Opiates (OPI)	2000 ng/mL
Methadone (MTD)	300 ng/mL
Norfentanyl (NFTY)	5 ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL
Propoxyphene (PPX)	300 ng/mL
Nortriptyline (TCA)	1000 ng/mL
Cannabinoids (THC)	50 ng/mL or 20 ng/mL
Tramadol (TRA)	100 ng/mL

SAFElife™ T-Dip Multi-Drug Urine Test Panel offers any combinations from 1 to 15 drugs of abuse tests but only one cutoff concentration under same drug condition will be included per device. It is for in vitro diagnostic use.

The tests may yield positive results for the prescription drugs Buprenorphine, Fentanyl, Nortriptyline, Oxazepam, Secobarbital, Oxycodone and Tramadol when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/Mass Spectrometry (LC/MS) is the recommended confirmatory method.

SAFElife™ T Dip Multi Drug Urine Test Panel Dx

SAFElife™ T Dip Multi Drug Urine Test Panel Dx is a competitive binding, lateral flow immunochromatographic assay for qualitative and simultaneous detection of 6 Monoacetylmorphine, Amphetamine, Secobarbital, Buprenorphine, Oxazepam, Cocaine, 2 ethylidene 1,5 dimethyl 3,3 diphenylpyrrolidine (EDDP), Fentanyl, Methylenedioxymethamphetamine, Methamphetamine, Morphine, Opiates, Methadone, Norfentanyl, Oxycodone, Phencyclidine, Propoxyphene, Nortriptyline, Cannabinoids and Tramadol in human urine with below cutoff concentrations and approximate detection time:

Drug (Identifier)	Calibrator	Cut-off
6-Monoacetylmorphine (6-MAM)	6-Monoacetylmorphine	10 ng/mL
Amphetamine (AMP500)	d-Amphetamine	500 ng/mL
Amphetamine (AMP1000)	d-Amphetamine	1000 ng/mL
Secobarbital (BAR)	Secobarbital	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Oxazepam (BZO)	Oxazepam	300 ng/mL

Cocaine (COC150)	Benzoyllecgonine	150 ng/mL
Cocaine (COC300)	Benzoyllecgonine	300 ng/mL
2-ethylidene-1,5 dimethyl-3, 3-diphenylpyrrolidine (EDDP)	2-ethylidene-1,5 dimethyl-3, 3-diphenylpyrrolidine	300 ng/mL
Fentanyl (FTY)	Fentanyl	1 ng/mL
Methylenedioxymethamphetamine (MDMA)	3,4-Methylenedioxymethamphetamine	500 ng/mL
Methamphetamine (MET500/mAMP500)	D (+)-Methamphetamine	500 ng/mL
Methamphetamine (MET1000/mAMP1000)	D (+)-Methamphetamine	1000 ng/mL
Morphine (MOP300)	Morphine	300 ng/mL
Opiates (OPI2000)	Morphine	2000 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Norfentanyl (NFTY)	Norfentanyl	5 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL
Propoxyphene (PPX)	d-Propoxyphene	300 ng/mL
Nortriptyline (TCA)	Nortriptyline	1000 ng/mL
Cannabinoids (THC20)	11-nor- Δ^9 -THC-9-COOH	20 ng/mL
Cannabinoids (THC50)	11-nor- Δ^9 -THC-9-COOH	50 ng/mL
Tramadol (TRA)	Tramadol	100 ng/mL

SAFElife™ T-Dip Multi-Drug Urine Test Panel Dx offers any combinations from 1 to 15 drugs of abuse tests with or without on-board adulteration/specimen validity test (SVT) but only one cutoff concentration under same drug condition will be included per device. It is for in vitro diagnostic use.

The tests may yield positive results for the prescription drugs Buprenorphine, Fentanyl, Nortriptyline, Oxazepam, Secobarbital, Oxycodone and Tramadol when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/Mass Spectrometry (LC/MS) is the recommended confirmatory method.

C Special Conditions for Use Statement(s):

OTC

D Special Instrument Requirements:

N/A

IV Device/System Characteristics:

A Device Description:

SAFElife™ T-Dip Multi-Drug Urine Test Panel and SAFElife™ T-Dip Multi-Drug Urine Test Panel Dx are immunochromatographic assays that use a lateral flow system for the qualitative detection of Amphetamine, Buprenorphine, Secobarbital, Oxazepam, Cocaine, 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Opiates, Methadone, Oxycodone, Phencyclidine, Propoxyphene, Nortriptyline, Cannabinoids, 6-Monoacetylmorphine, Fentanyl, Norfentanyl, and Tramadol in human urine. Each SAFElife™ T-Dip Multi-Drug Urine Test Panel or SAFElife™ T-Dip Multi-Drug Urine Test Panel Dx device consists of a test panel and a package insert. Each test panel is sealed with one desiccant in an aluminum pouch.

B Principle of Operation:

SAFElife™ T-Dip Multi-Drug Urine Test Panel and SAFElife™ T-Dip Multi Drug Urine Test Panel Dx are competitive immunochromatographic assays that are used to screen for the presence of drugs of abuse in urine. During testing, the urine specimen mixes with the drug specific monoclonal antibody conjugate, and flows across the membrane. When sample drug levels are zero or below the target cutoff, drug monoclonal antibody conjugate binds to the respective drug protein conjugate immobilized in the Test Region (T). This produces a colored band in the Test Region (T) that, regardless of its intensity, indicates a negative result. When drug concentration in the sample is at or above the target cutoff, the drug in the sample binds to the respective drug monoclonal antibody conjugate preventing the respective drug monoclonal antibody conjugate from binding to the respective drug protein conjugate immobilized in the Test Region (T). This prevents the development of a colored band in the Test Region (T), indicating a potentially positive result.

V Substantial Equivalence Information:

A Predicate Device Name(s):

Wondfo T-Dip® Multi-Drug Urine Test Panel
Wondfo T-Dip® Multi-Drug Urine Test Panel Rx

B Predicate 510(k) Number(s):

K202567

C Comparison with Predicate(s):

Device & Predicate Device(s):	K260065	Predicate (K202567)
Device Trade Name	SAFElife™ T-Dip Multi-Drug Urine Test Panel SAFElife™ T-Dip Multi-Drug Urine Test Panel Dx	Wondfo T-Dip® Multi-Drug Urine Test Panel Wondfo T-Dip® Multi-Drug Urine

		Test Panel Rx	
General Device Characteristic Similarities			
Intended Use/Indications For Use	For the qualitative determination of drugs of abuse in human urine	Same	
Methodology	Competitive binding, lateral flow immunochromatographic assay	Same	
Type of Test	Qualitative	Same	
Specimen Type	Human urine	Same	
Target Drug and Cut Off Values	Target Drugs	Cutoff	Same except that no FTY1, NFTY5, THC20, and TRA100
	Amphetamine (AMP)	1000 ng/mL or 500 ng/mL	
	Buprenorphine (BUP)	10 ng/mL	
	Secobarbital (BAR)	300 ng/mL	
	Oxazepam (BZO)	300 ng/mL	
	Cocaine (COC)	300 ng/mL or 150 ng/mL	
	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL	
	Methamphetamine (MET/mAMP)	1000 ng/mL or 500 ng/mL	
	Methylenedioxymethamphetamine (MDMA)	500 ng/mL	
	Morphine (MOP)	300 ng/mL	
	Opiates (OPI)	2000 ng/mL	
	Methadone (MTD)	300 ng/mL	
	Oxycodone (OXY)	100 ng/mL	
	Phencyclidine (PCP)	25 ng/mL	
	Propoxyphene (PPX)	300 ng/mL	
	Nortriptyline (TCA)	1000 ng/mL	
Cannabinoids (THC)	50 ng/mL		
6-Monoacetylmorphine	10 ng/mL		

	(6-MAM)	
	Fentanyl (FTY)	1 ng/mL
	Norfentanyl (NFTY)	5 ng/mL
	Cannabinoids (THC)	20 ng/mL
	Tramadol (TRA)	100 ng/mL
General Device Characteristic Differences		
Target Drug and Cut Off Values	See above	No FTY1, NFTY5, THC20, and TRA100

VI Standards/Guidance Documents Referenced:

N/A

VII Performance Characteristics (if/when applicable):

A Analytical Performance:

1. Precision/Reproducibility:

Precision studies were carried out for samples with concentrations of -100% cutoff, -75% cut off, -50% cutoff, -25% cutoff, cutoff, +25% cutoff, +50% cutoff, +75% cutoff and +100% cutoff for 6-Monoacetylmorphine (6-MAM10), Fentanyl (FTY1), Norfentanyl (NFTY5), Cannabinoids (THC20), and Tramadol (TRA100). Samples with concentration of -100% cutoff were drug-free urine samples. Other samples were prepared by spiking target drug in drug-free urine samples. Each drug concentration was confirmed by LC/MS or GC/MS. For each concentration, tests were performed two runs per day for 25 days using three lots of test panels. The results of a representative lot are summarized in the following tables. Precision results for AMP500/1000, BAR300, BUP10, BZO300, COC150/300, EDDP300, MDMA500, MET500/MET1000, MOP300, MTD300, OPI2000, OXY100, PCP25, PPX300, TCA1000 and THC50 were previously reviewed under K202567.

Drug	Lot Number	+100% cutoff	+75% cutoff	+50% cutoff	+25% cutoff	Cutoff	-25% cutoff	-50% cutoff	-75% cutoff	-100% cut-off
6AM 10	Lot 1	50-/0+	50-/0+	50-/0+	49-/1+	25-/25+	1-/49+	0-/50+	0-/50+	0-/50+
	Lot 2	50-/0+	50-/0+	50-/0+	49-/1+	26-/24+	2-/48+	0-/50+	0-/50+	0-/50+
	Lot 3	50-/0+	50-/0+	50-/0+	48-/2+	25-/25+	2-/48+	0-/50+	0-/50+	0-/50+
FTY1	Lot 1	50-/0+	50-/0+	50-/0+	49-/1+	24-/26+	0-/50+	0-/50+	0-/50+	0-/50+
	Lot 2	50-/0+	50-/0+	50-/0+	48-/2+	25-/25+	0-/50+	0-/50+	0-/50+	0-/50+
	Lot 3	50-/0+	50-/0+	50-/0+	49-/1+	27-/23+	0-/50+	0-/50+	0-/50+	0-/50+
	Lot 1	50-/0+	50-/0+	50-/0+	48-/2+	25-/25+	1-/49+	0-/50+	0-/50+	0-/50+

NFTY 5	Lot 2	50-/0+	50-/0+	50-/0+	49-/1+	24-/26+	1-/49+	0-/50+	0-/50+	0-/50+
	Lot 3	50-/0+	50-/0+	50-/0+	49-/1+	26-/24+	2-/48+	0-/50+	0-/50+	0-/50+
THC20	Lot 1	50-/0+	50-/0+	50-/0+	48-/2+	25-/25+	2-/48+	0-/50+	0-/50+	0-/50+
	Lot 2	50-/0+	50-/0+	50-/0+	49-/1+	27-/23+	1-/49+	0-/50+	0-/50+	0-/50+
	Lot 3	50-/0+	50-/0+	50-/0+	50-/0+	24-/26+	1-/49+	0-/50+	0-/50+	0-/50+
TRA100	Lot 1	50-/0+	50-/0+	50-/0+	49-/1+	25-/25+	1-/49+	0-/50+	0-/50+	0-/50+
	Lot 2	50-/0+	50-/0+	50-/0+	49-/1+	26-/24+	2-/48+	0-/50+	0-/50+	0-/50+
	Lot 3	50-/0+	50-/0+	50-/0+	48-/2+	25-/25+	2-/48+	0-/50+	0-/50+	0-/50+

2. Linearity:

Not applicable. This device is intended for qualitative use only.

3. Analytical Specificity/Interference:

Cross-reactivity/Interference of AMP500/1000, BAR300, BUP10, BZO300, COC150/300, EDDP300, MDMA500, MET500/MET1000, MOP300, MTD300, OPI2000, OXY100, PCP25, PPX300, TCA1000 and THC50 were previously reviewed under K202567. In this submission, 6-Monoacetylmorphine (6-MAM10), Fentanyl (FTY1), Norfentanyl (NFTY5), Cannabinoids (THC20), and Tramadol (TRA100) study data are summarized and presented below.

Cross-reactivity:

To test the specificity, drug metabolites and other components that are likely to cross-react in urine samples were spiked into drug-free urine. The relative cross-reactivity represents the minimum concentration necessary to yield a result similar to the cutoff level of the respective assay. Percent cross-reactivity, provided in the below table, was calculated as the concentration of analyte tested that yielded a positive result, divided by the cutoff concentration, multiplied by 100; compounds that did not yield a positive result at the highest concentration tested have relative cross reactivity results. The results obtained are summarized in tables below.

Drug/Cutoff	Compound	Minimum concentration required to obtain a positive result (ng/mL)	% Cross-reactivity
6-MAM	Heroin	60	16.7%
	Morphine	75000	0.01%
	Normorphine	150000	Not detected
	Nalorphine HCl	150000	Not detected
	Hydrocodone	150000	Not detected
	Hydromorphone	150000	Not detected
	Chlordiazepoxide	150000	Not detected
	Clobazam	150000	Not detected
	D-Amphetamine	150000	Not detected
	(±)-Amphetamine	150000	Not detected

	Levorphanol tartrate	150000	Not detected
	Codeine	150000	Not detected
	Ethylmorphine	150000	Not detected
	Morphine-3- β -D-glucuronide	150000	Not detected
	Norcodeine	150000	Not detected
	Oxycodone	150000	Not detected
	Oxymorphone	150000	Not detected
	Procaine hydrochloride	150000	Not detected
	Thebaine	150000	Not detected
	6-Acetylcodeine	150000	Not detected
	Buprenorphine	150000	Not detected
	Dihydrocodeine	150000	Not detected
	Dextromethorphan	150000	Not detected
	Imipramine hydrochloride	150000	Not detected
	Meperidine	150000	Not detected
	(\pm)-Methadone	150000	Not detected
	Mitragynine(kratom)	150000	Not detected
	Morphine-6- β -D-glucuronide	150000	Not detected
	Naloxone hydrochloride	150000	Not detected
	Naltrexone hydrochloride	150000	Not detected
	Naproxen	150000	Not detected
	Norbuprenorphine	150000	Not detected
	Norbuprenorphine glucuronide	150000	Not detected
	Noroxycodone HCL	150000	Not detected
	Noroxymorphone HCL	150000	Not detected
	(+)-Norpropoxyphene maleate	150000	Not detected
	Oxymorphone-3 β -D-glucuronide	150000	Not detected
	Tapentadol HCl	150000	Not detected
	Tramadol	150000	Not detected
FTY1	Acetyl fentanyl	16	6.25%
	Acrylfentanyl	1	100.00%
	ω -1-Hydroxyfentanyl	20,000	0.005%
	Isobutyryl fentanyl	1	100.00%
	Ocfentanil	2.3	43.48%
	Butyryl fentanyl	2	50.00%
	Furanyl fentanyl	1	100.00%
	Valeryl fentanyl	2.5	40.00%
	(\pm) β -hydroxythiofentanyl	2.5	40.00%
	4-Fluoro-isobutyrylfentanyl	3	33.33%

	Para-fluorobutyryl fentanyl	4	25.00%
	Para-fluoro fentanyl	2.5	40.00%
	(+)-3-cis-methyl fentanyl	50	2.00%
	Carfentanil	2	50.00%
	Sufentanil	15	6.67%
	Alfentanil	7500	0.01%
	Despropionyl fentanyl (4-ANPP)	2,000	0.05%
	Remifentanil	150000	Not detected
	Norfentanyl	150000	Not detected
	Acetyl norfentanyl	150000	Not detected
	Norcarfentanil	150000	Not detected
	Trazodone	25000	0.004%
NFTY5	Fentanyl	10	50%
	Acetyl fentanyl	150	3.3%
	Acetyl Norfentanyl	200	2.5%
	(±)-β-Hydroxythiofentanyl HCl	2500	0.2%
	Acryl Fentanyl	2500	0.2%
	Butyryl Fentanyl	5000	0.1%
	Furanyl Fentanyl	10000	0.05%
	Para-fluoro butyryl Fentanyl (P-FBF)	80000	0.006%
	Para-fluoro Fentanyl	40000	0.013%
	9-HydroxyRisperidone	10000	0.05%
	Alfentanil	20000	0.025%
	Isobutyryl Fentanyl	5000	0.1%
	Remifentanil	15000	0.03%
	Valeryl Fentanyl	20000	0.025%
	Thienyl Fentanyl	50	10%
	(+)-3-cis-methyl fentanyl	50	10%
	4-Fluoro-isobutyryl Fentanyl	30000	Not detected
	Despropionyl fentanyl (4-ANPP)	30000	Not detected
	MT-45 diHCL	150000	Not detected
	Ocfentanil	150000	Not detected
	Risperidone	150000	Not detected
	Sufentanil	150000	Not detected
	Carfentanil	150000	Not detected
	Labetalol Hydrochloride	150000	Not detected
	Trazodone	150000	Not detected
	U-47700	150000	Not detected
	ω-1-Hydroxyfentanyl	30000	Not detected
	6-Acetyl morphine	150000	Not detected

	(±)-Amphetamine	150000	Not detected
	Buprenorphine	150000	Not detected
	Buprenorphine-3β-D-glucuronide	150000	Not detected
	Codeine	150000	Not detected
	Dextromethorphan	150000	Not detected
	Dihydrocodeine	150000	Not detected
	EDDP	150000	Not detected
	EMDP	150000	Not detected
	Fluoxetine	150000	Not detected
	Heroin	150000	Not detected
	Hydrocodone	150000	Not detected
	Hydromorphone	150000	Not detected
	Ketamine	150000	Not detected
	Levorphanol tartrate	150000	Not detected
	Meperidine	150000	Not detected
	(±)-Methadone	150000	Not detected
	Morphine	150000	Not detected
	Morphine-3-β-D-glucuronide	150000	Not detected
	Naloxone hydrochloride	150000	Not detected
	Naltrexone hydrochloride	150000	Not detected
	Norbuprenorphine	150000	Not detected
	Norcodeine	150000	Not detected
	Norketamine	150000	Not detected
	Normeperidine	150000	Not detected
	Normorphine	150000	Not detected
	Noroxycodone	150000	Not detected
	Oxycodone	150000	Not detected
	Oxymorphone	150000	Not detected
	Pentazocine (Talwin)	150000	Not detected
	Pipamperone	150000	Not detected
	Tapentadol hydrochloride	150000	Not detected
	Thioridazine	150000	Not detected
	Tilidine	150000	Not detected
	Tramadol	150000	Not detected
	O-Desmethyl -cris-Tramadol	150000	Not detected
	N-Desmethyl -cris-Tramadol	150000	Not detected
	Norcarfentanil	150000	Not detected
THC20	(-)-11-Nor-Δ ⁹ -THC- 9-carboxylic acid glucuronide	20	100%
	(±)-11-Hydroxy-Δ ⁹ -THC	30	66.7%
	(-)-11-nor-9-carboxy-Δ ⁹ -THC	20	100%

	(-)- Δ^9 -THC	6000	0.3%
	(-)- Δ^8 -THC	4000	0.5%
	Cannabinol	8000	0.25%
	Cannabidiol	150000	Not detected
TRA100	n-Desmethyl -cris-Tramadol	400	25%
	o-Desmethyl -cris-Tramadol	1000	10%
	o-Desmethyl Venlafaxine	15000	Not detected
	Venlafaxine HCl	150000	Not detected

Interference:

Potential interference from compounds chemically dissimilar to the target drugs (6-Monoacetylmorphine (6-MAM10), Fentanyl (FTY1), Norfentanyl (NFTY5), Cannabinoids (THC20), and Tramadol (TRA100)) and from endogenous agents was performed by spiking the substances into pooled urine containing target drugs at near-cutoff concentrations (at +50% and -50% of cutoff). Unless otherwise indicated, substances were tested for potential interference at concentrations of 100 mg/mL. No interference was detected for following substances.

Acetaminophen	Effexor	Nimodipine
Acetophenetidin	Enalapril Maleate	Nitroglycerin
Acetylsalicylic Acid	Erythromycin	Norethindrone
Acyclovir	Esomeprazole Magnesium	N-Acetylprocainamide
Afrin	β -Estradiol	O-Hydroxyhippuric Acid
Albumin (100mg/dL)	1% ethanol	Olanzapine
Aminophylline	Fenofibrate	Omeprazole
Aminopyrine	Fenoprofen	Oxalic Acid
Amiodarone Hydrochloride	Fentanyl Citrate	Oxolinic Acid
Amlodipine Mesylate	Fluoxetine Hydrochloride	Oxymetazoline
Amoxicillin	Fluvoxamine	Ondansetran
Ampicillin	Furosemide	Paliperidone
Apomorphine	Gabapentin	Pantoprazole
Aripiprazole	Gentisic Acid	Papaverine
Aspartame	Glibenclamide	Paroxetine Hydrochloride
Atomoxetine	Gliclazide	Penfluridol
Atorvastatin Calcium	Glipizide	Penicillin V Potassium
Atropine	Glucose	Penicillin-G
Benzilic Acid	Haloperidol	Phenelzine
Benzoic Acid	Hemoglobin	Pioglitazone Hydrochloride
Bilirubin	Hydrochlorothiazide	Piracetam
Bupropion	Hydrocortisone	Pravastatin Sodium
Captopril	3-Hydroxytyramine	Prednisone
Carbamazepine	Isosorbide Dinitrate	Propylthiouracil
Cefradine	Isoxsuprine	Quetiapine Fumarate
Cephalexin	Ibuprofen	Quinine
Chloral Hydrate	Ketoconazole	Ranitidine
Chloramphenicol	Ketoprofen	Rifampicin

Chlorothiazide	Ketamine	Risperidone
Cholesterol	Kratom powder	Salicylic Acid
Ciprofloxacin Hydrochloride	Labetalol	Serotonin
Citalopram	Lamotrigine	Sertraline Hydrochloride
Clarithromycin	Levofloxacin Hydrochloride	Sildenafil Citrate
Clonidine	Levonorgestrel	Simvastatin
Clopidogrel Hydrogen Sulphate	Levothyroxine Sodium	Sodium Valproate
Clozapine	Lidocaine Hydrochloride	Spirolactone
Conjugated Estrogens	Lisinopril	Sulfamethazine
Cortisone	Lithium Carbonate	Sulindac
Creatinine	Liverite	Tetracycline
(-) Cotinine	Loperamide	Tetrahydrocortisone 3-acetate
chlorpheniramine	Loratadine	Tetrahydrocortisone 3- (β -D glucuronide)
D,L-Octopamine	Magnesium	Tetrahydrozoline
D,L-Propranolol	Meperidine	Thiamine
D,L-Tyrosine	Meprobamate	Thioridazine
Deoxycorticosterone	Metoprolol Tartrate	Topiramate
Dextromethorphan	Mifepristone	Tramadol Hydrochloride
Diclofenac	Mirtazapine	Trazodone Hydrochloride
Diffunisal	Montelukast Sodium	Triamterene
Digoxin	Mosapride Citrate	Trifluoperazine
Diphenhydramine	Minocycline	Trimethoprim
Dirithromycin	Nalidixic Acid	Uric Acid
Domperidone	Naproxen	Valproate
D-Pseudoephedrine	Niacinamide	Verapamil
Duloxetine	Nifedipine	Vitamin B2
Dicyclomine	Nikethamide	Vitamin C
Chloroquine	Ecgonine Methyl Ester	Promethazine

Interference by pH and specific gravity were also evaluated using pooled urine specimens containing target drugs at near-cutoff concentrations (at +50% and -50% of cutoff). The results demonstrated that pH levels of 4 to 9 and specific gravity levels of 1.000 to 1.035 do not affect the results of the assays.

4. Detection Limit and Assay Reportable Range:

Not applicable. This device is intended for qualitative use only.

5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):

The stability results support a shelf life of 24 months at 2-30°C based on the real-time stability study data.

6. Assay Cut-Off:

Characterization of how the device performs analytically around the claimed cutoff concentration appears in the precision section, A.1., above.

B Comparison Studies:

1. Method Comparison with Predicate Device:

Method comparison results of AMP500/1000, BAR300, BUP10, BZO300, COC150/300, EDDP300, MDMA500, MET500/MET1000, MOP300, MTD300, OPI2000, OXY100, PCP25, PPX300, TCA1000 and THC50 were previously reviewed under K202567.

The method comparison study for 6-Monoacetylmorphine (6-MAM10), Fentanyl (FTY1), Cannabinoids (THC20), Tramadol (TRA100) and Norfentanyl (NFTY5), was performed in-house by three operators with 80 unaltered urine samples. These samples were blind labeled and compared to LC/MS results. The results are shown in the table below.

Drug test	Test Cup Result		Drug-Free	Low Negative by LC-MS/MS (less than - 50%)	Near Cutoff Negative by LC-MS/MS (Between - 50% and the Cutoff)	Near Cutoff Positive by LC-MS/MS (Between the cutoff and +50%)	High Positive by LC-MS/MS (greater than +50%)
6AM10	Operator A	+	0	0	2	21	18
		-	13	14	11	1	0
	Operator B	+	0	0	2	21	18
		-	13	14	11	1	0
	Operator C	+	0	0	3	20	18
		-	13	14	10	2	0
FTY1	Operator A	+	0	0	2	20	18
		-	12	16	10	2	0
	Operator B	+	0	0	1	21	18
		-	12	16	11	1	0
	Operator C	+	0	0	3	20	18
		-	12	16	9	2	0
NFTY5	Operator A	+	0	0	3	22	16
		-	10	16	11	2	0
	Operator B	+	0	0	2	23	16
		-	10	16	12	1	0
	Operator C	+	0	0	2	22	16
		-	10	16	12	2	0
THC20	Operator A	+	0	0	2	29	10
		-	9	15	14	1	0

	Operator B	+	0	0	2	29	10
		-	9	15	14	1	0
	Operator C	+	0	0	2	29	10
		-	9	15	14	1	0
TRA100	Operator A	+	0	0	2	27	12
		-	10	18	10	1	0
	Operator B	+	0	0	1	27	12
		-	10	18	11	1	0
	Operator C	+	0	0	2	27	12
		-	10	18	10	1	0

Discordant results:

Drug	Operator	Sample Number	LC/MS Results (ng/mL)	Discordant Device Results
6-MAM 10	Operator A	SU25050057	8.193	Positive
	Operator A	SU25050075	9.192	Positive
	Operator A	SU25050041	10.825	Negative
	Operator B	SU25050045	7.867	Positive
	Operator B	SU25050057	8.193	Positive
	Operator B	SU25050074	10.359	Negative
	Operator C	SU25050043	7.861	Positive
	Operator C	SU25050045	7.867	Positive
	Operator C	SU25050075	9.192	Positive
	Operator C	SU25050029	10.863	Negative
	Operator C	SU25050015	10.942	Negative
FTY 1	Operator A	SU25050271	0.848	Positive
	Operator A	SU25050284	0.965	Positive
	Operator A	SU25050249	1.136	Negative
	Operator A	SU25050246	1.218	Negative
	Operator B	SU25050258	0.77	Positive
	Operator B	SU25050249	1.136	Negative
	Operator C	SU25050258	0.77	Positive
	Operator C	SU25050271	0.848	Positive
	Operator C	SU25050284	0.965	Positive
	Operator C	SU25050270	1.184	Negative
	Operator C	SU25050246	1.218	Negative
NFTY 5	Operator A	SU25060003	4.17	Positive
	Operator A	SU25060030	4.765	Positive
	Operator A	SU25060011	4.879	Positive
	Operator A	SU25060065	5.341	Negative
	Operator A	SU25060073	5.699	Negative
	Operator B	SU25060042	4.315	Positive
	Operator B	SU25060030	4.765	Positive
	Operator B	SU25060075	5.526	Negative
	Operator C	SU25060042	4.315	Positive

	Operator C	SU25060011	4.879	Positive
	Operator C	SU25060075	5.526	Negative
	Operator C	SU25060073	5.699	Negative
THC 20	Operator A	SU25050194	18.614	Positive
	Operator A	SU25050181	19.363	Positive
	Operator A	SU25050220	20.645	Negative
	Operator B	SU25050185	18.422	Positive
	Operator B	SU25050194	18.614	Positive
	Operator B	SU25050167	20.499	Negative
	Operator C	SU25050185	18.422	Positive
	Operator C	SU25050181	19.363	Positive
	Operator C	SU25050167	20.499	Negative
TRA 100	Operator A	SU25050130	98.136	Positive
	Operator A	SU25050098	98.772	Positive
	Operator A	SU25050144	108.094	Negative
	Operator B	SU25050130	98.136	Positive
	Operator B	SU25050120	102.471	Negative
	Operator C	SU25050135	83.391	Positive
	Operator C	SU25050098	98.772	Positive
	Operator C	SU25050120	102.471	Negative

2. Matrix Comparison:

Not applicable.

C Clinical Studies:

1. Clinical Sensitivity:

2. Clinical Specificity:

3. Clinical Cut-Off:

4. Other Clinical Supportive Data (When 1. and 2. Are Not Applicable):

A lay user study was performed involving a total of 140 participants from 3 sites. 68 males and 72 females tested one configuration of SAFELife™ T-Dip Multi-Drug Urine Test Panel, including 6-MAM10, THC20, TRA100, NFTY5, AMP1000, BAR300, BUP10, BZO300, COC300, FTY1, mAMP1000, MDMA500, MTD300, OPI2000, OXY100. Each participant was provided one package insert, one blind labeled test solution, and one test device. Test solutions were randomly assigned to participants, one for each. Following testing, users completed a study questionnaire to assess usability and user comprehension, and the results from this questionnaire were found to be acceptable. Participants aged 18 and over, with diverse educational backgrounds, were recruited. Results from the lay user testing are provided in the table below. Lay user results of AMP500, COC150, EDDP300, MET500,

MOP300, PCP25, PPX300, TCA1000, and THC50 were previously reviewed under K202567.

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
6-MAM	10	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
AMP	1000	Negative	20	20	20	20	1	0	0
		Positive	0	0	0	0	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	100%	95%	100%	100%
BAR	300	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
BUP	10	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
BZO	300	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
COC	150	Negative	20	20	20	19	0	0	0
		Positive	0	0	0	1	20	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	100%	100%	100%
FTY	1	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
mAMP	1000	Negative	20	20	20	20	1	0	0
		Positive	0	0	0	0	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	100%	95%	100%	100%
MDMA	500	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%
MTD	300	Negative	20	20	20	19	2	0	0

		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
NFTY	5	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
OPI	2000	Negative	20	20	20	20	1	0	0
		Positive	0	0	0	0	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	100%	95%	100%	100%
OXY	100	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
THC	20	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
TRA	100	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%

D Expected Values/Reference Range:

Not applicable.

VIII Proposed Labeling:

The labeling support the finding of substantial equivalence for this device.

IX Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.