



April 8, 2026

Shenzhen TPH Technology Co., Ltd.  
Dale Wang  
Regulatory Affairs Engineer  
Rm. 203, 2nd Floor, 29th Bldg., Lianchuang Technology Park  
# 21 Bulan Rd., Xialilang Community, Nanwan St. Longgang Dis  
Shenzhen, 518000  
CHINA

Re: K260033  
Trade/Device Name: Wearable Breast Pump (Model S21A)  
Regulation Number: 21 CFR 884.5160  
Regulation Name: Powered Breast Pump  
Regulatory Class: II  
Product Code: HGX  
Dated: November 22, 2025  
Received: January 6, 2026

Dear Dale Wang:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: The Center for Devices and Radiological Health (CDRH) does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, the Food and Drug Administration (FDA) may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality Management System Regulation (QMSR) (21 CFR Part 820), which includes, but is not limited to, ISO 13485 clause 7.3 (Design controls), ISO 13484 clause 8.3 (Nonconforming product), and ISO 13485 clause 8.5 (Corrective and preventative action). Please note that regardless of whether a change requires premarket review, the QMSR requires device manufacturers to review and approve changes to device design and production (ISO 13485 clause 7.3 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the Quality Management System Regulation (QMSR) (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the

Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

**Monica D. Garcia -S**

Monica D. Garcia, Ph.D.  
Assistant Director  
DHT3B: Division of Reproductive,  
Gynecology, and Urology Devices  
OHT3: Office of Gastrorenal, ObGyn,  
General Hospital, and Urology Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K260033

Device Name  
Wearable Breast Pump (Model S21A)

### Indications for Use (Describe)

The Wearable Breast Pump (Model S21A) is a powered breast pump intended to be used by lactating women to express and collect milk from their breasts. It is intended for a single user.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary – K260033

### 1. Submitter Information

Applicant: Shenzhen TPH Technology Co., Ltd.  
Address: Room 203, 2<sup>nd</sup> floor, 29<sup>th</sup> Building, Lianchuang Technology Park, No.21 Bulan Road,  
Xialilang Community, Nanwan Street, Longgang District, Shenzhen, China 518000  
Tel.: 86-19129492601

### 2. Correspondent Information

Contact: Dale WANG  
Title: Regulatory Affairs Engineer  
Company: Shenzhen TPH Technology Co., Ltd.  
Email: dale@tph-tech.com

3. Date prepared: April 7, 2026

### 4. Device Information

Device Name: Wearable Breast Pump (Model S21A)  
Common Name: Powered Breast Pump  
Regulation Number: 21 CFR 884.5160  
Regulation Name: Powered Breast Pump  
Product Code: HGX (Pump, Breast, Powered)  
Regulatory Class: Class II

### 5. Predicate Device Information

Device Name: Wearable Breast Pump (Model S21)  
510(k) Number: K250368  
Manufacturer: Shenzhen TPH Technology Co., Ltd.  
The predicate device has not been subject to a design-related recall.

### 6. Device Description

The Wearable Breast Pump (Model S21A) is a powered breast pump intended to be used by lactating women to

express and collect milk from their breasts, it is intended for a single user. The Wearable Breast Pump (Model S21A) is a breast pump powered by lithium battery, utilizing an embedded control program to manage all device functions. The main components of this pump include: pump, valve, control board, and milk collector. The user interface allows the user to switch from stimulation, expression, and auto modes and control the vacuum levels within those modes.

All available modes consist of 12 vacuum levels. The Wearable Breast Pump (Model S21A) is capable of providing vacuum levels from 40-150 mmHg with cycling rates from 67-109 cycles per minute in stimulation mode, vacuum levels from 100-245 mmHg with cycling rates from 25-93 cycles per minute in expression mode, and vacuum levels from 40-245 mmHg with cycling rates from 25-120 cycles per minute in auto mode. The Wearable Breast Pump (Model S21A) is charged with a 5 V DC adaptor and powered by an internal rechargeable lithium-ion polymer battery. The motor unit operates on embedded software. Software updates by end-users are not supported. The subject device is for repeated use by a single user in a home environment. The device is provided not sterile.

The motor unit operates on a rechargeable battery (3.7V, 1000mAh) and does not function when charging. The rechargeable battery can be charged from the external USB adapter if the motor unit is not in operation.

The breast pump expresses milk by creating a seal around the nipple using the flange and applying and releasing suction to the nipple. The milk is collected in a milk collection container, which can be used for storage. To prevent milk from flowing into the vacuum system, a backflow protection membrane physically separates the milk-contacting pathway from the vacuum system.

All other components (i.e., motor unit/housing) of the subject device are not in contact with the breast. All milk contacting components are compliant with 21 CFR 177.

## **7. Indications for Use**

The Wearable Breast Pump (Model S21A) is a powered breast pump intended to be used by lactating women to express and collect milk from their breasts. It is intended for a single user.

## **8. Comparison of Intended Use and Technological Characteristics with the Predicate Device**

The table below compares the intended use and technological characteristics of the subject and predicate device.

**Table 1: Comparator Table for Subject and Predicate Devices**

	<b>Wearable Breast Pump (Model S21A) K260033 Subject Device</b>	<b>Wearable Breast Pump (Model S21) K250368 Predicate Device</b>	<b>Comparison</b>
Product Name	Wearable Breast Pump (Model S21)	Wearable Breast Pump (Model S12)	N/A
Code	HGX	HGX	Same
Indications for use	The Wearable Breast Pump (Model S21A) is a powered breast pump intended to be used by lactating women to express and collect milk from their breasts. It is intended for a single user.	The Wearable Breast Pump (Model S21) is a powered breast pump intended to be used by lactating women to express and collect milk from their breasts. it is intended for a single user.	Same
Single User	Yes	Yes	Same
Single/Double pump	Single	Single	Same
Pump Type	Wearable	Wearable	Same
Media Separation (backflow protection)	Yes	Yes	Same
Cycling Control Mechanism	Microcontroller	Microcontroller	Same
Expression Pattern	2-Phase	2-Phase	Same
Power Supply	Li-Ion Battery	Li-Ion Battery	Same
Suction levels (Stimulation)	40-150 mmHg	40-120 mmHg	Different
Suction levels (Expression)	100-245 mmHg	120-245 mmHg	Different
Suction levels (Auto)	40-245 mmHg	40-245 mmHg	Same

Cycles per minute (Stimulation)	67-109 cpm	62-104 cpm	Different
Cycles per minute (Expression)	25-93 cpm	22-86 cpm	Different
Cycles per minute (Auto)	25-120 cpm	22-104 cpm	Different
Suction levels	12 Levels for expression, stimulation, auto each	12 Levels for expression, stimulation, auto each	Same
User Interface	On-Off switch, mode selection/long press, vacuum adjustment, LED display, night light button	On-Off switch, mode selection/long press, vacuum adjustment, LED display	Similar
Heating element	No	No	Same
Adjustable Suction Levels	Yes	Yes	Same
Mobile Application	No	No	Same
Design	Wearable pump with combined Milk Collector and Flange	Wearable pump with combined Milk Collector and Flange	Same
Materials	Milk collector/Linker: Polypropylene Flange/Valve/Diaphragm: Silicone Pump motor/outer housing: ABS	Milk collector/Linker: Polypropylene Flange/Valve/Diaphragm: Silicone Pump motor/outer housing: ABS	Same

The indications for use of the subject and predicate device are identical, and they have the same intended use (i.e., for collection of breast milk from the breasts of lactating women).

The subject and predicate devices have similar technological features, including wearable operation, power supply, materials, and user interface. However, as shown in the table above, there are technological differences between the subject and predicate device, including different overall vacuum/cycle specifications. The different technological characteristics of the subject device, as compared to the predicate device, do not raise different questions of safety and effectiveness.

## 9. Summary of Non-Clinical Performance Testing

### **Biocompatibility**

Biocompatibility information was leveraged from K250463 and was in accordance with the 2023 FDA guidance document *Use of International Standard ISO 10993-1, "Biological Evaluation of Medical Devices – Part 1: Evaluation and testing within a risk management process."* The following tests were conducted on the leveraged device:

- Cytotoxicity (ISO 10993-5:2009)
- Skin Sensitization (ISO 10993-10:2010)
- Skin Irritation (ISO 10993-10:2010)

The user-contacting materials were shown to be non-cytotoxic, non-irritating, and non-sensitizing.

### **Electrical Safety**

Testing was conducted in accordance with ANSI/AAMI ES60601- 1:2005/A2:2010 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance), IEC 62133-2:2017, Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems, and IEC 60601-1-11:2015 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

### **Electromagnetic Compatibility**

Testing was conducted in accordance with IEC 60601-1-2:2014 Medical Electrical Equipment - Part 1-2: *"General Requirements For Basic Safety And Essential Performance - Collateral Standard: Electromagnetic Compatibility - Requirements And Tests."*

### **Software**

Software was evaluated at the Basic Documentation level as recommended in the 2023 FDA guidance document *"Content of Premarket Submissions for Device Software Functions."*

### **Performance Testing**

Other performance testing was conducted to show that the device meets its design requirements and performs as intended. The performance tests include:

- Vacuum level verification testing at each mode/cycle demonstrated that the devices meet mode/cycle specifications.
- Backflow protection testing was conducted to verify liquid does not backflow into the tubing.
- Use life testing was conducted to demonstrate that the device maintains its specifications throughout its proposed use life.
- Battery performance testing was conducted to demonstrate that the battery remains functional during its stated battery use-life.
- Battery status indicator testing was conducted to demonstrate that the battery status indicator remains functional during its stated battery life.

## **10. Conclusion**

The results of the performance testing described above demonstrate that The Wearable Breast Pump (Model S21A) is as safe and effective as the predicate device and supports a determination of substantial equivalence.