



June 3, 2026

Vista Industries,, LLC
% Andrew Wang
Consultant
Shanghai SUNGO Management Consulting Co., Ltd.
14th Floor, Dongfang Bldg., 1500# Century Ave.
Shanghai, Shanghai 200122
China

Re: K260806

Trade/Device Name: Nitrile Disposable Examination Gloves
Regulation Number: 21 CFR 880.6250
Regulation Name: Non-Powdered Patient Examination Glove
Regulatory Class: Class I, reserved
Product Code: LZA
Dated: March 12, 2026
Received: March 12, 2026

Dear Andrew Wang:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality Management System Regulation (QMSR) (21 CFR Part 820), which includes, but is not limited to, ISO 13485 clause 7.3 (Design controls), ISO 13485 clause 8.3 (Nonconforming product), ISO 13485 clause 8.5.2 (Corrective action), and ISO 13485 clause 8.5.3 (Preventative action). Please note that regardless of whether a change requires premarket review, the QMSR requires device manufacturers to review and approve changes to device design and production (ISO 13485 clause 7.3 and ISO 13485 clause 7.5) and document changes and approvals in the Medical Device File (ISO 13485 clause 4.2.3).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the Quality Management System Regulation (QMSR) (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory->

[assistance/contact-us-division-industry-and-consumer-education-dice](#)) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

BIFENG QIAN -S

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4C: Division of Infection
Control Devices
OHT4: Office of Surgical and
Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K260806

Device Name

Nitrile Disposable Examination Gloves

Indications for Use (Describe)

Nitrile disposable examination gloves are non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

K260806

Document Prepared Date: 2026-05-18

1. Applicant Information

Vista Industries, LLC.

Address: 1560 DRUMHILL DR, HACIENDA HEIGHTS, CA 91745

Contact: James Yu

Email: medical@xingyugloves.com

Submission Correspondent:

Primary Contact: Andrew Wang

Shanghai SUNGO Management Consulting Co., Ltd.

14th Floor, Dongfang Building, 1500# Central Ave., Shanghai 200122, China

Tel: 0086-21-58817802

Email: yuanju.wang@sungoglobal.com

2. Subject Device

Device Name: Nitrile Disposable Examination Gloves

Models: D2000, D3000, D6600, D9000

Sizes: S, M, L, XL

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination

Glove Regulatory Class: Class I

Product Code: LZA

3. Predicate Device

510(k) Number: K240361

Product Name: Powder Free Green Nitrile Examination Gloves (XS/S/M/L/XL); Powder Free Pink Nitrile Examination Gloves (XS/S/M/L/XL); Powder Free Purple Nitrile Examination Gloves (XS/S/M/L/XL)

Applicant/Sponsor: Hebei Titans Hongsen Medical Technology Co., Ltd.

4. Device Description

Nitrile disposable examination gloves are non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. The gloves are pink, yellow, grey or green color, powder free, nitrile ambidextrous gloves. The gloves are offered in S, M, L or XL size, packed in a paper box.

The gloves are designed and manufactured in accordance with the ASTM D6319-19 standard.

5. Indication for Use

Nitrile disposable examination gloves are non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

6. Substantial Equivalent (SE) Comparison

Device		Proposed Device	Predicate Device	Result
510K #		K260803	K240361	
Product Name		Nitrile Disposable Examination Gloves	Powder Free Green/Pink/Purple Nitrile Examination Gloves	
Product Code		LZA	LZA	Same
Regulation Number		21 CFR 880.6250	21 CFR 880.6250	Same
Indications for use		Nitrile disposable examination gloves are non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.	Powder Free Green Nitrile Examination Gloves, Powder Free Pink Nitrile Examination Gloves and Powder Free Purple Nitrile Examination Gloves is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.	Same
Powder free		Yes	Yes	Same
Design feature		Ambidextrous	Ambidextrous	Same
Material		Nitrile	Nitrile	Same
Color		Pink, Yellow, Grey, Green	Green, Pink, Purple	Different 1
OTC use		Yes	Yes	Same
Sterility		Non-sterile	Non-sterile	Same
Use		Single use	Single use	Same
Label		Single-use indication, powder free, device color, device name, glove size and quantity, non-sterile.	Single-use indication, powder free, device color, device name, glove size and quantity, non-sterile.	Same
Size and Dimensions (mm)		Length(mm): S/M/L/XL: 240 ± 10mm Width(mm): S: 85±10mm; M: 96±10mm; L: 107±10mm; XL: 118±10mm	Length(mm): XS/S/M/L/XL: ≥230 Width(mm): XS: 70±10mm; S: 85±10mm; M: 95±10mm; L: 110±10mm; XL: 120±10mm	Different 2
Thickness (mm)		Middle finger: 0.11±0.02mm; Palm: 0.06±0.02mm	Finger: ≥0.05; Palm: ≥0.05	Different 3
Physical Properties	Before Aging	Tensile≥14MPa	Tensile Strength: 14MPa, min	Same
		Elongation≥500%	Ultimate Elongation: 500% min	Same
	After Aging	Tensile≥14MPa	Tensile Strength: 14MPa, min	Same
		Elongation≥400%	Ultimate Elongation: 400% min	Same
Freedom from Holes		Be free from holes when tested in Accordance with ASTM D5151	Be free from holes when tested in Accordance with ASTM D5151	Same

Device		Proposed Device	Predicate Device	Result
Powder Content		Meet the requirements of ASTM D6124 Less than 2mg per glove	Meet the requirements of ASTM D6124 Less than 2mg per glove	Same
Biocompatibility	Irritation ISO 10993-23	Under the conditions of the study, not an irritant. Comply with ISO 10993-23.	Under the conditions of the study, not an irritant. Comply with ISO 10993-23.	Same
	Sensitization ISO 10993-10	Under the conditions of the study, not a sensitizer. Comply with ISO 10993-10.	Under the conditions of the study, not a sensitizer. Comply with ISO 10993-10.	Same
	Cytotoxicity ISO 10993-11	Under the conditions of the study, the device did not show acute systemic toxicity in vivo. Comply with ISO 10993-11.	Under the conditions of the study, the device did not show acute systemic toxicity in vivo. Comply with ISO 10993-11.	Same

Different 1-Color

The color of the proposed device is different from the predicate device. However, the biocompatibility and specification performance tests have been performed on each of the proposed devices and the test results show no significant concerns. Therefore, the difference will not affect the safety and effectiveness of the proposed device.

Different 2-Size and dimensions

The size and dimensions of the proposed device is not exactly same as the predicate device. However, the size and dimension of the proposed device has been covered by ASTM D6319-19. The user can select appropriate model depended on size of user's hand. In addition, its dimension has been tested and met the requirement of ASTM D6319-19. Therefore, the difference will not affect the safety and effectiveness of the proposed device.

Different 3-Thickness

The thickness of the proposed device is not exactly same as the predicate device. However, the thickness of the proposed device has been tested and met the requirement of ASTM D6319-19. Therefore, the difference will not affect the safety and effectiveness of the proposed device.

7. Non Clinical Test

7.1 Product Material Safety

The following tests for the subject device were conducted to evaluated the biocompatibility of Nitrile Disposable Examination Gloves:

- ISO 10993-11: 2017, Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
- ISO 10993-10: 2021, Biological evaluation of medical devices - Part 10: Tests for skin sensitization
- ISO 10993-23: 2021, Biological evaluation of medical devices - Part 23: Tests for irritation

7.2 Performance of the products

The testing was done according to the following standards:

- ASTM D6319-19:2023, Standard Specification for Nitrile Examination Gloves for Medical Application
- ASTM D6124-06:2022, Standard Test Method for Residual Powder on Medical Gloves
- ASTM D5151-19:2023, Standard Test Method for Detection of Holes in Medical Gloves.

Test Method	Purpose	Acceptance Criteria	Results
ASTM D6319-19	Physical Dimensions Length	Min 220mm for size XS, S Min 230mm for size M- XXL Thickness (mm): Finger: ≥ 0.05 ; Palm: ≥ 0.05	Pass
ASTM D6319-19	Physical Dimensions Palm Width	XS: 70 ± 10 mm, S: 80 ± 10 mm, M: 95 ± 10 mm, L: 110 ± 10 mm, XL: 120 ± 10 mm, XXL: 130 ± 10 mm	Pass
ASTM D6319-19	Physical Dimensions Thickness	Finger: 0.05mm (min) Palm: 0.05mm (min)	Pass
ASTM D6319-19	Physical Properties	Tensile Strength (Min 14MPa) and Elongation (Before Aging 500% and after aging 400%) Min	Pass
ASTM D6319-19 ASTM D5151-19	Water leak test	AQL 2.5 (ISO 2859-1)	Pass
ASTM D6319-19 ASTM D6124-06	Powder Residue	Max 2mg/glove	Pass
ISO 10993-11	Acute systemic toxicity study	Under the conditions of the study, the device extract does not induce acute systemic toxicity response.	Pass
ISO 10993-10	Skin sensitization	Under the condition of the study, the device is not a sensitizer	Pass
ISO 10993-23	Irritation	Under the condition of the study, the device is not an irritant	Pass

8. Summary of Clinical Testing

Not Applicable

9. Conclusion

Based on the nonclinical tests performed, the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device, Disposable Medical Examination Nitrile Gloves, cleared under K240361.