

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY AND INSTRUMENT **COMBINATION** TEMPLATE**

A. 510(k) Number:

K173330

B. Purpose for Submission:

To obtain a Substantially Equivalent determination for the PanNAT Shiga Toxin-producing *E. coli* (STEC) Test

C. Measurand:

Conserved regions of the Shiga toxin genes (*stx1*, *stx2*) and the *E. coli* O157 O-antigen gene cluster (*fli*)

D. Type of Test:

Qualitative real-time Polymerase Chain Reaction (PCR)

E. Applicant:

Micronics, Inc.

F. Proprietary and Established Names:

PanNAT STEC Test

G. Regulatory Information:

1. Regulation section:

21 CFR 866.3990: Gastrointestinal microorganism multiplex nucleic acid-based assay

2. Classification:

Class II (Special Controls)

3. Product code:

PCI: Gastrointestinal Bacterial Panel Multiplex Nucleic Acid-Based Assay System

PCH: Gastrointestinal Pathogen Panel Multiplex Nucleic Acid-Based Assay System

OOI: Real Time Nucleic Acid Amplification System

4. Panel:

83 - Microbiology

H. Intended Use:

1. Intended use(s):

The Micronics PanNAT STEC Test is a qualitative, *in vitro* nucleic acid amplification-based test intended for the simultaneous detection and identification of the *stx1* and *stx2* Shiga toxin genes and the O-antigen gene cluster of *E. coli* O157 (*fcl*). Testing is performed in a unitized cartridge on the PanNAT System on soft to diarrheal unpreserved or Cary-Blair preserved stool specimens from individuals with signs and symptoms of gastrointestinal infection. The PanNAT STEC Test is intended for use in conjunction with clinical presentation, laboratory findings and epidemiological risk factors, as an aid in detection of specific Shiga-toxin expressing strains of *E. coli* (“STEC”) from patients with diarrhea.

The results of this test should not be used as the sole basis for diagnosis, treatment or other patient management decisions. Positive PanNAT STEC Test results do not rule out the potential for coinfection with other pathogens that are not detected by this device and may not be indicative of the sole or definitive cause of patient illness. Negative PanNAT STEC Test results in the context of clinical illness compatible with gastroenteritis may be due to infection by pathogens that are not detected by this test or non-infectious causes such as ulcerative colitis, irritable bowel syndrome, or Crohn’s disease.

2. Indication(s) for use:

Same as Intended Use.

3. Special conditions for use statement(s):

For prescription use only.

The PanNAT STEC Test may exhibit reduced sensitivity and/or failure detection failure with *stx1d*, *stx2b*, *stx2f* and *stx2g*. These subtypes of Shiga toxin are infrequently isolated from human specimens.

4. Special instrument requirements:

The PanNAT STEC Test is indicated for use with the PanNAT System.

I. Device Description:

The PanNAT STEC Test is a cartridge-based real-time PCR-based assay for the detection of the Shiga toxin genes *stx1* and *stx2* and the *Escherichia coli* O157-specific gene cluster, *fliC*. The assay is performed on the PanNAT System which comprises an instrument and software for automated sample processing, nucleic acid amplification and detection and result interpretation. The PanNAT STEC Test is for use on preserved (Cary-Blair Transport Medium) or unpreserved soft or diarrheal stool specimens from subjects with signs and symptoms of gastrointestinal infection.

The PanNAT STEC Test kit contains components to process and analyze up to 25 specimens and includes PanNAT Sample Transfer Packs (comprised of a swab, Sample Buffer Tube and Adaptor Cap) and PanNAT STEC Test Cartridges.

To perform the PanNAT STEC Test, a swab is used to transfer part of the stool specimen to a Sample Buffer Tube which is then fitted with an Adaptor Cap and attached to the inlet port of a PanNAT STEC Test Cartridge. The cartridge-tube assembly is then inserted into the PanNAT System for automated processing.

The PanNAT STEC Test cartridge contains all the components necessary to perform the assay and includes separate reaction mixtures for each target analyte (*stx1*, *stx2* and O157) that are dried in different wells. Each reaction mixture also contains primers and a detector probe for an endogenous human gene that is extracted from the sample and co-amplified with the target analyte (if present). There are three reaction mixtures per target analyte. The first three reaction mixtures (one per target) are rehydrated with the processed sample and analyzed to report the presence or absence of the target nucleic acids. The remaining wells are used as Positive and Negative Controls and receive none of the processed sample but instead are rehydrated with buffer. Each of the control wells must produce the expected results for the results from the test wells to be reported.

Detection of the amplification reactions occurs in real-time using quenched fluorescent probes and results are reported automatically.

J. Substantial Equivalence Information:

1. Predicate device name(s):

BD Max Enteric Bacterial Panel

2. Predicate 510(k) number(s):

K140111

3. Comparison with predicate:

Similarities		
Item	Device	Predicate
	PanNAT STEC Test (K173330)	BD MAX Enteric Bacterial Panel (K140111)
Regulation	21 CFR 866.3990	Same
Primary Product Code	PCI	Same
Device Class	II	Same
Intended Use	<p>The Micronics PanNAT STEC Test is a qualitative, <i>in vitro</i> nucleic acid amplification-based test intended for the simultaneous detection and identification of the <i>stx1</i> and <i>stx2</i> Shiga toxin genes and the O-antigen gene cluster of <i>E.coli</i> O157 (<i>fcl</i>). Testing is performed in a unitized cartridge on the PanNAT System on soft to diarrheal unpreserved or Cary-Blair preserved stool specimens from individuals with signs and symptoms of gastrointestinal infection. The PanNAT STEC Test is intended for use in conjunction with clinical presentation, laboratory findings and epidemiological risk factors, as an aid in detection of specific Shiga-toxin expressing strains of <i>E. coli</i> (“STEC”) from patients with diarrhea.</p> <p>The results of this test should not be used as the sole basis for diagnosis, treatment or other patient management decisions.</p>	<p>The BD MAX Enteric Bacterial Panel performed on the BD MAX System is an automated <i>in vitro</i> diagnostic test for the direct qualitative detection and differentiation of enteric bacterial pathogens. The BD MAX Enteric Bacterial Panel detects nucleic acids from:</p> <ul style="list-style-type: none"> • <i>Salmonella</i> spp. • <i>Campylobacter</i> spp. (<i>jejuni</i> and <i>coli</i>) • <i>Shigella</i> spp. / Enteroinvasive <i>E. coli</i> (EIEC) • Shiga toxin 1 (<i>stx1</i>) / Shiga toxin 2 (<i>stx2</i>) genes (found in Shiga toxin-producing <i>E. coli</i> [STEC]) as well as <i>Shigella dysenteriae</i>, which can possess a Shiga toxin gene (<i>stx</i>) that is identical to the <i>stx1</i> gene of STEC. <p>Testing is performed on unpreserved soft to diarrheal stool specimens or Cary-Blair preserved stool specimens from symptomatic patients with suspected acute</p>

Similarities		
Item	Device	Predicate
	PanNAT STEC Test (K173330)	BD MAX Enteric Bacterial Panel (K140111)
	<p>Positive PanNAT STEC Test results do not rule out the potential for coinfection with other pathogens that are not detected by this device and may not be indicative of the sole or definitive cause of patient illness. Negative PanNAT STEC Test results in the context of clinical illness compatible with gastroenteritis may be due to infection by pathogens that are not detected by this test or non-infectious causes such as ulcerative colitis, irritable bowel syndrome, or Crohn's disease.</p>	<p>gastroenteritis, enteritis or colitis. The test is performed directly on the specimen, utilizing real-time polymerase chain reaction (PCR) for the amplification of <i>spaO</i>, a <i>Campylobacter</i> specific <i>tuf</i> gene sequence, <i>ipaH</i> and <i>stx1/stx2</i>. The test utilizes fluorogenic sequence-specific hybridization probes for detection of the amplified DNA.</p> <p>This test is intended for use, in conjunction with clinical presentation, laboratory findings, and epidemiological information, as an aid in the differential diagnosis of <i>Salmonella</i>, <i>Shigella</i>/EIEC, <i>Campylobacter</i> and Shiga toxin-producing <i>E. coli</i> (STEC) infections. Results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions. Positive results do not rule out co-infection with other organisms that are not detected by this test, and may not be the sole or definitive cause of patient illness. Negative results in the setting of clinical illness compatible with gastroenteritis may be due to infection by pathogens that are not</p>

Similarities		
Item	Device	Predicate
	PanNAT STEC Test (K173330)	BD MAX Enteric Bacterial Panel (K140111)
		detected by this test or non-infectious causes such as ulcerative colitis, irritable bowel syndrome, or Crohn's disease.
Analytes	Shiga toxin genes <i>stx1</i> and <i>stx2</i>	Same
Specimen Type	Soft or diarrheal stool: unpreserved or preserved in Cary-Blair medium	Same
Assay Format	Real-time PCR	Same
Result Format	Qualitative	Same

Differences		
Item	Device	Predicate
	PanNAT STEC Test (K173330)	BD MAX Enteric Bacterial Panel (K140111)
Instrument System	PanNAT System	BD MAX
Test Format	Integrated cartridge for sample processing, DNA amplification/detection	Separate consumables for sample processing and amplification/detection
DNA Extraction	Silica membranes	Magnetic affinity beads
Liquid Movement	Pneumatic, microfluidic	Automated pipettor
Detection Chemistry	Hybridization probes	Hydrolysis probes
Process Control	Endogenous human DNA	Exogenous
Analytes	<ul style="list-style-type: none"> • Shiga toxin genes <i>stx1</i> and <i>stx2</i> (differentiated) • <i>E. coli</i> O157 <i>fli</i> gene cluster 	Specific DNA sequences from: <ul style="list-style-type: none"> • Shiga toxin genes <i>stx1</i> and <i>stx2</i> (not differentiated) • <i>Salmonella</i> spp. • <i>Campylobacter</i> spp. (<i>jejuni</i> and <i>coli</i>) • <i>Shigella</i> spp. / Enteroinvasive <i>E. coli</i>

Differences		
Item	Device	Predicate
	PanNAT STEC Test (K173330)	BD MAX Enteric Bacterial Panel (K140111)
		(EIEC)
Positive/Negative Controls	Integrated within test cartridge	Performed in separate test strips/amplification wells
Batch Size	1 per run	Up to 24 per run
Run Time	~1 hour 10 minutes	~3 hours

K. Standard/Guidance Document Referenced (if applicable):

1. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - 3rd Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
2. CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline - 2nd Edition*. CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
3. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - 2nd Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
4. ISO 14971. Medical device – Application of risk management to medical devices.
5. CLSI. *Molecular Diagnostic Methods for Infectious Diseases - 3rd Edition*. CLSI report MM03. Wayne, PA: Clinical and Laboratory Standards Institute; 2015.
6. IEC 61010-2-101 Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 2-101: Particular requirements for *In-Vitro* Diagnostic (IVD) Medical Equipment; 2002.
7. IEC 61010-2-010 Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use – Part 010: Particular requirements for laboratory equipment for the heating of materials; 2014.
8. IEC 61010-2-081 Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes; 2015.
9. Class II Special Controls Guideline: Gastrointestinal Microorganism Multiplex Nucleic Acid-Based Assays for Detection and Identification of Microorganisms and Toxin Genes from Human Stool Specimens; Guideline for Industry and FDA Staff.

L. Test Principle:

The PanNAT STEC Test is performed on soft or diarrheal stool specimens from patients suspected of gastrointestinal infection. The specimens may be preserved in Cary-Blair transport medium or unpreserved. To perform the test, a swab is used to transfer an aliquot of the stool specimen to a Sample Buffer Tube which is sealed with an Adaptor Cap and

inserted into the sample port of the PanNAT STEC Test cartridge. The cartridge is then loaded into the PanNAT System for automated processing.

Each test cartridge contains all the reagents necessary for nucleic acid extraction, amplification and detection for one sample. Processing of the sample within the cartridge is fully automated and fluid movement occurs via pneumatic pressure. Following bacterial lysis using chaotropic salts and detergent, DNA is captured on silica membranes and the eluate is used to rehydrate the amplification reagents.

Each cartridge contains nine amplification wells, three per target gene, that contain dried PCR reagents. Each well includes a pair of primers and a detector probe for an endogenous human gene that acts as a Process Control for nucleic acid recovery and amplification. One chamber corresponding to each target analyte (*stx1*, *stx2* and O157) is rehydrated with the extracted sample. The remaining chambers are rehydrated with buffer and serve as controls. The three Positive Control wells contain recombinant plasmid DNA bearing the target region for each of the target analytes as well as the Process Control. The three Negative Control wells contain no target or Process Control DNA. The expected results must be obtained with each Positive and Negative Control well for patient results to be reported. In addition, to report a negative result for a target analyte, the Process Control must be detected in the corresponding well, otherwise the result is reported as “Invalid.” The results are interpreted automatically by the instrument and may be viewed on-screen, printed or exported via USB (Universal Serial Bus) flash drive. A summary of the algorithm for result interpretation is provided in [Table 1](#).

Table 1. Result interpretation for the PanNAT STEC Test

Positive/Negative Control Wells	Sample Well		Interpretation ²
	Target Analyte ¹	Process Control	
Valid ³	Negative	Negative	Invalid Assay
	Negative	Positive	Valid
	Positive	Negative	Valid
	Positive	Positive	Valid
Invalid ⁴	Negative	Negative	Invalid Assay
	Negative	Positive	Invalid Assay
	Positive	Negative	Invalid Assay
	Positive	Positive	Invalid Assay

¹ *stx1*, *stx2* or O157

² For Positive/Negative results for the target analytes to be reported, *each* of the Sample Wells (*stx1*, *stx2* and O157) must produce a Valid result

³ Valid means that all 3 Positive and all 3 Negative Control wells generated the expected results

⁴ Invalid means that one or more Positive or Negative Control wells generated an unexpected result

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. *Precision/Reproducibility:*

The reproducibility of the PanNAT STEC Test between sites was evaluated in a study that was performed by two operators at each of three sites using a total of 13 PanNAT Systems. Over multiple days, each operator tested 15 replicates of a blinded 9-member panel of unpreserved stool samples containing different levels of the PanNAT STEC Test target analytes (3 sites X 2 operators X 9 panel members X 15 replicates = 810 PanNAT STEC Test results; 90 per panel member). The panel members used in the study are summarized in [Table 2](#). A different lot of reagents was used at each of the testing sites. A summary of the qualitative test results obtained with each panel member, stratified by site/reagent lot is shown in [Table 3](#) and an analysis of cycle numbers for each of the target analytes and the Process Control is provided in [Table 4](#).

Table 2. Panel members evaluated in the PanNAT STEC Test Reproducibility Study

Panel Member	<i>E. coli</i> Strain	Analyte	Target Level	CFU/mL ¹
1	None	None	Negative	0
2	TW14003	<i>stx1</i>	Low	1 x 10 ⁶
3			Medium	3 x 10 ⁶
4	TW08023	<i>stx2</i>	Low	1 x 10 ⁶
5			Medium	3 x 10 ⁶
6	CDC B6914-MS1 (ATCC 43888)	O157	Low	1 x 10 ⁶
7			Medium	3 x 10 ⁶
8	EDL933 (ATCC 43895)	<i>stx1/stx2/O157</i>	Low	1 x 10 ⁶
9			Medium	3 x 10 ⁶

¹ 1 x 10⁶ CFU/mL = 1X LoD; 3 x 10⁶ CFU/mL = 3X LoD ([Section M\(1\)\(d\)](#))

There was a total of 33/810 (4.1%) Invalid results during the study. In addition, run errors occurred with 4/810 (0.5%) samples. Therefore, a total of 773 valid PanNAT STEC Test results was included in the analysis of device reproducibility ([Tables 3](#) and [4](#)).

For each panel member, agreement with expected results was assessed collectively for all three analytes. For panel members #1-7 and #9, agreement with the expected results for all three analytes was ≥97.7%. For panel member #8, the Low Positive sample for the *stx1*⁺/*stx2*⁺/O157⁺ *E. coli* O157 strain EDL933, the agreement at Sites A and B was 93.3% (28/30) and 96.4% (27/28), respectively, whereas at Site C it was 76.7% (23/30). The seven (7) discordant results at Site C were investigated. Of these, two (2) were due to failure to detect *stx1* (*stx1* agreement 28/30, 93.3%), four (4) were due to failure to detect *stx2* (*stx2* agreement 26/30, 86.6%) and one (1) was due to failure to detect O157 (O157 agreement 29/30, 96.7%). In each case, the sample

was correctly reported positive for the other two analytes and the observed discordance in the composite results (for *stx1*, *stx2* and O157 combined) was therefore not due to a systematic failure for any single analyte. Overall, the reproducibility of the PanNAT STEC Test was determined to be acceptable.

Table 3. Summary of qualitative test results from the PanNAT STEC Test Reproducibility Study, stratified by site/reagent lot

Panel Member	Analyte(s)	Level	Agreement (%) ^{1,2}			
			Site A	Site B	Site C	Overall
1	None	Negative	28/28 (100)	29/29 (100)	26/26 (100)	84/84 (100)
2	<i>stx1</i>	Low	28/28 (100)	29/29 (100)	26/28 (92.9)	83/85 (97.6)
3		Medium	27/27 (100)	29/30 (96.7)	25/25 (100)	81/82 (98.8)
4	<i>stx2</i>	Low	27/27 (100)	29/29 (100)	28/30 (93.3)	84/86 (97.7)
5		Medium	28/28 (100)	29/29 (100)	30/30 (100)	87/87 (100)
6	O157	Low	30/30 (100)	29/29 (100)	28/29 (96.6)	87/88 (98.9)
7		Medium	30/30 (100)	28/28 (100)	29/29 (100)	87/87 (100)
8	<i>stx1/stx2/O157</i>	Low	28/30 (93.3)	27/28 (96.4)	23/30 ³ (76.7)	78/88 (88.6)
9		Medium	26/27 (96.3)	30/30 (100)	30/30 (100)	86/87 (98.9)

PanNAT Systems used: Site A: 9; Site B: 3; Site C: 6; 5 of the PanNAT Systems were used at more than one site (3 at Sites A and C; 2 at Sites A and B)

¹ In comparison to the expected results for the panel member for all three analytes combined

² Each site used a different reagent lot: site and reagent lot are therefore confounded

³ Of the 7 discordant results, 2 were due to failure to detect *stx1*, 4 were due to failure to detect *stx2* and 1 was due to failure to detect O157. In each case, the sample was reported positive for the other two analytes.

Table 4. Summary of mean cycle numbers for panel members in the PanNAT STEC Test Reproducibility Study

Panel Member	Analyte(s)	Level	Mean Cycle Number (Standard Deviation)					
			<i>stx1</i>	PrC	<i>stx2</i>	PrC	O157	PrC
1	None	Negative		28.9 (0.9)		28.8 (1.1)		28.5 (0.9)
2	<i>stx1</i>	Low	35.6 (0.7)			29.1 (1.0)		29.0 (0.9)
3		Medium	34.5 (1.0)			29.0 (0.8)		28.7 (0.9)
4	<i>stx2</i>	Low		29.1 (0.9)	35.3 (0.9)			28.7 (1.0)
5		Medium		29.1 (0.8)	33.7 (0.9)			28.9 (0.8)
6	O157	Low		29.1 (0.9)		29.0 (0.9)	35.5 (0.9)	
7		Medium		29.1 (1.1)		29.0 (1.0)	34.2 (0.7)	
8	<i>stx1/stx2/O157</i>	Low	36.1 (0.9)		35.3 (1.0)		35.4 (0.8)	
9		Medium	34.6 (0.8)		33.6 (1.0)		33.8 (0.9)	

PrC: Process Control

Note: Specimens with false negative results were excluded from the calculation of mean and standard deviation

The PanNAT STEC Test demonstrated acceptable reproducibility across sites, PanNAT Systems, operators, panel members and reagent lots.

b. Linearity/assay reportable range:

Not applicable.

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Specimen Stability

The stability of preserved and unpreserved stool specimens for use with the PanNAT STEC Test was evaluated analytically by testing contrived samples that were seeded with an enumerated suspension of *E. coli* O157 strain EDL933 (*stx1*⁺/*stx2*⁺/*fli*⁺) at levels close the LoD of the assay. The samples were stored at 2-8°C and tested at predefined intervals. Five assay replicates were tested at each time point. The results of the study support the claimed stability of preserved and unpreserved stool specimens at 2-8°C for up to 7 days.

The results of additional testing that was performed to establish the stability of frozen, unpreserved stool specimens are described in [Section M\(2\)\(b\)](#).

Process Control

The PanNAT STEC Test Process Control is comprised of an endogenous human gene that is co-extracted from the stool specimen and co-amplified along with the with the

target analytes (if present). The Process Control is designed to monitor for the adequacy of sample addition to the test cartridge as well as reagent and process integrity.

Integrated Positive/Negative Controls

Each PanNAT STEC Test cartridge includes integrated amplification wells for Positive and Negative Controls. These wells contain the primers and probes for detection of *stx1*, *stx2* and O157 target genes as well as the Process Control. In addition, the Positive Control wells contain plasmid DNA that carries the target sequences for each of the PanNAT STEC Test analytes; whereas the Negative Control wells contain no target DNA for any of the target analytes. Each of the Positive and Negative Control wells must produce the expected result in order for sample results to be reported. If one or more of the control wells fails the test sample is reported as “Invalid.” Refer to [Table 1](#), above, for additional information regarding the interpretation of PanNAT STEC Test results.

External Positive/Negative Controls

External Controls should be tested according to guidelines or requirements of local, provincial and/or federal regulations or accreditation organizations.

Commercially prepared External Positive and Negative Controls were tested throughout the prospective Clinical Study described in [Section M\(3\)\(a\)](#). The Positive External Control comprised a standardized suspension of *E. coli* ATCC 43895 strain EDL933 (5×10^6 CFU/mL) and a human cell line ($4.8\text{--}9.6 \times 10^5$ cells/mL). The negative External Control comprised diluent containing the human cell line alone. A total of 170 pairs of External Controls were tested over the course of the study. On initial testing, 6/170 (3.5%) Positive Controls and 12/170 (7.1%) Negative Controls produced Invalid results. Of the remainder, 162/164 (98.7%) Positive Controls and 158/158 (100%) Negative Controls produced the expected results.

As an alternative to commercially prepared controls, a strain of *E. coli* that is known to harbor the genes targeted by the PanNAT STEC Test (e.g., ATCC 43895 strain EDL933) may be used as a Positive Control. Another strain of *E. coli* that does not carry any of the target genes (e.g., *E. coli* K12) may be used as a Negative Control. Clinical specimens that are known to be positive or negative for the PanNAT STEC Test analytes may also be used as Positive and Negative External Controls, respectively.

d. Detection limit:

Limit of Detection

The Limit of Detection (LoD) of the PanNAT STEC Test was initially estimated using seven Shiga toxin positive/negative strains of O157 and non-O157 *E. coli* that were diluted preserved and unpreserved stool matrix at three different concentrations. The estimated LoD for each strain was the lowest concentration at which all of three replicates produced positive results. For two strains, the LoD was then confirmed by

testing additional replicates at the estimated LoD concentration. The LoD was defined as the concentration at which $\geq 95\%$ of replicates produced a positive result for both strains. The results of the study are summarized in **Tables 5** and **6** for unpreserved and preserved stool specimens, respectively. The confirmed LoD with unpreserved stool specimens was 1×10^6 CFU/mL and with preserved specimens in Cary-Blair medium it was 5×10^5 CFU/mL (equivalent to 1.75×10^6 CFU/mL of raw, undiluted stool).

Table 5. Summary of results from the PanNAT STEC Test LoD Study with unpreserved stool specimens

<i>E. coli</i> Strain	Genotype			Number Positive/Tested			
				LoD Titration			LOD Confirmation
	<i>stx1</i>	<i>stx2</i>	O157	1×10^5 CFU/mL	1×10^6 CFU/mL	1×10^7 CFU/mL	1×10^6 CFU/mL
EDL933	+	+	+	2/3 (66.7%)	3/3 (100%)	3/3 (100%)	27/27 ¹ (100%)
92-3265	+	-	+	0/3 (0%)	2/2 ¹ (100%)	3/3 (100%)	
B6914-MS1	-	-	+	3/3 (100%)	3/3 (100%)	3/3 (100%)	
92-3099	-	+	+	1/3 (33.3%)	3/3 (100%)	3/3 (100%)	
MT#2	-	+	-	1/3 (33.3%)	3/3 (100%)	3/3 (100%)	
M105-14	+	-	-	2/3 (66.7%)	3/3 (100%)	3/3 (100%)	
3007-85	+	+	-	1/3 (33.3%)	3/3 (100%)	2/2 ¹ (100%)	28/28 (100%)

¹ 1/28 assay replicates was reported as Invalid and no retest was performed

Table 6. Summary of results from the PanNAT STEC Test LoD Study with stool specimens preserved in Cary-Blair Transport Medium ¹

<i>E. coli</i> Strain	Genotype			Number Positive/Tested				
				LoD Titration			LoD Confirmation	
	<i>stx1</i>	<i>stx2</i>	O157	2.9 x 10 ⁴ CFU/mL	2.9 x 10 ⁵ CFU/mL	2.9 x 10 ⁶ CFU/mL	2.9 x 10 ⁵ CFU/mL	5 x 10 ⁵ CFU/mL
EDL933	+	+	+	0/3 (0%)	3/3 (100%)	2/2 ² (100%)	18/19 ^{3,5} (94.7%)	28/28 (100%)
92-3265	+	-	+	0/3 (0%)	3/3 (100%)	3/3 (100%)		
B6914-MS1	-	-	+	1/3 (33.3%)	3/3 (100%)	3/3 (100%)		
92-3099	-	+	+	0/3 (0%)	3/3 (100%)	3/3 (100%)		
MT#2	-	+	-	0/3 (0%)	3/3 (100%)	3/3 (100%)		
M105-14	+	-	-	1/2 ¹ (50.0%)	3/3 (100%)	3/3 (100%)		
3007-85	+	+	-	0/3 (0%)	3/3 (100%)	3/3 (100%)	17/19 ^{4,5} (89.5%)	28/28 (100%)

Note: Retesting of samples with Invalid or Incomplete results was not performed

¹ Preserved specimens were prepared by diluting 1 volume of raw stool in 2.5 volumes of Cary-Blair transport medium.

² 1/3 assay replicates was reported as Invalid

³ 1/20 assay replicates was reported as Incomplete

⁴ 1/20 replicates was reported as Invalid

⁵ In the initial confirmation study with preserved specimens, both strains yielded <95% proportion positive. The confirmation study was therefore repeated at a higher target level (5 x 10⁵ vs 2.9 x 10⁵ CFU/mL).

Inclusivity (Analytical Reactivity)

The inclusivity of the PanNAT STEC Test was evaluated by testing 24 well-characterized strains of *E. coli* including seven strains of serogroup O157 and 17 strains of Shiga toxin producing non-O157 serogroups (six of which non-typeable). An O-antigen deficient laboratory strain of *E. coli* (K12) was included as a Negative Control. Each strain was tested in triplicate in unpreserved stool matrix at a concentration equivalent to 3X LoD (3 x 10⁶ CFU/mL). The results of the study are summarized in [Table 7](#). Three strains failed to produce the expected results due to mismatches in the target regions for the PanNAT STEC Test *stx1/2* primers and/or probes. A Limitation was therefore added to the device labeling to indicate that false negative results may occur due to the presence of sequence variations in the regions of the *stx1* and *stx2* genes targeted by the PanNAT STEC Test.

Table 7. Strains of *E. coli* evaluated in the PanNAT STEC Test Inclusivity Study

<i>E. coli</i> Strain	Serogroup	Genotype (subtype) ¹			PanNAT STEC Test Correct
		<i>stx1</i>	<i>stx2</i>	O157	
ATCC 35150	O157:H7	+	+	+	3/3
ATCC 43888	O157:H7	-	-	+	3/3
ATCC 43889	O157:H7	-	+	+	3/3
ATCC 43890	O157:H7	+	-	+	3/3
ATCC 43894	O157:H7	+	+	+	3/3
ATCC 700376	O157:NM	+	-	+	3/3
ATCC 700377	O157:NM	-	+	+	3/3
ATCC BAA-2326	O104:H4	-	+	-	3/3
DEC142	O165:H25	-	+ (2c/d+)	-	3/3
DEC171	ONT:H25	+ (1a)	+	-	3/3
DEC99	O121:H19	-	+ (2a)	-	3/3
MSU TW01664	O145:H16	+	-	-	3/3
MSU TW06315	O111:NM	+	+	-	3/3
MSU TW08023	O121:H19	-	+	-	3/3
MSU TW14003	O45:H2	+	-	-	3/3
TW00971	O26:H11	+	-	-	3/3
TW05672	ONT:H24	-	+	-	3/3
TW05997	O103:NM	+	-	-	3/3
TW07731	O146:H21	+ (1c)	+	-	0/3 ²
TW07740	O146:H21	+ (1d)	-	-	1/3 ³
TW07754	O146:H21	+ (1d)	-	-	1/3 ³
TW07929	ONT:HNT	+	-	-	3/3
TW07930	ONT:HNM	+	-	-	3/3
TW08111	ONT:H49	-	+	-	3/3

¹ *stx1* and *stx2* subtypes are listed in parenthesis if known

² 3/3 replicates were reported positive for *stx1* but negative for *stx2*. Sequence analysis showed that the *stx2* variant carried by this strain most closely resembled *stx2b* which is not detected by the PanNAT STEC Test.

³ Sequence analysis showed the presence of mismatches in the *stx1* target region with the PanNAT STEC Test detector probe. Both TW07740 and TW07754 were initially characterized as subtype *stx1c* but further analysis suggested that they more closely resemble strains of subtype *stx1d* which is detected with reduced sensitivity by the PanNAT STEC Test.

Bioinformatic Analysis

The inclusivity of the PanNAT STEC Test primers and probes for the *stx1*, *stx2* and *fliC* gene targets was analyzed *in silico* using the Basic Local Alignment Search Tool (BLAST). The targeted regions these genes were generally well conserved, although the potential for false-negative PanNAT STEC Test results due to sequence variation in the *stx1* and *stx2* genes and with certain subtypes is noted as a Limitation in the device labeling.

e. Analytical specificity:

Cross-reactivity Study

The analytical specificity of the PanNAT STEC Test was evaluated by testing a panel of 49 organisms and viruses that may be found in stool specimens ([Table 8](#)). Testing was performed in triplicate by spiking unpreserved stool with each potentially cross-reactive species at $\geq 10^6$ CFU or organisms/mL for bacteria, yeast and parasites and $\geq 10^5$ TCID₅₀/mL for viruses. Negative PanNAT STEC Test results were obtained

with all replicates for all the organisms and viruses tested except for a strain of *Shigella dysenteriae* that carries the *stxA* gene, a closely related homologue of *E. coli stx1*, and which is therefore expected to produce a positive result for this analyte. The potential for the PanNAT STEC Test to give positive results for Shiga toxin genes with strains of *Shigella* or other Enterobacteriaceae is noted as a Limitation in the device labeling.

Table 8. Panel of microorganisms and viruses evaluated for potential cross-reaction and/or interference in the PanNAT STEC Test

Bacteria	
<i>Acinetobacter baumannii</i>	<i>Klebsiella pneumoniae</i> ¹
<i>Aeromonas hydrophila</i>	<i>Listeria monocytogenes</i>
<i>Bacillus cereus</i>	<i>Morganella morganii</i>
<i>Bacteroides fragilis</i>	<i>Plesiomonas shigelloides</i>
<i>Campylobacter jejuni</i>	<i>Proteus mirabilis</i>
<i>Citrobacter koseri</i>	<i>Pseudomonas aeruginosa</i>
<i>Clostridium difficile</i>	<i>Salmonella enterica</i>
<i>Clostridium parabotulinum</i>	<i>Serratia marcescens</i>
<i>Clostridium perfringens</i>	<i>Shigella dysenteriae (stxA)</i> ²
<i>Enterobacter cloacae</i>	<i>Shigella flexneri</i>
<i>Enterococcus faecalis</i>	<i>Shigella sonnei</i>
<i>Enterococcus faecium</i>	<i>Staphylococcus aureus</i> ¹
<i>Escherichia coli</i> K12	<i>Staphylococcus epidermidis</i>
<i>Escherichia coli</i> O111a (EAEC)	<i>Streptococcus agalactiae</i>
<i>Escherichia coli</i> O111a:2 (EPEC)	<i>Streptococcus dysgalactiae</i>
<i>Escherichia coli</i> O29:NM (EIEC)	<i>Vibrio</i> spp.
<i>Escherichia coli</i> O8:H19 (ETEC)	<i>Yersinia enterocolitica</i>
<i>Escherichia hermanii</i>	<i>Yersinia pseudotuberculosis</i>
<i>Hafnia alvei</i>	
Viruses	
Adenovirus 18	Echovirus 11
Adenovirus 40	Enterovirus 68
Coxsackie B4	Norovirus (recombinant)
Cytomegalovirus	Rotavirus
Parasites and Yeast	
<i>Candida albicans</i> ³	<i>Entamoeba histolytica</i> (lysate)
<i>Cryptosporidium parvum</i>	<i>Giardia lamblia</i>

EAEC: Enteroaggregative *Escherichia coli*; EPEC: Enteropathogenic *Escherichia coli*; Enteroinvasive *Escherichia coli*; ETEC: Enterotoxigenic *Escherichia coli*

¹ Cross-reactivity Study: 1/3 assay replicates reported as Incomplete; not retested

² Cross-reactivity Study: 3/3 replicates reported as positive for *stx1*, negative for *stx2* and O157. This result is expected.

³ Microbial Interference Study: 1/3 assay replicates reported as Incomplete; not retested (refer to [Section M\(1\)\(g\)](#))

Bioinformatic Analysis

In silico analysis was performed to determine the potential for cross-reaction of the PanNAT STEC Test primers and probes with non-target nucleic acids. The potential for cross-reaction with the *stxA* gene of *Shigella* spp. is noted in the device labeling. Otherwise, no significant homology/complementarity was observed that was predicted to produce false positive results.

Contamination Study

The potential for false-positive results with the PanNAT STEC Test due to run-to-run contamination was evaluated by testing an alternating series of *E. coli* O157 (*stx1*⁺/*stx2*⁺/O157⁺) “high positive” (10⁷ CFU/mL) and negative unprocessed stool samples in successive runs on three PanNAT Systems. One negative sample was reported as Invalid and not retested. The expected results were obtained for all the remaining “high positive” (15/15; 100%) and negative samples (14/14; 100%). These results are acceptable.

f. Assay cut-off:

The methodology for threshold calculation was optimized through analysis of >1800 contrived specimens and validated in the Clinical Studies described in [Section M\(3\)\(a\)](#). Individual thresholds for each cartridge well and optical channel are determined based on an analysis of the changes in fluorescence during the early cycles of PCR amplification. A cycle number is then assigned where the rate of fluorescence accumulation exceeds the applicable threshold for the well and optical channel, and the accumulation of fluorescence is determined to be consistent with amplification. The decision algorithm described in [Table 1](#) is then applied to determine the result status of the sample.

g. Assay interference:

Interfering Substances

The potential for interference with the PanNAT STEC Test was evaluated with 25 endogenous and exogenous substances that may be found in stool specimens ([Table 9](#)). The study was performed using two strains of Shiga toxin producing *E. coli* O157 (EDL931 and EDL933; both *stx1*⁺/*stx2*⁺/*fliC*⁺) at 1X and 3X LoD (1 x 10⁶ and 3 x 10⁶ CFU/mL) in unpreserved stool samples, as well as with unpreserved stool samples that were negative for each of the target analytes.

False negative results were obtained in the presence of three substances at the concentrations tested: 5% (v/v) Imodium AD, 10% (w/v) Gaviscon and 48 mg/mL stearic acid. In each case, only 1/3 replicates under the test condition produced a negative result and in two cases (with Imodium AD and stearic acid), the failures were at the 1X LoD target level. No false positive results were observed under any of the test conditions.

A total of 9 Invalid results were observed during the study in association with 8 of the potentially interfering substances that were tested. Although none of the samples with Invalid results was retested, there was no evidence of systemic failures with any individual substance.

Overall, the results of the Interfering Substances Study are considered acceptable.

Table 9. Substances evaluated for potential interference with the PanNAT STEC Test

Potential Source	Substance	Concentration in Unpreserved Stool
Endogenous	Whole Blood (EDTA)	40% v/v
	Mucin	35 mg/mL
	Cholesterol	48 mg/mL
	Human DNA	10 µg/mL
	Bile	250 mg/mL
Exogenous	Kaopectate	3.5% v/v
	Pepto Bismol	5% v/v
	Imodium AD ¹	5% v/v
	Petroleum Jelly	500 mg/mL
	Baby Wipes	50% v/v
	Hydrocortisone Cream	500 mg/mL
	Diaper Rash Cream	500 mg/mL
	Stearic Acid ²	48 mg/mL
	Palmitic Acid	48 mg/mL
	Gaviscon ³	100 mg/mL
	Milk of Magnesia	100 mg/mL
	Aleve	0.5 mg/mL
	Preparation H	95 mg/mL
	Mineral Oil	500 mg/mL
	Amoxicillin	50 mg/mL
	Nystatin	4.3 mg/mL
	Metronidazole	50 mg/mL
	Ciprofloxacin	12.5 mg/mL
	Sennica Laxative	10 mg/mL
	Diaper Eluent	25% v/v

EDTA: Ethylenediaminetetraacetic acid (anticoagulant)

¹ 1/3 replicates for *E. coli* O157 strain EDL933 at 1X LoD produced a negative result (*stx1*⁺/*stx2*⁺/O157)

² 1/3 replicates for *E. coli* O157 strain EDL931 at 1X LoD produced a negative result (*stx1*⁺/*stx2*⁺/O157)

³ 1/3 replicates for *E. coli* O157 strain EDL933 at 3X LoD produced a negative result (*stx1*⁺/*stx2*⁺/O157)

Microbial Interference

The potential for interference with the PanNAT STEC Test by organisms that may be present in stool specimens was investigated using the same list of species that was evaluated for potential cross-reactivity (refer to [Section M\(1\)\(e\)](#) and [Table 8](#)).

Testing was performed with each potentially interfering species in triplicate using unpreserved stool specimens in the presence of *E. coli* O157 strain EDL933 (ATCC 43895; *stx1*⁺/*stx2*⁺/O157⁺) at 1X LoD (10⁶ CFU/mL). The potentially interfering species were tested at 10⁶ CFU or organisms/mL for bacteria, yeast and parasites and 10⁵ PFU/mL for viruses. No interference was observed with any of the microorganisms or viruses tested, although the potential for false-positive results with species of *Shigella* that carry the *stxA* gene homolog of *stx1* is noted in the device labeling.

2. Comparison studies:

a. *Method comparison with predicate device:*

Not applicable.

b. *Matrix comparison:*

To support use of frozen, retrospective unpreserved specimens to establish clinical performance and in analytical testing, a study was conducted to compare the results obtained from testing of “fresh” and frozen stool samples at different target levels. Fresh (never frozen) unpreserved stool was spiked with *E. coli* O157 strain EDL933 (*stx1⁺/stx2⁺/O157⁺*) at different target levels, aliquoted and either tested on the day of preparation or after different periods of storage at -80°C. The results of the study are presented in [Table 10](#) and showed that freezing of specimens for up to 68 weeks at -80°C had no effect on PanNAT STEC Test results. These results are acceptable.

Table 10. Summary of results from testing of unpreserved “fresh” and frozen stool samples

Target Level (CFU/mL)	Agreement with Expected Results (%)			
	“Fresh”	Duration of Storage at -80°C		
		4 Weeks	10 Weeks	68 Weeks
Negative	57/57 ¹ (100)	60/60 (100)	60/60 (100)	60/61 ^{1, 2} (98.4)
1.5X LoD (1.5 x 10 ⁶)	18/18 (100)	18/18 (100)	18/18 (100)	18/18 (100)
5X LoD (5 x 10 ⁶)	20/20 ¹ (100)	21/21 (100)	21/21 (100)	21/21 (100)
100X LoD (1 x 10 ⁸)	19/19 ³ (100)	21/21 (100)	20/20 ¹ (100)	21/21 (100)

¹ 1 replicate reported as Invalid

² 1 false positive reported for *stx1*

³ 2 replicates reported as Invalid

3. Clinical studies:

a. *Clinical Sensitivity:*

Prospective Clinical Study

The performance of the PanNAT STEC Test was evaluated in a prospective Clinical Study that was conducted at six sites in the U.S. using leftover de-identified stool specimens obtained from subjects who presented with symptoms of acute diarrhea. The specimens were either unpreserved (“fresh”) or preserved in Cary-Blair transport medium. Each specimen was divided into three aliquots of ≥250µL each: #1: for analysis with the PanNAT STEC Test (stored at 2-8°C for up to 7 days); #2: for testing by the comparator culture method (stored at 2-8°C for up to three days or

longer at -70°C); #3: for analysis in the event of discordant results between the PanNAT STEC Test and the reference method.

The reference culture method comprised separate procedures for growth and detection of Shiga-toxin producing *E. coli* (STEC) and *E. coli* O157. For detection of STEC, 50µL of unpreserved or 175µL of preserved stool was inoculated into 5mL MacConkey enrichment broth which was incubated at 37±2°C for 16-24 hours before testing with an immunoassay for the presence of Shiga Toxin 1 (ST1) and Shiga Toxin 2 (ST2) according to the manufacturer's instructions. For detection of *E. coli* O157, an aliquot of the specimen was streaked onto a Sorbitol MacConkey agar plate that was incubated at 36±2°C for 18-24 hours. Clear or beige colonies that were presumptively identified as *E. coli* O157 were sub-cultured, serotyped by latex agglutination and, if positive for the O157 antigen, their identity was confirmed by biochemical analysis.

A total of 1331 specimens from unique subjects were initially enrolled in the prospective Clinical Study. Of these, 30 (29 preserved and one (1) unpreserved) were excluded from the analysis of performance due inadvertent data loss (21) or protocol deviation (9). Of the remaining 1301, 99 (7.6%) initially produced Invalid PanNAT STEC Test results. Per protocol, all 99 specimens were retested and 78/99 (78.8%) produced valid results for a final Invalid rate of 1.6% (21/1301). Of the 1280 specimens with valid PanNAT STEC Test results, 1184 (92.5%) were preserved in Cary-Blair Transport Medium, while 96 (7.5%) were unpreserved stool specimens.

The results from testing of prospectively collected stool specimens with the PanNAT STEC Test in comparison to the reference method are shown in **Tables [11](#), [12](#) and [13](#)**.

Table 11. Results from testing prospectively collected stool specimens preserved in Cary-Blair Transport Medium with the PanNAT STEC Test in comparison to the reference method

Detection of <i>stx1</i>		Reference Method		
		Positive	Negative	Total
PanNAT STEC Test	Positive	4	8 ¹	12
	Negative	0	1172	1172
	Total	4	1180	1184
Positive Percent Agreement		100% (4/4); 95% CI: 51.0-100%		
Negative Percent Agreement		99.3% (1172/1180); 95% CI: 98.7-99.7%		
Positive Predictive Value		33.3% (4/12)		
Negative Predictive Value		100% (1172/1172)		
Prevalence		0.3% (4/1184)		
Detection of <i>stx2</i>		Reference Method		
		Positive	Negative	Total
PanNAT STEC Test	Positive	5	7 ²	12
	Negative	0	1172	1172
	Total	5	1179	1184
Positive Percent Agreement		100% (5/5); 95% CI: 56.6-100%		
Negative Percent Agreement		99.4% (1172/1179); 95% CI: 98.8-99.7%		
Positive Predictive Value		41.7% (5/12)		
Negative Predictive Value		100% (1172/1172)		
Prevalence		0.4% (5/1184)		
Detection of <i>E. coli</i> O157		Reference Method		
		Positive	Negative	Total
PanNAT STEC Test	Positive	3 ³	2 ⁴	5
	Negative	0	1179	1179
	Total	3	1181	1184
Positive Percent Agreement		100% (3/3); 95% CI: 43.9-100%		
Negative Percent Agreement		99.8% (1179/1181); 95% CI: 99.4-100%		
Positive Predictive Value		60.0% (3/5)		
Negative Predictive Value		100% (1179/1179)		
Prevalence		0.3% (3/1184)		

95% CI: 95% score confidence interval

¹ 7/8 positive for *stx1* by PCR/bidirectional sequencing

² 3/7 positive for *stx2* by PCR/bidirectional sequencing

³ 3/3 positive for *stx2*; negative for *stx1*

⁴ 1/2 positive for *E. coli* O157 by PCR/bidirectional sequencing

Table 12. Results from testing prospectively collected unpreserved stool specimens with the PanNAT STEC Test in comparison to the reference method

Detection of <i>stx1</i>		Reference Method		
		Positive	Negative	Total
PanNAT STEC Test	Positive	0	0	0
	Negative	0	96	96
	Total	0	96	96
Positive Percent Agreement		Not applicable		
Negative Percent Agreement		100% (96/96); 95% CI: 96.2-100%		
Positive Predictive Value		Not applicable		
Negative Predictive Value		100% (96/96)		
Prevalence		0% (0/96)		
Detection of <i>stx2</i>		Reference Method		
		Positive	Negative	Total
PanNAT STEC Test	Positive	2	1 ¹	3
	Negative	0	93	93
	Total	2	94	96
Positive Percent Agreement		100% (2/2); 95% CI: 34.2-100%		
Negative Percent Agreement		98.9% (93/94); 95% CI: 94.2-99.8%		
Positive Predictive Value		66.7% (2/3)		
Negative Predictive Value		100% (93/93)		
Prevalence		2.1% (2/96)		
Detection of <i>E. coli</i> O157		Reference Method		
		Positive	Negative	Total
PanNAT STEC Test	Positive	2 ²	0	2
	Negative	0	94	94
	Total	2	94	96
Positive Percent Agreement		100% (2/2); 95% CI: 34.2-100%		
Negative Percent Agreement		100% (94/94); 95% CI: 96.1-100%		
Positive Predictive Value		100% (2/2)		
Negative Predictive Value		100% (94/94)		
Prevalence		2.1% (2/96)		

95% CI: 95% score confidence interval

¹ 1/1 positive for *stx2* by PCR/bidirectional sequencing

² 2/2 positive for *stx2*; negative for *stx1*

Table 13. Results from testing of prospectively collected stool specimens with the PanNAT STEC Test, stratified by target analyte/analyte combination

Analyte(s) ¹	Agreement	
	Preserved Specimens	Unpreserved Specimens
<i>stx1</i> Only	100% (4/4)	Not applicable
<i>stx2</i> Only	100% (2/2)	Not applicable
O157 Only	Not applicable	Not applicable
<i>stx1</i> + <i>stx2</i>	Not applicable	Not applicable
<i>stx1</i> + O157	Not applicable	Not applicable
<i>stx2</i> + O157	100% (3/3)	100% (2/2)
<i>stx1</i> + <i>stx2</i> + O157	Not applicable	Not applicable
Negative	98.7% (1160/1175) ²	98.9% (93/94) ³

¹ As determined by the reference method

² Of the 15 specimens with false positive results by the PanNAT STEC Test, 8 were positive for *stx1*, 5 for *stx2* and 2 for both *stx2* and O157

³ 1/94 negative specimens was positive by the PanNAT STEC Test for *stx2*

Analysis of Retrospective Specimens

Due to the low prevalence of the target analytes in the prospective Clinical Study, additional testing was performed with retrospective (frozen), unpreserved stool specimens that were characterized using a well validated PCR/bidirectional sequencing assay to confirm the presence/absence of the *stx1*, *stx2* and/or O157 targets. One hundred and fourteen (114) retrospective specimens were initially enrolled in the study of which eight (8) were excluded from the analysis of performance due to the absence of a result for the comparator PCR/bidirectional sequencing method (7) or because previous enrollment (1). On initial testing, 8/106 (7.5%) retrospective specimens produced Invalid results with the PanNAT STEC Test. Per protocol, all were retested and 5/8 (62.5%) produced valid results upon repeat, for a final invalid rate of 2.8% (3/106). Therefore, a total of 103 retrospective specimens produced evaluable results with both the comparator method and the PanNAT STEC Test. The results of the study are summarized in **Tables 14** and **15**.

Table 14. Results from testing retrospective (frozen), unpreserved specimens with the PanNAT STEC Test

Detection of <i>stx1</i>		PCR/Bidirectional Sequencing		
		Positive	Negative	Total
PanNAT STEC Test	Positive	23	2	25
	Negative	3 ¹	75	78
	Total	26	77	103
Positive Percent Agreement		88.5% (23/26); 95% CI: 71.0-96.0%		
Negative Percent Agreement		97.4% (75/77); 95% CI: 91.0-99.3%		
Detection of <i>stx2</i>		PCR/Bidirectional Sequencing		
		Positive	Negative	Total
PanNAT STEC Test	Positive	27	1	28
	Negative	3 ²	72	75
	Total	30	73	103
Positive Percent Agreement		90.0% (27/30); 95% CI: 74.4-96.5%		
Negative Percent Agreement		98.6% (72/73); 95% CI: 92.6-99.8%		
Detection of <i>E. coli</i> O157		PCR/Bidirectional Sequencing		
		Positive	Negative	Total
PanNAT STEC Test	Positive	26	1	27
	Negative	3 ³	72	75
	Total	29	73	102 ⁴
Positive Percent Agreement		89.7% (26/29); 95% CI: 73.6-96.4%		
Negative Percent Agreement		98.6% (72/73); 95% CI: 92.6-99.8%		

95% CI: 95% score confidence interval

¹ 3/3 specimens were negative by standard of care culture but positive by PCR/bidirectional sequencing

² 2/3 specimens were negative by standard of care culture but positive by PCR/bidirectional sequencing

³ 1/3 specimens was negative by standard of care culture but positive by PCR/bidirectional sequencing

⁴ Insufficient PCR product was obtained from one specimen to enable bidirectional sequencing for the O157 target

Table 15. Results from testing retrospective (frozen), unpreserved stool specimens with the PanNAT STEC Test, stratified by target analyte/analyte combination

Analyte(s) ¹	Agreement
<i>stx1</i> Only	25.0% (1/4)
<i>stx2</i> Only	0.0% (0/2)
O157 Only	Not applicable
<i>stx1</i> + <i>stx2</i>	Not applicable
<i>stx1</i> + O157	0.0% (0/1)
<i>stx2</i> + O157	85.7% (6/7)
<i>stx1</i> + <i>stx2</i> + O157	95.2% (20/21)
Negative	98.5% (67/68) ²

¹ As determined by PCR/bidirectional sequencing

² 1/68 negative specimens was positive by the PanNAT STEC Test for *stx1*, *stx2* and O157

Analysis of Contrived Specimens

To supplement the prospective Clinical Study and analysis of retrospective specimens, additional testing was also performed with 16 contrived (spiked) unpreserved stool specimens, each of which was spiked with an enumerated stock of a different strain of Shiga-toxin producing *E. coli* or *E. coli* O157:H7 to a final concentration of 5×10^6 CFU/mL of stool (5X LoD). The contrived specimens were frozen and tested at the clinical sites, interspersed with retrospective specimens (above). No Invalid results were obtained with the contrived specimens (Invalid rate 0/16 = 0%). The results of the study are summarized in **Tables 16** and **17**.

Table 16. Results from testing frozen, contrived specimens with the PanNAT STEC Test

Detection of <i>stx1</i>		Expected Result		
		Positive	Negative	Total
PanNAT STEC Test	Positive	12	0	12
	Negative	0	4	4
	Total	12	4	16
Positive Percent Agreement		100% (12/12); 95% CI: 75.8-100%		
Negative Percent Agreement		100% (4/4); 95% CI: 51.0-100%		
Detection of <i>stx2</i>		Expected Result		
		Positive	Negative	Total
PanNAT STEC Test	Positive	8	0	8
	Negative	0	8	8
	Total	8	8	16
Positive Percent Agreement		100% (8/8); 95% CI: 67.6-100%		
Negative Percent Agreement		100% (8/8); 95% CI: 67.6-100%		
Detection of <i>E. coli</i> O157		Expected Result		
		Positive	Negative	Total
PanNAT STEC Test	Positive	2	0	2
	Negative	0	14	14
	Total	2	14	16
Positive Percent Agreement		100% (2/2); 95% CI: 34.2-100%		
Negative Percent Agreement		100% (14/14); 95% CI: 78.5-100%		

95% CI: 95% score confidence interval

Table 17. Results from testing contrived, unpreserved stool specimens with the PanNAT STEC Test, stratified by target analyte/analyte combination

Analytes ¹	Positive Agreement
<i>stx1</i> Only	100% (7/7)
<i>stx2</i> Only	100% (3/3)
O157 Only	100% (1/1)
<i>stx1</i> + <i>stx2</i>	100% (4/4)
<i>stx1</i> + O157	Not applicable
<i>stx2</i> + O157	Not applicable
<i>stx1</i> + <i>stx2</i> + O157	100% (1/1)

¹ Expected results based on genetic determinants known to be present within each strain

b. Clinical specificity:

Refer to [Section M\(3\)\(a\)](#), above.

c. Other clinical supportive data (when a. and b. are not applicable):

Not applicable.

4. Clinical cut-off:

Not applicable

4. Expected values/Reference range:

The performance of the PanNAT STEC Test was evaluated in a prospective Clinical Study conducted at six (6) sites in the US ([Section M\(3\)\(a\)](#)). The overall prevalence of *stx1*, *stx2* and *E. coli* O157 in stool specimens was 0.9%, 1.2% and 0.5%, respectively as determined by the PanNAT assay and 0.3%, 0.5% and 0.4% respectively as determined by the reference culture method. In [Table 18](#), the prevalence of each analyte as determined by the PanNAT assay is stratified by the age of the subjects.

Table 18. Prevalence of *stx1*, *stx2* and *E. coli* O157 positive subjects observed in the prospective PanNAT STEC Test Clinical Study, stratified by age of the subjects

Age	Number	PanNAT STEC Test Positive (% Prevalence)		
		<i>stx1</i>	<i>stx2</i>	<i>E. coli</i> O157
≤5 years	150	1 (0.7)	5 (3.3)	4 (2.7)
6-21 years	191	4 (2.1)	3 (1.6)	1 (0.5)
22-59 years	495	5 (1.1)	6 (1.2)	2 (0.4)
≥60 years	444	2 (0.5)	1 (0.2)	0 (0)
Total	1280	12 (0.9)	15 (1.2)	7 (0.5)

N. Instrument Name:

PanNAT System

O. System Descriptions:

1. Modes of Operation:

Does the applicant’s device contain the ability to transmit data to a computer, webserver, or mobile device?

Yes _____ or No _____

Does the applicant’s device transmit data to a computer, webserver, or mobile device using wireless transmission?

Yes _____ or No _____

2. Software:

FDA has reviewed applicant’s Hazard Analysis and software development processes for this line of product types:

Yes _____ or No _____

3. Specimen Identification:

Specimens are identified by scanning or typing in the applicable patient identification code.

4. Specimen Sampling and Handling:

Using the swab provided in the PanNAT STEC Sample transfer pack, the operator transfers a portion of the stool specimen to a Sample Buffer Tube, seals the tube using an Adaptor Cap and attaches the cap to the inlet port of the PanNAT STEC Test cartridge. The operator then inserts the cartridge with attached Sample Buffer Tube into the

PanNAT System for automated processing.

5. Calibration:

The PanNAT System is factory calibrated and does not require user calibration. A system check using Positive and Negative Controls is recommended if the instrument is moved after initial installation.

6. Quality Control:

Refer to [Section M\(1\)\(c\)](#).

P. Other Supportive Instrument Performance Characteristics Data Not Covered In The “Performance Characteristics” Section above:

Not applicable.

Q. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of [21 CFR Parts 801](#) and [809](#), as applicable and the special controls for this device type under [21 CFR 866.3990](#).

R. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.