



BRUKER DALTONICS, INC  
MARKUS KOSTRZEWA  
VICE PRESIDENT CLINICAL MASS SPECTROMETRY R&D  
FAHRENHEILSTRASSE-4  
BREMEN 28359  
DE

March 27, 2015

Re: K142677  
Trade/Device Name: MALDI Biotyper CA System  
Regulation Number: 21 CFR 866.3361  
Regulation Name: Mass spectrometer system for clinical use for the identification of  
microorganisms  
Regulatory Class: II  
Product Code: PEX  
Dated: February 27, 2015  
Received: March 2, 2015

Dear Dr. Kostrzewa:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the

electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Uwe Scherf -S** for

Sally Hojvat, M.Sc., PhD.

Director

Division of Microbiology Devices

Office of In Vitro Diagnostics

and Radiological Health

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K142677

Device Name

MALDI Biotyper CA System

Indications for Use (Describe)

The Bruker Daltonics, Inc. MALDI Biotyper CA System is a mass spectrometer system using matrix-assisted laser desorption/ionization – time-of-flight (MALDI-TOF) for the identification of microorganisms cultured from human specimens.

The MALDI Biotyper CA System is a qualitative in vitro diagnostic device indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of bacterial and yeast infections.

### BACTERIA

Achromobacter xylosoxidans  
Acinetobacter baumannii complex [4]  
Acinetobacter haemolyticus  
Acinetobacter johnsonii  
Acinetobacter junii  
Acinetobacter lwoffii  
Acinetobacter radioresistens  
Acinetobacter ursingii  
Actinomyces meyeri  
Actinomyces neuii  
Actinomyces odontolyticus  
Actinomyces oris  
Aerococcus urinae  
Aerococcus viridans  
Aeromonas salmonicida  
Aeromonas sp[7]  
Alcaligenes faecalis  
Anaerococcus vaginalis  
Bacteroides fragilis  
Bacteroides ovatus group  
Bacteroides thetaiotaomicron group  
Bacteroides uniformis  
Bacteroides vulgatus group  
Bordetella group[3]  
Bordetella hinzii  
Brevibacterium casei  
Brevundimonas diminuta group  
Burkholderia cepacia complex [13]  
Burkholderia gladioli  
Burkholderia multivorans  
Campylobacter coli  
Campylobacter jejuni  
Campylobacter ureolyticus  
Capnocytophaga ochracea  
Capnocytophaga sputigena

---

Chryseobacterium gleum  
Chryseobacterium indologenes  
Citrobacter amalonaticus complex  
Citrobacter freundii complex  
Citrobacter koseri  
Clostridium difficile  
Clostridium perfringens  
Corynebacterium amycolatum  
Corynebacterium aurimucosum group  
Corynebacterium bovis  
Corynebacterium diphtheriae  
Corynebacterium glucuronolyticum  
Corynebacterium jeikeium  
Corynebacterium kroppenstedtii  
Corynebacterium macginleyi  
Corynebacterium minutissimum  
Corynebacterium propinquum  
Corynebacterium pseudodiphtheriticum  
Corynebacterium riegeli  
Corynebacterium striatum group  
Corynebacterium tuberculostearicum  
Corynebacterium ulcerans  
Corynebacterium urealyticum  
Corynebacterium xerosis  
Cronobacter sakazakii group  
Cupriavidus pauculus group  
Delftia acidovorans group  
Dermacoccus nishinomiyaensis  
Edwardsiella tarda  
Eikenella corrodens  
Elizabethkingia meningoseptica group  
Enterobacter aerogenes  
Enterobacter amnigenus  
Enterobacter cloacae complex  
Enterococcus avium group  
Enterococcus casseliflavus  
Enterococcus faecalis  
Enterococcus faecium  
Enterococcus gallinarum  
Enterococcus hirae  
Escherichia coli  
Finegoldia magna  
Fusobacterium canifelinum  
Fusobacterium necrophorum  
Fusobacterium nucleatum  
Gardnerella vaginalis  
Gemella haemolysans  
Gemella sanguinis  
Granulicatella adiacens  
Haemophilus haemolyticus  
Haemophilus influenzae  
Haemophilus parahaemolyticus group  
Haemophilus parainfluenzae  
Hafnia alvei

---

Kingella kingae  
Klebsiella oxytoca / Raoultella ornithinolytica  
Klebsiella pneumoniae  
Kocuria kristinae  
Kytococcus sedentarius  
Lactococcus garvieae  
Lactococcus lactis  
Leuconostoc mesenteroides  
Macrococcus caseolyticus  
Micrococcus luteus  
Moraxella sg Branhamella catarrhalis  
Moraxella sg Moraxella nonliquefaciens  
Moraxella sg Moraxella osloensis  
Morganella morganii  
Myroides odoratimimus  
Myroides odoratus  
Oligella ureolytica  
Oligella urethralis  
Pantoea agglomerans  
Parabacteroides distasonis  
Pasteurella multocida  
Pediococcus pentosaceus  
Peptoniphilus harei group  
Peptostreptococcus anaerobius  
Plesiomonas shigelloides  
Porphyromonas gingivalis  
Prevotella bivia  
Prevotella buccae  
Prevotella denticola  
Prevotella intermedia  
Prevotella melaninogenica  
Propionibacterium acnes  
Proteus mirabilis  
Proteus vulgaris group  
Providencia rettgeri  
Providencia stuartii  
Pseudomonas aeruginosa  
Pseudomonas fluorescens group  
Pseudomonas oryzihabitans  
Pseudomonas putida group  
Pseudomonas stutzeri  
Rhizobium radiobacter  
Rothia aerea  
Rothia dentocariosa  
Rothia mucilaginosa  
Salmonella sp  
Serratia liquefaciens  
Serratia marcescens  
Serratia plymuthica  
Serratia rubidaea  
Staphylococcus aureus  
Staphylococcus auricularis  
Staphylococcus capitis  
Staphylococcus caprae

---

Staphylococcus carnosus  
Staphylococcus cohnii  
Staphylococcus epidermidis  
Staphylococcus equorum  
Staphylococcus felis  
Staphylococcus haemolyticus  
Staphylococcus hominis  
Staphylococcus lugdunensis  
Staphylococcus pasteurii  
Staphylococcus pettenkoferi  
Staphylococcus pseudintermedius  
Staphylococcus saccharolyticus  
Staphylococcus saprophyticus  
Staphylococcus schleiferi  
Staphylococcus simulans  
Staphylococcus vitulinus  
Staphylococcus warneri  
Stenotrophomonas maltophilia  
Streptococcus agalactiae  
Streptococcus anginosus  
Streptococcus constellatus  
Streptococcus dysgalactiae  
Streptococcus gallolyticus  
Streptococcus gordonii  
Streptococcus intermedius  
Streptococcus lutetiensis  
Streptococcus mitis / oralis group  
Streptococcus mutans  
Streptococcus pneumoniae  
Streptococcus pyogenes  
Streptococcus salivarius  
Sutterella wadsworthensis  
Vibrio parahaemolyticus  
Vibrio vulnificus  
Yersinia enterocolitica  
Yersinia pseudotuberculosis

#### YEAST

Candida albicans  
Candida boidinii  
Candida dubliniensis  
Candida duobushaemulonii  
Candida famata  
Candida glabrata  
Candida guilliermondii  
Candida haemulonis  
Candida inconspicua  
Candida kefyr  
Candida krusei  
Candida lambica  
Candida lipolytica  
Candida lusitanae  
Candida metapsilosis  
Candida norvegensis

---

Candida orthopsilosis  
Candida parapsilosis  
Candida pararugosa  
Candida pelliculosa  
Candida tropicalis  
Candida valida  
Cryptococcus gattii  
Cryptococcus neoformans var grubii  
Cryptococcus neoformans var neoformans  
Geotrichum candidum  
Geotrichum capitatum  
Kloeckera apiculata  
Pichia ohmeri  
Saccharomyces cerevisiae  
Trichosporon asahii

---

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



## 510(k) SUMMARY

**Date of Summary:** March 27, 2015

**Product Name** MBT-CA System

**Sponsor** Bruker Daltonics, Inc.  
40 Manning Road,  
Billerica, MA 01821

**Correspondent** Bruker Daltonik GmbH  
Markus Kostrzewa, Vice President Clinical Mass Spectrometry R&D  
Fahrenheitstrasse 4  
Bremen, 28359  
Phone: +49 (0)421-2205 1258 Fax: +49 (0)421-2205 205  
Email: Markus.Kostrzewa@bruker.com

**Device Identification**

*Trade or Proprietary Name:* MALDI Biotyper CA System

*Common or Usual Name:* System, mass spectrometry, maldi tof, microorganism identification, cultured isolates

*Product Code:* PEX

*Regulation Section:* 21 CFR 866.3361

*Device Class:* Class II

*Panel:* Microbiology



## Substantial Equivalency

The Bruker Daltonics, Inc. MBT-CA System is substantially equivalent to the bioMérieux Vitek MS (K124067) and the Bruker Daltonics, Inc. MBT-CA System (K130831). Table 1 compares the characteristics of the MBT-CA System (New Device) and the Vitek MS (Predicate Device).

**Table 1. Substantial Equivalency Table**

<i>Similarities</i>			
Characteristic	NEW DEVICE Bruker Daltonics, Inc. MBT-CA System (TBD)	PRIMARY PREDICATE DEVICE Vitek® MS (K124067)	PREDICATE DEVICE Bruker Daltonics, Inc. MBT-CA System (K130831)
Product Codes	PEX	PEX	PEX
Intended use	<p>The Bruker Daltonics, Inc. MALDI Biotyper CA System is a mass spectrometer system using matrix-assisted laser desorption/ionization – time-of-flight (MALDI-TOF) for the identification of microorganisms cultured from human specimens.</p> <p>The MALDI Biotyper CA System is a qualitative <i>in vitro</i> diagnostic device indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of bacterial and yeast infections.</p>	<p>The Vitek® MS is a mass spectrometer system using matrix-assisted laser desorption/ionization-time-of-flight (MALDI-TOF) for the identification of microorganisms cultured from human specimen.</p> <p>The VITEK MS is a qualitative <i>in vitro</i> diagnostic device indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of bacterial and yeast infections.</p>	See “Differences”
Sample type	<p>Isolated colony from any patient sample source.</p> <p>Acceptable media:</p> <ul style="list-style-type: none"> <li>• Columbia blood agar with 5% sheep blood</li> <li>• Trypticase soy agar with 5% sheep Blood</li> <li>• Chocolate agar</li> <li>• MacConkey Agar</li> <li>• Columbia CNA agar with 5% sheep blood</li> <li>• Brucella Agar with 5% horse blood</li> <li>• CDC anaerobe Agar with 5% sheep blood</li> <li>• CDC anaerobe 5% sheep blood Agar with phenylethyl alcohol</li> <li>• CDC anaerobe laked sheep blood Agar with kanamycin and vancomycin</li> <li>• Bacteroides bile esculin Agar with amikacin</li> <li>• Clostridium difficile Agar with 7% sheep blood</li> <li>• Sabouraud-Dextrose Agar</li> <li>• Brain Heart Infusion Agar</li> <li>• Campylobacter Agar with 5 Antimicrobics and 10% Sheep Blood</li> <li>• Bordet Gengou Agar with 15% sheep blood</li> </ul>	<p>Isolated colony from any patient sample source.</p> <p>Acceptable media:</p> <ul style="list-style-type: none"> <li>• Columbia blood agar with 5% sheep blood</li> <li>• Trypticase soy agar with 5% sheep Blood</li> <li>• Chocolate polyvitex agar</li> <li>• Campylosel agar</li> <li>• MacConkey Agar</li> <li>• Modified Sabouraud dextrose Agar</li> <li>• ChromID CPS</li> </ul>	<p>Isolated colony from any patient sample source.</p> <p>Acceptable media:</p> <ul style="list-style-type: none"> <li>• Columbia blood agar with 5% sheep blood</li> <li>• Trypticase soy agar with 5% sheep Blood</li> <li>• Chocolate agar</li> <li>• MacConkey Agar</li> </ul>
Type of Test	Automated Mass Spectrometry System	Automated Mass Spectrometry System	Automated Mass Spectrometry System



<i>Similarities</i>			
Characteristic	NEW DEVICE Bruker Daltonics, Inc. MBT-CA System (TBD)	PRIMARY PREDICATE DEVICE Vitek® MS (K124067)	PREDICATE DEVICE Bruker Daltonics, Inc. MBT-CA System (K130831)
Matrix	α-Cyano-4-hydroxycinnamic acid	α-Cyano-4-hydroxycinnamic acid	α-Cyano-4-hydroxycinnamic acid
Method of Testing	Bacteria & Yeast: Direct testing  If after initial analysis the log(score) is reported at <2.00, organisms may be processed using the Extraction (Ext) procedure or extended Direct Transfer (eDT, 70% aqueous formic acid) procedure. If eDT procedure still yields log(score) <2.00, organism may be processed via Ext procedure.	Bacteria & Yeast: Direct testing	Bacteria: Direct testing  If after initial analysis the log(score) is reported at <2.00, organisms are processed using the Extraction procedure.
Result Reporting	Organism identification is reported with high confidence if the log(score) is ≥2.00.  An organism identification is reported with low confidence if the log(score) is between 1.70 and <2.00.	A single identification is displayed, with a confidence value from 60.0 to 99.9, when one significant organism or organism group is retained.  “Low-discrimination” identifications are displayed when more than one but not more than four significant organisms or organism groups are retained.  When more than four organisms or organism groups are found, or when no match is found, the organism is considered unidentified.	Organism identification is reported with high confidence if the log(score) is ≥2.00.  An organism identification is reported with low confidence if the log(score) is between 1.70 and <2.00.
Matching Algorithm	Calculates matches by comparing a new spectrum against each single reference entry of a reference database.	Uses a proprietary process called "mass binning." In this process, the spectrum between 3,000 and 17,000 Daltons are divided into 1300 pre-defined intervals called “bins”. Next, an algorithm based on supervised machine learning known as the “Advanced Spectrum Classifier”, is used to determine how informative each bin was in differentiating that species from all other species in the database.	Calculates matches by comparing a new spectrum against each single reference entry of a reference database.
Recorded mass range	2,000 - 20,000 m/z	2,000 - 20,000 m/z	2,000 - 20,000 m/z



<i>Similarities</i>			
Characteristic	NEW DEVICE Bruker Daltonics, Inc. MBT-CA System (TBD)	PRIMARY PREDICATE DEVICE Vitek® MS (K124067)	PREDICATE DEVICE Bruker Daltonics, Inc. MBT-CA System (K130831)
Calibration	Bruker US IVD Bacterial Test Standard (BTS)	See "Differences"	Bruker US IVD Bacterial Test Standard (BTS)
MALDI Target Plate	US IVD 48 Spot Target <ul style="list-style-type: none"> <li>• 48 positions reusable steel targets</li> </ul>	See "Differences"	US IVD 48 Spot Target <ul style="list-style-type: none"> <li>• 48 positions reusable steel targets</li> </ul>
MALDI-TOF MS instruments	Bruker microflex (benchtop)	See "Differences"	Bruker microflex (benchtop)
Database	MALDI Biotyper for Clinical Applications (MBT-CA)	See "Differences"	MALDI Biotyper for Clinical Applications (MBT-CA)

<i>Differences</i>			
Characteristic	NEW DEVICE Bruker Daltonics, Inc. MBT-CA System (TBD)	PREDICATE DEVICE Vitek® MS (K124067)	PREDICATE DEVICE Bruker Daltonics, Inc. MBT-CA System (K130831)
Culture Age	Bacteria and yeasts growth should be between 18h to 48h (+12h storage at RT) Specific requirements: <ul style="list-style-type: none"> <li>• <i>Bordetella</i>: Incubation on BG agar should not be longer than 24h (+12h storage at RT).</li> <li>• <i>Campylobacter</i>: Incubation can be prolonged to 72h (+12h storage at RT).</li> <li>• <i>Streptococcus pneumoniae</i>: Incubation should not be longer than 24h (+12h storage at RT) due to possible autolysis.</li> </ul>	Bacteria and yeast growth should be between 24 to 72 hours.	Bacteria growth should be between 18h to 36h
Calibration	Bruker US IVD Bacterial Test Standard (BTS)	<i>E. coli</i> ATCC 8739	See "Similarities"
MALDI Target Plate	US IVD 48 Spot Target <ul style="list-style-type: none"> <li>• 48 positions reusable steel targets</li> </ul>	VITEK MS-DS Target Slides <ul style="list-style-type: none"> <li>• 48 positions disposable plastic targets</li> </ul>	See "Similarities"
MALDI-TOF MS instruments	Bruker microflex (benchtop)	Shimadzu AXIMA® Assurance MS (floor standing)	See "Similarities"
Database	MALDI Biotyper for Clinical Applications (MBT-CA)	VITEK® MS V2.0 Knowledge Base	See "Similarities"



<i>Differences</i>			
Characteristic	NEW DEVICE Bruker Daltonics, Inc. MBT-CA System (TBD)	PREDICATE DEVICE Vitek® MS (K124067)	PREDICATE DEVICE Bruker Daltonics, Inc. MBT-CA System (K130831)
Intended Use	<p>The Bruker Daltonics, Inc. MALDI Biotyper CA System is a mass spectrometer system using matrix-assisted laser desorption/ionization – time-of-flight (MALDI-TOF) for the identification of microorganisms cultured from human specimens.</p> <p>The MALDI Biotyper CA System is a qualitative <i>in vitro</i> diagnostic device indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of bacterial and yeast infections.</p>	See “Similarities”	<p>The Bruker Daltonics, Inc. MALDI Biotyper CA System is a qualitative <i>in vitro</i> diagnostic mass spectrometer system for the identification of Gram-negative bacterial colonies cultured from human specimens using matrix-assisted laser desorption/ionization – time-of-flight (MALDI-TOF) mass spectrometry technology.</p> <p>The MALDI Biotyper CA System is indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of Gram-negative bacterial infections.</p>

These differences do not affect substantial equivalence of the MBT-CA System, Vitek® MS system and MBT-CA System (K130831). All systems are mass spectrometers using matrix-assisted laser desorption/ionization – time-of-flight (MALDI-TOF) for the identification of microorganisms cultured from human specimens. The differences noted above do not impact the intended use and do not raise questions as to the safety and effectiveness of the test (new) device.

## Intended Use

The Bruker Daltonics, Inc. MALDI Biotyper CA System is a mass spectrometer system using matrix-assisted laser desorption/ionization – time-of-flight (MALDI-TOF) for the identification of microorganisms cultured from human specimens.

The MALDI Biotyper CA System is a qualitative *in vitro* diagnostic device indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of bacterial and yeast infections.

The following organisms are claimed:

### Bacteria:

Achromobacter xylosoxidans	Cupriavidus pauculus group	Propionibacterium acnes
Acinetobacter haemolyticus	Delftia acidovorans group	Proteus mirabilis
Acinetobacter johnsonii	Dermacoccus nishinomiyaensis	Proteus vulgaris group
Acinetobacter junii	Edwardsiella tarda	Providencia rettgeri
Acinetobacter lwoffii	Eikenella corrodens	Providencia stuartii
Acinetobacter radioresistens	Elizabethkingia meningoseptica group	Pseudomonas aeruginosa
Acinetobacter ursingii	Enterobacter aerogenes	Pseudomonas fluorescens group
Acinetobacter baumannii complex [4]	Enterobacter amnigenus	Pseudomonas oryzihabitans
Actinomyces meyeri	Enterobacter cloacae complex	Pseudomonas putida group
Actinomyces neuii	Enterococcus casseliflavus	Pseudomonas stutzeri
Actinomyces odontolyticus	Enterococcus faecalis	Rhizobium radiobacter
Actinomyces oris	Enterococcus faecium	Rothia aeria
Aerococcus urinae	Enterococcus gallinarum	Rothia dentocariosa
Aerococcus viridans	Enterococcus hirae	Rothia mucilaginosa
Aeromonas salmonicida	Enterococcus avium group	Salmonella sp
Aeromonas sp[7]	Escherichia coli	Serratia liquefaciens
Alcaligenes faecalis	Finegoldia magna	Serratia marcescens
Anaerococcus vaginalis	Fusobacterium canifelinum	Serratia plymuthica
Bacteroides fragilis	Fusobacterium necrophorum	Serratia rubidaea
Bacteroides uniformis	Fusobacterium nucleatum	Staphylococcus aureus
Bacteroides ovatus group	Gardnerella vaginalis	Staphylococcus auricularis
Bacteroides thetaiotaomicron group	Gemella haemolysans	Staphylococcus capitis
Bacteroides vulgatus group	Gemella sanguinis	Staphylococcus caprae
Bordetella group[3]	Granulicatella adiacens	Staphylococcus carnosus
Bordetella hinzii	Haemophilus haemolyticus	Staphylococcus cohnii
Brevibacterium casei	Haemophilus influenzae	Staphylococcus epidermidis
Brevundimonas diminuta group	Haemophilus parainfluenzae	Staphylococcus equorum
Burkholderia gladioli	Haemophilus parahaemolyticus group	Staphylococcus felis
Burkholderia multivorans	Hafnia alvei	Staphylococcus haemolyticus
Burkholderia cepacia complex [13]	Kingella kingae	Staphylococcus hominis
Campylobacter coli	Klebsiella pneumoniae	Staphylococcus lugdunensis
Campylobacter jejuni	Klebsiella oxytoca / Raoultella ornithinolytica	Staphylococcus pasteurii
Campylobacter ureolyticus	Kocuria kristinae	Staphylococcus pettenkoferi
Capnocytophaga ochracea	Kytococcus sedentarius	Staphylococcus pseudintermedius
Capnocytophaga sputigena	Lactococcus garvieae	Staphylococcus saccharolyticus

Chryseobacterium gleum	Lactococcus lactis	Staphylococcus saprophyticus
Chryseobacterium indologenes	Leuconostoc mesenteroides	Staphylococcus schleiferi
Citrobacter amalonaticus complex	Macrocococcus caseolyticus	Staphylococcus simulans
Citrobacter koseri	Micrococcus luteus	Staphylococcus vitulinus
Citrobacter freundii complex	Moraxella sg Branhamella catarrhalis	Staphylococcus warneri
Clostridium difficile	Moraxella sg Moraxella nonliquefaciens	Stenotrophomonas maltophilia
Clostridium perfringens	Moraxella sg Moraxella osloensis	Streptococcus agalactiae
Corynebacterium amycolatum	Morganella morgani	Streptococcus anginosus
Corynebacterium bovis	Myroides odoratimimus	Streptococcus constellatus
Corynebacterium diphtheriae	Myroides odoratus	Streptococcus dysgalactiae
Corynebacterium glucuronolyticum	Oligella ureolytica	Streptococcus gallolyticus
Corynebacterium jeikeium	Oligella urethralis	Streptococcus gordonii
Corynebacterium kroppenstedtii	Pantoea agglomerans	Streptococcus intermedius
Corynebacterium macginleyi	Parabacteroides distasonis	Streptococcus lutetiensis
Corynebacterium minutissimum	Pasteurella multocida	Streptococcus mutans
Corynebacterium propinquum	Pediococcus pentosaceus	Streptococcus pneumoniae
Corynebacterium pseudodiphtheriticum	Peptoniphilus harei group	Streptococcus pyogenes
Corynebacterium riegelii	Peptostreptococcus anaerobius	Streptococcus salivarius
Corynebacterium tuberculostearicum	Plesiomonas shigelloides	Streptococcus mitis / oralis group
Corynebacterium ulcerans	Porphyromonas gingivalis	Sutterella wadsworthensis
Corynebacterium urealyticum	Prevotella bivia	Vibrio parahaemolyticus
Corynebacterium xerosis	Prevotella buccae	Vibrio vulnificus
Corynebacterium aurimucosum group	Prevotella denticola	Yersinia enterocolitica
Corynebacterium striatum group	Prevotella intermedia	Yersinia pseudotuberculosis
Cronobacter sakazakii group	Prevotella melaninogenica	

### Yeast:

Candida albicans	Candida parapsilosis
Candida boidinii	Candida pararugosa
Candida dubliniensis	Candida pelliculosa
Candida duobushaemulonii	Candida tropicalis
Candida glabrata	Candida valida
Candida famata	Cryptococcus gattii
Candida guilliermondii	Cryptococcus neoformans_var_grubii
Candida haemulonis	Cryptococcus neoformans_var_neoformans
Candida inconspicua	Geotrichum candidum
Candida kefyr	Geotrichum capitatum
Candida krusei	Kloeckera apiculata
Candida lambica	Pichia ohmeri
Candida lipolytica	Saccharomyces cerevisiae
Candida lusitanae	Trichosporon asahii
Candida metapsilosis	
Candida norvegensis	
Candida orthopsilosis	

## Methodology

Biochemical methods are currently the most commonly used methods for the identification of microorganisms. Organisms are tested against a range of reagents and organism identification is based on a microorganism's reaction to these reagents.

The MBT-CA System uses a different methodology for organism identification based on unique protein patterns of the microorganisms obtained from mass spectrometry. The test organism's spectrum (a pattern of mass peaks) is compared with a reference spectra library (database). Using biostatistical analysis, a probability ranking of the organism identification is generated. The probability ranking is represented as a log(score) between 0.00 and 3.00. Organism identification is reported with high confidence if the log(score) is  $\geq 2.00$ . An organism identification is reported with low confidence if the log(score) is between 1.70 and  $< 2.00$ .

Some MBT-CA identifications displayed are non-clinically validated organisms. In the interest of public health, these organisms are displayed but are grayed out in the MBT-CA report as a means of directing the required additional laboratory testing. These results are not reported; identifications must be confirmed using alternate laboratory methods. Results for non-clinically validated organisms cannot be transmitted from the MBT-CA to the laboratory information system.

Organisms to be identified with the MBT-CA System should be isolated for purity on appropriate isolation media.

*Direct Transfer (DT):* An individual colony from a subculture plate is transferred to a selected position on an US IVD 48 Spot Target plate (target) and overlaid with US IVD HCCA portioned (matrix). The standard solvent (50% acetonitrile / 47.5% H<sub>2</sub>O / 2.5% trifluoroacetic acid) in the matrix solution extracts proteins (mainly ribosomal proteins, which are present in high concentration) from the microorganisms. When dried matrix crystallizes, the inoculated target is ready to be analyzed on the MBT-CA System. If after initial analysis the log(score) is reported as  $< 2.00$ , organisms can be processed using the extended Direct Transfer (eDT) procedure or the Extraction (Ext) procedure and analysis repeated. If eDT is employed and log(score) is reported as  $< 2.00$ , reanalysis via the Extraction procedure may be used.

### *extended Direct Transfer (eDT):*

If DT analysis yields a (log(score)  $< 2.00$ ) result, an individual colony from a subculture plate may be transferred to a selected position on a target and overlaid with 70% aqueous formic acid solution. The target is air-dried and then matrix is overlaid. When dried matrix crystallizes, the inoculated target is ready to be analyzed on the MBT-CA System. If a high confidence result is not achieved (log(score) is reported at  $< 2.00$ ), organisms can be processed using the Extraction procedure and analysis repeated.

*Extraction procedure (Ext):* If after initial analysis and eDT procedure the log(score) is reported at <2.00, organisms are processed using the Extraction procedure and analysis repeated. For this purpose, isolated colonies from the subculture plate are extracted in accordance with MBT-CA System user manual. Afterwards they are transferred to the target and treated as described above.

*MALDI-TOF Analysis:*

Samples are analyzed using MALDI (matrix-assisted laser desorption/ionization) TOF (time-of-flight) mass spectrometry. The matrix transfers protons onto the extracted proteins and absorbs UV light. A laser in the MALDI-TOF mass spectrometer irradiates the matrix sample composite, causing evaporation and release of positively charged intact proteins and peptides ("soft" ionization technique). These ions are electrostatically accelerated over a short distance and arrive in the flight tube at a mass-dependent speed. As different proteins/peptides have different masses, ions arrive at the detector at different times (time-of-flight). The system measures the time (in the nanosecond range) between pulsed acceleration and the corresponding detector signal, the speed is converted into an exact molecular mass. The mass-to-charge ratio of an ion is proportional to the square of its drift time.

Highly abundant microbial proteins (mainly ribosomal proteins) result in a mass spectrum with characteristic mass and intensity distribution. It is species-specific for many bacteria and is interpreted as a molecular fingerprint to identify the test organism.

Data acquisition is controlled with MBT-CA Software. The spectrum of the unknown organism is first transformed into a peak list. Using a biostatistical algorithm, this peak list is compared to the reference peak lists of organisms in the reference library (database) and a log(score) is generated. A higher log(score) indicates a higher degree of similarity to the organism in the reference library. Organism identification is reported with high confidence if the log(score) is  $\geq 2.00$ . An organism identification is reported with low confidence if the log(score) is between 1.70 and <2.00.

The log(score) ranges, defined in the MBT-CA System, are indicative of the probability of organism identification. Results should be reviewed by a trained microbiologist and final organism identification should be based on all relevant information available. This information includes but is not limited to: Gram staining, colony morphology, growth characteristics, sample matrix, etc.

### **Performance Data**

#### **Precision/Repeatability:**

Validation of the whole MALDI Biotyper CA System was performed on six (6) working days with two (2) runs/day following manufacturer's instructions for use. Ten (10) test organisms were tested in triplicate via Direct Transfer (DT) and extended Direct Transfer (eDT) in each run. If a replicate yielded a log(score) <2.00, the test organism was repeated in triplicate via Extraction. The study also tested multiple sources of system variability including two (2) test operators, two

(2) microflex LT/SH instruments and two (2) target plates. Overall results from the precision/repeatability study are presented below.

**Table 2: Overall Precision per Test Organism**

Test organism	# samples measured	# samples $\geq 2.0$ ID (DT)	# samples $\geq 2.0$ ID (eDT)	# samples $\geq 2.0$ ID (DT/eDT+Ext)
<i>Brevibacterium casei</i>	36	36 (100%)	36 (100%)	<b>36 (100%)</b>
<i>Enterococcus faecalis</i>	36	34 (94%)	36 (100%)	<b>36 (100%)</b>
<i>Micrococcus luteus</i>	36	21 (58%)	36 (100%)	<b>36 (100%)</b>
<i>Staphylococcus aureus</i>	36	36 (100%)	36 (100%)	<b>36 (100%)</b>
<i>Staphylococcus epidermidis</i>	36	36 (100%)	36 (100%)	<b>36 (100%)</b>
<i>Streptococcus agalactiae</i>	36	34 (94%)	36 (100%)	<b>36 (100%)</b>
<i>Candida albicans</i>	36	18 (50%)	30 (83%)	<b>36 (100%)</b>
<i>Candida parapsilosis</i>	36	6 (17%)	32 (89%)	<b>36 (100%)</b>
<i>Candida tropicalis</i>	36	34 (94%)	35 (97%)	<b>36 (100%)</b>
<i>Saccharomyces cerevisiae</i>	36	20 (56%)	27 (75%)	<b>36 (100%)</b>

Based upon the data presented, the study confirms repeatability and precision of the MALDI Biotyper CA System independent from:

- System Operators
- microflex LT/SH instruments
- Target plates

Limit of Detection/Dynamic Range:

The Limit of Detection study was conducted to estimate the dynamic range (in terms of sample amount) of Gram-positive bacteria and yeasts to be identified on the MALDI Biotyper CA System. Six (6) frequently occurring clinically relevant test organisms [three (3) Gram-positive and three (3) yeast] were chosen for this study. [NOTE: Due to the nature of yeast organisms, dynamic range studies using known yeast concentration was not feasible for the Direct Transfer procedure].

Turbidity measurements of stock suspensions containing microbial material were performed at a wavelength of 600 nm. To determine the amount of cfu/ $\mu$ L the stock suspensions of each test-organism were diluted in a series of 1:10 dilutions resulting in a final dilution of  $10^7$  (Gram-positive bacteria) and  $10^6$  (yeasts). 10  $\mu$ L from the final diluted test-suspensions were transferred to TSA isolation media plates and incubated for 18-24h at  $(37\pm 2)^\circ\text{C}$  for Gram-positive bacteria and at  $(29\pm 2)^\circ\text{C}$  for yeasts, respectively. To account for random errors, the determination of each suspension's concentration in cfu/ $\mu$ L containing microbial material was done in triplicate. All suspensions were tested in replicates of eight (8) via each testing methodology (DT, eDT, Ext). Study results concluded that the estimated dynamic range for the Direct, extended Direct and Extraction procedure are as follows:

Test Organism	DT		eDT		EXT	
	Lower limit [cfu/μL]	Upper limit [cfu/μL]	Lower limit [cfu/μL]	Upper limit [cfu/μL]	Lower limit [cfu/μL]	Upper limit [cfu/μL]
<i>Enterococcus faecalis</i>	1.2 x 10 <sup>6</sup>	6.0 x 10 <sup>7</sup>	3.6 x 10 <sup>6</sup>	1.8 x 10 <sup>8</sup>	3.6 x 10 <sup>6</sup>	1.8 x 10 <sup>8</sup>
<i>Enterococcus faecium</i>	4.5 x 10 <sup>7</sup>	4.5 x 10 <sup>7</sup>	2.1 x 10 <sup>6</sup>	1.1 x 10 <sup>8</sup>	2.1 x 10 <sup>6</sup>	1.1 x 10 <sup>8</sup>
<i>Staphylococcus aureus</i>	3.5 x 10 <sup>5</sup>	1.8 x 10 <sup>8</sup>	4.1 x 10 <sup>4</sup>	2.1 x 10 <sup>8</sup>	4.1 x 10 <sup>5</sup>	2.1 x 10 <sup>8</sup>
<i>Candida albicans</i>	N/A	N/A	2.0 x 10 <sup>5</sup>	2.0 x 10 <sup>6</sup>	2.0 x 10 <sup>6</sup>	1.0 x 10 <sup>7</sup>
<i>Candida parapsilosis</i>	N/A	N/A	2.5 x 10 <sup>5</sup>	2.5 x 10 <sup>6</sup>	2.5 x 10 <sup>6</sup>	1.3 x 10 <sup>7</sup>
<i>Saccharomyces cerevisiae</i>	N/A	N/A	1.5 x 10 <sup>5</sup>	1.5 x 10 <sup>6</sup>	1.5 x 10 <sup>5</sup>	7.5 x 10 <sup>6</sup>

### Media and Colony Stability

With the inclusion of Gram-negative microaerophilic, Gram-negative anaerobic, Gram-positive aerobic and anaerobic and yeast organisms, a study on the following media was conducted to confirm acceptability of the recommended agar and the stability of the colony for up to 12 hours at room temperature after initial plate incubation prior to analysis.

- Chocolate Agar (CHOC)
- Columbia CNA agar with 5% sheep blood (CNA)
- Brucella Agar with 5% horse blood (BRU)
- CDC anaerobe Agar with 5% sheep blood (CDC)
- CDC anaerobe 5% sheep blood Agar with phenylethyl alcohol (CDC/PEA)
- CDC anaerobe laked sheep blood Agar with kanamycin and vancomycin (CDC/LKV)
- Bacteroides bile esculin Agar with amikacin (BBE)
- Clostridium difficile Agar with 7% sheep blood (CDA)
- Trypticase Soy Agar with 5% sheep blood (TSA)
- Sabouraud-Dextrose Agar (SDA)
- Brain-Heart Infusion Agar (BHI)
- Campylobacter Agar with 5 Antimicrobics and 10% sheep blood (CAMPY BAP)
- Bordet Gengou Agar with 15% sheep blood (BGA)
- Columbia Blood Agar with 5% sheep blood (CBA)

Testing was conducted using three (3) Gram-positive bacteria, three (3) yeasts, five (5) anaerobic bacteria, two (2) *Campylobacter* and three (3) *Bordetella* species at varying incubation time points in replicates of eight (8). After initial testing, isolates were further tested at room temperature after twelve (12) hours post-incubation.

The study results confirmed the acceptability of all culture media tested with the following parameters:

- Bacteria and yeasts growth should be between 18h to 48h (+12h storage at room temperature (RT)).

Specific requirements:

- *Bordetella*: Incubation on BG agar should not be longer than 24h (+12h storage at RT).
- *Campylobacter*: Incubation on CAMPY BAP can be prolonged to 72h (+12h storage at RT).
- *Streptococcus pneumoniae*: Incubation should not be longer than 24h (+12h storage at RT) due to possible autolysis.

#### Organism Stability Prior to MBT-CA Analysis

This study was conducted to assess Gram-positive and yeast isolate stability on the target plate prior to matrix overlay via Direct Transfer (DT), extended Direct transfer (eDT) and Extraction (Ext) procedure. Two (2) Gram-positive organisms were inoculated eight times and overlaid with matrix at five (5) different time points. Five (5) yeasts organisms were inoculated eight times and overlaid with matrix at (4) different time points. Extracts of Gram-positive bacteria and yeasts were stored at room temperature and inoculated eight times and overlaid with matrix at five (5) different time points. All testing was performed in duplicate.

The study results confirmed that Gram-positive bacteria and yeast organisms are stable on the target plate for up to 60 minutes and 30 minutes respectively prior to matrix addition. Extracts of Gram-positive bacteria and yeasts are stable at room temperature for up to 24 hours and 4 hours respectively.

#### Sample Stability Overlaid with Matrix

This study was conducted to assess test organism stability overlaid with matrix after inoculation on the target plate. For this study, six (6) organisms were tested (three (3) Gram-positive and three (3) yeast organisms). All organisms were subcultured and aging experiments were conducted at two (2) temperatures and two (2) different relative humidity conditions to stress the plate. Plates were inoculated with the test organism via DT, eDT and Ext procedure. Plates were then read immediately (0h) and then incubated at each test condition and analyzed at three (3) additional time points (4 hours, 8 hours and 24 hours).

The study results confirmed that organisms overlaid with matrix on the target plate are stable for up to 24 hours when stored at room temperature.

Reproducibility:

The reproducibility study for Gram-positive aerobic bacteria, Gram-negative microaerophilic bacteria, Gram-positive anaerobic bacteria, Gram-negative anaerobic bacteria and yeast organisms was carried out to confirm day-to-day reproducibility and precision of the MALDI Biotyper CA System at different clinical study sites. The study was conducted for five (5) days with two (2) runs (two (2) operators) each day per clinical site. The sources of variability tested were:

- \* Two (2) operators/each clinical study site
- \* Three (3) clinical study sites
- \* At least four (4) target plates/each clinical study site
- \* Four (4) microflex LT/SH instruments

Ten (10) well-characterized organisms were chosen for this study and tested in duplicate via Direct Transfer and extended Direct Transfer procedure in accordance with product instructions. When the DT and/or eDT log(score) was <2.00, per product instructions, the test organism was tested following Extraction procedure in duplicate.

**Table 3: Overall Reproducibility Panel Testing per Test Organism using  $\geq 2.0$  MBT-CA log(scores)**

Blinded Test Organism	Reproducibility Panel	$\geq 2.0$ ID (DT)	$\geq 2.0$ ID (eDT)	$\geq 2.0$ ID (DT+eDT+Ext)
<i>Enterococcus faecalis</i>	REPRO-1	60/60 (100%)	60/60 (100%)	<b>60/60 (100%)</b>
<i>Staphylococcus epidermidis</i>	REPRO-2	58/60 (97%)	59/60 (98%)	<b>60/60 (100%)</b>
<i>Streptococcus agalactiae</i>	REPRO-3	58/60 (97%)	55/60 (92%)	<b>60/60 (100%)</b>
<i>Bacteroides fragilis</i>	REPRO-4	60/60 (100%)	60/60 (100%)	<b>60/60 (100%)</b>
<i>Fusobacterium necrophorum</i>	REPRO-5	56/60 (93%)	53/60 (88%)	<b>58/60 (97%)</b>
<i>Clostridium perfringens</i>	REPRO-6	53/60 (88%)	58/60 (97%)	<b>59/60 (98%)</b>
<i>Propionibacterium acnes</i>	REPRO-7	53/60 (88%)	49/60 (82%)	<b>60/60 (100%)</b>
<i>Candida albicans</i>	REPRO-8	33/60 (55%)	41/60 (68%)	<b>60/60 (100%)</b>
<i>Saccharomyces cerevisiae</i>	REPRO-9	5/60 (8%)	0/60 (0%)	<b>31/60 (52%)</b>
<i>Cryptococcus neoformans var grubii</i>	REPRO-10	31/60 (52%)	44/60 (73%)	<b>53/60 (88%)</b>
<b>TOTAL</b>		<b>467/600 (78%)</b>	<b>479/600 (80%)</b>	<b>561/600 (94%)</b>

94% of test organisms were correctly identified with a log(score)  $\geq 2.00$  result. In addition, no isolates were falsely identified. Thus, data confirm reproducibility and precision of the whole MALDI Biotyper CA System independent from:

- Clinical Site
- System operator
- microflex LT/SH instrument
- Target plate

Challenge Panel:

A panel of 55 organisms (24 Gram-positive aerobic bacteria, 1 Gram-negative microaerophilic bacterium, 4 Gram-negative anaerobic bacteria, 6 Gram-positive anaerobic bacteria, 20 yeasts) was tested at four (4) study sites. Fifty-three (53) of the organisms included in the panel were selected from stored organisms tested during the clinical study. Two (2) were selected from strain collections. The study reference laboratory, prepared the panel. Organism identifications were blinded to test sites. Each site tested the challenge panel member via Direct Transfer and extended Direct Transfer procedure in accordance with product instructions. If DT and/or eDT result yielded a log(score)  $< 2.00$ , the organism was retested using the Extraction procedure.

**Table 4: Challenge Panel Study Summary**

Test procedure	Site A $\geq 2.0$ ID	Site B $\geq 2.0$ ID	Site C $\geq 2.0$ ID	Site D $\geq 2.0$ ID	TOTAL $\geq 2.0$ ID
DT method	46/55 (84%)	47/55 (85%)	35/55 (64%)*	36/55 (65%)	164/220 (75%)
eDT method	47/55 (85%)	50/55 (91%)	44/55 (80%)*	43/55 (78%)	184/220 (84%)
Ext method	54/55 (98%)	53/55 (96%)	49/55 (89%)	37/55 (67%)	193/220 (88%)
<b>MBT-CA workflow</b>	<b>54/55 (98%)</b>	<b>55/55 (100%)</b>	<b>49/55 (89%)</b>	<b>46/55 (84%)</b>	<b>204/220 (93%)</b>

\* One sample was incorrectly identified due to a mixed culture.

93% of test organisms were correctly identified with a log(score)  $\geq 2.00$  result applying MBT-CA workflow. Testing of the challenge panel confirms intra laboratory performance of the MALDI Biotyper CA System.

Method Comparison:

To demonstrate performance of the MALDI Biotyper CA (MBT-CA) System, a method comparison study was performed at six (6) US clinical test sites and in-house laboratory. 4,395 (generating 4,399 data points) fresh and stored organisms were tested on the MALDI Biotyper CA System in accordance to manufacturer’s instructions for use. All organisms included in the study were sub-cultured for purity. Testing on the MBT-CA System was done from a fresh isolated colony. Due to the rarity of some organisms, replicates of these rarer species were tested by multiple testing sites to generate additional data to support performance of the MBT-CA System. Results from the 3,802 replicate testing results were analyzed separately from the method comparison isolates.

All organisms included in the study were also sub-cultured on to an agar slant or appropriate media for isolation and shipped to the study interim reference laboratory. The interim reference laboratory stored all organisms included in the study and sent all organisms to the sequencing reference laboratory for sequencing and protein sequencing when requested.

The following Gram-negative, Gram-positive and yeast isolates are included in the reference library (please refer to K130831 for previously claimed Gram-negative organisms)

Table 5: Claimed Organisms

Organisms	Organisms
<i>Acinetobacter haemolyticus</i>	<i>Candida guilliermondii</i>
<i>Acinetobacter johnsonii</i>	<i>Candida haemulonis</i>
<i>Acinetobacter junii</i>	<i>Candida inconspicua</i>
<i>Actinomyces meyeri</i>	<i>Candida kefyr</i>
<i>Actinomyces neuii</i>	<i>Candida krusei</i>
<i>Actinomyces odontolyticus</i>	<i>Candida lambica</i>
<i>Actinomyces oris</i>	<i>Candida lipolytica</i>
<i>Aerococcus urinae</i>	<i>Candida lusitaniae</i>
<i>Aerococcus viridans</i>	<i>Candida metapsilosis</i>
<i>Aeromonas salmonicida</i>	<i>Candida norvegensis</i>
<i>Anaerococcus vaginalis</i>	<i>Candida orthopsilosis</i>
<i>Bacteroides fragilis</i>	<i>Candida parapsilosis</i>
<i>Bacteroides ovatus</i> group	<i>Candida pararugosa</i>
<i>Bacteroides thetaiotaomicron</i> group	<i>Candida pelliculosa</i>
<i>Bacteroides uniformis</i>	<i>Candida tropicalis</i>
<i>Bacteroides vulgatus</i> group	<i>Candida valida</i>
<i>Bordetella</i> group[3]	<i>Capnocytophaga ochracea</i>
<i>Bordetella hinzii</i>	<i>Capnocytophaga sputigena</i>
<i>Brevibacterium casei</i>	<i>Chryseobacterium gleum</i>
<i>Brevundimonas diminuta</i> group	<i>Chryseobacterium indologenes</i>
<i>Campylobacter coli</i>	<i>Clostridium difficile</i>
<i>Campylobacter jejuni</i>	<i>Clostridium perfringens</i>
<i>Campylobacter ureolyticus</i>	<i>Corynebacterium amycolatum</i>
<i>Candida albicans</i>	<i>Corynebacterium aurimucosum</i> group
<i>Candida boidinii</i>	<i>Corynebacterium bovis</i>
<i>Candida dubliniensis</i>	<i>Corynebacterium diphtheriae</i>
<i>Candida duobushaemulonii</i>	<i>Corynebacterium glucuronolyticum</i>
<i>Candida famata</i>	<i>Corynebacterium jeikeium</i>
<i>Candida glabrata</i>	<i>Corynebacterium kroppenstedtii</i>

Organisms
<i>Corynebacterium macginleyi</i>
<i>Corynebacterium minutissimum</i>
<i>Corynebacterium propinquum</i>
<i>Corynebacterium pseudodiphtheriticum</i>
<i>Corynebacterium riegellii</i>
<i>Corynebacterium striatum</i> group
<i>Corynebacterium tuberculostearicum</i>
<i>Corynebacterium ulcerans</i>
<i>Corynebacterium urealyticum</i>
<i>Corynebacterium xerosis</i>
<i>Cronobacter sakazakii</i> group
<i>Cryptococcus gattii</i>
<i>Cryptococcus neoformans</i> var <i>grubii</i>
<i>Cryptococcus neoformans</i> var <i>neoformans</i>
<i>Cupriavidus pauculus</i> group
<i>Delftia acidovorans</i> group
<i>Dermacoccus nishinomiyaensis</i>
<i>Edwardsiella tarda</i>
<i>Elizabethkingia meningoseptica</i> group
<i>Enterobacter amnigenus</i>
<i>Enterococcus avium</i> group
<i>Enterococcus casseliflavus</i>
<i>Enterococcus faecalis</i>
<i>Enterococcus faecium</i>
<i>Enterococcus gallinarum</i>
<i>Enterococcus hirae</i>
<i>Fingoldia magna</i>
<i>Fusobacterium canifelinum</i>
<i>Fusobacterium necrophorum</i>
<i>Fusobacterium nucleatum</i>
<i>Gardnerella vaginalis</i>
<i>Gemella haemolysans</i>
<i>Gemella sanguinis</i>
<i>Geotrichum candidum</i>
<i>Geotrichum capitatum</i>
<i>Granulicatella adiacens</i>
<i>Haemophilus haemolyticus</i>
<i>Haemophilus influenzae</i>
<i>Haemophilus parahaemolyticus</i> group
<i>Kingella kingae</i>
<i>KloECKera apiculata</i>
<i>Kocuria kristinae</i>

Organisms
<i>Kytococcus sedentarius</i>
<i>Lactococcus garvieae</i>
<i>Lactococcus lactis</i>
<i>Leuconostoc mesenteroides</i>
<i>Macrococcus caseolyticus</i>
<i>Micrococcus luteus</i>
<i>Moraxella</i> sg <i>Moraxella nonliquefaciens</i>
<i>Myroides odoratimimus</i>
<i>Myroides odoratus</i>
<i>Oligella ureolytica</i>
<i>Oligella urethralis</i>
<i>Parabacteroides distasonis</i>
<i>Pediococcus pentosaceus</i>
<i>Peptoniphilus harei</i> group
<i>Peptostreptococcus anaerobius</i>
<i>Pichia ohmeri</i>
<i>Plesiomonas shigelloides</i>
<i>Porphyromonas gingivalis</i>
<i>Prevotella bivia</i>
<i>Prevotella buccae</i>
<i>Prevotella denticola</i>
<i>Prevotella intermedia</i>
<i>Prevotella melaninogenica</i>
<i>Propionibacterium acnes</i>
<i>Pseudomonas oryzihabitans</i>
<i>Pseudomonas stutzeri</i>
<i>Rhizobium radiobacter</i>
<i>Rothia aeria</i>
<i>Rothia dentocariosa</i>
<i>Rothia mucilaginoso</i>
<i>Saccharomyces cerevisiae</i>
<i>Serratia plymuthica</i>
<i>Serratia rubidaea</i>
<i>Staphylococcus aureus</i>
<i>Staphylococcus auricularis</i>
<i>Staphylococcus capitis</i>
<i>Staphylococcus caprae</i>
<i>Staphylococcus carnosus</i>
<i>Staphylococcus cohnii</i>
<i>Staphylococcus epidermidis</i>
<i>Staphylococcus equorum</i>
<i>Staphylococcus felis</i>

Organisms
<i>Staphylococcus haemolyticus</i>
<i>Staphylococcus hominis</i>
<i>Staphylococcus lugdunensis</i>
<i>Staphylococcus pasteurii</i>
<i>Staphylococcus pettenkoferi</i>
<i>Staphylococcus pseudintermedius</i>
<i>Staphylococcus saccharolyticus</i>
<i>Staphylococcus saprophyticus</i>
<i>Staphylococcus schleiferi</i>
<i>Staphylococcus simulans</i>
<i>Staphylococcus vitulinus</i>
<i>Staphylococcus warneri</i>
<i>Streptococcus agalactiae</i>
<i>Streptococcus anginosus</i>
<i>Streptococcus constellatus</i>
<i>Streptococcus dysgalactiae</i>
<i>Streptococcus gallolyticus</i>

Organisms
<i>Streptococcus gordonii</i>
<i>Streptococcus intermedius</i>
<i>Streptococcus lutetiensis</i>
<i>Streptococcus mitis / oralis</i> group
<i>Streptococcus mutans</i>
<i>Streptococcus pneumoniae</i>
<i>Streptococcus pyogenes</i>
<i>Streptococcus salivarius</i>
<i>Sutterella wadsworthensis</i>
<i>Trichosporon asahii</i>
<i>Vibrio parahaemolyticus</i>
<i>Vibrio vulnificus</i>

Tables 6 - 8 below show the overall isolate performance.

Table 6: Overall Isolate Performance - claim 2

Overall Performance - claim 2				
MBT-CA RESULT	REFERENCE ALGORITHM			Total
	high resolution species	low resolution species / genus	Negative	
Organism ID $\geq 2.0$ (High Confidence)	3817	392	18	4227
Organism ID ( $\geq 1.7$ ; $< 2.0$ ) (Low Confidence)	107	13	9	129
- INCORRECT MBT-CA ID ( $\geq 1.7$ ) - NO ID ( $< 1.7$ )	42	1	n/a	43
Total	3966	406	27	4399

Positive		Negative
high resolution	high & low resolution	
96.27%	99.02%	n/a

Table 7: Overall Bacteria Performance

Overall Performance BACTERIA				
MBT-CA RESULT	REFERENCE ALGORITHM			
	high resolution species	low resolution species / genus	Negative	Total
Organism ID $\geq 2.0$ (High Confidence)	3079	389	17	3485
Organism ID ( $\geq 1.7$ ; $< 2.0$ ) (Low Confidence)	52	12	9	73
- INCORRECT MBT-CA ID ( $\geq 1.7$ ) - NO ID ( $< 1.7$ )	25	1	n/a	26
Total	3156	402	26	3584

Positive		Negative
high resolution	high & low resolution	
97.47%	99.27%	n/a

Table 8: Overall Yeast Performance

Overall Performance YEAST				
MBT-CA RESULT	REFERENCE ALGORITHM			
	high resolution species	low resolution species / genus	Negative	Total
Organism ID $\geq 2.0$ (High Confidence)	738	3	1	742
Organism ID ( $\geq 1.7$ ; $< 2.0$ ) (Low Confidence)	55	1	0	56
- INCORRECT MBT-CA ID ( $\geq 1.7$ ) - NO ID ( $< 1.7$ )	17	0	0	17
Total	810	4	1	815

Positive		Negative
high resolution	high & low resolution	
91.03%	97.91%	0.00%

### Statement of Safety and Efficacy

The data presented clearly demonstrate the safety and efficacy of the Bruker Daltonics, Inc. MBT-CA System as compared to the reference algorithm, when the instructions for use are followed.