



June 28, 2019

iCubate, Inc.
% Fran White
President
MDC Associates, LLC
180 Cabot Street
Beverly, Massachusetts 01915

Re: K190341

Trade/Device Name: iC-GN iC-Cassette for use on the iC-System

Regulation Number: 21 CFR 866.3365

Regulation Name: Multiplex nucleic acid assay for identification of microorganisms and resistance markers from positive blood cultures

Regulatory Class: Class II

Product Code: PEN

Dated: February 11, 2019

Received: February 14, 2019

Dear Fran White:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for

Uwe Scherf, M.Sc., Ph.D.
Director
Division of Microbiology Devices
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

510(k) SUMMARY

Date of Summary: June 22, 2019

Product Name: iC-GN Assay™ for use on the iC-System™

Sponsor:

iCubate, Inc.
601 Genome Way
Huntsville, AL 35806

Correspondent:

MDC Associates, Inc.
Fran White, President
180 Cabot Street
Beverly, MA 01915
Phone: (978) 705 5011
Fax: (866) 540 3448
Email: regulatory@mdcassoc.com

Common Name:

Gram-Negative Bacteria and Associated Resistance Markers

Regulation Number:

866.3365

Classification:

PEN, Class II

Substantial Equivalency

Characteristic	iCubate, Inc. iC-GN Assay™ for use on the iC-System™ (New Device)	Nanosphere, Inc. Verigene® Gram Negative Blood Culture Nucleic Acid Test (GC-GN) K132843 (Predicate Device)																		
<i>Similarities</i>																				
<p>Intended Use</p>	<p>The iCubate, Inc. iC-GN Assay™ for use on the iC-System™ is a qualitative, multiplexed, <i>in vitro</i> diagnostic test for the detection and identification of potentially pathogenic gram negative bacteria, which may cause bloodstream infection (BSI). The iC-GN Assay™ is performed directly on positive blood cultures, confirmed by Gram stain to contain gram negative bacilli. Cultures demonstrating mixed Gram stain results should not be tested on the assay. The iC-GN Assay™ is validated for use with select <i>BACTEC</i>™, <i>Bact/ALERT</i>® and <i>VersaTREK</i>® blood culture bottles. The iC-GN Assay™ is indicated for use in conjunction with other clinical and laboratory findings, such as culture, to aid in the diagnosis of bacterial bloodstream infections; however, it is not used to monitor bloodstream infections.</p> <p>The iC-GN Assay™ detects target DNA and identifies the following:</p> <table border="1" data-bbox="521 1171 964 1665"> <thead> <tr> <th data-bbox="521 1171 732 1234">Bacterial Genera and Species</th> <th data-bbox="732 1171 964 1234">Resistance Markers</th> </tr> </thead> <tbody> <tr> <td data-bbox="521 1234 732 1329"><i>Acinetobacter baumannii</i> complex</td> <td data-bbox="732 1234 964 1329">KPC (bla_{KPC})-associated with resistance to carbapenems</td> </tr> <tr> <td data-bbox="521 1329 732 1392"><i>Enterobacter cloacae</i> complex</td> <td data-bbox="732 1329 964 1392">NDM (bla_{NDM})-associated with resistance to carbapenems</td> </tr> <tr> <td data-bbox="521 1392 732 1413"><i>Escherichia coli</i></td> <td data-bbox="732 1392 964 1413"></td> </tr> <tr> <td data-bbox="521 1413 732 1434"><i>Klebsiella oxytoca</i></td> <td data-bbox="732 1413 964 1434"></td> </tr> <tr> <td data-bbox="521 1434 732 1455"><i>Klebsiella pneumoniae</i></td> <td data-bbox="732 1434 964 1455">CTX-M group</td> </tr> <tr> <td data-bbox="521 1455 732 1476"><i>Pseudomonas aeruginosa</i></td> <td data-bbox="732 1455 964 1476">1(bla_{CTX-M} group 1)-associated with</td> </tr> <tr> <td data-bbox="521 1476 732 1497"><i>Proteus</i> species</td> <td data-bbox="732 1476 964 1497">resistance to</td> </tr> <tr> <td data-bbox="521 1497 732 1518"><i>Serratia marcescens</i></td> <td data-bbox="732 1497 964 1518">extended spectrum beta-lactams</td> </tr> </tbody> </table> <p>In mixed growth, the iC-GN Assay™ does not specifically attribute detection of KPC, NDM, or CTX-M group 1 to a specific genera or species.</p> <p>Sub-culturing of positive blood cultures is necessary to recover organisms for</p>	Bacterial Genera and Species	Resistance Markers	<i>Acinetobacter baumannii</i> complex	KPC (bla _{KPC})-associated with resistance to carbapenems	<i>Enterobacter cloacae</i> complex	NDM (bla _{NDM})-associated with resistance to carbapenems	<i>Escherichia coli</i>		<i>Klebsiella oxytoca</i>		<i>Klebsiella pneumoniae</i>	CTX-M group	<i>Pseudomonas aeruginosa</i>	1(bla _{CTX-M} group 1)-associated with	<i>Proteus</i> species	resistance to	<i>Serratia marcescens</i>	extended spectrum beta-lactams	<p>The Verigene® Gram Negative Blood Culture Nucleic Acid Test (BC-GN), performed using the sample-to-results Verigene System, is a qualitative multiplexed <i>in vitro</i> diagnostic test for the simultaneous detection and identification of selected gram-negative bacteria and resistance markers. BC-GN is performed directly on blood culture media using blood culture bottles identified as positive by a continuous monitor blood culture system and which contain gram-negative bacteria as determined by Gram stain. BC-GN detects and identifies the following:</p> <p><i>Acinetobacter</i> spp. <i>Citrobacter</i> spp. <i>Enterobacter</i> spp. <i>Proteus</i> spp. <i>Escherichia coli</i> <i>Klebsiella pneumoniae</i> <i>Klebsiella oxytoca</i> <i>Pseudomonas aeruginosa</i></p> <p>BC-GN is indicated for use in conjunction with other clinical and laboratory findings to aid in the diagnosis of bacterial bloodstream infections; however, is not to be used to monitor these infections. Sub-culturing of positive blood cultures is necessary to recover organisms for susceptibility testing, identification of organisms not detected by BC-GN, differentiation of mixed growth, association of antimicrobial resistance marker genes to a specific organism, or for epidemiological typing.</p>
Bacterial Genera and Species	Resistance Markers																			
<i>Acinetobacter baumannii</i> complex	KPC (bla _{KPC})-associated with resistance to carbapenems																			
<i>Enterobacter cloacae</i> complex	NDM (bla _{NDM})-associated with resistance to carbapenems																			
<i>Escherichia coli</i>																				
<i>Klebsiella oxytoca</i>																				
<i>Klebsiella pneumoniae</i>	CTX-M group																			
<i>Pseudomonas aeruginosa</i>	1(bla _{CTX-M} group 1)-associated with																			
<i>Proteus</i> species	resistance to																			
<i>Serratia marcescens</i>	extended spectrum beta-lactams																			

Characteristic	iCubate, Inc. iC-GN Assay™ for use on the iC-System™ (New Device)	Nanosphere, Inc. Verigene® Gram Negative Blood Culture Nucleic Acid Test (GC-GN) K132843 (Predicate Device)
	susceptibility testing, identification of organisms not detected by the iC-GN Assay™, differentiation of mixed growth, association of antimicrobial resistance marker genes to a specific organism, or for epidemiological typing.	
Sample Type	Positive Blood Culture	Positive Blood Culture
<i>Differences</i>		
INSTRUMENT REQUIREMENTS	iC-System™	Verigene System
TEST PRINCIPLE	ARM-PCR	Gold nanoparticle probe-based PCR
COMPATIBLE BLOOD CULTURE BOTTLES	BD BACTEC Standard/10 Aerobic/F BD BACTEC Standard/10 Anaerobic/F BD BACTEC Plus Aerobic/F BD BACTEC Plus Anaerobic/F BD BACTEC Lytic/10 Anaerobic/F BacT/Alert SA Standard Aerobic BacT/Alert SN Standard Anaerobic BacT/Alert FA Aerobic FAN BacT/Alert FN Anaerobic FAN BacT/Alert FA Plus Aerobic BacT/Alert FN Plus Anaerobic VersaTREK REDOX 1 VersaTREK REDOX 2	BACTEC™ Plus Aerobic/F BacT/ALERT FA FAN
THROUGHPUT	Four (4) samples/iC-Processor™	One (1) Sample/Processor

Intended Use

The iCubate, Inc. iC-GN Assay™ for use on the iC-System™ is a qualitative, multiplexed, *in vitro* diagnostic test for the detection and identification of potentially pathogenic gram negative bacteria, which may cause bloodstream infection (BSI). The iC-GN Assay™ is performed directly on positive blood cultures, confirmed by Gram stain to contain gram negative bacilli. Cultures demonstrating mixed Gram stain results should not be tested on the assay. The iC-GN Assay™ is validated for use with select *BACTEC*™, *BacT/ALERT*® and *VersaTREK*® blood culture bottles. The iC-GN Assay™ is indicated for use in conjunction with other clinical and laboratory findings, such as culture, to aid in the diagnosis of bacterial bloodstream infections; however, it is not used to monitor bloodstream infections.

The iC-GN Assay™ detects target DNA and identifies the following:

Bacterial Genera and Species	Resistance Markers
<i>Acinetobacter baumannii</i> complex	KPC (<i>bla_{KPC}</i>)- associated with resistance to carbapenems NDM (<i>bla_{NDM}</i>)- associated with resistance to carbapenems CTX-M group 1(<i>bla_{CTX-M}</i> group 1)- associated with resistance to extended spectrum beta-lactams
<i>Enterobacter cloacae</i> complex	
<i>Escherichia coli</i>	
<i>Klebsiella oxytoca</i>	
<i>Klebsiella pneumoniae</i>	
<i>Pseudomonas aeruginosa</i>	
<i>Proteus</i> species	
<i>Serratia marcescens</i>	

In mixed growth, the iC-GN Assay™ does not specifically attribute detection of KPC, NDM, or CTX-M group 1 to a specific genera or species.

Sub-culturing of positive blood cultures is necessary to recover organisms for susceptibility testing, identification of organisms not detected by the iC-GN Assay™, differentiation of mixed growth, association of antimicrobial resistance marker genes to a specific organism, or for epidemiological typing.

Limitations

For prescription use only.

Please refer to the iC-GN Assay™ labeling for a more complete list of warnings, precautions and contraindications.

Methodology

The iC-GN Assay™ utilizes polymerase chain reaction (PCR) for the multiplex amplification of specific targets and detects the amplified targets with microarray hybridization. Targets are detected directly from patient positive blood cultures confirmed by Gram stain to contain gram negative bacilli. The iC-GN Assay utilizes proprietary ARM-PCR (Amplicon Rescued Multiplex PCR) technology allowing for multiple targets to be amplified in one reaction. Testing is done in a self-contained, automated, disposable cassette using the iCubate™ processor (iC-Processor™). After the reaction is complete, the cassette is read on the

iCubate® reader (iC-Reader™). Results from the iC-Reader™ are interpreted by iC-Report™ software and a final report is displayed on the iMac® computer.

To operate, the user opens the iC-Cassette™ cap and pipettes an aliquot of the diluted positive blood culture sample into the sample/PCR well in the bottom well plate of the cassette. Once inoculated, the cassette cap is closed, and all extraction, amplification and detection processes are completed in the cassette, a closed system. Extraction, amplification and detection sequences are defined by an assay script controlled by the iC-Processor™.

The processing script is defined within a barcode label positioned on the top of each iC-Cassette™ which communicates with the iC-Processor™. To access and pierce the foil-sealed reagent wells located in the bottom well plate of the cassette, the processor manipulates the cassette to move the cassette pipette horizontally and vertically. The script directs the transfer of reagents between the wells in the bottom well plate and finally to the array within the cassette. The iC-Processor™ is capable of processing four (4) iC-Cassettes™ with random access.

Once processing is complete, the cassette is manually transferred from the iC-Processor™ to the iC-Reader™ where the microarray within the cassette is read. The iC-Reader™ is capable of reading up to four (4) iC-Cassettes™ at one time. The results are interpreted via the iC-Report™ software and displayed for the user on the iMac®. Raw data and result interpretations are stored within the iMac®; raw data is accessible to iCubate® service personnel only and not to the end user.

When finished with a loaded iC-GN Cassette™, it should be disposed as biohazardous waste.

Performance Data

For ease of reference, the following table defines iC-GN target organisms and common acronyms used in the study descriptions.

TABLE 1: iC-GN Assay Targets	
Target	Acronym
<i>Acinetobacter baumannii</i> complex	ABX
<i>Enterobacter cloacae</i> complex	ECX
<i>Escherichia coli</i>	EC
<i>Klebsiella oxytoca</i>	KO
<i>Klebsiella pneumoniae</i>	KPN
<i>Proteus mirabilis</i>	PM
<i>Pseudomonas aeruginosa</i>	PA
<i>Serratia marcescens</i>	SM
KPC carbapenemase resistance marker	KPC
NDM carbapenemase resistance marker	NDM
CTX-M group 1 extended spectrum β-lactamase resistance marker	CTXM

Reproducibility

To confirm the site-to-site, operator-to-operator, system-to-system, and lot-to-lot reproducibility of the iC-GN Assay, a representative panel of target organisms and one non-target organism were evaluated at two clinically relevant concentrations: initial bottle positivity and eight hours beyond initial bottle positivity. Organisms were grown to the appropriate concentrations in BD BACTEC Plus Aerobic blood culture bottles with human blood added on the BD BACTEC System. Testing was performed by two independent operators at each of three sites, two external and one internal. Each operator tested the eighteen-organism panel in triplicate across five, non-consecutive days. Testing was performed on six iC-GN Cassette lots and multiple iC-Systems. Performance is based on all expected targets detected and no false positive targets detected. Table 2 below summarizes Reproducibility results stratified by iC-GN target and concentration. Overall Reproducibility performance was 99.3%, confirming that iC-GN Assay performance is reproducible across sites, operators, systems and lots.

Target/Concentration	Overall Performance	Overall Performance % [95% CI]	False Negatives	False Positives	PC Check Failures	System Failures
A. baumannii complex Bottle Ring	90/90	100.0 [95.91-100.0]	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
A. baumannii complex Bottle Ring + 8 hours	87/90	96.7 [90.65-98.86]	0/90 (0.00%)	3/90 (3.33%)	0/90 (0.00%)	0/90 (0.00%)
E. cloacae complex Bottle Ring	86/88	97.7 [92.09-99.37]	1/88 (1.14%)	1/88 (1.14%)	2/90 (2.22%)	0/90 (0.00%)
E. cloacae complex Bottle Ring + 8 hours	90/90	100.0 [95.91-100.0]	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
E. coli Bottle Ring	90/90	100.0 [95.91-100.0]	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
E. coli Bottle Ring + 8 hours	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	0/90 (0.00%)	1/90 (1.11%)
K. oxytoca Bottle Ring	89/90	98.9 [93.97-99.80]	0/90 (0.00%)	1/90 (1.11%)	0/90 (0.00%)	0/90 (0.00%)
K. oxytoca Bottle Ring + 8 hours	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	1/90 (1.11%)	0/90 (0.00%)
K. pneumoniae Bottle Ring	90/90	100.0 [95.91-100.0]	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
K. pneumoniae Bottle Ring + 8 hours	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	1/90 (1.11%)	0/90 (0.00%)
Proteus species Bottle Ring	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	1/90 (1.11%)	0/90 (0.00%)
Proteus species Bottle Ring + 8 hours	88/88	100.0 [95.92-100.0]	0/88 (0.00%)	0/88 (0.00%)	0/90 (0.00%)	2/90 (2.22%)

TABLE 2: iC-GN Assay Reproducibility Performance by Target						
Target/Concentration	Overall Performance	Overall Performance % [95% CI]	False Negatives	False Positives	PC Check Failures	System Failures
<i>P. aeruginosa</i> Bottle Ring	88/89	98.9 [93.91-99.80]	1/89 (1.12%)	0/89 (0.00%)	1/90 (1.11%)	0/90 (0.00%)
<i>P. aeruginosa</i> Bottle Ring + 8 hours	89/90	98.9 [93.97-99.80]	1/90 (1.11%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
<i>S. marcescens</i> Bottle Ring	87/89	97.8 [92.17-99.38]	0/89 (0.00%)	2/89 (2.25%)	1/90 (1.11%)	0/90 (0.00%)
<i>S. marcescens</i> Bottle Ring + 8 hours	87/89	97.8 [92.17-99.38]	0/89 (0.00%)	2/89 (2.25%)	1/90 (1.11%)	0/90 (0.00%)
CTX-M group 1 Bottle Ring	90/90	100.0 [95.91-100.0]	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
CTX-M group 1 Bottle Ring + 8 hours	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	0/90 (0.00%)	1/90 (1.11%)
KPC Bottle Ring	90/90	100.0 [95.91-100.0]	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)
KPC Bottle Ring + 8 hours	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	1/90 (1.11%)	0/90 (0.00%)
NDM Bottle Ring	89/89	100.0 [95.86-100.0]	0/89 (0.00%)	0/89 (0.00%)	1/90 (1.11%)	0/90 (0.00%)
NDM Bottle Ring + 8 hours	89/90	98.9 [93.97-99.80]	1/90 (1.11%)	0/90 (0.00%)	0/90 (0.00%)	0/90 (0.00%)

Limit of Detection (LoD)

A study was performed to determine the limit of detection for each iC-GN Assay target, defined as the lowest concentration (CFU/mL) of analyte that can be detected approximately 95% of the time. For the eleven targets detected by the iC-GN Assay, a panel of twenty-seven representative strains were evaluated, a minimum of three per target. For complex and genus level targets, at least two representative species were evaluated. LoD testing was conducted in two phases, the first to narrow the range for LoD analysis. In phase II, the approximated 95% performance point determined in phase I was confirmed by testing a minimum of twenty replicates on each of three unique cassette lots. Plating and subsequent colony counts were used to determine organism concentrations. The final limit of detection for each target, provided in Table 3 below, was defined as the concentration that produced a positive result $\geq 95\%$ but $< 100\%$ of the time.

TABLE 3: iC-GN Assay LoD Results			
Target	Strain	Concentration (CFU/mL)	Defined Target LoD (CFU/mL)
<i>A. baumannii</i> complex	307-0294	5.3×10^5	$5.3 \times 10^5 - 5.2 \times 10^6$
	CDC-83	5.2×10^6	
	ATCC 23055	9.0×10^5	

TABLE 3: iC-GN Assay LoD Results			
Target	Strain	Concentration (CFU/mL)	Defined Target LoD (CFU/mL)
<i>E. cloacae</i> complex	Z101	5.0×10^6	$4.9 \times 10^5 - 5.5 \times 10^6$
	CDC-164	5.5×10^6	
	ATCC 700323	4.9×10^5	
<i>E. coli</i>	ATCC 43895	7.7×10^5	$7.7 \times 10^5 - 8.4 \times 10^5$
	ATCC BAA-2326	7.9×10^5	
	CDC-55	8.4×10^5	
<i>K. oxytoca</i>	Z115	6.2×10^5	$5.4 \times 10^5 - 1.1 \times 10^6$
	ATCC 13182	5.4×10^5	
	CDC-147	1.1×10^6	
<i>K. pneumoniae</i>	ATCC 35657	1.9×10^6	$6.0 \times 10^5 - 4.2 \times 10^6$
	CDC-40	3.6×10^6	
	CDC-42	1.9×10^6	
	KPC-2	4.2×10^6	
	LACNY 11	6.0×10^5	
<i>Proteus</i> species	Z050	1.1×10^6	$6.9 \times 10^5 - 1.1 \times 10^6$
	CDC-59	9.9×10^5	
	Z028	7.6×10^5	
	Z129	6.9×10^5	
<i>P. aeruginosa</i>	Z139	1.2×10^6	$5.0 \times 10^5 - 1.2 \times 10^6$
	CDC-231	5.0×10^5	
	CDC-250	6.9×10^5	
<i>S. marcescens</i>	ATCC 43297	7.2×10^5	$6.4 \times 10^5 - 8.1 \times 10^5$
	ATCC 21212	8.1×10^5	
	CDC-91	6.4×10^5	
CTX-M group 1	ATCC BAA-2326 (CTX-M-15)	7.9×10^5	$7.9 \times 10^5 - 2.3 \times 10^6$
	CDC-40 (CTX-M-15)	2.3×10^6	
	CDC-42 (CTX-M-15)	1.9×10^6	
KPC	CDC-147 (KPC-3)	2.3×10^6	$1.5 \times 10^5 - 4.2 \times 10^6$
	KPC-2	4.2×10^6	
	CDC-231 (KPC-5)	1.5×10^5	
NDM	CDC-83 (NDM-1)	5.2×10^6	$3.3 \times 10^5 - 5.2 \times 10^6$
	CDC-55 (NDM-1)	4.0×10^6	
	CDC-250 (NDM-1)	3.3×10^5	

Bottle Ring

A study was performed to establish the levels of each iC-GN target organism at two clinically relevant concentrations: initial bottle positivity (bottle “ring”) and eight hours beyond initial positivity. Twenty-seven representative organisms were evaluated, a minimum of three per iC-GN target. Organisms were grown in BD BACTEC Plus Aerobic blood culture bottles with human blood added on the BD BACTEC System. Within two hours of initial bottle positivity, the bottles were removed for plating and subsequent colony counts to determine organism concentrations. The bottles were then returned to the incubator and approximately eight hours after initial bottle positivity, the bottles were again removed for plating and subsequent colony counts to determine organism concentrations. Three bottles were grown for each

strain, and the average concentrations at initial bottle positivity and eight hours beyond initial bottle positivity are provided in Table 4 below. The concentrations at initial bottle positivity, representative of the lowest levels that may be observed in a clinical setting, are above the limits of detection determined for each strain.

TABLE 4: iC-GN Target Organism Concentrations at Bottle “Ring”			
Organism	Strain ID	Initial Bottle Positivity Average Concentration (CFU/mL)	Bottle Positivity + 8 Average Concentration (CFU/mL)
<i>Acinetobacter baumannii</i>	307-0294	4.24×10^8	8.27×10^8
<i>Acinetobacter baumannii</i>	CDC-83	3.39×10^8	7.23×10^8
<i>Acinetobacter calcoaceticus</i>	ATCC 23055	6.78×10^7	2.93×10^8
<i>Enterobacter cloacae</i>	Z101	2.17×10^8	1.97×10^9
<i>Enterobacter cloacae</i>	CDC-164	5.62×10^8	2.31×10^9
<i>Enterobacter hormaechei</i>	ATCC 700323	4.36×10^8	2.75×10^9
<i>Escherichia coli</i>	ATCC 43895	1.50×10^8	9.48×10^8
<i>Escherichia coli</i>	ATCC BAA-2326	6.23×10^8	1.52×10^9
<i>Escherichia coli</i>	CDC-55	4.93×10^8	1.51×10^9
<i>Klebsiella oxytoca</i>	Z115	5.32×10^8	2.07×10^9
<i>Klebsiella oxytoca</i>	ATCC 13182	4.16×10^8	4.52×10^9
<i>Klebsiella oxytoca</i>	CDC-147	9.67×10^8	1.31×10^9
<i>Klebsiella pneumoniae</i>	ATCC 35657	9.78×10^8	1.08×10^9
<i>Klebsiella pneumoniae</i>	CDC-40	2.16×10^8	1.36×10^9
<i>Klebsiella pneumoniae</i>	CDC-42	2.55×10^8	1.10×10^9
<i>Klebsiella pneumoniae</i>	KPC-2	7.70×10^8	1.66×10^9
<i>Klebsiella pneumoniae</i>	LACNY 11	5.43×10^7	1.67×10^9
<i>Proteus mirabilis</i>	Z050	1.71×10^8	7.40×10^8
<i>Proteus mirabilis</i>	CDC-59	7.37×10^7	8.10×10^8
<i>Proteus penneri</i>	Z028	8.88×10^7	4.33×10^8
<i>Proteus vulgaris</i>	Z129	4.37×10^7	5.00×10^8
<i>Pseudomonas aeruginosa</i>	Z139	9.18×10^7	1.37×10^{10}
<i>Pseudomonas aeruginosa</i>	CDC-231	3.26×10^8	7.98×10^8
<i>Pseudomonas aeruginosa</i>	CDC-250	1.64×10^8	8.97×10^8
<i>Serratia marcescens</i>	ATCC 43297	8.55×10^8	2.03×10^9
<i>Serratia marcescens</i>	ATCC 21212	1.07×10^8	8.83×10^8
<i>Serratia marcescens</i>	CDC-91	7.28×10^8	1.67×10^9

Blood Culture Bottle Equivalency

Commonly used blood culture bottle (BCB) media types were evaluated to demonstrate that variability in BCB media composition does not interfere with iC-GN Assay performance. Twenty-seven (27) representative iC-GN target organisms plus one non-target organism were tested in thirteen (13) BCB media types. Target organisms were tested near LoD

concentrations (2-3×LoD). Each strain was tested in triplicate in each BCB media type. Target performance is based on all expected targets detected and no false positive targets detected. Non-target performance is based on all expected negative results. In the event of a false negative result, the strain was retested in replicates of ten. In the event of a false positive result or other failure, the strain was retested in triplicate. The results of iC-GN BCB equivalency testing are summarized in Table 5 below. Performance in all bottle types met the acceptance criteria of ≥ 95% performance; all bottle types are validated for use with the iC-GN Assay.

TABLE 5: iC-GN Assay BCB Equivalency Results					
BCB Media Type	Overall Performance (%)	False Negatives (%)	False Positives (%)	PC Check Failures (%)	System Failures (%)
BACTEC Standard Aerobic	93/94 (98.9%)	1/94 (1.1%)	0/94 (0.0%)	3/97 (3.1%)	0/97 (0.0%)
BACTEC Standard Anaerobic	85/86 (98.8%)	0/86 (0.0%)	1/86 (1.2%)	0/87 (0.0%)	1/87 (1.1%)
BACTEC Plus Aerobic	93/94 (98.9%)	1/94 (1.1%)	0/94 (0.0%)	2/97 (2.1%)	1/97 (1.1%)
BACTEC Plus Anaerobic	95/96 (98.6%)	1/96 (1.0%)	0/96 (0.0%)	2/100 (2.0%)	2/100 (2.0%)
BACTEC Lytic/10 Anaerobic	81/81 (100.0%)	0/81 (0.0%)	0/81 (0.0%)	0/81 (0.0%)	0/81 (0.0%)
BACT/ALERT SA Standard Aerobic	97/99 (98.0%)	1/99 (1.0%)	1/99 (1.0%)	4/103 (3.9%)	0/103 (0.0%)
BACT/ALERT SN Standard Anaerobic	87/88 (98.9%)	0/88 (0.0%)	1/88 (1.1%)	2/90 (2.2%)	0/90 (0.0%)
BACT/ALERT FA Aerobic FAN	94/96 (97.9%)	0/96 (0.0%)	2/96 (2.1%)	1/97 (1.0%)	0/97 (0.0%)
BACT/ALERT FN Anaerobic FAN	92/94 (97.9%)	0/94 (0.0%)	2/94 (2.1%)	2/97 (2.1%)	1/97 (1.0%)
BACT/ALERT FA Plus Aerobic	94/95 (98.9%)	1/95 (1.1%)	0/95 (0.0%)	1/97 (1.1%)	1/97 (1.1%)
BACT/ALERT FN Plus Anaerobic	87/87 (100.0%)	0/87 (0.0%)	0/87 (0.0%)	2/90 (2.2%)	1/90 (1.1%)
VersaTREK REDOX 1	81/81 (100.0%)	0/81 (0.0%)	0/81 (0.0%)	0/81 (0.0%)	0/81 (0.0%)
VersaTREK REDOX 1	92/93 (98.9%)	1/93 (1.1%)	0/93 (0.0%)	1/94 (1.1%)	0/94 (0.0%)

An increased rate of false positive *Proteus* results was observed in some lots of BD BACTEC blood culture bottles. The high rate of false positive results observed prompted an

investigation by the manufacturer, BD Life Sciences. The false positives are due to the presence of nucleic acids or non-viable organisms present in the culture media at concentrations near or above the target’s limit of detection. While the observed contamination was resolved at the time of publication, positive *Proteus* results observed in BD BACTEC media types should be confirmed using alternative methods.

Inclusivity

To demonstrate the inclusivity of the iC-GN Assay, eighty-two (82) representative strains were evaluated, a minimum of ten strains for each target analyte. Strains were tested at the lowest level of bottle positivity, considered within two hours of bottle “ring.” Organisms were grown in BD BACTEC Plus Aerobic blood culture bottles with human blood added on the BD BACTEC System. Each strain was tested in triplicate. Performance is based on all expected targets detected and no false positive targets detected. In the event of a false negative result, the strain was retested in replicates of ten. In the event of a false positive result or other failure, the strain was retested in triplicate. The results of iC-GN Inclusivity testing are summarized in Table 6 below. Two strains were not detected by the iC-GN Assay: *Acinetobacter calcoaceticus* ATCC 31926 was not detected as *A. baumannii* complex and *Enterobacter kobei* ATCC BAA-260 was not detected as *E. cloacae* complex. An *in silico* analysis was also performed, and the predicted reactivity of each resistance marker detected by the iC-GN Assay is summarized in Tables 7-9 below.

TABLE 6: iC-GN Assay Inclusivity Results			
Organism	Strain	Targets	Performance
<i>Acinetobacter baumannii</i>	ATCC 19606	ABX	5/5
<i>Acinetobacter baumannii</i>	NCIMB 12457	ABX	3/3
<i>Acinetobacter baumannii</i>	CDC-36	ABX	3/3
<i>Acinetobacter baumannii</i>	CDC-37	ABX, NDM-1	3/3
<i>Acinetobacter baumannii</i>	CDC-45	ABX	3/3
<i>Acinetobacter baumannii</i>	CDC-52	ABX	3/3
<i>Acinetobacter baumannii</i>	CDC-56	ABX	3/3
<i>Acinetobacter baumannii</i>	CDC-88	ABX, NDM-1	3/3
<i>Acinetobacter baumannii</i>	CDC-101	ABX	3/3
<i>Acinetobacter calcoaceticus</i>	ATCC 14987	ABX	3/3
<i>Acinetobacter calcoaceticus</i>	ATCC 31926	ABX	2/11 ¹
<i>Enterobacter cloacae</i>	ATCC BAA-1143	ECX	3/3
<i>Enterobacter cloacae</i>	ATCC BAA-2341	ECX, KPC	3/3
<i>Enterobacter cloacae</i>	NCTC 10005	ECX	14/16 ²
<i>Enterobacter cloacae</i>	NCTC 13464	ECX	3/3
<i>Enterobacter cloacae</i>	CDC-32	ECX, KPC-3	3/3
<i>Enterobacter cloacae</i>	CDC-38	ECX, CTX-M-15, NDM-1	3/3
<i>Enterobacter cloacae</i>	CDC-65	ECX	3/3
<i>Enterobacter cloacae</i>	CDC-163	ECX, CTX-M-15, KPC-2	3/3
<i>Enterobacter asburiae</i>	ATCC 35923	ECX	3/3
<i>Enterobacter hormaechei</i>	ATCC 49162	ECX	3/3

TABLE 6: iC-GN Assay Inclusivity Results			
Organism	Strain	Targets	Performance
<i>Enterobacter kobei</i>	ATCC BAA-260	ECX	0/13 ³
<i>Escherichia coli</i>	ATCC 10536	EC	3/3
<i>Escherichia coli</i>	ATCC BAA-2469	EC, NDM-1	3/3
<i>Escherichia coli</i>	NCTC 9001	EC	3/3
<i>Escherichia coli</i>	NCTC 10538	EC	5/5
<i>Escherichia coli</i>	NCTC 13476	EC	3/3
<i>Escherichia coli</i>	CDC-48	EC, CTX-M-15, NDM-1	3/3
<i>Escherichia coli</i>	CDC-61	EC, KPC-3	3/3
<i>Escherichia coli</i>	CDC-104	EC, KPC-4	7/8 ⁴
<i>Escherichia coli</i>	CDC-119	EC, CTX-M-15, NDM-1	3/3
<i>Escherichia coli</i>	CDC-162	EC, CTX-M-15, NDM-7	3/3
<i>Klebsiella oxytoca</i>	ATCC 8724	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 43086	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 43165	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 43863	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 49134	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 49334	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 51817	KO	3/3
<i>Klebsiella oxytoca</i>	ATCC 700324	KO	3/3
<i>Klebsiella oxytoca</i>	NCTC 11686	KO	3/3
<i>Klebsiella oxytoca</i>	CDC-71	KO	3/3
<i>Klebsiella pneumoniae</i>	ATCC-13882	KPN	3/3
<i>Klebsiella pneumoniae</i>	ATCC BAA-1705	KPN, KPC-2	3/3
<i>Klebsiella pneumoniae</i>	NCTC 9633	KPN	3/3
<i>Klebsiella pneumoniae</i>	NCTC 13438	KPN, KPC-3	3/3
<i>Klebsiella pneumoniae</i>	NCTC 13443	KPN, CTX-M-15, NDM-1	3/3
<i>Klebsiella pneumoniae</i>	CDC-44	KPN, CTX-M-15	3/3
<i>Klebsiella pneumoniae</i>	CDC-46	KPN, CTX-M-15	5/5
<i>Klebsiella pneumoniae</i>	CDC-49	KPN, CTX-M-15, NDM-1	3/3
<i>Klebsiella pneumoniae</i>	CDC-66	KPN, CTX-M-15	3/3
<i>Klebsiella pneumoniae subsp. ozaenae</i>	ATCC 11296	KPN	3/3
<i>Proteus mirabilis</i>	ATCC 7002	Proteus	3/3
<i>Proteus mirabilis</i>	ATCC 21100	Proteus	3/3
<i>Proteus mirabilis</i>	ATCC 43071	Proteus	3/3
<i>Proteus mirabilis</i>	NCIMB 13283	Proteus	3/3
<i>Proteus mirabilis</i>	CDC-155	Proteus, KPC-6	3/3
<i>Proteus mirabilis</i>	CDC-156	Proteus, KPC-2	3/3
<i>Proteus mirabilis</i>	CDC-159	Proteus, NDM-1	3/3
<i>Proteus penneri</i>	ATCC 33519	Proteus	3/3
<i>Proteus vulgaris</i>	ATCC 9484	Proteus	3/3
<i>Proteus vulgaris</i>	ATCC 29905	Proteus	3/3
<i>Pseudomonas aeruginosa</i>	ATCC 10145	PA	3/3
<i>Pseudomonas aeruginosa</i>	ATCC 19429	PA	3/3

TABLE 6: iC-GN Assay Inclusivity Results			
Organism	Strain	Targets	Performance
<i>Pseudomonas aeruginosa</i>	ATCC BAA-1744	PA	3/3
<i>Pseudomonas aeruginosa</i>	CDC-54	PA	3/3
<i>Pseudomonas aeruginosa</i>	CDC-64	PA	3/3
<i>Pseudomonas aeruginosa</i>	CDC-90	PA, KPC-5	3/3
<i>Pseudomonas aeruginosa</i>	CDC-94	PA	3/3
<i>Pseudomonas aeruginosa</i>	CDC-105	PA	3/3
<i>Pseudomonas aeruginosa</i>	CDC-108	PA	3/3
<i>Pseudomonas aeruginosa</i>	CDC-246	PA, NDM-1	5/5
<i>Serratia marcescens</i>	ATCC 8100	SM	3/3
<i>Serratia marcescens</i>	ATCC 13880	SM	3/3
<i>Serratia marcescens</i>	ATCC 14041	SM	3/3
<i>Serratia marcescens</i>	ATCC 14756	SM	3/3
<i>Serratia marcescens</i>	ATCC 29634	SM	3/3
<i>Serratia marcescens</i>	ATCC 29635	SM	3/3
<i>Serratia marcescens</i>	ATCC 43861	SM	3/3
<i>Serratia marcescens</i>	ATCC 43862	SM	3/3
<i>Serratia marcescens</i>	NCTC 9743	SM	3/3
<i>Serratia marcescens</i>	CDC-99	SM	3/3

- 1) 2/2 false negative ABX in initial testing. 7/9 false negative ABX in repeat testing. See limitation.
- 2) 1/3 false positive ABX in initial testing. 1/3 false positive ABX in repeat testing. Strain repeated in replicates of 10, 10/10 repeats passed.
- 3) 3/3 false negative ECX in initial testing. 10/10 false negative ECX in repeat testing. See limitation.
- 4) 1/3 processor error in initial testing. 1/3 false positive KPN in repeat testing. Strain repeated in triplicate, 3/3 repeats passed.

TABLE 7: Predicted (<i>in silico</i>) Reactivity for CTX-M group 1			
Associated Target Organism	Variant Detected	Associated Target Organism	Variant Detected
<i>Acinetobacter baumannii</i> complex	CTX-M-3	<i>Klebsiella oxytoca</i>	CTX-M-3
	CTX-M-15		CTX-M-15
<i>Enterobacter cloacae</i> complex	CTX-M-1		CTX-M-35
	CTX-M-3		CTX-M-36
	CTX-M-15		CTX-M-162
	CTX-M-22		CTX-M-1
	CTX-M-37	CTX-M-3	
	CTX-M-55	CTX-M-15	
	CTX-M-167	CTX-M-22	
	CTX-M-177	CTX-M-28	
	CTX-M-187	CTX-M-32	
	CTX-M-224	CTX-M-54	
<i>Escherichia coli</i>	CTX-M-1	<i>Klebsiella pneumoniae</i>	CTX-M-55
	CTX-M-2		CTX-M-71
	CTX-M-3		CTX-M-72
	CTX-M-4		CTX-M-118
	CTX-M-5		CTX-M-124
	CTX-M-6		CTX-M-129
	CTX-M-7		CTX-M-130

TABLE 7: Predicted (*in silico*) Reactivity for CTX-M group 1

<i>Associated Target Organism</i>	<i>Variant Detected</i>	<i>Associated Target Organism</i>	<i>Variant Detected</i>
	CTX-M-8		CTX-M-133
	CTX-M-9		CTX-M-135
	CTX-M-10		CTX-M-138
	CTX-M-11		CTX-M-139
	CTX-M-12		CTX-M-173
	CTX-M-15		CTX-M-176
	CTX-M-28		CTX-M-183
	CTX-M-29		CTX-M-188
	CTX-M-32		CTX-M-197
	CTX-M-33		CTX-M-204
	CTX-M-36		CTX-M-208
	CTX-M-42		CTX-M-210
	CTX-M-55		CTX-M-220
	CTX-M-58		CTX-M-15
	CTX-M-69		CTX-M-66
	CTX-M-71		CTX-M-116
	CTX-M-79	<i>Proteus species</i>	CTX-M-136
	CTX-M-82		CTX-M-164
	CTX-M-90		CTX-M-167
	CTX-M-101		CTX-M-212
	CTX-M-102		CTX-M-1
	CTX-M-103	<i>Pseudomonas aeruginosa</i>	CTX-M-15
	CTX-M-109		CTX-M-32
	CTX-M-117		CTX-M-3
	CTX-M-120		CTX-M-15
	CTX-M-125	<i>Serratia marcescens</i>	CTX-M-22
	CTX-M-127		CTX-M-55
	CTX-M-128		CTX-M-221
	CTX-M-131		
	CTX-M-132		
	CTX-M-134		
	CTX-M-137		
	CTX-M-138		
	CTX-M-139		
	CTX-M-140		
	CTX-M-141		
	CTX-M-142		
	CTX-M-143		
	CTX-M-146		
	CTX-M-158		
	CTX-M-163		
	CTX-M-166		
	CTX-M-167		
	CTX-M-170		
	CTX-M-172		
	CTX-M-175		
	CTX-M-178		
	CTX-M-179		
	CTX-M-180		
	CTX-M-181		
	CTX-M-182		
	CTX-M-184		

TABLE 7: Predicted (<i>in silico</i>) Reactivity for CTX-M group 1			
<i>Associated Target Organism</i>	<i>Variant Detected</i>	<i>Associated Target Organism</i>	<i>Variant Detected</i>
	CTX-M-186		
	CTX-M-188		
	CTX-M-189		
	CTX-M-193		
	CTX-M-194		
	CTX-M-202		
	CTX-M-203		
	CTX-M-207		
	CTX-M-211		
	CTX-M-216		
	CTX-M-218		
	CTX-M-222		
	CTX-M-226		

TABLE 8: Predicted (<i>in silico</i>) Reactivity for KPC			
<i>Associated Target Organism</i>	<i>Variant Detected</i>	<i>Associated Target Organism</i>	<i>Variant Detected</i>
<i>Acinetobacter baumannii</i> complex	KPC-2		KPC-1
	KPC-3		KPC-2
	KPC-10		KPC-3
<i>Enterobacter cloacae</i> complex	KPC-1		KPC-4
	KPC-2		KPC-5
	KPC-3		KPC-6
	KPC-4		KPC-7
	KPC-13		KPC-8
	KPC-18		KPC-11
	KPC-47		KPC-14
<i>Escherichia coli</i>	KPC-2		KPC-15
	KPC-3		KPC-16
	KPC-12		KPC-17
	KPC-18	KPC-19	
	KPC-20	KPC-22	
	KPC-21	KPC-23	
	KPC-28	KPC-25	
<i>Klebsiella oxytoca</i>	KPC-2	KPC-26	
	KPC-3	KPC-27	
<i>Proteus</i> species	KPC-1	KPC-29	
	KPC-2	KPC-30	
<i>Pseudomonas aeruginosa</i>	KPC-2	KPC-31	
	KPC-5	KPC-32	
<i>Serratia marcescens</i>	KPC-2	KPC-33	
		KPC-34	
		KPC-35	
		KPC-36	
		KPC-37	
		KPC-38	
		KPC-39	
		KPC-42	
		KPC-43	
		KPC-59	

TABLE 9: Predicted (<i>in silico</i>) Reactivity for NDM			
Associated Target Organism	Variant Detected	Associated Target Organism	Variant Detected
<i>Acinetobacter baumannii</i> complex	NDM-1	<i>Klebsiella oxytoca</i>	NDM-1
	NDM-2		NDM-3
	NDM-3		NDM-4
	NDM-4	<i>Klebsiella pneumoniae</i>	NDM-1
	NDM-5		NDM-3
	NDM-7		NDM-4
	NDM-14		NDM-5
<i>Enterobacter cloacae</i> complex	NDM-1		NDM-6
	NDM-4		NDM-7
	NDM-5		NDM-9
	NDM-7		NDM-10
	NDM-22		NDM-16
<i>Escherichia coli</i>	NDM-1		
	NDM-2		NDM-28
	NDM-3	<i>Proteus</i> species	NDM-1
	NDM-4	<i>Pseudomonas aeruginosa</i>	NDM-1
	NDM-5		NDM-5
	NDM-6	<i>Serratia marcescens</i>	NDM-1
	NDM-7		NDM-4
	NDM-9		NDM-12
	NDM-11		
	NDM-12		
	NDM-13		
	NDM-15		
	NDM-16		
	NDM-17		
	NDM-18		
	NDM-19		
	NDM-20		
	NDM-21		
	NDM-27		

Exclusivity

To demonstrate the exclusivity of the iC-GN Assay, a comprehensive panel of non-target organisms that may be encountered in positive blood cultures was evaluated. A total of 114 strains were tested including organisms phylogenetically related to iC-GN target organisms as well as common blood culture contaminants. Potential cross-reactivity was evaluated by testing exclusivity panel organisms at the highest possible concentrations, considered eight hours beyond initial bottle positivity or the equivalent. Organisms were grown in BD BACTEC Plus Aerobic blood culture bottles with human blood added on the BD BACTEC System. Each strain was tested in triplicate. Performance is based on the observation of all expected negative results. In the event of a false positive result or other failure, the organism was retested in replicates of three (3) or ten (10). Exclusivity results are presented in Table 10 below. Three (3) strains demonstrated reproducible cross-reactivity with iC-GN Assay targets: *Acinetobacter haemolyticus* cross-reacted with *Acinetobacter baumannii* complex, *Klebsiella*

variicola cross-reacted with *Klebsiella pneumoniae*, and *Serratia odorifera* cross-reacted with *Serratia marcescens*.

TABLE 10: iC-GN Assay Exclusivity Results			
Organism	Strain	Concentration (CFU/mL)	Performance
<i>Acinetobacter haemolyticus</i>	ATCC 19002	7.20×10^8	0/3 ¹
<i>Acinetobacter lwoffii</i>	Z141	2.45×10^8	3/3
<i>Acinetobacter radioresistens</i>	ATCC 43998	5.20×10^7	3/3
<i>Acinetobacter schindleri</i>	ATCC BAA618	3.50×10^8	3/3
<i>Acinetobacter ursingii</i>	ATCC BAA617	3.80×10^8	3/3
<i>Aerococcus viridans</i>	Z219	2.24×10^7	3/3
<i>Aeromonas hydrophila</i>	Z161	8.10×10^8	3/3
<i>Alcaligenes faecalis</i>	Z218	9.70×10^8	3/3
<i>Aspergillus niger</i>	Z105	1.62×10^8	3/3
<i>Bacillus cereus</i>	Z091	ND	3/3
<i>Bacteroides fragilis</i>	Z029	8.40×10^9	3/3
<i>Brevundimonas vesicularis</i>	ATCC 11426	3.80×10^8	5/5
<i>Burkholderia cepacia</i>	ATCC 25416	5.40×10^8	3/3
<i>Campylobacter coli</i>	Z293	3.90×10^8	3/3
<i>Campylobacter jejuni</i>	Z086	4.60×10^8	11/13 ²
<i>Candida albicans</i>	Z006	ND	3/3
<i>Candida glabrata</i>	Z007	3.20×10^7	3/3
<i>Candida krusei</i>	Z009	1.90×10^7	3/3
<i>Candida parapsilosis</i>	Z011	9.00×10^6	3/3
<i>Candida tropicalis</i>	Z012	3.50×10^7	4/4
<i>Cedecea davisae</i>	ATCC 33431	6.20×10^8	3/3
<i>Citrobacter amalonaticus</i>	Z051	8.4×10^8	3/3
<i>Citrobacter braakii</i>	ATCC 51113	4.90×10^8	3/3
<i>Citrobacter freundii</i>	Z064	2.25×10^8	3/3
<i>Citrobacter koseri</i>	Z039	1.14×10^9	3/3
<i>Citrobacter sedlakii</i>	ATCC 51115	9.80×10^8	2/2
<i>Clostridium difficile (NAP-1 toxigenic)</i>	NAP1	4.87×10^7	4/4
<i>Clostridium difficile (non-toxigenic)</i>	Z228	5.93×10^7	3/3
<i>Clostridium novyi*</i>	Z179	1.14×10^7	5/5
<i>Corynebacterium amycolatum</i>	Z284	9.26×10^8	3/3
<i>Corynebacterium genitalium</i>	Z328	1.35×10^8	3/3
<i>Corynebacterium jeikeium</i>	Z232	8.50×10^8	4/4
<i>Corynebacterium striatum</i>	MCW000	2.07×10^9	5/6 ³
<i>Cronobacter muytjensii</i>	ATCC 51329	2.79×10^8	3/3
<i>Cronobacter sakazakii</i>	ATCC 29544	6.90×10^8	3/3
<i>Cryptococcus neoformans</i>	Serotype A	2.15×10^8	3/3
<i>Edwardsiella tarda</i>	Z183	8.70×10^7	4/5 ⁴
<i>Enterobacter aerogenes</i>	Z052	1.77×10^9	5/5
<i>Enterobacter amnigenus</i>	ATCC 51816	7.50×10^8	3/3
<i>Enterococcus avium</i>	Z171	2.58×10^8	5/6 ⁵

TABLE 10: iC-GN Assay Exclusivity Results			
Organism	Strain	Concentration (CFU/mL)	Performance
<i>Enterococcus casseliflavus</i>	Z002	2.44×10^9	4/4
<i>Enterococcus cecorum</i>	Z208	1.03×10^9	5/6 ⁶
<i>Enterococcus faecalis</i>	ATCC 51299	2.13×10^9	3/3
<i>Enterococcus faecium</i>	ATCC 700221	7.20×10^8	3/3
<i>Enterococcus gallinarum</i>	Z209	1.35×10^9	3/3
<i>Enterococcus hirae</i>	Z193	2.37×10^8	3/3
<i>Enterococcus raffinosus</i>	ATCC 49427	5.40×10^8	3/3
<i>Escherichia fergusonii</i>	ATCC 35469	8.70×10^8	3/3
<i>Escherichia hermannii</i>	Z184	1.01×10^9	5/5
<i>Escherichia vulneris</i>	ATCC 33821	7.50×10^8	3/3
<i>Fusobacterium varium</i>	Z361	2.49×10^9	3/3
<i>Hafnia alvei</i>	ATCC 51815	1.37×10^9	3/3
<i>Haemophilus influenzae</i>	ATCC 10211	3.09×10^9	3/3
<i>Haemophilus parainfluenzae</i>	ATCC 9796	1.33×10^8	3/3
<i>Klebsiella variicola</i>	ATCC 31488	4.40×10^8	0/3 ⁷
<i>Kluyvera ascorbata</i> (KPC+)	CDC-0144	1.40×10^9	3/3
<i>Kocuria kristinae</i>	Z250	7.20×10^7	3/3
<i>Kytococcus schroeteri</i>	ATCC BAA2410	1.50×10^7	3/3
<i>Lactobacillus acidophilus</i>	Z048	6.00×10^8	3/3
<i>Lactobacillus plantarum</i>	17-5	5.30×10^8	3/3
<i>Lactobacillus reuteri</i>	Z333	5.80×10^7	5/5
<i>Lactococcus lactis</i>	Z169	9.30×10^7	3/3
<i>Leclercia adecarboxylata</i>	ATCC 23216	1.01×10^9	3/3
<i>Leminorella grimontii</i>	Z364	4.00×10^9	3/3
<i>Leuconostoc mesenteroides</i>	Z197	4.00×10^7	5/5
<i>Listeria monocytogenes</i>	ATCC 19115	2.03×10^9	3/3
<i>Micrococcus luteus</i>	Z100	1.80×10^8	3/3
<i>Moraxella catarrhalis</i>	ATCC 25238	1.27×10^9	3/3
<i>Morganella morganii</i>	ATCC 25830	1.23×10^9	3/3
<i>Neisseria gonorrhoeae</i>	ATCC 19424	ND	3/3
<i>Neisseria lactamica</i>	ATCC 23970	2.90×10^8	3/3
<i>Neisseria meningitidis</i>	Serotype A	2.55×10^8	5/5
<i>Neisseria mucosa</i>	ATCC 49233	5.80×10^8	4/4
<i>Neisseria sicca</i>	ATCC 9913	1.43×10^8	3/3
<i>Pantaea agglomerans</i>	ATCC 27155	2.00×10^6	3/3
<i>Pasturella multocida</i>	ATCC 12945	2.84×10^9	2/2
<i>Pediococcus pentosaceus</i>	Z226	1.91×10^8	3/3
<i>Planococcus citreus</i>	ATCC 14404	1.95×10^8	3/3
<i>Pluralibacter gergoviae</i>	ATCC 33028	1.27×10^9	3/3
<i>Propionibacterium acnes</i>	Z144	7.90×10^8	5/5
<i>Providencia alcalifaciens</i>	Z292	3.10×10^9	3/3
<i>Providencia rettgeri</i>	Z370	2.20×10^9	3/3
<i>Providencia stuartii</i>	Z213	1.70×10^9	3/3

TABLE 10: iC-GN Assay Exclusivity Results			
Organism	Strain	Concentration (CFU/mL)	Performance
<i>Pseudomonas fluorescens</i>	ATCC 13525	2.43×10^8	3/3
<i>Pseudomonas luteola</i>	ATCC 43273	1.09×10^8	3/3
<i>Pseudomonas mendocina</i>	ATCC 25411	1.23×10^9	3/3
<i>Pseudomonas nitroreducens</i>	ATCC 33634	5.30×10^8	3/3
<i>Pseudomonas oryzihabitans</i>	ATCC 43272	1.70×10^7	4/5 ⁸
<i>Pseudomonas putida</i>	Z030	3.30×10^8	3/3
<i>Pseudomonas stutzeri</i>	ATCC 17588	6.20×10^8	3/3
<i>Raoultella planitcola</i>	ATCC 33558	1.25×10^9	3/3
<i>Rothia mucilaginosa</i>	Z033	5.50×10^7	3/3
<i>Salmonella enterica</i>	ATCC BAA1715	2.23×10^9	3/3
<i>Serratia fonticola</i>	ATCC 29844	1.18×10^9	3/3
<i>Serratia liquefaciens</i>	ATCC 27592	1.24×10^9	11/12 ⁹
<i>Serratia odorifera</i>	ATCC 33077	2.19×10^9	11/13 ¹⁰
<i>Serratia rubidaea</i>	ATCC 19278	1.54×10^8	5/5
<i>Staphylococcus aureus</i>	ATCC 700699	4.30×10^7	3/3
<i>Staphylococcus capitis</i>	Z192	2.13×10^8	3/3
<i>Staphylococcus epidermidis</i>	ATCC 700566	5.90×10^7	3/3
<i>Staphylococcus haemolyticus</i>	Z067	2.70×10^7	5/5
<i>Staphylococcus hominis</i>	Z031	9.90×10^7	3/3
<i>Staphylococcus intermedius</i>	Z112	3.30×10^7	3/3
<i>Staphylococcus lugdunensis</i>	Z097	2.43×10^8	3/3
<i>Staphylococcus schleiferi</i>	Z294	2.52×10^9	3/3
<i>Stenotrophomonas maltophilia</i>	ATCC BAA84	1.74×10^9	5/5
<i>Streptococcus agalactiae</i>	Z019	5.80×10^8	3/3
<i>Streptococcus anginosus</i>	Z179	9.50×10^8	3/3
<i>Streptococcus bovis</i>	Z167	8.00×10^8	3/3
<i>Streptococcus dysgalactiae</i>	Z068	2.65×10^8	3/3
<i>Streptococcus intermedius</i>	Z126	1.40×10^7	5/6 ¹¹
<i>Streptococcus pneumoniae</i>	ATCC 6301	3.80×10^8	5/6 ¹²
<i>Streptococcus pyogenes</i>	Z018	4.80×10^7	3/3
<i>Veillonella parvula</i>	Z379	6.70×10^9	5/6 ¹³

- 1) 3/3 false positive *A. baumannii* complex in initial testing. See limitation.
- 2) 2/3 false positive *E. coli* in initial testing. 10/10 repeats negative.
- 3) 1/3 false positive *S. marcescens* in initial testing. 3/3 repeats negative.
- 4) 1/3 positive control check failure in initial testing. 1/3 false positive *S. marcescens* in repeat testing.
- 5) 1/3 false positive *S. marcescens* in initial testing. 3/3 repeats negative.
- 6) 1/3 false positive *E. coli* in initial testing. 3/3 repeats negative.
- 7) 3/3 false positive *K. pneumoniae* in initial testing. See limitation.
- 8) 1/3 false positive *S. marcescens* in initial testing. 2/2 repeats negative.
- 9) 1/3 false positive *S. marcescens* in initial testing. 9/9 repeats negative.
- 10) 1/3 false positive *S. marcescens* in initial testing. 1/10 false positive *S. marcescens* in repeat testing. See limitation.
- 11) 1/3 false positive *S. marcescens* in initial testing. 3/3 repeats negative.

- 12) 1/3 false positive *S. marcescens* in initial testing. 3/3 repeats negative.
 13) 1/3 false positive *E. coli* in initial testing. 3/3 repeats negative.

Microbial Interference

Potential microbial interference was evaluated by testing high concentrations of gram negative exclusivity organisms in combination with low concentrations of iC-GN target organisms. A total of sixty (60) gram negative exclusivity strains were tested at the highest possible concentrations, considered eight hours beyond initial bottle positivity or the equivalent. Eight (8) representative iC-GN target organisms were tested at concentrations below the lowest levels of bottle positivity. Each organism combination was tested in triplicate. Performance was based on all expected targets detected and no false positive targets detected. In the event of a false negative result, the combination was retested in replicates of ten (10). In the event of a false positive result or other failure, the combination was retested in replicates of three (3) or ten (10). Microbial interference results are presented in Table 11 below.

TABLE 11: iC-GN Assay Microbial Interference Results											
Organism	ABX	ECX	EC	KO	KPN	PM	PA	SM	CTX-M-15	KPC-2	NDM-1
<i>A. lwoffii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>A. radioresistens</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>A. schindleri</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>A. ursingii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>A. hydrophila</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>A. faecalis</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>B. fragilis</i>	3/3	3/3	3/3	3/3	3/3	4/5 ¹	3/3	3/3	3/3	3/3	3/3
<i>B. vesicularis</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>B. cepacia</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. coli</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. jejuni</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. davisae</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. amalonaticus</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. braakii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. freundii</i>	3/3	3/3	3/3	3/3	3/3	2/2	3/3	3/3	3/3	3/3	3/3
<i>C. koseri</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. sedlakii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. muytjensii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>C. sakazakii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>E. tarda</i>	3/3	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	5/5
<i>E. aerogenes</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>E. amnigenus</i>	3/3	3/3	5/5	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3
<i>E. fergusonii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>E. hermannii</i>	3/3	3/3	5/5	3/3	13/13	3/3	3/3	3/3	5/5	12/13 ²	3/3

TABLE 11: iC-GN Assay Microbial Interference Results											
Organism	ABX	ECX	EC	KO	KPN	PM	PA	SM	CTX-M-15	KPC-2	NDM-1
<i>E. vulneris</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>F. varium</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>H. alvei</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>H. influenzae</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>H. parainfluenzae</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>K. ascorbata</i>	10/13 ³	3/3	11/12 ⁴	3/3	3/3	3/3	3/3	3/3	10/12 ⁴	3/3	3/3
<i>L. adecarboxylata</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>L. grimontii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>M. catarrhalis</i>	12/13 ⁵	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>M. morgani</i>	3/3	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	5/5
<i>N. gonorrhoeae</i>	5/5	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>N. lactamica</i>	3/3	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	5/5
<i>N. meningitidis</i>	3/3	3/3	5/5	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3
<i>N. mucosa</i>	5/5	3/3	5/5	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3
<i>N. sicca</i>	5/6 ⁶	3/3	3/3	5/5	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. agglomerans</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. multocida</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. gergoviae</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. alcalifaciens</i>	5/5	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. rettgeri</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. stuartii</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. fluorescens</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. luteola</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. mendocina</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3
<i>P. nitroreducens</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>P. oryzihabitans</i>	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	3/3	5/5	3/3
<i>P. putida</i>	3/3	3/3	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	5/5
<i>P. stutzeri</i>	3/3	5/5	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	5/5
<i>R. planitcola</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>S. enterica</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>S. fonticola</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>S. liquefaciens</i>	3/3	5/5	3/3	3/3	12/13 ⁷	3/3	3/3	3/3	3/3	13/13	3/3
<i>S. odorifera</i>	3/3	3/3	3/3	12/13 ⁸	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>S. rubidaea</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>S. maltophilia</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
<i>V. parvula</i>	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3

- 1) 1/3 array registration error in initial testing. 1/3 false positive *E. cloacae* complex in repeat testing.
- 2) 1/3 false negative KPC in initial testing. 10/10 repeats passed.
- 3) 3/3 false negative *A. baumannii* complex in initial testing; concentration was determined to be below the target limit of detection. 10/10 repeats passed.

- 4) 1/2 false negative CTX-M in initial testing. 1/10 false positive *K. pneumoniae* in repeat testing.
- 5) 1/3 false negative *A. baumannii* complex in initial testing. 10/10 repeats passed.
- 6) 1/3 false positive *E. coli* in initial testing. 3/3 repeats passed.
- 7) 1/3 false negative *K. pneumoniae* in initial testing. 10/10 repeats passed.
- 8) 1/3 false positive *S. marcescens* in initial testing. 10/10 repeats passed.

Competitive Inhibition

iC-GN Assay performance was evaluated with combinations of target analytes that may be found in mixed positive blood cultures. One target organism was prepared at the lowest level of bottle positivity, considered within two hours of bottle “ring”, while the second target organism was prepared at the highest possible concentration, considered eight hours after initial bottle positivity. All organisms were grown in BD BACTEC Plus Aerobic blood cultures bottles with human blood added on the BD BACTEC System. The organisms were combined at a ratio of one part “low” to four parts “high”. Each low concentration organism was tested in combination with each high concentration organism in triplicate. Performance was based on all expected targets detected. In the event of a false negative result, the organism combination was retested in replicates of ten (10) at the same “low” and “high” organism ratio. In the event of a reproducible false negative result, the organism combination was retested in replicates of ten (10) at a ratio of one part “low” to one part “high.” Competitive inhibition results are presented in Table 12 below. All high concentration iC-GN targets were detected. Due to competitive inhibition, low concentration targets were not detected in 1.7% of tests (3/178). When iC-GN target organisms were present at similar concentrations, all targets were detected.

TABLE 12: iC-GN Assay Competitive Inhibition Performance					
		Target Performance			
Low Organism	High Organism	Low Organism	Low Marker	High Organism	High Marker
<i>A. baumannii</i>	ECX	3/3	NA	3/3	NA
	EC (CTX-M-15+)	3/3	NA	3/3	3/3
	KO	3/3	NA	3/3	NA
	KPN (KPC-2+)	3/3	NA	3/3	3/3
	PM	3/3	NA	3/3	NA
	PA (NDM-1+)	3/3	NA	3/3	3/3
	SM	3/3	NA	3/3	NA
<i>E. cloacae</i>	ABX	3/3	NA	3/3	1 FP KPC
	EC (CTX-M-15+)	3/3	NA	3/3	3/3
	KO	3/3	NA	3/3	NA
	KPN (KPC-2+)	3/3	NA	3/3	3/3
	PM	3/3	NA	3/3	NA
	PA (NDM-1+)	3/3	NA	3/3	3/3
	SM	3/3	NA	3/3	NA
<i>E. coli</i>	ABX	3/3	3/3	3/3	NA

TABLE 12: iC-GN Assay Competitive Inhibition Performance					
		Target Performance			
Low Organism	High Organism	Low Organism	Low Marker	High Organism	High Marker
(CTX-M-15+)	ECX	3/3	3/3	3/3	NA
	KO	3/3	3/3	3/3	NA
	KPN (KPC-2+)	3/3	3/3	3/3	3/3
	PM	3/3	3/3	3/3	NA
	PA (NDM-1+)	3/3	3/3	3/3	3/3
	SM	3/3	3/3	3/3	NA
<i>K. oxytoca</i>	ABX	3/3	NA	3/3	NA
	ECX	3/3	NA	3/3	NA
	EC (CTX-M-15+)	3/3	NA	3/3	3/3
	KPN (KPC-2+)	3/3	NA	3/3	3/3
	PM	3/3	NA	3/3	NA
	PA (NDM-1+)	3/3	NA	3/3	3/3
	SM	3/3	3/3	3/3	NA
<i>K. pneumoniae</i> (KPC-2+)	ABX	3/3	3/3	3/3	NA
	ECX	3/3	3/3	3/3	NA
	EC (CTX-M-15+)	3/3	3/3	3/3	3/3
	KO	3/3	3/3	3/3	NA
	PM	3/3	3/3	3/3	NA
	PA (NDM-1+)	3/3	3/3	3/3	3/3
	SM	3/3	3/3	3/3	NA
<i>P. mirabilis</i>	ABX	3/3	NA	3/3	NA
	ECX	3/3	NA	3/3	NA
	EC (CTX-M-15+)	3/3	NA	3/3	3/3
	KO	5/5	NA	5/5	NA
	KPN (KPC-2+)	3/3	NA	3/3	3/3
	PA (NDM-1+)	3/3	NA	3/3	3/3
	SM	3/3	NA	3/3	NA
<i>P. aeruginosa</i> (NDM-1+)	ABX	3/3	3/3	3/3	NA
	ECX	3/3	3/3	3/3	NA
	EC (CTX-M-15+)	12/13	11/13	13/13	13/13
	EC (CTX-M-15+) 1:1	10/10	10/10	10/10	10/10
	KO	3/3	3/3	3/3	NA
	KPN (KPC-2+)	3/3	3/3	3/3	3/3
	PM	3/3	3/3	3/3	NA
	SM	3/3	3/3	3/3	NA
<i>S. marcescens</i>	ABX	3/3	NA	3/3	NA
	ECX	3/3	NA	3/3	NA
	EC (CTX-M-15+)	3/3	NA	3/3	3/3
	KO	3/3	NA	3/3	NA

TABLE 12: iC-GN Assay Competitive Inhibition Performance					
		Target Performance			
Low Organism	High Organism	Low Organism	Low Marker	High Organism	High Marker
	KPN (KPC-2+)	3/3	NA	3/3	3/3
	PM	3/3	NA	3/3	NA
	PA (NDM-1+)	3/3	NA	3/3	3/3

Interfering Substances

iC-GN Assay performance was evaluated in the presence of potentially inhibiting substances that may be encountered in blood and blood culture media. Eight representative target organisms plus one non-target organism were evaluated. Organisms were tested at the lowest levels of bottle positivity, considered within two hours of bottle “ring.” Potential interferents were tested at concentrations exceeding the highest concentrations that may be encountered in blood and blood culture media (Table 12). Target performance is based on all expected targets detected and no false positive targets detected. Non-target performance is based on all negative results. In the event of a false negative result, the organism/interferent combination was retested in replicates of ten (10). In the event of a false positive result or other failure, the organism/interferent combination was retested in triplicate. If the discordant result was observed in repeat testing, the combination was retested at a decreased inhibitor concentration. Interference results are presented in Table 13 below. Interference testing was performed in BD BACTEC Plus Aerobic blood culture bottle media, which has a sodium polyanetholesulfonate (SPS) concentration of 0.05% w/v. Additional SPS at a concentration greater than 0.05% w/v was found to interfere with the performance of some iC-GN Assay targets, resulting in increased false negative results and positive control check failures.

TABLE 13: Interfering Substances Test Panel		
Interference Compound	Clinically Relevant Concentration	Test Concentration
Hemoglobin	1-2 g/L	10 g/L
Conjugated Bilirubin	0.1-0.4 mg/dL	10 mg/dL
Unconjugated Bilirubin	0.1-0.8 mg/dL	10 mg/dL
Protein (γ-globulin + albumin)	0.7-1.7 g/dL	4 g/dL
Triglyceride	300-500 mg/dL	1500 mg/dL
Human Genomic DNA	NA	1 × 10 ⁶ cells/mL
Sodium Polyanetholesulfonate (SPS)	0.02-0.05% w/v	0.1% w/v
Cefepime	16 µg/mL	80 µg/mL

TABLE 14: iC-GN Assay Interfering Substances Performance												
Interference Compound	Target Performance											
	ABX	ECX	EC	KO	KPN	PM	PA	SM	SE	KPC-2	NDM-1	CTX-M-15
Vancomycin	3/3	3/3	3/3	3/3	5/5	3/3	3/3	3/3	3/3	5/5	3/3	3/3

- 1) 2/3 false negative *A. baumannii* complex in initial testing. 7/9 false negative *A. baumannii* complex in repeat testing.
- 2) 3/3 false negative *P. aeruginosa* in initial testing. 7/9 false negative *P. aeruginosa* in repeat testing.
- 3) 4/9 false negative NDM in repeat testing.
- 4) 1/3 false positive *K. pneumoniae* in initial testing. 3/3 repeats passed.

Method Comparison

A method comparison study was performed at five (5) geographically dispersed clinical sites. Sites tested 1002 leftover de-identified specimens from anaerobic and aerobic blood culture bottles flagged as positive by their respective continuous monitoring blood culture system. Three of the commonly used blood culture systems were included in the study: Thermo Fisher VersaTREK, BD BACTEC and BioMerieux BacT/ALERT.

Patient positive blood cultures confirmed by Gram stain to be positive for gram negative bacilli were enrolled in the study. Any positive blood cultures showing an initial mixed Gram stain were not enrolled or were subsequently withdrawn from the study dataset.

Final performance of the iC-GN Assay organism targets was compared to reference culture followed by MALDI identification per the study protocol. Final performance of the iC-GN Assay resistance marker targets was compared to PCR amplification followed by confirmatory bi-directional sequencing. Phenotypic antimicrobial susceptibility testing (AST) was also performed on all specimens to identify additional samples which required sequencing. Discordant samples were also sequenced.

To supplement performance of observed lower prevalence organisms, 170 contrived samples were prepared using verified strains. Contrived samples were prepared at iCubate using BD BACTEC Plus Aerobic Blood Culture Bottles with 10mL of human blood added (in accordance with BACTEC instructions). Organisms were spiked into bottles at a concentration of 5-30 CFU/bottle and incubated until bottles were flagged as positive. Aliquots of samples were frozen and provided to the sites (frozen) for testing.

Of the 1107 positive blood culture specimens enrolled in the study, a total of 105 specimens were excluded/withdrawn from the study and all subsequent performance analyses. Of the 1002 specimens remaining, 976 were fresh prospective specimens and 26 (2.6%) were frozen prospective specimens.

The total specimens excluded from the iC-GN Assay Method Comparison Study (n=105) are listed by site and reason for exclusion in the table below. The most common reasons for exclusion included incomplete reference testing and repeat iC-GN errors.

Table 15: Withdrawn Summary

Site Code	Unresolved iC-GN Error	Incomplete Reference Method	Outside Fresh Stability Window	Didn't Meet Inclusion Criteria	Total Withdrawn
LAC	2	26	0	0	28
IU	0	10	0	0	10
MCW	2	23	0	2	27
TC	4	25	5	0	34
TGH	1	3	1	1	6
Total	9*	87	6	3	105

*Please Note: the nine (9) samples that were excluded from performance analysis due to Unresolved iC-GN are included in the calculation of instrument errors; please refer to **Table 17: No-Calls**

Gender Demographics

Gender was reported when available for all clinical samples collected for the study. of the 1002 clinical samples included in performance analysis 47.3% were males and 52.6% were female; gender was not provided for 1/1002. The table below summarizes this data.

Table 16: Gender Stratification

Site	MALES		FEMALES		Not Provided		Clinical Samples
	#	%	#	%	#	%	
LAC	156	43.9%	199	56.1%	0	-	355
IU	41	43.6%	53	56.4%	0	-	94
MCW	109	58.0%	78	41.5%	1	0.5%	188
TC	71	45.5%	85	54.5%	0	-	156
TGH	97	46.4%	112	53.6%	0	-	209
Total	474	47.3%	527	52.6%	1	0.1%	1002

Error Rate

Throughout the course of the study, an initial error rate of 2.9 % (34/1181) was observed. Reasons for error included the following: *Positive controls check failure (27), Array registration error (6), and Processor/System error (1)*. When an error was observed, repeat testing was performed with the iC-GN Assay per the protocol. Upon repeat testing, the error rate was reduced to 0.8% (9/1181).

Table 17: No-Calls

Internal Positive Control Failure		Instrument Errors		Total Non-Reportable Rate	
Initial #fail/#total (95% CI)	Final #fail/#total (95% CI)	Initial #fail/#total (95% CI)	Final #fail/#total (95% CI)	Initial #fail/#total (95% CI)	Final #fail/#total (95% CI)
2.3%	0.8%	0.6%	-	2.9%	0.8%
27/1181	9/1181	7/1181	0/1181	34/1181	9/1181
[1.6-3.3%]	[0.4-1.4%]	[0.3-1.2%]	[0.0-0.3%]	[2.1-4.0%]	[0.4-1.4%]

When performance of the iC-GN Assay was compared to reference culture followed by MALDI identification or PCR/bi-directional sequencing, there was no apparent difference in performance noted between the five study sites or between the three blood culture systems. Performance for all positive bottle types/systems combined is presented in the tables below for detection of the iC-GN Assay targets as compared to culture and MALDI or PCR/bi-directional sequencing. Results are stratified by prospectively tested fresh specimens, prospectively collected/retrospectively tested frozen specimens and contrived specimens.

Table 18: iC-GN Assay Performance: *Acinetobacter baumannii* complex (*ppa*)

Specimen Type	N=	Percent Agreement		Comparator Method	
		Positive (95% CI)	Negative (95% CI)		
Prospective	Fresh	976	100% 7/7 (64.6-100)	99.9% 968-969** (99.4-100)	Culture & MALDI
	Frozen	26	- 0/0 -	100% 26/26 (87.1-100)	
	TOTAL	1002	100% 7/7 (64.6-100)	99.9% 994/995 (99.4-100)	
Contrived	170		100% 45/45 (92.1-100)	100% 125/125 (97.0-100)	

**1/1 false positive observed was negative for *A. baumannii* complex by PCR/bi-directional sequencing

Table 19: iC-GN Assay Performance: *Enterobacter cloacae* complex (*ramA*)

Specimen Type	N=	Percent Agreement		Comparator Method	
		Positive (95% CI)	Negative (95% CI)		
Prospective	Fresh	976	94.5% 52/55* (85.1-98.1)	100% 921/921 (99.6-100)	Culture & MALDI
	Frozen	26	100% 5/5 (56.6-100)	100% 21/21 (84.5-100)	
	TOTAL	1002	95.0% 57/60 (86.3-98.3)	100% 942/942 (99.6-100)	
Contrived	170		100% 17/17 (81.6-100)	100% 153/153 (97.6-100)	

*1/3 false negatives observed was negative for *E. cloacae* complex by PCR/bi-directional sequencing; 2/3 were positive for *E. cloacae* complex by PCR/bi-directional sequencing

Table 20: iC-GN Assay Performance: *Escherichia coli* (*uidA*)

Specimen Type	N=	Percent Agreement		Comparator Method	
		Positive (95% CI)	Negative (95% CI)		
Prospective	Fresh	976	98.4% 480/488* (96.8-99.2)	100% 488/488 (99.2-100)	Culture & MALDI
	Frozen	26	100% 6/6 (61.0-100)	100% 20/20 (83.9-100)	
	TOTAL	1002	98.4% 486/494 (96.8-99.2)	100% 508/508 (99.2-100)	
Contrived	170		100% 15/15 (79.6-100)	100% 155/155 (97.6-100)	

*4/8 false negatives observed were negative for *E. coli* by PCR/bi-directional sequencing; 3/8 were positive for *E. coli* by PCR/bi-directional sequencing; 1/8 was not available for sequencing

Table 21: iC-GN Assay Performance: *Klebsiella oxytoca* (*pehX*)

Specimen Type	N=	Percent Agreement		Comparator Method	
		Positive (95% CI)	Negative (95% CI)		
Prospective	Fresh	976	95.8% 23/24* (79.8-99.3)	99.7% 949/952** (99.1-99.9)	Culture & MALDI
	Frozen	26	- 0/0 -	100% 26/26 (87.1-100)	
	TOTAL	1002	95.8% 23/24 (79.8-99.3)	99.7% 975/978 (99.1-99.9)	
Contrived	170		100% 30/30 (88.6-100)	100% 140/140 (97.3-100)	

*1/1 false negative observed was negative for *K. oxytoca* by PCR/bi-directional sequencing
 **3/3 false positives observed were negative for *K. oxytoca* by PCR/bi-directional sequencing

Table 22: iC-GN Assay Performance: *Klebsiella pneumoniae* (*parC*)

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	976	96.8% 150/155* (92.7-98.6)	99.3% 815/821** (98.4-99.7)	Culture & MALDI
	Frozen	26	100% 3/3 (43.9-100)	100% 23/23 (85.7-100)	
	TOTAL	1002	96.8% 153/158 (92.8-98.6)	99.3% 838/844 (98.4-99.7)	
Contrived		170	100% 21/21 (84.5-100)	99.3% 148/149 (96.3-99.9)	

*3/5 false negatives observed were negative for *K. pneumoniae* by PCR/bi-directional sequencing; 2/3 were positive for *K. pneumoniae* by PCR/bi-directional sequencing

**6/6 false positives observed were negative for *K. pneumoniae* by PCR/bi-directional sequencing

Nineteen (19) samples were excluded from *Proteus mirabilis* performance analysis due to confirmed *Proteus* contamination within the BD BACTEC Bottles leaving a total of 983 evaluable specimens.

Table 23: iC-GN Assay Performance: *Proteus mirabilis* (*rpoB*)

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	957	97.4% 37/38* (86.5-99.5)	99.5% 914/919** (98.7-99.8)	Culture & MALDI
	Frozen	26	100% 9/9 (70.1-100)	100% 17/17 (81.6-100)	
	TOTAL	983	97.9% 46/47 (88.9-99.6)	99.5% 931/936 (98.8-99.8)	
Contrived		170	100% 12/12 (75.8-100)	100% 158/158 (97.6-100)	

*1/1 false negative observed was positive for *P. mirabilis* by PCR/bi-directional sequencing

**3/5 false positives observed were negative for *P. mirabilis* by PCR/bi-directional sequencing; 2/5 were not available for sequencing

Table 24: iC-GN Assay Performance: *Pseudomonas aeruginosa* (*algD*)

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	976	95.1% 78/82* (88.1-98.1)	99.8% 892/894** (99.2-99.9)	Culture & MALDI
	Frozen	26	100% 1/1 (20.7-100)	100% 25/25 (86.7-100)	
	TOTAL	1002	95.2% 79/83 (88.3-98.1)	99.8% 917/919 (99.2-99.9)	
Contrived		170	100% 10/10 (72.2-100)	100% 160/160 (97.7-100)	

*4/4 false negatives observed were positive for *P. aeruginosa* by PCR/bi-directional sequencing

**2/2 false positives observed were negative for *P. aeruginosa* by PCR/bi-directional sequencing

Table 25: iC-GN Assay Performance: *Serratia marcescens* (*gyrB*)

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	976	100% 29/29 (88.3-100)	99.6% 943/947** (98.9-99.8)	Culture & MALDI
	Frozen	26	- 0/0 -	100% 26/26 (87.1-100)	
	TOTAL	1002	100% 29/29 (88.3-100)	99.6% 969/973 (98.9-99.8)	
Contrived		170	100% 20/20 (83.9-100)	99.3% 149/150 (96.3-99.9)	

**1/4 false positives observed was positive for *S. marcescens* by PCR/bi-directional sequencing; 3/4 were negative for *S. marcescens* by PCR/bi-directional sequencing

Table 26: iC-GN Assay Performance: CTX-M

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	976	97.0% 64/66 (89.6-99.2)	99.9% 909/910 (99.4-100)	PCR/Bi-directional sequencing
	Frozen	26	100% 1/1 (20.7-100)	100% 25/25 (86.7-100)	
	TOTAL	1002	97.0% 65/67 (89.8-99.2)	99.9% 934/935 (99.4-100)	
Contrived		170	100% 15/15 (79.6-100)	100% 155/155 (97.6-100)	

Table 27: iC-GN Assay Performance: KPC

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	976	100% 1/1 (20.7-100)	99.9% 974/975 (99.4-100)	PCR/Bi-directional sequencing
	Frozen	26	- 0/0 -	100% 26/26 (87.1-100)	
	TOTAL	1002	100% 1/1 (20.7-100)	99.9% 1000/1001 (99.4-100)	
Contrived		170	100% 50/50 (92.9-100)	99.2% 119/120 (95.4-99.9)	

Table 28: iC-GN Assay Performance: NDM

Specimen Type		N=	Percent Agreement		Comparator Method
			Positive (95% CI)	Negative (95% CI)	
Prospective	Fresh	976	- 0/0 -	100% 976/976 (99.6-100)	PCR/Bi-directional sequencing
	Frozen	26	- 0/0 -	100% 26/26 (87.1-100)	
	TOTAL	1002	- 0/0 -	100% 1002/1002 (99.6-100)	
Contrived		170	100% 50/50 (92.9-100)	100% 120/120 (96.9-100)	

Analysis of Mixed Culture Results:

In the method comparison study, there were thirty (30) mixed culture specimens that were detected by the iC-GN Assay, culture and MALDI, or both. The tables below list the mixed target combinations detected by iC-GN and the comparator method in the clinical study. There were twelve (12) discrepant mixed samples for which iC-GN detected a target that was not detected by the comparator assay. There were four (4) discrepant mixed samples for which the comparator assay detected targets that were not detected by iC-GN. Due to competitive

inhibition, target organisms present at low concentrations may not be detected by the iC-GN Assay when a second target organism is present at higher concentrations.

TABLE 29: Multiple Organism Detections by iC-GN as Compared to Culture/MALDI							
Multiple Detections by iC-GN					Total Targets Detected by iC-GN	No of Discrepant Targets	Discrepant Results (Targets Not Detected by culture/MALDI)
Site	ID	Target 1	Target 2	Target 3			
LAC	1102	<i>E. coli</i>	<i>K. pneumoniae</i>		2	0	
LAC	1118	<i>E. coli</i>	<i>K. pneumoniae</i>		2	0	
LAC	1141	<i>E. coli</i>	<i>K. pneumoniae</i>		2	0	
LAC	1220	<i>E. cloacae</i> complex	<i>E. coli</i>		2	0	
LAC	1236	<i>E. coli</i>	<i>K. oxytoca</i>	<i>K. pneumoniae</i>	3	1	<i>K. oxytoca</i>
LAC	1285	<i>K. pneumoniae</i>	<i>S. marcescens</i>		2	1	<i>S. marcescens</i>
LAC	1307	<i>E. coli</i>	<i>K. oxytoca</i>		2	0	
LAC	1378	<i>E. coli</i>	<i>K. pneumoniae</i>		2	1	<i>K. pneumoniae</i>
LAC	1382	<i>K. oxytoca</i>	<i>K. pneumoniae</i>		2	1	<i>K. oxytoca</i>
MCW	2023	<i>E. coli</i>	<i>K. pneumoniae</i>		2	0	
MCW	2032	<i>E. cloacae</i> complex	<i>K. oxytoca</i>		2	0	
MCW	2038	<i>K. oxytoca</i>	<i>K. pneumoniae</i>		2	0	
MCW	2041	<i>E. coli</i>	<i>P. mirabilis</i>		2	0	
MCW	2104	<i>E. coli</i>	<i>S. marcescens</i>		2	1	<i>S. marcescens</i>
MCW	2193	<i>K. pneumoniae</i>	<i>S. marcescens</i>		2	2	<i>K. pneumoniae</i> , <i>S. marcescens</i>
TC	3015	<i>K. oxytoca</i>	<i>P. aeruginosa</i>		2	2	<i>K. oxytoca</i> , <i>P. aeruginosa</i>
TC	3096	<i>E. coli</i>	<i>K. pneumoniae</i>		2	1	<i>K. pneumoniae</i>
TC	3131	<i>E. cloacae</i> complex	<i>K. pneumoniae</i>		2	0	
TC	3183	<i>K. pneumoniae</i>	<i>S. marcescens</i>		2	1	<i>S. marcescens</i>
TGH	4031	<i>E. coli</i>	<i>P. mirabilis</i>		2	1	<i>P. mirabilis</i>
TGH	4037	<i>E. coli</i>	<i>P. aeruginosa</i>		2	0	
TGH	4124	<i>E. cloacae</i> complex	<i>P. aeruginosa</i>		2	1	<i>P. aeruginosa</i>
TGH	4132	<i>E. cloacae</i> complex	<i>K. pneumoniae</i>		2	0	
IU	5025	<i>A. baumannii</i> complex	<i>K. pneumoniae</i>		2	1	<i>A. baumannii</i> complex
IU	5031	<i>E. coli</i>	<i>K. pneumoniae</i>		2	0	
IU	5042	<i>E. cloacae</i> complex	<i>K. pneumoniae</i>		2	0	

TABLE 29: Multiple Organism Detections by Culture/MALDI as Compared to iC-GN						
Multiple Detections by culture/MALDI				Total Targets Detected by Culture	Discrepant Targets	Discrepant Targets (Targets Not Detected by iC-GN)
Site	ID	Target 1	Target 2			
LAC	1102	<i>K. pneumoniae</i>	<i>E. coli</i>	2	0	
LAC	1118	<i>K. pneumoniae</i>	<i>E. coli</i>	2	0	
LAC	1141	<i>E. coli</i>	<i>K. pneumoniae</i>	2	0	
LAC	1220	<i>E. cloacae</i> complex	<i>E. coli</i>	2	0	
LAC	1236	<i>K. pneumoniae</i>	<i>E. coli</i>	2	0	
LAC	1268	<i>P. aeruginosa</i>	<i>P. mirabilis</i>	2	1	<i>P. aeruginosa</i>
LAC	1307	<i>E. coli</i>	<i>K. oxytoca</i>	2	0	
LAC	1338	<i>E. coli</i>	<i>K. pneumoniae</i>	2	1	<i>K. pneumoniae</i>
MCW	2023	<i>K. pneumoniae</i>	<i>E. coli</i>	2	0	
MCW	2032	<i>K. oxytoca</i>	<i>E. cloacae</i> complex	2	0	
MCW	2038	<i>K. pneumoniae</i>	<i>K. oxytoca</i>	2	0	
MCW	2041	<i>P. mirabilis</i>	<i>E. coli</i>	2	0	
TC	3006	<i>E. coli</i>	<i>K. pneumoniae</i>	2	1	<i>K. pneumoniae</i>
TC	3131	<i>E. cloacae</i> complex	<i>K. pneumoniae</i>	2	0	
TGH	4007	<i>E. coli</i>	<i>P. aeruginosa</i>	2	1	<i>P. aeruginosa</i>
TGH	4037	<i>E. coli</i>	<i>P. aeruginosa</i>	2	0	
TGH	4132	<i>K. pneumoniae</i>	<i>E. cloacae</i> complex	2	0	
IU	5031	<i>E. coli</i>	<i>K. pneumoniae</i>	2	0	
IU	5042	<i>K. pneumoniae</i>	<i>E. cloacae</i> complex	2	0	

Expected Values:

A total of 1002 prospectively collected fresh and frozen blood culture specimens were obtained from five geographically dispersed clinical sites. The number and percentage of positive cases (positivity rate) determined by the iC-GN Assay stratified by U.S. state for each of the organisms and resistance markers detected by the assay are presented below. Overall, the iC-GN Assay detected at least one organism in 89% (901/1002) prospectively collected specimens and at least one resistance marker in 6.8% (68/1002) prospectively collected specimens. Expected values are presented in the table below.

TABLE 30: Positivity by the iC-GN Assay as Observed in the Clinical Study							
Organism	U.S. State	NY	WI	NM	FL	IN	TOTAL
	TOTAL n	355	188	156	209	94	1002
<i>Acinetobacter baumannii</i> complex	POSITIVE n	1	1	0	3	3	8
	% Positivity	0.3%	0.5%	0.0%	1.4%	3.2%	0.8%
	POSITIVE n	27	11	5	13	1	57

TABLE 30: Positivity by the iC-GN Assay as Observed in the Clinical Study							
Organism	U.S. State	NY	WI	NM	FL	IN	TOTAL
	TOTAL n	355	188	156	209	94	1002
<i>Enterobacter cloacae complex</i>	% Positivity	7.6%	5.9%	3.2%	6.2%	1.1%	5.7%
<i>Escherichia coli</i>	POSITIVE n	182	74	94	88	48	486
	% Positivity	51.3%	39.4%	60.3%	42.1%	51.1%	48.5%
<i>Klebsiella oxytoca</i>	POSITIVE n	9	7	2	4	4	26
	% Positivity	2.5%	3.7%	1.3%	1.9%	4.3%	2.6%
<i>Klebsiella pneumoniae</i>	POSITIVE n	51	32	28	33	15	159
	% Positivity	14.4%	15.4%	17.3%	15.8%	16.0%	15.9%
<i>Proteus mirabilis</i>	POSITIVE n	21	9	1	13	7	51
	% Positivity	5.9%	4.8%	0.6%	6.2%	7.4%	5.1%
<i>Pseudomonas aeruginosa</i>	POSITIVE n	24	19	5	25	8	81
	% Positivity	6.8%	10.1%	3.2%	12.0%	8.5%	8.1%
<i>Serratia marcescens</i>	POSITIVE n	9	8	2	8	6	33
	% Positivity	2.5%	4.3%	1.3%	3.8%	6.4%	3.3%
Resistance Marker	TOTAL n	355	188	156	209	94	1002
<i>KPC</i>	POSITIVE n	0	2	0	0	0	2
	% Positivity	0.0%	1.1%	0.0%	0.0%	0.0%	0.2%
<i>NDM</i>	POSITIVE n	0	0	0	0	0	0
	% Positivity	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
<i>CTX-M</i>	POSITIVE n	12	22	7	15	10	66
	% Positivity	3.4%	11.7%	4.5%	7.2%	10.6%	6.6%

Statement of Safety and Effectiveness

The data presented clearly demonstrates the safety and efficacy of the iC-GN Assay™ for use on the iC-System as compared to the reference method when the product Instructions for Use are followed.