## 510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY DEVICE ONLY TEMPLATE

# A. 510(k) Number:

k032299

**B.** Analyte:

Ampicillin 0.5 – 32 µg/mL Gram-Negative AST

### C. Type of Test:

Antimicrobial Susceptibility Test (Quantitative) colorimetric oxidation-reduction, growth-based

# **D.** Applicant:

Becton, Dickinson & Company

**E. Proprietary and Established Names:** BD Phoenix<sup>TM</sup> Automated Microbiology System – Ampicillin Gram Negative

# F. Regulatory Information:

- <u>Regulation section:</u>
  21 CFR 866.1645 Fully Automated Short-Term Incubation Cycle Antimicrobial
- 2. <u>Classification:</u> Class II
- 3. <u>Product Code:</u> LON
- 4. <u>Panel:</u> 83

### G. Intended Use:

1. Intended use(s):

BD Phoenix<sup>TM</sup> Automated Microbiology System:

The BD Phoenix<sup>TM</sup> Automated Microbiology System is intended for in vitro quantitative determination of antimicrobial susceptibility by minimal inhibitory concentration of gram-negative aerobic and facultative anaerobic bacteria belonging to the family *Enterobacteriaceae* and non – *Enterobacteriaceae* and gram-positive bacteria belonging to the genera *Staphylococcus* and *Enterococcus*.

The BD Phoenix<sup>™</sup> GN Panel:

The BD Phoenix<sup>TM</sup> Automated Microbiology System is intended for the *in vitro* rapid identification (ID) and quantitative determination of antimicrobial susceptibility by minimal inhibitory concentration (MIC) of gram-negative aerobic and facultative anaerobic bacteria belonging to the family *Enterobacteriaceae* and non – *Enterobacteriaceae*.

2. Indication(s) for use:

This submission is for the addition of the antibiotic ampicillin at concentrations of  $0.5 - 32 \ \mu g/mL$  to the gram negative susceptibility panel for testing *Enterobacteriaceae*.

- 3. <u>Special condition for use statement(s):</u> Not applicable
- 4. <u>Special instrument Requirements:</u> Not applicable

# H. Device Description:

The BD Phoenix<sup>™</sup> Automated Microbiology System includes instrumentation and software, sealed and self-inoculating molded polystyrene trays with 136 micro-wells containing dried reagents, and specific inoculum broth formulations for ID and AST Indicator. The organism to be tested must be a pure culture and be preliminarily identified as gram positive or gram negative. Colonies are then suspended in broth, and equated to a 0.5 McFarland with the recommendation to use the BD CrystalSpec<sup>™</sup> Nephelometer. A further dilution is made into an AST broth, which contains an AST indicator, prior to inoculating the panel. The AST broth is a cation-adjusted formulation of Mueller-Hinton broth containing 0.01% Tween 80. After adding the indicator solution to the AST inoculum the color is blue and after inoculation and incubation goes to pink to colorless as reduction in the panel well proceeds. Inoculated panels are barcode scanned and loaded into the BD Phoenix<sup>TM</sup> Automated Microbiology System instrument where the panels are continuously incubated at 35°C. The AST has a final inoculum of  $5 \times 10^5$ CFU/ml. The instrument incubates, reads and records the results of the biochemical substrates and antimicrobial agents and interprets the reactions to give an ID of the isolate and MIC value and category interpretation of the antimicrobial agents. Organisms growing in the presence of a given antimicrobic agent reduce the indicator, signaling organism growth and resistance to the antimicrobic agent. Organisms killed or inhibited by a given antimicrobic do not cause reduction of the indicator and therefore do not produce a color change. Additional interpretation is done using software driven "EXPERT" System using rules derived from the NCCLS documentation.

Readings are taken every 20 minutes with an ID result available between 2-12 hours and an AST result available between 4-16 hours. This is only an autoread result; there are no manual readings possible.

# I. Substantial Equivalence Information:

- 1. <u>Predicate device name(s):</u> VITEK® System
- 2. <u>Predicate K number(s):</u> N50510
- 3. <u>Comparison with predicate:</u>

Similarities								
Item	Device	Predicate						
1.	Isolated colonies from	Isolated colonies from						
	culture used	culture used						
2.	Inoculum density equated to	Inoculum density equated to						
	0.5 McFarland standard	0.5 McFarland standard						
3.	Report results as minimum	Report results as minimum						
	inhibitory concentration	inhibitory concentration						
	(MIC) and categorical	(MIC) and categorical						
	interpretation (SIR)	interpretation (SIR)						
4.	<16 hours	<16 hours						
Differences								
Item	Device	Predicate						
1.	Results are determined from	Results are determined from						
	serial twofold dilutions of	extrapolation of doubling						
	antimicrobial agents	dilutions						
2.	Automated growth based	Automated growth based						
	enhanced by use of a redox	with detection using an						
	indicator (colorimetric	attenuation of light						
	oxidation-reduction) to	measured by an optical						
	detect organism growth.	scanner.						

# J. Standard/Guidance Document Referenced (if applicable):

"Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems; Guidance for Industry and FDA"; NCCLS M7 (M100-S13) "Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard."

# K. Test Principle:

The system employs conventional, colorimetric, fluorogenic and chromogenic substrates to identify the genus and species of the isolate. The AST portion of the BD Phoenix<sup>TM</sup> Automated Microbiology System is a broth based microdilution method that utilizes a redox indicator (colorimetric oxidation-reduction) to enhance detection of organism growth. The MIC is determined by comparing growth in wells containing serial two-fold dilutions of an antibiotic to the growth in "growth control wells" which contain no antibiotic.

### L. Performance Characteristics (if/when applicable):

- 1. <u>Analytical performance:</u>
  - a. Precision/Reproducibility:

Reproducibility within sites was determined using the QC isolates for >95% reproducibility. Between sites was performed at three sites for three days in triplicate for >95% reproducibility on 16 isolates.

- *b. Linearity/assay reportable range:* Not applicable
- c. Traceability (controls, calibrators, or method):

The recommended QC isolate was tested a sufficient number of times with acceptable results with the reference method. The Phoenix results demonstrate that the system can produce QC results in the recommended range.

ORGANISM	conc.	Reference	Phoenix
E. coli	2	3	383
ATCC 25922	4	318	1
Expected Range:	8	63	3
$2-8 \ \mu g/mL$	>32	1	

Inoculum density control: The organism suspension density of the ID broth was equivalent to a 0.5 McFarland standard using the BBL<sup>TM</sup> CrystalSpec<sup>TM</sup> Nephelometer which was verified each day of testing. Internal data was used to demonstrate that the use of the BBL<sup>TM</sup> CrystalSpec<sup>TM</sup> Nephelometer would produce reproducible results. Five different instruments were used.

- *d. Detection limit:*
- Not applicable e. Analytical specificity:
- Not applicable
- f. Assay cut-off: Not applicable
- 2. Comparison studies:
  - *a.* Method comparison with predicate device: The NCCLS recommended broth dilution reference panel was prepared according to the NCCLS recommendation. Clinical testing was performed at six sites. The testing included both fresh clinical isolates and stock isolates along with a challenge set with known results. The test device had a growth rate of >99%. A comparison was provided to the reference method with the following agreement.

	EA	EA	EA	Eval	Eval	Eval	CA	CA	#R	min	maj	vmj
	Tot	Ν	%	EA Tot	EA N	EA %	Ν	%			_	_
Clinical	802	788	98.3	342	331	96.8	790	98.5	398	7	3	2
Challenge	36	34	94.4	5	5	100.0	34	94.4	29	1	0	1
Combined	838	822	98.1	347	336	96.8	824	98.3	427	8	3	3

EA-Essential Agreement CA-Category Agreement R-resistant isolates **maj**-major discrepancies **vmj**-very major discrepancies **min**- minor discrepancies

Essential agreement (EA) is when the BD Phoenix<sup>TM</sup> panels agree with the reference test panel results exactly or within one doubling

dilution of the reference method. Category agreement (CA) is when the BD Phoenix<sup>TM</sup> panel result interpretation agrees exactly with the reference panel result interpretation.

- *b. Matrix comparison:* Not applicable
- 3. Clinical studies:
  - *a. Clinical sensitivity:* Not applicable
  - *b. Clinical specificity:* Not applicable
  - *c. Other clinical supportive data (when a and b are not applicable):* Not applicable
- 4. <u>Clinical cut-off:</u> Not applicable
- 5. Expected values/Reference range:

 $\leq 8(S), 16(I), \geq 32(R)$ 

The expected value range, interpretative criteria and QC are the same as recommended in NCCLS. All values will be included in the package insert.

## M. Conclusion:

This demonstrates acceptable performance as described in the FDA guidance document, "Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems; Guidance for Industry and FDA" and therefore the testing of ampicillin on the BD Phoenix<sup>TM</sup> Automated Microbiology System is substantially equivalent to other commercial devices such as bioMerieux VITEK® AST panels.