



**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY  
ASSAY ONLY**

**I Background Information:**

**A 510(k) Number**

K211748

**B Applicant**

Selux Diagnostics, Inc

**C Proprietary and Established Names**

Selux AST System; Model AST Gen 1.0

**D Regulatory Information**

Product Code(s)	Classification	Regulation Section	Panel
LON	Class II	21 CFR 866.1645 - Fully Automated Short-Term Incubation Cycle Antimicrobial Susceptibility System	MI - Microbiology
LTT	Class II	21 CFR 866.1640 - Antimicrobial susceptibility test powder	MI - Microbiology
LTW	Class II	21 CFR 866.1640 - Antimicrobial susceptibility test powder	MI - Microbiology

**II Submission/Device Overview:**

**A Purpose for Submission:**

To obtain substantial equivalence determination for Gram-negative organisms tested with the Selux AST System to determine the minimum inhibitory concentration of specific antimicrobials with specific Gram-negative organisms.

**B Measurand:**

**Antimicrobial**

Amikacin  
Amoxicillin-Clavulanate  
Ampicillin  
Ampicillin-Sulbactam

**Reportable Range**

≤0.12 to ≥256 µg/mL  
≤0.5 to ≥128 µg/mL  
≤0.25 to ≥128 µg/mL  
≤0.5 to ≥128 µg/mL

<b>Antimicrobial</b>	<b>Reportable Range</b>
Aztreonam	≤0.03 to ≥128 µg/mL
Cefazolin	≤0.12 to ≥128 µg/mL
Cefepime	≤0.25 to ≥128 µg/mL
Cefoxitin	≤1 to ≥128 µg/mL
Ceftazidime	≤0.25 to ≥64 µg/mL
Ceftazidime-Avibactam	≤0.12 to ≥64 µg/mL
Ceftriaxone	≤0.25 to ≥32 µg/mL
Ciprofloxacin	≤0.03 to ≥16 µg/mL
Eravacycline	≤0.016 to ≥4 µg/mL
Ertapenem	≤0.03 to ≥16 µg/mL
Gentamicin	≤0.06 to ≥64 µg/mL
Imipenem-Relebactam	≤0.03 to ≥128 µg/mL
Levofloxacin	≤0.06 to ≥32 µg/mL
Meropenem	≤0.12 to ≥64 µg/mL
Meropenem-Vaborbactam	≤0.06 to ≥64 µg/mL
Minocycline	≤0.25 to ≥64 µg/mL
Piperacillin-Tazobactam	≤0.25 to ≥512 µg/mL
Tobramycin	≤0.12 to ≥128 µg/mL
Trimethoprim-Sulfamethoxazole	≤0.12 to ≥32 µg/mL

### C Type of Test:

Quantitative antimicrobial susceptibility test (AST) system that utilizes colorimetric, oxidation-reduction and growth-based strategies to determine the minimum inhibitory concentration (MIC) of specific antimicrobial-organism combinations.

## III Intended Use/Indications for Use:

### A Intended Use(s):

The Selux AST System is intended to be used for the automated quantitative or qualitative susceptibility testing for most clinically significant aerobic microorganisms. The Selux AST System does not provide organism identification.

### B Indication(s) for Use:

The Selux Gram-Negative Panel is intended for use with the Selux AST System as an *in vitro* test to determine the susceptibility of isolated colonies of specific gram-negative bacilli to specific antimicrobial agents when used as instructed.

The Selux Gram-Negative Panel is a quantitative test for the following antimicrobial agents with the specific organisms identified below:

- Amikacin: *Pseudomonas aeruginosa*
- Amoxicillin-Clavulanate: *Escherichia coli*, *Klebsiella* species (including *K. oxytoca*, *K. pneumoniae*), *Proteus mirabilis*
- Ampicillin: *Escherichia coli*, *Proteus mirabilis*
- Ampicillin-Sulbactam: *Acinetobacter baumannii* complex, *Escherichia coli*, *Klebsiella* species (including *K. oxytoca*, *K. pneumoniae*), *Proteus mirabilis*, *Proteus vulgaris*
- Aztreonam: *Escherichia coli*

- Cefazolin: *Escherichia coli*, *Klebsiella pneumoniae*
- Cefepime: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*
- Cefoxitin: *Escherichia coli*, *Klebsiella* species (including *K. oxytoca*, *K. pneumoniae*), *Morganella morganii*
- Ceftazidime: *Citrobacter* species (including *C. freundii* complex, *C. koseri*), *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella* species (including *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*), *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*
- Ceftazidime-Avibactam: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Serratia marcescens*
- Ceftriaxone: *Citrobacter freundii* complex, *Citrobacter koseri*, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Proteus mirabilis*
- Ciprofloxacin: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*
- Eravacycline: *Citrobacter freundii* complex, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella oxytoca*
- Ertapenem: *Citrobacter freundii* complex, *Citrobacter koseri*, *Escherichia coli*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*
- Gentamicin: *Citrobacter* species (including *C. freundii* complex, *C. koseri*), *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella* species (including *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*), *Proteus* species (including *P. mirabilis*, *P. vulgaris*), *Pseudomonas aeruginosa*, *Serratia marcescens*
- Imipenem-Relebactam: *Citrobacter freundii* complex, *Citrobacter koseri*, *Escherichia coli*, *Klebsiella oxytoca*, *Pseudomonas aeruginosa*
- Levofloxacin: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*
- Meropenem: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Serratia marcescens*
- Meropenem-Vaborbactam: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morganii*, *Serratia marcescens*
- Minocycline: *Escherichia coli*, *Klebsiella* species (including *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*)
- Piperacillin-Tazobactam: *Citrobacter koseri*, *Escherichia coli*, *Klebsiella pneumoniae*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*
- Tobramycin: *Pseudomonas aeruginosa*
- Trimethoprim-Sulfamethoxazole: *Enterobacter cloacae* complex, *Klebsiella* species (including *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*)

### C Special Conditions for Use Statement(s):

- Rx - For Prescription Use Only

The following limitations were added to the device labeling based on performance demonstrated in the current submission:

- The ability of the Selux AST system to detect resistance in the following antimicrobial/organism combinations is unknown because an insufficient number of resistant isolates were available at the time of comparative testing:
  - Ampicillin-Sulbactam: *Proteus vulgaris*
  - Cefepime: *Klebsiella aerogenes*, *Morganella morganii*, *Proteus vulgaris*, *Serratia marcescens*
  - Cefoxitin: *Klebsiella oxytoca*
  - Ceftazidime-avibactam: *Citrobacter koseri*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*
  - Ciprofloxacin: *Citrobacter koseri*, *Proteus vulgaris*
  - Eravacycline: *Escherichia coli*, *Klebsiella oxytoca*
  - Ertapenem: *Citrobacter koseri*, *Klebsiella oxytoca*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*
  - Gentamicin: *Klebsiella aerogenes*, *Morganella morganii*, *Proteus vulgaris*, *Serratia marcescens*
  - Imipenem-Relebactam: *Citrobacter koseri*, *Klebsiella aerogenes*, *Klebsiella oxytoca*
  - Levofloxacin: *Citrobacter koseri*
  - Meropenem: *Citrobacter koseri*, *Klebsiella oxytoca*, *Morganella morganii*, *Proteus vulgaris*
  - Meropenem-vaborbactam: *Citrobacter koseri*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Morganella morganii*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*
  - Minocycline: *Klebsiella oxytoca*
  - Piperacillin-tazobactam: *Citrobacter koseri*, *Proteus mirabilis*, *Proteus vulgaris*
- Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):
  - Amikacin: *Acinetobacter baumannii*
  - Ampicillin: *Proteus mirabilis* when the Selux MIC is 2 µg/mL due to one very major error
  - Ampicillin-Sulbactam: *Proteus vulgaris* when the Selux MIC is 32 µg/mL due to 3 major errors
  - Aztreonam: *Pseudomonas aeruginosa*
  - Cefazolin: *Proteus mirabilis*
  - Cefepime: *Klebsiella aerogenes* when the Selux MIC is  $\geq 64$  µg/mL due to 2 major errors, *Pseudomonas aeruginosa*
  - Cefoxitin: *Proteus mirabilis*, *Proteus vulgaris*
  - Ceftazidime: *Acinetobacter baumannii*, *Pseudomonas aeruginosa*
  - Ceftriaxone: *Enterobacter cloacae*, *Serratia marcescens*
  - Ciprofloxacin: *Serratia marcescens*
  - Eravacycline: *Klebsiella pneumoniae*
  - Ertapenem: *Citrobacter freundii* complex when the Selux MIC is 2, 4, or 8 µg/mL due to 3 major errors, *Enterobacter cloacae*, *Klebsiella aerogenes*
  - Gentamicin: *Pseudomonas aeruginosa* when the Selux MIC is  $\leq 0.06$  µg/mL

- Imipenem: *Acinetobacter baumannii*
- Imipenem-Relebactam: *Klebsiella oxytoca* when the Selux MIC is 4 µg/mL due to one major error, *Klebsiella aerogenes*, *Klebsiella pneumoniae*
- Levofloxacin: *Citrobacter freundii* complex when the Selux MIC is 2 µg/mL due to 2 major errors, *Klebsiella aerogenes* when the Selux MIC is 2 µg/mL due to one major error, *Pseudomonas aeruginosa*
- Meropenem: *Acinetobacter baumannii*, *Citrobacter freundii* complex when the Selux MIC is 16 or 32 µg/mL due to 3 major errors, *Serratia marcescens* when the Selux MIC is 16 or 32 µg/mL due to 4 major errors
- Meropenem-Vaborbactam: *Proteus mirabilis*, *Proteus vulgaris*
- Minocycline: *Acinetobacter baumannii*, *Escherichia coli*, *Klebsiella pneumoniae* when the Selux MIC is 4 µg/mL due to 2 very major errors
- Piperacillin-Tazobactam: *Acinetobacter baumannii*, *Citrobacter freundii*, *Enterobacter cloacae*, *Klebsiella aerogenes*, *Pseudomonas aeruginosa*, *Serratia marcescens*
- Trimethoprim-Sulfamethoxazole: *Citrobacter freundii*, *Escherichia coli*

#### D Special Instrument Requirements:

Selux AST System, Software version 1, subversion 8.168

### IV Device/System Characteristics:

#### A Device Description:

The Selux AST System is an antimicrobial susceptibility test (AST) system that consists of a Sample Prep Station, an Inoculator, an Analyzer, a computer workstation, and the reagents and consumables required to perform AST testing. The system is operated via software that guides users through the manual sample preparation process and operates the automated Inoculator and Analyzer. The software includes an algorithm that enables the system to determine the susceptibility of an organism to the variety of antimicrobials provided in the Selux AST panels. The system is designed so that only Gram stain information is required to initiate testing to select the gram-negative or gram-positive antimicrobial panel. While complete system testing can be performed without species-level identification, this information is required for the system to report susceptibility results. Species identification information can be either manually input to the Selux system or automatically downloaded from the laboratory information system (LIS) at any time, once the sample ID is entered into the LIS.

The Selux AST Systems provides the following consumables: AST panels, panel lids, tubes and caps for sample prep, troughs for inoculator, McFarland Standards for densitometer performance checks, inoculator pipette tips, inoculator waste bins, inoculator growth media, inoculator cleaning solution, analyzer reagent packs, analyzer liquid waste disinfectant tablets, analyzer absorbent waste pads, and an analyzer solid waste bag.

The Selux AST System is divided into three primary AST steps performed at three different stations: workbench (for sample preparation), panel inoculator (for inoculation), and AST analyzer (for antimicrobial susceptibility testing), as described in detail in the Device Description section of [K211759](#).

The Selux AST System with the Gram-Negative Panel can determine the MIC of various antimicrobials when tested against specific organisms **Table 1**.

**Table 1.** Reportable MIC Ranges and Organism-Specific Breakpoints for Antimicrobials included in the Selux AST System, Gram-Negative Panel

Antimicrobial	Organism Group	Selux AST System Reportable Range (µg/mL)	FDA-Recognized/ Approved Breakpoints (µg/mL)		
			S	I	R
Amikacin	<i>Pseudomonas aeruginosa</i>	≤0.12 to ≥256	≤16	32	≥64
Amoxicillin-Clavulanate	Enterobacterales	≤0.5 to ≥128	≤8	16	≥32
Ampicillin	Enterobacterales	≤0.25 to ≥128	≤8	16	≥32
Ampicillin-Sulbactam	<i>Acinetobacter baumannii</i> complex Enterobacterales	≤0.5 to ≥128	≤8	16	≥32
Aztreonam	Enterobacterales	≤0.03 to ≥128	≤4	8	≥16
Cefazolin	Enterobacterales	≤0.12 to ≥128	≤1	2	≥4
Cefepime	Enterobacterales	≤0.25 to ≥128	≤2	4-8	≥16
Cefoxitin	Enterobacterales	≤1 to ≥128	≤4	8	≥16
Ceftazidime	Enterobacterales	≤0.25 to ≥64	≤4	8	≥16
Ceftazidime-Avibactam	Enterobacterales <i>Pseudomonas aeruginosa</i>	≤0.12 to ≥64	≤8	-	≥16
Ceftriaxone	Enterobacterales	≤0.25 to ≥32	≤1	2	≥4
Ciprofloxacin	Enterobacterales	≤0.03 to ≥16	≤0.25	0.5	≥1
	<i>Pseudomonas aeruginosa</i>		≤0.5	1	≥2
Eravacycline	Enterobacterales	≤0.016 to ≥4	≤0.5	-	-
Ertapenem	Enterobacterales	≤0.03 to ≥16	≤0.5	1	≥2
Gentamicin	Enterobacterales	≤0.06 to ≥64	≤4	8	≥16
	<i>Pseudomonas aeruginosa</i>				
Imipenem-Relebactam	Enterobacterales	≤0.03 to ≥128	≤1	2	≥4
	<i>Pseudomonas aeruginosa</i>		≤2	4	≥8
Levofloxacin	Enterobacterales	≤0.06 to ≥32	≤0.5	1	≥2
Meropenem	Enterobacterales	≤0.12 to ≥64	≤1	2	≥4
	<i>Pseudomonas aeruginosa</i>		≤2	4	≥8
Meropenem-Vaborbactam	Enterobacterales	≤0.06 to ≥64	≤4	8	≥16
Minocycline	Enterobacterales	≤0.25 to ≥64	≤4	8	≥16
Piperacillin-Tazobactam	Enterobacterales	≤0.25 to ≥512	≤16	32-64	≥128
Tobramycin	<i>Pseudomonas aeruginosa</i>	≤0.12 to ≥128	≤4	8	≥16
Trimethoprim-Sulfamethoxazole	Enterobacterales	≤0.12 to ≥32	≤2	-	≥4

S, Susceptible; I, Intermediate; R, Resistant; -, no breakpoint (interpretive criterion) recognized

### B Principle of Operation:

The Selux AST test requires that the Gram type of the organism be known prior to testing as the information is necessary to select the proper AST panel to use (i.e., gram-positive or gram-negative). The organism identification (ID) is not needed for Selux AST processing to be performed; however, the organism ID is necessary for a final result to be reported.

The Selux AST System performs antimicrobial susceptibility testing similar to the gold-standard broth microdilution method. To get an accurate reading of microbial growth, the sufficient growth assay monitors growth in dedicated AST panel wells that contain organisms and cation-adjusted Mueller-Hinton Broth but no antimicrobials or probes. Sufficient growth assay wells are monitored by fluorescence to those wells which the standard metabolism-based viability assay pair resazurin/methylene blue have been added and/or by optical absorbance. Two probe-based assays commence across all wells in the panel after the threshold in the sufficient growth well has been met. The “viability assay” uses a fluorescent metabolic probe as an indicator of cellular activity and the “surface area assay” uses a fluorescent cationic probe that binds microorganism surfaces as a proxy for cellular surface area. The viability and surface area assays are performed in each AST panel well, providing two complementary datasets for each well. These data are input to an MIC-determining algorithm that provides results when organism IDs are available.

Refer to the [K211759](#) Decision Summary for additional information on the Selux AST System.

**V Substantial Equivalence Information:**

**A Predicate Device Name(s):**

BD Phoenix Automated Microbiology System - GN Ceftaroline (0.0156-4 µg/mL)

**B Predicate 510(k) Number(s):**

K190905

**C Comparison with Predicate(s):**

Device & Predicate Device(s):	Device: <u>K211748</u>	Predicate: <u>K190905</u>
Device Trade Name	Selux AST System	BD Phoenix Automated Microbiology System – GN Ceftaroline (0.0156-4 µg/mL)
<b>General Device Characteristic Similarities</b>		
Intended Use/Indications For Use	<p>The Selux AST System is intended to be used for the automated quantitative or qualitative susceptibility testing for most clinically significant aerobic microorganisms. The Selux AST System does not provide organism identification.</p> <p>The Selux Gram-Negative Panel is intended for use with the Selux AST System as an <i>in vitro</i> test to determine the susceptibility of isolated colonies of specific gram-negative bacilli to specific antimicrobial agents when used as instructed.*</p>	<p>The BD Phoenix Automated Microbiology System is intended for <i>in vitro</i> quantitative determination of antimicrobial susceptibility by minimal inhibitory concentration (MIC) of most Gram-negative aerobic and facultative anaerobic bacteria isolates from pure culture for <i>Enterobacteriaceae</i> and Non-<i>Enterobacteriaceae</i> and most Gram-positive bacteria isolates from pure culture belonging to the genera <i>Staphylococcus</i>, <i>Enterococcus</i> and <i>Streptococcus</i>.</p>
Source of Microorganisms	Bacterial colonies isolated from culture	Same
Technology	Automated growth-based detection using metabolic indicators to detect organism growth	Automated growth based enhanced by use of a redox indicator (colorimetric oxidation-reduction) to detect organism growth
Panel Type	Gram-negative	Same

<b>Device &amp; Predicate Device(s):</b>	<b>Device: <u>K211748</u></b>	<b>Predicate: <u>K190905</u></b>
Read Method	Automated	Same
Inoculation Method	Automated	Same
Result Reported	Report results as minimum inhibitory concentration (MIC) and categorical interpretation (S, I, R, NS)	Report results as minimum inhibitory concentration (MIC) and categorical interpretation (S, I, R)
Indicated Species	<i>Acinetobacter baumannii</i> complex, specific Enterobacterales species, and <i>Pseudomonas aeruginosa</i>	Specific Enterobacterales species
<b>General Device Characteristic Differences</b>		
Antimicrobial Agent and Reporting Range	Amikacin ≤0.12 to ≥256 µg/mL Amoxicillin-Clavulanate ≤0.5 to ≥128 µg/mL Ampicillin ≤0.25 to ≥128 µg/mL Ampicillin-Sulbactam ≤0.5 to ≥128 µg/mL Aztreonam ≤0.03 to ≥128 µg/mL Cefazolin ≤0.12 to ≥128 µg/mL Cefepime ≤0.25 to ≥128 µg/mL Cefoxitin ≤1 to ≥128 µg/mL Ceftazidime ≤0.25 to ≥64 µg/mL Ceftazidime-Avibactam ≤0.12 to ≥64 µg/mL Ceftriaxone ≤0.25 to ≥32 µg/mL Ciprofloxacin ≤0.03 to ≥16 µg/mL Eravacycline ≤0.016 to ≥4 µg/mL Ertapenem ≤0.03 to ≥16 µg/mL Gentamicin ≤0.06 to ≥64 µg/mL Imipenem ≤0.016 to ≥64 µg/mL Imipenem-Relebactam ≤0.03 to ≥128 µg/mL Levofloxacin ≤0.06 to ≥32 µg/mL Meropenem ≤0.12 to ≥64 µg/mL Meropenem-Vaborbactam ≤0.06 to ≥64 µg/mL Minocycline ≤0.25 to ≥64 µg/mL Piperacillin-Tazobactam ≤0.25 to ≥512 µg/mL Tobramycin ≤0.12 to ≥128 µg/mL Trimethoprim-Sulfamethoxazole ≤0.12 to ≥32 µg/mL	Ceftaroline 0.0156-4 µg/mL
IVD Functions	AST	ID and AST
Instrument	Selux AST System	BD Phoenix Automated Microbiology System

\* The Selux Gram-Negative Panel is a quantitative test for the antimicrobial agents with the specific organisms listed in section III.B (Indications for Use Statement).

## VI Standards/Guidance Documents Referenced:

- FDA Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems; Guidance for Industry and FDA (Issued August 28, 2009)
- CLSI M100-M30, “Performance Standards for Antimicrobial Susceptibility Testing”; Thirtieth Edition (January 2020)

## VII Performance Characteristics (if/when applicable):

### A Analytical Performance:

#### 1. Precision/Reproducibility:



Reproducibility testing for the Selux AST System with the Gram-Negative Panel was conducted at 3 testing sites (2 external and 1 internal site). Panel members generally consisted of species indicated for use with each respective antimicrobial. To accommodate the numerous antimicrobial/organism combinations being tested concurrently, up to three non-indicated species, within an indicated genus or family, were considered acceptable for testing. Reproducibility was determined from the total number (and percent) of results that fell within one dilution (+/- one doubling dilution) of the modal MIC result for quantitative assays divided by the total number of results. Reproducibility was evaluated between sites (inter-site) and within sites (intra-site). Both best-case (assumes that off-scale results are within one dilution of the mode) and worst-case (assumes that off-scale results are more than one dilution of the mode) performance was determined for each antimicrobial, as outlined in the AST Special Controls Guidance.

In the initial study, inter-site reproducibility was evaluated at three sites by testing at least 25 isolates with on-scale MIC values for each antimicrobial, for a minimum of 75 results per antimicrobial (25 isolates x 3 sites = 75 results/antimicrobial). In general, best-case inter-site reproducibility was acceptable ( $\geq 95\%$ ) and a sufficient number of results obtained ( $\geq 75$  results). However, antimicrobials with inter-site worst-case reproducibility of  $< 89\%$  or an insufficient number of results generated (i.e.,  $< 75$ ) in the initial study were further evaluated in a supplemental study in which testing was performed with three instruments. Data from both studies are collated and summarized in **Table 2**. Performance is summarized for each antimicrobial tested with all organisms as well as indicated organisms only. Inter-site reproducibility was determined to be acceptable.

**Table 2.** Inter-site Reproducibility of Selux AST System (Gram Negative Panel)

Antimicrobial	All organisms (combined)		Indicated organisms only	
	Best-case (%)	Worst-case (%)	Best-case (%)	Worst-case (%)
Amikacin	98/104 (94.2%)	98/104 (94.2%)	92/98 (93.9%)	92/98 (93.9%)
Amoxicillin-Clavulanate	77/78 (98.7%)	76/78 (97.4%)	71/72 (98.6%)	70/72 (97.2%)
Ampicillin	72/75 (96%)	72/75 (96%)	68/69 (98.6%)	68/69 (98.6%)
Ampicillin-Sulbactam	74/75 (98.7%)	74/75 (98.7%)	69/69 (100%)	69/69 (100%)
Aztreonam	77/78 (98.7%)	75/78 (96.2%)	77/78 (98.7%)	75/78 (96.2%)
Cefazolin <sup>a</sup>	77/81 (95.1%)	76/81 (93.8%)	68/72 (94.4%)	67/72 (93.1%)
Cefepime	77/78 (98.7%)	76/78 (97.4%)	77/78 (98.7%)	76/78 (97.4%)
Cefoxitin	73/75 (97.3%)	72/75 (96.0%)	70/72 (97.2%)	69/72 (95.8%)
Ceftazidime	77/81 (95.1%)	76/81 (93.8%)	77/81 (95.1%)	76/81 (93.8%)
Ceftazidime-Avibactam	140/144 (97.2%)	139/144 (96.5%)	140/144 (97.2%)	139/144 (96.5%)
Ceftriaxone	140/141 (99.3%)	140/141 (99.3%)	140/141 (99.3%)	140/141 (99.3%)
Ciprofloxacin <sup>a</sup>	74/75 (98.7%)	71/75 (94.7%)	74/75 (98.7%)	71/75 (94.7%)
Eravacycline	78/78 (100%)	78/78 (100%)	78/78 (100%)	78/78 (100%)
Ertapenem	143/145 (98.6%)	143/145 (98.6%)	143/145 (98.6%)	143/145 (98.6%)
Gentamicin <sup>a</sup>	77/80 (96.3%)	77/80 (96.3%)	77/80 (96.3%)	77/80 (96.3%)
Imipenem	145/149 (97.3%)	140/149 (94.0%)	145/149 (97.3%)	140/149 (94.0%)
Imipenem-Relebactam	71/75 (94.7%)	71/75 (94.7%)	71/75 (94.7%)	71/75 (94.7%)
Levofloxacin	76/78 (97.4%)	76/78 (97.4%)	73/75 (97.3%)	73/75 (97.3%)
Meropenem	75/78 (96.2%)	73/78 (93.6%)	70/72 (97.2%)	68/72 (94.4%)
Meropenem-Vaborbactam	142/145 (97.9%)	132/145 (91.0%)	142/145 (97.9%)	132/145 (91.0%)

Antimicrobial	All organisms (combined)		Indicated organisms only	
	Best-case (%)	Worst-case (%)	Best-case (%)	Worst-case (%)
Minocycline <sup>a</sup>	73/75 (97.3%)	73/75 (97.3%)	64/66 (97.3%)	64/66 (97.3%)
Piperacillin-Tazobactam <sup>a</sup>	72/75 (96.0%)	69/75 (92.0%)	63/66 (95.5%)	60/66 (90.9%)
Tobramycin	74/78 (94.9%)	73/78 (93.6%)	74/78 (94.9%)	73/78 (93.6%)
Trimethoprim-Sulfamethoxazole	143/148 (96.6%)	142/148 (95.9%)	143/148 (96.6%)	142/148 (95.9%)

<sup>a</sup> In instances where the original study had lower than the required number, supplemental testing was conducted and data collated with original data.

Intra-site reproducibility was evaluated by testing a minimum of five isolates in triplicate on three days at one internal site for each antimicrobial to generate a minimum of 45 results per antimicrobial (5 isolates x 3 replicates x 3 days = 45 results/antimicrobial). In general, best-case intra-site reproducibility was acceptable ( $\geq 95\%$ ) and a sufficient number of results obtained ( $\geq 45$  results). However, antimicrobials with worst-case reproducibility  $< 89\%$  or an insufficient number of results generated (i.e.,  $< 45$ ) in the initial study were further evaluated in a supplemental study. Data from both studies are collated and summarized in **Table 3**. Performance is summarized for each antimicrobial tested with all organisms as well as indicated organisms only. Intra-site reproducibility was determined to be acceptable.

**Table 3.** Intra-site Reproducibility of Selux AST System (Gram Negative Panel)

Antimicrobial	All organisms (combined)		Indicated organisms only	
	Best-case (%)	Worst-case (%)	Best-case (%)	Worst-case (%)
Amoxicillin-Clavulanate	63/63 (100%)	63/63 (100%)	63/63 (100%)	63/63 (100%)
Ampicillin	45/45 (100%)	45/45 (100%)	45/45 (100%)	45/45 (100%)
Ampicillin-Sulbactam	72/72 (100%)	72/72 (100%)	72/72 (100%)	72/72 (100%)
Aztreonam	71/74 (95.9%)	68/74 (91.9%)	71/74 (95.9%)	68/74 (91.9%)
Cefazolin	61/63 (96.8%)	61/63 (96.8%)	61/63 (96.8%)	61/63 (96.8%)
Cefepime	45/47 (95.7%)	45/47 (95.7%)	45/47 (95.7%)	45/47 (95.7%)
Cefoxitin	54/54 (100%)	54/54 (100%)	54/54 (100%)	54/54 (100%)
Ceftazidime	92/93 (98.9%)	91/93 (97.8%)	92/93 (98.9%)	91/93 (97.8%)
Ceftazidime-Avibactam	47/47 (100%)	43/47 (91.5%)	47/47 (100%)	43/47 (91.5%)
Eravacycline	103/103 (100%)	99/103 (96.1%)	103/103 (100%)	99/103 (96.1%)
Gentamicin	116/121 (95.9%)	116/121 (95.9%)	116/121 (95.9%)	116/121 (95.9%)
Imipenem-Relebactam	72/76 (94.7%)	72/76 (94.7%)	72/76 (94.7%)	72/76 (94.7%)
Levofloxacin	179/181 (98.9%)	174/181 (96.1%)	179/181 (98.9%)	174/181 (96.1%)
Meropenem	57/58 (98.3%)	56/58 (96.6%)	57/58 (98.3%)	56/58 (96.6%)
Meropenem-Vaborbactam	45/47 (95.7%)	43/47 (91.5%)	45/47 (95.7%)	43/47 (91.5%)
Piperacillin-Tazobactam	54/56 (96.4%)	54/56 (96.4%)	54/56 (96.4%)	54/56 (96.4%)
Trimethoprim-Sulfamethoxazole	62/65 (95.4%)	61/65 (93.8%)	62/65 (95.4%)	61/65 (93.8%)

<sup>a</sup> In instances where the original study had lower than the required number, supplemental testing was conducted and data collated with original data.

2. Linearity:  
Not applicable
3. Analytical Specificity/Interference:  
Not Applicable
4. Assay Reportable Range:  
Not applicable
5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):

**Quality Control Testing.** The Selux AST System software contains QC Sets (A, B, C) such that when the Gram-Negative Panel is selected, the system will indicate the appropriate QC strains to be prepared for testing. CLSI-recommended QC strains and expected QC ranges for use with the Gram-Negative Panel are summarized in **Table 4**, which also includes the designation of the QC set as appropriate for each antimicrobial. Notably, QC Set A comprises gram-positive ATCC QC organisms which is acceptable to use for QC testing of the Gram-Negative Panel provided that QC is run concurrently with Set B and/or C, which comprise gram-negative ATCC QC organisms. This is necessary because the metabolic assay in the Selux AST System utilizes an additional biochemical reagent, an electron transfer agent, exclusively with gram-negative organisms. The following QC reporting scheme is included in the device labeling (and is supported by software information) to ensure that gram-negative organisms will be processed correctly by the Selux AST System:

- *Gram-Negative Panel QC failure when testing gram-positive QC organisms:* clinical results should not be reported for the antimicrobial(s) in which QC failed but may be reported for the antimicrobials for which QC passed.
- *Gram-Negative Panel QC failure when testing gram-negative QC organisms:*
  - If all gram-negative organisms fail QC: all QC results (including those when testing gram-positive organisms) should be considered failed and clinical results should not be reported for any antimicrobial.
  - If one (or fewer than all) gram-negative organisms fail QC: clinical results should not be reported for the antimicrobial(s) in which QC failed but may be reported for the antimicrobials for which QC passed (including those when testing gram-positive organisms).

**Table 4.** CLSI-Recommended QC Organisms for Gram-Negative Panel

Antimicrobial	Abbreviation	QC Set	QC Strain	Expected QC Range (µg/mL)
Amikacin	AMK	A	ATCC 29213 <sup>a</sup>	1-4
Amoxicillin-Clavulanate	AMC	B	ATCC 700603	4-16
Ampicillin	AMP	C	ATCC 25922	2-8
Ampicillin-Sulbactam	SAM	B	ATCC 700603	8-32
Aztreonam	ATM	C	ATCC 25922	0.06-0.25
Cefazolin	CFZ	C	ATCC 25922	1-4
Cefepime	FEP	A	ATCC 29213 <sup>a</sup>	1-4
Cefoxitin	FOX	C	ATCC 25922	2-8
Ceftazidime	CAZ	A	ATCC 29213 <sup>a</sup>	4-16
Ceftazidime-Avibactam	CZA	B	ATCC 700603	0.25-2

Antimicrobial	Abbreviation	QC Set	QC Strain	Expected QC Range (µg/mL)
Ceftriaxone	CRO	A	ATCC 29213 <sup>a</sup>	1-8
Ciprofloxacin	CIP	A	ATCC 29213 <sup>a</sup>	0.12-0.5
Eravacycline	ERV	C	ATCC 25922	0.016-0.12
Ertapenem	ETP	A	ATCC 29213 <sup>a</sup>	0.06-0.25
Gentamicin	GEN	A	ATCC 29212 <sup>a</sup>	4-16
Imipenem	IMP	B	ATCC 700603	0.06-0.5
Imipenem-Relebactam	IMR	C	ATCC BAA-2814	0.06-0.5
Levofloxacin	LVX	A	ATCC 29212 <sup>a</sup>	0.25-2
Meropenem	MEM	A	ATCC 29212 <sup>a</sup>	2-8
Meropenem-Vaborbactam	MEV	C	ATCC BAA-2814	0.12-0.5
Minocycline	MIN	A	ATCC 29212 <sup>a</sup>	1-4
Piperacillin-Tazobactam	TZP	B	ATCC 700603	8-32
Tobramycin	TOB	A	ATCC 29212 <sup>a</sup>	8-32
Trimethoprim-Sulfamethoxazole	SXT	C	ATCC 25922	0.12-0.5 <sup>b</sup>

<sup>a</sup> Set A comprises gram-positive ATCC QC organisms which is acceptable to use for QC testing of the Gram-Negative Comprehensive panel provided that QC is run concurrently with Set B and/or C, which comprise gram-negative ATCC QC organisms. This is necessary because the metabolic assay in the Selux AST System utilizes an additional biochemical reagent, an electron transfer agent, exclusively with gram-negative organisms. This ensures that gram-negative organisms will be processed correctly by the Selux AST System.

<sup>b</sup> Expected QC range for SXT in CLSI M100 is ≤0.5; Selux AST reporting range for SXT is ≤0.12 to ≥32 µg/mL. Selux MIC results of ≤0.5 µg/mL were considered acceptable.

Quality control testing was performed each day that testing was conducted. CLSI-recommended QC strains for each antimicrobial (described in **Table 4**) were tested a sufficient number of times (i.e., at least 20 times/site) at each testing site using the Selux AST System. QC testing of the broth microdilution reference method was performed at a single, central site. In some cases, the site used CLSI-recommended QC organisms that were different from QC organisms tested with the Selux AST System. This is noted in **Table 5**.

QC expected ranges and results for the Selux AST System with Gram-Negative panel are summarized in **Table 5**. For all antimicrobials, greater than 95% of results were within the expected range, which is acceptable.

*Due to an administrative error, reference data for minocycline and tobramycin was not previously included in Table 5. This data is included in the updated table for both minocycline and tobramycin testing of E. coli ATCC 25922.*

**Table 5.** QC Results for the Selux AST System Assay with Gram Negative Panel

Antimicrobial	QC Strain	Expected Range (µg/mL)	No. in Range (%)	
			Reference	Selux
Amikacin	<i>S. aureus</i> ATCC 29213	1-4	56/56 (100)	218/219 (99.5)
Amoxicillin	<i>K. pneumoniae</i> ATCC 700603 <sup>a</sup>	≥128 <sup>c</sup>	47/47 (100)	219/221 (99.1)
Amoxicillin-Clavulanate	<i>K. pneumoniae</i> ATCC 700603	4-16	49/50 (98.0)	211/221 (95.5)
Ampicillin	<i>E. coli</i> ATCC 25922	2-8	47/48 (97.9)	210/210 (100)
	<i>K. pneumoniae</i> ATCC 700603 <sup>a</sup>	≥128 <sup>c</sup>	50/50 (100)	219/221 (99.1)
Ampicillin-Sulbactam	<i>K. pneumoniae</i> ATCC 700603	8-32	46/47 (97.9)	221/221 (100)
Aztreonam	<i>E. coli</i> ATCC 25922	0.06-0.25	56/57 (98.2)	210/210 (100)

Antimicrobial	QC Strain	Expected Range (µg/mL)	No. in Range (%)	
			Reference	Selux
Cefazolin	<i>E. coli</i> ATCC 25922	1-4	53/53 (100)	210/210 (100)
Cefepime	<i>S. aureus</i> ATCC 29213	1-4	54/55 (98.2)	218/219 (99.5)
Cefoxitin	<i>E. coli</i> ATCC 25922	2-8	45/46 (97.8)	209/210 (99.5)
Ceftazidime	<i>S. aureus</i> ATCC 29213	4-16	NT	217/219 (99.1)
	<i>E. coli</i> ATCC 25922 <sup>b</sup>	0.06-0.5	56/57 (98.2)	NT
	<i>K. pneumoniae</i> ATCC 700603 <sup>a</sup>	16-64	66/66 (100)	218/221 (98.6)
Ceftazidime-Avibactam	<i>K. pneumoniae</i> ATCC 700603	0.25-2	59/59 (100)	221/221 (100)
Ceftriaxone	<i>S. aureus</i> ATCC 29213	1-8	45/49 (97.8)	218/219 (99.5)
Ciprofloxacin	<i>S. aureus</i> ATCC 29213	0.12-0.5	57/58 (98.3)	218/219 (99.5)
Eravacycline	<i>E. coli</i> ATCC 25922	0.016-0.12	49/49 (100)	210/210 (100)
Ertapenem	<i>S. aureus</i> ATCC 29213	0.06-0.25	51/51 (100)	217/219 (99.1)
Gentamicin	<i>E. faecalis</i> ATCC 29212	4-16	51/53 (96.2)	216/217 (99.5)
Imipenem	<i>K. pneumoniae</i> ATCC 700603	0.06-0.5	64/64 (100)	221/221 (100)
	<i>K. pneumoniae</i> ATCC BAA-2814 <sup>a</sup>	16-64	65/65 (100)	210/210 (100)
Imipenem-Relebactam	<i>K. pneumoniae</i> ATCC BAA-2814	0.06-0.5	65/65 (100)	207/210 (98.6)
Levofloxacin	<i>E. faecalis</i> ATCC 29212	0.25-2	38/40 (95.0)	89/90 (98.9)
Meropenem	<i>E. faecalis</i> ATCC 29212	2-8	55/56 (98.2)	217/217 (100)
Meropenem-Vaborbactam	<i>K. pneumoniae</i> ATCC BAA-2814	0.12-0.5	48/49 (98.0)	204/210 (97.1)
Minocycline	<i>E. faecalis</i> ATCC 29212	1-4	NT	216/217 (99.5)
	<i>E. coli</i> ATCC 25922 <sup>b</sup>	0.25-1	43/44 (97.7)	NT
Piperacillin-Tazobactam	<i>K. pneumoniae</i> ATCC 700603	8-32	66/66 (100)	221/221 (100)
Tobramycin	<i>E. faecalis</i> ATCC 29212	8-32	NT	211/217 (97.2)
	<i>E. coli</i> ATCC 25922 <sup>b</sup>	0.25-1	52/53 (98.1)	NT
Trimethoprim-Sulfamethoxazole	<i>E. coli</i> ATCC 25922	0.12-0.5 <sup>b</sup>	44/100 (100)	208/210 (99.0)

NT: Not tested

<sup>a</sup> Tested to confirm the integrity of the QC strain for testing with the beta-lactam/beta-lactam-inhibitor combination antimicrobial.

<sup>b</sup> BMD reference method testing was performed at a single, central site which utilized their standard, CLSI-recommended QC organisms.

<sup>c</sup> Expected QC range in CLSI M100 is >128; highest concentration of Amoxicillin and Ampicillin on the Selux AST result panel is 128 µg/mL. Selux MIC results of ≥128 µg/mL were considered acceptable.

**Inoculum Density Check:** Refer to the [K211759](#) decision summary for description and results since studies were performed concurrently.

**Device Failure:** Refer to the [K211759](#) decision summary for description and results since studies were performed concurrently.

**Growth Failure Rate:** There were 13 growth failures (13/1401 tests = 0.93%) when using the Gram-Negative Panel on the Selux AST System. All samples that generated Selux results also generated reference method results.

**Purity Check:** Purity plates were prepared from the inoculum suspensions of every sample tested. AST results were only reported for pure isolates; data generated from plates that generated multiple colony morphologies was excluded from analyses.

## 6. Detection Limit:

Not Applicable

7. Assay Cut-Off:  
Not Applicable

## **B Comparison Studies:**

1. Method Comparison with Predicate Device:

Clinical performance testing on the Selux AST System was initially performed at three U.S. test sites (2 external, 1 internal). For instances in which testing was required to supplement existing data from the original study and support specific claims, testing was performed on three Selux AST System instruments at two testing sites.

Performance was evaluated using contemporary and stock frozen clinical isolates as well as challenge isolates, which were selected for their resistance profiles. Clinical isolates of non-fastidious bacteria were recovered from cultures of clinical specimens (e.g., blood, stool, urine, respiratory, wound, aspirates) from diverse geographic locations across the U.S.

Contemporary frozen isolates were defined as isolates that had been collected and frozen within six months of testing while stock frozen isolates were tested six or more months after collection.

A total of 1401 clinical (426 contemporary and 975 stock) and 222 challenge isolates from 12 Enterobacterales species, *Acinetobacter baumannii* complex, and *Pseudomonas aeruginosa* were tested to evaluate the Selux AST System performance for 24 antimicrobials using the Selux Gram-Negative panel. Depending on the spectrum of activity, breakpoints and the claimed organisms (species/group) for each antimicrobial on the panel, the number of datapoints for the various antimicrobial-organisms tested varied and ranged from 165 (e.g. *P. aeruginosa*/Imipenem-Relebactam) to 977 (e.g. Enterobacterales/Ertapenem). Selux results were compared to the modal value of triplicate broth microdilution reference results performed at an independent reference laboratory.

Performance was determined generally based on criteria outlined in the Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems including essential agreement (EA), categorical agreement (CA), and categorical errors (minor, major and very major errors). EA was calculated as the percentage of Selux MIC results that were within plus or minus one serial two-fold dilution of the reference result. CA was calculated as the percentage of Selux interpretive results (S/I/R) that were identical to the interpretive results of the reference result. EA of evaluable results (on-scale Selux and reference results or results in which an off-scale result was at least two doubling dilutions from the on-scale result) were also calculated. Performance was considered acceptable if the EA and CA were  $\geq 90\%$ , major error rate was  $\leq 3\%$ , and very major error rate was  $\leq 2\%$ .

A trending analysis using combined clinical and challenge isolate results was also conducted to evaluate antimicrobial-organism combinations for which Selux MIC results were determined to be one or more doubling dilutions lower or higher than the reference result. MIC results that were off-scale for both the reference and Selux were not considered in the trending analysis. Antimicrobial-organism combinations for which the difference between the percentage of isolates with higher or lower MIC values was  $\geq 30\%$  with a statistically

significant confidence interval were considered to have evidence of trending and is addressed in the labeling.

A high-level summary of the Selux AST System performance using the Gram-Negative Panel is described below for each antimicrobial and indicated species. Complete details and results including EA, CA and error rate analyses are summarized in **Table 6** and trending analyses are summarized in **Table 7**.

**Amikacin.** A total of 111 *A. baumannii* isolates were evaluated with amikacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 55.0% and CA of 86.5%, which is not acceptable. There were 6 minor, 8 major ( $8/80 = 10\%$ ) and 1 very major ( $1/25 = 4\%$ ) errors, which is not acceptable. Due to the unacceptable performance for *A. baumannii*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *A. baumannii* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Amikacin: *Acinetobacter baumannii*

A total of 165 *P. aeruginosa* isolates were evaluated with amikacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 90.9% and CA of 98.2%. There were 3 minor, 0 major and 0 very major errors. Overall, performance is acceptable.

Analysis of trending indicated that MIC values for *P. aeruginosa* tended to be one doubling dilution lower than the reference MIC value. The following statement is included in the device labeling as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Amikacin: *Pseudomonas aeruginosa*

**Amoxicillin-Clavulanate.** A total of 593 Enterobacterales isolates (indicated species: 179 *E. coli*, 48 *K. oxytoca*, 142 *K. pneumoniae*, 88 *P. mirabilis*; non-indicated species: 74 *C. koseri*, 62 *P. vulgaris*) were evaluated with amoxicillin-clavulanate. The combined results from clinical and challenge isolate testing demonstrated an EA of 93.3% and CA of 92.8%. There were 37 minor, 5 major and 1 very major errors. When evaluating results by individual species, *E. coli* had a CA of 88.7% due to 17 minor, 1 major and 0 very major errors which was determined to be acceptable since the majority of categorical errors were minor and the EA of evaluable results was good (96.5%). *P. mirabilis* had one very major error ( $1/8 = 16.7\%$ ), which was considered a random error due to the limited number of resistant isolates tested.

Analysis of trending indicated that MIC values for *K. pneumoniae* tended to be one doubling dilution lower than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Amoxicillin-Clavulanate: *Klebsiella pneumoniae*

**Ampicillin.** A total of 254 Enterobacterales isolates (165 *E. coli*, 89 *P. mirabilis*) were evaluated with ampicillin. The combined results from clinical and challenge isolate testing demonstrated an EA of 94.9% and CA of 98.8%. There was 1 minor, 1 major and 1 very major error. When evaluating results by individual species, *P. mirabilis* had 1 major error (1/60 = 1.7%) and 1 very major error (1/26 = 3.8%) at a Selux MIC value of 2 µg/mL which represents 4.5% (4/89) of results with this drug/organism combination. The following limitation is included in the device labeling to address the very major error:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Ampicillin: *Proteus mirabilis* when the Selux MIC is 2 µg/mL due to one very major error

Analysis of trending indicated that MIC values for *E. coli* tended to be one doubling dilution lower than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Ampicillin: *Escherichia coli*

**Ampicillin-Sulbactam.** A total of 123 *A. baumannii* isolates were evaluated with ampicillin-sulbactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 91.9% and CA of 92.7%. There were 6 minor, 2 major and 1 very major error. Overall, performance is acceptable.

A total of 574 Enterobacterales isolates (204 *E. coli*, 29 *K. oxytoca*, 116 *K. pneumoniae*, 75 *M. morgani*, 88 *P. mirabilis*, 62 *P. vulgaris*) were evaluated with ampicillin-sulbactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 96.9% and CA of 83.5%. There were 90 minor, 5 major and 0 very major errors. When evaluating results by individual species, *K. oxytoca* and *P. mirabilis* each had CA <90% which was considered acceptable since all of the categorical errors were minor and the EA of evaluable results was good (>95%). *P. vulgaris* had a CA of 45.2% due to 31 minor errors as well as 3 major errors (3/43 = 7.0%) at a Selux MIC value of 32 µg/mL which represents 9.7% (6/62) of results with this drug/organism combination. *M. morgani* had a CA of 72% due to 19 minor and 2 major errors (2/4 = 50%), which is not acceptable. Due to the unacceptable performance for *M. morgani*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling address the *P. vulgaris* major errors and restrict reporting of *M. morgani* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):



- Ampicillin-Sulbactam: *Morganella morganii*; *Proteus vulgaris* when the Selux MIC is 32 µg/mL due to 3 major errors

A limitation statement is included in the device labeling to address the lack of testing with resistant *P. vulgaris* isolates.

Analysis of trending indicated that MIC values for *K. oxytoca* tended to be one doubling dilution lower than the reference MIC value while MIC values for *P. vulgaris* tended to be one doubling dilution higher than the reference MIC value. The following statements are included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Ampicillin-Sulbactam: *Klebsiella oxytoca*

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Ampicillin-Sulbactam: *Proteus vulgaris*

**Aztreonam.** A total of 183 Enterobacterales isolates (183 *E. coli*) were evaluated with aztreonam. The combined results from clinical and challenge isolate testing demonstrated an EA of 97.8% and CA of 97.3%. There were 5 minor, 0 major and 0 very major errors. Overall, performance is acceptable.

A total of 219 *P. aeruginosa* isolates were evaluated with aztreonam. The combined results from clinical and challenge isolate testing demonstrated an EA of 91.3% and CA of 87.7%. There were 20 minor, 3 major (3/158 = 1.9%) and 4 very major (4/48 = 8.3%) errors. Due to the unacceptable performance for *P. aeruginosa*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *P. aeruginosa* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Aztreonam: *Pseudomonas aeruginosa*

**Cefazolin.** A total of 412 Enterobacterales isolates (indicated species: 177 *E. coli*, 115 *K. pneumoniae*, 66 *P. mirabilis*; non-indicated species: 54 *C. koseri*) were evaluated with cefazolin. The combined results from clinical and challenge isolate testing demonstrated an EA of 94.2% and CA of 74.8%. There were 102 minor, 1 major and 1 very major errors. When evaluating results by individual species, *E. coli* and *K. pneumoniae* each had CA <90% which was considered acceptable since most of the categorical errors were minor and the EA of evaluable results was good (>90%). *P. mirabilis* had an EA of 86.4% and an EA of evaluable results of 84.8% (50/59), which is not acceptable. Due to the unacceptable performance for *P. mirabilis*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *P. mirabilis* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Cefazolin: *Proteus mirabilis*

Analysis of trending indicated that MIC values for *C. koseri*, *E. coli* and *K. pneumoniae* tended to be one doubling dilution lower than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Cefazolin: *Citrobacter koseri*, *Escherichia coli*, *Klebsiella pneumoniae*

**Cefepime.** A total of 956 Enterobacterales isolates (77 *C. freundii* complex, 74 *C. koseri*, 77 *E. cloacae* complex, 170 *E. coli*, 52 *K. aerogenes*, 29 *K. oxytoca*, 111 *K. pneumoniae*, 75 *M. morgani*, 88 *P. mirabilis*, 100 *P. vulgaris*, 103 *S. marcescens*) were evaluated with cefepime. The combined results from clinical and challenge isolate testing demonstrated an EA of 94.8% and CA of 96.0%. There were 26 minor, 10 major and 2 very major errors. When evaluating results by individual species, *M. morgani* had 1 very major error (1/1 = 100%) and *P. vulgaris* had 1 very major error (1/2 = 50%), which were considered random errors due to the insufficient number of resistant isolates tested. *K. aerogenes* had 2 major errors (2/51 = 3.9%) at Selux MIC values of  $\geq 64$   $\mu\text{g/mL}$  which represents 5.8% (3/52) of results with this drug/organism combination. The following limitation is included in the device labeling to address the major errors:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Cefepime: *Klebsiella aerogenes* when the Selux MIC is to  $\geq 64$   $\mu\text{g/mL}$  due to two major errors

A total of 215 *P. aeruginosa* isolates were evaluated with cefepime. The combined results from clinical and challenge isolate testing demonstrated an EA of 87.4% and CA of 86.8%. There were 0 minor, 13 major (13/178 = 7.3%) and 7 very major (7/37 = 18.9%) errors. Due to the lack of an intermediate interpretive criterion, further analysis of the errors was performed and adjustments were made by considering the MIC values of the errors compared to the reference MIC value. Nine of the 13 major errors and 6 of the 7 very major errors had an MIC value that was in essential agreement with the reference MIC value. Therefore, the adjusted major error rate is 2.2% (4/178), which is acceptable. However, the adjusted very major error rate is 2.7% (1/37) which is not acceptable. Due to the unacceptable performance for *P. aeruginosa*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *P. aeruginosa* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Cefepime: *Pseudomonas aeruginosa*

A limitation statement is included in the device labeling to address the lack of testing with resistant *K. aerogenes*, *M. morgani*, *P. vulgaris* and *S. marcescens* isolates.

Analysis of trending indicated that MIC values for *C. freundii*, *C. koseri*, *E. cloacae*, *E. coli*, *K. aerogenes*, *K. oxytoca*, *M. morgani*, *P. mirabilis*, *P. vulgaris*, *S. marcescens* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Cefepime: *Citrobacter freundii*, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Morganella morgani*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*

**Cefoxitin.** A total of 637 Enterobacterales isolates (indicated species: 168 *E. coli*, 29 *K. oxytoca*, 133 *K. pneumoniae*, 64 *M. morgani*, 89 *P. mirabilis*, 100 *P. vulgaris*; non-indicated species: 54 *C. koseri*) were evaluated with cefoxitin. The combined results from clinical and challenge isolate testing demonstrated an EA of 92.3% and CA of 80.2%. There were 114 minor, 9 major and 3 very major errors. When evaluating results by individual species, *K. oxytoca* and *M. morgani* had CA of <90% which was considered acceptable since most of the categorical errors were minor and the EA of evaluable results was good (>90%). *P. mirabilis* had a CA of 88.8% due to 6 minor, 2 major (2/75 = 2.7%) and 2 very major (2/7 = 28.6%) errors, which is not acceptable. *P. vulgaris* had an EA of 86% and a CA of 53% due to 46 minor errors and 1 major error (1/53 = 1.9%), which is not acceptable. Due to the unacceptable performance for *P. mirabilis* and *P. vulgaris*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *P. mirabilis* and *P. vulgaris* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Cefoxitin: *Proteus mirabilis*, *Proteus vulgaris*

A limitation statement is included in the device labeling to address the lack of testing with resistant *K. oxytoca* and *P. vulgaris* isolates.

There was no evidence of trending observed with any indicated species with cefoxitin.

**Ceftazidime.** A total of 127 *A. baumannii* isolates were evaluated with ceftazidime. The combined results from clinical and challenge isolate testing demonstrated an EA of 86.6% and CA of 94.5%. There were 7 minor, 0 major and 0 very major error. In addition, the EA of evaluable results was 81.5% (75/92). Due to the unacceptable performance for *A. baumannii*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *A. baumannii* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Ceftazidime: *Acinetobacter baumannii*

A total of 787 Enterobacterales isolates (64 *C. freundii*, 55 *C. koseri*, 84 *E. cloacae*, 186 *E. coli*, 36 *K. aerogenes*, 29 *K. oxytoca*, 145 *K. pneumoniae*, 66 *P. mirabilis*, 55 *P. vulgaris*, 67 *S. marcescens*) were evaluated with ceftazidime. The combined results from clinical and challenge isolate testing demonstrated an EA of 95.6% and CA of 96.4%. There were 26 minor, 1 major and 1 very major errors. Overall, performance is acceptable.

A total of 224 *P. aeruginosa* isolates were evaluated with ceftazidime. The combined results from clinical and challenge isolate testing demonstrated an EA of 88.0% and CA of 93.3%. There were 9 major (9/178 = 5.1%) and 6 very major (6/46 = 13%) errors. Due to the lack of an intermediate interpretive criterion, further analysis of the errors was performed and adjustments were made by considering the MIC values of the errors compared to the reference MIC value. One of the 9 major errors and 1 of the 6 very major errors had an MIC value that was in essential agreement with the reference MIC value. Therefore, the adjusted major error rate is 4.5% (8/177) and the adjusted very major error rate is 10.9% (5/46), which is not acceptable. Due to the unacceptable performance for *P. aeruginosa*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *P. aeruginosa* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Ceftazidime: *Pseudomonas aeruginosa*

Analysis of trending indicated that MIC values for *C. freundii* complex, *C. koseri*, *E. cloacae*, *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*, *P. mirabilis*, *P. vulgaris*, and *S. marcescens* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Ceftazidime: *Citrobacter freundii* complex, *Citrobacter koseri*, *Enterobacter cloacae*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*

**Ceftazidime-avibactam.** A total of 815 Enterobacterales isolates (61 *C. freundii*, 54 *C. koseri*, 75 *E. cloacae*, 170 *E. coli*, 52 *K. aerogenes*, 28 *K. oxytoca*, 126 *K. pneumoniae*, 62 *M. morgani*, 66 *P. mirabilis*, 55 *P. vulgaris*, 66 *S. marcescens*) were evaluated with ceftazidime-avibactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 97.3% and CA of 99.4%. There were 5 major and 0 very major errors. Overall, performance is acceptable.

A total of 163 *P. aeruginosa* isolates were evaluated with ceftazidime-avibactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 89.6% and CA of 96.9%. There were 4 major errors (4/153 = 2.6%) and 1 very major error (1/10 = 10%). Due to the lack of an intermediate interpretive criterion, further analysis of the errors was performed and adjustments were made by considering the MIC values of the errors compared to the reference MIC value. Two of the 4 major errors and 1 of the 1 very major errors had an MIC value that was in essential agreement with the reference MIC value.

Therefore, the adjusted major error rate is 1.3% (2/153) and the adjusted very major error rate is 0% (0/10), which is acceptable.

A limitation statement is included in the device labeling to address the lack of testing with resistant *C. koseri*, *K. aerogenes*, *K. oxytoca*, *M. morgani*, *P. mirabilis*, *P. vulgaris* and *S. marcescens* isolates.

Analysis of trending indicated that MIC values for *C. freundii*, *C. koseri*, *E. coli*, *K. aerogenes*, *K. oxytoca*, *M. morgani*, *P. mirabilis*, and *P. vulgaris* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Ceftazidime-avibactam: *Citrobacter freundii*, *Citrobacter koseri*, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Morganella morgani*, *Proteus mirabilis*, *Proteus vulgaris*

**Ceftriaxone.** A total of 722 Enterobacterales isolates (64 *C. freundii*, 54 *C. koseri*, 84 *E. cloacae*, 175 *E. coli*, 35 *K. aerogenes*, 29 *K. oxytoca*, 109 *K. pneumoniae*, 67 *P. mirabilis*, 105 *S. marcescens*) were evaluated with ceftriaxone. The combined results from clinical and challenge isolate testing demonstrated an EA of 97.5% and CA of 98.5%. There were 3 minor, 4 major and 4 very major errors. When evaluating results by individual species, *E. cloacae* had 3 very major errors (3/34 = 8.8%) and *S. marcescens* had 1 very major error (1/12 = 8.3%), which are not acceptable. Due to the unacceptable performance for *E. cloacae* and *S. marcescens*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *E. cloacae* and *S. marcescens* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Ceftriaxone: *Enterobacter cloacae*, *Serratia marcescens*

Analysis of trending indicated that MIC values for *C. freundii*, *C. koseri*, *E. coli*, *K. aerogenes*, *K. oxytoca*, and *P. mirabilis* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Ceftriaxone: *Citrobacter freundii*, *Citrobacter koseri*, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Proteus mirabilis*

**Ciprofloxacin.** A total of 818 Enterobacterales isolates (62 *C. freundii*, 55 *C. koseri*, 76 *E. cloacae*, 170 *E. coli*, 32 *K. aerogenes*, 29 *K. oxytoca*, 105 *K. pneumoniae*, 63 *M. morgani*, 67 *P. mirabilis*, 55 *P. vulgaris*, 104 *S. marcescens*) were evaluated with ciprofloxacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 96.7% and CA of 96.2%. There were 26 minor, 3 major and 2 very major errors. When evaluating

results by individual species, *S. marcescens* had 2 very major errors (2/12 = 16.7%), which is not acceptable. Due to the unacceptable performance for *S. marcescens*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *S. marcescens* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Ciprofloxacin: *Serratia marcescens*

A total of 169 *P. aeruginosa* isolates were evaluated with ciprofloxacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 95.9% and CA of 94.7%. There were 5 minor, 4 major and 0 very major errors. Overall, performance is acceptable.

A limitation statement is included in the device labeling to address the lack of testing with resistant *C. koseri* and *P. vulgaris* isolates.

Analysis of trending indicated that MIC values for *P. aeruginosa* tended to be one doubling dilution lower than the reference MIC value while MIC values for *C. freundii*, *C. koseri*, *E. cloacae*, *E. coli*, *K. aerogenes*, *K. oxytoca*, *M. morgani*, and *P. vulgaris* tended to be one doubling dilution higher than the reference MIC value. The following statements are included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Ciprofloxacin: *Pseudomonas aeruginosa*

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Ciprofloxacin: *Citrobacter freundii*, *Citrobacter koseri*, *Enterobacter cloacae*, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Morganella morgani*, *Proteus vulgaris*

**Eravacycline.** A total of 487 Enterobacterales isolates (62 *C. freundii*, 78 *E. cloacae*, 167 *E. coli*, 47 *K. oxytoca*, 133 *K. pneumoniae*) were evaluated with eravacycline. The combined results from clinical and challenge isolate testing demonstrated an EA of 98.2% and CA of 98.6%. There were 4 major and 3 very major errors. When evaluating results by individual species, *E. cloacae* had 3 major errors (3/59 = 5.1%), *K. oxytoca* had 1 very major error (1/1 = 100%), and *K. pneumoniae* had 2 very major errors (2/30 = 6.7%). Due to the lack of an intermediate interpretive criterion (eravacycline only has a “susceptible” category), further analysis of the errors was performed and adjustments were made by considering the MIC values of the errors compared to the reference MIC value. Two of the 3 *E. cloacae* major errors and 1 of the 1 *K. oxytoca* very major errors had an MIC value that was in essential agreement with the reference MIC value. Therefore, the adjusted major error rate for *E. cloacae* is 1.7% (1/59) and the adjusted very major error rate for *K. oxytoca* is 0% (0/1), which is acceptable. None of the 2 *K. pneumoniae* very major errors were in essential agreement resulting in an unacceptable very major error rate. Due to the unacceptable

performance for *K. pneumoniae*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *K. pneumoniae* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Eravacycline: *Klebsiella pneumoniae*

A limitation statement is included in the device labeling to address the lack of testing with resistant *E. coli* and *K. oxytoca* isolates.

Analysis of trending indicated that MIC values for *C. freundii* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Eravacycline: *Citrobacter freundii*

**Ertapenem.** A total of 977 Enterobacterales isolates (56 *C. freundii*, 74 *C. koseri*, 95 *E. cloacae*, 210 *E. coli*, 55 *K. aerogenes*, 28 *K. oxytoca*, 141 *K. pneumoniae*, 62 *M. morgani*, 88 *P. mirabilis*, 100 *P. vulgaris*, 68 *S. marcescens*) were evaluated with ertapenem. The combined results from clinical and challenge isolate testing demonstrated an EA of 94.7% and CA of 97.4%. There were 11 minor, 14 major and 0 very major errors. When evaluating results by individual species, *C. freundii* had 3 major errors (3/53 = 5.7%) at Selux MIC values of 2, 4 and 8 µg/mL which represent 5.4% (3/56) of results from this drug/organism combination. *E. cloacae* had an EA of 80%, a CA of 88.4%, and 6 major errors (6/70 = 8.6%), which is not acceptable. *K. aerogenes* had an EA of 80% and 2 major errors (2/50 = 4%), which is not acceptable. Due to the unacceptable performance for *E. cloacae* and *K. aerogenes*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation is included in the device labeling to address the *C. freundii* major errors and restrict reporting of *E. cloacae* and *K. aerogenes* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Ertapenem: *Citrobacter freundii* when the Selux MIC is 2, 4 or 8 µg/mL due to 3 major errors, *Enterobacter cloacae*, *Klebsiella aerogenes*

A limitation statement is included in the device labeling to address the lack of testing with resistant *C. koseri*, *K. oxytoca*, *M. morgani*, *P. mirabilis*, and *P. vulgaris* isolates.

Analysis of trending indicated that MIC values for *M. morgani* tended to be one doubling dilution lower than the reference MIC value while MIC values for *C. freundii*, *C. koseri*, *E. coli*, *P. mirabilis*, *P. vulgaris* and *S. marcescens* tended to be one doubling dilution higher than the reference MIC value. The following statements are included as footnotes to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Ertapenem: *Morganella morganii*

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Ertapenem: *Citrobacter freundii*, *Citrobacter koseri*, *Escherichia coli*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*

**Gentamicin.** A total of 817 Enterobacterales isolates (indicated species: 56 *C. freundii*, 74 *C. koseri*, 95 *E. cloacae*, 210 *E. coli*, 55 *K. aerogenes*, 28 *K. oxytoca*, 141 *K. pneumoniae*, 88 *P. mirabilis*, 100 *P. vulgaris*, 68 *S. marcescens*; non-indicated species: 62 *M. morganii*) were evaluated with gentamicin. The combined results from clinical and challenge isolate testing demonstrated an EA of 96.0% and CA of 98.8%. There were 7 minor, 1 major and 2 very major errors. When evaluating results by individual species, *M. morganii* had 1 very major error (1/8 = 12.5%) and *P. mirabilis* had 1 very major error (1/6 = 16.7%), which were considered random errors due to the limited number of resistant isolates tested.

A total of 218 *P. aeruginosa* were evaluated with gentamicin. The combined results from clinical and challenge isolate testing demonstrated an EA of 88.5% and CA of 97.3%. There were 6 minor, 0 major and 0 very major errors. Restriction of reporting MIC values of  $\leq 0.06$   $\mu\text{g/mL}$  (lowest reportable MIC value that represents 4.1% (9/218) of results from this drug/organism combination) improved the EA to 90.4%, which is acceptable. The following limitation statement is included in the device labeling to restrict reporting of *P. aeruginosa* due to the low EA:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Gentamicin: *Pseudomonas aeruginosa* when the Selux MIC is  $\leq 0.06$   $\mu\text{g/mL}$

A limitation statement is included in the device labeling to address the lack of testing with resistant *K. aerogenes*, *M. morganii*, *P. vulgaris*, and *S. marcescens* isolates.

Analysis of trending indicated that MIC values for *C. freundii*, *E. cloacae* complex, *E. coli*, *K. oxytoca*, and *K. pneumoniae* tended to be one doubling dilution lower than the reference MIC value while MIC values for *S. marcescens* tended to be one doubling dilution higher than the reference MIC value. The following statements are included as footnotes to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Gentamicin: *Citrobacter freundii*, *Enterobacter cloacae*, *Escherichia coli*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Gentamicin: *Serratia marcescens*



**Imipenem.** A total of 126 *A. baumannii* isolates were evaluated with imipenem. The combined results from clinical and challenge isolate testing demonstrated an EA of 88.9% and CA of 98.4%. There were 0 minor, 0 major and 2 very major errors (2/62 = 3.2%), which is not acceptable. Due to the unacceptable performance for *A. baumannii*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *A. baumannii* due to the unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Imipenem: *Acinetobacter baumannii*

**Imipenem-relebactam.** A total of 471 Enterobacterales isolates (55 *C. freundii*, 55 *C. koseri*, 172 *E. coli*, 31 *K. aerogenes*, 27 *K. oxytoca*, 131 *K. pneumoniae*) were evaluated with imipenem-relebactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 93.4% and CA of 97.5%. There were 9 minor, 3 major and 0 very major errors. When evaluating results by individual species, *K. oxytoca* had 1 major error (1/27 = 3.7%) at a Selux MIC value of 4 µg/mL which represents 3.7% (1/27) of results from this drug/organism combination. *K. aerogenes* and *K. pneumoniae* had an EA <90%, which is not acceptable. Due to the unacceptable performance for *K. aerogenes* and *K. pneumoniae*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation is included in the device labeling to address the *K. oxytoca* major error and restrict reporting of *K. aerogenes* and *K. pneumoniae* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Imipenem-relebactam: *Klebsiella oxytoca* when the Selux MIC is 4 µg/mL due to one major error, *Klebsiella pneumoniae*, *Klebsiella aerogenes*

A total of 165 *P. aeruginosa* isolates were evaluated with imipenem-relebactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 96.4% and CA of 98.2%. There were 3 minor, 0 major and 0 very major errors. Overall, performance was acceptable.

A limitation statement is included in the device labeling to address the lack of testing with resistant *C. koseri*, *K. aerogenes*, and *K. oxytoca* isolates.

Analysis of trending indicated that MIC values for *C. freundii* tended to be one doubling dilution lower than the reference MIC value while MIC values for *K. oxytoca* tended to be one doubling dilution higher than the reference MIC value. The following statements are included as footnotes to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Imipenem-relebactam: *Citrobacter freundii*

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Imipenem-relebactam: *Klebsiella oxytoca*

**Levofloxacin.** A total of 784 Enterobacterales isolates (63 *C. freundii*, 54 *C. koseri*, 75 *E. cloacae*, 171 *E. coli*, 32 *K. aerogenes*, 27 *K. oxytoca*, 110 *K. pneumoniae*, 62 *M. morgani*, 67 *P. mirabilis*, 55 *P. vulgaris*, 68 *S. marcescens*) were evaluated with levofloxacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 96.1% and CA of 94.9%. There were 33 minor, 6 major and 1 very major errors. When evaluating results by individual species, *K. aerogenes* had 1 major error (1/30 = 3.3%) at a Selux MIC value of 2 µg/mL, which represents 3.1% (1/32) of results from this drug/organism combination. *C. freundii* had 2 major errors (2/55 = 3.6%) at a Selux MIC value of 2 µg/mL which represent 4.8% (3/63) of results from this drug/organism combination, as well as 1 very major error (1/7 = 14.3%), which was considered a random error due to the limited number of resistant isolates tested. The following limitation is included in the device labeling to address the major errors:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Levofloxacin: *Citrobacter freundii* when the Selux MIC is 2 µg/mL due to two major errors, *Klebsiella aerogenes* when the Selux MIC is 2 µg/mL due to one major error

A total of 171 *P. aeruginosa* isolates were evaluated with levofloxacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 94.2% and CA of 90.1%. There were 14 minor, 1 major and 2 very major errors (2/42 = 4.8%), which is not acceptable. Due to the unacceptable performance for *P. aeruginosa*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting of *P. aeruginosa* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Levofloxacin: *Pseudomonas aeruginosa*

A limitation statement is included in the device labeling to address the lack of testing with resistant *C. koseri* isolates.

Analysis of trending indicated that MIC values for *C. freundii*, *C. koseri*, *E. cloacae*, *E. coli*, *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*, *M. morgani*, *P. mirabilis*, and *P. vulgaris* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Levofloxacin: *C. freundii*, *C. koseri*, *E. cloacae*, *E. coli*, *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*, *M. morgani*, *P. mirabilis*, *P. vulgaris*

**Meropenem.** A total of 126 *A. baumannii* isolates were evaluated with meropenem. The combined results from clinical and challenge isolate testing demonstrated an EA of 86.5%

and CA of 98.4%. There were 0 minor, 2 major (2/63 = 3.2%) and 0 very major errors. Due to the unacceptable performance for *A. baumannii*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *A. baumannii* due to the unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Meropenem: *Acinetobacter baumannii*

A total of 865 Enterobacterales isolates (indicated species: 78 *C. freundii*, 74 *C. koseri*, 74 *E. cloacae*, 166 *E. coli*, 27 *K. oxytoca*, 126 *K. pneumoniae*, 62 *M. morgani*, 66 *P. mirabilis*, 55 *P. vulgaris*, 105 *S. marcescens*; non-indicated species: 32 *K. aerogenes*) were evaluated with meropenem. The combined results from clinical and challenge isolate testing demonstrated an EA of 96.3% and CA of 98.0%. There were 4 minor, 12 major and 1 very major errors. When evaluating results by individual species, *C. freundii* has 3 major errors (3/74 = 4.1%) at Selux MIC values of 16 and 32 µg/mL, which represent 6.4% (5/78) of results from this drug/organism combination. *S. marcescens* had 4 major errors (4/101 = 4.0%) at Selux MIC value of 16 and 32 µg/mL, which represent 6.7% (7/105) of results from this drug/organism combination, as well as 1 very major error (1/4 = 25%) which was considered a random error due to the limited number of resistant isolates tested. The following limitations are included in the device labeling to address the major errors:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Meropenem: *Citrobacter freundii* complex when the Selux MIC is 16 or 32 µg/mL due to 3 major errors, *Serratia marcescens* when the Selux MIC is 16 or 32 µg/mL due to 4 major errors

A total of 175 *P. aeruginosa* isolates were evaluated with levofloxacin. The combined results from clinical and challenge isolate testing demonstrated an EA of 93.1% and CA of 96.0%. There were 7 minor, 0 major and 0 very major errors. Overall, performance was acceptable.

A limitation is included in the device labeling to address the lack of testing with resistant *C. koseri*, *K. oxytoca*, *M. morgani*, *P. vulgaris* isolates.

Analysis of trending indicated that MIC values for *C. freundii*, *C. koseri*, *E. cloacae*, *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*, *M. morgani*, *P. mirabilis*, *P. vulgaris*, and *S. marcescens* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Meropenem: *Citrobacter freundii*, *Citrobacter koseri*, *Enterobacter cloacae*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morgani*, *Proteus mirabilis*, *Proteus vulgaris*, *Serratia marcescens*

**Meropenem-vaborbactam.** A total of 865 Enterobacterales isolates (indicated species: 64 *C. freundii*, 54 *C. koseri*, 70 *E. cloacae*, 167 *E. coli*, 50 *K. aerogenes*, 47 *K. oxytoca*, 113 *K. pneumoniae*, 64 *M. morgani*, 66 *P. mirabilis*, 65 *S. marcescens*; non-indicated species: 100 *P. vulgaris*) were evaluated with meropenem-vaborbactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 91.7% and CA of 96.9%. There were 10 minor, 16 major and 1 very major errors. When evaluating results by individual species, *P. mirabilis* had an EA of 62%, CA of 97.0%, 2 minor, 0 major and 0 very major errors. *P. vulgaris* had an EA of 62%, CA of 81%, 15 major errors (15/100 = 15%). Due to the unacceptable performance for *P. mirabilis* and *P. vulgaris*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation is included in the device labeling to restrict reporting due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Meropenem-vaborbactam: *Proteus mirabilis*; *Proteus vulgaris*

A limitation is included in the device labeling to address the lack of testing with resistant *C. koseri*, *K. aerogenes*, *K. oxytoca*, *M. morgani*, *P. mirabilis*, *P. vulgaris*, and *S. marcescens* isolates.

Analysis of trending indicated that MIC values for *C. freundii*, *C. koseri*, *E. cloacae* complex, *E. coli*, *K. aerogenes*, *K. oxytoca*, *K. pneumoniae*, *M. morgani*, and *S. marcescens* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Meropenem-vaborbactam: *Citrobacter freundii*, *Citrobacter koseri*, *Enterobacter cloacae* complex, *Escherichia coli*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Morganella morgani*, *Serratia marcescens*

**Minocycline.** A total of 199 *A. baumannii* isolates were evaluated with minocycline. The combined results from clinical and challenge isolate testing demonstrated an EA of 88.9% and CA of 91.0%. There were 13 minor, 5 major (5/150 = 3.3%) and 0 very major errors. Due to the unacceptable performance for *A. baumannii*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *A. baumannii* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Minocycline: *Acinetobacter baumannii*

A total of 404 Enterobacterales isolates (198 *E. coli*, 53 *K. aerogenes*, 27 *K. oxytoca*, 126 *K. pneumoniae*) were evaluated with minocycline. The combined results from clinical and challenge isolate testing demonstrated an EA of 89.6% and CA of 94.3%. There were 16 minor, 3 major and 4 very major errors. When evaluating results by individual species, *E. coli* had an EA of 87.4% as well as 8 minor, 2 major (2/182 = 1.1%) and 1 very major (1/7 =

14.3%) errors, which was considered a random error due to the limited number of resistant isolates tested. *K. aerogenes* had 1 very major error ( $1/3 = 33.3\%$ ), which was considered a random error due to the limited number of resistant isolates tested. *K. pneumoniae* had an EA of 88.9% and 2 very major errors ( $2/33 = 6.1\%$ ) at a Selux MIC of 4. Restricting *K. pneumoniae* reporting at a Selux MIC of 4 µg/mL, which represents 6.3% (8/126) of results from this drug/organism combination addresses both very major errors and improves EA to 89.8%. Due to the unacceptable performance for *E. coli*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *E. coli* due to unacceptable performance and address the *K. pneumoniae* very major errors:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Minocycline: *Escherichia coli*, *Klebsiella pneumoniae* when the Selux MIC is 4 µg/mL due to two very major errors

A limitation is included in the device labeling to address the lack of testing with resistant *K. oxytoca* isolates.

There was no evidence of trending observed with Enterobacterales with minocycline.

**Piperacillin-tazobactam.** A total of 113 *A. baumannii* isolates were evaluated with piperacillin-tazobactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 81.4% and CA of 91.2%. There were 9 minor errors, 0 major errors and 1 very major error ( $1/58 = 1.7\%$ ). Due to the unacceptable performance for *A. baumannii*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *A. baumannii* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Piperacillin-tazobactam: *Acinetobacter baumannii*

A total of 926 Enterobacterales isolates (indicated species: 75 *C. koseri*, 176 *E. coli*, 142 *K. pneumoniae*, 75 *M. morgani*, 88 *P. mirabilis*, 55 *P. vulgaris*, 60 *S. marcescens*; non-indicated species: 78 *C. freundii*, 96 *E. cloacae*, 53 *K. aerogenes*, 28 *K. oxytoca*) were evaluated with piperacillin-tazobactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 90.2% and CA of 94.8%. There were 35 minor, 9 major and 4 very major errors. When evaluating results by individual species, *C. freundii* had an EA of 79.5%, CA of 84.6% as well as 11 minor, 0 major, and 1 very major ( $1/11 = 9.1\%$ ) errors. *E. cloacae* had an EA and CA of 86.5% as well as 9 minor and 4 major ( $4/61 = 6.6\%$ ) errors. *K. aerogenes* had an EA of 86.8% as well as 2 minor and 2 major ( $2/42 = 4.8\%$ ) errors. *S. marcescens* had an EA of 66.7% as well as 1 minor and 1 very major (1/2) error, which was considered a random error due to the limited number of resistant isolates tested. Due to the unacceptable performance for reporting of *C. freundii*, *E. cloacae*, *K. aerogenes* and *S. marcescens*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to

restrict reporting of *C. freundii*, *E. cloacae*, *K. aerogenes* and *S. marcescens* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Piperacillin-tazobactam: *Citrobacter freundii*, *Enterobacter cloacae*, *Klebsiella aerogenes*, *Serratia marcescens*

A total of 217 *P. aeruginosa* isolates were evaluated with piperacillin-tazobactam. The combined results from clinical and challenge isolate testing demonstrated an EA of 80.7% and CA of 83.9%. There were 27 minor errors, 2 major errors (2/167 = 1.2%), and 6 very major errors (6/28 = 21.4%). Due to the unacceptable performance for *P. aeruginosa*, this drug/organism combination is not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *P. aeruginosa* due to unacceptable performance:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combinations:

- Piperacillin-tazobactam: *Pseudomonas aeruginosa*

A limitation is included in the device labeling to address the lack of testing with resistant *C. koseri*, *P. mirabilis*, and *P. vulgaris* isolates.

Analysis of trending indicated that MIC values for *K. oxytoca* and *K. pneumoniae* tended to be one doubling dilution lower than the reference MIC value while MIC values for *M. morgani*, *P. mirabilis* and *P. vulgaris* tended to be one doubling dilution higher than the reference MIC value. The following statements are included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution lower than the reference MIC value:

- Piperacillin-tazobactam: *Klebsiella oxytoca*, *Klebsiella pneumoniae*

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:

- Piperacillin-tazobactam: *Morganella morgani*, *Proteus mirabilis*, *Proteus vulgaris*

**Tobramycin.** A total of 166 *P. aeruginosa* isolates were evaluated with tobramycin. The combined results from clinical and challenge isolate testing demonstrated an EA of 93.4% and CA of 97.6%. There were 3 minor errors, 1 major error (1/150 = 0.7%), and 0 very major errors (6/28 = 21.4%). Overall, performance was acceptable.

There was no evidence of trending observed with *P. aeruginosa* with tobramycin.

**Trimethoprim-Sulfamethoxazole.** A total of 526 Enterobacterales isolates (indicated species: 76 *E. cloacae*, 203 *E. coli*, 32 *K. aerogenes*, 29 *K. oxytoca*, 109 *K. pneumoniae*; non-indicated species: 77 *C. freundii*) were evaluated with trimethoprim-sulfamethoxazole. The combined

results from clinical and challenge isolate testing demonstrated an EA of 97.0% and CA of 98.5%. There were 0 minor, 2 major and 6 very major errors. When evaluating results by individual species, *C. freundii* had 1 very major error (1/17 = 5.9%), *E. coli* had 4 very major errors (4/67 = 6.0%), *K. pneumoniae* had 2 major (2/61 = 3.3%) and 1 very major (1/28 = 3.6%) errors. Due to the lack of an intermediate interpretive criterion, further analysis of the errors was performed and adjustments were made by considering the MIC values of the errors compared to the reference MIC value. One of the 2 major errors and 1 of the 1 very major errors from *K. pneumoniae* had an MIC value that was in essential agreement with the reference MIC value. Therefore, *K. pneumoniae* has an adjusted major error rate of 1.6% (1/61) and an adjusted very major error rate of 0% (0/28), which is acceptable. Due to the unacceptable performance for *C. freundii* and *E. coli*, these drug/organism combinations are not indicated for use with the Selux AST System. The following limitation statement is included in the device labeling to restrict reporting of *C. freundii* and *E. coli*:

Perform an alternative method of testing prior to reporting results for the following antibiotic/organism combination(s):

- Trimethoprim-Sulfamethoxazole: *Citrobacter freundii*, *Escherichia coli*

Analysis of trending indicated that MIC values for *E. cloacae*, *K. aerogenes*, *K. oxytoca*, and *K. pneumoniae* tended to be one doubling dilution higher than the reference MIC value. The following statement is included as a footnote to the AST performance table:

Selux MIC values for the following antimicrobial/organism combinations tended to be one doubling dilution higher than the reference MIC value:  
Trimethoprim-Sulfamethoxazole: *Enterobacter cloacae*, *Klebsiella aerogenes*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*

**Table 6.** Selux AST System – Gram-Negative Panel Performance

	Tot	No. EA	EA %	Eval EA Tot	No. Eval EA	Eval EA %	No. CA	CA %	No. R or NS	No. S	min	maj	vmj
aAmikacin - <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 16 (S), 32 (I), 64 (R)]													
Clinical	86	47	54.65	77	38	49.35	78	90.7	14	69	3	4	1
Challenge	25	14	56	16	5	31.25	18	72	11	11	3	4	0
Combined	111	61	54.95	93	43	46.24	96	86.49	25	80	6	8	1
Amikacin - <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 16 (S), 32 (I), 64 (R)]													
Clinical	159	145	91.19	158	144	91.14	156	98.11	3	154	3	0	0
Challenge	6	5	83.33	5	4	80	6	100	4	2	0	0	0
Combined	165	150	90.91	163	148	90.8	162	98.18	7	156	3	0	0
Amoxicillin-Clavulanate – Enterobacterales [Breakpoints (µg/mL): 8 (S), 16 (I), 32 (R)]													
Clinical	502	463	92.23	454	415	91.41	460	91.63	40	412	37	4	1
Challenge	91	90	98.9	49	48	97.96	90	98.9	71	19	0	1	0
Combined	593	553	93.25	503	463	92.05	550	92.75	111	431	37	5	1
Ampicillin: Enterobacterales [Breakpoints (µg/mL): 8 (S), 16 (I), 32 (R)]													
Clinical	247	234	94.74	138	125	90.58	244	98.79	111	135	1	1	1
Challenge	7	7	100	0	0	NA	7	100	7	0	0	0	0
Combined	254	241	94.88	138	125	90.58	251	98.82	118	135	1	1	1
Ampicillin-Sulbactam: <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 8 (S), 16 (I), 32 (R)]													
Clinical	121	111	91.74	111	101	90.99	112	92.56	50	65	6	2	1
Challenge	2	2	100	2	2	100	2	100	2	0	0	0	0
Combined	123	113	91.87	113	103	91.15	114	92.68	52	65	6	2	1
Ampicillin-Sulbactam: Enterobacterales [Breakpoints (µg/mL): 8 (S), 16 (I), 32 (R)]													

	Tot	No. EA	EA %	Eval EA Tot	No. Eval EA	Eval EA %	No. CA	CA %	No. R or NS	No. S	min	maj	vmj
Clinical	518	502	96.91	475	459	96.63	433	83.59	163	279	81	4	0
Challenge	56	54	96.43	26	24	92.31	46	82.14	39	11	9	1	0
Combined	574	556	96.86	501	483	96.41	479	83.45	202	290	90	5	0
<sup>a</sup> Aztreonam: Enterobacterales [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	159	156	98.11	144	141	97.92	154	96.86	22	129	5	0	0
Challenge	24	23	95.83	10	9	90	24	100	24	0	0	0	0
Combined	183	179	97.81	154	150	97.4	178	97.27	46	129	5	0	0
Aztreonam: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 8 (S), 16 (I), 32 (R)]													
Clinical	208	189	90.87	202	183	90.59	181	87.02	37	158	20	3	4
Challenge	11	11	100	4	4	100	11	100	11	0	0	0	0
Combined	219	200	91.32	206	187	90.78	192	87.67	48	158	20	3	4
Cefazolin: Enterobacterales [Breakpoints (µg/mL): 1 (S), 2 (I), 4 (R)]													
Clinical	368	344	93.48	290	266	91.72	264	71.74	174	99	102	1	1
Challenge	44	44	100	0	0	NA	44	100	44	0	0	0	0
Combined	412	388	94.17	290	266	91.72	308	74.76	218	99	102	1	1
Cefepime: Enterobacterales [Breakpoints (µg/mL): 2 (S), 4-8 (I), 16 (R)]													
Clinical	905	857	94.7	127	79	62.2	867	95.8	78	810	26	10	2
Challenge	51	49	96.08	17	15	88.24	51	100	49	2	0	0	0
Combined	956	906	94.77	144	94	65.28	918	96.03	127	812	26	10	2
<sup>a</sup> Cefepime: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 8 (S), - (I), 16 (R)]													
Clinical	208	182	87.5	203	177	87.19	188	90.38	30	178	0	13	7
Challenge	7	6	85.71	2	1	50	7	100	7	0	0	0	0
Combined	215	188	87.44	205	178	86.83	195	90.7	37	178	0	13	7
Cefoxitin: Enterobacterales [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	600	551	91.83	576	527	91.49	474	79	71	394	114	9	3
Challenge	37	37	100	15	15	100	37	100	35	2	0	0	0
Combined	637	588	92.31	591	542	91.71	511	80.22	106	396	114	9	3
<sup>a</sup> Ceftazidime: <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 8 (S), 16 (I), 32 (R)]													
Clinical	121	105	86.78	90	74	82.22	114	94.21	52	61	7	0	0
Challenge	6	5	83.33	2	1	50	6	100	5	1	0	0	0
Combined	127	110	86.61	92	75	81.52	120	94.49	57	62	7	0	0
Ceftazidime: Enterobacterales [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	670	635	94.78	140	105	75	642	95.82	95	563	26	1	1
Challenge	117	117	100	11	11	100	117	100	113	4	0	0	0
Combined	787	752	95.55	151	116	76.82	759	96.44	208	567	26	1	1
<sup>a</sup> Ceftazidime: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 8 (S), - (I), 16 (R)]													
Clinical	211	185	87.68	205	179	87.32	196	92.89	34	177	0	9	6
Challenge	13	12	92.31	8	7	87.5	13	100	12	1	0	0	0
Combined	224	197	87.95	213	186	87.32	209	93.3	46	178	0	9	6
Ceftazidime-avibactam: Enterobacterales [Breakpoints (µg/mL): 8 (S), - (I), 16 (R)]													
Clinical	772	751	97.28	211	190	90.05	767	99.35	2	770	0	5	0
Challenge	43	42	97.67	11	10	90.91	43	100	30	13	0	0	0
Combined	815	793	97.3	222	200	90.09	810	99.39	32	783	0	5	0
Ceftazidime-avibactam: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 8 (S), - (I), 16 (R)]													
Clinical	157	140	89.17	153	136	88.89	153	97.45	5	152	0	4	0
Challenge	6	6	100	2	2	100	5	83.33	5	1	0	0	1
Combined	163	146	89.57	155	138	89.03	158	96.93	10	153	0	4	1
Ceftriaxone: Enterobacterales [Breakpoints (µg/mL): 1 (S), 2 (I), 4 (R)]													
Clinical	666	648	97.3	68	50	73.53	655	98.35	134	530	3	4	4
Challenge	56	56	100	1	1	100	56	100	54	1	0	0	0
Combined	722	704	97.51	69	51	73.91	711	98.48	188	531	3	4	4
Ciprofloxacin: Enterobacterales [Breakpoints (µg/mL): 0.25 (S), 0.5 (I), 1 (R)]													
Clinical	770	743	96.49	150	123	82	739	95.97	130	622	26	3	2



	Tot	No. EA	EA %	Eval EA Tot	No. Eval EA	Eval EA %	No. CA	CA %	No. R or NS	No. S	min	maj	vmj
Challenge	48	48	100	6	6	100	48	100	47	1	0	0	0
Combined	818	791	96.7	156	129	82.69	787	96.21	177	623	26	3	2
Ciprofloxacin: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 0.5 (S), 1 (I), 2 (R)]													
Clinical	159	152	95.6	147	140	95.24	150	94.34	23	132	5	4	0
Challenge	10	10	100	1	1	100	10	100	10	0	0	0	0
Combined	169	162	95.86	148	141	95.27	160	94.67	33	132	5	4	0
Eravacycline: Enterobacterales [Breakpoints (µg/mL): 0.5 (S), - (I), - (R)]													
Clinical	441	432	97.96	436	427	97.94	434	98.41	23	418	0	4	3
Challenge	46	46	100	39	39	100	46	100	31	15	0	0	0
Combined	487	478	98.15	475	466	98.11	480	98.56	54	433	0	4	3
Ertapenem: Enterobacterales [Breakpoints (µg/mL): 0.5 (S), 1 (I), 2 (R)]													
Clinical	880	832	94.55	123	75	60.98	857	97.39	31	842	9	14	0
Challenge	97	93	95.88	14	10	71.43	95	97.94	90	7	2	0	0
Combined	977	925	94.68	137	85	62.04	952	97.44	121	849	11	14	0
Gentamicin: Enterobacterales [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	745	713	95.7	704	672	95.45	735	98.66	60	679	7	1	2
Challenge	72	71	98.61	20	19	95	72	100	61	11	0	0	0
Combined	817	784	95.96	724	691	95.44	807	98.78	121	690	7	1	2
Gentamicin: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	208	183	87.98	200	175	87.5	202	97.12	7	195	6	0	0
Challenge	10	10	100	2	2	100	10	100	8	2	0	0	0
Combined	218	193	88.53	202	177	87.62	212	97.25	15	197	6	0	0
<sup>a</sup> Imipenem: <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 2 (S), 4 (I), 8 (R)]													
Clinical	121	107	88.43	77	63	81.82	119	98.35	57	64	0	0	2
Challenge	5	5	100	1	1	100	5	100	5	0	0	0	0
Combined	126	112	88.89	78	64	82.05	124	98.41	62	64	0	0	2
Imipenem-Relebactam: Enterobacterales [Breakpoints (µg/mL): 1 (S), 2 (I), 4 (R)]													
Clinical	433	406	93.76	427	400	93.68	422	97.46	4	428	8	3	0
Challenge	38	34	89.47	37	33	89.19	37	97.37	29	8	1	0	0
Combined	471	440	93.42	464	433	93.32	459	97.45	33	436	9	3	0
Imipenem-Relebactam: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 2 (S), 4 (I), 8 (R)]													
Clinical	158	152	96.2	158	152	96.2	155	98.1	3	153	3	0	0
Challenge	7	7	100	4	4	100	7	100	3	4	0	0	0
Combined	165	159	96.36	162	156	96.3	162	98.18	6	157	3	0	0
Levofloxacin: Enterobacterales [Breakpoints (µg/mL): 0.5 (S), 1 (I), 2 (R)]													
Clinical	732	701	95.77	235	204	86.81	692	94.54	112	600	33	6	1
Challenge	52	52	100	7	7	100	52	100	49	3	0	0	0
Combined	784	753	96.05	242	211	87.19	744	94.9	161	603	33	6	1
<sup>a</sup> Levofloxacin: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 1 (S), 2 (I), 4 (R)]													
Clinical	159	149	93.71	148	138	93.24	142	89.31	30	117	14	1	2
Challenge	12	12	100	3	3	100	12	100	12	0	0	0	0
Combined	171	161	94.15	151	141	93.38	154	90.06	42	117	14	1	2
<sup>a</sup> Meropenem: <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 2 (S), 4 (I), 8 (R)]													
Clinical	121	104	85.95	42	25	59.52	119	98.35	58	63	0	2	0
Challenge	5	5	100	0	0	NA	5	100	5	0	0	0	0
Combined	126	109	86.51	42	25	59.52	124	98.41	63	63	0	2	0
Meropenem: Enterobacterales [Breakpoints (µg/mL): 1 (S), 2 (I), 4 (R)]													
Clinical	805	776	96.4	54	25	46.3	789	98.01	25	777	4	12	0
Challenge	60	57	95	25	22	88	59	98.33	59	1	0	0	1
Combined	865	833	96.3	79	47	59.49	848	98.03	84	778	4	12	1
Meropenem: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 2 (S), 4 (I), 8 (R)]													
Clinical	159	148	93.08	131	120	91.6	152	95.6	23	131	7	0	0
Challenge	16	15	93.75	11	10	90.91	16	100	15	1	0	0	0

	Tot	No. EA	EA %	Eval EA Tot	No. Eval EA	Eval EA %	No. CA	CA %	No. R or NS	No. S	min	maj	vmj
Combined	175	163	93.14	142	130	91.55	168	96	38	132	7	0	0
Meropenem-vaborbactam: Enterobacterales [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	817	748	91.55	152	83	54.61	791	96.82	4	810	9	16	1
Challenge	43	41	95.35	12	10	83.33	42	97.67	37	6	1	0	0
Combined	860	789	91.74	164	93	56.71	833	96.86	41	816	10	16	1
<sup>a</sup> Minocycline: <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	196	174	88.78	96	74	77.08	178	90.82	38	148	13	5	0
Challenge	3	3	100	3	3	100	3	100	1	2	0	0	0
Combined	199	177	88.94	99	77	77.78	181	90.95	39	150	13	5	0
Minocycline: Enterobacterales [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	381	341	89.5	339	299	88.2	358	93.96	30	334	16	3	4
Challenge	23	21	91.3	15	13	86.67	23	100	13	9	0	0	0
Combined	404	362	89.6	354	312	88.14	381	94.31	43	343	16	3	4
<sup>a</sup> Piperacillin-Tazobactam: <i>Acinetobacter baumannii</i> [Breakpoints (µg/mL): 16(S), 32-64 (I), 128 (R)]													
Clinical	108	87	80.56	31	10	32.26	98	90.74	53	51	9	0	1
Challenge	5	5	100	1	1	100	5	100	5	0	0	0	0
Combined	113	92	81.42	32	11	34.38	103	91.15	58	51	9	0	1
Piperacillin-Tazobactam: Enterobacterales [Breakpoints (µg/mL): 16(S), 32-64 (I), 128 (R)]													
Clinical	815	728	89.33	641	554	86.43	768	94.23	56	730	35	8	4
Challenge	111	107	96.4	36	32	88.89	110	99.1	92	19	0	1	0
Combined	926	835	90.17	677	586	86.56	878	94.82	148	749	35	9	4
<sup>a</sup> Piperacillin-Tazobactam: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 16(S), 32-64 (I), 128 (R)]													
Clinical	211	170	80.57	208	167	80.29	177	83.89	22	167	26	2	6
Challenge	6	5	83.33	5	4	80	5	83.33	6	0	1	0	0
Combined	217	175	80.65	213	171	80.28	182	83.87	28	167	27	2	6
Tobramycin: <i>Pseudomonas aeruginosa</i> [Breakpoints (µg/mL): 4 (S), 8 (I), 16 (R)]													
Clinical	158	147	93.04	152	141	92.76	154	97.47	6	150	3	1	0
Challenge	8	8	100	2	2	100	8	100	8	0	0	0	0
Combined	166	155	93.37	154	143	92.86	162	97.59	14	150	3	1	0
Trimethoprim Sulfamethoxazole: Enterobacterales [Breakpoints (µg/mL): 2 (S), - (I), 4 (R)]													
Clinical	480	464	96.67	52	36	69.23	472	98.33	116	364	0	2	6
Challenge	46	46	100	2	2	100	46	100	43	3	0	0	0
Combined	526	510	96.96	54	38	70.37	518	98.48	159	367	0	2	6

<sup>a</sup> Due to unacceptable performance, do not report Selux MIC results for these drug/organism combinations.

EA – Essential Agreement  
CA – Category Agreement  
EVAL – Evaluable isolates

R – Resistant isolates  
NS – Non-susceptible isolates  
S – Susceptible isolates

min – minor errors  
maj – major errors  
vmj – very major errors

**Table 7. Selux AST System – Gram-Negative Panel Trending**

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
Amikacin	<i>Pseudomonas aeruginosa</i>	164	76 (46.34)	78	10 (6.1)	-40% (-48%, -31%)	Yes
Amoxicillin-Clavulanate	<i>Escherichia coli</i>	179	64 (35.75)	97	18 (10.06)	-26% (-34%, -17%)	No
Amoxicillin-Clavulanate	<i>Klebsiella oxytoca</i>	48	21 (43.75)	18	9 (18.75)	-25% (-41%, -6%)	No
Amoxicillin-Clavulanate	<i>Klebsiella pneumoniae</i>	142	89 (62.68)	46	7 (4.93)	-58% (-66%, -48%)	Yes
Amoxicillin-Clavulanate	<i>Proteus mirabilis</i>	88	11 (12.5)	42	35 (39.77)	27% (14%, 39%)	No
Ampicillin	<i>Escherichia coli</i>	84	64	18	2	-74%	Yes

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
			(76.19)		(2.38)	(-82%, -62%)	
Ampicillin	<i>Proteus mirabilis</i>	81	25 (30.86)	42	14 (17.28)	-14% (-26%, 0%)	No
Ampicillin-sulbactam	<i>Acinetobacter baumannii</i> (complex)	123	32 (26.02)	60	31 (25.2)	-1% (-12%, 10%)	No
Ampicillin-sulbactam	<i>Escherichia coli</i>	198	53 (26.77)	118	27 (13.64)	-13% (-21%, -5%)	No
Ampicillin-sulbactam	<i>Klebsiella oxytoca</i>	29	17 (58.62)	8	4 (13.79)	-45% (-63%, -20%)	Yes
Ampicillin-sulbactam	<i>Klebsiella pneumoniae</i>	116	50 (43.1)	35	31 (26.72)	-16% (-28%, -4%)	No
Ampicillin-sulbactam	<i>Morganella morganii</i>	75	9 (12)	47	19 (25.33)	13% (1%, 26%)	No
Ampicillin-sulbactam	<i>Proteus mirabilis</i>	88	11 (12.5)	58	19 (21.59)	9% (-2%, 20%)	No
Ampicillin-sulbactam	<i>Proteus vulgaris</i>	62	2 (3.23)	27	33 (53.23)	50% (35%, 62%)	Yes
Aztreonam	<i>Escherichia coli</i>	179	73 (40.78)	81	25 (13.97)	-27% (-35%, -18%)	No
Cefazolin	<i>Escherichia coli</i>	171	81 (47.37)	69	21 (12.28)	-35% (-44%, -26%)	Yes
Cefazolin	<i>Klebsiella pneumoniae</i>	113	67 (59.29)	29	17 (15.04)	-44% (-54%, -32%)	Yes
Cefepime	<i>Citrobacter freundii</i> (complex)	77	11 (14.29)	7	59 (76.62)	62% (48%, 72%)	Yes
Cefepime	<i>Citrobacter koseri</i>	72	0 (0)	0	72 (100)	100% (93%, 100%)	Yes
Cefepime	<i>Enterobacter cloacae</i> (complex)	77	13 (16.88)	9	55 (71.43)	55% (40%, 66%)	Yes
Cefepime	<i>Escherichia coli</i>	158	21 (13.29)	8	129 (81.65)	68% (59%, 75%)	Yes
Cefepime	<i>Klebsiella aerogenes</i>	41	2 (4.88)	0	39 (95.12)	90% (74%, 95%)	Yes
Cefepime	<i>Klebsiella oxytoca</i>	20	5 (25)	1	14 (70)	45% (14%, 66%)	Yes
Cefepime	<i>Klebsiella pneumoniae</i>	111	38 (34.23)	13	60 (54.05)	20% (7%, 32%)	No
Cefepime	<i>Morganella morganii</i>	17	1 (5.88)	0	16 (94.12)	88% (58%, 95%)	Yes
Cefepime	<i>Proteus mirabilis</i>	86	3 (3.49)	2	81 (94.19)	91% (81%, 95%)	Yes
Cefepime	<i>Proteus vulgaris</i>	90	4 (4.44)	4	82 (91.11)	87% (77%, 92%)	Yes
Cefepime	<i>Serratia marcescens</i>	98	3 (3.06)	2	93 (94.9)	92% (83%, 95%)	Yes
Cefoxitin	<i>Escherichia coli</i>	166	41 (24.7)	87	38 (22.89)	-2% (-11%, 7%)	No
Cefoxitin	<i>Klebsiella oxytoca</i>	29	6 (20.69)	16	7 (24.14)	3% (-18%, 24%)	No
Cefoxitin	<i>Klebsiella pneumoniae</i>	109	27 (24.77)	62	20 (18.35)	-6% (-17%, 5%)	No
Cefoxitin	<i>Morganella morganii</i>	59	10 (16.95)	39	10 (16.95)	0% (-14%, 14%)	No

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
Ceftazidime	<i>Citrobacter freundii</i> (complex)	53	12 (22.64)	12	29 (54.72)	32% (14%, 48%)	Yes
Ceftazidime	<i>Citrobacter koseri</i>	45	1 (2.22)	2	42 (93.33)	91% (77%, 96%)	Yes
Ceftazidime	<i>Enterobacter cloacae</i> (complex)	52	9 (17.31)	10	33 (63.46)	46% (28%, 60%)	Yes
Ceftazidime	<i>Escherichia coli</i>	74	31 (41.89)	9	34 (45.95)	4% (-12%, 20%)	No
Ceftazidime	<i>Klebsiella aerogenes</i>	29	2 (6.9)	7	20 (68.97)	62% (38%, 77%)	Yes
Ceftazidime	<i>Klebsiella oxytoca</i>	23	2 (8.7)	3	18 (78.26)	70% (42%, 83%)	Yes
Ceftazidime	<i>Klebsiella pneumoniae</i>	73	12 (16.44)	9	52 (71.23)	55% (40%, 66%)	Yes
Ceftazidime	<i>Proteus mirabilis</i>	63	1 (1.59)	4	58 (92.06)	90% (79%, 95%)	Yes
Ceftazidime	<i>Proteus vulgaris</i>	50	2 (4)	1	47 (94)	90% (76%, 95%)	Yes
Ceftazidime	<i>Serratia marcescens</i>	56	3 (5.36)	9	44 (78.57)	73% (58%, 83%)	Yes
Ceftazidime-avibactam	<i>Citrobacter freundii</i> (complex)	61	11 (18.03)	15	35 (57.38)	39% (22%, 53%)	Yes
Ceftazidime-avibactam	<i>Citrobacter koseri</i>	54	0 (0)	6	48 (88.89)	89% (76%, 95%)	Yes
Ceftazidime-avibactam	<i>Enterobacter cloacae</i> (complex)	75	19 (25.33)	27	29 (38.67)	13% (-2%, 27%)	No
Ceftazidime-avibactam	<i>Escherichia coli</i>	136	15 (11.03)	20	101 (74.26)	63% (53%, 71%)	Yes
Ceftazidime-avibactam	<i>Klebsiella aerogenes</i>	52	2 (3.85)	11	39 (75)	71% (55%, 81%)	Yes
Ceftazidime-avibactam	<i>Klebsiella oxytoca</i>	28	3 (10.71)	4	21 (75)	64% (40%, 78%)	Yes
Ceftazidime-avibactam	<i>Klebsiella pneumoniae</i>	73	23 (31.51)	19	31 (42.47)	11% (-5%, 26%)	No
Ceftazidime-avibactam	<i>Morganella morganii</i>	62	1 (1.61)	3	58 (93.55)	92% (81%, 96%)	Yes
Ceftazidime-avibactam	<i>Proteus mirabilis</i>	65	0 (0)	2	63 (96.92)	97% (88%, 99%)	Yes
Ceftazidime-avibactam	<i>Proteus vulgaris</i>	55	0 (0)	2	53 (96.36)	96% (85%, 99%)	Yes
Ceftazidime-avibactam	<i>Pseudomonas aeruginosa</i>	158	70 (44.3)	63	25 (15.82)	-28% (-38%, -19%)	No
Ceftazidime-avibactam	<i>Serratia marcescens</i>	66	10 (15.15)	33	23 (34.85)	20% (5%, 33%)	No
Ceftriaxone	<i>Citrobacter freundii</i> (complex)	49	8 (16.33)	7	34 (69.39)	53% (34%, 67%)	Yes
Ceftriaxone	<i>Citrobacter koseri</i>	27	1 (3.7)	0	26 (96.3)	93% (72%, 97%)	Yes
Ceftriaxone	<i>Escherichia coli</i>	44	7 (15.91)	1	36 (81.82)	66% (47%, 78%)	Yes
Ceftriaxone	<i>Klebsiella aerogenes</i>	34	6 (17.65)	3	25 (73.53)	56% (33%, 71%)	Yes
Ceftriaxone	<i>Klebsiella oxytoca</i>	29	5 (17.24)	1	23 (79.31)	62% (37%, 77%)	Yes

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
Ceftriaxone	<i>Klebsiella pneumoniae</i>	107	52 (48.6)	2	53 (49.53)	1% (-12%, 14%)	No
Ceftriaxone	<i>Proteus mirabilis</i>	37	5 (13.51)	0	32 (86.49)	73% (53%, 84%)	Yes
Ciprofloxacin	<i>Citrobacter freundii</i> (complex)	57	14 (24.56)	7	36 (63.16)	39% (21%, 53%)	Yes
Ciprofloxacin	<i>Citrobacter koseri</i>	36	1 (2.78)	1	34 (94.44)	92% (75%, 96%)	Yes
Ciprofloxacin	<i>Enterobacter cloacae</i> (complex)	68	14 (20.59)	5	49 (72.06)	51% (35%, 64%)	Yes
Ciprofloxacin	<i>Escherichia coli</i>	170	51 (30)	17	102 (60)	30% (20%, 40%)	Yes
Ciprofloxacin	<i>Klebsiella aerogenes</i>	30	0 (0)	2	28 (93.33)	93% (75%, 98%)	Yes
Ciprofloxacin	<i>Klebsiella oxytoca</i>	28	1 (3.57)	6	21 (75)	71% (48%, 84%)	Yes
Ciprofloxacin	<i>Klebsiella pneumoniae</i>	23	7 (30.43)	5	11 (47.83)	17% (-10%, 42%)	No
Ciprofloxacin	<i>Morganella morganii</i>	58	3 (5.17)	6	49 (84.48)	79% (65%, 87%)	Yes
Ciprofloxacin	<i>Proteus mirabilis</i>	24	9 (37.5)	5	10 (41.67)	4% (-22%, 30%)	No
Ciprofloxacin	<i>Proteus vulgaris</i>	17	3 (17.65)	3	11 (64.71)	47% (14%, 68%)	Yes
Ciprofloxacin	<i>Pseudomonas aeruginosa</i>	151	62 (41.06)	74	15 (9.93)	-31% (-40%, -22%)	Yes
Eravacycline	<i>Citrobacter freundii</i> (complex)	62	2 (3.23)	32	28 (45.16)	42% (28%, 54%)	Yes
Eravacycline	<i>Enterobacter cloacae</i> (complex)	76	13 (17.11)	47	16 (21.05)	4% (-9%, 16%)	No
Eravacycline	<i>Escherichia coli</i>	167	40 (23.95)	99	28 (16.77)	-7% (-16%, 1%)	No
Eravacycline	<i>Klebsiella oxytoca</i>	47	4 (8.51)	30	13 (27.66)	19% (3%, 34%)	No
Ertapenem	<i>Citrobacter freundii</i> (complex)	24	3 (12.5)	5	16 (66.67)	54% (27%, 72%)	Yes
Ertapenem	<i>Citrobacter koseri</i>	72	0 (0)	0	72 (100)	100% (93%, 100%)	Yes
Ertapenem	<i>Escherichia coli</i>	210	11 (5.24)	11	188 (89.52)	84% (78%, 88%)	Yes
Ertapenem	<i>Klebsiella oxytoca</i>	3	1 (33.33)	1	1 (33.33)	0% (-53%, 53%)	No
Ertapenem	<i>Klebsiella pneumoniae</i>	23	10 (43.48)	6	7 (30.43)	-13% (-38%, 14%)	No
Ertapenem	<i>Morganella morganii</i>	10	7 (70)	2	1 (10)	-60% (-81%, -17%)	Yes
Ertapenem	<i>Proteus mirabilis</i>	85	2 (2.35)	0	83 (97.65)	95% (87%, 98%)	Yes
Ertapenem	<i>Proteus vulgaris</i>	98	0 (0)	1	97 (98.98)	99% (93%, 100%)	Yes
Ertapenem	<i>Serratia marcescens</i>	65	6 (9.23)	11	48 (73.85)	65% (49%, 75%)	Yes
Gentamicin	<i>Citrobacter freundii</i> (complex)	64	41 (64.06)	22	1 (1.56)	-62% (-73%, -49%)	Yes

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
Gentamicin	<i>Citrobacter koseri</i>	54	23 (42.59)	23	8 (14.81)	-28% (-43%, -11%)	No
Gentamicin	<i>Enterobacter cloacae</i> (complex)	74	36 (48.65)	31	7 (9.46)	-39% (-51%, -25%)	Yes
Gentamicin	<i>Escherichia coli</i>	173	84 (48.55)	78	11 (6.36)	-42% (-50%, -34%)	Yes
Gentamicin	<i>Klebsiella aerogenes</i>	31	8 (25.81)	18	5 (16.13)	-10% (-29%, 11%)	No
Gentamicin	<i>Klebsiella oxytoca</i>	29	14 (48.28)	15	0 (0)	-48% (-66%, -28%)	Yes
Gentamicin	<i>Klebsiella pneumoniae</i>	126	58 (46.03)	62	6 (4.76)	-41% (-50%, -31%)	Yes
Gentamicin	<i>Proteus mirabilis</i>	64	12 (18.75)	29	23 (35.94)	17% (2%, 32%)	No
Gentamicin	<i>Proteus vulgaris</i>	55	13 (23.64)	27	15 (27.27)	4% (-13%, 20%)	No
Gentamicin	<i>Pseudomonas aeruginosa</i>	211	68 (32.23)	79	64 (30.33)	-2% (-11%, 7%)	No
Gentamicin	<i>Serratia marcescens</i>	68	3 (4.41)	23	42 (61.76)	57% (43%, 68%)	Yes
Imipenem-Relebactam	<i>Citrobacter freundii</i> (complex)	54	28 (51.85)	19	7 (12.96)	-39% (-53%, -22%)	Yes
Imipenem-Relebactam	<i>Citrobacter koseri</i>	55	18 (32.73)	24	13 (23.64)	-9% (-25%, 8%)	No
Imipenem-Relebactam	<i>Escherichia coli</i>	172	40 (23.26)	101	31 (18.02)	-5% (-14%, 3%)	No
Imipenem-Relebactam	<i>Klebsiella oxytoca</i>	27	0 (0)	19	8 (29.63)	30% (11%, 48%)	Yes
Imipenem-Relebactam	<i>Pseudomonas aeruginosa</i>	162	22 (13.58)	119	21 (12.96)	-1% (-8%, 7%)	No
Levofloxacin	<i>Citrobacter freundii</i> (complex)	60	6 (10)	12	42 (70)	60% (44%, 71%)	Yes
Levofloxacin	<i>Citrobacter koseri</i>	54	0 (0)	0	54 (100)	100% (91%, 100%)	Yes
Levofloxacin	<i>Enterobacter cloacae</i> (complex)	68	0 (0)	5	63 (92.65)	93% (82%, 97%)	Yes
Levofloxacin	<i>Escherichia coli</i>	169	13 (7.69)	31	125 (73.96)	66% (58%, 73%)	Yes
Levofloxacin	<i>Klebsiella aerogenes</i>	32	1 (3.12)	1	30 (93.75)	91% (72%, 96%)	Yes
Levofloxacin	<i>Klebsiella oxytoca</i>	27	1 (3.7)	4	22 (81.48)	78% (54%, 89%)	Yes
Levofloxacin	<i>Klebsiella pneumoniae</i>	84	9 (10.71)	13	62 (73.81)	63% (50%, 73%)	Yes
Levofloxacin	<i>Morganella morganii</i>	61	4 (6.56)	3	54 (88.52)	82% (68%, 89%)	Yes
Levofloxacin	<i>Proteus mirabilis</i>	65	3 (4.62)	9	53 (81.54)	77% (63%, 85%)	Yes
Levofloxacin	<i>Proteus vulgaris</i>	55	1 (1.82)	6	48 (87.27)	85% (72%, 92%)	Yes
Levofloxacin	<i>Serratia marcescens</i>	68	15 (22.06)	34	19 (27.94)	6% (-9%, 20%)	No
Meropenem	<i>Citrobacter freundii</i> (complex)	78	2 (2.56)	1	75 (96.15)	94% (84%, 97%)	Yes

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
Meropenem	<i>Citrobacter koseri</i>	72	0 (0)	1	71 (98.61)	99% (91%, 100%)	Yes
Meropenem	<i>Enterobacter cloacae</i> (complex)	74	5 (6.76)	6	63 (85.14)	78% (66%, 86%)	Yes
Meropenem	<i>Escherichia coli</i>	5	2 (40)	0	3 (60)	20% (-32%, 60%)	No
Meropenem	<i>Klebsiella oxytoca</i>	27	0 (0)	0	27 (100)	100% (82%, 100%)	Yes
Meropenem	<i>Klebsiella pneumoniae</i>	106	14 (13.21)	15	77 (72.64)	59% (47%, 69%)	Yes
Meropenem	<i>Morganella morganii</i>	62	1 (1.61)	1	60 (96.77)	95% (85%, 98%)	Yes
Meropenem	<i>Proteus mirabilis</i>	66	1 (1.52)	3	62 (93.94)	92% (82%, 96%)	Yes
Meropenem	<i>Proteus vulgaris</i>	51	0 (0)	0	51 (100)	100% (90%, 100%)	Yes
Meropenem	<i>Pseudomonas aeruginosa</i>	167	52 (31.14)	77	38 (22.75)	-8% (-18%, 1%)	No
Meropenem	<i>Serratia marcescens</i>	102	3 (2.94)	2	97 (95.1)	92% (84%, 96%)	Yes
Meropenem-vaborbactam	<i>Citrobacter freundii</i> (complex)	57	0 (0)	1	56 (98.25)	98% (88%, 100%)	Yes
Meropenem-vaborbactam	<i>Citrobacter koseri</i>	48	0 (0)	0	48 (100)	100% (90%, 100%)	Yes
Meropenem-vaborbactam	<i>Enterobacter cloacae</i> (complex)	52	0 (0)	1	51 (98.08)	98% (87%, 100%)	Yes
Meropenem-vaborbactam	<i>Escherichia coli</i>	29	4 (13.79)	0	25 (86.21)	72% (49%, 84%)	Yes
Meropenem-vaborbactam	<i>Klebsiella aerogenes</i>	50	0 (0)	1	49 (98)	98% (87%, 100%)	Yes
Meropenem-vaborbactam	<i>Klebsiella oxytoca</i>	47	0 (0)	0	47 (100)	100% (89%, 100%)	Yes
Meropenem-vaborbactam	<i>Klebsiella pneumoniae</i>	100	4 (4)	6	90 (90)	86% (77%, 91%)	Yes
Meropenem-vaborbactam	<i>Morganella morganii</i>	64	0 (0)	35	29 (45.31)	45% (32%, 57%)	Yes
Meropenem-vaborbactam	<i>Serratia marcescens</i>	60	0 (0)	1	59 (98.33)	98% (89%, 100%)	Yes
Minocycline	<i>Escherichia coli</i>	198	43 (21.72)	83	72 (36.36)	15% (6%, 23%)	No
Minocycline	<i>Klebsiella aerogenes</i>	52	17 (32.69)	30	5 (9.62)	-23% (-38%, -7%)	No
Minocycline	<i>Klebsiella oxytoca</i>	27	3 (11.11)	15	9 (33.33)	22% (0%, 42%)	No
Minocycline	<i>Klebsiella pneumoniae</i>	126	38 (30.16)	66	22 (17.46)	-13% (-23%, -2%)	No
Piperacillin-tazobactam	<i>Citrobacter koseri</i>	75	15 (20)	49	11 (14.67)	-5% (-17%, 7%)	No
Piperacillin-tazobactam	<i>Escherichia coli</i>	175	47 (26.86)	91	37 (21.14)	-6% (-15%, 3%)	No
Piperacillin-tazobactam	<i>Klebsiella oxytoca</i>	27	11 (40.74)	15	1 (3.7)	-37% (-56%, -15%)	Yes
Piperacillin-tazobactam	<i>Klebsiella pneumoniae</i>	132	80 (60.61)	37	15 (11.36)	-49% (-58%, -38%)	Yes

Drug	Organism Name	Total On Scale for Trending	≥ 1 Dilution Lower # (%)	Exact # (%)	≥ 1 Dilution Higher # (%)	Percent Difference (95% CI)	Significant Trending Noted <sup>a</sup>
Piperacillin-tazobactam	<i>Morganella morganii</i>	64	12 (18.75)	16	36 (56.25)	38% (21%, 51%)	Yes
Piperacillin-tazobactam	<i>Proteus mirabilis</i>	70	15 (21.43)	12	43 (61.43)	40% (24%, 53%)	Yes
Piperacillin-tazobactam	<i>Proteus vulgaris</i>	50	0 (0)	16	34 (68)	68% (52%, 79%)	Yes
Tobramycin	<i>Pseudomonas aeruginosa</i>	162	25 (15.43)	91	46 (28.4)	13% (4%, 22%)	No
Trimethoprim Sulfamethoxazole	<i>Enterobacter cloacae</i> (complex)	72	17 (23.61)	11	44 (61.11)	38% (22%, 51%)	Yes
Trimethoprim Sulfamethoxazole	<i>Klebsiella aerogenes</i>	32	2 (6.25)	1	29 (90.62)	84% (64%, 92%)	Yes
Trimethoprim Sulfamethoxazole	<i>Klebsiella oxytoca</i>	22	0 (0)	0	22 (100)	100% (79%, 100%)	Yes
Trimethoprim Sulfamethoxazole	<i>Klebsiella pneumoniae</i>	66	6 (9.09)	5	55 (83.33)	74% (60%, 83%)	Yes

<sup>a</sup>Significant trending as determined by a statistically significant confidence interval.

**Characterized Resistance Mechanisms.** The molecular characteristics of the challenge isolates evaluated in the clinical study, provided by the FDA-CDC AR Isolate Bank, are provided in **Table 8**.



**Table 8. Molecular characteristics of resistance markers in challenge isolates.**

Antimicrobial	Drug-Class Specific	Truncated Porins
Amikacin	aac(3)-IIa, aac(3)-Id, aac(6)-II, aadA6, ant(2'')-Ia, aph(3')-Ic, aph(3')-VIa, armA, strA, strB, aac(3)-Ia	
Amoxicillin-Clavulanate	CMY-2, CMY-4, CMY-42, CMY-6, CTX-M-15, CTX-M-2, KPC-2, KPC-3, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-232, OXA-9, SHV-1, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-30, SHV-83, TEM-1A, TEM-1B, VIM-1, VIM-27	OmpC2, OmpK35, OmpK36, OmpF
Ampicillin	CMY-42, CTX-M-15, KPC-3, NDM-5, OXA-1, OXA-9, SHV-12, TEM-1A, TEM-1B	
Ampicillin-Sulbactam	CMY-4, CMY-6, CMY-42, CTX-M-2, CTX-M-15, KPC-3, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-2, OXA-23, OXA-232, OXA-69, OXA-9, SHV-11, SHV-12, SHV-83, TEM-1A, TEM-1B, TEM-1D, VIM-27	OmpF, OmpK35, OmpK36
Aztreonam	carB-2, CMY-2, CMY-4, CMY-42, CMY-6, CTX-M-15, KPC-2, KPC-3, NDM-1, NDM-5, OXA-1, OXA-2, OXA-10, OXA-50, OXA-9, SHV-12, TEM-1A, TEM-1B, VEB-1, VEB-5	OmpF
Cefazolin	CMY-2, CMY-4, CMY-42, CMY-6, CTX-M-15, CTX-M-2, KPC-2, NDM-1, OXA-1, OXA-10, OXA-181, OXA-232, OXA-9, SHV-1, SHV-11, SHV-12, SHV-154, SHV-26, SHV-83, TEM-1A, TEM-1B, VEB-5	OmpF, OmpK35, OmpK36
Cefepime	ACT-7, ACT-16, carB-2, CMY-42, CMY-6, CTX-M-15, CTX-M-2, IMP-1, KPC-2, KPC-6, NDM-1, OXA-1, OXA-10, OXA-181, OXA-2, OXA-232, OXA-50, OXA-9, PAO, SHV-1, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-83, TEM-1A, TEM-1B, VEB-1, VEB-5, VIM-1	OmpC2, OmpK35, OmpK36, OmpF
Cefoxitin	CMY-4, CMY-42, CMY-6, CTX-M-15, CTX-M-2, KPC-2, NDM-1, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-232, OXA-9, SHV-1, SHV-100, SHV-11, SHV-12, SHV-28, SHV-83, TEM-1A, TEM-1B	OmpC2, OmpK35, OmpK36, OmpF
Ceftazidime	ACT-7, ACT-16, ADC-25, CMY-2, CMY-4, CMY-42, CMY-6, CTX-M-15, CTX-M-2, IMP-1, IMP-8, KPC-2, KPC-3, KPC-6, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-23, OXA-232, OXA-24, OXA-50, OXA-65, OXA-66, OXA-69, OXA-71, OXA-9, OXA-94, PAO, SHV-1, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-30, SHV-5, SHV-83, SME-3, TEM-1A, TEM-1B, TEM-1D, VEB-1, VEB-5, VIM-1, VIM-2, VIM-27	OmpC2, OmpF, OmpK35, OmpK36
Ceftazidime-Avibactam	ACT-7, CMY-4, CMY-42, CMY-6, CTX-M-15, IMP-1, IMP-8, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-101, OXA-181, OXA-2, OXA-232, OXA-50, OXA-9, PAO, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-30, SME-3, TEM-1A, TEM-1B, VEB-1, VIM-1, VIM-2, VIM-27	OmpC2, OmpF, OmpK35
Ceftriaxone	ACT-7, CMY-2, CMY-4, CMY-42, CTX-M-15, CTX-M-2, KPC-2, NDM-1, NDM-7, OXA-1, OXA-10, OXA-181, OXA-232, OXA-9, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-83, SME-3, TEM-1A, TEM-1B, VEB-5, VIM-27	OmpC2, OmpK35, OmpK36, OmpF
Ciprofloxacin	QnrB1, QnrB7, QnrS1, oqxA, oqxB, QnrB34	OmpC2, OmpK35, OmpK36, OmpF
Eravacycline	tet(A), tet(D), tet(G), tet(R)	OmpC2, OmpK35
Ertapenem	ACT-16, ACT-7, CMY-2, CMY-4, CMY-42, CMY-6, CTX-M-15, CTX-M-2, IMP-8, KPC-2, KPC-6, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-232, OXA-9, SHV-1, SHV-100, SHV-11, SHV-12, SHV-154, SHV-26, SHV-28, SHV-83, SME-3, TEM-1A, TEM-1B, VEB-5, VIM-27	OmpC2, OmpK35, OmpK36, OmpF
Gentamicin	aac(3)-IIa, aac(3)-IId, aac(3)-Id, aac(6)-IIa, aac(6)-IIc, aac(6)-IIan, aac(6)-Ib-cr, aac(6)-II, aadA1, aadA2, aadA5, aadA6, aadB, ant(2'')-Ia, aph(3')-Ia, aph(3')-Ib, aph(3')-VI, aph(3')-XV, aph(6)-Id, armA, rmtB, rmtC, rmtG, strA, strB, aac(6)-Ib	OmpC2, OmpK35, OmpK36, OmpF
Imipenem	ADC-25, OXA-23, OXA-24, OXA-65, OXA-66, OXA-69, OXA-71, SHV-5, TEM-1D	
Imipenem-Relebactam	BCR1, CMY-4, CMY-42, CMY-6, CTX-M-15, IMP-1, KPC-2, KPC-3, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-232, OXA-50, OXA-9, PAO, SHV-1, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-30, TEM-1A, TEM-1B, VIM-1, VIM-2, VIM-27	OmpC2, OmpK35, OmpK36, OmpF
Levofloxacin	QnrB1, QnrB34, QnrB7, QnrS1, oqxA, oqxB, norA	OmpC2, OmpF, OmpK35, OmpK36
Meropenem	ACT-7, ADC-25, BCR1, carB-2, CMY-4, CMY-42, CMY-6, CTX-M-15, IMP-1, KPC-2, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-101, OXA-181, OXA-2, OXA-23, OXA-232, OXA-24, OXA-50, OXA-65, OXA-66, OXA-69, OXA-71, OXA-9, OXA-94, PAO, SHV-1, SHV-100, SHV-11, SHV-12, SHV-28, SHV-5, SME-3, TEM-1A, TEM-1B, TEM-1D, VEB-1, VIM-2, VIM-27	OmpC2, OmpF, OmpK35
Meropenem-Vaborbactam	CMY-4, CMY-42, CMY-6, CTX-M-15, KPC-2, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-232, OXA-9, SHV-1, SHV-100, SHV-11, SHV-12, SHV-28, TEM-1A, TEM-1B, VIM-27	OmpC2, OmpK35, OmpK36, OmpF
Minocycline	tet(A), tet(A), tet(D), tet(K), tet(R), tet(38)	OmpK35
Piperacillin-Tazobactam	ACT-7, ACT-16, ADC-25, carB-2, CMY-2, CMY-4, CMY-42, CMY-6, CTX-M-15, KPC-2, KPC-3, KPC-6, NDM-1, NDM-5, NDM-7, OXA-1, OXA-10, OXA-181, OXA-2, OXA-23, OXA-232, OXA-24, OXA-50, OXA-65, OXA-69, OXA-71, OXA-9, OXA-94, SHV-1, SHV-100, SHV-11, SHV-12, SHV-26, SHV-28, SHV-5, TEM-1A, TEM-1B, TEM-1D, VIM-27	OmpC2, OmpF, OmpK35
Tobramycin	aac(6)-II, aadA1, aadA6, aadB, ant(2'')-Ia, strA, strB, aac(3)-Id	
Trimethoprim-Sulfamethoxazole	dfrA1, dfrA12, dfrA14, dfrA15, dfrA17, dfrA29, dfrB5	OmpC2, OmpK35, OmpK36, OmpF

**Non-indicated species.** As required under 511A(b)(2)(C)(ii)(I) of the Federal Food, Drug and Cosmetic Act, the following statement is included in the Precautions section of the device labeling to address testing and reporting of non-indicated species:

Per the FDA-Recognized Susceptibility Test Interpretive Criteria Website, the safety and efficacy of antimicrobial drugs, for which antimicrobial susceptibility is tested by this AST device, may or may not have been established in adequate and well-controlled clinical trials for treating infections due to microorganisms outside of those found in the indications and usage in the drug label. The clinical significance of susceptibility information in those instances is unknown. The approved labeling for specific antimicrobial drugs provides the uses for which the antimicrobial drug is approved.

- Matrix Comparison:  
Not applicable.

**C Clinical Studies:**

1. Clinical Sensitivity:  
Not applicable.
2. Clinical Specificity:  
Not applicable.
3. Other Clinical Supportive Data (When 1. and 2. Are Not Applicable):  
Not applicable.

**D Clinical Cut-Off:**

Not applicable.

**E Expected Values/Reference Range:**

The FDA-recognized/approved susceptibility test interpretive criteria (breakpoints) for the antimicrobials evaluated with the Selux AST System are listed in the table below.

**Table 9.** FDA-Approved or Recognized Interpretive Criteria<sup>1</sup>

Antimicrobial	Organism	Minimum Inhibitory Concentration (µg/mL)		
		S	I	R
Amikacin	<i>Pseudomonas aeruginosa</i>	≤16	32	≥64
Amoxicillin-Clavulanate	Enterobacterales	≤8	16	≥32
Ampicillin	Enterobacterales	≤8	16	≥32
Ampicillin-sulbactam	<i>Acinetobacter</i> spp., Enterobacterales	≤8	16	≥32
Aztreonam	Enterobacterales	≤4	8	≥16
Cefazolin	Enterobacterales	≤1	2	≥4
Cefepime	Enterobacterales	≤2	4-8	≥16
Cefoxitin	Enterobacterales	≤4	8	≥16
Ceftazidime	Enterobacterales	≤4	8	≥16
Ceftazidime-Avibactam	Enterobacterales, <i>Pseudomonas aeruginosa</i>	≤8	-	≥16
Ceftriaxone	Enterobacterales	≤1	2	≥4
Ciprofloxacin	Enterobacterales	≤0.25	0.5	≥1
Ciprofloxacin	<i>Pseudomonas aeruginosa</i>	≤0.5	1	≥2
Eravacycline	Enterobacterales	≤0.5	-	-
Ertapenem	Enterobacterales	≤0.5	1	≥2
Gentamicin	Enterobacterales, <i>Pseudomonas aeruginosa</i>	≤4	8	≥16
Imipenem-Relebactam	Enterobacterales	≤1	2	≥4
Imipenem-Relebactam	<i>Pseudomonas aeruginosa</i>	≤2	4	≥8

Antimicrobial	Organism	Minimum Inhibitory Concentration (µg/mL)		
		S	I	R
Levofloxacin	Enterobacterales	≤0.5	1	≥2
Meropenem	Enterobacterales	≤1	2	≥4
Meropenem	<i>Pseudomonas aeruginosa</i>	≤2	4	≥8
Meropenem-Vaborbactam	Enterobacterales	≤4	8	≥16
Minocycline	Enterobacterales	≤4	8	≥16
Piperacillin-tazobactam	<i>Pseudomonas aeruginosa</i>	≤16	32-64	≥128
Tobramycin	<i>Pseudomonas aeruginosa</i>	≤4	8	≥16
Trimethoprim-Sulfamethoxazole	Enterobacterales	≤2	-	≥4

S = Susceptible; I = Intermediate; R = Resistant; - = no interpretive criterion recognized

<sup>1</sup> According to the FDA-Recognized Antimicrobial Susceptibility Test Interpretive Criteria Website

<https://www.fda.gov/Drugs/DevelopmentApprovalProcess/DevelopmentResources/ucm410971.htm>

### VIII Proposed Labeling:

The labeling supports the finding of substantial equivalence for this device.

### IX Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.

To support the implementation of changes to FDA-recognized susceptibility test interpretive criteria (i.e., breakpoints), this submission included a breakpoint change protocol that was reviewed and accepted by FDA. This protocol addresses future revisions to device labeling in response to breakpoint changes that are recognized on the FDA STIC webpage (<https://www.fda.gov/drugs/development-resources/fda-recognized-antimicrobial-susceptibility-test-interpretive-criteria>). The protocol outlined the specific procedures and acceptance criteria that Selux intends to use to evaluate the Selux AST System when revised breakpoints for the following antimicrobials are published on the FDA STIC webpage: Amikacin, Amoxicillin-Clavulanate, Ampicillin, Ampicillin-sulbactam, Aztreonam, Cefazolin, Cefepime, Cefoxitin, Ceftazidime, Ceftazidime-Avibactam, Ceftriaxone, Ciprofloxacin, Eravacycline, Ertapenem, Gentamicin, Imipenem-Relebactam, Levofloxacin, Meropenem, Meropenem-Vaborbactam, Minocycline, Piperacillin-Tazobactam, Tobramycin, Trimethoprim-Sulfamethoxazole.

The breakpoint change protocol included with the submission indicated that if specific criteria are met, Selux will update the device label to include (1) the new breakpoints, (2) an updated performance section after re-evaluation of data in this premarket notification with the new breakpoints, and (3) any new limitations as determined by their evaluation.