



APR 19 1996

* 510(k) SUMMARY *

Date Prepared: December 4, 1995

Contact Person: Eric S. Hoy, Ph.D.

Name of Device:

- Trade Name - RELISA® Jo-1 Antibody Test System
- Common Name - Jo-1 Antibody Test System
- Classification Name - Extractable Antinuclear Antibody (21 CFR 866.5100)

Legally marketed device with which this device has been shown to be equivalent:

RELISA® ENA Antibody Screening Tests System, K935129

Description:

This is an enzyme immunoassay for the detection of antibodies to nuclear antigen Jo-1 in human serum.

Intended Use:

This test system is for in vitro diagnostic use for the detection of antibodies to nuclear antigen Jo-1 in human serum.

Summary of Technological Characteristics Compared to the Predicate Device:

This device is identical to the predicate device with the following exceptions:

- a) The predicate device has six different autoantigens coated on individual microwells; the present device has only Jo-1 autoantigen coated on the microwells.
- b) The predicate device includes a procedure control well on each strip of microwells, the present device includes a calibrator serum in the kit.

Description of Laboratory Data That Indicate Substantial Equivalence:

For direct determination of relative sensitivity and specificity, we used the Immuno Concepts RELISA® Screening Assay (K935129) as a reference method. The data obtained in this comparison are shown in the following Table.

Table 1. Detection of antibodies to the Jo-1 autoantigen.

		Immuno Concepts RELISA® Screening Assay		
		Positive	Borderline	Negative
Immuno Concepts RELISA® Jo-1	Positive	21	1	0
	Borderline	2	4	1
	Negative	0	0	140

If we assume that "borderline" results are actually positive, these data yield the following statistics:
relative sensitivity, 100.0%; relative specificity, 99.3%; and overall agreement, 99.4%

In accordance with 21 CFR 807.92(b)(3), we conclude from these data that the present device is substantially equivalent to the predicate device.