

K961462

510(k) SUMMARY

JUN 27 1996

<b>Subject:</b>	Premarket Notification, 510(k) CEDIA® Theophylline Assay application: BM/Hitachi 704 Safety and Effectiveness Summary
<b>Manufacturer:</b>	A Boehringer Mannheim Corporation 2400 Bisso Lane P.O. Box 4117 Concord CA 94524-4117 Phone (510) 674 - 0667 Fax 510 674 1850
<b>Contact:</b>	Betsy Soares-Maddox, Manager of Regulatory Affairs and Quality Assurance
<b>Date:</b>	14 June 1996
<b>Proprietary Name</b>	CEDIA® Theophylline Assay
<b>Common Name</b>	Homogeneous Enzyme Immunoassay for the Determination of Theophylline Levels in Serum and Plasma.
<b>Classification Name</b>	Theophylline Test System
<b>Predicate Device</b>	CEDIA® Theophylline Assay: BM/Hitachi 911
<b>Device Description</b>	CEDIA® Technology  The CEDIA® Theophylline Assay is an in-vitro homogeneous enzyme immunoassay used for the measurement of theophylline in serum and plasma. It is based on competitive binding concepts employing theophylline labeled enzymatic fragments ( $\beta$ -galactosidase) competing with sample theophylline for the theophylline -specific antibody.

## 510(k) SUMMARY cont.

**Device Description cont.** Using recombinant DNA techniques, the  $\beta$ -galactosidase molecule has been split into two totally inactive polypeptide subunits called enzyme acceptor and enzyme donor. Theophylline has been covalently linked to the enzyme donor in a manner that does not prevent spontaneous reassociation of the subunits to yield active  $\beta$ -galactosidase enzyme. Theophylline -specific antibody, by binding to the Theophylline derivative on the enzyme donor will inhibit enzyme reassociation, thereby regulating the level of  $\beta$ -galactosidase formed. The amount of enzyme formed is proportional to the amount of theophylline as monitored by the hydrolysis of the substrate chlorophenol red- $\beta$ -D-galactopyranoside (CPRG).

**Intended Use** The CEDIA Theophylline Assay is a homogeneous enzyme immunoassay for the in vitro assay of theophylline in human serum and plasma. Measurements are used in the diagnosis and treatment of theophylline overdose and in monitoring levels of theophylline to ensure proper therapy.

**Statement of Similarities and Differences** The following table outlines the similarities and differences between the CEDIA Theophylline Assay on the BM/Hitachi 911 to the BM/Hitachi 704.

Parameter	BM/Hitachi 911	BM/Hitachi 704
<b>Methodology</b>	Homogeneous Enzyme Immunoassay	Homogeneous Enzyme Immunoassay
<b>Intended Use</b>	Quantitative Determination of Theophylline in human serum and plasma	Quantitative Determination of Theophylline in human serum and plasma
<b>Detection Mechanism</b>	Spectrophotometer at 570 nm	Spectrophotometer at 570 nm
<b>Reaction Volumes</b>		
<b>Sample</b>	3 $\mu$ L	4 $\mu$ L
<b>R 1</b>	130 $\mu$ L	185 $\mu$ L
<b>R 2</b>	130 $\mu$ L	170 $\mu$ L

# 510(k) SUMMARY cont.

Statement of Similarities and Differences, cont.						
Parameter	BM/Hitachi 911			BM/Hitachi 704		
Reagents	<p>Enzyme Acceptor lyophilized with buffer salts, bulking agent, detergent and preservative.</p> <p>Enzyme Donor lyophilized with substrate, stabilizer and preservative.</p> <p>Enzyme Acceptor Reconstitution Buffer with primary antibody, buffer salts, monoclonal anti-Theophylline antibody, stabilizer and preservative.</p> <p>Enzyme Donor Reconstitution Buffer with buffer salts and preservative.</p>			Same		
Sensitivity (LDD)	0.8 µg/mL			0.8 µg/mL		
Precision	Dose, µg/mL: NCCLS modified			Dose, µg/mL: NCCLS modified		
Control Level	Low	Mid	High	Low	Mid	High
Within-Run	5.1	15.1	29.3	4.7	15.1	29.2
%CV	3.3	1.9	1.3	5.5	2.8	2.1
Total	5.1	15.1	29.3	4.7	15.1	29.2
%CV	5.1	2.4	2.0	6.4	3.2	2.5

# 510(k) SUMMARY cont.

## Statement of Similarities and Differences, cont.

Parameter	BM/Hitachi 911	BM/Hitachi 704
<b>Method Comparison</b>		
<b>Versus:</b>	Fluorescence Polarization Immunoassay	BM/Hitachi 911
<b>Slope</b>	1.01	1.09
<b>Intercept</b>	-0.38	-0.60
<b>Correlation</b>	0.997	0.996

## Performance Characteristics

Within-run and total precision were analyzed and the following results were obtained:

### Within-run

	<u>Concentration Level</u>		
	Low	Mid	High
Mean, µg/mL	4.7	15.1	29.2
SD, µg/mL	0.26	0.42	0.62
CV, %	5.5	2.8	2.1
N	120	120	120

### Total Precision

	<u>Concentration Level</u>		
	Low	Mid	High
Mean, µg/mL	4.7	15.1	29.2
SD, µg/mL	0.30	0.48	0.72
CV, %	6.4	3.2	2.5
N	120	120	120

### Method Comparison:

A total of 126 serum samples having theophylline values throughout the assay range were tested with new CEDIA Theophylline Assay on the BM/Hitachi 911 and the BM/Hitachi 704, and yielded the following results:

Number of Observations	Slope	Intercept	Correlation Coefficient
126	1.09	-0.60	0.996

The performance information establishes the basis for substantial equivalence to the predicate device.