

K961487



Diagnostics

JUN 12 1996

510(k) Summary

-
- Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.
-
- 1. Submitter name, address, contact** Boehringer Mannheim Corporation
2400 Bisso Lane
P.O. Box 4117
Concord, CA 94524-4117
(510) 674 - 0690 extension 8415

Contact Person: Mary Koning

Date Prepared: April 9, 1996
-
- 2. Device name** Proprietary name: Elecsys® HCG Assay

Common name: Electrochemiluminescence assay for the determination of human chorionic gonadotropin (HCG).

Classification name: System, Test, Human Chorionic Gonadotropin
-
- 3. Predicate device** We claim substantial equivalence to the Enzymun® HCG Assay(K896901).
-
- 4. Device Description** The Elecsys® test principle is based on competition principle. Total duration of assay: 18 minutes (37° C).
•1st incubation (9 minutes): Sample (15 µL), a biotinylated HCG specific antibody (75 µL), and a specific anti-HCG antibody labeled with a ruthenium complex (75 µL).
•2nd incubation (9 minutes): After addition of streptavidin-coated microparticles (35 µL), the still-free binding sites of the labeled antibody become occupied, with formation of an antibody-hapten complex. The entire complex is bound to the solid phase via interaction of biotin and streptavidin.
-

Continued on next page



510(k) Summary, Continued

**4.
Device
Description**

•The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier (0.4 second read frame).
•Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent bar code.

**5.
Intended use**

Immunoassay for the in vitro quantitative determination of human chorionic gonadotropin in human serum and plasma.

**6.
Comparison
to predicate
device**

The Boehringer Mannheim Elecsys® HCG Assay is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed Enzymun® HCG Assay (K896901).

The following table compares the Elecsys® HCG Assay with the predicate device, Enzymun® HCG Assay. Specific data on the performance of the test have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device is provided in attachment 6.

Similarities:

- Intended Use: Immunoassay for the in vitro quantitative determination of human chorionic gonadotropin (HCG)
- Sample type: Serum and plasma
- Antibody: Polyclonal Sheep anti-HCG antibodies
- Solid phase binding principle: Streptavidin/Biotin
- Assay standardization: World Health Organization Standard (WHO) 75/537 (1st IRP).

Continued on next page

**510(k) Summary, Continued****6.
Comparison
to predicate
device cont.****Differences:**

Feature	Elecsys® HCG	Enzymun® HCG
Reaction test principle	Electrochemiluminescence	ELISA/1-step sandwich assay using streptavidin technology
Instrument required	Elecsys® 2010	ES 300
Calibration Stability	A calibration is recommended every 7 days if kits is not consumed; 4 weeks with same reagent lot if reagent is consumed within 7 days.	Calibration required every run

Performance Characteristics:

Feature	Elecsys® HCG			Enzymun® HCG		
Precision	Modified NCCLS (mIU/mL):			Modified NCCLS (µg/dL):		
Level	<u>Serum</u>	<u>Control 1</u>	<u>Control 2</u>	<u>Low</u>	<u>Mid</u>	<u>High</u>
N	60	60	60	58	58	53
Within-Run	24.80	35.39	854.30	19.9	181.5	593.8
%CV	4.5	3.3	2.7	4.6	2.9	3.1
Total	24.80	35.39	854.30	19.9	181.5	593.8
%CV	5.8	3.9	4.5	5.7	4.6	4.1
Lower Detection Limit	0.5 mIU/mL			1.5 mIU/mL		

Continued on next page



510(k) Summary, Continued

6.
Comparison
to predicate
device, (cont.)

Performance Characteristics:

Feature	Elecsys® HCG	Enzymun® HCG																								
Linearity	0.5 - 1,000 mIU/mL (with a deviation from a linear line of $\pm 10\%$)	1.5 - 600 mIU/mL (with a deviation from a linear line of $\pm 10\%$)																								
Method Comparison	<p>Vs Enzymun-Test® HCG</p> <p><u>Least Squares</u></p> <p>$y = 1.35x - 9.21$ $r = 0.989$ $SEE = 17.50$ $N = 64$</p> <p><u>Passing/Bablok</u></p> <p>$y = 1.29x - 4.05$ $r = 0.989$ $SEE = 17.50$ $N = 64$</p>	<p>Vs Enzymun-Test® HCG</p> <p><u>Least Squares</u></p> <p>$y = 1.047x - 4.92$ $r = 0.996$ $SEE = 13.056$ $N = 49$</p>																								
Interfering substances	No interference at:	No interference at:																								
Bilirubin	25 mg/dL	51.7 mg/dL																								
Hemoglobin	1 g/dL	1 g/dL																								
Lipemia	1500 mg/dL	1250 mg/dL																								
Biotin	30 ng/mL	200 ng/mL																								
Specificity	<table><tr><th></th><th>Level tested (mIU/mL)</th><th>% Cross-reactivity</th></tr><tr><td>LH</td><td>1000</td><td>0.07</td></tr><tr><td>FSH</td><td>1000</td><td>0.09</td></tr><tr><td>TSH</td><td>2500</td><td>0.000</td></tr></table>		Level tested (mIU/mL)	% Cross-reactivity	LH	1000	0.07	FSH	1000	0.09	TSH	2500	0.000	<table><tr><th></th><th>Level tested (mIU/mL)</th><th>% Cross-reactivity</th></tr><tr><td></td><td>400</td><td>0.15</td></tr><tr><td></td><td>400</td><td>0.28</td></tr><tr><td></td><td>2000</td><td>5.0</td></tr></table>		Level tested (mIU/mL)	% Cross-reactivity		400	0.15		400	0.28		2000	5.0
	Level tested (mIU/mL)	% Cross-reactivity																								
LH	1000	0.07																								
FSH	1000	0.09																								
TSH	2500	0.000																								
	Level tested (mIU/mL)	% Cross-reactivity																								
	400	0.15																								
	400	0.28																								
	2000	5.0																								