

**Micro Diagnostics Corporation**

K961803

JUL - 5 1996

June 21, 1996

PAGE A

Food and Drug Administration
Center for Devices and Radiological Health
Office of Device Evaluation
Document Mail Center (HFZ-401)
9200 Corporate Blvd
Rockville MD 20850

RE: K961803, Submission 510(k) Summary for Spuncrit™

Dear Sir or Madam:

A summary of the safety and effectiveness information upon which
substantial equivalence is based as follows:

CLASSIFICATION:

Class II Hematocrit Centrifuge

PROPRIETARY NAME:

Spuncrit™ (Model DRC-40)

**PRODUCTS UPON WHICH
SUBSTANTIAL EQUIVALENCE
IS BASED:**

o Compur M1100 Minicentrifuge
FDA No. K850391
Date: May 10, 1985

**SIMILARITIES TO PREDICATE DEVICES CURRENTLY ON THE US
MARKET****A) Predicate's Name and Manufacturer Address**Hematocrit

- o Compur M1100 Hematocrit System
Bayer Diagnostic + Electronic GmbH
Steinerstrasse 15
8000 Muchen 70
Germany
C/O Hans J. Brouwers and Associates
Box 823
Mt Prospect, IL 60056
FDA Reference No. K850391

PAGE B

A1) Comparison of Spuncrit™ to Compur (Predicate Device)

MICRO DIAGNOSTICS
SPUNCRIPT™COMPUR
M1100

- o Micro Hematocrit Centrifuge
- o Operations Manual
- o Battery Charger
- o 100 Micro Capillary/
Micro Capillary
Holders (heparinized)

- o Micro Hematocrit Centrifuge
- o Operation Manual
- o Battery Replacement
- o 100 Capillaries
(heparinized)

A2) The general features of the Spuncrit™ (Model DRC-40) as compared to the Compur M1100 Minicentrifuge.

SPUNCRIPT™COMPUR M1100

- o General Laboratory Device
- o Methodology:
Centrifugal
- o Measures Hematocrit
- o Battery Operated
- o Portable, Hand-Held
- Safety Features:
 - (1) Cover Lock
 - (2) Low Battery Charge Indicator
 - (3) Motor failure indicator
 - o Optically read
(digital readout)

- o General Laboratory Device
- o Methodology:
Centrifugal
- o Measure Hematocrit
- o Battery Operated
- o Portable, Hand-held
- Safety Features:
 - (1) Cover Lock
 - (2) Low Battery Charge Indicator
 - (3) No such feature
 - o Manual, visual read

PAGE C

- o Uses plastic disposable
- o Range 14% to 57%
- o Operates at 10,000 RPM
- o RCF= 4260
- o Rotor Radius= 1.5 inches

- o Uses glass disposable
- o Range 10% to 80%
- o Operates at 11,500 RPM
- o RCF=5396
- o Rotor Radius= 1.5 inches

A3) The safety features of the Spuncrit™ as compared to the Compur M1100 Minicentrifuge.

SPUNCRIT™

COVER LOCK

Mechanical lock switch which prevents cover from being opened while centrifuge is in operation

MECHANICAL BRAKE

In the event the cover is opened prior to the Spuncrit completing the 150 second cycle time the rotor automatically is braked, stopping its rotation within a few seconds

LOW BATTERY CHARGE

Digital display reads out "C" indicating to operator to charge batteries

INSTRUMENT MALFUNCTION

In the event of an electronic failure within the instrument the digital readout will not display "8's", as the DRC-40 is running.

COMPUR M1100

COVER LOCK

Mechanical lock switch which prevents cover from being opened while centrifuge is in operation

MECHANICAL BRAKE

No such feature

LOW BATTERY CHARGE

An illuminated red LED indicates batteries need to be replaced

INSTRUMENT MALFUNCTION

No such feature

PAGE D

- A4) Comparative data (Performance/Effectiveness Information) between the Spuncrit™ (Model DRC-40) and Compur's M1100 Minicentrifuge.

RESULTS OBSERVED WHEN OPERATING THE SPUNCRITCOMPARISON DATA

An "In Service" to explain the operation of the Spuncrit was performed in each of the locations. Approximately 10 minutes were utilized for training of persons operating the Spuncrit. The comparative study at each location included a total of 45 whole blood samples encompassing three levels ranging from below the normal range to the upper level of the normal range. The samples were split such that each location ran 15 samples at each level.


Spuncrit operating in a Laboratory:

- o Typical laboratory with technician running the Spuncrit

N = 45

Range = 14% Hct to 49% Hct

R = .99



Each data point represents an average of 16 whole blood samples.


Spuncrit operating in Physician Office - #1

- o Physician Office employs six Physicians with nurse running Spuncrit.

N = 45

Range = 14% Hct to 49% Hct

R = .99



Each data point represents an average of 16 whole blood samples.


Spuncrit operating in Physician Office - #2

- o Physician Office employs twelve physicians with a nurse running Spuncrit

N = 45

Range = 19% Hct to 53% Hct

R = .99



Each data point represents an average of 16 whole blood samples.

PAGE E


Spuncrit operating in
an Alternate Site:

- o Blood Bank with a
Phlebotomist running
the Spuncrit.

N = 45

Range = 14% Hct to 49% Hct

R = .99



Each data point represents an average of 15
whole blood samples.

PRECISION DATA

The precision study was performed within the physician office and another alternate site; a blood bank. Assayed whole blood controls were used in the precision study. Three levels of controls were used. A total of fifteen (15) samples were processed per each level for a total of 45 samples. The range varied from below the normal range to slightly above the normal range for the hematocrit. The study was repeated utilizing assayed whole blood controls from the same lot number seven (7) days after the first study.

WITHIN RUN PRECISION DATA

PHYSICIAN OFFICE # 1

Day One (Hematocrit)

Level 1	Level 2	Level 3	
18.1	30.8	46.3	X
.26	.94	1.03	SD
.01	.03	.02	Cv
15	15	15	N

Day Seven (Hematocrit)

Level 1	Level 2	Level 3
16.8	32.4	46.4
.77	.63	.63
.04	.02	.02
15	15	15

PHYSICIAN OFFICE # 2

Day One (Hematocrit)

Level 1	Level 2	Level 3	
17.9	33.5	46.1	X
.70	.52	.74	SD
.04	.02	.02	Cv
15	15	15	N

Day Seven (Hematocrit)

Level 1	Level 2	Level 3
17.4	32.7	47.7
.64	.59	.59
.04	.02	.01
15	15	15

ALTERNATE SITE (BLOOD BANK) # 3

Day One (Hematocrit)

Level 1	Level 2	Level 3	
18.6	33.7	47.3	X
.51	.88	1.11	SD
.03	.02	.02	Cv
15	15	15	N

Day Seven (Hematocrit)

Level 1	Level 2	Level 3
18.8	33.4	47.5
.68	.51	.52
.04	.02	.01
15	15	15

PAGE G

BETWEEN RUN PRECISION DATA

Day One and Day 7 (Hematocrit)

Physician Office #1	\bar{X}	SD	Cv	N
Level 1	17.4	.92	.05	2
Level 2	31.6	1.13	.04	2
Level 3	46.4	.07	<.01	2

Physician Office #2	\bar{X}	SD	Cv	N
Level 1	17.6	.35	.02	2
Level 2	33.1	.56	.02	2
Level 3	46.9	1.13	.02	2

Alternate Site Blood Bank # 3	\bar{X}	SD	Cv	N
Level 1	18.7	.14	.01	2
Level 2	33.6	.21	.01	2
Level 3	47.4	.14	<.01	2

Level 1,2,3 whole blood controls were Para 4 from Steck Laboratories Lot Numbers 61410065, 61410066, and 61410067 respectively. Hematocrit values(%) for Streck controls are 16 +3 or -3, 32 +4 or -4, and 47 +5 or -5 respectfully.

NOTE:

An additional precision study was performed in a typical laboratory with 10 whole blood samples run at each level, a total of five (5) levels were used. The range varied from the lower end of the Spuncrit to the highest end of the Spuncrit. this study was repeated utilizing the same whole blood samples 13 days after the first study.

Whole blood samples (5 levels, A through E) were run on Day 1 and Day 13. These samples were prepared from the same lot, run on two different dates (Day 1 and Day 13).

PAGE . H

WITHIN RUN PRECISION DATA

LABORATORY LOCATION

	<u>DAY ONE (Hematocrit)</u>				<u>DAY 13 (Hematocrit)</u>			
	\bar{X}	SD	Cv	N	\bar{X}	SD	Cv	N
LEVEL A	14.1	.32	.02	10	14.6	.52	.04	10
LEVEL B	29.3	.48	.02	10	30.3	.67	.02	10
LEVEL C	37.3	.48	.02	10	36.8	.92	.02	10
LEVEL D	53.6	.58	.01	10	52.5	1.51	.02	10
LEVEL E	56.3	.68	.01	10	56.6	.52	.01	10

BETWEEN RUN PRECISION DATA

LABORATORY LOCATION

	<u>DAY ONE and DAY THIRTEEN (Hematocrit)</u>			
	\bar{X}	SD	Cv	N
LEVEL A	14.4	.35	.02	2
LEVEL B	29.8	.71	.02	2
LEVEL C	37.1	.35	.02	2
LEVEL D	53.1	.78	.02	2
LEVEL E	56.4	.21	<.01	2

CONCLUSION: Based on the above information it is believed that the two products are substantially equivalent and the results provided by the Spuncrit provide the same clinical utility as the Compur M1100 Centrifuge for Hematocrit.

PAGE I

A6) Limitations of Procedure

- o Hematocrit
- o Estimated Hemoglobin

There are no known interfering substances that will interfere with the Hematocrit determination performed by the Spuncrit. However, severely hemolyzed blood should not be used in determining the Hematocrit by the Spuncrit.

The range of the Spuncrit for Hematocrit is 14% to 57%. The estimated Hemoglobin range is from 12.2 g/dL to 18.2 g/dL. The estimated Hemoglobin only estimates Hemoglobin as a function of the PCV. The estimated Hemoglobin is not representative of anemic problems. As a result the Spuncrit will only estimate the Hemoglobin for normal range only. Genetic Anemias or anemia caused by chronic disease states will not be estimated.