

GAYMAR®

GAYMAR INDUSTRIES, INC.

APR 30 1997

K962788

Premarket Notification [510(K)] Summary

Submitter: Peter Scott

Phone: 716/662-2551, ext. 676

FAX: 716/662-8763

**Date Summary
Prepared:** January 30, 1997

Device Name: CLINI-DYNE Rotational Therapy System

Common Name: Rotation/Low Air Loss System

Classification Name: Alternating Pressure Air Flotation Mattress per
21 CFR, Section 880.5550

Predicate Device: CLA 1400 System

Intended use of Device:

The CLINI-DYNE rotational/low air loss system is a patient rotation, pressure relieving support system used for the treatment and prevention of decubitus ulcers, urinary tract blockage, to aid circulation and for complications associated with immobility.

Description: The CLINI-DYNE Rotational Therapy System consists of an electromechanical inflation pump control unit which is connected to a mattress with cells (air filled bladders) through hose assemblies connected to the inflation pump. The system affords programmed patient position changes by gently turning the patient from supine to left side or right side and back to supine position when staff is unable to do so or needs assistance.

Substantial Equivalence:

The following tables, CLINI-DYNE substantial equivalence matrix and the safety testing comparison table, summarize the technological characteristics and the nonclinical performance data upon which the substantial equivalence submission was made to the Food and Drug Administration.

CLINI-DYNE SUBSTANTIAL EQUIVALENCE MATRIX

7/11/96

CLA-1400 (Predicate)	CLD-1000
--------------------------------	-----------------

PERFORMANCE		
Patient Support Surface	Overlay	Mattress
Air Source	Blower	Electro-mechanical reciprocating pump
Pressure control mechanism	Microprocessor controlled pressure relief valve	Microprocessor controlled pressure relief valve
Voltage	120v 60Hz Nominal	120v 60Hz Nominal
Ambient	60-90°F	60-90°F
Tissue interface pressure measurements	Pressure relief	Pressure relief
Rotation angle	0-45°	0-40° **
Safety		
Maximum allowable pressure controlled by	Mechanical relief valve (55-60 mmHg)	Mechanical relief valve (56 mmHg)
CPR Deflate mechanism	Quick-disconnect feature of hose from pump	Remove panel at mattress
Physical characteristics		
Construction (Patient support surface)	Air filled cells used as an overlay	Air filled cells contained in a foam crib used as a mattress
Function & Intended uses		
Type of therapy provided	Treatment of the diseases of immobility and treatment of pressure ulcers	Treatment of the diseases of immobility and treatment of pressure ulcers
(Primary) Target population	Acute care and alternate care settings	Acute care and alternate care settings
Life Expectancy (pump)	Reusable*	Reusable*
Life Expectancy (cell)	Reusable*	Reusable*
Product Labeling	Refer to TAB C - LABELING	Refer to TAB B - LABELING
Approvals	UL 544	UL 2601-1 (submitted)

* Cleanable using suggested disinfectants and cleaning agents.

** Dependent on physical characteristics of the patient. (40° minimum when set at maximum angle setting for each patient.)

COMPARISON OF MAJOR CHARACTERISTICS

Characteristic	CLA-1400 (PREDICATE)	CLD-1000
<u>Alerts</u>		
Audible Alerts	yes	yes
Alert Silence (minutes)	30	no
Head of Bed Elevation Alert	yes	no*
Out of Rotation Alert	yes	no*
Siderail UP/Down Alert	yes	no*
Visual Alerts	yes	yes

<u>Rotation</u>		
Automatic Left Only	yes	yes
Automatic Left to Right	yes	yes
Automatic Right Only	yes	yes
Center Only (no rotation)	yes	yes
Manual Left Only	yes	yes
Manual Right Only	yes	yes
Var.Angle, Min.range (degrees)	0-45	0-45
Variable Rate	yes	yes
Rotation Cycle Times (minutes)	90, 30, 10, 2	60, 45, 30, 15

<u>Cushion</u>		
Placed onto Hospital Bed	yes	no (mattress replacement)
CPR Deflate	yes	yes
Low Air Loss	yes	yes
Pressure Relieving	yes	yes
Support Surface Inflation Control	yes	yes
Disposable Covers	yes	no**
Disposable Support Surface	yes	no**
Footboard/Foot Support Cushion	yes	yes

<u>Air Control Pump</u>		
120 Volt AC	yes	yes
60 Hz	yes	yes
Amperage	1-5	1
Control Panel Lock-out	yes	yes
Microprocessor Control	yes	yes
U.L. Listed	yes	yes
Weight (lbs.)	25	19

* See Substantial Equivalence

**Reusable mattress

5/6/96

VI. SAFETY SUMMARY

Safety Testing Comparison Table

	CLA 1400	CLD 1000
Biocompatibility		
a) primary skin irritation	pass	pass
b) delayed contact sensitization	pass	pass
c) cytotoxicity	pass	pass
 Tissue Interface Testing	 pass	 pass
 Flammability	 pass	 pass
 UL Approvals	 UL 544	 UL 2601-1 (pending)
 Alerts		
a) hi/low pressure	yes	yes
b) out of range	yes	yes