

K 965062

JUN - 4 1997

510 (K) Notification  
Siemens SC9000/ SC9015 MULTIGas<sup>™</sup> and MULTIGas+<sup>™</sup> Modules

**510(k) SUMMARY**  
as required per 807.92(c)

2: Submitters Name, Address:

Siemens Medical Systems, Inc.  
16 Electronics Avenue  
Danvers, MA 01923  
Tel: (508) 750-7500  
Fax: (508) 777-3398  
Official Correspondent: David Simard, Director Quality Assurance and Regulatory Affairs  
Contact person for this submission: Jacqueline E. M. Emery  
Date submission was prepared: December 9, 1996

3: Trade Name, Common Name and Classification Name:

A. Trade Name: Siemens SC9000/ SC9015 Series MULTIGas<sup>™</sup> Module  
Siemens SC9000/SC9015 Series MULTIGas+<sup>™</sup> Module

B. Common Name, Classification Number, Class and Regulation Number:

| Common Name                                    | Classification Number | Class | Regulation Number |
|--|-----------------------|-------|-------------------|
| Cardiac monitor                                | 74DRT                 | II    | 21 CFR 870.2300   |
| Arrhythmia detector<br>& Alarm System          | 74DSI                 | III   | 21 CFR 870.1025   |
| Breathing frequency<br>monitor                 | 73BZQ                 | II    | 21 CFR 868.2375   |
| Pulse rate monitor                             | 74BWS                 | II    | 21 CFR 870.2300   |
| Non-indwelling blood<br>pressure monitor       | 74DXN                 | II    | 21 CFR 870.1130   |
| Clinical electronic thermometer                | 80BWX                 | II    | 21 CFR 880.2910   |
| Pulse Oximeter                                 | 74DQA                 | II    | 21 CFR 870.2700   |
| Cardiac Output Monitor                         | 74KFN                 | II    | 21 CFR 870.1435   |
| end-tidal Carbon-Dioxide Monitor               | 73CCK                 | II    | 21 CFR 868.1400   |
| Indwelling Blood Pressure Monitor              | 74CAA                 | II    | 21 CFR 870.1110   |
| Indwelling blood<br>pressure monitor           | 74CAA                 | II    | 21 CFR 870.1110   |
| Heart Rate Monitor, Neonatal                   | 74FLO                 | II    | 21 CFR 870.2300   |
| Ventilatory Effort Monitor<br>(Apnea Detector) | 73FLS                 | II    | 21 CFR 868.2375   |
| Monitor Blood Pressure, Neonatal,<br>Invasive  | 74FLP                 | II    | 21 CFR 870.1110   |

**Siemens Medical Systems, Inc.**

Electromedical Group

16 Electronics Avenue  
Danvers, MA 01923  
USA

Tel: (508) 750-7500  
Fax: (508) 777-3398  
Telex: 511958 (Siemensm SD)

## 510 (K) Notification

### Siemens SC9000/ SC9015 MULTIGas™ and MULTIGas+™ Modules

|  |       |    |               |
|--|-------|----|---------------|
| Multi Gas Monitor                      |       |    |               |
| Analyzer, Gas, Carbon Dioxide, Gaseous | 73CCK | II | 21CFR868.1400 |
| Analyzer, Gas, Enflurane Gaseous       | 73CBQ | II | 21CFR868.1500 |
| Analyzer, Gas, Halothane, Gaseous      | 73CBS | II | 21CFR868.1620 |
| Analyzer, Gas, Nitrous Oxide, Gaseous  | 73CBR | II | 21CFR868.1700 |
| Analyzer, Gas, Oxygen, Gaseous         | 73CCL | II | 21CFR868.1720 |

#### 4: Predicate Device Identification:

The MULTIGas modules use the same hardware, provided by the same OEM manufacturer, as the Hewlett-Packard model M1026A Anesthesia Gas Module Gas Module, cleared under 510(K) number K951127.

The SC9000/9015 Bedside Monitor System (510 (K) number K946306) provides the display and user interface capabilities for the MULTIGas™ modules.

#### 5. Device Description

There are two multi gas modules being submitted as part of this 510(K) Notification: the MULTIGas™ (MGM) and the MULTIGas+™ (MGM+). Both modules are free standing units that perform sidestream measurements of respiratory gases (CO<sub>2</sub>, N<sub>2</sub>O, and O<sub>2</sub>) and anesthetic gases. Both modules automatically identify and report measurement data to the SC9000/SC9015 Bedside Monitor for display.

The MULTIGas™ and MULTIGas+™ modules differ only in the way they measure O<sub>2</sub>. The basic MULTIGas Module measures O<sub>2</sub> using a galvanic cell, and calculates average inspiratory values for O<sub>2</sub> (labeled iO<sub>2</sub>). The MULTIGas+ incorporates a faster-responding paramagnetic sensor that provides both inspired and expired O<sub>2</sub> measurements (iO<sub>2</sub> and etO<sub>2</sub>).

The enhanced software (version VB1.1) is compatible with previously sold versions of the monitor. A retrofit will be offered to the owners of units with previous software versions.

#### 6. Intended Use:

The intended use of the Siemens SC9000/ SC9015 MULTIGas™ Module and MULTIGas™+ Module is to measure inspiratory and expiratory carbon dioxide, inspiratory and expiratory oxygen, inspiratory and expiratory Nitrous Oxide and anesthetic agents. This device will produce visual and audible alarms if any of these parameters vary beyond preset limits and produce recordings. This device will connect to either the Siemens SIRENET or Infinity (Olympus) network.

---

### Siemens Medical Systems, Inc.

Electromedical Group

16 Electronics Avenue  
Danvers, MA 01923  
USA

Tel: (508) 750-7500  
Fax: (508) 777-3398  
Telex: 511958 (Siemensm SD)

510 (K) Notification  
Siemens SC9000/ SC9015 MULTIGas<sup>™</sup> and MULTIGas+<sup>™</sup> Modules

---

7. Table of device similarities and differences to predicate device

|                            | Substantial Equivalence<br>Hewlett Packard<br>Model M1026A Anesthesia Gas<br>Module | SC 9000/9015 Series MULTIGAS <sup>™</sup><br>and MULTIGAS+ <sup>™</sup> Modules |
|----------------------------|---|---|
| Manufacturer               | Hewlett-Packard   | Siemens Medical Systems -<br>Electromedical Group                               |
| 510K Number                | K951127   | New   |
| Intended Population        | Not stated in literature  | Adult and Pediatric   |
| Module/Stand-Alone Monitor | Module: Communicates only<br>with HP Anesthesia Component<br>Monitoring System      | Module: Communicates only<br>with Siemens SC9000/9015<br>Bedside Monitors       |
| Displayed Parameters       | CO <sub>2</sub> , N <sub>2</sub> O, Respiration Rate,<br>Anesthetic Agents          | Same  |
| Principle of Operation     | Non-dispersive Infra Red  | Same  |
| Measuring Methods          | sidestream  | Same  |
| Waveform Display           | Inspired and expired<br>concentrations (in %) plus<br>waveform and trend data       | Same  |
| Dimensions HxWxD (mm/in)   | 80x370x439 / 3.5x14.5x17.3  | 146x183x451 / 5.7x7.2/17.8  |
| Weight kg /lb              | 8.2 / 18  | 7.3 / 16.0  |

8. Assesment of non-clinical performance data for equivalence: Both the Siemens MULTIGas<sup>™</sup> Modules and the HP M1026A Anesthesia Gas module employ, internally, the Andros Inc Model 4700 MGM Multi Gas Module as the means for measuring the various gases.

9. Assesment of clinical performance data for equivalence: Not applicable

10. Biocompatibility: Not applicable

11. Sterilization: Not applicable

---

**Siemens Medical Systems, Inc.**

Electromedical Group

16 Electronics Avenue  
Danvers, MA 01923  
USA

Tel: (508) 750-7500  
Fax: (508) 777-3398  
Telex: 511958 (Siemens SD)



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUN - 4 1997

Ms. Jacqueline E. M. Emery  
Siemens Medical Systems, Inc.  
16 Electronics Avenue  
Danvers, Massachusetts 01923

Re: K965062  
Siemens SC9000/SC9015 MULTIGas™ and MULTIGas+™ Modules  
Regulatory Class: III (three)  
Product Code: 74 DSI  
Dated: March 25, 1997  
Received: March 27, 1997

Dear Ms. Emery:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Jacqueline E. M. Emery

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Thomas J. Callahan, Ph.D.  
Director  
Division of Cardiovascular,  
Respiratory, and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indicated Use Statement**

The Siemens SC9000/SC9015 MULTIGas<sup>™</sup> and MULTIGas+<sup>™</sup> Modules are indicated for adult and pediatric patient populations in an environment where patient care is provided by healthcare professionals (Physician, Nurse, Technician) when the professional determines that these devices are required to measure any one or more of the following parameters:

- Respiration rate
- inspired and expired Carbon Dioxide (CO<sub>2</sub>)
- inspired and expired Oxygen (MULTIGas+ only)
- average inspired Oxygen (MULTIGas only)
- inspired and expired gas concentrations of Enflurane, Halothane, Isoflurane, Desflurane, Sevoflurane, and Nitrous Oxide

**MRI Compatibility Statement:**

The Siemens SC9000/9015 MULTIGas<sup>™</sup> and MULTIGas+<sup>™</sup> modules are not intended for use in an MRI magnetic field

Prescription Use

per 21 CFR 801.109 ✓

or

  
PCRND Sign-Off

(Division Sign-Off)  
Division of Cardiovascular, Respiratory,  
and Neurological Devices

510(k) Number over-the-counter use

---

**Company Confidential**

**Siemens Medical Systems, Inc.**

Electromedical Group

16 Electronics Avenue  
Danvers, MA 01923  
USA

Tel: (508) 750-7500  
Fax: (508) 777-3398  
Telex: 511958 (Siemensm SD)