

FEB 12 1997

K970365

**510(k) SUMMARY FOR PASSY-MUIR OXYGEN (O₂) ADAPTOR FOR USE WITH THE
PASSY-MUIR LOW-PROFILE TRACHEOSTOMY AND VENTILATOR SPEAKING VALVES
(PMV 2000 AND 2001)**

Contact: Patricia E. Passy Phone: (714) 833-8255
Passy-Muir, Inc. Facsimile: (714) 833-8299
4521 Campus Drive, Suite 273
Irvine, CA 92715

Establishment Registration Number: 2024841

Date Summary Prepared: 1/29/97

Name of the Device: Passy-Muir Oxygen (O₂) Adapter
Product Number PMA 2000

Common or Usual Name: Oxygen (O₂) Adapter; Accessory O₂ Adapter; Supplemental
O₂ Adapter; Speaking Valve O₂ Adapter; Airway
Connector; O₂ Port; Supplemental O₂ Port

Classification Name: TBD: Accessory O₂ Adapter for the Tracheostomy Speaking Valves

Possible Product Classification Codes:

1. Airway Connector (Extension): Classification Number: 73BZA

21 CFR 868.5810 - "An airway connector is a device
intended to connect a breathing gas source to a
tracheal tube, Tracheostomy tube or mask."
Formerly a Class II device.

2. Oxygen Mask: Classification Number: 73BYG

21 CFR 868.5680 - "An oxygen mask is a device placed over
a patient's nose, mouth, or Tracheostomy to
administer oxygen or aerosols." Formerly a Class II
device.

Predicate Devices:

1. Trach Ring Adapter (Marquest Medical): Product # 178548, 15mm
I.D. x 22 mm OD with side port for supplemental oxygen supply to
trach or ET tube

2. Aerosol Tee Connector (Hudson RCI): Product # 1648, 22mm OD
horizontal ports with a 15mm ID vertical port

3. Tracheostomy Mask (Hudson RCI): Product # 1075, Adult mask
for tracheostomy and laryngectomy aerosol therapy

Description of Device:

The Passy-Muir O₂ Adapter is a one-piece, injection-molded plastic
device designed to attach to the outer surface of either one of the Passy-
Muir Low-Profile Tracheostomy and Ventilator Speaking Valves (PMV
2000 or PMV 2001) and provide a means for connecting a supplemental
oxygen supply tube.

Intended Use:

The P-M O₂ Adapter is designed for use with the companion PMV 2000 (clear) and PMV 2001 (purple) Low Profile Tracheostomy and Ventilator Speaking Valves. The P-M O₂ Adapter allows for improved mobility of those patients requiring a tracheostomy tube, speaking valve, and low flow supplemental oxygen.

Principles of Operation:

The P-M O₂ Adapter is designed to snap around the body of the low-profile speaking valve, held in place by friction, the mechanical strength of two gripper flanges, and the shape of the device (which hugs the outer profile of the valve). The O₂ Adapter is equipped with an integral, tapered tube fitting that provides a friction fit with the inside diameter of flexible, supplemental oxygen tubing (not supplied by Passy-Muir). The low-pressure oxygen (flowing through the flexible oxygen supply tubing) passes first through the tube fitting of the O₂ Adapter, then into a duct-like channel created by the juxtaposed walls of the Adapter and the outer surface of the Speaking Valve, exiting *at the front edge* of the valve, where it is either inhaled (along with ambient air) or allowed to dissipate into the atmosphere. This is an inherently *open* system, functioning in the manner of a Tracheostomy mask or an aerosol "T" piece, commonly used in weaning tracheostomized patients of off ventilators.

Summary of the Basis for the Finding of Substantial Equivalence:

The safety and effectiveness of supplemental oxygen devices is based on a long history of O₂ delivery systems.

Supplemental O₂ delivery systems have been used with tracheostomized and ventilator dependent patients since before the passage of the Medical Device Amendments of 1976. The Passy-Muir Oxygen (O₂) Adapter and the predicate devices have the same general intended use as a delivery system for low levels of supplemental oxygen. It is similar to the predicate devices as it attaches directly to (or in immediate proximity of) the tracheostomy speaking valve and delivers O₂ *in front* of speaking valve diaphragm, usually in an *open system* configuration where O₂ delivery is continuous.