

## 510(k) Summary

K981799  
July 22, 1998

**Submitter's Name/Address**

Abbott Laboratories  
1920 Hurd Drive  
Irving, Texas 75038

**Contact Person**

Mark Littlefield  
Section Manager MS 1-8  
ADD Regulatory Affairs  
(972) 518-6062  
Fax (972) 753-3367

**Date of Preparation of this Summary:**

May 20, 1998

**Device Trade or Proprietary Name:**

Crea

**Device Common/Usual Name or Classification Name:**

Creatinine

**Classification Number/Class:**

75C6X/Class II

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: \_\_\_\_\_

**Test Description:**

Creatinine is an *in vitro* diagnostic assay for the quantitative determination of creatinine in human serum, plasma, or urine. The Creatinine assay is a clinical chemistry assay in which the creatinine in the sample, at an alkaline pH, reacts with picrate to form a creatinine-picrate complex. The rate of increase in absorbance at 500 nm due to the formation of this complex is directly proportional to the amount of creatinine in the sample.

**Substantial Equivalence:**

The Creatinine assay is substantially equivalent to the Boehringer Mannheim® Creatinine assay on the Hitachi® 717 Analyzer (K812095) for the serum, plasma, or urine applications.

32

0000011

These assays yield similar Performance Characteristics.

**Similarities to Boehringer Mannheim:**

- Both assays are *in vitro* clinical chemistry methods.
- Both assays can be used for the quantitative determination of creatinine.
- Both assays yield similar clinical results.

**Differences to Boehringer Mannheim:**

- There is a minor difference between the assay range.

**Intended Use:**

The Creatinine assay is used for the quantitation of creatinine in human serum, plasma, or urine.

**Performance Characteristics:**

Comparative performance studies were conducted using the AEROSSET™ System.

The Creatinine assay method comparison yielded acceptable correlation with the Boehringer Mannheim Creatinine assay on the Hitachi 717 Analyzer for the serum and urine applications. For the serum application, the correlation coefficient = 0.9997, slope = 0.873, and Y-intercept = 0.078 mg/dL. For the urine applications, the correlation coefficient = 0.9947, slope = 1.056, and Y-intercept = -2.410 mg/dL.

Precision studies were conducted using the Creatinine assay. Within-run, between-run, and between-day studies were performed using two levels of control material. For the serum application, the total %CV for Level 1/Panel 101 is 3.0% and Level 2/Panel 102 is 2.6%. For the urine application, the total %CV for Level 1/Panel 201 is 1.6% and Level 2/Panel 202 is 2.1%. The Creatinine assay is linear up to 38.7 mg/dL for the serum application, and 757.0 mg/dL for the urine application. The limit of quantitation (sensitivity) of the Creatinine assay is 0.10 mg/dL. These data demonstrate that the performance of the Creatinine assay is

83

substantially equivalent to the performance of the Boehringer Mannheim Creatinine assay on the Hitachi 717 Analyzer for the serum and urine applications.

**Conclusion:**

The Creatinine assay is substantially equivalent to the Boehringer Mannheim Creatinine assay on the Hitachi 717 Analyzer for the serum and urine applications as demonstrated by results obtained in the studies

37



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUL 22 1998

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Mark Littlefield  
Section Manager, Regulatory Affairs  
Abbott Laboratories  
1920 Hurd Drive  
Irving, Texas 75038

Re: K981799  
Crea  
Regulatory Class: II  
Product Code: CGX  
Dated: May 20, 1998  
Received: May 21, 1998

Dear Mr. Littlefield:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing-Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health


Enclosure

510(k) Number (if known): \_\_\_\_\_

Device Name: Creatinine

Indications For Use:

The Creatinine assay is used for the quantitation of creatinine in human serum, plasma, or urine. Creatinine measurements are used in the diagnosis and treatment of renal diseases and in monitoring renal dialysis.

  
\_\_\_\_\_  
(Signature Sign-Off)  
\_\_\_\_\_  
Clinical Laboratory Devices  
(510(k) Number) K981799

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)  
Prescription Use ☒ OR Over-The-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109)

(Optional Format 1-2-96)