

**510(k) PREMARKET NOTIFICATION  
ACE® LDL-C Reagent****SUMMARY OF SAFETY AND EFFECTIVENESS**

**In lieu of a 510(k) statement under 513(i) of the Act, this Summary of Safety and Effectiveness is provided as a 510(k) summary for disclosure to any other persons/companies without specific written authorization from Schiapparelli Biosystems, Inc.**

**Submitter**

Schiapparelli Biosystems, Inc.  
368 Passaic Avenue  
Fairfield, NJ 07004  
Phone: (973) 882-8630

**Contact Person**

Steven Dalessio  
Manager, Quality Assurance/Regulatory Affairs  
Phone: (973) 882-8630

**Device Names**

**Proprietary Name:** ACE® LDL-C Reagent  
LDL-C Calibrator  
LDL-C Controls

**Common Name:** Homogeneous assay for low density lipoprotein cholesterol

**Classification Name:** Low density lipoprotein cholesterol test  
Calibrator, Primary  
Low density lipoprotein cholesterol control

**Predicate Device:** Genzyme N-geneous® LDL Cholesterol Reagent  
[510(k) Number K971573]  
Genzyme N-geneous® LDL Cholesterol Calibrator  
Genzyme LDL Control Set

**Device Description**

The ACE® LDL-C Reagent contains two reagents. An aliquot of serum is added to the first reagent, which contains a unique detergent that selectively solubilizes the non LDL lipoproteins. Enzymes also present in the first reagent consume the cholesterol in a non color forming reaction. The second reagent contains another detergent that releases the remaining LDL lipoproteins. The enzyme reaction with LDL cholesterol, in the presence of a chromogenic coupler, produces color that is directly proportional to the amount of LDL cholesterol in the sample.

**Intended Use of the Device**

ACE® LDL-C Reagent is intended for use in the quantitative determination of low density lipoprotein cholesterol in human serum.

**510(k) PREMARKET NOTIFICATION**  
**ACE® LDL-C Reagent**

**SUMMARY OF SAFETY AND EFFECTIVENESS**

**COMPARATIVE FEATURES OF PREDICATE AND PROPOSED DEVICES**

<b>PARAMETER</b>	<b>PREDICATE DEVICE</b>	<b>PROPOSED DEVICE</b>
<b>Trade Name</b>	Genzyme N-geneous® LDL Cholesterol Reagent	ACE® LDL-C Reagent
<b>Reference No.</b>	K971573	TBD
<b>Analyte</b>	LDL cholesterol	LDL cholesterol
<b>Intended Use</b>	Quantitative determination of LDL cholesterol	Quantitative determination of LDL cholesterol
<b>Methodology</b>	Homogeneous, Direct	Homogeneous, Direct
<b>Reagents</b>		
<b>Reagent 1 Volume</b>	Liquid; Detergent, Enzymes 300 µL	Liquid; Detergent, Enzymes 300 µL
<b>Reagent 2 Volume</b>	Liquid; Detergent, Chromogenic coupler 100 µL	Liquid; Detergent, Chromogenic coupler 100 µL
<b>Specimen Type</b>		
<b>Volume</b>	Serum and plasma 3 µL	Serum 3 µL
<b>Assay System</b>		
<b>Reagent 1 + Sample</b>	Incubate 300 sec and Read	Incubate 300 sec and Read
<b>Reagent 2</b>	Read at 300 sec	Read at 300 sec
<b>Temperature</b>	37 °C	37 °C
<b>Detection Method</b>		
<b>Type</b>	Spectrophotometric	Spectrophotometric
<b>Wavelength, nm</b>	Bichromatic: 546/660	Bichromatic: 544/692

**510(k) PREMARKET NOTIFICATION**  
**ACE® LDL-C Reagent**

**SUMMARY OF SAFETY AND EFFECTIVENESS**

**PERFORMANCE ASSESSMENT**

Non-clinical test results submitted in the premarket notification include within-run and between-run precision and method correlation. Following is a data summary.

PARAMETER	PREDICATE DEVICE	PROPOSED DEVICE
<b>Performance Summary</b> <b>Assay Range</b> <b>Precision</b> <b>Within Run</b> <b>Between Run</b>	6.6 mg/dL to 992 mg/dL  ≤ 0.73 %CV ≤ 2.27 %CV	3 mg/dL to 800 mg/dL  ≤ 2.6 %CV ≤ 3.2 %CV
<b>Correlation vs</b> <b>Slope</b> <b>Intercept</b> <b>r</b> <b>n</b>	<b>Reference method (Ultracentrifugation)</b> 0.95 +3.02 0.96 54	<b>Hitachi 911</b> 1.111 -15.5 0.9747 70
<b>Correlation vs</b> <b>Slope</b> <b>Intercept</b> <b>r</b> <b>n</b>	<b>Immunoseparation method</b> 0.94 +4.46 0.97 92	<b>Friedewald calculation</b> 1.09 -17.8 0.9728 66

Based on these data, the Schiapparelli Biosystems ACE® LDL-C Reagent is substantially equivalent to the predicate device Genzyme N-geneous® LDL Cholesterol Reagent. On this basis, the reagent is determined to be safe and effective for its intended use. Performance details are included in the reagent product labeling.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

JUL 13 1999

Mr. Steven Dalessio  
Manager, Quality Assurance/Regulatory Affairs  
Schiapparelli Biosystems, Inc.  
368 Passaic Avenue  
Fairfield, New Jersey 07004

Re: K991733  
Trade Name: ACE® LDL-C Reagent  
Regulatory Class: I reserved (Product codes: MRR, JJX)  
II (Product code: JIS)  
Dated: May 20, 1999  
Received: May 21, 1999

Dear Mr. Dalessio:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.

A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification"(21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, flowing style.

Steven I. Gutman, M.D, M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K 991733

Device Name: ACE® LDL-C Reagent

Indications For Use:

ACE® LDL-C Reagent is intended for the quantitative determination of LDL cholesterol in serum using the ACE® clinical chemistry system.

LDL-C Calibrator is intended for the calibration of the ACE® LDL-C Assay.

LDL-C Controls are intended to monitor the performance of the ACE® LDL-C Assay.

The measurement of LDL cholesterol is a factor in the pathogenesis of atherosclerosis and coronary artery disease.

Jean Cooper  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K 991733

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)