

SEP 8 1999

K992472

**FDA Premarket Notification [510(k)] Summary**

Preparation Date: 21 July 1999

Contact: Gerald S. Palecki, Director, Quality and Regulatory

Device Name: Proprietary: *SureGuide* CO<sub>2</sub> Laser Beam Delivery SystemCommon: CO<sub>2</sub> Laser WaveguidePredicate Devices: CO<sub>2</sub> Laser Waveguide, reference K943543 & K960475.**Device Description:**

The *SureGuide* CO<sub>2</sub> Laser Beam Delivery System is an accessory for CO<sub>2</sub> laser systems. Its two primary components are a laser system interface adapter and a fiber cable assembly consisting of a hollow silica fiber having an internal coating that reflects and propagates CO<sub>2</sub> laser energy and a flexible protective outer cover.

The *SureGuide* CO<sub>2</sub> Fiber Cable has FSMA 905 series fiber optic connectors on each end and may be used with CO<sub>2</sub> laser systems that provide such connectors between the laser system and various beam delivery accessories.

**Intended Use:**

The *SureGuide* is intended for use with CO<sub>2</sub> laser systems for general and plastic surgery procedures, neurosurgery, ophthalmology, oral surgery, oto-rhino-laryngology, podiatry, gynecology, and urology procedures for incision, excision, vaporization, ablation, coagulation, and cauterization of soft tissue.

Refer to the laser system Directions for Use manual for specific indications for Use.

**Technological Characteristics Compared to Predicate Devices:**

The *SureGuide* CO<sub>2</sub> Laser Beam Delivery System has the same technological characteristics and materials as the Medical Optics, Inc. flexible fiber cable assembly and it is functionally equivalent to the Luxar fiber assembly. Clinicon Corporation purchased the products of Medical Optics, Inc. and holds a license from Rutgers University to manufacture the *SureGuide* flexible cable assembly.

The laser adapter design is derived from the Medical Optics laser system flexible cable assembly interface adapter and has a section designed to fit the Luxar LX-20 mast interface. Beam alignment and focus optics are built into the interface adapter to allow correction of inherent beam alignment variations between laser systems and provides optimization of laser beam coupling into the *SureGuide*.

**Non-Clinical Tests:**

The *SureGuide* performance characteristics have been evaluated through testing and analysis of laser power loss and beam quality from the guide when the laser is energized, compared to similar devices cleared for marketing in the past. The performance of the *SureGuide* is comparable.

**Conclusions Drawn from Tests and Analysis:**

The predicted energy transmission levels and beam quality meet criteria established through evaluation of the *SureGuide* on various commercially available medical CO<sub>2</sub> laser systems.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 8 1999

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Gerald S. Palecki  
Director, Quality and Regulatory  
Clinicon Inc.  
5825 Avenida Encinas  
Carlsbad, California 92008

Re: K992472  
Trade Name: SUREGUIDE CO<sub>2</sub> Laser Beam Delivery System  
Regulatory Class: II  
Product Code: GEX  
Dated: July 21, 1999  
Received: July 26, 1999

Dear Mr. Palecki:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 – Mr. Gerald S. Palecki

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



for Celia M. Witten, Ph.D., M.D.

Director  
Division of General and  
Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**SureGuide** CO<sub>2</sub> Laser Beam Delivery SystemPage 1 of 1510(k) Number: K 992472Device Name: SureGuide CO<sub>2</sub> Laser Beam Delivery System

## Indications for Use:

The *SureGuide* CO<sub>2</sub> Laser Beam Delivery System is indicated for use in general and plastic surgery, neurosurgery, ophthalmology, oral surgery, otorhino-laryngology, podiatry, gynecology, and urology procedures for incision, excision, vaporization, ablation, coagulation, and cauterization of soft tissue.


The specific *SureGuide* CO<sub>2</sub> Laser Beam Delivery System indications are dependent upon the cleared indications for use of the laser system and laser system accessories to which it is attached.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_

\_\_\_\_\_

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of General Restorative Devices

510(k) Number

K992472

(Optional Format 3-10-98)

Prescription Use X  
(Per 21 CFR 801.109)

2-2