

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY TEMPLATE**

A. 510(k) Number:

k123518

B. Purpose for Submission:

New device

C. Measurand:

Carbamazepine

D. Type of Test:

Quantitative immunoassay

E. Applicant:

Microgenics Corporation

F. Proprietary and Established Names:

Abbott Carbamazepine Assay

G. Regulatory Information:

1. Regulation section:

21 CFR 862.3645, Neuroleptic drugs radioreceptor assay test system

2. Classification:

Class II

3. Product code:

KLT, Enzyme Immunoassay, Carbamazepine

4. Panel:

Toxicology (91)

H. Intended Use:

1. Intended use(s):

See indication(s) for use below

2. Indication(s) for use:

The Abbott Carbamazepine Assay is used for the in vitro quantitative measurement of carbamazepine in human serum or plasma on the ARCHITECT cSystems. The measurements obtained are used in monitoring levels of carbamazepine to help ensure appropriate therapy.

3. Special conditions for use statement(s):

For prescription use only

4. Special instrument requirements:

ARCHITECT c8000 System

I. Device Description:

The Carbamazepine Assay kit is supplied ready-to-use in liquid form, for storage at 2 to 8°C. Each Carbamazepine Assay kit contains three bottles of Antibody Reagent (R1; 3 x 27mL), three bottles of Microparticle Reagent (R2; 3 x 9mL), and the package insert. Each kit is sufficient for 300 tests. The R1 antibody reagent contains <1.0% anti-carbamazepine monoclonal antibody (mouse) in Bis-Tris buffer and <0.09% NaN₃ as preservative. The R2 microparticle reagent contains <1.0% carbamazepine-coated microparticles in Tris buffer and <0.09% NaN₃ as preservative.

J. Substantial Equivalence Information:

1. Predicate device name(s):

Abbott Aeroset® Carbamazepine Assay

2. Predicate 510(k) number(s):

K993028

3. Comparison with predicate:

Item	Candidate device: Abbott Carbamazepine Assay (k123518)	Predicate device: Abbott Aeroset® Carbamazepine Assay (k993028)
Intended Use	The Abbott carbamazepine assay is used for the in vitro quantitative measurement of carbamazepine in human serum or plasma on the ARCHITECT cSystems. The measurements obtained are used in monitoring levels of carbamazepine to help ensure appropriate therapy.	Same

Methodology	Homogenous particle-enhanced turbidimetric inhibition immunoassay (PETINIA).	Homogeneous enzyme immunoassay based on competition for antibody binding sites between the analyte drug in the specimen and exogenous drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH).
Sample Matrix	Human serum or human plasma (including Li heparin, Na heparin, K EDTA, Na citrate, Na fluoride/K oxalate)	Same
Reagent	Liquid ready-to-use	Same
Analyzer	Architect cSystems	Abbott Aeroset analyzer
Claimed Assay range	1.9 to 20.0 µg/mL	0.5 to 20.0 µg/mL
Calibrators	Liquid ready-to-use, six levels (0, 2, 4, 8, 12 and 20 µg/mL) Not provided with the kit.	Liquid ready-to-use, six levels (0, 2, 4, 8, 12 and 20 µg/mL) Provided with the kit.

K. Standard/Guidance Document Referenced (if applicable):

CLSI EP5-A2; Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition

CLSI EP9-A2; Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition

CLSI EP17-A2; Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition

L. Test Principle:

The Abbott carbamazepine assay is a homogenous particle-enhanced turbidimetric inhibition immunoassay (PETINIA) used for the analysis of carbamazepine in serum or plasma. The assay is based on competition between drug in the sample and drug coated onto a microparticle of antibody binding sites of the carbamazepine antibody reagent. The carbamazepine-coated microparticle reagent is rapidly agglutinated in the presence of the anti-carbamazepine antibody reagent and in the absence of any competing drug in the sample. The rate of absorbance change is measured photometrically, and is directly proportional to the rate of agglutination of the particles. When a sample containing carbamazepine is added, the agglutination reaction is partially inhibited, slowing down the rate of absorbance change. A concentration-dependent classic agglutination inhibition curve can be obtained, with maximum rate of agglutination at the lowest carbamazepine concentration and the

lowest agglutination rate at the highest carbamazepine concentration.

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

All studies were performed on the Architect c8000.

a. Precision/Reproducibility:

Precision was evaluated following the CLSI EP5-A2 guidance. Pooled patient samples at five carbamazepine concentrations and negative serum pools spiked with carbamazepine stock solution (1 mg/mL) at five carbamazepine concentrations were tested. Testing was performed at two locations (an internal and an external site) on an Architect c8000 clinical chemistry analyzer using one lot each of the reagents and the calibrators. Each carbamazepine sample was assayed in replicates of two, twice a day for 20 days with a total of 40 runs and n=80. Each run per day was separated by at least two hours. The mean, within-run and total-run SD, and with-run and total-run %CV were calculated. The results obtained are as shown in the table below:

Sample	Mean (µg/mL)	Within-run		Total	
		SD	%CV	SD	%CV
Internal Site					
Spike 1	2.25	0.05	2.1	0.09	3.9
Spike 2	3.49	0.07	2.1	0.1	2.8
Spike 3	8.18	0.13	1.6	0.2	2.4
Spike 4	15.27	0.24	1.5	0.48	3.1
Spike 5	17.19	0.31	1.8	0.33	1.9
Patient 1	1.93	0.06	3.0	0.12	6.3
Patient 2	3.99	0.07	1.6	0.1	2.6
Patient 3	8.15	0.08	1.0	0.18	2.1
Patient 4	12.05	0.18	1.5	0.24	2.0
Patient 5	18.09	0.31	1.7	0.39	2.2
External Site					
Spike 1	2.32	0.06	2.7	0.1	4.1
Spike 2	3.17	0.07	2.2	0.09	2.9
Spike 3	7.94	0.09	1.1	0.12	1.6
Spike 4	13.46	0.32	2.4	0.38	2.8
Spike 5	15.15	0.16	1.1	0.25	1.7

b. Linearity/assay reportable range:

The claimed assay range is 1.9 to 20 µg/mL. To evaluate linearity, a series of carbamazepine concentrations in serum were prepared by spiking into negative serum and diluting the high sample to ten concentrations spanning (and exceeding) the assay range. Concentrations included carbamazepine ranged from 1.24 to 20.6 µg/mL. Four replicates were measured for each

concentration and the mean of measured concentration was compared with its expected concentration. The expected concentration was based on dilution factors and the high sample concentration. The concentration of the high sample was determined by chromatographic testing traceable to purified carbamazepine.

The linear regression results are:

Study Site	Slope	Intercept	R ²
Manufacturer Site	1.073	0.029	0.9990
External Site	0.983	0.065	0.9984

Recoveries relative to the target concentrations were within $\pm 10\%$ or within ± 0.4 ug/mL.

This evaluation supports linearity across the claimed measuring range of 1.9 to 20.0 ug/mL.

Spike/Recovery

Purified carbamazepine was gravimetrically spiked into negative human serum at 2, 8 and 18 μ g/mL. Twenty one replicates for each sample were tested in the assay. The average measured concentration was compared to its expected concentration for each spiked sample, as determined by a chromatographic method. The percent recoveries were less than $\pm 10\%$ for the mid and high samples (8 and 18 μ g/mL) and were within 0.4 μ g/mL for samples $< 4\mu$ g/mL.

To validate labeling recommendations for dilution of high samples, high concentration carbamazepine samples were prepared by adding reference material to serum. Samples were manually diluted into Calibrator 1 or saline. Three high sample concentrations (20, 25, and 30 ug/mL) were tested with three replicates and two separate 4-fold dilutions each for a total of 18 measurements diluted into saline and 18 measurements diluted into Calibrator 1. Recoveries ranged from 101-107% of expected. The automated dilution into saline was similarly tested and recoveries were within 95-101% of expected.

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Calibrators were previously cleared under k120936. See k120936 for traceability, stability, and expected values information on calibrators.

d. *Detection limit:*

The sponsor performed a study following CLSI EP17 guidelines to determine the LoB and LoQ of the Abbott Carbamazepine assay on the Architect c8000 analyzer. The LoB was calculated to be 0.1018 μ g/mL.

Precision and bias at the LoQ were evaluated by adding reference material to carbamazepine-free serum and carbamazepine-free plasma to attain concentrations ranging from 0.5 to 3.2 ug/mL. Four replicates per day were measured over 17 days with 2 lots of reagent. At the claimed lower limit of 1.9 ug/mL the upper 95% confidence limit of the CV was within <7%, and bias relative to the known concentration was within ± 0.4 ug/mL.

e. Analytical specificity:

Interference by metabolites, drugs, and endogenous compounds was tested by spiking the potential interferents into human serum containing carbamazepine at concentrations of 4, 8, and 12 ug/mL. Corresponding control samples were prepared by adding the equal volume of solvent that was in the stock solution of the compound tested. Samples were measured in duplicate by the carbamazepine assay. The compounds listed below did not interfere at the concentrations shown below in the table. Recoveries were all within $\pm 6\%$ of expected.

Compound tested for interference	Concentration ($\mu\text{g/mL}$)
5-(p-Hydroxyphenyl)-5-phenylhydantoin	1000
Acetaminophen	200
N-acetylcysteine	150
Acetylsalicylic acid	1000
Amitriptyline	2
Amobarbital	50
Ampicillin-Na	100
Ascorbic acid	30
Cefoxin	2500
Cetirizine dihydrochloride	3
Chlordiazepoxide	30
Chlorpromazine	100
Clonazepam	12
Cyclosporine	5
Desipramine	3
Diazepam	25
Ethosuximide	1000
Ethotoin	50
Glutethimide	50
Hydroxyzine Dihydrochloride	1
Ibuprofen	500
Imipramine	6

besilate (hydroquinonesulfonic acid potassium salt)	200
Levodopa (3,4-Dihydroxy-L-phenylalanine)	20
Mephenytoin	150
Methsuximide	50
Methyldopa sesquihydrate ((-)-3-(3,4-Dihydroxyphenyl)-2-methyl-L-alanine sesquihydrate)	20
Metronidazole	200
Nortriptyline	1
Phenobarbital	500
Phenothiazine	200
Phenylbutazoneeeee	16
Phenytoin	1000
Primidone	1000
Probenecid	500
Promethazine	100
Probenecid	60
Secobarbital	50
Tetracycline	50
Theophylline	100
Valproic Acid	1000

The following compounds were tested for cross-reactivity at the following three serum carbamazepine concentrations. Results are shown in the table:

Compound	Percent Bias/Recovery			% -Cross-reactivity		
	4 ug/mL CBMZP	8 ug/mL CBMZP	12 ug/mL CBMZP	4 ug/mL CBMZP	8 ug/mL CBMZP	12 ug/mL CBMZP
Eslicarbazepine	4.3%	5%	1%	9.0%	20%	6%
Carbamazepine-10,11-epoxide	15.2%	.9	8.2	31.5%	35%	549.8
10-Hydroxycarbamazepine	6.2%	6.3	-1.3%	12.8%	24%	-8%
Oxcarbazepine	3.5%	3.9	8.2	7.2%	15%	49.8%

Potentially interfering endogenous compounds were spiked at the highest possible

concentrations expected from the intended use population into human serum pool containing approximately 4, 8, and 12 µg/mL of carbamazepine. The potential interference of Rheumatoid Factor (RF) and triglyceride was assessed using RF patient specimens and two triglyceride specimens naturally high in these substances and then spiking with carbamazepine at approximately 8µg/mL. Corresponding control samples were prepared by supplementing the human serum pool with the equal volume of solvent that was in the stock solution of the compound tested or by spiking carbamazepine into negative human serum. Each prepared sample and its control solution were measured in replicates of three. The following compounds at listed concentrations cause minimal interference to the assay. Carbamazepine recoveries observed were within ±3% of expected, and no trends were observed.

Compound	Concentration tested
Bilirubin, unconjugated	60.0 mg/dL
Bilirubin, conjugated	30.0 mg/dL
Hemoglobin	800.0 mg/dL
Triglycerides	1000.0 mg/dL
Cholesterol	500.0 mg/dL
Human serum albumin (HSA)	7.5 g/dL
Gamma globulin (IgG)	12.0 g/dL
HAMA	400.0 ng/mL
Rheumatoid factor (RF)	919.0 IU/mL

Recoveries for the potential interferents above at all carbamazepine concentrations were within +/- 8% of expected values and no trends were observed.

f. Assay cut-off:

Not applicable

2. Comparison studies:

a. Method comparison with predicate device:

Two correlation studies were performed using a modified CLSI guideline EP9-A2 to evaluate the accuracy of the Abbott's Carbamazepine Assay by comparing patient sample measurements to the predicate Abbott Aeroset Carbamazepine Assay. In the first study, 103 serum samples with carbamazepine concentrations ranging from 2.2 to 18.2 µg/mL were tested by the candidate carbamazepine assay (on-test) and Abbott Aeroset Carbamazepine assay. In the second study, a different set of 105 serum samples ranging from approximately 2-15 µg/mL were evaluated with the candidate carbamazepine assay (on-test) and Abbott Aeroset Carbamazepine assay.

Data analysis was performed using Passing-Bablok regression model. Results are summarized below:

Comparative Methods	N	Slope (95% CI)	Intercept (95% CI)	Correlation Coeff. (R)
Manufacturer Site				
On-test vs. Abbott Aeroset assay	103	0.905 (0.857 – 0.960)	0.564 (0.202 – 1.00)	0.9675
External Laboratory				
On-test vs. Abbott Aeroset assay	105	0.956 (0.923 – 0.989)	0.013 (-0.208 – 0.283)	0.9827

b. Matrix comparison:

Matrix comparison studies were performed to evaluate the performance of the following: serum in plastic, serum separator tubes (SST) in plastic, plasma with sodium fluoride/potassium oxalate in plastic, plasma with sodium heparin in plastic and glass, plasma with lithium heparin in plastic with or without gel, plasma with K3 EDTA in glass and plastic, plasma with K2 EDTA in plastic, and sodium citrate in plastic and glass, compared to the control matrix (serum in glass) in the candidate carbamazepine assay.

Samples were spiked with carbamazepine at multiple concentrations spanning the assay range. Each spiked sample was measured in duplicate by the Carbamazepine Assay. The results of samples prepared in the evaluating matrix were compared to those of serum in glass by determining the slope, intercept, and correlation coefficient, and recovery for each sample relative to expected.

Matrix Comparison	N	Range $\mu\text{g/mL}$	Slope	Intercept	Correlation Coeff. (R)
X: Serum in glass	24	1.60-19.10	0.991	0.030	0.9982
Y: Serum in plastic					
X: Serum in glass	24	1.60-19.10	1.003	-0.086	0.9986
Y: SST in plastic					
X: Serum in glass	23	1.35-19.85	0.997	-0.066	0.9987
Y: Plasma w/ NaF / K Oxalate					
X: Serum in glass	25	1.51-18.53	1.007	0.029	0.9986
Y: Plasma w/ Na Heparin in glass					
X: Serum in glass	25	1.51-18.53	0.969	0.196	0.9709
Y: Plasma w/ Na Heparin in plastic					
X: Serum in glass	24	1.81-19.80	0.988	0.095	0.9986
Y: Plasma w/ Li Heparin w/gel in plastic					
X: Serum in glass	24	1.81-19.80	0.993	-0.115	0.9992

Y: Plasma w/ Li Heparin w/o gel plastic					
X: Serum in glass	23	1.36-18.70	1.003	-0.082	0.9982
Y: Plasma w/ K3 EDTA in glass					
X: Serum in glass	23	1.52-19.15	1.010	-0.087	0.9988
Y: Plasma w/ K3 EDTA in plastic					
X: Serum in glass	23	1.36-18.70	0.988	0.115	0.9993
Y: Plasma w/ K2 EDTA in plastic					
X: Serum in glass	24	1.38-18.61	1.042	-0.302	0.9985
Y: Plasma w/ Na Citrate in glass					
X: Serum in glass	24	1.38-18.61	1.047	-0.275	0.9988
Y: Plasma w/ Na Citrate in plastic					

Recoveries were largely within $\pm 4\%$ of the expected value (based on serum in glass) and were all within $\pm 8\%$.

3. Clinical studies:

a. *Clinical Sensitivity:*

Not applicable

b. *Clinical specificity:*

Not applicable

c. Other clinical supportive data (when a. and b. are not applicable):

Not applicable

4. Clinical cut-off:

Not applicable

5. Expected values/Reference range:

The package insert includes the following:

Therapeutic effect is usually achieved in the serum and plasma carbamazepine concentration range of 4.0 to 12.0 $\mu\text{g/mL}$ (17.0 to 51.0 $\mu\text{mol/L}$). Peak carbamazepine concentrations above 15.0 $\mu\text{g/mL}$ (63.0 $\mu\text{mol/L}$) are often associated with toxicity.^{1,2}

In addition, the sponsor has included the following in the labeling: For effective treatment, some patients may require serum levels of carbamazepine outside these ranges. Therefore, the expected range is provided only as a guide, and individual patient results should be interpreted in light of other clinical signs and symptoms.

¹Levy, R.H., Wilensky, A.J., Freil, P.N. Other antiepileptic drugs. In: Evans, W.E., Schentag, J.J., Jusko, W.J., editors. Applied Pharmacokinetics: Principles of Therapeutic Drug Monitoring, 2nd ed. Spokane, W.A.: Applied Therapeutics; 1986: 540-569.

²Goldman, L., Ausiello, D., editors. Cecil Medicine, 23rd ed. Philadelphia, PA: Elsevier Saunders; 2008:2994.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.