

JUN 29 1998

K 980887

510(k) Summary

Introduction

According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

1. Submitter name, address, contact

Boehringer Mannheim Corporation
4300 Hacienda Drive
P.O. Box 9002
Pleasanton, CA 94566-0900
(510) 730-8215

Contact Person: Patricia M. Klimley
Date Prepared: March 6, 1998

2. Device name

Proprietary name: Elecsys® CEA Assay

Common name: Electrochemiluminescence assay for the determination of Carcinoembryonic antigen (CEA).

Classification name: Kit, Test , Carcinoembryonic antigen

3. Predicate device

The Boehringer Mannheim Elecsys® CEA on Elecsys® 1010 is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed Elecsys® CEA on Elecsys® 2010.

4. Device Description

The Elecsys® test principle is based on sandwich principle. Total duration of assay: 18 minutes (37° C).

- 1st incubation (9 minutes): Sample (30 µL), biotinylated monoclonal CEA-specific antibody (60 µL), and a monoclonal CEA-specific antibody labeled with a ruthenium complex (60 µL) react to form a sandwich complex.

- 2nd incubation (9 minutes): After addition of streptavidin-coated microparticles (50 µL), the complex is bound to the solid phase via interaction of biotin and streptavidin.

Continued on next page

510(k) Summary, Continued

-
- | | |
|--|---|
| 4.
Device
Description
(con't) | <ul style="list-style-type: none">•The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier (0.4 second read frame).•Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent bar code. |
|--|---|
-
- | | |
|----------------------------------|--|
| 5.
Intended use | <p>Immunoassay for the in vitro quantitative determination of carcinoembryonic antigen (CEA) in human serum and plasma. The Elecsys CEA assay is further indicated for serial measurement of CEA to aid in the management of cancer patients.</p> <p>The electrochemiluminescence immunoassay "ECLIA" is intended for use on the Boehringer Mannheim Elecsys 1010 and 2010 immunoassay analyzers.</p> <p>Summary</p> |
|----------------------------------|--|

**6.
Comparison to
predicate device**

The Boehringer Mannheim Elecsys® CEA Assay has been approved for use on the Elecsys 2010 immunoassay analyzer (K964368). The application of the Elecsys® CEA Assay on the Elecsys 1010 immunoassay analyzer is substantially equivalent to the same assay (Elecsys CEA Assay) on the Elecsys 2010.

The following table compares the Elecsys® CEA Assay on Elecsys® 1010 with the predicate device, Elecsys® CEA Assay on Elecsys® 2010 . Specific data on the performance of this test for both the Elecsys 1010 and 2010 have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device in attachment 6 will be replaced upon the clearance of this premarket notification submission with the combined Elecsys 2010 and 1010 insert (attachment 5).

Similarities:

- Intended Use: Immunoassay for the in vitro quantitative determination of Carcinoembryonic Antigen (CEA). The assay is further indicated for the serial measurement of CEA to aid in the management of cancer patients.
 - Assay range: 0-1000 ng/mL
 - Assay methodology: Sandwich immunoassay
 - Kit (cat. No.) also cleared for use on the Elecsys 2010 (K964368)
 - Sample and reagent volumes
 - Package insert
 - Performance specifications
-

510(k) Summary, Continued

6. Comparison to predicate device cont.

Differences:

Feature	Elecsys® 1010	Elecsys® 2010
Instrument required	Elecsys 1010	Elecsys 2010
Instrument Type	Batch	Random access
Reagent Storage Temp (C)	Ambient Temperature	20° C

Performance Characteristics:

Feature	Elecsys® 1010		Elecsys® 2010		
Precision	Modified NCCLS (ng/mL):		Modified NCCLS (ng/mL):		
Level	<u>Control 1</u>	<u>Control 2</u>	<u>Control 1</u>	<u>Control 2</u>	<u>Pool 1</u>
N	60	60	60	60	60
Within-Run	4.44	32.96	4.9	34.1	2.2
%CV	2.36	1.71	2.5	1.7	5.0
Total	4.44	35.51	4.9	34.1	2.2
%CV	3.07	2.09	3.6	3.0	5.4
	Modified NCCLS (ng/mL):		Modified NCCLS (ng/mL):		
	<u>Pool 2</u>	<u>Pool 3</u>	<u>Pool 2</u>	<u>Pool 3</u>	
N	60	60	60	60	
Within-Run	11.86	113.64	19.6	528	
%CV	1.57	2.89	1.6	1.3	
Total	11.86	113.64	19.6	528	
%CV	2.83	3.97	2.3	2.0	

Continued on next page

510(k) Summary, Continued

Performance Characteristics:

Feature	Elecsys® 1010	Elecsys® 2010
Lower Detection Limit	0.2 ng/mL	0.2 ng/mL
Linearity	0.2 - 1000 ng/mL (with a deviation from a linear line of $\pm 10\%$)	0.2 - 1000 ng/mL (with a deviation from a linear line of $\pm 10\%$)
Method Comparison	<p>Vs Elecsys 2010</p> <p><u>Least Squares</u> $y = 1.006x - 0.64$ $r = 0.995$ $N = 117$</p> <p><u>Passing/Bablok</u> $y = 0.958x + 0.14$ $r = 0.995$ $N = 117$</p>	
Hook Effect	No Hook Effect up to 200,000 ng/mL CEA	No Hook Effect up to 200,000 ng/mL CEA



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUN 29 1998

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Ms. Patricia M. Klimley
Manager, Regulatory Affairs
Boehringer Mannheim Corporation
4300 Hacienda Drive
P.O. Box 9002
Pleasanton, California 94566-0900

Re: K980887/S1
Trade Name: Elecsys® CEA Assay
Regulatory Class: II
Product Code: DHX
Dated: May 4, 1998
Received: May 5, 1998

Dear Ms. Klimley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal Laws or Regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770)488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>"

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical
Laboratory Devices
Office of Device Evaluation
Center for Devices and
Radiological Health


Enclosure

510(k) Number (if known): K980887

Device Name: Elecsys® CEA Assay

Intended use

Immunoassay for the in vitro quantitative determination of carcinoembryonic antigen in human serum and plasma. The Elecsys CEA assay is further indicated for serial measurement of CEA to aid in the management of cancer patients. The electrochemiluminescence immunoassay "ECLIA" is intended for use on the Boehringer Mannheim Elecsys 1010 and 2010 immunoassay analyzers.



(Division Sign Off)
Division of In Vitro Diagnostic Laboratory Devices
510(k) Number K980887

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

(Optional Format 1-2-96)