

OCT 25 1999

Page 1 of 3

510(k) SUMMARY

K991853

PREPARED BY: INTERNATIONAL DISTRIBUTORS OF
ELECTRONICS FOR MEDICINE, INC.
(IDEM)
4814 East Second Street
Benicia, CA 94510

CONTACT PERSON: Donna Ward, President

TELEPHONE: 707-746-6334

DATE ON WHICH THE SUMMARY
WAS PREPARED: May 28, 1999

NAME IF DEVICE: Interacoustics Model EP15
ABR Stand Alone Unit

COMMON NAME: Evoked Response Auditory Stimulator

PREDICATE DEVICE: ICS Medical Corp. Chartr EP System

DESCRIPTION OF DEVICE:

The Interacoustics EP15 ABR Stand Alone Unit is a diagnostic tool for computer-based testing. Data is collected in a Windows 95/98 database format that allows for easy inspection of results. The EP15 Stand Alone Unit produces a sound stimulus for use in evoked response measurements or electroencephalogram activation.

Comparison of the Interacoustics Model EP15 and the ICS Medical Chartr EP System:

Indication for use – Identical for both units.

Similarities and differences:

EP15 ABR	Chartr EP
Display Description: 12.1" Wide Angle Industrial TFT	15" SVGA Monitor (others available)
Software Features: Windows 95/98 Operating System	Same
AMPLIFIERS: Channels: 2 Gain: 80 dB	2 (optionally up to 8) X50 - X5000000
PREAMPLIFIERS: Frequency response: 100 to 8000 Hz CMR Ratio: >120dB at 50/60 Hz. Noise: 5.5nV/Hz, 0.3µV RMS (0– 3 kHz); 5.5nV/Hz, 0.5µV RMS (100-8000 Hz)	Same >100dB at 50/60 Hz <1.5 µV RMS (0.1 Hz – 5 kHz)
Points Per Trace: 450 (30 kHz sample)	600
ANALOG FEATURES: Low pass: 1000, 1500, 2000, 3000, 5000 Hz, none (12000Hz). 33 taps FIR filter. High pass: None (5 Hz), 50, 100, 150, 300 500 Hz. 6dB/octave.	15 Hz to 25 kHz, 12 dB/octave DC, 0.002 – 1 kHz, 6 dB/octave
STIMULATORS: Stimuli: Click, Tone Burst, Blackman, Gaussian, Hanning, Hamming, Bartlett, Rectangle, Manual (rise/plateau/fall) Parameters for tonal stimuli; frequency, intensity, rise/fall time, plateau duration, envelope shape: Programmable Rate: 7.1, 13.1, 21.1, 23.1, 29.1, 31.1, 35.1, 41.1, 45.1, 49.1 Hz, ext. trigger input Intensity: 20 dB to 130 dB peSPL; 10 dB to 100 dB HL in 1 dB step. Masking: White noise; masking ear is opposite side of stimulus Transducers: EAR 3A Insert Phones-ABR	Click, Tone Burst Programmable 0.1 – 100/sec -12 dB to +128 dB SPL White noise, programmable intensity TDH-49 earphones, insert earphones, bone conduction transducer
SAFETY CHARACTERISTICS: Optical isolation: Yes Built-in isolation transformer: Yes	Yes Yes

(CONTINUED – COMPARISON)	
EP15 ABR	Chartr EP
Controlled parameters: Stimuli Rate; Number of Stimuli; Stimuli Polarity; Click; Tone Burst (frequency, number of sin waves, window); Stimulus Intensity; Number of Curve Pr. Intensity; Intensity (ascend, descend); Soft Attenuator; Stimulus Ear (right, left, simultaneously); Masking Level; Preliminary Filter Setting (low, high pass filter); Recording Onset; Automatic Next Intensity (wave repro level setting); Ext. Trigger Output Duration; Rejection System Rejection Level; Gain (manual, automatic); Display Options (Invert curves on screen; Origin line; Latency Norm. Report Templates; After Filtering; Print out; Manual Stimulus to Familiarization; Talk Forward; Talk Back Monitor.	Stimulated Ear, Masked Ear, Stimulus Intensity, Masking Stimulus, Stimulus Transducer, Stimulus Type, Stimulus Polarity, Stimulus Characteristics, Number of Sweeps Acquires, Stimulus Presentation Rate, Sweep Time, Number of Channels, Amplifier Gain, Filter Characteristics, Inclusion of Notch Filter, Inclusion of Artifact Rejection
Data Collection: Impedance Test; Waveform Buffer (A/B, contra, ipsi-contra, A-B=Noise); Curve (Hide, Fixate, Merge, Delete); Show Online EEG, Store Waveforms in unlimited Storage database	Tests Impedance of Patient Electrode Connections, Display Waveform Buffers During Examination, Displays on-going EEG Activity, Stores Waveforms, Stores the Waveform Presentation
Dimensions: 14" x 10" x 15"	22.5" x 17.5" x 7.8"
Weight: 26.5 lbs.	58 lbs.
Power Supply: Input volts: 90 to 250 VAC Universal Input Switch Mode; Safety: VDE750, EN60601-1, IEC601, IEC1010, UL544, CSA 22.2	200 Watts, 50/60 Hz, 120 or 240 volts
Keyboard: 101-key IBM standard	101 – key IBM standard
Ancillary Functions: Help system. Exports data of one patient to diskette.	Help system. Exports data of one patient to diskette. Backup procedure automatically implemented.

SAFETY AND EFFECTIVENESS:

The Interacoustics Model EP15 ABR Stand Alone Unit is in compliance with the following performance and safety standards:

VDE750; EN 60601-01 (General Safety) Class I, Type BF; EN 60601-1-1 (Safety of Systems); EN 60601-1-2 (EMC); EN 60601-2-26 (Electroencephalographs); EN 60645-3 (Auditory test signals)



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

APR - 9 2012

Ms. Donna Ward
President
IDEM INTERNATIONAL DISTRIBUTORS OF ELECTRONICS FOR MEDICINE, INC.
4814 East Second Street
Benicia, California 94510

Re: K991853

Trade/Device Name: Interacoustics Model EP15 ABR Stand Alone Unit
Regulation Number: 21 CFR 882.1900
Regulation Name: Evoked response auditory stimulator
Regulatory Class: II
Product Code: GWJ
Dated (Date on orig SE ltr): August 5, 1999
Received (Date on orig SE ltr): August 20, 1999

Dear Ms. Ward:

This letter corrects our substantially equivalent letter of October 25, 1999.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



 Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic, Neurological,
and Ear, Nose and Throat Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K 991853Device Name: Interacoustics Model EP15 ABR Stand Alone Unit

Indications For Use:

The Interacoustics EP15 ABR Stand Alone Unit is a diagnostic tool for computer-based testing. Data is collected in a Windows 98 database format that allows for easy inspection of results. The EP15 Stand Alone Unit produces a sound stimulus for use in evoked response measurements or electroencephalogram activation.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

[Signature]
(Division Sign-Off)
Division of General Restorative Devices K991853
510(k) Number K991853

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

(Optional Format 1-2-96)