



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

November 21, 2016

Beijing Choice Electronic Technology Co., Ltd.
Lei Chen
Quality Director
No.9 Shuangyuan Road, Badach Hi-tech Zone
Shijingshan District
Beijing, 100041 CN

Re: K160508

Trade/Device Name: Electronic Pulse Stimulator
Regulation Number: 21 CFR 882.5890
Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief
Regulatory Class: Class II
Product Code: NUH
Dated: February 22, 2016
Received: February 24, 2016

Dear Lei Chen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in

the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

William J.
Heetderks -A

Digitally signed by William J. Heetderks -A
DN: c=US, o=U.S. Government, ou=HHS, ou=NIH,
ou=People,
0.9.2342.19200300.100.1.1=0010149848,
cn=William J. Heetderks -A
Date: 2016.11.21 21:22:10 -05'00'

for Carlos L. Peña, PhD, MS
Director
Division of Neurological
and Physical Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Section II Indications for Use Statement

Indications for Use

510(k) Number (if known): _____

Device Name: Electronic Pulse Stimulator

Model: MDTS100

Indications for Use:

The Electronic Pulse Stimulator MDTS100 is to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household and work activities.

Prescription Use _____

AND/OR

Over-The-Counter Use ✓

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

510(k) Summary

This summary of 510(k) is being submitted in accordance with the requirements of 21 CFR 807.92.

3.1 Submitter Information

Manufacturer Name:

Establishment Registration Number: 3005569927
Beijing Choice Electronic Technology Co., Ltd.
Room 4104, No. A12 Yuquan Road Haidian District
100143 Beijing, P.R. China

Contact Person:

Mr. Lei Chen
Beijing Choice Electronic Technology Co., Ltd.
North Building 3F, No. 9 Shuangyuan Road,
Badachu Hi-tech Zone, Shijingshan District
Beijing China 100041

Phone: +86-10-88798300 Ext 6020

Fax: 215-4052545

Email: cc@choicemed.com

Date prepared: October 31, 2016

3.2 Proposed Device Information

Device Common Name: Transcutaneous Electrical Nerve Stimulator, OTC

Device Trade/Proprietary Name: Electronic Pulse Stimulator

Model: MDTS100

Classification Name: Transcutaneous Electrical Nerve Stimulator for pain relief

Regulation Number: 21 CFR 882.5890

Product Code: NUH

Class: II

Use: Over -The-Counter Use

Premarket Notification 510(k) Submission—510(k) Summary

Panel: Neurology

3.3 Predicate Device

510(k) Number: K122744

Common Name: Transcutaneous electrical nerve stimulator OTC

Device Trade/Proprietary Name: Prospera OTC TENSE Electronic Pulse Massager

Model: PL029

Classification Name: Transcutaneous Electrical Nerve Stimulator for pain relief

Product Code: NUH, NGX

Regulation Number: 21 CFR 882.5890

Device Class: II

Panel: Neurology

Use: Over-The-Counter Use

Manufacturer: Prospera Corporation

Intended Use:

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities(arm), and lower extremities(leg) due to strain from exercise or normal household work activities.

3.4 Device Description

The MDTS100 Electronic Pulse Stimulator is an electrically powered device intended for over the counter use and used to apply an electrical current to electrodes on a user's skin to relieve pain. The proposed device is intended for at home use in delivering electric pulses to tired and sore muscles.

The proposed device consists of controlling host device, output connecting cable and skin electrodes. The proposed device has many safety measures such as output load detection, timing indication and so on. It is adopted battery as the power supply. The output waveform is adopted "square waveform" as the basic wave. It simulates the various output modes by various frequency combination. The waveform is mainly for single rectangular. The output waveform intensity is continuous adjustable. The working time is controlled automatic.

The MDTS100 is adopted the method that inputting the specific low frequency pulse

Premarket Notification 510(k) Submission—510(k) Summary

current into the human body through the human skin to treat the pain.

The MDTS100 Electronic Pulse Stimulator has seven stimulation modes. These stimulation modes can be selected by pressing the MODE button on the front of the unit.

The MDTS100 comes with 2 independent output channels by pressing the CHA button.

Channel A&B are selected by default. Pressing the ◀ & ▶ buttons will decrease & increase the intensity of stimulation. LCD screen display the left remaining stimulation time in the lower left corner. The default time is 20 minutes.

The proposed device is not for life-supporting or life-sustaining, not for implant.

The device is not sterile and the transducers are reusable and do not need sterilization and re-sterilization.

The device is for over the counter use.

The device does not contain drug or biological products.

The device is software-driven and software validation is provided in *software*.

3.5 Indication for Use

The Electronic Pulse Stimulator MDTS100 is to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household and work activities.

3.6 Comparison list of technical characteristics

Table 3-1 Basic Unit Characteristics Comparison Table between the Proposed Device (MDTS100) and Predicate Device

Parameter		Proposed Device	Predicate Device
510(k) number		(to be assigned)	K122744
Device Name and Model		MDTS100	PL029
Manufacturer		Beijing Choice Electronic Technology Co., Ltd.	Prospera
Power Source		DC 3V, 2 AAA batteries	3V Battery
-Method of Line Current Isolation		N/A for DC current	Battery Supply N/A
-patient Leakage Current		--	---
-Normal Condition(μA)		$\leq 10\mu\text{A}$	$2.0\mu\text{A}$
-Single Fault Condition(μA)		$\leq 50\mu\text{A}$	$3.0\mu\text{A}$
Average DC current through electrodes when device is on but no pulses are being applied(μA)		$\leq 10\mu\text{A}$	$0\mu\text{A}$
Number of Output Modes		7	8
Number of Output Channels	Synchronous or Alternating?	Alternating	Alternating
	Method of Channel Isolation	No	By software
Regulated Current or Regulated Voltage?		Regulated Voltage	Voltage control
Software/Firmware/Microprocessor Control?		Yes	Yes

Premarket Notification 510(k) Submission—510(k) Summary

Automatic Overload Trip?		No	No
Automatic No-Load Trip?		Yes	No
Automatic Shut Off?		Yes	Yes
User Override Control?		Yes	Yes
Indicator Display:	On/OFF Status?	Yes	Yes
	Low Battery?	Yes	No
	Voltage/Current Level?	Yes	No
Timer Range (minutes)		20 minutes	5minutes, 10 minutes
Compliance with Voluntary Standards?		Yes	Yes
Compliance with 21 CFR 898?		Yes	Yes
Weight (lbs., oz.)		0 lb., 2.2 oz. (Battery Excluded)	0 lb, 8.21oz.
Dimensions (in.) [W x H x D]		2.17x 5.12x0.87	2.24x7.87x0.83
Housing Materials and Construction		Enclosure: ABS	Enclosure: ABS

Table 3-2 Output Specification Comparison Table between the Proposed Device (MDTS100) and Predicate Device

Parameter	Proposed Device	Predicate Device
Device Name and Model	MDTS100	PL029
Waveform	pulsed monophasic	Monophasic
Shape	Rectangular	Rectangular
Maximum Output Voltage (volts) (+/- <u>20</u> %)	<u>150V</u> @500 Ω	<u>49.6V</u> @500 Ω
	<u>160V</u> @2k Ω	<u>99.2V</u> @2k Ω
	<u>165V</u> @10k Ω	<u>114V</u> @10k Ω
Maximum Output Current (specify units) (+/- <u>20</u> %)	<u>300mA</u> @500 Ω	<u>18mA</u> @500 Ω
	<u>80mA</u> @2k Ω	<u>3.2mA</u> @2k Ω

Premarket Notification 510(k) Submission—510(k) Summary

		<u>16.5mA @10k Ω</u>	<u>0.6mA @10k Ω</u>
Duration of primary (depolarizing) phase (μ sec)		0 μ sec	40msec
Pulse Duration (μ sec)		50~140 μ sec	120~6800 μ sec
Frequency (Hz) [or Rate (pps)]		0.9Hz~82Hz	<u>0.5~86Hz</u>
Net Charge (microcoulombs (μ C) per pulse) @500 Ω		<u>42μC</u>	<u>18000μC</u>
Maximum Phase Charge, (μ C) @500 Ω		<u>42μC</u>	<u>23μC</u>
Maximum Current Density, (mA/cm ² , r.m.s.)		3.3mA/cm ² @1.57"×1.57"Electrode Pad	1.4mA/cm ²
Maximum Average Current (average absolute value), mA		2.8mA	---
Maximum Average Power Density, (W/cm ²), (using smallest electrode conductive surface rate)		0.02W/cm ² @1.57"×1.57"Electrode Pad	0.23 W/cm ²
Burst Mode	(a) Pulses per burst	197	1
	(b) Burst per second	0.185	0-7
	(c) Burst duration (seconds)	2.4s	1
	(d) Duty Cycle: Line(b) x Line(c)	0.44	7
ON TIME (seconds)		≤1second	40ms
OFF TIME (seconds)		≤1second	18ms
Additional Features		---	---

Premarket Notification 510(k) Submission—510(k) Summary

Conclusion: The data in Table 3-1 and Table 3-2 indicates that the technical characteristics, features, specifications and intended use of the two devices are very similar. Both systems are designed to relief pain in the muscle areas. Outputs are almost the identical, each can produce square electronic pulses and positive-going waveform output. Both devices have selection modes to allow users to select preset programs and adjust intensity level by pressing the keypad on the main unit. Both devices are transiting pulses by the use of Electrode Pads. Technological characteristics, features, specifications, materials and intended uses of the new device are substantially equivalent to the quoted predicated device. The parameters such as waveform, frequency, pulse width, output current and voltage, although not exactly similar to this particular predicate, they do not raise any new questions of safety or effectiveness. The differences that exist between the devices are insignificant in terms of safety and effectiveness.

3.7 Functional and Safety Testing:

Non-clinical tests were performed on the subject matter device in order to validate the design and to assure conformance with the following voluntary design standards in connection with medical device electrical safety, electromagnetic compatibility, and particular requirements for the safety of nerve and muscle stimulator:

EMC and electrical Safety

IEC 60601-1:2005 Medical electrical equipment-Part 1: General requirements for basic safety and essential performance

IEC 60601-1-2:2007 Medical electrical equipment-Part 1-2: General requirements for basic safety and essential performance-Collateral Standard: Electromagnetic compatibility-Requirements and tests.

IEC 60601-1-11: 2010 Medical electrical equipment-Part 1-11: General requirements for basic safety and essential performance-Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

IEC 60601-2-10: 2012 Medical electrical equipment-Part 2-10: Particular requirements for the safety of nerve and muscle stimulator.

Biocompatibility

The electrode pads are the mainly body contacting parts. Electrode pads, as accessories of the subject matter device, are adopted the Self-adhesive Electrode manufactured by Wuxi Jiajian Medical Instrument Co., Ltd. and have been cleared by FDA in June 23, 2009 as K090198. The patient contacting electrodes comply with ISO 10993 and are considered Class A $\leq$ 24 hours contact.

Software

In addition to the compliance of voluntary standards, the software verification has been carries out according to the FDA Guidance for the Content of Premarket Submission for Software Contained in Medical Devices.

Cleaning

The cleaning instructions as described in the Instruction Manual have been tested to be sufficient. The test result was presented in the system test report in Performance Testing.

3.8 Conclusion

The Electronic Pulse Stimulator MDTS100 has the similar intended use and technological characteristics as the predicate device. Moreover, bench testing and safety report documentation demonstrate that the submitted device could maintain the same safety and

Premarket Notification 510(k) Submission—510(k) Summary

effectiveness as that of the predicate device. In other words, the differences do not affect the intended use and do not raise any new questions of safety or effectiveness or alter the fundamental scientific technology of the device. Thus, the proposed device MDTS100 is substantially equivalent to the predicate device.