



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

November 3, 2016

L&K Biomed Co., Ltd.
Yerim An
RA Specialist
#201, 202 16-25, Dongbaekjungang-ro 16 Beon-gil
Giheung-gu, Yongin-si, Gyeonggi-do 17015
KOREA

Re: K162801

Trade/Device Name: CastleLoc Spinal Fixation System
Regulation Number: 21 CFR 888.3070
Regulation Name: Pedicle Screw Spinal System
Regulatory Class: Class II
Product Code: MNH, MNI
Dated: October 4, 2016
Received: October 5, 2016

Dear Yerim An:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in

the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

Vincent J. Devlin -S

for

Mark N. Melkerson

Director

Division of Orthopedic Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K162801

Device Name

CastleLoc Spinal Fixation System

Indications for Use (Describe)

The CastleLoc Spinal Fixation System is a pedicle screw system indicated for the treatment of severe Spondylolisthesis (Grade 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft having implants attached to the lumbar and sacral spine (L3 to sacrum) with removal of the implants after the attainment of a solid fusion. In addition, the CastleLoc Spinal Fixation System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine: degenerative Spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor and failed previous fusion (pseudarthrosis).

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY

The following 510(k) summary is being submitted as required by 21 CFR 807.92(a):

1. **Submitter:** L&K BIOMED Co., Ltd.
 #201, 202 16-25, Dongbaekjungang-ro 16 beon-gil
 Giheung-gu, Yongin-si, Gyeonggi-do, 17015, Korea
 Phone. 82-2-6717-1985
 e-mail: yerim2706@lnkbiomed.com
- Contact Person:** Yerim An
- Date prepared** Nov, 02, 2016

2. Device Identification

Trade Name	CastleLoc Spinal Fixation System
Classification	Class II Pedicule Screw Spinal System Orthopaedic 21 CFR 888.3070 Product Code: MNH, MNI

3. Purpose of 510(k)

The L&K BIOMED Co. Ltd., here by submits this special 510(k) for adding new type of screws and adding new brand name.

4. Predicate or legally marketed devices which are substantially equivalent

- Primary Predicate Device: K103085 VENUS BASIC Spinal Fixation System (LnK Basic Spinal Fixation System)

5. Description of the Device

The CastleLoc Spinal Fixation System is a top-loading multiple component, posterior spinal fixation system which consists of pedicle screws, rods, set screws, and a transverse (cross) linking mechanism.

The CastleLoc Spinal Fixation System will allow surgeons to build a spinal implant construct to stabilize and promote spinal fusion. CastleLoc Spinal Fixation System implants components are supplied non-sterile are single use and are fabricated from titanium alloy (Ti-6Al-4V ELI) that conforms to ASTM F136.

The purpose of this submission is to add components of the type2 pedicle screw in CastleLoc Spinal Fixation System. Various sizes of these implants are available.

6. Indication for Use

The CastleLoc Spinal Fixation System is a pedicle screw system indicated for the treatment of severe Spondylolisthesis (Grade 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft having implants attached to the lumbar and sacral spine (L3 to sacrum) with removal of the implants after the attainment of a solid fusion. In addition, the CastleLoc Spinal Fixation System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine: degenerative Spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor and failed previous fusion (pseudarthrosis).

7. Comparison of the technological characteristics of the subject and predicate devices

The CastleLoc Spinal Fixation System is considered substantially equivalent to legally marketed devices LnK Basic Spinal Fixation System. They are similar in design, material, scientific technologies and indications for use. The only different thing is shape of bone screw. We have evaluated equivalency of bone screw shape via worst case justification.

8. Performance Testing

Additional components of this system are not worst case. To verify worst case, we have performed mechanical test followed below standard.

ASTM F 543

9. Conclusion

The CastleLoc Spinal Fixation System is substantially equivalent to legally marketed predicates.