



January 5, 2018

Zhejiang kindly Medical Devices Co., Ltd
% Diana Hong
General Manager
Mid-Link Consulting Co., Ltd
P.O. Box 120-119
Shanghai, 200120 China

Re: K171518

Trade/Device Name: Epidural Anesthesia Needles, Spinal Anesthesia Needles, Combined Anesthesia Needles

Regulation Number: 21 CFR 868.5150

Regulation Name: Anesthesia Conduction Needle

Regulatory Class: Class II

Product Code: BSP

Dated: December 11, 2017

Received: December 13, 2017

Dear Diana Hong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good

manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Tara A. Ryan -S
2018.01.05 11:04:19 -05'00'

for

Tina Kiang, Ph.D.

Acting Director

Division of Anesthesiology,

General Hospital, Respiratory,

Infection Control, and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K171518

Device Name

Spinal Anesthesia Needles, Epidural Anesthesia Needles, Combined Anesthesia Needles

Indications for Use (Describe)

The Epidural Anesthesia Needles are intended to be used for injection into the epidural space/or placing the epidural catheter into the epidural space.

The Spinal Anesthesia Needles are intended to be used for injection of local anesthetic agent into the subarachnoid cavity for pain management.

The Combined Anesthesia Needles are intended for injection of local anesthetics into the spinal and epidural spaces of a patient to provide regional anesthesia. The administration of the spinal anesthesia allows rapid anesthesia onset and the placement of an epidural catheter allows for bolus injections or continuous infusion of local anesthetics or other drugs into the epidural space.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Exhibit #2 510(k) Summary

This 510(k) Summary is being submitted in accordance with requirements of Title 21, CFR Section 807.92.

The assigned 510(k) Number: K171518

1. Date of Preparation: 01/02/2017
2. Sponsor Identification

Zhejiang kindly medical devices Co., Ltd.

No.758, 5th Binhai Road, Binhai Industrial Park, Longwan District,
325025 Wenzhou, Zhejiang Province, PRC.

Establishment Registration Number: Not yet registered

Contact Person: Zenghua Zhang

Position: General Manager

Tel: +86-577-86872005

Fax: +86-577-86374972

Email: zhangzh@kdlchina.com

3. Designated Submission Correspondent

Ms. Diana Hong (Primary Contact Person)

Ms. Ying Xu (Alternative Contact Person)

Mid-Link Consulting Co., Ltd

P.O. Box 120-119, Shanghai, 200120, China

Tel: +86-21-22815850,

Fax: 240-238-7587

Email: info@mid-link.net

4. Identification of Proposed Device

Trade Name: Spinal Anesthesia Needles

Epidural Anesthesia Needles

Combined Anesthesia Needles

Common Name: Anesthesia conduction needle

Regulatory Information

Classification Name: Anesthesia conduction needle

Classification: II;

Product Code: BSP;

Regulation Number: 21CFR 868.5150

Review Panel: Anesthesiology

Intended Use Statement:

The Epidural Anesthesia Needles are intended to be used for injection into the epidural space/or placing the epidural catheter into the epidural space.

The Spinal Anesthesia Needles are intended to be used for injection of local anesthetic agent into the subarachnoid cavity for pain management.

The Combined Anesthesia Needles are intended for injection of local anesthetics into the spinal and epidural spaces of a patient to provide regional anesthesia. The administration of the spinal anesthesia allows rapid anesthesia onset and the placement of an epidural catheter allows for bolus injections or continuous infusion of local anesthetics or other drugs into the epidural space.

Device Description

The proposed device, Epidural Anesthesia Needles are available in a series combination of needle size and length. The needle tubing is stabilized during puncture with use of an inner stylet. The stylet is withdrawn after the epidural anesthesia needle has reached its anatomical site for neuraxial anesthesia. Then an epidural catheter is introduced into epidural cavities for convenience of injecting anesthetic continuously.

The proposed device, Spinal Anesthesia Needles are available in a series combination of needle size and length. The spinal anesthesia needle has a tightly fitting removable stylet that completely occludes the lumen to avoid block. The stylet is withdrawn after the spinal anesthesia needle has penetrated into the subarachnoid space for injecting anesthetic. The needles are available in Sprotte and Quincke two type.

The proposed device, Combination Anesthesia Needle are instruments for a spinal (subarachnoid) injection of anesthetics, followed by the placement of an epidural catheter to allow modification of the spinal analgesia if necessary, or continuous infusion of local anesthetics into the epidural space for subsequent pain relief if required. The needles are a matched set. The epidural needle and the spinal needle are locked by conical fitting that enable the spinal needle to be locked into position once the dura has been pierced so that it is secured to the epidural needle to prevent accidental displacement.

5. Identification of Predicate Devices

Predicate Device 1

510(k) Number: K142553

Product Name: Uniever Disposable Epidural Anesthesia Needle

Predicate Device 2

510(k) Number: K141126

Product Name: UNIEVER Disposable Spinal Anesthesia Needle

Predicate Device 3

510(k) Number: K993619

Product Name: CSEcure™ Combined Spinal/Epidural Anesthesia System with Lock

6. Non-Clinical Test Conclusion

Non clinical tests were conducted to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The test results demonstrated that the proposed device complies with the following standards:

- ISO 9626:1991 AMD 1:2001 Stainless Steel Needle Tubing for the Manufacture of Medical Devices.
- ISO 7864:1993 Sterile Hypodermic Needles for Single Use
- ISO 594-1:1986 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment- Part 1: General requirements
- ISO 594-2:1998 Conical fittings with 6% (Luer) taper for syringes, needles and certain other medical equipment- Part 2: Lock fittings
- ISO 10993-7:2008 Biological evaluation of medical devices- Part 7: Ethylene oxide sterilization residuals;
- ASTM F88/88M-09 Standard test method for seal strength of flexible barrier materials;
- ASTM F1140/1140M-13 Standard test methods for internal pressurization failure resistance of

unrestrained packages

- USP38-NF33 <85> Bacterial Endotoxins Test
- ISO 10993-5:2009 Biological evaluation of medical devices-Part 5: Tests for in vitro cytotoxicity.
- ISO 10993-10:2010 Biological evaluation of medical devices- Part 10: Test for irritation and delayed-type hypersensitivity
- ISO 10993-11:2006 Biological evaluation of medical devices- Part 11: Tests for systemic toxicity
- ISO 10993-4: 2002/Amd 1: 2006 Biological evaluation of medical devices- Part 4: Selection of tests for interactions with blood

7. Clinical Test Conclusion

No clinical study is included in this submission.

8. Substantially Equivalent (SE) Comparison

Table 1 Comparison of Technology Characteristics for Epidural Anesthesia Needles

Item	Proposed Device	Predicate Device 1 K142553
Product Code	BSP	Same
Regulation Number	21 CFR 868.5150	Same
Class	II	Same
Intended Use	The Epidural Anesthesia Needles are intended to be used for injection into the epidural space/or placing the epidural catheter into the epidural space.	The Uniever Disposable Epidural Anesthesia Needle is intended to be used for injection into the epidural space/or placing the epidural catheter into the epidural space.
Configuration	Stylet Hub Stylet Needle Tube Needle Hub Needle Hub Insert Protective Cap of Needle	Stylet Hub Stylet Needle Tube Needle Hub Needle Hub Insert Protective Cap of Needle
Needle Gauge	14G~22G	14G~25G
Needle Length	65mm, 70mm, 80mm, 90mm, 100mm, 110mm, 120mm, 150mm	30~160mm
Sterile	EO Sterilized 10 ⁻⁶	Same
Single Use	Single Use	Same
Needle Tip	Tuohy	Huber, Hustead and Crawford
Label/Labeling	Conform with 21 CFR 801	Same
Biocompatibility		
Cytotoxicity	No Cytotoxicity	Conform with ISO 10993
Intracutaneous Reactivity	No Intracutaneous Reactivity	
Skin Sensitization	No Sensitization	
Acute Systemic Toxicity	No Systemic Toxicity	
Hemolysis	No Hemolysis	
Pyrogen	No Pyrogen	

Table 2 Comparison of Technology Characteristics for Spinal Anesthesia Needles

Item	Proposed Device	Predicate Device 2 K141126
Product Code	BSP	Same
Regulation Number	21 CFR 868.5150	Same
Class	II	Same
Intended Use	The Spinal Anesthesia Needles are intended to be used for injection of local anesthetic agent into the subarachnoid cavity for pain management.	UNIEVER Disposable Spinal Anesthesia Needle is intended to be used for injection of local anesthetic agent into the subarachnoid cavity for pain management.
Configuration	Stylet Hub Stylet Needle Tube Needle Hub Needle Hub Insert Protective Cap of Needle Cannula Needle Tube	Stylet Hub Stylet Needle Tube Needle Hub Needle Hub Insert Protective Cap of Needle Cannula Needle Tube
Needle Gauge	18G~27G	18G~29G
Needle Length	40mm, 50mm, 65mm, 70mm, 80mm, 90mm, 120mm	30~150mm
Sterile	EO Sterilized	Same
	10 ⁻⁶	
Needle Tip	Quincke and Pencile	Sharp and blunt
Single Use	Single Use	Same
Label/Labeling	Conform with 21 CFR 801	Same
Biocompatibility		
Cytotoxicity	No Cytotoxicity	Conform with ISO 10993
Intracutaneous Reactivity	No Intracutaneous Reactivity	
Skin Sensitization	No Sensitization	
Acute Systemic Toxicity	No Systemic Toxicity	
Hemolysis	No Hemolysis	
Pyrogen	No Pyrogen	

Table 3 Comparison of Technology Characteristics for Combined Anesthesia Needles

Item	Proposed Device		Predicate Device 3 K993619	
Product Code	BSP		Same	
Regulation Number	21 CFR 868.5150		Same	
Class	II		Same	
Intended Use	The Combined Anesthesia Needles are intended for injection of local anesthetics into the spinal and epidural spaces of a patient to provide regional anesthesia. The administration of the spinal anesthesia allows rapid anesthesia onset and the placement of an epidural catheter allows for bolus injections or continuous infusion of local anesthetics or other drugs into the epidural space.		The combined spinal epidural needle kit is intended for injection of local anesthetics into the spinal and epidural spaces of a patient to provide regional anesthesia. The administration of the spinal anesthesia allows rapid anesthesia onset and the placement of an epidural catheter allows for bolus injections or continuous infusion of local anesthetics or other drugs into the epidural space.	
Configuration	Stylet Hub Stylet Needle Tube Needle Hub Needle Hub Insert Protective Cap of Needle Cannula Needle Tube		Stylet Hub Stylet Needle Tube Needle Hub Needle Hub Insert Protective Cap of Needle Cannula Needle Tube	
Needle Gauge and Length	Epidural Needle	Spinal Needle	Epidural Needle	Spinal Needle
	14G~18G	22G~27G	17G, 18G	26G, 27G
	60~120mm	110~150mm	90mm	132mm
Sterile	EO Sterilized		Same	
	10 ⁻⁶			
Single Use	Single Use		Same	
Label/Labeling	Conform with 21 CFR 801		Same	
Biocompatibility				
Cytotoxicity	No Cytotoxicity		Conform with ISO 10993	
Intracutaneous Reactivity	No Intracutaneous Reactivity			
Skin Sensitization	No Sensitization			
Acute Systemic Toxicity	No Systemic Toxicity			
Hemolysis	No Hemolysis			
Pyrogen	No Pyrogen			

From above comparison table, the proposed devices have different combination of needle gauge and needle length compared with predicate devices. However, the different specification will be selected by physician per patient's condition. Therefore, this difference is not considered to affect substantially equivalence.

9. Substantially Equivalent (SE) Conclusion

Based on the comparison and analysis above, the proposed devices are determined to be Substantially Equivalent (SE) to the predicate devices.