



May 31, 2018

Orantech Inc.  
% Mei Tan  
RA Consultant  
Chonconn Medical Device Consulting Co. LTD.  
22A, Hai Jing Square, No.18 Taizi Road, Nanshan District,  
Shenzhen, 51800  
China

Re: K171828  
Trade/Device Name: ETCO2 Sensor, Model CTM-RP01  
Regulation Number: 21 CFR 868.1400  
Regulation Name: Carbon Dioxide Gas Analyzer  
Regulatory Class: Class II  
Product Code: CCK  
Dated: April 30, 2018  
Received: April 30, 2018

Dear Mei Tan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Geeta K.  
Pamidimukkala -S

for Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K171828

Device Name

EtCO2 Sensor, Model CTM-RP01

Indications for Use (Describe)

The intended use of the ETCO2 sensor is to provide carbon dioxide monitoring to a host monitoring system during anesthesia / recovery, in the intensive care unit (ICU), and in Emergency Medicine/Transport or Respiratory care.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(K) Summary K171828

### 1. Prepared Date: 2018/5/30

### 2. Submitter Information

|                |                                                                                                                  |
|----------------|------------------------------------------------------------------------------------------------------------------|
| Name           | Orantech Inc.                                                                                                    |
| Address        | Zone#A,4F,1st Bld,7th Industrial Zone,Yulv Community,<br>GongMing, Guangming New District, Shenzhen,China 518106 |
| Contact person | Mr. Hsin Xiong                                                                                                   |
| Tel            | 0755-23699939                                                                                                    |
| Fax            | 0755-23243594                                                                                                    |
| Email          | Yunxi.xiong@orantech.cn                                                                                          |

### 3. Submission Correspondent

|                |                                                                                                                             |
|----------------|-----------------------------------------------------------------------------------------------------------------------------|
| Name           | Chonconn Medical Device Consulting Co.LTD.                                                                                  |
| Address        | 22A,HaiJing Square,No.18 Taizi Road,Nanshan District,<br>Shenzhen,Guangdong,P.R.China,                                      |
| Contact person | Mei(RA consultant)                                                                                                          |
| E-mail         | <a href="mailto:FDA@chonconn.com">FDA@chonconn.com</a> / <a href="mailto:Xinmeitan@chonconn.com">Xinmeitan@chonconn.com</a> |
| Tel            | +86-755-33941160                                                                                                            |

### 4. Proposed Device Information

|                       |                                            |
|-----------------------|--------------------------------------------|
| Trade Name            | ETCO <sub>2</sub> Sensor                   |
| Model                 | CTM-RP01                                   |
| Common name           | Carbon Dioxide (CO <sub>2</sub> ) Analyzer |
| Regulatory class      | II                                         |
| Production regulation | 21 CFR §868.1400                           |
| Product code          | CCK                                        |
| Panel                 | Anesthesiology                             |

### 5. Predicate Device Information

| 510(K)No. | Trade Name/model                   | Submitter                   |
|-----------|------------------------------------|-----------------------------|
| K042601   | Capnostat 5 CO <sub>2</sub> sensor | Respironics Novamatrix, LLC |

### 6. Device description

The ETCO<sub>2</sub> Sensor is comprised of three main components, a sensor with photo detector and light emitter, cable and connector, The CO<sub>2</sub> Sensor incorporates an infrared light source, of specified wavelength, and an infrared detector. The photo ETCO<sub>2</sub> Sensor-510(k) Submission

detector and light emitter end of the CO<sub>2</sub> sensor is connected to an airway adapter. The airway adapter is connected between the patient airway and the respirator. As the patient completes an expiratory breath the sensor measures the CO<sub>2</sub> levels in the expiratory breath and sends that data to the compatible CO<sub>2</sub> monitor. The monitor then reads the data and converts the data so it can be displayed on screen.

The ETCO<sub>2</sub> sensor are described as follows:

| Model    | Description      | Compatible Monitor                              |
|----------|------------------|-------------------------------------------------|
| CTM-RP01 | Mainstream/ 3.0M | General Meditech INC. Patient Monitor model G3C |

## 7. Intended use & Indications for Use

The intended use of the ETCO<sub>2</sub> sensor is to provide carbon dioxide monitoring to a host monitoring system during anesthesia / recovery, in the intensive care unit (ICU), and in Emergency Medicine/Transport or Respiratory care.

## 8. Comparison to predicate device

The following is provided as a summary of how the technological characteristics of the device compare to the predicate device

Table1 Comparison between CTM-RP01 sensor and Capnostat 5 CO<sub>2</sub> sensor

| Comparison item                    | Subject Device                                                                                                                                                                                                                                | Predicate Device                                                                                                                                                                                                                                       | Comments |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| 510(K) number                      | Pending application                                                                                                                                                                                                                           | K042601                                                                                                                                                                                                                                                | --       |
| Device                             | CTM-RP01                                                                                                                                                                                                                                      | Capnostat 5 CO <sub>2</sub> sensor                                                                                                                                                                                                                     | --       |
| Intended use & Indications for Use | The intended use of the ETCO <sub>2</sub> sensor is to provide carbon dioxide monitoring to a host monitoring system during anesthesia / recovery, in the intensive care unit (ICU), and in Emergency Medicine/Transport or Respiratory care. | The intended use of the Capnostat 5 CO <sub>2</sub> sensor is to provide carbon dioxide monitoring to a host monitoring system during anesthesia / recovery, in the intensive care unit (ICU), and in Emergency Medicine/Transport or Respiratory care | Same     |
| Environment of Use                 | Hospitals, clinics, and other healthcare environments                                                                                                                                                                                         | Hospitals, clinics, and other healthcare environments                                                                                                                                                                                                  | Same     |
| Measurement Technique              | Non-dispersive infrared (NDIR) Infrared spectroscopy                                                                                                                                                                                          | Non-dispersive infrared (NDIR) Infrared spectroscopy                                                                                                                                                                                                   | Same     |
| CO <sub>2</sub> measurement Range  | 0 to 150 mmHg, 0 to 19.7%                                                                                                                                                                                                                     | 0 to 150 mmHg, 0 to 19.7%                                                                                                                                                                                                                              | Same     |
| CO <sub>2</sub> Resolution         | 0.1 mmHg 0 to 69 mmHg<br>0.25 mmHg 70 to 150 mmHg                                                                                                                                                                                             | 0.1 mmHg 0 to 69 mmHg<br>0.25 mmHg 70 to 150 mmHg                                                                                                                                                                                                      | Same     |

|                           |                                                                                                                                                                                                            |                                                                                                                                                                                                            |         |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|
| CO <sub>2</sub> Accuracy  | 0 - 40 mmHg ± 2 mmHg<br>41 - 70 mmHg ± 5% of reading<br>71 - 100 mmHg ± 8% of reading<br>101 - 150 mmHg ± 10% of reading<br>Above 80 breath per minute ± 12% of reading<br>* NOTE: Gas temperature at 25°C | 0 - 40 mmHg ± 2 mmHg<br>41 - 70 mmHg ± 5% of reading<br>71 - 100 mmHg ± 8% of reading<br>101 - 150 mmHg ± 10% of reading<br>Above 80 breath per minute ± 12% of reading<br>* NOTE: Gas temperature at 25°C | Same    |
| Respiration Rate Range    | 0 to 150 breaths per minute (bpm)                                                                                                                                                                          | 0 to 150 breaths per minute (bpm)                                                                                                                                                                          | Same    |
| Respiration Rate Accuracy | ± 1 breaths per minute (bpm)                                                                                                                                                                               | ± 1 breaths per minute (bpm)                                                                                                                                                                               | Same    |
| Voltage Requirements      | 5.0 VDC ±5%                                                                                                                                                                                                | 5.0 VDC ±5%                                                                                                                                                                                                | Same    |
| Interconnection           | Lemo Redel 8-pin plastic                                                                                                                                                                                   | Lemo Redel 8-pin plastic                                                                                                                                                                                   | Same    |
| Temperature and Humidity  | Operating: 0 to 40°C, 15 to 85% RH, non-condensing<br>Storage: -10 to 40°C, <90% RH, non-condensing                                                                                                        | Operating: 0 to 45°C, 10 to 90% RH, non-condensing<br>Storage: -40 to 70°C, <90% RH, non-condensing                                                                                                        | Similar |
| Water Resistance          | IPX2 – Splash-proof (sensor head only)                                                                                                                                                                     | IPX4 – Splash-proof (sensor head only)                                                                                                                                                                     | Similar |
| Compliance                | IEC 60601-1(safety)<br>IEC 60601-1-2(EMC)<br>ISO 80601-2-55(Performance)                                                                                                                                   | IEC 60601-1(safety)<br>IEC 60601-1-2(EMC)<br>ISO80601-2-55(Performance)                                                                                                                                    | Same    |

From the comparison form above, the subject and predicate devices are based on the following same technological elements:

- Same measurement techniques
- Same intended use & indications for Use
- Same CO<sub>2</sub> Resolution & Accuracy
- Same respiration rate range & Accuracy
- Same voltage Requirements
- Both comply to the same FDA recognized standards

The following slightly technological differences exist between the subject and predicate devices:

- Similar temperature and humidity range
- Water Resistance of predicate device is higher than Orantech's

These differences do not raise different questions of safety and effectiveness.

## **9. Clinical test data**

Not applicable.

## **10. Non-clinical test data**

1) Industry Standards for Electrical Safety, EMC and Essential Performance  
Testing of the ECO<sub>2</sub> Sensor has been completed to verify compliance with FDA recognized national and international standards for electrical safety and performance for medical devices, and particular requirements applicable to this device including:

- IEC 60601-1:2012 Basic safety and essential performance
- IEC 60601-1-2:2007 and 2014 EMC
- ISO 80601-2-55:2011 Respiratory gas monitors

### 2) Biocompatibility

The ETCO<sub>2</sub> Sensor does not contain any direct or indirect patient contacting components. Patient contacting accessories are previously cleared by FDA. Biocompatibility test is not applicable for subject device.

### 3) Software Verification

Software verification and validation testing were conducted and documentation was provided as recommended by FDA Software Guidance.

## **11. Conclusion**

The results for testing demonstrate that the ETCO<sub>2</sub> Sensor is substantially equivalent to the above listed predicate device.